

GE HealthCare

Renewal Parts and Tooling for 1.5T and 3.0T MRI Magnets

1.5T R, CR, RD and HM Series Magnets

3.0T W, WB and UA Series Magnets



Operating Documentation

DOC1472639

Revision 03

Important Information, Language

ПРЕДУПРЕЖДЕНИЕ

(BG)

Това упътване за работа е налично само на английски език.

- Ако доставчикът на услугата на клиента изиска друг език, задължение на клиента е да осигури превод.
- Не използвайте оборудването, преди да сте се консултирали и разбрали упътването за работа.
- Неспазването на това предупреждение може да доведе до нараняване на доставчика на услугата, оператора или пациента в резултат на токов удар, механична или друга опасност.

警告

(ZH-CN)

本维修手册仅提供英文版本。

- 如果客户的维修服务人员需要非英文版本，则客户需自行提供翻译服务。
- 未详细阅读和完全理解本维修手册之前，不得进行维修。
- 忽略本警告可能对维修服务人员、操作人员或患者造成电击、机械伤害或其他形式的伤害。

警告

(ZH-HK)

本服務手冊僅提供英文版本。

- 倘若客戶的服務供應商需要英文以外之服務手冊，客戶有責任提供翻譯服務。
- 除非已參閱本服務手冊及明白其內容，否則切勿嘗試維修設備。
- 不遵從本警告或會令服務供應商、網絡供應商或病人受到觸電、機械性或其他危險。

警告

(ZH-TW)

本維修手冊僅有英文版。

- 若客戶的維修廠商需要英文版以外的語言，應由客戶自行提供翻譯服務。
- 請勿試圖維修本設備，除非 您已查閱並瞭解本維修手冊。
- 若未留意本警告，可能導致維修廠商、操作員或病患因觸電、機械或其他危險而受傷。

UPOZORENJE

(HR)

Ovaj servisni priručnik dostupan je na engleskom jeziku.

- Ako davatelj usluge klijenta treba neki drugi jezik, klijent je dužan osigurati prijevod.
- Ne pokušavajte servisirati opremu ako niste u potpunosti pročitali i razumjeli ovaj servisni priručnik.
- Zanimarite li ovo upozorenje, može doći do ozljede davatelja usluge, operatera ili pacijenta uslijed strujnog udara, mehaničkih ili drugih rizika.

GE Healthcare**Renewal Parts and Tooling****VÝSTRAHA**

(CS)

Tento provozní návod existuje pouze v anglickém jazyce.

- V případě, že externí služba zákazníkům potřebuje návod v jiném jazyce, je zajištění překladu do odpovídajícího jazyka úkolem zákazníka.
- Nesnažte se o údržbu tohoto zařízení, aniž byste si přečetli tento provozní návod a pochopili jeho obsah.
- V případě nedodržování této výstrahy může dojít k poranění pracovníka prodejního servisu, obslužného personálu nebo pacientů vlivem elektrického proudu, respektive vlivem mechanických či jiných rizik.

ADVARSEL

(DA)

Denne servicemanual findes kun på engelsk.

- Hvis en kundes tekniker har brug for et andet sprog end engelsk, er det kundens ansvar at sørge for oversættelse.
- Forsøg ikke at servicere udstyret uden at læse og forstå denne servicemanual.
- Manglende overholdelse af denne advarsel kan medføre skade på grund af elektrisk stød, mekanisk eller anden fare for teknikeren, operatøren eller patienten.

WAARSCHUWING

(NL)

Deze onderhoudshandleiding is enkel in het Engels verkrijgbaar.

- Als het onderhoudspersoneel een andere taal vereist, dan is de klant verantwoordelijk voor de vertaling ervan.
- Probeer de apparatuur niet te onderhouden alvorens deze onderhoudshandleiding werd geraadpleegd en begrepen is.
- Indien deze waarschuwing niet wordt opgevolgd, zou het onderhoudspersoneel, de operator of een patiënt gewond kunnen raken als gevolg van een elektrische schok, mechanische of andere gevaren.

WARNING

(EN)

This service manual is available in English only.

- If a customer's service provider requires a language other than English, it is the customer's responsibility to provide translation services.
- Do not attempt to service the equipment unless this service manual has been consulted and is understood.
- Failure to heed this warning may result in injury to the service provider, operator or patient from electric shock, mechanical or other hazards.

GE Healthcare**Renewal Parts and Tooling****HOIATUS****(ET)**

See teenindusjuhend on saadaval ainult inglise keeles.

- Kui klienditeeninduse osutaja nõuab juhendit inglise keelest erinevas keeles, vastutab klient tõlketeenuse osutamise eest.
- Ärge üritage seadmeid teenindada enne eelnevalt käesoleva teenindusjuhendiga tutvumist ja sellest aru saamist.
- Käesoleva hoiatuse eiramine võib põhjustada teenuseosutaja, operaatori või patsiendi vigastamist elektrilöögi, mehaanilise või muu ohu tagajärjel.

VAROITUS**(FI)**

Tämä huolto-ohje on saatavilla vain englanniksi.

- Jos asiakkaan huoltohenkilöstö vaatii muuta kuin englanninkielistä materiaalia, tarvittavan käännöksen hankkiminen on asiakkaan vastuulla.
- Älä yritä korjata laitteistoa ennen kuin olet varmasti lukenut ja ymmärtänyt tämän huolto-ohjeen.
- Mikäli tätä varoitusta ei noudateta, seurauksena voi olla huoltohenkilöstön, laitteiston käyttäjän tai potilaan vahingoittuminen sähköiskun, mekaanisen vian tai muun vaaratilanteen vuoksi.

ATTENTION**(FR)**

Ce manuel d'installation et de maintenance est disponible uniquement en anglais.

- Si le technicien d'un client a besoin de ce manuel dans une langue autre que l'anglais, il incombe au client de le faire traduire.
- Ne pas tenter d'intervenir sur les équipements tant que ce manuel d'installation et de maintenance n'a pas été consulté et compris.
- Le non-respect de cet avertissement peut entraîner chez le technicien, l'opérateur ou le patient des blessures dues à des dangers électriques, mécaniques ou autres.

WARNUNG**(DE)**

Diese Serviceanleitung existiert nur in englischer Sprache.

- Falls ein fremder Kundendienst eine andere Sprache benötigt, ist es Aufgabe des Kunden für eine entsprechende Übersetzung zu sorgen.
- Versuchen Sie nicht diese Anlage zu warten, ohne diese Serviceanleitung gelesen und verstanden zu haben.
- Wird diese Warnung nicht beachtet, so kann es zu Verletzungen des Kundendiensttechnikers, des Bedieners oder des Patienten durch Stromschläge, mechanische oder sonstige Gefahren kommen.

GE Healthcare**Renewal Parts and Tooling****ΠΡΟΕΙΔΟΠΟΙΗΣΗ**

(EL)

Το παρόν εγχειρίδιο σέρβις διατίθεται μόνο στα αγγλικά.

- Εάν ο τεχνικός σέρβις ενός πελάτη απαιτεί το παρόν εγχειρίδιο σε γλώσσα εκτός των αγγλικών, αποτελεί ευθύνη του πελάτη να παρέχει τις υπηρεσίες μετάφρασης.
- Μην επιχειρήσετε την εκτέλεση εργασιών σέρβις στον εξοπλισμό αν δεν έχετε συμβουλευτεί και κατανοήσει το παρόν εγχειρίδιο σέρβις.
- Αν δεν προσέξετε την προειδοποίηση αυτή, ενδέχεται να προκληθεί τραυματισμός στον τεχνικό σέρβις, στο χειριστή ή στον ασθενή από ηλεκτροπληξία, μηχανικούς ή άλλους κινδύνους.

FIGYELMEZTETÉS

(HU)

Ezen karbantartási kézikönyv kizárólag angol nyelven érhető el.

- Ha a vevő szolgáltatója angoltól eltérő nyelvre tart igényt, akkor a vevő felelőssége a fordítás elkészítése.
- Ne próbálja elkezdni használni a berendezést, amíg a karbantartási kézikönyvben leírtakat nem értelmezték.
- Ezen figyelmeztetés figyelmen kívül hagyása a szolgáltató, működtető vagy a beteg áramütés, mechanikai vagy egyéb veszélyhelyzet miatti sérülését eredményezheti.

AÐVÖRUN

(IS)

Þessi þjónustuhandbók er aðeins fánleg á ensku.

- Ef að þjónustuveitandi viðskiptamanns þarfnast annas tungumáls en ensku, er það skylda viðskiptamanns að skaffa tungumálþjónustu.
- Reynið ekki að afgreiða tækið nema að þessi þjónustuhandbók hefur verið skoðuð og skilin.
- Brot á sinna þessari aðvörun getur leitt til meiðsla á þjónustuveitanda, stjórnanda eða sjúklings frá raflosti, vélrænu eða öðrum áhættum.

AVVERTENZA

(IT)

Il presente manuale di manutenzione è disponibile soltanto in lingua inglese.

- Se un addetto alla manutenzione richiede il manuale in una lingua diversa, il cliente è tenuto a provvedere direttamente alla traduzione.
- Procedere alla manutenzione dell'apparecchiatura solo dopo aver consultato il presente manuale ed averne compreso il contenuto.
- Il mancato rispetto della presente avvertenza potrebbe causare lesioni all'addetto alla manutenzione, all'operatore o ai pazienti provocate da scosse elettriche, urti meccanici altri rischi.

GE Healthcare**Renewal Parts and Tooling****警告**

(JA)

このサービスマニュアルには英語版しかありません。

- サービスを担当される業者が英語以外の言語を要求される場合、翻訳作業はその業者の責任で行うものとさせていただきます。
- このサービスマニュアルを熟読し理解せずに、装置のサービスを行わないでください。
- この警告に従わない場合、サービスを担当される方、操作員あるいは患者さんが、感電や機械的又はその他の危険により負傷する可能性があります。

경고

(KO)

본 서비스 매뉴얼은 영어로만 이용하실 수 있습니다.

- 고객의 서비스 제공자가 영어 이외의 언어를 요구할 경우, 번역 서비스를 제공하는 것은 고객의 책임입니다.
- 본 서비스 매뉴얼을 참조하여 숙지하지 않은 이상 해당 장비를 수리하려고 시도하지 마십시오.
- 본 경고 사항에 유의하지 않으면 전기 쇼크, 기계적 위험, 또는 기타 위험으로 인해 서비스 제공자, 사용자 또는 환자에게 부상을 입힐 수 있습니다.

BRĪDINĀJUMS

(LV)

Šī apkopes rokasgrāmata ir pieejama tikai angļu valodā.

- Ja klienta apkopes sniedzējam nepieciešama informācija citā valodā, klienta pienākums ir nodrošināt tulkojumu.
- Neveiciet aprikojuma apkopi bez apkopes rokasgrāmatas izlasīšanas un saprašanas.
- Šī brīdinājuma neievērošanas rezultātā var rasties elektriskās strāvas trieciena, mehānisku vai citu faktoru izraisītu traumu risks apkopes sniedzējam, operatoram vai pacientam.

ĮSPĖJIMAS

(LT)

Šis eksploatavimo vadovas yra tik anglų kalba.

- Jei kliento paslaugų tiekėjas reikalauja vadovo kita kalba – ne anglų, suteikti vertimo paslaugas privalo klientas.
- Nemėginkite atlikti įrangos techninės priežiūros, jei neperskaitėte ar nesupratote šio eksploatavimo vadovo.
- Jei nepaisysite šio įspėjimo, galimi paslaugų tiekėjo, operatoriaus ar paciento sužalojimai dėl elektros šoko, mechaninių ar kitų pavojų.

GE Healthcare**Renewal Parts and Tooling****ADVARSEL**

(NO)

Denne servicehåndboken finnes bare på engelsk.

- Hvis kundens serviceleverandør har bruk for et annet språk, er det kundens ansvar å sørge for oversettelse.
- Ikke forsøk å reparere utstyret uten at denne servicehåndboken er lest og forstått.
- Manglende hensyn til denne advarselen kan føre til at serviceleverandøren, operatøren eller pasienten skades på grunn av elektrisk støt, mekaniske eller andre farer.

OSTRZEŻENIE

(PL)

Niniejszy podręcznik serwisowy dostępny jest jedynie w języku angielskim.

- Jeśli serwisant klienta wymaga języka innego niż angielski, zapewnienie usługi tłumaczenia jest obowiązkiem klienta.
- Nie próbować serwisować urządzenia bez zapoznania się z niniejszym podręcznikiem serwisowym i zrozumienia go.
- Niezastosowanie się do tego ostrzeżenia może doprowadzić do obrażeń serwisanta, operatora lub pacjenta w wyniku porażenia prądem elektrycznym, zagrożenia mechanicznego bądź innego.

ATENÇÃO

(PT-BR)

Este manual de assistência técnica encontra-se disponível unicamente em inglês.

- Se outro serviço de assistência técnica solicitar a tradução deste manual, caberá ao cliente fornecer os serviços de tradução.
- Não tente reparar o equipamento sem ter consultado e compreendido este manual de assistência técnica.
- A não observância deste aviso pode ocasionar ferimentos no técnico, operador ou paciente decorrentes de choques elétricos, mecânicos ou outros.

ATENÇÃO

(PT-PT)

Este manual de assistência técnica só se encontra disponível em inglês.

- Se qualquer outro serviço de assistência técnica solicitar este manual noutro idioma, é da responsabilidade do cliente fornecer os serviços de tradução.
- Não tente reparar o equipamento sem ter consultado e compreendido este manual de assistência técnica.
- O não cumprimento deste aviso pode colocar em perigo a segurança do técnico, do operador ou do paciente devido a choques eléctricos, mecânicos ou outros.

GE Healthcare**Renewal Parts and Tooling****ATENȚIE****(RO)**

Acest manual de service este disponibil doar în limba engleză.

- Dacă un furnizor de servicii pentru clienți necesită o altă limbă decât cea engleză, este de datoria clientului să furnizeze o traducere.
- Nu încercați să reparați echipamentul decât ulterior consultării și înțelegerii acestui manual de service.
- Ignorarea acestui avertisment ar putea duce la rănirea depanatorului, operatorului sau pacientului în urma pericolelor de electrocutare, mecanice sau de altă natură.

ОСТОРОЖНО!**(RU)**

Данное руководство по техническому обслуживанию представлено только на английском языке.

- Если сервисному персоналу клиента необходимо руководство не на английском, а на каком-то другом языке, клиенту следует самостоятельно обеспечить перевод.
- Перед техническим обслуживанием оборудования обязательно обратитесь к данному руководству и поймите изложенные в нем сведения.
- Несоблюдение требований данного предупреждения может привести к тому, что специалист по техобслуживанию, оператор или пациент получит удар электрическим током, механическую травму или другое повреждение.

UPOZORENJE**(SR)**

Ovo servisno uputstvo je dostupno samo na engleskom jeziku.

- Ako klijentov serviser zahteva neki drugi jezik, klijent je dužan da obezbedi prevodilačke usluge.
- Ne pokušavajte da opravite uređaj ako niste pročitali i razumeli ovo servisno uputstvo.
- Zanemarivanje ovog upozorenja može dovesti do povređivanja servisera, rukovaoca ili pacijenta usled strujnog udara ili mehaničkih i drugih opasnosti.

UPOZORNENIE**(SK)**

Tento návod na obsluhu je k dispozícii len v angličtine.

- Ak zákazníkovi poskytovateľ služieb vyžaduje iný jazyk ako angličtinu, poskytnutie prekladateľských služieb je zodpovednosťou zákazníka.
- Nepokúšajte sa o obsluhu zariadenia, kým si neprečítate návod na obsluhu a neporozumiete mu.
- Zanedbanie tohto upozornenia môže spôsobiť zranenie poskytovateľa služieb, obsluhujúcej osoby alebo pacienta elektrickým prúdom, mechanické alebo iné ohrozenie.

GE Healthcare**Renewal Parts and Tooling****ATENCION**

(ES)

Este manual de servicio sólo existe en inglés.

- Si el encargado de mantenimiento de un cliente necesita un idioma que no sea el inglés, el cliente deberá encargarse de la traducción del manual.
- No se deberá dar servicio técnico al equipo, sin haber consultado y comprendido este manual de servicio.
- La no observancia del presente aviso puede dar lugar a que el proveedor de servicios, el operador o el paciente sufran lesiones provocadas por causas eléctricas, mecánicas de otra naturaleza.

VARNING

(SV)

Den här servicehandboken finns bara tillgänglig på engelska.

- Om en kunds servicetekniker har behov av ett annat språk än engelska, ansvarar kunden för att tillhandahålla översättningstjänster.
- Försök inte utföra service på utrustningen om du inte har läst och förstår den här servicehandboken.
- Om du inte tar hänsyn till den här varningen kan det resultera i skador på serviceteknikern, operatören eller patienten till följd av elektriska stötar, mekaniska faror eller andra faror.

OPOZORILO

(SL)

Ta servisni priročnik je na voljo samo v angleškem jeziku.

- Če ponudnik storitve stranke potrebuje priročnik v drugem jeziku, mora stranka zagotoviti prevod.
- Ne poskušajte servisirati opreme, če tega priročnika niste v celoti prebrali in razumeli.
- Če tega opozorila ne upoštevate, se lahko zaradi električnega udara, mehanskih ali drugih nevarnosti poškoduje ponudnik storitev, operater ali bolnik.

DİKKAT

(TR)

Bu servis kılavuzunun sadece ingilizcesi mevcuttur.

- Eğer müşteri teknisyeni bu kılavuzu ingilizce dışında bir başka lisandan talep ederse, bunu tercüme ettirmek müşteriye düşer.
- Servis kılavuzunu okuyup anlamadan ekipmanlara müdahale etmeyiniz.
- Bu uyarıya uyulmaması, elektrik, mekanik veya diğer tehlikelerden dolayı teknisyen, operatör veya hastanın yaralanmasına yol açabilir.

GE Healthcare

Renewal Parts and Tooling

Revision History

REV	Affected Page(s)	Description of Content Changed	Reason for Change or Change Control Number	Author	Effective date of Document	Periodic review date
03	Multiple throughout the document	Newly introduced parts and tools added. Discontinued parts and tools removed. Corrected typos. Articles re-arranged	ECR2165186; ECR2166059; ECR2050875; ECR2174634; ECR2199883; ECR2150117; PCN188426; ECR2172749; ECR2196551; ECR2208446; ECR2175019; ECR2196967; ECR2206658; ECR2177535; ECR2174843; ECR2167958; ECR2222028; eNPI022066; ECR2261509; ECO2253937; ECR2103322; DOC2183415; DOC1862868; SN DOC2315815; eNPI023302; eNPI019594; ECR2213185; MR2204860; DOC2269356; Direction 5166012; ECR2178759; Direction 5166013; eNPI022869; ECO2032098; Direction 5265188; ECR2282159; ECR2174907; 5265188TPH; ECO2189951; ECR2093027; DOC2315174; ECR2167278; ECR2292067; CD32ZZ-056; CD32ZZ-060; ECO2171940; ECO2172733; ECO2175817 ECO2238819; ECO2019447; ECR2151173; eNPI020346; ECR2306382; ECR2244580; ECR2305932; ECR2092643; DOC2375398; ECR2139622; ECO2120342; ECR2163373; MR2242025; MR2234676; ECR2110409; MR2307344; ECO2185979; ECO2177535; ECO2108477; ECR2102049; MR2238685; ECR2292299; ECR2239692; ECO2114860; ECR2289055; ECR2376996; ECR2291781; ECR2259976; MR2267590; ECR2248101; MR2234185; MR2292177; ECR2287728; ECR2205927; MR2281844; ECO2096570; MR2263378; MR2285577; MR2274613; ECR2317963; ECR2321794; ECO2185914; SN DOC2383693	Nick Logan Alexey Prokofyev	My Workshop Release	Does not Apply to Service Manuals
Yes	No	NA	Enter Validation Information or Rationale if "NO" or "NA" for Validation or Training			
☒	☐		Typo or Format Change	According to applicable change above		
☐	☒	☐	Regulatory Filing			
☒	☐	☐	Validation	According to applicable change above		
☐	☒	☐	Training	No change to existing techniques.		
☒	☐	☐	Communication	Upon release.		
☒	☐	☐	Labeling, (list)	Document is a Service Technical Publication (classified as Labeling)		
☐	☒	☐	Other Documents, (list)			
☐	☒	☐	Other (describe)			

Table of Contents

Important Information, Language.....	2
Revision History	10
Table of Contents.....	11
1. Safety and Introduction	19
1.1. MR Service Safety Manual, Direction 5452735	19
1.2. Direction Application	19
1.3. Documentation Organization	19
1.4. Damage in Transportation.....	19
2. Handheld Meters	20
2.1. FRU Assemblies – Handheld Meters	20
2.1.1. 46-294052G1, Revision 0: Water Flow Meter Kit.....	20
2.2. Tool Kits – Handheld Meters	20
2.2.1. 46-301477G2, Revision 1: Lakeshore Thermometer Kit.....	20
2.2.2. 5119507, Revision 2: Cryogenic Temperature Meter [MTM101]	20
2.2.3. 362622, Revision 0: Portable Helium Meter [LHM 100]	20
2.2.4. 2270571, Revision 0: Portable Liquid Helium Level Meter with Battery Charger, Cable and Case 20	
2.2.5. 2183710, Revision 0: Magnet Specialty Tool Kit [RuO Temperature Meter]	21
2.2.6. 5886219, Revision 1: Magnet Temperature Monitor Kit [LakeShore 218]	21
3. Magnet Rundown Unit Items (MRUs) & Magnet Monitor (MagMon) Kits	22
3.1. Magnet Rundown Unit Items (MRUs)	22
3.1.1. 5180313, Revision 7: Kit, Magnet Rundown Unit.....	22
3.1.2. 5180187, Revision 3: Kit, Remote Magnet Rundown Unit	22
3.2. Magnet Monitor (MagMon) Kits	23
3.2.1. 5771373, Revision 5: Magnet Monitor 4	23
3.2.2. 5267626, Revision 3: Collector Kit, Mag Mon	23
3.2.3. 5180301, Revision 1: Remote Alarm Kit, Magnet Monitor 3	23
3.2.4. 2222520-3, Revision 3: Collector Kit, Magnet Monitor 3	23
3.2.5. 2222520-4, Revision 7: MagMon Mobile Collector Kit.....	23
3.2.6. 5193795, Revision 6: Kit, Magnet Monitor Cables, DVMR.....	24

GE Healthcare**Renewal Parts and Tooling**

3.2.7.	5548703, Revision 1: Kizuna Magnet Monitor Cable Kit	24
3.2.8.	5782696, Revision 2: RIO Mag Mon Cable Collector	24
3.2.9.	2350809, Revision 3: Core Kit, Magnet Monitor 3	25
3.2.10.	2350809-3, Revision 6: Magnet Monitor Non-LCC Upgrade HW Kit.....	25
3.2.11.	5113506, Revision 2: Rebuild Kit for Water Flow Meter 2333825-2.....	25
4.	Cryogenics	26
	Coldhead and Compressor.....	26
4.1.	FRU Assemblies – Coldhead and Compressor	26
4.1.1.	2188440-2, Revision 4: Compressor, Water Cooled, Sumitomo CSW-71D.....	26
4.1.2.	5314887, Revision 2: Compressor, Water Cooled, Sumitomo F-50H.....	26
4.1.3.	2166592, Revision 3: Compressor, Air Cooled, Sumitomo CSA-71	26
4.1.4.	5428192, Revision 3: Kit Compressor F-70H, SHIG #267446A	26
4.1.5.	5812867, Revision 1: F-50SH compressor service FRU.....	26
4.1.6.	5368565, Revision 1: Service Spares, Sumitomo F-50LV Compressor	26
4.1.7.	5840878, Revision 2: Sumitomo FA-50SH High Voltage Air-cooled Helium Compressor	26
4.1.8.	5840879, Revision 2: Sumitomo FA-50SL Low Voltage Air-cooled Helium Compressor.....	27
4.1.9.	2384329, Revision 3: Compressor, Outdoor, CNA-61D, Sumitomo	27
4.1.10.	2384329-2, Revision 2: Compressor, Indoor, CNA-61D, Sumitomo	27
4.1.11.	5491924, Revision 1: Compressor, Outdoor, CAN-61C3-C, Sumitomo	27
4.1.12.	5491924-2, Revision 1: Compressor, Indoor, CNA-61C3-E, Sumitomo	27
4.1.13.	5438959, Revision 2: Kit, A3 FRU, Sumitomo Coldhead and Gasket.....	27
4.1.14.	5814003, Revision 1: Kit A4 FRU, Sumitomo Coldhead and Gasket.....	27
4.1.15.	5814003-R, Revision 1: A4-R FRU Kit, Sumitomo Coldhead and Gasket.....	27
4.1.16.	2171620, Revision 11: Kit, Cryocooler Replacement Gaskets.....	28
4.1.17.	46-318696G1, Revision 0: Water Tee Assembly.....	28
4.1.18.	2194003-2, Revision 1: Kit, A2 FRU, Sumitomo Coldhead and Gasket.....	28
4.2.	FRU Parts – Coldhead and Compressor	29
4.3.	Tool Kits – Coldhead and Compressor.....	30
4.3.1.	46-281088G4, Revision 1: Coldhead/Compressor Maintenance Kit	30
4.3.2.	2300306-3, Revision 1: TOOL KIT, UNIVERSAL COLDHEAD REPLACE WITH INDIUM REMOVAL TOOL	31
4.3.3.	5180672, Revision 3: Coldsleeve Heater Kit	31
4.3.4.	5179792, Revision 1: Cable, CE Approved Coldsleeve Heater Extension.....	31
4.3.5.	2232674, Revision 0: Fixed Site Coldhead Change Out Kit.....	32

GE Healthcare**Renewal Parts and Tooling**

4.3.6.	2380982, Revision 0: Kit, Coldhead Insertion/Removal Tool	32
4.3.7.	2380985, Revision 0: LCC Mobile Coldhead Tool Kit	32
4.3.8.	5395677, Revision 1: Kit, Coldhead Alignment Tool Clamp	33
4.3.9.	2385069, Revision 2: 3.0T Coldhead Transport Kit	33
4.3.10.	46-294804G1, Revision AA: Nonmagnetic Coldhead Wrenches Kit.....	33
4.3.11.	46-318688P1, Revision AA: Cryocooler HE Line Oil Contamination Checker	33
4.3.12.	5102120, Revision 1: Charging Tool, Field Service, Precision Cryogenics	33
4.3.13.	5262081, Revision 1: Pressure Relief Tool, Male, -8 Coupling	33
4.3.14.	5272260, Revision 1: Pressure Relief Tool, Female, -8 Coupling.....	33
4.3.15.	2214919, Revision 0: Stand, Coldhead, LCC	34
4.3.16.	46-318784G2, Revision 1: Shield Cooler Test Kit.....	34
4.3.17.	5503375, Revision 1: Repair Kit, Cold Sleeve Pump Out Port	34
Thermal Stability.....		35
4.4.	FRU Assemblies – Thermal Stability	35
4.4.1.	5177897-3, Revision 2: Kit, Pressure Transducer Replacement.....	35
4.5.	Tool Kits – Thermal Stability	35
4.5.1.	5263845, Revision 1: De-Icing Tool Kit	35
4.5.2.	2374573, Revision 2, TOOL KIT, DEICEING	35
Helium Fill		36
4.6.	FRU Assemblies – Helium Fill.....	36
4.6.1	46-260555G2, Revision 0: Transfer Line Extension Tip for 500L Dewar	36
4.6.2	46-294513G1, Revision 4: Liquid Helium Transfer Line Kit Collector	36
4.6.3	46-271135P1, Revision AA: HOSE ASSEMBLY, 40 FT, DEWAR TO GAS CYLINDER	37
4.7.	Tool Kits –Cool Down and Top Off.....	37
4.7.1	46-306734G1, Revision AB: Helium Tank Regulator and Hose Kit	37
4.7.2	46-294705G1, Revision 5: Liquid Helium Transfer Line Kit	37
4.7.3	2166440, Revision 0: Transfer Line, Helium	37
4.7.4	5508132, Revision 3: Stinger Kit, Flexible Bottom-Fill.....	37
Ramping and Shimming.....		38
5. Ramping		38
5.1	Tool Kits - Ramping	38
5.1.1	46-260703G7, Revision 1: Magnet Ramping Equipment Kit	38
5.1.2	2236233 , Revision 0: HOLD DOWN TOOL HARDWARE KIT.....	38
5.1.3	2279142, Revision A: LCC Hold-Down Tool Assembly	38

GE Healthcare**Renewal Parts and Tooling**

5.1.4	Ramp and Shim Power Supplies	38
5.1.5	2141701, Revision 1: B0 Power Supply Kit	39
5.1.6	2353394-3, Revision 3: 1,000 Amp Ramp Cable Kit with Ground	39
5.1.7	2217315, Revision 0: Cable, Backup Ramp Down [Aux Heater Cable]	39
5.1.8	46-318314G2, Revision AA: Ramp Cable Holder	39
5.1.9	2171719, Revision 1: Ramp and Shim Cable Kit for Conquest	39
5.1.10	46-251865G5, Revision 0: Metrolab Teslameter Kit	40
5.1.11	5849221, Revision 1: METROLAB PT2026 TESLAMETER KIT	40
Shimming		41
5.2	FRU Assemblies – Shimming.....	41
5.2.1	5263916, Revision 3: Kit, Field Passive Shim, 3.0T HDv	41
5.2.2	5411648, Revision 1: LCC 1.5T Passive Shim Kit	41
5.2.3	5428096, Revision 5: MR750w Field Passive Shimming Kit.....	41
5.2.4	5342951, Revision 8: 1.5T DVw Field Passive Shimming Kit.....	42
5.2.5	2173354, Revision 10: SHIM LEAD ASM,RECTACTABLE,LCC.....	42
5.2.6	2385710, Revision 6: KIT, FIELD PASSIVE SHIM, G3.....	43
5.3	Tool Kits - Shimming	44
5.3.1	2135558, Revision 1: Shim Cable Kit	44
5.3.2	2386028, Revision 1: 3.0T Field Service Shim Camera Kit	44
5.3.3	2386028-2, Revision 1: Kit, Field Service Shim Camera, Dual 1.5T and 3.0T.....	44
5.3.4	2386028-3, Revision 1: Kit, Field Service Shim Camera, Dual 3.0T and 1.5T DVw	45
5.3.5	5266042-4, Revision 2: Kit, Universal Field Mapping Fixture, Non-Magnetic, Extended.....	45
5.3.6	5266042-5, Revision 1: Kit, Universal Field Mapping Fixture, Non-Magnetic Extended.....	46
5.3.7	5431812, Revision 2: Kit, STET Transport	46
5.3.8	5431812-2, Revision 1: KIT, STET TRANSPORT RIO AND 750	47
5.3.9	5779462, Revision 2: USB to RS-232 Conversion Kit	47
6. Gradient Related Items		48
6.1	FRU Assemblies – BRM and BRM2 gradient coils.....	48
6.1.1	2144093-5, Revision 7: BRM-D Coil, Gray Cart, Cradle	48
6.1.2	5788960, Revision 1: BRM2 WITH GRAY CART, FRU	48
6.1.3	5758309, Revision 4: SHIM TRAY ASSEMBLY, BRM2	48
6.1.4	2214166-4, Revision 20, BRM-F LCC GRAD ISO INSTALL PART LIST	49
6.1.5	2214166-5, Revision 16, BRM-F LCC GRAD ISO MOBILE INSTALL	50
6.1.6	5110336, Revision 1, Kit, BRM Gradient Cable Lead Board.....	51

GE Healthcare**Renewal Parts and Tooling**

6.2	BRM Gradient Miscellaneous	51
6.2.1	5182681-2, Revision 1, Gradient Cable Clamp Block Kit - BRM Injection Molded	51
6.2.2	5182681, Revision 1, Gradient Cable Clamp Block Kit - Injection Molded	51
6.3	FRU Assemblies – TRM gradient coil	52
6.3.1	2144093-9, Revision 5: 3T 1.5T TRM Coil, Blue Cart, Cradle	52
6.3.2	2299498, Revision 3: TRM Cooling Circuit Repair Kit	52
6.3.3	5836225, Revision 1: XRMB, TRM and BRM/BRM2 shipping brackets Collector	52
6.4	FRU Assemblies – XRMB gradient coil	53
6.4.1	2144093-14, Revision 3: Collector, XRMB FRU Gradient Coil with HOS, Gray Cart with Cradle	53
6.4.2	5335389, Revision 1: Collector, 1.5T Magnet, XRMB Gradient Mounting	53
6.4.3	5508396, Revision 1: Collector, XRMB Busbar Assembly Replacement	53
6.4.4	5147904, Revision 7: Assembly, Passive Shim Tray, XRMB Gradient	54
6.4.5	5534976, Revision 1: MR750 Safety Cover FRU Collector	54
6.4.6	5407946, Revision 1: Collector, XRMB Coolant Manifold Replacement	54
6.5	FRU Assemblies – XRMw gradient coil	55
6.5.1	5407761, Revision 6: XRMW FRU Gradient Coil Asm	55
6.5.2	5508397, Revision 4: Collector, DVw Busbar XZ Assembly Replacement	55
6.5.3	5508398, Revision 4: Collector, DVw Busbar Y Assembly Replacement	56
6.5.4	5342776, Revision 7: XRMW Gradient Coil Mounting Kit FRU	56
6.5.5	5508399, Revision 4: Collector, MR750W Busbar XZ Assembly Replacement	56
6.5.6	5508401, Revision 2: Collector, MR750w Busbar Y Assembly Replacement	57
6.5.7	5845904, Revision 2: 450w M12 Busbar XZ Assembly Replacement	57
6.5.8	5845905, Revision 1: 450w M12 Busbar Y Assembly Replacement	57
6.5.9	5845907, Revision 1: 750w M12 Busbar XZ Assembly Replacement	58
6.5.10	5845908, Revision 1: 750w M12 Busbar Y Assembly Replacement	58
6.5.11	5407963, Revision 4: XRMw Gradient Coil Mounting Kit, FRU, MR750w	58
6.5.12	5431091, Revision 3: Shim Tray Locking Ring, XRMw, FRU	58
6.5.13	5534973, Revision 1: MR450w Safety Cover FRU Collector	59
6.5.14	5534975, Revision 2: MR750w Safety Cover FRU Collector	59
6.5.15	5538535, Revision 1: Collector FRU, Clamp Replacement, MR750w	59
6.5.16	5407756, Revision 1: XRMW Coolant Manifold FRU	59
6.5.17	5850819, Revision 1: XRMW and VRMW Shipping Collector, FRU	59
6.5.18	5620700, Revision 2: collector, Kizuna 3T busbar Y assembly replacement	60
6.5.19	5620701, Revision 2: collector, Kizuna 3T busbar XZ assembly replacement	60

GE Healthcare**Renewal Parts and Tooling**

6.5.20	5877116, Revision 5: Kizuna 3T Terminal Cover Kit	60
6.5.21	5483563, Revision 1: PETMR New FG Catalogs	60
6.5.22	5788902, Revision 2: Gradient termination Lug kit.....	60
6.5.23	5861176, Revision 1, DV Gradient Cable Lug Kit - 1 AWG (Pink Lugs)	61
6.5.24	5861177, Revision 1, DV Gradient Cable Lug Kit - 2/0 AWG (Orange Lugs)	61
6.6	FRU Assemblies – VRMw gradient coil	62
6.6.1	5534741, Revision 5: VRMW FRU Gradient Coil Asm.....	62
6.6.2	5473880, Revision 2: HM Magnet, Molded Shim Tray Assembly, VRMw and XRMw.....	62
6.6.3	5473879, Revision 2: UA Magnet, Molded Shim Tray Assembly, VRMw and XRMw.....	62
6.6.4	5542955, Revision 1: Collector, VRMW Coolant Manifold FRU	62
6.6.5	5850819, Revision 1: XRMW and VRMW Shipping Collector, FRU	63
6.7	FRU Assemblies – HRMw gradient coil	64
6.7.1	5768102, Revision 2: Collector, HRMw Gradient Coil FRU.....	64
6.7.2	5863473, Revision 4: COLLECTOR, HRMW GRADIENT BRAID-WOUND HOS FRU	64
6.7.3	5927118, Revision 2: HRMW ATLANTIS FRU + PRO DRIVE.....	64
6.7.4	5768103, Revision 1: Collector, HRMw Gradient Coil Mounting Kit, UA, FRU	65
6.7.5	5768105, Revision 2: HRMw Busbar, UA, FRU	65
6.7.6	5778991, Revision 2: HRMw Busbar, G3, FRU.....	66
6.8	FRU Assemblies – IRMw GRADIENT COIL	67
6.9	FRU Assemblies – Artist EVO GRADIENT COIL.....	67
6.10	Miscellaneous parts for gradients	67
6.10.1	5803053, Revision 1: Nordlock Washer Kit, FRU.....	67
6.10.2	5426726, Revision 2: IBC Temperature Sensor Replacement Kit.....	67
6.10.3	5414896, Revision 3, Hardware FRU Kit, FRMW Gradient Asm	67
6.10.4	2241521, Revision 0: GRADIENT CABLE SPLICE KIT	68
6.10.5	5849225, Revision 1: M12 Gradient Asm Hardware FRU Kit.....	68
6.10.6	5867176, Revision 1: DV Gradient Cable Lug Kit - 1 AWG (Pink Lugs).....	68
6.10.7	5867177, Revision 1: DV Gradient Cable Lug Kit - 2/0 AWG (Orange Lugs)	68
6.11	Tool Kits – Gradient	69
6.11.1	2164744-9, Revision 2: Gradient Insertion Lifting Tool for BRM, TRM, XRMB and XRMw	69
6.11.2	Torque wrenches kits	70
6.11.3	Short Rm, Blue Coil Cart, Cradle	70
7.	Miscellaneous Items	71
7.1	Non-Magnetic/Non-Ferrous Hand Tools	71

GE Healthcare**Renewal Parts and Tooling**

7.1.1	2385097, Revision 0: FE Basic Nonmagnetic Tool Kit	71
7.1.2	46-320273G4, Revision 0: Nonmagnetic Tool Kit	72
7.1.3	46-320718G1, Revision 0: Nonmagnetic Metric Tool Kit	73
7.1.4	46-320719G1, Revision 0: Nonmagnetic Inch Tool Kit	74
7.2	FRU Assemblies, Parts and Kits – Miscellaneous.....	75
7.2.1	46-271137G1, Revision 2: Cryogen Safety Equipment Kit.....	75
7.2.2	Miscellaneous Safety Items	75
7.2.3	46-260888G4, Revision 1: LCC Magnet Leveling Kit	75
7.2.4	5819317, Revision 1: MAGNET LEVELING SHIM KIT (PM and AR magnet series)	75
7.2.5	5926419, Revision 3: Lifting Rail FRU Kit: LCC/450/450W magnet (R, RD, CR, HM)	76
7.2.6	5926420, Revision 3: Lifting Rail FRU Kit: G3/750w magnet (W, WB, UA).....	76
7.2.7	5926421, Revision 3: Lifting Rail FRU Kit: IPM/Starter magnet (PM, AR).....	76
7.2.8	5926709, Revision 3: Hardware Collector Kit, LCC, LCCW & DVW	76
7.2.9	5926710, Revision 3: Hardware Collector Kit, 3TLC and Corona.....	77
7.2.10	5926711, Revision 3: Hardware Collector Kit, PM and AR	77
7.2.11	5366120, Revision 4: Stud Replacement Kit, DVw	77
7.2.12	5339638, Revision 2: IO Angle Bracket Asm.....	78
7.2.13.	46-318565G1, Revision 4: Vacuum Tool Kit [Emergency Vacuum Port Operator]	78
7.2.14.	Miscellaneous Vacuum Maintenance Equipment.....	78
7.3	FRU Parts common – Miscellaneous	79
7.3.1	FRU Parts for R Series Magnets - Miscellaneous.....	79
7.3.2	FRU Parts for HM/UA Series Magnets– Miscellaneous.....	79
7.3.3	FRU Parts for W/WB Series Magnets - Miscellaneous	80
7.3.4	FRU Parts for RD Series Magnets - Miscellaneous	80
7.3.5	5265652, Revision 4: Kit, Field Spares, HDv	80
7.3.6	46-294744G5, Revision 10: Field spare Kit LCC	81
7.4	Helium Vent Kits	82
7.4.1	46-318057G6, Revision 5: LCC Vent Kit Assembly.....	82
7.4.2	46-318057G2, Revision 0: Vent Kit Assembly, Mobile Sites.....	82
7.5	Tool Kits – Miscellaneous	83
7.5.1	2230707-2, Revision 0: Service Platform, Wide Open, Mobile	83
7.5.2	46-294872G3, Revision 2: Sav-Con Kit.....	83
7.5.3	5158913, Revision 1: Latching Lemo Tool Kit	83
7.5.4	5105189, Revision 1: Kit, Instrumentation Connector Reinsertion.....	83

GE Healthcare

Renewal Parts and Tooling

7.5.5	2262769, Revision 1: Plumbing Asm, LCC Replacement	84
7.5.6	5555227, Revision 12: Assembly Plumbing, PED Certified	85
7.5.7	5175534, Revision 3: Plumbing Assembly, 3.0T, SST	86
7.5.8	2293050, Revision 0: GAUSS METER W/PROBE, CABLE, CASE	87
7.6	MCATs – Miscellaneous.....	87
7.6.1	M1060SR, Revision 0: 2.5M Ceiling Ht Kit (M1060CH for China)	87
7.6.2	M7000GM, Revision 2: Ceiling Height Kit-passive.....	87

Chapter 1 – SAFETY and INTRODUCTION

1. Safety and Introduction

1.1. MR Service Safety Manual, Direction 5452735

- 1.1.1. Prior to and while performing any actions contained within this document, one must read, review and follow all guidelines as defined in the MR Service Safety Manual, Direction 5452735. The most current version of the MR Service Safety Manual should be available at all times while performing any actions contained within this document.
- 1.1.2. The latest release of the MR Service Safety Manual can be obtained through SIMS or through your GE HealthCare Field Service Representative.

1.2. Direction Application

- 1.2.1. This document contains information pertaining to Magnet and Gradient related Service Tools and Field Replaceable Units associated with Zero Boil Off Magnets (R, CR, RD, W, WB, HM and UA Series).
- 1.2.2. GPRS Status “Yes” indicates the item is available for order via the FE ordering system. The immediate item availability is not guaranteed by this status. It depends from stocking and supply chain processes.
- 1.2.3. GPRS Status “No” indicates the item is not available for order via the FE ordering system. If item is needed and GPRS Status is “No,” a Priority Purchase Order or upper level unit order would be required.
- 1.2.4. Tools, not available from GPRS, must be purchased using applicable tool provisioning process
- 1.2.5. The parts, marked as “substitute” will be used in the kit once vendor or manufacturing will run out of primary part.

1.3. Documentation Organization

- 1.3.1. Documents referenced within this manual are available through GE HealthCare Customer Documentation Library portal or through your GE Healthcare Field Service Representative.

1.4. Damage in Transportation

- 1.4.1. All packages should be closely examined at time of delivery.
- 1.4.2. If damage is apparent, have notation “damage in shipment” written on all copies of the freight or express bill before delivery is accepted or “signed for” by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 days period.
 - ◆ Call 1–800–548–3366 and use option 8.
 - ◆ Contact your local service coordinator for more information on this process.

GE Healthcare**Renewal Parts and Tooling****Chapter 2 – HANDHELD METERS****2. Handheld Meters****2.1. FRU Assemblies – Handheld Meters****2.1.1. 46-294052G1, Revision 0: Water Flow Meter Kit**

Item No	Part Number	Qty	Description	GPRS
	46-294052G1		WATER FLOWMETER KIT	Yes
1	46-294052P1	1	FLOWMETER, IN-LINE, FLOW RANGE .5-4.0 GPM, MEDIA: WATER/1.0 HEDLAND PART NO. 613-004	No
2	46-294052P2	1	HEX NIPPLE, BRASS, .5 MNPT X .5 MNPT CAJON NO. B-8-HN	No
3	46-294052P3	1	HOSE CONNECTOR, BRASS, .5 MNPT X .5 ID HOSE CAJON NO. B-8-HC-1-8	No
4	46-252065P19		TEFLON TAPE .50 WIDE	Yes

2.2. Tool Kits – Handheld Meters**2.2.1. 46-301477G2, Revision 1: Lakeshore Thermometer Kit**

Item No	Part Number	Qty	Description	GPRS
	46-301477G2		KIT, LAKESHORE THERMOMETER	No
1	46-301453P6	1	BLOW MOLDED, 14 X 10-3/4 X 5-3/4	No
2	46-301618P1	1	TOP - CONVOLUTING PATTERN	No
3	46-301618P2	1	BOTTOM - 2 LAYERS	No
4	46-301478P1	1	TOOL OLD THERMOMETER, DIGITAL, 3W, 90-125VAC FOUR DIGIT LED DISPLAY, (8) CHANNELS 1.4-475K W/O PROBE, 23-473K W/PROBE. SUPPLIED WITH POWER CORD AND 36-PIN ~D~ CONNECTOR. (DWG REV 3/92, JFK - OK).	No
5	46-301619P1	1	INTERCONNECTING TO BALZER DIODES	No
6	46-301620P1	1	INTERFACE TO GE MAGNET	Yes
7	2139962	1	LABEL, LAKESHORE THERM. INTERFAC KIT	No
8	46-193027P782	1	X-RAY RATING PLATE	No
9	2129919	1	CABLE	Yes

2.2.2. 5119507, Revision 2: Cryogenic Temperature Meter [MTM101]

Item No	Part Number	Qty	Description	GPRS
	5119507		CRYOGENIC TEMPERATURE METER (MTM 101)	No

2.2.3. 362622, Revision 0: Portable Helium Meter [LHM 100]

Item No	Part Number	Qty	Description	GPRS
	2362622		PORTABLE HELIUM METER	No

2.2.4. 2270571, Revision 0: Portable Liquid Helium Level Meter with Battery Charger, Cable and Case

Item No	Part Number	Qty	Description	GPRS
	2270571		PORTABLE LIQUID HELIUM LEVEL METER WITH BATTERY CHARGER, CABLE AND CASE. AMI #150A-120-N-18.5-IN-PCT-C WITH CASE	No

GE Healthcare**Renewal Parts and Tooling****2.2.5. 2183710, Revision 0: Magnet Specialty Tool Kit [RuO Temperature Meter]**

Item No	Part Number	Qty	Description	GPRS
	2183710		MAGNET SPECIALTY TOOLS	No
1	2171219	1	MONITOR, RUO TEMPERATURE SI NO. 9650-I-7-7 (RO600GE-10)	No
2	2167169	1	KIT, BATTERY CHARGER	Yes
3	2174214	1	PLATE, PURGE COVER	Yes
4	2183760	1	COVER PLATE AND O-RING	No
5	2183709	1	CASE, FOAM AND LABEL ASSEMBLY	No

2.2.6. 5886219, Revision 1: Magnet Temperature Monitor Kit [LakeShore 218]

Item No	Part Number	Qty	Description	GPRS
	5886219		Magnet Temperature Monitor Kit	No

GE Healthcare**Renewal Parts and Tooling****Chapter 3 –MAGNET RUNDOWN UNIT & MAGNET MONITOR ITEMS****3. Magnet Rundown Unit Items (MRUs) & Magnet Monitor (MagMon) Kits****3.1. Magnet Rundown Unit Items (MRUs)****3.1.1. 5180313, Revision 7: Kit, Magnet Rundown Unit**

Item №	Part Number	Qty	Description	GPRS
	5180313		KIT, MOD MRU2006 PRIMARY MRU	No
1	5196918-2	1	MAGNET RUNDOWN UNIT, PRIMARY, MRU2006	Yes
2	2128699	1	606 MAGNET END TO ERU MRU	Yes
3	5199742	1	PARTS PACK, MRU2006, POWER RESISTOR-FUSE-HOLE PLUG,	No

3.1.2. 5180187, Revision 3: Kit, Remote Magnet Rundown Unit

Item №	Part Number	Qty	Description	GPRS
	5180187		KIT, MOD MRU2006 REMOTE MRU	No
1	5196919	1	REMOTE MAGNET RUNDOWN UNIT, MRU2006	Yes
2	5196920	1	CABLE, MRU2006 PRIMARY TO PEN PANEL	Yes
3	5196921	1	CABLE, MRU2006 REMOTE TO PEN PANEL	Yes
4	5128562	1	LABEL, MRU DANGER	No

3.1.3. 5730424, Revision 2: MRU Tester Assembly

Item №	Part Number	Qty	Description	GPRS
	5730424		MRU Tester Assembly	Yes

3.1.4. 5199490, Revision 5: Battery Pack Asm, MRU2006

Item №	Part Number	Qty	Description	GPRS
	5199490		Battery Pack Asm, MRU2006	Yes

3.1.5. 5363631, Revision 2: Kit, Spare Fuse for MRU

Item №	Part Number	Qty	Description	GPRS
	5363631		Kit, Spare Fuse for MRU	Yes
1	5196919-2	5	Fuse, F 1.25A, 250V, 5x20mm, Littelfuse or Bussman, F3 and F4, mru2006	No
2	5196919-3	5	Fuse, T 2.5A, 250V, 5x20mm, Littelfuse or Bussman, F2, mru2006	No
3	5196919-9	5	Fuse, GMD-400-R, Bussmann, 120V, 400MA, 5x20mm,	No
4	5196919-10	5	Fuse, GMD-200-R, Bussmann, 220VAC, 200MA, 5x20mm,	No

GE Healthcare**Renewal Parts and Tooling****3.2. Magnet Monitor (MagMon) Kits****3.2.1. 5771373, Revision 5: Magnet Monitor 4**

Item №	Part Number	Qty	Description	GPRS
	5771373		Magnet Monitor 4	Yes

3.2.2. 5267626, Revision 3: Collector Kit, Mag Mon

Item №	Part Number	Qty	Description	GPRS
	5267626		COLLECTOR KIT, HDV MAGNET MON 3	No
1	5771373	1	Magnet Monitor 4*	Yes
2	5193795	1	DVMR MAGNET MONITOR CABLE KIT	No
3	2219341	1	RUTHENIUM OXIDE SENSOR AMPLIFIER	Yes
4	5807125	1	CABLE, RUN 850 MON J-2 to COMPRESSOR RS232	Yes

* MagMon4 FRU replaces MagMon3 FRU 2394952 and 2394952-100

3.2.3. 5180301, Revision 1: Remote Alarm Kit, Magnet Monitor 3

Item №	Part Number	Qty	Description	GPRS
	5180301		REMOTE ALARM KIT, MAGMON3	No
1	2266628 or 2266628-100	1	MAGNET MONITOR REMOTE ALARM BOX (2266628) or MAGNET MONITOR REMOTE ALARM BOX - ROHS (2266628-100)	Yes
2	2263160	1	(915) MSM1A2-J1 TO FJ6	Yes
3	2263159	1	(916) MSM1-J7 TO FJ1, FJ5 AND FJ6	No

3.2.4. 2222520-3, Revision 3: Collector Kit, Magnet Monitor 3

Item №	Part Number	Qty	Description	GPRS
	2222520-3		MAGMON3 COLLECTOR KIT, FIXED, TRANSPORTABLE, AND RELOCATABLE	No
1	2394952 or 2394952-100	1	MAGNET MONITOR (2394952) or MAGNET MONITOR 3 – ROHS (2394952-100)	No*
2	2350809-2	1	CORE CABLE KIT, MAGMON3	No
3	2219341	1	RUTHENIUM OXIDE SENSOR AMPLIFIER	Yes

* Replaced with 5771373

3.2.5. 2222520-4, Revision 7: MagMon Mobile Collector Kit

Item №	Part Number	Qty	Description	GPRS
	2222520-4		MagMon Mobile Collector Kit	No
1	5771373	1	Magnet Monitor 4	Yes
2	2350809-105	1	KIT, HD MAGNET CORE CABLE, MAG MON 3 - ROHS	No
3	2219341	1	RUTHENIUM OXIDE SENSOR AMPLIFIER	Yes
5	2276094	1	UNINTERRUPTIBLE POWER SUPPLY	Yes
6	5180301	1	REMOTE ALARM KIT, MAGMON3	No
7	5224349	1	ENCLOSURE CABLE KIT	No
8	5807125	1	CABLE, RUN 850 MON J-2 to COMPRESSOR RS232	Yes

GE Healthcare**Renewal Parts and Tooling****3.2.6. 5193795, Revision 6: Kit, Magnet Monitor Cables, DVMR**

Item №	Part Number	Qty	Description	GPRS
	5193795		DVMR MAGNET MONITOR CABLE KIT	No
1	5357236	1	RUN M3023 SPW-J117 TO MAG-FJ2,14000MM	No
2	5351557	1	RUN E3023 SPW-J117 TO MON-J8,17000MM	No
3	5351627	1	RUN E3022,MON-J7 TO SPW-J116, INSTRUMENTATION BOX	No
4	5352885	1	RUN E3024, MON-J9 TO MON-FJ4 AND MON-FJ3, FLOW MONITORING	No
5	5352519	1	RUN E3026, HEC-FL-TEMP TO MON-FJ3, COMPRESSOR INTERFACE	No
6	5351365	1	RUN E3020 PEN-J13 TO MON-J10,17000MM	No
7	5351806	1	RUN E3025, CRY TO MON-FJ4, FLOW MONITOR-SUMITOMO	No
8	5351746	1	RUN M3022, SPW-J116 TO MAG-P403, 14000 MM	No
9	5351608	1	RUN E3034 MON-J5 TO NETWORK, 30000MM	No
10	5115594-3	1	POWER CORD IEC 60320 C13 TO NEMA 5-15P 72 INCHES LONG	No
11	5352964	1	CABLE, MAGNET ROOM INTERFACE	No

3.2.7. 5548703, Revision 1: Kizuna Magnet Monitor Cable Kit

Item №	Part Number	Qty	Description	GPRS
	5548703		Kizuna Magnet Monitor Cable Kit	No
1	5351557	1	Run E3023 SPW-J117 to MON-J8,17000mm	No
2	5351627	1	Run E3022,MON-J7 to SPW-J116, Instrumentation Box	No
3	5508124	1	Run 826 FJ3 TO MSM1-A1-J9	No
4	5508128	1	Run 827 FJ3 TO FJ4 and MS5A5	Yes
5	5508127	1	Run 833 FJ4 TO MS5-A1-A6-JR	Yes
6	2333825	1	WATER FLOW TURBINE_TEMP SENSOR ASM, SS BARB, SS PIPE, SS COVER	Yes
7	46-265067P1	6	I/F PANEL CABLE SCREWLOCK CONNECTOR	Yes
8	5351365	1	Run E3020 Pen-J13 to Mon-J10,17000mm	No
9	5351608	1	Run E3034 MON-J5 to Network, 30000mm	No
10	5115594-3	1	POWER CORD IEC 60320 C13 to NEMA 5-15P 72 inches long	No

3.2.8. 5782696, Revision 2: RIO Mag Mon Cable Collector

Item №	Part Number	Qty	Description	GPRS
	5782696		RIO Mag Mon Cable Collector	No
1	5357236	1	Run M3023 SPW-J117 to Mag-FJ2,14000mm	No
2	5351557	1	Run E3023 SPW-J117 to MON-J8,17000mm	No
3	5351627	1	Run E3022,MON-J7 to SPW-J116, Instrumentation Box	No
4	5508124	1	Run 826 FJ3 TO MSM1-A1-J9	No
5	5508128	1	Run 827 FJ3 TO FJ4 and MS5A5	Yes
6	5351365-2	1	RUN E3020 PEN-J13 OR ISC IO J13	No
7	5508127	1	Run 833 FJ4 TO MS5-A1-A6-JR	Yes
8	5351746	1	Run M3022, SPW-J116 to MAG-P403, 14000 mm	No
9	5351608	1	Run E3034 MON-J5 to Network, 30000mm	No
10	5115594-3	1	POWER CORD IEC 60320 C13 to NEMA 5-15P 72 inches long	No
11	5352964	1	CABLE, MAGNET ROOM INTERFACE	No

GE Healthcare

Renewal Parts and Tooling

3.2.9. 2350809, Revision 3: Core Kit, Magnet Monitor 3

Item №	Part Number	Qty	Description	GPRS
	2350809		CORE KIT, MAGMON3	Yes
1	46-317359G975	1	RUN 824TR PP1-J10 TO FJ1	Yes
2	46-328000G976	1	RUN 829 MS1-FJ2 TO PP1-J48	No
3	46-317359G977	1	RUN 825TR PP1-J48 TO MSM1-A1-J8	No
4	46-328000G980	1	RUN 831 FJ1 TO MSM1-A1-J7	Yes
5	46-328000G978	1	RUN 826 FJ3 TO MSM1-A1-J9	No
6	46-328578P1	1	RUN 827 FJ3 TO FJ4 AND MS5A5	Yes
7	46-328000G979	1	RUN 823 MSM1-J10 TO MR2-A11-J24	No
9	46-328000G981	1	RUN 833 FJ4 TO MS5-A1-A6-JR	Yes
10	2333825	1	WATER FLOW (TURBINE)/TEMP SENSOR ASM	Yes
11	46-265067P1	10	I/F PANEL CABLE SCREWLOCK CONNECTOR	Yes
12	46-328000G974	1	RUN 828 MS1-A3-A1-J403 TO PP1-J10	Yes
13	2222797	1	830 CABLE, MAGNET ROOM INTERFACE	Yes
14	46-296816P1	1	ADAPTER, BNC TEE, 50 OHM AMPHENOL #31-200 - ALL CENTER CONDUCTORS ARE SILVER PLATED. (REV TK, 11/91). HAS REXOLITE INSULATION	Yes
15	2276568	1	(940) MSM3-RS-232 TO MSM1-J1	Yes

3.2.10. 2350809-3, Revision 6: Magnet Monitor Non-LCC Upgrade HW Kit

Item №	Part Number	Qty	Description	GPRS
	2350809-3		Magnet Monitor Non-LCC Upgrade HW Kit	Yes
1	46-317359G975	1	RUN 824TR PP1-J10 TO FJ1	Yes
2	46-328000G976	1	RUN 829 MS1-FJ2 TO PP1-J48	No
3	46-317359G977	1	RUN 825TR PP1-J48 TO MSM1-A1-J8	No
4	46-328000G980	1	RUN 831 FJ1 TO MSM1-A1-J7	Yes
5	46-328000G978	1	RUN 826 FJ3 TO MSM1-A1-J9	No
6	46-328578P1	1	RUN 827 FJ3 TO FJ4 AND MS5A5	Yes
7	46-328000G979	1	RUN 823 MSM1-J10 TO MR2-A11-J24	No
8	2333825	1	WATER FLOW (TURBINE)/TEMP SENSOR ASM	Yes
9	46-328000G974	1	RUN 828 MS1-A3-A1-J403 TO PP1-J10	Yes
9	2108726	1	MAG ROOM INTERFACE	Yes
10	2109440	2	SWAGLOK PORT CONNECTOR	No
11	2204485	1	CABLE, SYSCAB INTERFACE	Yes
12	5192503	1	Pressure Transmitter	Yes
13	2369032	1	SHIELD, PRESSURE TRANSDUCER RF	No
14	46-252294P1	1	1/4 TEE	Yes
15	46-260793P11	2	SCREW,#4-40 X .25LG, ROUND HD, BRASS	No
16	46-294167P98	1	Nipple, Pipe, Brass, .125 NPT Sch.40 x 1.50 Long	No
17	5180922	1	0.25 inch OD Union Elbow	No
18	46-252557P2	1	.25 NPT, BRASS	Yes

3.2.11. 5113506, Revision 2: Rebuild Kit for Water Flow Meter 2333825-2.

Item №	Part Number	Qty	Description	GPRS
	5113506		Rebuild Kit for Water Flow Meter 2333825-2	Yes

GE Healthcare

Renewal Parts and Tooling

Chapter 4 – CRYOGENICS

4. Cryogenics

Coldhead and Compressor

4.1. FRU Assemblies – Coldhead and Compressor

4.1.1. 2188440-2, Revision 4: Compressor, Water Cooled, Sumitomo CSW-71D

Item №	Part Number	Qty	Description	GPRS
	2188440-2		COMPRESSOR, SHIELD COOLER, SUMITOMO CSW-71D, HV WATER COOLED, 460V,60HZ 380V,50HZ, THREE PHASE*	No

* Replaced with 5314887

4.1.2. 5314887, Revision 2: Compressor, Water Cooled, Sumitomo F-50H

Item №	Part Number	Qty	Description	GPRS
	5314887		SERVICE SPARE, SUMITOMO F-50H COMPRESSOR	Yes
1	5197842	1	COMPRESSOR, WATER COOLED, AC 380-415V AT 50HZ OR AC 460-480V AT 60HZ, 3 PHASE, SUMITOMO F-50H	No
2	5314775	1	CRATE FOR SUMITOMO	No

4.1.3. 2166592, Revision 3: Compressor, Air Cooled, Sumitomo CSA-71

Item №	Part Number	Qty	Description	GPRS
	2166592 or 5497024		COMPRESSOR, SUMITOMO, CSA-71	Yes

4.1.4. 5428192, Revision 3: Kit Compressor F-70H, SHIG #267446A

Item №	Part Number	Qty	Description	GPRS
	5428192		Kit, Compressor F-70H, HV, SHIG #267446A	Yes
1	5428182	1	Compressor, F-70H, SHI-APD, High Volt Variant, Water-Cooled, UL 480V #8 Gas Couplings, SHIG #267446A	Yes
2	5428185	1	Power Cord, HV, 90 Degree Down, Side Entry, NEMA L 16-20P, Plug SHIG #267195C10GE	No
3	5428187	1	Cable, Cold HD, F-70H Compressor to RF Filter, 6MM Terminal, 15M LG, SHIG #2677447C50GE	No
4	5428182-100	1	Compressor Manual, F70H Sumitomo, SHIG @267449A	No
5	5199282	1	Cable, Sumitomo HC-10 Compressor, US Cable	Yes
6	5449002	1	Instruction, F-70H Compressor Installation	No

4.1.5. 5812867, Revision 1: F-50SH compressor service FRU

Item №	Part Number	Qty	Description	GPRS
	5812867		F-50SH compressor service FRU	Yes
1	5772873	1	Sumitomo F50SH High-Voltage Water-Cooled Helium Compressor	No
2	5314775	1	Crate for Sumitomo	No

4.1.6. 5368565, Revision 1: Service Spares, Sumitomo F-50LV Compressor

Item №	Part Number	Qty	Description	GPRS
	5368565		Service Spares, Sumitomo F-50LV Compressor*	Yes
1	5366042	1	Compressor, Water Cooled, Sumitomo F-50, Low Voltage	No
2	5314775	1	Crate for Sumitomo	No

*Replaces obsolete 2188184 - COMPRESSOR, SUMITOMO, WATERCOOLED (CSW-71C).

4.1.7. 5840878, Revision 2: Sumitomo FA-50SH High Voltage Air-cooled Helium Compressor

Item №	Part Number	Qty	Description	GPRS
	5840878		Sumitomo FA-50SH High Voltage Air-cooled Helium Compressor	Yes

*Replacement CSA-71A with FA-50SH requires also replacement input power cable with 5849910.

GE Healthcare**Renewal Parts and Tooling****4.1.8. 5840879, Revision 2: Sumitomo FA-50SL Low Voltage Air-cooled Helium Compressor**

Item №	Part Number	Qty	Description	GPRS
	5840879		Sumitomo FA-50SL Low Voltage Air-cooled Helium Compressor	Yes

*Replacement CSA-71A with FA-50SL requires also replacement input power cable with 5849911.

4.1.9. 2384329, Revision 3: Compressor, Outdoor, CNA-61D, Sumitomo

Item №	Part Number	Qty	Description	GPRS
	2384329	1	Compressor, Outdoor, CNA-61D, Sumitomo	Yes

4.1.10. 2384329-2, Revision 2: Compressor, Indoor, CNA-61D, Sumitomo

Item №	Part Number	Qty	Description	GPRS
	2384329-2	1	Compressor, Indoor, CNA-61D, Sumitomo	Yes

4.1.11. 5491924, Revision 1: Compressor, Outdoor, CAN-61C3-C, Sumitomo

Item №	Part Number	Qty	Description	GPRS
	5491924		Compressor, Outdoor, CAN-61C3-C, Sumitomo	Yes

*Can be used in place of 5160751, however, NOT COMPATIBLE with indoor unit 5160751-2 until FMI 67817 installed. No FMI required if used with indoor unit 5491924-2.

4.1.12. 5491924-2, Revision 1: Compressor, Indoor, CNA-61C3-E, Sumitomo

Item №	Part Number	Qty	Description	GPRS
	5491924-2		Compressor, Indoor, CNA-61C3-E, Sumitomo	Yes

4.1.13. 5438959, Revision 2: Kit, A3 FRU, Sumitomo Coldhead and Gasket

Item №	Part Number	Qty	Description	GPRS
	5438959		A3 FRU KIT, SUMITOMO COLDHEAD AND GASKET	Yes
1	5445412	1	Cold Head Asm, A3_R	No
2	2171620	1	KIT, CRYOCOOLER REPLACEMENT GASKETS	Yes

4.1.14. 5814003, Revision 1: Kit A4 FRU, Sumitomo Coldhead and Gasket

Item №	Part Number	Qty	Description	GPRS
	5814003		A4 FRU KIT, SUMITOMO COLDHEAD AND GASKET	Yes
1	5815457	1	Assembly, Sumitomo Coldhead RDE 412A4-RE	No
1*	5720168	1	Assembly, Sumitomo Cold Head RDE 412A4	No
2	2171620	1	KIT, CRYOCOOLER REPLACEMENT GASKETS	Yes

*The coldhead 5720168 is a substitute for 5815457 within 5814003 kit

4.1.15. 5814003-R, Revision 1: A4-R FRU Kit, Sumitomo Coldhead and Gasket

Item №	Part Number	Qty	Description	GPRS
	5814003-R		A4-R FRU KIT, SUMITOMO COLDHEAD AND GASKET	Yes
1	5812683	1	Assembly, Sumitomo Cold Head RDE-412A4, R	No
1*	5812685	1	Assembly, Sumitomo Cold Head A3 to A4 Upgrade, R	No
1*	5815457	1	Assembly, Sumitomo Coldhead RDE 412A4-RE	No
2	2171620	1	KIT, CRYOCOOLER REPLACEMENT GASKETS	Yes

*The coldheads 5815457 and 5812685 are substitute for 5812683 within 5814003-R (5814003-RA) kit

GE Healthcare

Renewal Parts and Tooling

4.1.16. 2171620, Revision 11: Kit, Cryocooler Replacement Gaskets

Item No	Part Number	Qty	Description	GPRS
	2171620		KIT, CRYOCOOLER REPLACEMENT GASKETS	Yes
1	46-281247P1	1	4.987 ID, X .103 THK, PARKER#2-159	Yes
2	46-281241P1	1	INDIUM GASKET, 3.37 OD	Yes
3	2171621	1	KIT, INDIUM GASKET	Yes
4	2188185	1	FEELEER GAUGE, .003 THK.	Yes
5	2188185-2	1	FEELEER GAUGE, .005 THK.	Yes
6	46-252065P200	1	Gloves, ambidextrous cotton or cotton blend, no coating, color: white, size: Mens large. Magid Jersey Master pn T2705 or equivalent.	No
7	2290272	1	KIT,CHARCOAL PACKET	Yes
8	2391747	1	FILTER,DISPOSABLE IN-LINE DESSICCANT	Yes
9	2173358	48	SPRING, DISC, 12.5 OD 6.2 ID .70 THK TYPE 302 SST.	Yes
10	46-252065P141	2	AEROQUIP O-RING, FLAT, RUBBER	Yes
11	46-252065P157	1	AEROQUIP GASKET, CATALOG PART NUMBER 22008-4, .25 IN SIZE, FLAT, RUBBER	No
12	5159706	2	SCREW HEX HD STST M6X45L	No

4.1.17. 46-318696G1, Revision 0: Water Tee Assembly

Item No	Part Number	Qty	Description	GPRS
	46-318696G1		Water Tee Assembly	Yes
1	46-294167P9	1	TUBING CLEAR PVC, 1/2 ID X 27/32 OD 1 FT. LONG WORKING PRESSURE 150 PSI @70 DEGREES. MCMaster-CARR CATALOG PART NUMBER 5238K758 OR EQUAL.	No
2	46-294167P10	1	VALVE PVC BALL, 1/2 PIPE SIZE, MAX PRESSURE W.O.G., 150 PSI. MCMaster- CARR CATALOG PART NUMBER 4953K12 OR EQUAL.	No
3	46-294167P11	1	FITTING TEE MALE BRANCH, 1/2 ID X 1/2 ID X 1/2 NPT NYLON HOSE FITTING. MCMaster-CARR CATALOG PART NUMBER 5375K39 OR EQUAL.	No
4	46-294167P12	1	FITTING, NIPPLE MALE, 1/2 ID X 1/2 NPT NYLON HOSE FITTING. MCMaster- CARR CATALOG PART NUMBER 5372K25 OR EQUAL.	No
5	46-294167P13	3	CLAMP HOSE, 1/2 X 1 1/4 X 9/16 WIDE SWIVEL LOCK WORM DRIVE. MCMaster- CARR CATALOG PART NUMBER 5322K14 OR EQUAL.	No
6	46-294167P16	1	QUICKSNAP EXPANDABLE PLUG FOR 1/2~ I.D. TUBE.	No

4.1.18. 2194003-2, Revision 1: Kit, A2 FRU, Sumitomo Coldhead and Gasket

Item No	Part Number	Qty	Description	GPRS
	2194003-2		KIT, SUMITOMO COLDHEAD AND GASKET*	No*
1	2218465-3	1	COLDHEAD ASM,HIGH EFFNCY, SUMITOMO	No
2	2171620	1	KIT, CRYOCOOLER REPLACEMENT GASKETS	Yes

* This kit is not available any longer, listed for reference only. Replaced with 5814003

GE Healthcare

Renewal Parts and Tooling

4.2. FRU Parts – Coldhead and Compressor

Item №	Part Number	Qty	Description	GPRS
	5264644		Adsorber, Sumitomo F-50 ^{a)}	Yes
	5836219		Temperature Sensor Replacement Kit, Sumitomo F-50H/L ^{h)}	Yes
	5823886		Fuse Kit, Sumitomo F-50SH Compressor	Yes
	2154502-2		621 GAS LINE, FLEXIBLE SUPPLY 65.6 FT. OR 20 M LONG	Yes
	2154502		GAS LINE, FLEXIBLE (SUPPLY) 19.7 FT. OR 6 M LONG	Yes
	2154505-2		622 GAS LINE, FLEXIBLE RETURN 65.6 FT. OR 20 M LONG	Yes
	2154505		GAS LINE, FLEXIBLE (RETURN) 19.7 FT. OR 6 M LONG	Yes
	2273967		ADAPTER, ELBOW, SUMITOMO FLEX LINE (Supply/Yellow) ^{b)}	Yes
	2273967-2		ADAPTER, ELBOW, SUMITOMO FLEX LINE (Return/Green) ^{b)}	Yes
	5480609-4		Compressor Quiet Line Kit 12m ^{c)}	Yes
	5481535		Kit, Quiet Extension Line ^{d)}	Yes
	5412653		adapter, elbow, sumitomo flex line (Supply/Yellow) ^{e)}	Yes
	5412653-2		adapter, elbow, sumitomo flex line (Return/Green) ^{e)}	Yes
	2259837		CABLE ASM, COLDHEAD CONTROL, LCC ^{f)}	Yes
	2200835		CABLE, INPUT POWER, SUMITOMO CSW-71D HIGH VOLTAGE, WATER COOLED COMPRESS.	Yes
	2274831		CABLE, INPUT POWER, 5M LG, SUMITOMO CSW-71D HV WATER COOLED COMPRESSOR	Yes
	5145460-38		FUSE, COLDSLEEVE HEATER TOOL	Yes
	5440278		150 W CARTRIDGE HEATER	Yes
	2251611		RACEWAY, FLEXIBLE CORRUGATED, ORANGE	Yes
	5480609		COMPRESSOR QUIET LINE KIT ^{g)}	Yes
	2191112-13		Fuse, 2A, SOC ST4N1-2A, SHI part RE71WT1679	Yes
	2191112-3		3A CLASS G FUSE, SUMITOMO CRYOCOOLER SUMITOMO P/N RE71WT0601	Yes
	2191112-2		2A CLASS G FUSE, SUMITOMO CRYOCOOLER SUMITOMO P/N RE71WT0602	Yes
	5507283		Water Line, Rotating Adaptor 90	Yes
	5819315		HOSE NIPPLE FRU KIT	Yes
	5335002		Flow Control Adapter, NPT	Yes
	46-260802P1		.50 WIDE ALLOY #1145-0 SOFT ALUMINUM FOIL TAPE	Yes
	46-260802P6		TAPE, ALUMINUM FOIL TAPE, .125 WIDE	Yes
	46-252282P1		TAPE, PERMACEL TYPE P-212, .006 X.50 GLASS CLOTH TAPE WITH SILSCON ADHESIVE	Yes
	5849910		CABLE, INPUT POWER, SUMITOMO FA-50SH, HIGH VOLTAGE, AIR COOLED COMPRESSOR.	Yes
	5849911		CABLE, INPUT POWER, SUMITOMO FA-50SL, LOW VOLTAGE, AIR COOLED COMPRESSOR.	Yes
	5761761		F50 COMPRESSOR REMOTE PORT ADAPTOR	Yes
	2155316		CABLE, COMPRESSOR TO PEN PANEL, OBO 19.7 FT OR 6 M LONG (10K compressors only)	Yes
	2172239		CABLE, COLDHEAD TO PEN PANEL, OBO 19.69 FT OR 6.00 M LONG	Yes
	2224726		METRIC SCREW, HEX HD, SST, PART THRD	Yes
	5878762	1	SPARE KIT, TEMPERATURE STRIP (Pack of 10)	Yes
	46-306717G1	1	NON-MAGNETIC GAS CYLINDER CART WITH LABEL	Yes
	46-258150P1	1	NON-MAGNETIC GAS CYLINDER CART	Yes

a) Adsorber 5264644 is compatible with CSW-71C, CSW-71D, CSA-71, F-50H, F-50L, F-50SH, FA-50SH, FA-50SL compressors

b) Long radius smooth 90-degree flexline adapters

c) 5480609-4 consists of 12 m Supply and Return gas lines

d) 5481535 consists of 6 m Supply and Return M-F extension gas lines

e) Short radius sharp 90-degree flexline adapters

f) 2259837 contains RUN 623 p/n 2155316-2 and RUN 624 p/n 2218292

g) 5480609 consists of 20 m Supply (Run 621) and Return (Run 622) Quiet flex lines

h) Contains 60C sensor

GE Healthcare**Renewal Parts and Tooling****4.3. Tool Kits – Coldhead and Compressor****4.3.1. 46-281088G4, Revision 1: Coldhead/Compressor Maintenance Kit**

Item №	Part Number	Qty	Description	GPRS
	46-281088G4		COLDHEAD/COMPRESSOR MAINTENANCE KIT.	No
2	46-294150P16	2	1 5/8 IN. OPEN END WRENCH LEYBOLD P/N 722-81-021	Yes
3	46-294150P15	1	1 3/8 IN. OPEN END WRENCH LEYBOLD P/N 722-81-020	Yes
4	46-294150P14	2	1 3/16 IN. OPEN END WRENCH LEYBOLD P/N 722-81-019	Yes
5	46-294150P13	1	1 IN. OPEN END WRENCH LEYBOLD P/N 722-81-018	Yes
8	46-294003P1	1	OIL CHARGING HOSE	Yes
9	46-294009P1	1	HIGH PRESSURE, HELIUM	Yes
11	46-294002P1	1	CHARGING - AEROQUIP PART NO. 2807-4 OR EQUIVALENT	Yes
12	46-294006G1	1	#12 AEROQUIP DISCHARGE TOOL ASSY. PART #5401-S14-12.	Yes
13	46-294007G1	1	ADAPTER ASSEMBLY BACKFILL/PUMPDOWN	Yes
14	46-294010G1	1	PURGE COVER	Yes
15	46-260267P3	1	SEAL OFF VALVE OPERATOR WITH KF16 FLANGE	Yes
18	46-294000G1	1	#12 AEROQUIP MALE HALF, AEROQUIP PART NO. 5401-S17-12. INC TUBE ADAPT	Yes
19	46-281999G1	1	#8 AEROQUIP MALE HALF, AEROQUIP PART NO. 5401-S17-8. INC. TUBE ADAPT	Yes
20	46-252065P63	1	WELDERS~ BRUSH, MCMASTER-CAR CAT. NO. 7092T14.	Yes
21	46-294005G1	1	#4 AEROQUIP	Yes
24	46-294026P1	1	QF16 ALUMINUM HINGED, SARGENT WELCH PART NO. AC51100NW016 OR EQUIVALENT LEYBOLD #18341 AND MDC #701000	Yes
25	46-294030G1	1	QF16 CENTERING, SARGENT WELCH PART NO. AC51200NW016 OR EQUIVALENT LEYBOLD PART #88346 AND MDC #710000	Yes
26	46-294939G50	1	CASE, FOAM AND LABEL FOR SHIELD COOLER TOOL KIT.	No
27	46-294937P1	1	#4 MALE COUPLING WITH TUBE ADAPTER AEROQUIP PART #5401-S17-4.	Yes
28	46-294936P1	1	#8 FEMALE COUPLING WITH TUBE ADAPTER AEROQUIP PART #5401-S14-8.	Yes
29	46-317904P1	2	GAS DISCHARGE (PURGE) TOOL #8 (1/2~) AEROQUIP MALE COUPLING TO 1/8~ MALE ADAPTER WITH VALVE AND 0-400 PSI GAUGE.	No
30	46-317904P2	2	GAS DISCHARGE (PURGE) TOOL #12 (3/4~) AEROQUIP MALE COUPLING TO 1/8~ MALE ADAPTER WITH VALVE AND 0-400 PSE GAUGE	No
31	46-281989G1	1	COLD HEAD PORT COVER	No
32	46-318911P1	1	BALZERS COMPRESSOR TO COLDHEAD POWER CABLE.	No
33	46-318910P1	1	LEYBOLD COMPRESSOR TO COLDHEAD POWER CABLE.	Yes
34	2100316	1	LARGE ADAPTER FITTING THAT ALLOWS CONNECTION OF AEROQUIP FITTING TO HELIUM LINE.	No
35	2100317	1	SMALL ADAPTER FITTING THAT ALLOWS CONNECTION OF AEROQUIP FITTING TO HELIUM LINE.	No
37	2102181	1	SHIELD COOLER TOOL KIT LABEL UPGRADE	No
38	2132124	1	COLDHEAD CHANGEOUT TOOL	No
39	2157411	1	#8 AEROQUIP ADAPTER FITTING	Yes

GE Healthcare**Renewal Parts and Tooling****4.3.2. 2300306-3, Revision 1: TOOL KIT, UNIVERSAL COLDHEAD REPLACE WITH INDIUM REMOVAL TOOL**

Item No	Part Number	Qty	Description	GPRS
	2300306-3		TOOL KIT, UNIVERSAL COLDHEAD REPLACE WITH INDIUM REMOVAL TOOL	No
1	2306020-2	1	CASE AND LABEL, C.H. REPLACEMENT KIT WITH INDIUM REMOVAL TOOL	No
2	2287707-2	1	ADAPTER, COLDHEAD SERVICE,	No
3	2267197	1	COLDHEAD SLEEVE PUMPOUT LINE	Yes
4	2287716	1	INSTRUMENT, DIGITAL VACUUM	Yes
5	46-294030G1	7	QF16 CENTERING, SARGENT WELCH PART NO. AC51200NW016 OR EQUIVALENT LEYBOLD PART #88346 AND MDC #710000	Yes
6	46-294026P1	7	QF16 ALUMINUM HINGED, SARGENT WELCH PART NO. AC51100NW016 OR EQUIVALENT LEYBOLD #18341 AND MDC #701000	Yes
7	46-294028P1	2	BLANK-OFF, SARGENT WELCH PART NO. AC51050NW016 OR EQUIVALENT LEYBOLD #88346 AND MDC #712000	Yes
8	46-260267P3	1	SEAL OFF VALVE OPERATOR WITH KF16 FLANGE	Yes
9	5180672	1	FIELD SERVICE KIT, COLDSLEEVE HEATER	No
10	2303728	1	LEXAN COVER PLATE, FSDE	Yes
11	2303724	1	NITROGEN GAS LINE, FSDE	Yes
12	2303731	1	REGULATOR / FLOWMETER WITH QUICK-CON	Yes
13	46-320472P5	1	OPEN END 8MM X 10MM	Yes
14	2236606	1	PLUG,DOUBLE O-RING,BRASS	Yes
15	46-252065P24	1	VACUUM GREASE,DOW CORNING	Yes
16	46-294041P5	1	OIL, QUART SIZE, DUOSEAL SWSII, SARGENT-WELCH PART NO. 1407K-11 (WHEN ORDERING SPECIFY QUART SIZE AND QUANTITY)	Yes
18	46-260342P2	1	.612 I.D. X .818 O.D. ~O~RING	Yes
19	2321289	1	FLEX LINE COUPLING, COLDHEAD KIT	Yes
20	2365037	1	TOOL, INDIUM REMOVAL AND CASE	Yes
21	5111335	1	GAS HOSE LINE, .25 MNPT TO MALE QUICK CONNECT STEM, 25 FT	Yes
22	5111361	1	HELIUM VENT LINE, .25 FEMALE TO MALE QUICK CONNECT STEM, 5 FT	Yes
23	2305056PSP	1	Vacuum pump	No

4.3.3. 5180672, Revision 3: Coldsleeve Heater Kit

Item No	Part Number	Qty	Description	GPRS
	5180672		FIELD SERVICE KIT, COLDSLEEVE HEATER	No
2	5180672-2	1	TOOL BOX, FIELD SERVICE KIT, COLDSLEEVE HEATER	No
3	5180672-3	2	SPRING CLAMP, FIELD SERVICE KIT, COLDSLEEVE HEATER	No
4	5145460	1	TOOL, COLDSLEEVE HEATER	No
5	5179642	1	HOSE, COLDSLEEVE HEATER SHIM LEAD CONNECTOR	No
6	5440278	Available to Order	150 W CARTRIDGE HEATER	Yes
7	5145460-38	1	FUSE, COLDSLEEVE HEATER TOOL	Yes

4.3.4. 5179792, Revision 1: Cable, CE Approved Coldsleeve Heater Extension

Item No	Part Number	Qty	Description	GPRS
	5179792		CABLE, CE APPROVED COLDSLEEVE HEATER EXTENSION	No

GE Healthcare**Renewal Parts and Tooling****4.3.5. 2232674, Revision 0: Fixed Site Coldhead Change Out Kit**

Item №	Part Number	Qty	Description	GPRS
	2232674		COLDHEAD INSERTION AND PURGE COVER TOOL KIT ASSEMBLY	No
1	2232677	1	CASE FOR COLDHEAD INSERTION AND PURGE COVER TOOL KIT - PELICAN MODEL NUMBER 1600	No
2	2232673	1	COLDHEAD INSERTION AND PURGE COVER TOOLS KIT STICK-ON LABEL	No
3	2228876	1	FIXED SITE SHI COLDHEAD INSERTION TOOL	Yes
4	2232560	1	LEXAN PURGE COVER PLATE ASSEMBLY	Yes

4.3.6. 2380982, Revision 0: Kit, Coldhead Insertion/Removal Tool

Item №	Part Number	Qty	Description	GPRS
	2380982		KIT, COLDHEAD INSERTION/REMOVAL TOOL	No
1	2262749	1	CASE FOR COLDHEAD INSERTION / REMOVAL TOOL KIT	No
2	2262685	1	COLDHEAD INSERTION/REMOVAL TOOL KIT STICK-ON LABEL	No
3	2379887	1	INSERTION FIX., COLDHEAD, HFO	No
4	2254615	1	COLDHEAD HOIST	No
5	2275776	1	PURGER COVER PLATE ASSEMBLY	Yes
6	2279643	0	TOOL, REI ASCENDER, NMR - INCLUDED ON BOM 2275779 AND 2262686 AS A RENEWAL PART REFERENCE ONLY.	No
7	2279647	0	TOOL, HARKEN PULLEY, NMR - INCLUDED ON BOM 2275779 AND 2262686 AS A RENEWAL PART REFERENCE ONLY.	No
8	2270738	0	NYLON ROPE, .50~ DIA. X 25~ LONG INCLUDED ON BOM 2275779 AND 2262686 AS A RENEWAL PART REFERENCE ONLY.	No
9	2279644	0	TOOL, LOCKING PIN, NMR - INCLUDED ON BOM 2262686 AS A RENEWAL PART REFERENCE ONLY.	No
10	2275600	0	CAM-MATIC, HARKEN - INCLUDED ON BOM 2262686 AS A RENEWAL PART REFERENCE ONLY.	No

4.3.7. 2380985, Revision 0: LCC Mobile Coldhead Tool Kit

Item №	Part Number	Qty	Description	GPRS
	2380985		KIT, LCC MOBILE COLDHEAD TOOL KIT	No
1	2275777	1	CASE FOR LCC MOBILE COLDHEAD TOOL	No
2	2275778	1	LABEL FOR LCC MOBILE COLDHEAD TOOL KIT CASE	No
3	2379888	1	INSERTION FIX.,COLDHEAD,MOBILE MAG	No
4	2258802	1	MOBILE COLDHEAD TOOL	No
5	2275776	1	PURGER COVER PLATE ASSEMBLY	Yes
6	2279647	0	TOOL, HARKEN PULLEY, NMR - INCLUDED ON BOM 2275779 AND 2262686 AS A RENEWAL PART REFERENCE ONLY.	No
7	2279646	0	TOOL, REI CARABINER, NMR - INCLUDED ON BOM 2275779 AS A RENEWAL PART REFERENCE ONLY.	No
8	2279643	0	TOOL, REI ASCENDER, NMR - INCLUDED ON BOM 2275779 AND 2262686 AS A RENEWAL PART REFERENCE ONLY.	No
9	2270738	0	NYLON ROPE, .50~ DIA. X 25~ LONG INCLUDED ON BOM 2275779 AND 2262686 AS A RENEWAL PART REFERENCE ONLY.	No

GE Healthcare

Renewal Parts and Tooling

4.3.8. 5395677, Revision 1: Kit, Coldhead Alignment Tool Clamp

Item No	Part Number	Qty	Description	GPRS
	5395677		KIT, COLDHEAD ALIGNMENT TOOL CLAMP	No
1	2381228-15	2	SCREW, SS HEX SOCKET CAP, M6 X 70	No
2	5116429	1	BLOCK, STRAP BANDING	No
3	5391106	1	T-BOLT BAND CLAMP	No
4	2390368-4	2	SCREW, SS BUTT HD CAP, .25-20 X .75	No
5	46-252635P7**	2	.25 NOM., PLAIN WASHER	Yes
6	2109878-3*	6	WASHER PLAIN - NORMAL 6.4 MM 12 MM	No

*Replaced with the part 2109879-2

**Plain washer, used on 8 coldhead mounting bolts

4.3.9. 2385069, Revision 2: 3.0T Coldhead Transport Kit

Item No	Part Number	Qty	Description	GPRS
	2385069		KIT, 3.0T COLDHEAD TRANSPORT	No
1	2385070	1	CASE, 3.0T COLDHEAD TRANSPORT KIT	No
2	2385071	1	INSTRUCTIONS, PACKING, 3.0T COLDHEAD TRANSPORT KIT	No
3	2348452	1	TOOL, COLDHEAD INSERTION/REMOVAL	No
4	2375543	1	TOOL, COLDHEAD TRANSPORT, G3	No
5	2384270	1	WRENCH, T-HANDLE,10MM, TITANIUM	No
6	2384264	1	WRENCH, COMBINATION,10MM,TITANIUM	Yes

4.3.10. 46-294804G1, Revision AA: Nonmagnetic Coldhead Wrenches Kit

Item No	Part Number	Qty	Description	GPRS
	46-294804G1		NON-MAGNETIC AEROQUIP WRENCHES FOR COLDHEAD.	No
2	46-294800G1	1	CASE ASSEMBLY FOR BRONZE WRENCH KIT	No
3	46-294805P1	2	BRONZE, OPEN END, 1 5/8~	Yes
4	46-294805P2	2	BRONZE, OPEN END, 1 3/8~	Yes
5	46-294805P3	2	BRONZE, OPEN END, 1 3/16~	Yes
6	46-294805P4	1	BRONZE, OPEN END, 1 1/8~	Yes
7	46-294805P5	1	BRONZE, OPEN END, 1~	Yes

4.3.11. 46-318688P1, Revision AA: Cryocooler HE Line Oil Contamination Checker

Item No	Part Number	Qty	Description	GPRS
	46-318688P1		CRYOCOOLER HE LINE OIL CONTAMINATION CHECKER.	Yes

4.3.12. 5102120, Revision 1: Charging Tool, Field Service, Precision Cryogenics

Item No	Part Number	Qty	Description	GPRS
	5102120		CHARGING TOOL, FIELD SERVICE, PRECISION CRYOGENICS PART # 090203-GE	Yes

4.3.13. 5262081, Revision 1: Pressure Relief Tool, Male, -8 Coupling

Item No	Part Number	Qty	Description	GPRS
	5262081		PRESSURE RELIEF TOOL, MALE, -8 COUPLING	Yes

4.3.14. 5272260, Revision 1: Pressure Relief Tool, Female, -8 Coupling

Item No	Part Number	Qty	Description	GPRS
	5272260		PRESSURE RELIEF TOOL, FEMALE, -8 COUPLING	Yes

GE Healthcare**Renewal Parts and Tooling**

4.3.15. 2214919, Revision 0: Stand, Coldhead, LCC

Item No	Part Number	Qty	Description	GPRS
	2214919		STAND,COLDHEAD,LCC	No

4.3.16. 46-318784G2, Revision 1: Shield Cooler Test Kit

Item No	Part Number	Qty	Description	GPRS
	46-318784G2		SHIELD COOLER TEST KIT ASSEMBLY.	No
2	2100644	1	CASE ASSEMBLY FOR SHIELD COOLER TEST KIT	No
3	46-318649P1	1	DISSOLVED SOLIDS TESTER	Yes
4	46-318648P1	1	pH METER	No
5	2100976	1	UV LIGHT LABEL ASSEMBLY	No
7	46-318829P1	1	GOGGLES FOR UV LIGHT PROTECTION.	No
8	2100660	3	SPECIMEN CONTAINER	No
9	2100661	1	SOLUTION FOR SOLIDS	No
10	2100662	1	SOLUTION FOR PH	No
13	2101986	1	~SAMPLE~ LABEL	No
14	2101987	1	~SOLIDS~ LABEL	No
15	2101988	1	~pH~ LABEL	No
16	2101994	1	SHIELD COOLER TEST KIT (46-318784G1) OPERATING INSTRUCTIONS.	No

4.3.17. 5503375, Revision 1: Repair Kit, Cold Sleeve Pump Out Port

Item No	Part Number	Qty	Description	GPRS
	5503375		Repair Kit, Cold Sleeve Pump Out Port	Yes
1	2109873-21	1	SCREW PAN HEAD Z M4X40 A2-70	No
2	2381273-5	1	NUT, SS, HEX, M4	No
3	46-252322P10	1	WASHER, BRASS, .50 PLAIN	Yes
4	46-252635P22	1	#10 X .738 OD, PLAIN STAIN STEEL 304	No
5	2236606	1	PLUG,DOUBLE O-RING,BRASS	Yes

GE Healthcare

Renewal Parts and Tooling

Thermal Stability

4.4. FRU Assemblies – Thermal Stability

4.4.1. 5177897-3, Revision 2: Kit, Pressure Transducer Replacement

Item №	Part Number	Qty	Description	GPRS
	5177897-3		KIT, PRESSURE TRANSDUCER REPLACEMENT, WIKA*	Yes
1	5177897-4	1	INSTRUCTIONS, PRESSURE TRANSDUCER REPLACEMENT KIT, WIKA	No
2	5192503	1	PRESSURE TRANSMITTER	Yes
3	46-252557P2	1	.25 NPT, BRASS	Yes
4	5212372	1	FRACTIONAL TUBE ADAPTER	No
5	5180922	1	0.25 INCH OD UNION ELBOW	No
6	2109440	1	SWAGLOK PORT CONNECTOR	No
7	46-281174P2	2	NUT TUBE .25 OD, SWAGelok #B-402-1	Yes
8	46-281173P4	2	TUBE .25 OD, SWAGelok #B-404-1	Yes
9	46-281060P4	2	TUBE .25 OD, SWAGelok #B-403-1	Yes
10	46-252065P19	1	TEFLON TAPE .50 WIDE	Yes

* 0.25" NPT transducer for conventional plumbings 2281696; 2281696-2; 5139639; 5175534

4.5. Tool Kits – Thermal Stability

4.5.1. 5263845, Revision 1: De-Icing Tool Kit

Item №	Part Number	Qty	Description	GPRS
	5263845		KIT, DEICING TOOL	No
1	5264712	1	CASE, DEICING TOOL, PELICAN 1620	No
2	5266224	1	HOSE ASSEMBLY, 20 FT	No
3	5270565	1	HOSE ASSEMBLY, 3 FT	No
4	5248857	1	INSULATED FILL PORT CORING TOOL	No
5	5248855	1	INSULATING RAMPING PORT TOOL	No
6	2280035	1	TOOL, ICE COREING	Yes
7	5266208	1	GAS WAND VALVE	No
8	5248860	1	GAS WAND TUBE	No
9	5248856	1	INSULATING ICE CHIPPING TOOL	No
10	5248765	1	TOOL, FILL PORT CORING	No
11	5248852	1	ICE CHIPPING TOOL, LOW CEILING HEIGHT	No
12	5248766	1	RAMPING PORT CORING TOOL, LOW CEILING HEIGHT	No
13	5248851	1	TOOL, FILL PORT CORING, LOW CEILING HEIGHT	No
14	5270768	1	RECONDENSOR DEICING VALVE	No
15	5268231	1	VACUUM VALVE FOR RECONDENSOR	No
16	5266816	1	REGULATOR, 75 PSI W QUICK CONNECT	No
17	5268307	1	FLOWMETER W QUICK CONNECT	No
18	5271429	1	REGULATOR, 15 PSI W QUICK CONNECT.	No

4.5.2. 2374573, Revision 2, TOOL KIT, DEICEING

Item №	Part Number	Qty	Description	GPRS
	2374573		TOOL KIT, DEICEING	No
1	2303292	2	CAP ASM, VIEWER	No
2	2278972	1	TOOL, ICE COREING	No
3	2280035	1	TOOL, ICE COREING	No
4	2281483	1	VALVE ASM, DEICEING TOOL	No
5	2347725	1	TUBE, DEICEING PURGE	No
6	46-294765G1	1	LEXAN TURRET COVER	Yes
7	5101724	1	CASE, DEICING TOOL WITH LABEL	No
8	5101708	1	TOOL DATA SHEET, DEICING KIT 2374573	No

GE Healthcare**Renewal Parts and Tooling****Helium Fill****4.6. FRU Assemblies – Helium Fill****4.6.1 46-260555G2, Revision 0: Transfer Line Extension Tip for 500L Dewer**

Item №	Part Number	Qty	Description	GPRS
	46-260555G2		TRANSFER LINE EXTENSION TIP FOR 500L DEWER	Yes
1	46-260555P2	1	.500 DIA. X .500 LONG STAINLESS STEEL ROUND	No
2	46-260555P3	1	.250 DIA X .035 WALL X 10.00 LONG STAINLESS STEEL TUBING	No

4.6.2 46-294513G1, Revision 4: Liquid Helium Transfer Line Kit Collector

Item №	Part Number	Qty	Description	GPRS
	46-294513G1		LIQUID HELIUM TRANSFER LINE KIT COLLECTOR	No
1	46-294511P1	1	HELIUM DEWAR STINGER ASSEMBLY WITH 59.25~ LONG STINGER FOR A 250 LITER DEWAR.	No
2	46-294511P2	1	HELIUM DEWAR STINGER ASSEMBLY WITH 64.25~ LONG STINGER FOR A 500 LITER DEWAR.	No
3	46-294512P1	1	FLEXIBLE LIQUID HELIUM TRANSFER LINE X 12 FT. LONG WITH .50~ X 6.00~ MALE BAYONET ON BOTH ENDS.	Yes
4	46-294512P2	1	FLEXIBLE LIQUID HELIUM TRANSFER LINE X 8 FT. LONG WITH .50~ X 6.00~ MALE BAYONET ON BOTH ENDS.	No
5	46-294512P3	1	LIQUID HELIUM CRYOSTAT STINGER ASSY. WITH 7.00~ LONG STINGER WITH THREADED END AND TEFLON TIP.	Yes
6	46-294512P4	1	LIQ He CRYOSTAT STINGER ASSY. WITH 14.75~ LONG STINGER.	No
7	46-294512P5	1	LIQ He MALE BAYONET REDUCER WITH .50~ X 3.00~ LONG BAYONET.	Yes
8	46-294512P6	1	LIQUID HELIUM MALE BAYONET COUPLING NUT.	Yes
9	46-294512P7	1	LIQUID HELIUM MALE BAYONET NUT RETAINER RING.	Yes
10	46-294512P8	1	LIQUID HELIUM MALE BAYONET O-RING.	Yes
11	46-294512P9	2	LIQ He .50~ VAC PUMP-OUT PORT, RELIEF VALVE O-RING (LARGE)	Yes
12	46-294512P10	1	LIQ He .50~ VAC PUMP-OUT PORT, RELIEF VALVE O-RING (MED).	Yes
13	46-294512P11	2	LIQ He .50~ VAC PUMP-OUT PORT; RELIEF VALVE O-RING (SMALL).	Yes
14	46-294512P12	1	16.00~ LONG CRYOSTAT STINGER EXTENSION WITH THREADED END AND TEFLON TIP.	Yes
15	46-294512P13	1	PURGE/PRECOOL ADAPTER FOR COOLDOWN OF MAGNET HELIUM RESERVOIR.	No
16	46-294512P14	1	FOAM FILLED BLACK PLASTIC CASE 26.00~ X 14.50~ X 6.63~.	Yes
17	46-294512P15	1	KNURLED BRASS BLANKING CAP WITH RUBBER SEAL - 1.50~ I.D.	Yes
18	46-294512P16	1	NITROGEN FILL LINE ADAPTER - .960-14 NGO RIGHT HANDLE MALE THREAD TO .50~ ODT X 45 DEG. FLARE MALE FITTING.	Yes
19	46-294512P17	1	HELIUM GAS LINE ADAPTER - .50~ NPT MALE TO .50~ ODT X 45 DEG. FLARE FEMALE FITTING.	Yes
20	46-294512P18	1	VACUUM PUMP-OUT PORT CAP.	No
21	46-294512P19	1	VACUUM PUMP-OUT PORT VALVE.	No
22	46-294512P20	1	PVC BAYONET PROTECTORS - SET OF TWO FOR FS-8-2B6 AND FS-12-2B6 LIQUID HELIUM TRANSFER LINES.	No
23	46-294512P21	1	TEFLON TIP FOR CRYOSTAT STINGER AND EXTENSION - SET OF TWO.	No
24	46-294512P22	1	VALVE STEM CAP KIT.	Yes
25	46-252805P2	1	10 FT. LONG NITROGEN TRANSFER LINE WITH .50~ X 45 DEG. FLARE FITTINGS ON BOTH ENDS.	No
26	46-252805P3	1	15 FT. LONG NITROGEN TRANSFER LINE WITH .50~ X 45 DEG. FLARE FITTINGS ON BOTH ENDS.	No
27	46-294512P25	1	LIQ He CRYOSTAT STINGER ASSY. WITH 17.00~ LONG STINGER.	No
28	5264435	2	STINGER	No

GE Healthcare

Renewal Parts and Tooling

4.6.3 46-271135P1, Revision AA: HOSE ASSEMBLY, 40 FT, DEWAR TO GAS CYLINDER

Item №	Part Number	Qty	Description	GPRS
	46-271135P1	1	HOSE ASSEMBLY, 40 FT, DEWAR TO GAS	Yes

4.7. Tool Kits –Cool Down and Top Off

4.7.1 46-306734G1, Revision AB: Helium Tank Regulator and Hose Kit

Item №	Part Number	Qty	Description	GPRS
	46-306734G1		NON-MAGNETIC GAS CYLINDER REGULATOR/ HOSE ASSEMBLY CONSISTING OF REGULA- TOR, HOSE, CASE, FOAM AND LABEL	No
1	46-258151P1	1	NON-MAGNETIC GAS CYLINDER REGULATOR	No
2	46-271135P16	1	20 FT DEWAR TO GAS CYLINDER HOSE ASM	No
3	46-301453P8	1	BLOW MOLDED, 17 X 13 X 7-3/8	No
4	46-306735P1	1	FOAM FOR CASE TOP, GAS REG/HOSE KIT	No
5	46-306735P2	1	FOAM FOR CASE BOTTOM, GAS REG/HOSE KIT	No
6	46-306736P1	1	LABEL FOR GAS CYLINDER REG/HOSE KIT	No

4.7.2 46-294705G1, Revision 5: Liquid Helium Transfer Line Kit

Item №	Part Number	Qty	Description	GPRS
	46-294705G1		LIQUID HELIUM TRANSFER LINE STINGER/EXTENSION/ADAPTOR KIT.	No
1	46-294512P3	1	LIQUID HELIUM CRYOSTAT STINGER ASSY. WITH 7.00~ LONG STINGER WITH THREADED END AND TEFLON TIP.	Yes
2	46-294512P4	1	LIQ He CRYOSTAT STINGER ASSY. WITH 14.75~ LONG STINGER.	No
3	46-294512P5	1	LIQ He MALE BAYONET REDUCER WITH .50~ X 3.00~ LONG BAYONET.	Yes
4	46-294512P12	1	16.00~ LONG CRYOSTAT STINGER EXTENSION WITH THREADED END AND TEFLON TIP.	Yes
5	46-294512P13	1	PURGE/PRECOOL ADAPTER FOR COOLDOWN OF MAGNET HELIUM RESERVOIR.	No
6	46-294512P14	1	FOAM FILLED BLACK PLASTIC CASE 26.00~ X 14.50~ X 6.63~.	Yes
7	46-294512P15	1	KNURLED BRASS BLANKING CAP WITH RUBBER SEAL - 1.50~ I.D.	Yes
8	46-294512P16	1	NITROGEN FILL LINE ADAPTER - .960-14 NGO RIGHT HANDLE MALE THREAD TO .50~ ODT X 45 DEG. FLARE MALE FITTING.	Yes
9	46-294512P17	1	HELIUM GAS LINE ADAPTER - .50~ NPT MALE TO .50~ ODT X 45 DEG. FLARE FEMALE FITTING.	Yes
11	46-294512P25	1	LIQ He CRYOSTAT STINGER ASSY. WITH 17.00~ LONG STINGER.	No
12	5264435	2	STINGER	No

4.7.3 2166440, Revision 0: Transfer Line, Helium

Item №	Part Number	Qty	Description	GPRS
	2166440		TRANSFER LINE, HELIUM	No

4.7.4 5508132, Revision 3: Stinger Kit, Flexible Bottom-Fill

Item №	Part Number	Qty	Description	GPRS
	5508132		Stinger Kit, Flexible Bottom-Fill	No
1	5508130	1	Stinger, Flexible Bottom-Fill	No
2	5508131	1	Guide Tube, Flexible Stinger	No
3	46-294512P15	1	KNURLED BRASS BLANKING CAP WITH RUBBER SEAL - 1.50~ I.D. AIR PRODUCTS PART # BBC-1.5	Yes
4	46-318619P1	1	NUT EXTENDE	Yes
5	46-260389P4	1	O-RING, SILICONE, PARKER 2-205	Yes
6	46-260342P9	1	ORING	Yes
7	5535563	1	Pelican 1700 Long Case	No
8	5535574	1	Label, Flexible Stinger Asm	No

Note: Extended nut 46-318619P1, O-ring 46-260389P4 and Retaining ring 46-260342P9 are used for cryostat stingers 46-294512P3, 46-294512P4 and 46-294512P25

GE Healthcare**Renewal Parts and Tooling****Chapter 5 – RAMPING AND SHIMMING****Ramping and Shimming****5. Ramping****5.1 Tool Kits - Ramping****5.1.1 46-260703G7, Revision 1: Magnet Ramping Equipment Kit**

Item №	Part Number	Qty	Description	GPRS
	46-260703G7		MAGNET RAMPING EQUIPMENT	No
1	46-260817G3	2	LEAD, EXTENSION POWER MAIN	Yes
2	46-252366P4	3	SCREW BOLT	Yes
3	46-252320P19	6	NUT, HEX, .375-16, BRASS	Yes
4	46-252322P8	12	WASHER	Yes
5	46-281046P35	3	.375-16 X 1.25 LG, HEX HD, BRASS	Yes
6	46-252351P2	4	# 10 Brass RD HD Screw	Yes
7	46-252351P51	4	SCREW, RD HD #10-32 X .438 BR	Yes
8	46-281256P1	20	CONTACT BANDS	Yes
9	46-294104P1	4	NUT, VACUUM QUICK COUPLING	Yes
10	46-294105P1	4	RING, QUICK COUPLING RETAINING	Yes
11	46-260389P1	8	SILICONE, PARKER#2-116	Yes
12	2142687	1	RAMP LEAD HOLD DOWN - SIGNA V	Yes
13	46-294236G1	1	CASE FOR SERVICE RAMPING KIT	No
14	46-294204G1	2	LEAD EXTENSION POWER ASSEMBLY	Yes
15	2152359	1	FIXTURE ASM, RAMP LEAD PLATE	Yes
16	2185743	1	TOOL, RAMP LEAD HOLD DOWN, CONQUEST	Yes
17	2193922	1	FIXTURE, RAMP LEAD, LCC	Yes
18	5323611	1	RAMP LEAD HOLD TOOL, HDW	No

5.1.2 2236233 , Revision 0: HOLD DOWN TOOL HARDWARE KIT

Item №	Part Number	Qty	Description	GPRS
	2236233		HOLD DOWN TOOL HARDWARE KIT	Yes
1	46-281003P5	2	STUD, .25-20UNC X 7.12 LG, BRASS	No
2	46-252320P15	6	NUT, HEX, BRASS .250-20	No
3	46-252322P7	2	WASHER PLAIN, TYPE A NAR, BRASS .31	No

5.1.3 2279142, Revision A: LCC Hold-Down Tool Assembly

Item №	Part Number	Qty	Description	GPRS
	2279142		LCC RAMP LEAD HOLD DOWN TOOL ASSY.	No
1	2276755	1	FIXTURE, RAMP LEAD, SHORT	Yes
2	46-281003P5	2	STUD, .25-20UNC X 7.12 LG, BRASS	No
3	46-252320P15	6	NUT, HEX, BRASS .250-20	No
4	46-252322P7	2	WASHER PLAIN, TYPE A NAR, BRASS .31	No

5.1.4 Ramp and Shim Power Supplies

Item №	Part Number	Qty	Description	GPRS
	46-260776G4		1,000 AMP MAIN COIL POWER SUPPLY	No
	46-260777G3		SHIMMING POWER SUPPLY	No

GE Healthcare**Renewal Parts and Tooling****5.1.5 2141701, Revision 1: B0 Power Supply Kit**

Item №	Part Number	Qty	Description	GPRS
	2141701		B0 POWER SUPPLY KIT.	No
1	2141700	1	DESKTOP UNREGULATED POWER SUPPLY. 12 VDC WITH 6.00~ POWER CORD. JEROME INDUSTRIES PART NUMBER ILDK-108A.	No
2	2141700-2	1	STORAGE CONTAINER. UNDERWATER KINETICS MODEL ~409 DRY BOX~ PER DRAWING NUMBER 100391.	No
3	2140361	1	LABEL, B0 QUENCH POWER SUPPLY KIT	No
4	46-294151P39	1	ARMSTRONG A-12, 5 GRAM TWIN PAK	No
5	5327239	1	WARNING LABEL FERROUS MATERIAL HDV	Yes

5.1.6 2353394-3, Revision 3: 1,000 Amp Ramp Cable Kit with Ground

Item №	Part Number	Qty	Description	GPRS
	2353394-3		KIT, 1000AMP RAMP CART AND CABLE, WGND	No
1	2135434	1	RAMP CABLE CART WITH FOAM INSERTS IN SIDES, TOP AND BOTTOM.	No
2	46-260723G1	3	(601) POWER ASSEMBLY POSITIVE 4/0	Yes
3	46-260723G2	3	(602) POWER ASSEMBLY NEGATIVE 4/0	No
4	46-260724G1	1	(604) SWITCH HEATERS	Yes
5	46-281667G1	1	SWITCH HEATERS RESONA C and D	Yes
6	46-281235P1	1	VOLTAGE SENSE LEAD, 40 FEET	Yes
7	5402320	1	LABEL, RAMP KIT UPGRADE	No
8	46-294167P20	1	BLACK/YELLOW CHECKERED TAPE	No
9	2351842	1	CABLE, POWER SUPPLY GROUNDING, HFO	No
10	2351843	1	CABLE, POWER SUPPLY JUMPER, HFO	No
11	46-320460G1	1	PDU OUTLET TO MAG/SHIM POWER SUPPLY 5-#8 COND. W/30A PLUG; RECEPT. 50FT	Yes
12	5395581	1	CABLE, RAMP AXIAL HEATER, DVW	Yes

5.1.7 2217315, Revision 0: Cable, Backup Ramp Down [Aux Heater Cable]

Item №	Part Number	Qty	Description	GPRS
	2217315		CABLE, BACKUP RAMPDOWN	Yes

5.1.8 46-318314G2, Revision AA: Ramp Cable Holder

Item №	Part Number	Qty	Description	GPRS
	46-318314G2		RAMP CABLE HOLDER	No
2	46-260723G1	2	(601) POWER ASSEMBLY POSITIVE 4/0	Yes
3	46-260723G2	2	(602) POWER ASSEMBLY NEGATIVE 4/0	No
4	46-318373G1	2	CABLE, POWER	No
5	46-318314G1	1	RAMP CABLE HOLDER	No

5.1.9 2171719, Revision 1: Ramp and Shim Cable Kit for Conquest

Item №	Part Number	Qty	Description	GPRS
	2171719		RAMP AND SHIM CABLE KIT FOR CONQUEST	No
1	46-260724G1	1	(604) SWITCH HEATERS	Yes
2	46-281235P1	1	VOLTAGE SENSE LEAD, 40 FEET	Yes

GE Healthcare**Renewal Parts and Tooling**

5.1.10 46-251865G5, Revision 0: Metrolab Teslameter Kit

Item No	Part Number	Qty	Description	GPRS
	46-251865G5		METROLAB TESLAMETER KIT	No
1	46-255839P2	1	TESLAMETER MODEL 2025	No
2	46-255839P104	1	#5 PROBE WITHOUT AMPLIFIER	No
3	46-255839P121	1	CABLE FOR METER TO PRE-AMP	Yes
4	46-255839P200	1	ORDER 46-251865G50 FOR CASE ASSEMBLY FOR TESLAMETER KIT	No
5	46-255839P103	1	#4 PROBE WITHOUT AMPLIFIER	No
6	46-306832P1	1	FOAM FOR CASE BOTTOM, TESLAMETER KIT	No
7	46-306813P4	1	METROLAB TESLAMETER KIT CASE LABEL	No
8	46-306836G1	1	HOLDER TO SECURE PROBE MAPPING FIXTURE CARRIAGE (TESLAMETER)	No
9	46-255839P102	1	#3 PROBE WITHOUT AMPLIFIER	No
10	2295525-5	1	PROBE, MINIATURE, METROLAB, NMR PART # 1062/6/10M	No
11	2378859	1	TOOL DATA SHEET	No

5.1.11 5849221, Revision 1: METROLAB PT2026 TESLAMETER KIT

Item No	Part Number	Qty	Description	GPRS
	5849221		METROLAB PT2026 TESLAMETER KIT	No
1	5835888	1	Metrolab PT2026 Tesla Meter	No
2	5850517	1	MFC2046 Transit Case	No

GE Healthcare**Renewal Parts and Tooling****Shimming****5.2 FRU Assemblies – Shimming****5.2.1 5263916, Revision 3: Kit, Field Passive Shim, 3.0T HDv**

Item №	Part Number	Qty	Description	GPRS
	5263916		KIT, 3.0T HDV FIELD PASSIVE SHIM	Yes
1	5159201-3	60	PASSIVE SHIM, HALF WIDTH, .020 INCH THICK, XRMB GRADIENT	No
2	5159201-2	40	PASSIVE SHIM, .020 INCH THICK, XRMB GRADIENT	No
3	5159201	30	PASSIVE SHIM, .031 INCH THICK, XRMB GRADIENT	No
4	2379197-4	30	SHIM, 2.0MM, HALF WIDTH, M4 THRU	No
5	2379197	30	SHIM, 2.0MM, M4 THRU	No
6	2379197-5	30	SHIM, CARBON STEEL, 3MM THK, M4 THRU	No
7	5212343	100	SQUARE NUT, PASSIVE SHIM, XRMB GRADIENT	No
8	5214479	120	SCREW, PASSIVE SHIM, XRMB GRADIENT	No

5.2.2 5411648, Revision 1: LCC 1.5T Passive Shim Kit

Item №	Part Number	Qty	Description	GPRS
	5411648		LCC SHIM KIT	Yes
1	2190204	300	PASSIVE SHIM, .001 THK. - CONQUEST	No
2	2190204-2	300	PASSIVE SHIM, .003 THK. - CONQUEST	No
3	2190204-3	300	PASSIVE SHIM, .010 THK. - CONQUEST	No
4	5133586	200	DUCT TAPE STRIP, WHITE	No
5	46-294151P24	12	EPOXY 3M-DP420 - IN DUO PACK DISPENSER. DATE CODE REQUIREMENT: MATERIAL MUST BE USED ONE YEAR FROMDATE OF RECEIPT.	No
6	5265376	1	DISPENSER, EPX PLUS II, APPLICATOR	No
7	46-252065P154	12	MIXING NOZZLE 9742, 3M SCOTCH-WELD EPX APPLICATOR, 50 ML, PN 62974299356	No
8	2158613	12	TAPE SCALE FOR PASS. SHIMS, CONQUEST	No

5.2.3 5428096, Revision 5: MR750w Field Passive Shimming Kit

Item №	Part Number	Qty	Description	GPRS
	5428096		KIT, FIELD PASSIVE SHIM, MR750W	Yes
1	5430878-3	500	SHIM, SILICONE STEEL, 0.35MM, DIA 4.5 THRU	No
4	5430878-4	240	SHIM, SILICONE STEEL, HALF 0.35MM, DIA 4.5 THRU	No
5	5430878	80	SHIM, PASSIVE, SILICONE STEEL, 2.1MM, M4 THRU	No
6	5430878-5	80	SHIM, SILICONE STEEL, 3.5MM, M4 THRU	No
7	5212343	300	SQUARE NUT, PASSIVE SHIM, XRMB GRADIENT	No
8	5214479	280	SCREW, PASSIVE SHIM, XRMB GRADIENT	No
9	5430878-6	100	SHIM, SILICONE STEEL, 3.5MM, DIA 4.5 THRU CS DIA 8.0	No
10	5430878-7	200	SHIM, SILICONE STEEL, 4.9MM, M4 THRU	No
12	5430878-2	160	SHIM, SILICONE STEEL, 4.9MM, DIA 4.5 THRU CS DIA 8.0	No
14	2385706-5	150	SCREW, FLAT HEAD, BRASS, M4 X 10	No
15	46-252283P81	40	PANDUIT CABLE TIE PLT31-C	No

GE Healthcare**Renewal Parts and Tooling****5.2.4 5342951, Revision 8: 1.5T DVW Field Passive Shimming Kit**

Item №	Part Number	Qty	Description	GPRS
	5342951		KIT, FIELD PASSIVE SHIM, 1.5T DVW	Yes
1	5430878-3	250	SHIM, SILICONE STEEL, 0.35MM, DIA 4.5 THRU	No
2	5430878-4	120	SHIM, SILICONE STEEL, HALF 0.35MM, DIA 4.5 THRU	No
3	5430878	40	SHIM, PASSIVE, SILICONE STEEL, 2.1MM, M4 THRU	No
4	5430878-5	40	SHIM, SILICONE STEEL, 3.5MM, M4 THRU	No
5	5212343	150	SQUARE NUT, PASSIVE SHIM, XRMB GRADIENT	No
6	5214479	180	SCREW, PASSIVE SHIM, XRMB GRADIENT	No
7	5430878-6	50	SHIM, SILICONE STEEL, 3.5MM, DIA 4.5 THRU CS DIA 8.0	No
8	5430878-7	50	SHIM, SILICONE STEEL, 4.9MM, M4 THRU	No
9	5430878-2	40	SHIM, SILICONE STEEL, 4.9MM, DIA 4.5 THRU CS DIA 8.0	No
10	2385706-5	75	SCREW, FLAT HEAD, BRASS, M4 X 10	No

5.2.5 2173354, Revision 10: SHIM LEAD ASM,RECTACTABLE,LCC

Item №	Part Number	Qty	Description	GPRS
	2173354		SHIM LEAD ASM,RECTACTABLE,LCC	Yes
1	2180942	1	COVER,SHIM LEAD	No
2	2180821	1	HOUSING ASSEMBLY,CONNECTOR	No
3	2180818	1	BAFFLE ASM,LCC	Yes
4	2181116	1	PIN	No
5	2180588	1	ADAPTER,CONNECTOR HOUSING	No
6	46-260714P3	1	ITT CANNON 46 PIN, FK-46-21C-1	No
7	2180682	1	PLATE, SHIM FLANGE INSULATOR	No
8	46-281234P3	24	#18 AWG CARTRIDGE BRASS HEAVY FORMVAR COATED	No
9	46-294357P1	14	#32 AWG, PHOSPHOR BRONZE WIRE HEAVY FORMVAR COATED, 4 STRAND	No
10	2181023	1	VACUUM COUPLING	No
11	5124955	1	Backing Ring, Coupling	No
12	5126342	1	Cap,Vacuum Coupling	No
13	46-252880P3	2	RECEPTACLE RG FIXED LEMO CAT. EGG2B310CLL	No
14	46-260648P2	1	No Description provided*	No
15	46-260789P18	8	SCREW #6-32 X .50LG, HEX SOC HD STAINLESS STEEL 304	Yes
16	2184199	2	CAP,VACUUM COUPLING	Yes
17	2326325	1	CAP ASM, EXTENDED	Yes
18	46-260435P1	6	ALUMINUM ALLOY 2024-0 ADEC PT. #759-3 OR EQUIV. ADEL PART # 759-3 OR EQUIV.	Yes
19	46-260325P1	3	5 INCH LONG SEGMENTS #6 NICKEL PLATED BRASS CHAIN, ROUND, .125 DIA. MCMASTER CARR #3606TG4 OR EQUIV.	Yes
20	46-260298P2	3	DRIVE SCREW, #6 X .25 LENGTH, 302 OR 410 ST. STL.	Yes
21	46-294151P5	0	CREST or LORD EPOXY 3170 DATE CODE REQUIREMENTS: MATERIAL MUST BE USED ONE YEAR FROM DATE OF RECEIPT	No
22	46-260789P18	2	SCREW #6-32 X .50LG, HEX SOC HD STAINLESS STEEL 304	No
23	46-252644P6	0	EMERSON-CUMMINGS EPOXY STYCAST 2850FT	No
24	46-252065P69	0	AMERICAN MATRIX CORP., SILICON CARBIDE ,WHISKERS.	No
25	46-294151P4	0	ARMSTRONG A-12 EPOXY	Yes
26	46-294022P1	1	EPOXY RETAINER	No
27	2180811	1	PLATE,CONNECTOR LABEL	No
28	2182830-2	2	FLANGE NUT,.25-20, STAIN STL 302, LG	No
29	2180886	1	ADAPTER,BAFFLE ASM ATTACHMENT	No
30	46-294022P2	0	EPOXY RETAINER	No
31	2184241	2	SLEEVE,VACUUM COUPLING	No
32	2184241-2	1	SLEEVE,VACUUM COUPLING	Yes
33	46-260389P3	2	O-RING, NITRILE BUNA-N #2-116	Yes
34	46-294151P52	0.01	adhesive, eccobond 286 ab general purpose	No

* Part of 46-260648 BULKHEAD CONNECTOR, 54 PIN 0.5T METASTABLE

GE Healthcare**Renewal Parts and Tooling**

5.2.6 2385710, Revision 6: KIT, FIELD PASSIVE SHIM, G3

Item No	Part Number	Qty	Description	GPRS
	2385710*		KIT, FIELD PASSIVE SHIM, G3	No
1	2379197	80	SHIM, 2.0MM, M4 THRU	No
2	2379197-2	80	SHIM, 3.0MM, DIA 4.5 CS DIA 9.0 THRU	No
3	2379197-5	80	Shim,Carbon steel, 3mm thk, M4 thru	No
4	2379197-4	80	SHIM, 2.0MM, HALF WIDTH, M4 THRU	No
5	2379093-2	50	SHIM COVER, 4.5 DIA THRU CS 9.0 DIA	No
6	2385706-4	100	SCREW, FLAT HEAD, BRASS, M4 X 8	No
7	5198786	1	shield, 3t-trm aluminized mylar	Yes

* This passive shim kit required for warm bore shim of W/WB magnet series. It IS NOT available as FRU but can be requested from Manufacturing. Contact GE HealthCare MR technical support for assist with the order.

GE Healthcare

Renewal Parts and Tooling

5.3 Tool Kits - Shimming

5.3.1 2135558, Revision 1: Shim Cable Kit

Item №	Part Number	Qty	Description	GPRS
	2135558		SHIM CABLE CASE WITH CABLES	No
1	5178030	1	SHIM CABLE CASE FROM AMERIPACK.	No
2	46-260726G1	1	RAMP AND SHIM CABLE ASSEMBLY	No
3	46-260724G1	1	(604) SWITCH HEATERS	Yes
4	2135362	1	SHIM LEAD ADAPTER - SIGNA V	Yes
5	46-320460G1	1	PDU OUTLET TO MAG/SHIM POWER SUPPLY 5-#8 COND. W/30A PLUG; RECEPT. 50FT	Yes

5.3.2 2386028, Revision 1: 3.0T Field Service Shim Camera Kit

Item №	Part Number	Qty	Description	GPRS
	2386028		KIT, FIELD SERVICE SHIM CAMERA, 3.0T	No
1	2276277	1	MAGNETIC FIELD CAMERA, NMR MAGNET	No
2	2276277-13	1	PROBE ARRAY, SHIM CAMERA, HM 32-500, 3.0T	No
3	2276277-3	1	POWER SUPPLY, SHIM CAMERA, NMR	No
4	2276277-4	1	REMOTE BOX, SHIM CAMERA, NMR	No
5	2276277-5	1	CABLE, RS232, SHIM CAMERA, NMR	No
6	2276277-6	1	CABLE, MAIN UNIT TO PS, SHIM CAMERA	No
7	2276277-7	1	MANUAL, METROLAB SHIM CAMERA, NMR	No
8	2276277-8	1	CABLE, PS FOR U.S., METROLAB, NMR	No
9	2276277-9	1	CABLE, PS FOR EUROPE, METROLAB, NMR	No
10	2281993	1	TRANSIT CASE, MAGNETIC FIELD CAMERA SHIM KIT	No
11	46-317773P1	1	MAGNETIC MATERIALS (RED)	No
12	5102665	1	LABEL, 3.0T SHIM CAMERA CASE	No
13	46-328349P1	8	PLASTIC ARROW CLIP (.248 BODY X .54 LG)	No
14	5102783	1	SOFTWARE, SHIM TOOL	No
15	2378859	1	TOOL DATA SHEET	No

5.3.3 2386028-2, Revision 1: Kit, Field Service Shim Camera, Dual 1.5T and 3.0T

Item №	Part Number	Qty	Description	GPRS
	2386028-2		KIT, FIELD SERVICE SHIM CAMERA, DUAL 1.5T AND 3.0T	No
1	2276277	1	MAGNETIC FIELD CAMERA, NMR MAGNET	No
2	2276277-13	1	Probe Array, Shim Camera, HM 32-500, 3.0T	No
3	2276277-3	1	POWER SUPPLY, SHIM CAMERA, NMR	No
4	2276277-4	1	REMOTE BOX, SHIM CAMERA, NMR	No
5	2276277-5	1	CABLE, RS232, SHIM CAMERA, NMR	No
6	2276277-6	1	CABLE, MAIN UNIT TO PS, SHIM CAMERA	No
7	2276277-7	1	MANUAL, METROLAB SHIM CAMERA, NMR	No
8	2276277-8	1	CABLE, PS FOR U.S., METROLAB, NMR	No
9	2276277-9	1	CABLE, PS FOR EUROPE, METROLAB, NMR	No
10	2281993	1	TRANSIT CASE, MAGNETIC FIELD CAMERA SHIM KIT, GMW PART NO MFC-TC	No
11	46-317773P1	1	MAGNETIC MATERIALS (RED)	No
12	5102665-2	1	LABEL, DUAL 1.5T AND 3.0T SHIM CAMERA KIT	No
13	46-328349P1	8	PLASTIC ARROW CLIP (.248 BODY X .54 LG)	No
14	5102783	1	SOFTWARE, SHIM TOOL	No
15	2378859	1	TOOL DATA SHEET	No
16	2276277-12	1	PROBE-ARRAY, SHIM CAMERA, HM 32/450 FOR 1.5T MAGNETS*	No

* This is 45 cm DSV camera used for 1.5T LCC magnets (R and CR series) warm bore shimming. Not to be used for Artist EVO upgrade.

GE Healthcare

Renewal Parts and Tooling

5.3.4 2386028-3, Revision 1: Kit, Field Service Shim Camera, Dual 3.0T and 1.5T DVW

Item №	Part Number	Qty	Description	GPRS
	2386028-3		KIT, FIELD SERVICE SHIM CAMERA, DUAL 3.0T AND 1.5T DVW	No
1	2276277	1	MAGNETIC FIELD CAMERA, NMR MAGNET	No
2	2276277-13	1	Probe Array, Shim Camera, HM 32-500, 3.0T	No
3	2276277-3	1	POWER SUPPLY, SHIM CAMERA, NMR	No
4	2276277-4	1	REMOTE BOX, SHIM CAMERA, NMR	No
5	2276277-5	1	CABLE, RS232, SHIM CAMERA, NMR	No
6	2276277-6	1	CABLE, MAIN UNIT TO PS, SHIM CAMERA	No
7	2276277-7	1	MANUAL, METROLAB SHIM CAMERA, NMR	No
8	2276277-8	1	CABLE, PS FOR U.S., METROLAB, NMR	No
9	2276277-9	1	CABLE, PS FOR EUROPE, METROLAB, NMR	No
10	2281993	1	TRANSIT CASE, MAGNETIC FIELD CAMERA SHIM KIT, GMW PART NO MFC-TC	No
11	46-317773P1	1	MAGNETIC MATERIALS (RED)	No
12	5343428	1	LABEL, 3.0T AND 1.5T DVW SHIM CAMERA KIT	No
13	46-328349P1	8	PLASTIC ARROW CLIP (.248 BODY X .54 LG)	No
14	5102783	1	SOFTWARE, SHIM TOOL	No
15	2378859	1	TOOL DATA SHEET	No
16	2276277-14	1	PROBE ARRAY, SHIM CAMERA, HM 32-500, 1.5T	No

5.3.5 5266042-4, Revision 2: Kit, Universal Field Mapping Fixture, Non-Magnetic, Extended

Item №	Part Number	Qty	Description	GPRS
	5266042-4		KIT, UNIVERSAL FIELD MAPPING FIXTURE, NON-MAGNETIC, EXTENDED	No
1	5122854	1	ASSEMBLY, UNIVERSAL FIELD CAMERA MAPPING FIXTURE	No
2	2383579	4	TOOL, G3 GRADIENT JACK, FIELD SHIM	No
3	5101699	4	SCREW, M10 X 75 HEX HEAD CAP, SPHERICAL POINT, STAINLESS STEEL UNS S31600	No
4	46-294167P81	1	SCREWDRIVER, MCMASTER CARR PN 5699A14, TORQUE INDICATOR, 0-25 IN-LBS	No
5	5102367	1	CASE, G3 PASSIVE SHIM FIELD SERVICE KIT	No
6	2379661	1	TOOL, SHIM TRAY REMOVAL HANDLE	No
7	2378859	1	TOOL DATA SHEET	No
8	5405325	1	LABEL, UNIVERSAL FIELD SERVICE MAPPING KIT	No
9	46-317773P1	1	MAGNETIC MATERIALS (RED)	No
10	5103265	1	CASE, PLASTIC, PLANO MOLDING MODEL 3620 OR EQUIVALENT	No
11	5103265-4	3	CASE, PLASTIC, PLANO MOLDING PART NUMBER 3450-18 OR EQUIVALENT	No
12	5116737	1	TOOL, RAMP LEAD CLAMP	Yes
13	2137787	1	MRT HANDHELD SHIM DRAWER PULLER	No
14	5195587	1	LEFT ARM, PATIENT END, HDV MAPPING FIXTURE	No
15	5195589	1	MAPPING FIXTURE SPIDER ASSEMBLY	No
16	5195586	1	RIGHT ARM, PATIENT END, HDV MAPPING FIXTURE	No
17	5172460	1	TOOL, PASSIVE SHIM TRAY EXTRACTION, XRMB GRADIENT	No
18	5325761	1	TOOL, SHIM EXTRACTION, XRMW	No
19	5338453	2	HEX SPACER MAPPING FIXTURE	No
20	5337880	1	TOOL, MAPPING FIXTURE SHAFT, DVW	No
21	5146169	4	M10X24MM WASHER	No

GE Healthcare**Renewal Parts and Tooling****5.3.6 5266042-5, Revision 1: Kit, Universal Field Mapping Fixture, Non-Magnetic Extended**

Item No	Part Number	Qty	Description	GPRS
	5266042-5		Kit, Universal Field Mapping Fixture, Non-Magnetic Extended	No
1	5122854	1	assembly, universal field camera mapping fixture	No
2	2383579	4	TOOL, G3 GRADIENT JACK, FIELD SHIM	No
3	5101699	4	screw, m10 x 75 hex head cap, spherical point, stainless steel uns s31600	No
4	46-294167P81	1	screwdriver, McMaster Carr PN 5699A14, torque indicator, 0-25 in-lbs	No
5	5102367	1	case, g3 passive shim feild service kit	No
6	2379661	1	TOOL, SHIM TRAY REMOVAL HANDLE	No
7	2378859	1	TOOL DATA SHEET	No
8	5782137	1	LABEL, UNIVERSAL FIELD SERVICE MAPPING KIT	No
9	46-317773P1	1	MAGNETIC MATERIALS (RED)	No
10	5103265	1	case, plastic, plano molding model 3620 or equivalent	No
11	5103265-4	3	case, plastic, plano molding part number 3450-18 or equivalent	No
12	5116737	1	tool, ramp lead clamp	Yes
13	2137787	1	MRT HANDHELD SHIM DRAWER PULLER	No
14	5195587	1	Left Arm, Patient End, HDv Mapping Fixture	No
15	5195589	1	Mapping Fixture Spider Assembly	No
16	5195586	1	Right Arm, Patient End, HDv Mapping Fixture	No
17	5172460	1	tool, passive shim tray extraction, xrm b gradient	No
18	5325761	1	tool, shim extraction, xrm w	No
19	5338453	2	Hex Spacer Mapping Fixture	No
20	5337880	1	Tool, Mapping Fixture Shaft, DVw	No
21	5146169	4	M10x24mm Washer	No
22	5736661	1	extension, universal field camera mapping fixture	Yes

5.3.7 5431812, Revision 2: Kit, STET Transport

Item No	Part Number	Qty	Description	GPRS
	5431812		KIT, STET TRANSPORT	No
1	5431810	1	MANUAL, SHIM TRAY EXTRACTION TOOL	No
2	5431813	1	CASE, LONG, STET TRANSPORT KIT	No
3	5431814	1	CASE, SHORT, STET TRANSPORT KIT	No
4	5419765	1	TOOL, SHIM TRAY EXTRACTION	No
5	5423147	1	TOOL, MANUAL SHIM TRAY REMOVAL, MR750W	No
6	5432260	1	KIT, STET ACCESSORY, MR750w	No
7	5434295	1	KIT, STET SUPPORT TOOL	No
8	5439170	2	TRANSIT CASE PACKING INSTRUCTION, STET TRANSPORT KIT	No
9	5442391	1	PREVENTIVE MAINTENANCE CHECKLIST INSTRUCTION, STET TRANSPORT KIT	No

GE Healthcare**Renewal Parts and Tooling****5.3.8 5431812-2, Revision 1: KIT, STET TRANSPORT RIO AND 750**

Item No	Part Number	Qty	Description	GPRS
	5431812-2		KIT, STET TRANSPORT RIO AND 750	No
1	5431810-2	1	MANUAL, SHIM TRAY EXTRACTION TOOL	No
2	5431813-2	1	CASE, LONG, STET TRANSPORT KIT, RIO AND 750	No
3	5431814-2	1	CASE, SHORT, STET TRANSPORT KIT, RIO AND 750	No
4	5423147	1	tool, manual shim tray removal, MR750w	No
5	5432260	1	KIT, STET ACCESSORY, MR750w	No
6	5434295	1	KIT, STET SUPPORT TOOL	No
7	5439170-2	2	Transit Case Packing Instruction, STET Transport Kit, Rio and 750	No
8	5442391-2	1	Preventive Maintenance Checklist Instruction, STET Transport Kit, Rio and 750	No
9	5419765-2	1	Tool, Shim Tray Extraction Rio and 750	No

5.3.9 5779462, Revision 2: USB to RS-232 Conversion Kit

Item No	Part Number	Qty	Description	GPRS
	5779462		USB to RS-232 Conversion Kit	No
1	5778077	1	USB to Serial Adapter Cable	Yes
2	5779341	1	Null Modem Cable, DB9F to DB25M, 10 ft	Yes
3	5779451	1	Gender Changer, Female_Female	Yes

GE Healthcare**Renewal Parts and Tooling****Chapter 6 – GRADIENT RELATED ITEMS****6. Gradient Related Items****6.1 FRU Assemblies – BRM and BRM2 gradient coils****6.1.1 2144093-5, Revision 7: BRM-D Coil, Gray Cart, Cradle**

Item №	Part Number	Qty	Description	GPRS
	2144093-5		BRM-D COIL, GRAY CART, CRADLE*	Yes
1	2134810	1	BRM CART FOR A STANDARD ROOM	No
2	2124499-7	1	GEMSPRO FRU BRM- GRAD CL W SHIELD	No
3	2134810-3	1	ALUMINUM CRADLE FOR BRM CART	No
4	46-208758P5	100	SELF LOCKING CABLE TIE	Yes
5	5312618	1	Foam	No
6	5321374	1	label, magnet shipping temperature	No

* Replaced with BRM2 5788960

6.1.2 5788960, Revision 1: BRM2 WITH GRAY CART, FRU

Item №	Part Number	Qty	Description	GPRS
	5788960		BRM2 WITH GRAY CART, FRU	Yes
1	5788959	1	BRM2 GRADIENT COIL, TESTED	No
2	2134810	1	BRM CART FOR A STANDARD ROOM	No
3	2134810-3	1	ALUMINUM CRADLE FOR BRM CART	No
4	46-208758P5	100	SELF LOCKING CABLE TIE	Yes
5	5312618	1	FOAM	No
6	5321374	1	label, magnet shipping temperature	No
7	5758309	36	SHIM TRAY ASSEMBLY, BRM2	Yes

6.1.3 5758309, Revision 4: SHIM TRAY ASSEMBLY, BRM2

Item №	Part Number	Qty	Description	GPRS
	5758309		SHIM TRAY ASSEMBLY, BRM2	Yes
1	5762243	1	passive shim tray, BRM2	No
2	5782319	1	RETAINER KNOB, SHIM TRAY, BRM2	No
3	5306007-3	1	setscrew, m8 x 12 flat slotted flat point, short glass filled nylon	No

Each BRM2 contains 36 shim trays 5758309.

GE Healthcare

Renewal Parts and Tooling

6.1.4 2214166-4, Revision 20, BRM-F LCC GRAD ISO INSTALL PART LIST

Item No	Part Number	Qty	Description	GPRS
	2214166-4		BRM-F LCC GRAD ISO INSTALL PART LIST	Yes
1	2379255	1	FOAM, GRAY TEXTURED SURFACE, 116 LG.	Yes
2	2109867-24	6	SCREW HEX SOCKET CAP M10X30 A2-70	No
3	2109867-26	2	SCREW HEX SOCKET CAP M10X40 A2-70	No
4	2370629	1	LEAD BOARD, BRM CABLES, LCC	No
5	2370909	1	COVER, BRM LEAD BOARD, LCC	No
6	2213168-3	1	Front Lower Support, BRM	Yes
7	2213167-3	1	Rear Lower Support, BRM	Yes
8	46-294151P2	0.1	LOCTITE 242 TYPE 11 GRADE N DATE CODE REQUIREMENTS: MATERIAL M	No
9	2240708	2	SUPPORT MOUNTING	No
10	46-294151P8	0.5	BOSTICK NEVER-SEEZ,REGULAR GRADE,LUBRICATING AND ANTI-SEIZE COMPOUND.	Yes
11	2181231	4	FILLER PAD CONQUEST	Yes
12	2184009	14	WASHER PLAIN - NORMAL 10.5 MM 18 MM	No
13	2174921-25	8	SCREW SOCKET HEAD CAP SCREW, M10 x L40, LOW HEAD, FULL THREAD	No
14	2109866-20	4	SCREW HEXAGON HEAD 10 MM 25 MM	Yes
15	2109873-26	8	SCREW PAN HEAD 6 MM 20 MM	Yes
16	2235604	1	ENDBELL PLATE MOUNTING SPACER	Yes
17	2238907	1	SPACER	No
18	2238908	2	SPACER SHIM 1.5MM THK	Yes
19	2238908-2	6	SPACER SHIM 2.4MM THK	Yes
20	2241176	2	RUBBER ISOLATION PAD	Yes
21	2209823-3	1	SCREW CHEESE HEAD 10MM 20MM NYLON	No
22	2209823-4	2	SCREW CHEESE HEAD 10MM 25MM NYLON	No
23	2209823-5	1	SCREW CHEESE HEAD 10MM 35MM NYLON	Yes
24	2238907-2	1	SPACER	No
25	2293535	4	SET SCREW- MTG.PAD	Yes
26	2212715-3	4	LCC BRM Mounting Pad	Yes
27	46-294151P11	0.01	LOCTITE 271, LOCTITE NO. 27121, 10 ML. BOTTLE. DATE CODE REQUIREMENTS: MATERIAL MUST BE USED WITHIN ONE YEAR FROM DATE OF MANUFACTURE	No
28	2325862	4	METRIC HEX NUT, NYLON	Yes
29	2109866-50	1	SCREW HEXAGON HEAD	No
30	46-294167P73	1	FLAT WASHER, METRIC, NYLON, M10 SCREW, 10.5MM ID, 22MM OD, 2.5MM THK MCMaster #95610A200	No
31	5139350	2	Edge Protector Pad	No
32	2370892	6	G10 FIBERGLASS FLAT WASHER	Yes
33	46-294151P48	0.2	GE FAST CURE RTV 5818 TRANSLUCENT SILICONE RUBBER ADHESIVE SEALANT	No

GE Healthcare**Renewal Parts and Tooling****6.1.5 2214166-5, Revision 16, BRM-F LCC GRAD ISO MOBILE INSTALL**

Item No	Part Number	Qty	Description	GPRS
	2214166-5		BRM-F LCC GRAD ISO MOBILE INSTALL	No
1	2379255	1	FOAM, GRAY TEXTURED SURFACE, 116 LG.	Yes
2	2109867-24	6	SCREW HEX SOCKET CAP M10X30 A2-70	No
3	2109867-26	2	SCREW HEX SOCKET CAP M10X40 A2-70	No
4	2213168-3	1	Front Lower Support, BRM	Yes
5	2213167-3	1	Rear Lower Support, BRM	Yes
6	46-294151P2	0.1	LOCTITE 242 TYPE 11 GRADE N DATE CODE REQUIREMENTS: MATERIAL M	No
7	2240708	2	SUPPORT MOUNTING	No
8	46-294151P8	0.5	BOSTICK NEVER-SEEZ, REGULAR GRADE,LUBRICATING AND ANTI-SEIZE COMPOUND.	Yes
9	2181231	4	FILLER PAD CONQUEST	Yes
10	2157414	0.2	3M TARTAN FILAMENT TAPE 48MM X 55M	No
11	2184009	16	WASHER PLAIN - NORMAL 10.5 MM 18 MM	No
12	2174921-25	8	SCREW SOCKET HEAD CAP SCREW, M10 x L40, LOW HEAD, FULL THREAD	No
13	2214162	1	GRADIENT CABLE CLAMP	Yes
14	2109866-24	2	SCREW HEXAGON HEAD 10 MM 50 MM	Yes
15	2235604	1	ENDBELL PLATE MOUNTING SPACER	Yes
16	2239148-2	1	SPACER BLOCK AND NUT STRIP BLOCK	Yes
17	46-208758P3	12	SELF LOCKING CABLE TIE, 222 N TENSILE STRENGTH, 157 MM LENGTH, 38 MM BUNDLE DIAMETER, 4.8 MM WIDTH, 1.3 MM THICK, SERVICE TEMP. RANGE -40 to 85 DEG C	Yes
18	2239904-2	6	PIPE INSULATION	No
19	2238907	1	SPACER	No
20	2238908	2	SPACER SHIM 1.5MM THK	Yes
21	2238908-2	6	SPACER SHIM 2.4MM THK	Yes
22	2238908-3	1	SPACER SHIM 16MM THK	Yes
23	2241176	2	RUBBER ISOLATION PAD	Yes
24	2209823-3	1	SCREW CHEESE HEAD 10MM 20MM NYLON	No
25	2209823-4	3	SCREW CHEESE HEAD 10MM 25MM NYLON	No
26	2209823-5	1	SCREW CHEESE HEAD 10MM 35MM NYLON	Yes
27	2238907-2	1	SPACER	No
28	2293535	4	SET SCREW- MTG.PAD	Yes
29	2212715-3	4	LCC BRM Mounting Pad	Yes
30	46-294151P11	0.01	LOCTITE 271, LOCTITE NO. 27121, 10 ML. BOTTLE. DATE CODE REQUIREMENTS: MATERIAL MUST BE USED WITHIN ONE YEAR FROM DATE OF MANUFACTURE	No
31	2325862	4	METRIC HEX NUT, NYLON	Yes
32	2109867-21	1	SCREW HEX SOCKET CAP M10X16 A2-70	No
33	46-294167P73	1	FLAT WASHER, METRIC, NYLON, M10 SCREW, 10.5MM ID, 22MM OD, 2.5MM THK MCMaster #95610A200	No
34	5139350	2	Edge Protector Pad	No
35	2370892	6	G10 FIBERGLASS FLAT WASHER	Yes
36	46-294151P48	0.2	GE FAST CURE RTV 5818 TRANSLUCENT SILICONE RUBBER ADHESIVE SEALANT	No

GE Healthcare**Renewal Parts and Tooling****6.1.6 5110336, Revision 1, Kit, BRM Gradient Cable Lead Board**

Item No	Part Number	Qty	Description	GPRS
	5110336	1	Kit, BRM Gradient Cable Lead Board	Yes
1	2370629	1	LEAD BOARD, BRM CABLES, LCC	No
2	2370909	1	COVER, BRM LEAD BOARD, LCC	No
3	2377933	1	CLAMP BASE, HOSE CABLE, BRM LD BRD	No
4	2109866-20	4	SCREW HEXAGON HEAD 10 MM 25 MM	No
5	2109866-22	4	CAP SCREW HEXAGON HEAD 10 MM X 1.5 X 35 MM LONG PER DIN 933 MATERIAL: A2 STAINLESS STEEL FINISH: PLAIN	No
6	2109866-23	6	CAP SCREW HEXAGON HEAD, 10 MM X 1.5 X 40 MM, LONG PER DIN 933, MATERIAL:A2 STAINLESS STEEL, FINISH: PLAIN	No
7	2109873-26	8	SCREW PAN HEAD 6 MM 20 MM	No
8	46-294151P8	1	BOSTICK NEVER-SEEZ, REGULAR GRADE, LUBRICATING AND ANTI-SEIZE COMPOUND.	Yes
9	2227841	2	FOAM GASKET, HOSE CABLE CLAMP	No
10	2109866-21	6	SCREW HEXAGON HEAD 10 MM 30 MM	No
11	2239902-2	1	CONFOR FOAM	Yes
6	2239902-3	1	CONFOR FOAM	Yes
7	2228286	2	FOAM SEAL CABLE HOSE CLAMP	No
8	2228287	2	FOAM BLOCK, CABLE HOSE CLAMP	No
9	2185175	8	FOAM, CLOSED CELL, BRM CONQUEST	No
10	46-208758P5	100	SELF LOCKING CABLE TIE	Yes
11	2404642	1	GRADIENT CABLE LEAD BD UPGRADE FOR BRM	No

6.2 BRM Gradient Miscellaneous

Item No	Part Number	Qty	Description	GPRS
	5142289		Gradient Cable Lead Board (BRM)	Yes
	5142288-2		Gradient Cable Lead Board Cover	Yes
	5177681		Bracket Turtle Arm Front BRM Molded	Yes
	5177458		Bracket Turtle Arm Rear BRM Molded	Yes
	5177394		Bracket Turtle Arm Molded	Yes
	5142289-2		Gradient Cable Lead Board Molded Ring	Yes

6.2.1 5182681-2, Revision 1, Gradient Cable Clamp Block Kit - BRM Injection Molded

Item No	Part Number	Qty	Description	GPRS
	5182681-2	1	Gradient Cable Clamp Block Kit - BRM Injection Molded	Yes
1	5179382	6	Clamp Block, 3 cable.	No
2	5179303	6	Cover Cable Clamp, 3 Cable.	No
3	2109874-22	36	SCREW COUNTERSUNK HEA 6 MM 25 MM	No
4	2109866-17	2	SCREW HEXAGON HEAD 6 MM 50 MM	Yes
5	2109879-2	2	WASHER PLAIN - LARGE 6.4 MM, 18MM OD	Yes

6.2.2 5182681, Revision 1, Gradient Cable Clamp Block Kit - Injection Molded

Item No	Part Number	Qty	Description	GPRS
	5182681	1	Gradient Cable Clamp Block Kit - Injection Molded	Yes
1	5179382	6	Clamp Block, 3 cable.	No
2	5179382-2	6	Clamp block, tru holes	No
3	5179303	6	Cover Cable Clamp, 3 Cable.	No
4	2109874-22	36	SCREW COUNTERSUNK HEA 6 MM 25 MM	No
5	2109874-25	36	SCREW COUNTERSUNK HEA 6 MM 40 MM	No
6	2109866-17	2	SCREW HEXAGON HEAD 6 MM 50 MM	Yes
7	2109879-2	2	WASHER PLAIN - LARGE 6.4 MM, 18MM OD	Yes

GE Healthcare

Renewal Parts and Tooling

6.3 FRU Assemblies – TRM gradient coil

6.3.1 2144093-9, Revision 5: 3T 1.5T TRM Coil, Blue Cart, Cradle

Item №	Part Number	Qty	Description	GPRS
	2144093-9		3T 1.5T TRM Coil, Blue Cart, Cradle	Yes
1	2285211-2	1	TWIN GRAD. COIL WITH RF SHIELD, 3.0T	No
2	2134810-2	1	BRM/CRM/TRM MODIFIED CART--SHORT ROOM	No
3	2134810-3	1	ALUMINUM CRADLE FOR BRM CART	No
4	2346042	1	CRADLE EXTENSION	No
5	46-208758P5	5	SELF LOCKING CABLE TIE	Yes
6	2294981	2	POLYETHYLENE TEREPHTHALATE POLY FILM	No
7	2143771-2	0.01	2 INCH TEFLON TAPE	Yes
8	2144944	2	COIL SHIPPING SUPPORT BARS UFI	No
9	2109866-20	4	SCREW HEXAGON HEAD 10 MM 25 MM	Yes
10	46-230119P1	0.01	SINGLE-FACED, CORRUGATED FIBREBOARD. 3/16~ A-FLUTE STYLE, 250 FT X 36 IN. SINGLE-FACED AND CORRUGATED, A-FLUTE STYLE. (DWG REV 4/92, JFK - OK).	No
11	2212493	0.01	FILAMENT TAPE, PREMIUM GRADE	No
12	5479338	1	Collector for FRU accessories	No
13	5479349	4	CAP for Cooling Pipes	No

6.3.2 2299498, Revision 3: TRM Cooling Circuit Repair Kit

Item №	Part Number	Qty	Description	GPRS
	2299498		TRM Cooling Circuit Repair Kit	Yes
1	2297574	1	TRM Exterior Coolant Circuit ASM	No
2	2138069	1	Product Identifier	No
3	2318411	1	Polyethylene Siphon Pump	No
4	2293612-5	4	Male Conn for ½ Hose to 3/8 MNPT (White)	No

6.3.3 5836225, Revision 1: XRMB, TRM and BRM/BRM2 shipping brackets Collector

Item №	Part Number	Qty	Description	GPRS
	5836225		TRM* Return Shipping Collector	Yes
1	5836224	2	TRM, BRM, XRMB Shipping Bracket	No
2	2144944	2	COIL SHIPPING SUPPORT BARS UFI	No
3	5489254	4	Bolts M10x1.5 thread stainless steel, 25mm length, 17mm Hexagon head	No
4	46-208758P5	10	SELF LOCKING CABLE TIE	Yes
5	2381267-9	4	WASHER, SS SPRING LOCK, M10	No
6	2109879-3	4	WASHER PLAIN - LARGE 10.5 MM 30 MM	Yes
7	2001-M12-10	2	Large Flat Washer, ISO7093, M12, Stainless Steel 140HV, RoHS Fastener	No
8	1006-M10C050-31	6	HEXAGON HEAD SCREW, ISO 4017, M10-1.5 X 50MM LONG, A2-70 STAINLESS STEEL, RoHS FASTENER	No
9	1005-M12C090-31	2	Hexagon Head Bolt, ISO 4014, M12-1.75, X90mm Long, A2-70 Stainless Steel, Passivate, RoHS Fastener	No
10	3000-M12C-40	2	HEXAGON NUT STYLE 1, ISO 4032, M12-1.75 BRASS RoHS FASTENER	No

*Fits 60cm systems gradient coils: TRM, XRMB, BRM/BRM2

GE Healthcare**Renewal Parts and Tooling****6.4 FRU Assemblies – XRMB gradient coil****6.4.1 2144093-14, Revision 3: Collector, XRMB FRU Gradient Coil with HOS, Gray Cart with Cradle**

Item №	Part Number	Qty	Description	GPRS
	2144093-14		COLLECTOR, XRMB FRU GRADIENT COIL WITH HOS, GRAY CART WITH CRADLE	Yes
1	2134810	1	BRM CART FOR A STANDARD ROOM	No
2	2134810-3	1	ALUMINUM CRADLE FOR BRM CART	No
3	5342577-2	1	COLLECTOR, BAY TESTED XRMB GRADIENT W HOS FRU, XU9537 EPOXY	No
4	5308479	1	TOOL, LEADBOARD INSTALLATION GAGE	No
5	5321374	4	LABEL, MAGNET SHIPPING TEMPERATURE	No
6	5313540	1	FOAM SEAL REAR GRADIENT COIL	No
7	5214524-6	20	WASHER, M10 LARGE OD NORDLOCK, A4 STAINLESS	No
8	46-294151P11	1	LOCKTITE 271, LOCKTITE NO. 27121, 10 ML. BOTTLE. DATE CODE REQUIREMENTS: MATERIAL MUST BE USED WITHIN ONE YEAR FROM DATE OF MANUFACTURE.	No
9	5197020-100	1	TERMINATION CAP_ROHS	Yes
10	5127507	2	0.8125 INCH TO 1.5 INCH 316 SST WORM DRIVE HOSE CLAMP	No
11	5189441-2	2	WIRE, RF SHIELD GROUND, 762MM LG., XRMB GRADIENT, 1.5T	No
12	5321374	1	label, magnet shipping temperature	No

6.4.2 5335389, Revision 1: Collector, 1.5T Magnet, XRMB Gradient Mounting

Item №	Part Number	Qty	Description	GPRS
	5335389		COLLECTOR, 1.5T MAGNET, XRMB GRADIENT MOUNTING	Yes
1	5170098	4	POSITIONING BLOCK, LCC AXI-RADIAL, XRMB GRADIENT	No
2	5170015	2	POSITIONING BLOCK, LCC RADIAL, XRMB GRADIENT	No
3	5170035	4	MOUNTING WEDGE, LCC OUTWARD, XRMB GRADIENT	No
4	5170312	5	MOUNTING WEDGE, LCC INWARD, XRMB GRADIENT	No
5	5326293-4	15	SCREW, M10X30 HEX HEAD CAP, BRASS UNS C27000, PER ASTM F468M	No
6	5146169	15	M10X24MM WASHER	No
7	46-294151P11	1	LOCKTITE 271, LOCKTITE NO. 27121, 10 ML. BOTTLE. DATE CODE REQUIREMENTS: MATERIAL MUST BE USED WITHIN ONE YEAR FROM DATE OF MANUFACTURE.	No

6.4.3 5508396, Revision 1: Collector, XRMB Busbar Assembly Replacement

Item №	Part Number	Qty	Description	GPRS
	5508396		COLLECTOR, XRMB BUSBAR ASSEMBLY REPLACEMENT	Yes
1	5508430	1	MOLDED CABLE BUSBAR ASSEMBLY, XRMB GRADIENT, MR750	No
2	46-294167P96	8	BUSHING, LEADBOARD VIBRATION DAMPING, XRMB GRADIENT	No
3	5268439	12	SHIM, POWER CONNECTION LEADBOARD AXIAL ADJUSTMENT, XRMB GRADIENT	No
4	5214524-6	20	WASHER, M10 LARGE OD NORDLOCK, A4 STAINLESS	No
5	2341885	22	NUT M10 FULL BRASS	No
6	5321557	1	CRATE, CABLE LEADBOARD ASSEMBLY, XRMB GRADIENT	No
7	5308479	1	TOOL, LEADBOARD INSTALLATION GAGE	No
8	2274863-2	6	FLANGE NUT, M10 X 1.5, A2, PER DIN6923	No
9	5415261-2	1	LOCTITE 243 10 ML	No
10	5214523-15	6	SCREW, M10 X 60 CUP POINT HEX SOCKET SET WITH LOCKING COMPOUND, A2 OR A4 STAINLESS, DIN 916	No

GE Healthcare**Renewal Parts and Tooling****6.4.4 5147904, Revision 7: Assembly, Passive Shim Tray, XRMb Gradient**

Item №	Part Number	Qty	Description	GPRS
	5147904		ASSEMBLY, PASSIVE SHIM TRAY, XRMb GRADIENT	Yes
1	5147899	1	TRAY, PASSIVE SHIM, XRMb GRADIENT	No
2	5147900	1	FITTING, SHIM TRAY RETAINING, XRMb GRADIENT	No
3	5306007-3	1	SETSCREW, M8 X 12 FLAT SLOTTED FLAT POINT, SHORT GLASS FILLED NYLON	No

6.4.5 5534976, Revision 1: MR750 Safety Cover FRU Collector

Item №	Part Number	Qty	Description	GPRS
	5534976		MR750 SAFETY COVER FRU COLLECTOR	Yes
1	5506393	1	COVER ASSEMBLY, MR750	No
2	5506138	1	MOUNTING BRACKET, SAFETY COVER, MR750	No
3	1000-M6C016-31	2	HEXAGON SOCKET HEAD CAP SCREW, M6-1, X16MM LONG, A2-70 STAINLESS STEEL, PASSIVATE, ROHS FASTENER	No
4	1000-M8C012-31	2	HEXAGON SOCKET HEAD CAP SCREW, M8-1.25, X12MM LONG, A2-70 STAINLESS STEEL, PASSIVATE, ROHS FASTENER	No
5	1000-M10C050-31	2	HEXAGON SOCKET HEAD CAP SCREW, M10-1.5, X50MM LONG, A2-70 STAINLESS STEEL, PASSIVATE, ROHS FASTENER	No
6	2000-M8-30	2	FLAT WASHER, ISO7089, M8, NYLON, ROHS FASTENER	No
7	5146169	2	M10X24MM WASHER	No
8	5394712	2	WASHER M6, GLASS LAMINATE	No

6.4.6 5407946, Revision 1: Collector, XRMb Coolant Manifold Replacement

Item №	Part Number	Qty	Description	GPRS
	5407946		Collector, XRMb Coolant Manifold Replacement	Yes
1	5401890	1	Supply Manifold Assembly, XRMb	No
2	5401890-2	1	Supply Manifold Assembly, XRMb	No
3	46-252065P155	25	hose clamp, 9mm wide x 12-20mm clamping range, stainless steel, McMaster pn 5388K16	No
4	46-294167P91	1	Tape, Electrical, Rubber Splicing, .75 Inch x 30 FT. Roll, Scotch 130C	No
5	5401891	1	Return Manifold Assembly, XRMb	No
6	5401891-2	1	Return Manifold Assembly, XRMb	No
7	5394470	1	Hose, Push-Lok	No
8	46-294167P103	1	hose fitting, kynar barbed coupling, 1 inch hose id	No
9	5338553	2	S.S. hose clamp - 1.125-1.75 inches - 0.5 inch band. IDEAL part number: 6720M Material: 316 Stainless Steel	No
10	5189441-2	2	Wire, RF Shield Ground, 762mm Lg., XRMb Gradient, 1.5T	No
11	5127507	2	0.8125 inch to 1.5 inch 316 SST Worm Drive Hose Clamp	No

GE Healthcare

Renewal Parts and Tooling

6.5 FRU Assemblies – XRMw gradient coil

6.5.1 5407761, Revision 6: XRMW FRU Gradient Coil Asm

Item No	Part Number	Qty	Description	GPRS
	5407761		XRMW FRU GRADIENT COIL ASM	Yes
1	2134810	1	BRM CART FOR A STANDARD ROOM	No
2	2134810-3	1	ALUMINUM CRADLE FOR BRM CART	No
3	5351301	1	BAY TESTED XRMW GRADIENT	No
4	5262712	1	PARKER NIPPLE WITH VALVE 0.5 MNPT . SEALING STYLE OF COUPLING. FIRST USED ON DVMR XGA,XPS, 3T RF AMP. USED WITH GE PN 5271369	No
5	5368487	1	PUSH ON BARB NPTF TO 0.5 INCH	No
6	5355126	3	PUSH-LOK HOSE 0.5 INCH	No
7	46-252065P155	2	HOSE CLAMP, 9MM WIDE X 12-20MM CLAMPING RANGE, STAINLESS STEEL, NORMA PRODUCTS INC PN 0126 6704 016	No
8	5368493	2	PLUG, 0.5 INCH BARBED END, PVC	No
9	5214524-6	20	WASHER, M10 LARGE OD NORDLOCK, A4 STAINLESS	No
10	46-294151P11	1	LOCTITE 271, LOCTITE NO. 27121, 10 ML. BOTTLE. DATE CODE REQUIREMENTS: MATERIAL MUST BE USED WITHIN ONE YEAR FROM DATE OF MANUFACTURE.	No
11	5357473	2	XRMW GRADIENT COIL SHIPPING BRACKET	No
12	46-260950P49	2	.500-13UNC X 2.00 LG., HEX HD, STAINLESS STEEL	No
13	5406803	2	450W SECONDARY ANTI WHISTLE SEAL	Yes
14	5350784	2	ACOUSTIC BARRIER END SERVICE	No
15	46-252065P151	.12	GORILLA TAPE, 2 INCH WIDE X 35 YARD ROLL, GORILLA GLUE COMPANY PN 60035-6	No
16	46-252283P81	50	PANDUIT CABLE TIE PLT31-C	No
17	46-294167P91	1	TAPE, ELECTRICAL, RUBBER SPLICING, .75 INCH X 30 FT. ROLL, SCOTCH 130C	No
18	5473879	4	SHIM TRAY ASSEMBLY, MOLDED, EXTRACTOR KNOB, VRMW*	Yes
19	5765900-50	6	HHCS M12X30 A4-80 SMARTBOLT	No
20	5214524-8	12	washer, m12 large od nordlock, a4 stainless	No

* 5473879 tray is backward compatible with legacy UA tray 5408218 - Shim Tray Asm, XRMW

6.5.2 5508397, Revision 4: Collector, DVw Busbar XZ Assembly Replacement

Item No	Part Number	Qty	Description	GPRS
	5508397		COLLECTOR, DVW BUSBAR XZ ASSEMBLY REPLACEMENT*	Yes
1	5508422	1	CABLE BUSBAR, X AND Z BUSBAR ASSEMBLY, MACHINED BOARD, DVW	No
2	2341885	10	NUT M10 FULL BRASS	No
3	5268439	12	SHIM, POWER CONNECTION LEADBOARD AXIAL ADJUSTMENT, XRMW GRADIENT	No
4	46-294167P96	6	BUSHING, LEADBOARD VIBRATION DAMPING, XRMW GRADIENT	No
5	5146169	6	M10X24MM WASHER	No
6	5214524-6	4	WASHER, M10 LARGE OD NORDLOCK, A4 STAINLESS	No
7	5368483	1	SHIPPING CRATE, XZ BUSBAR, DVW	No
8	2274863-2	4	FLANGE NUT, M10 X 1.5, A2, PER DIN6923	No
9	5415261-2	1	LOCTITE 243 10 ML	No
10	5214523-15	4	SCREW, M10 X 60 CUP POINT HEX SOCKET SET WITH LOCKING COMPOUND, A2 OR A4 STAINLESS, DIN 916	No
11	5344896	4	spacer, gradient to busbar, dvw	No

* Compatible with XRMw Rev. 16 and below

GE Healthcare

Renewal Parts and Tooling

6.5.3 5508398, Revision 4: Collector, DVW Busbar Y Assembly Replacement

Item №	Part Number	Qty	Description	GPRS
	5508398		COLLECTOR, DVW BUSBAR Y ASSEMBLY REPLACEMENT*	Yes
1	5508424	1	CABLE BUSBAR, Y ASSEMBLY, MACHINED BOARD, DVW	No
2	2341885	12	NUT M10 FULL BRASS	No
3	5268439	20	SHIM, POWER CONNECTION LEADBOARD AXIAL ADJUSTMENT, XRGB GRADIENT	No
4	46-294167P96	10	BUSHING, LEADBOARD VIBRATION DAMPING, XRGB GRADIENT	No
5	5146169	10	M10X24MM WASHER	No
6	5214524-6	2	WASHER, M10 LARGE OD NORDLOCK, A4 STAINLESS	No
7	5368484	1	SHIPPING CRATE, Y BUSBAR, DVW	No
8	2274863-2	2	FLANGE NUT, M10 X 1.5, A2, PER DIN6923	No
9	5415261-2	1	LOCTITE 243 10 ML	No
10	5214523-15	2	SCREW, M10 X 60 CUP POINT HEX SOCKET SET WITH LOCKING COMPOUND, A2 OR A4 STAINLESS, DIN 916	No
11	5344896	2	spacer, gradient to busbar, dvw	No

* Compatible with XRMw Rev. 16 and below

6.5.4 5342776, Revision 7: XRMW Gradient Coil Mounting Kit FRU

Item №	Part Number	Qty	Description	GPRS
	5342776		XRMW GRADIENT COIL MOUNTING KIT FRU	Yes
1	5408938	8	MOUNTING WEDGE, XRMW GRADIENT	No
2	5352691	2	XRMW GRADIENT MOUNTING STOP BLOCK	No
3	5473683-2	1	XRMW GRADIENT MOUNTING PAD	No
4	5352697	16	XRMW GRADIENT MOUNTING BLOCK SHIM	No
5	5363536-20	4	SCREW, M10X30 HEX HD CAP, SILICON BRONZE UNS C65100, PER ASTM F468M	No
6	5326293-5	8	SCREW, M10X35 HEX HEAD CAP, BRASS UNS C27000, PER ASTM F468M	No
7	46-294151P11	1	LOCTITE 271, LOCTITE NO. 27121, 10 ML. BOTTLE. DATE CODE REQUIREMENTS: MATERIAL MUST BE USED WITHIN ONE YEAR FROM DATE OF MANUFACTURE.	No
8	5146169	14	M10X24MM WASHER	No
9	46-294151P8	.01	BOSTIK NEVER-SEEZE LUBRICATING AND ANTI-SEEZE COMPOUND. HELICOIL SYSTEMS GROUP - DANBURY, CT OR EQUIVALENT. PART 49501 1 OZ TUBE DATE CODE REQUIREMENTS: MATERIAL MUST BE USED WITHIN TEN YEARS FROM DATE OF RECEIPT.	Yes
10	5389085	2	XRMW GRADIENT STOP BLOCK LH	No
11	5389090-4	4	SCREW M10x30 FLAT HEX SOCKET HEAD STAIN ST	No
13	5566941	8	INWARD WEDGE BOTTOM	No
14	2104944-7	2	RF SHIELD RETAINER – 80.65	No

6.5.5 5508399, Revision 4: Collector, MR750W Busbar XZ Assembly Replacement

Item №	Part Number	Qty	Description	GPRS
	5508399		COLLECTOR, MR750W BUSBAR XZ ASSEMBLY REPLACEMENT*	Yes
1	5508425	1	ASSEMBLY, X AND Z CABLE BUSBAR, MR750W	No
2	3000-M10C-40	6	HEXAGON NUT, ISO4032, M10-1.5, BRASS, ROHS FASTENER	No
4	46-294167P96	6	BUSHING, LEADBOARD VIBRATION DAMPING, XRGB GRADIENT	No
5	5146169	6	M10X24MM WASHER	No
6	5214524-6	4	WASHER, M10 LARGE OD NORDLOCK, A4 STAINLESS	No
7	5368483	1	SHIPPING CRATE, XZ BUSBAR, DVW	No
8	2274863-2	4	FLANGE NUT, M10 X 1.5, A2, PER DIN6923	No
9	5415261-2	1	LOCTITE 243 10 ML	No
10	5214523-15	4	SCREW, M10 X 60 CUP POINT HEX SOCKET SET WITH LOCKING COMPOUND, A2 OR A4 STAINLESS, DIN 916	No
11	5538449	4	NYLON-INSERT HEX LOCKNUT, .250-20, BRASS	No

* Compatible with XRMw Rev. 16 and below

GE Healthcare

Renewal Parts and Tooling

6.5.6 5508401, Revision 2: Collector, MR750w Busbar Y Assembly Replacement

Item №	Part Number	Qty	Description	GPRS
	5508401		COLLECTOR, MR750W BUSBAR Y ASSEMBLY REPLACEMENT*	Yes
1	5508426	1	ASSEMBLY, Y CABLE BUSBAR, MR750W	No
2	3000-M10C-40	10	HEXAGON NUT, ISO4032, M10-1.5, BRASS, ROHS FASTENER	No
4	46-294167P96	10	BUSHING, LEADBOARD VIBRATION DAMPING, XRGB GRADIENT	No
5	5146169	10	M10X24MM WASHER	No
6	5214524-6	2	WASHER, M10 LARGE OD NORDLOCK, A4 STAINLESS	No
7	5368484	1	SHIPPING CRATE, Y BUSBAR, DVW	No
8	2274863-2	2	FLANGE NUT, M10 X 1.5, A2, PER DIN6923	No
9	5415261-2	1	LOCTITE 243 10 ML	No
10	5214523-15	2	SCREW, M10 X 60 CUP POINT HEX SOCKET SET WITH LOCKING COMPOUND, A2 OR A4 STAINLESS, DIN 916	No
11	5538449	2	NYLON-INSERT HEX LOCKNUT, .250-20, BRASS	No

* Compatible with XRMw Rev. 16 and below

6.5.7 5845904, Revision 2: 450w M12 Busbar XZ Assembly Replacement

Item №	Part Number	Qty	Description	GPRS
	5845904		450w M12 Busbar XZ Assembly Replacement*	Yes
1	5830725	1	CABLE BUSBAR, X AND Z BUSBAR ASSEMBLY, MACHINED BOARD, 450W	No
2	3000-M10C-40	6	Hexagon Nut, ISO4032, M10-1.5, Brass, RoHS Fastener	No
3	5268439	12	shim, power connection leadboard axial adjustment, xrgb gradient	No
4	46-294167P96	6	bushing, leadboard vibration damping, xrgb gradient	No
5	5146169	6	M10x24mm Washer	No
6	5214524-8	6	washer, m12 large od nordlock, a4 stainless	No
7	5765900-50	4	HHCS M12X30 A4-80 SMARTBOLT	No
8	5415261-2	1	Loctite 243 10 ml	No
9	5849267	4	M12 Spacer, Gradient to Busbar, DVw	No
10	5368483	1	Shipping Crate, XZ Busbar, DVw	No

* Compatible with XRMw Rev. 17 and above

** Not compatible with Artist IPM. MR750w busbars to be used on PM magnets with XRMw

6.5.8 5845905, Revision 1: 450w M12 Busbar Y Assembly Replacement

Item №	Part Number	Qty	Description	GPRS
	5845905		450w M12 Busbar Y Assembly Replacement*	Yes
1	5830724	1	CABLE BUSBAR, Y ASSEMBLY, MACHINED BOARD, 450W	No
2	3000-M10C-40	10	Hexagon Nut, ISO4032, M10-1.5, Brass, RoHS Fastener	No
3	5268439	20	shim, power connection leadboard axial adjustment, xrgb gradient	No
4	46-294167P96	10	bushing, leadboard vibration damping, xrgb gradient	No
5	5146169	10	M10x24mm Washer	No
6	5214524-8	2	washer, m12 large od nordlock, a4 stainless	No
7	5765900-50	2	HHCS M12X30 A4-80 SMARTBOLT	No
8	5415261-2	1	Loctite 243 10 ml	No
9	5849267	2	M12 Spacer, Gradient to Busbar, DVw	No
10	5368484	1	Shipping Crate, Y Busbar, DVw	No

* Compatible with XRMw Rev. 17 and above

** Not compatible with Artist IPM. MR750w busbars to be used on PM magnets with XRMw

GE Healthcare

Renewal Parts and Tooling

6.5.9 5845907, Revision 1: 750w M12 Busbar XZ Assembly Replacement

Item №	Part Number	Qty	Description	GPRS
	5845907		750w M12 Busbar XZ Assembly Replacement*	Yes
1	5831202	1	ASSEMBLY, X AND Z CABLE BUSBAR, 750W	No
2	3000-M10C-40	6	Hexagon Nut, ISO4032, M10-1.5, Brass, RoHS Fastener	No
3	46-294167P96	6	bushing, leadboard vibration damping, xrb gradient	No
4	5146169	6	M10x24mm Washer	No
5	5214524-8	6	washer, m12 large od nordlock, a4 stainless	No
6	5765900-50	4	HHCS M12X30 A4-80 SMARTBOLT	No
7	5415261-2	1	Loctite 243 10 ml	No
8	5538449	4	nylon-insert hex locknut, .250-20, brass	No
9	5368483	1	Shipping Crate, XZ Busbar, DVw	No

* Compatible with XRMw Rev. 17 and above

** Compatible with Artist IPM with XRMw Rev.17 and above

6.5.10 5845908, Revision 1: 750w M12 Busbar Y Assembly Replacement

Item №	Part Number	Qty	Description	GPRS
	5845908		750w M12 Busbar Y Assembly Replacement*	Yes
1	5831201	1	ASSEMBLY, Y CABLE BUSBAR, 750W	No
2	3000-M10C-40	10	Hexagon Nut, ISO4032, M10-1.5, Brass, RoHS Fastener	No
3	46-294167P96	10	bushing, leadboard vibration damping, xrb gradient	No
4	5146169	10	M10x24mm Washer	No
5	5214524-8	2	washer, m12 large od nordlock, a4 stainless	No
6	5765900-50	2	HHCS M12X30 A4-80 SMARTBOLT	No
7	5415261-2	1	Loctite 243 10 ml	No
8	5538449	2	nylon-insert hex locknut, .250-20, brass	No
9	5368484	1	Shipping Crate, Y Busbar, DVw	No

* Compatible with XRMw Rev. 17 and above

** Compatible with Artist IPM with XRMw Rev.17 and above

6.5.11 5407963, Revision 4: XRMw Gradient Coil Mounting Kit, FRU, MR750w

Item №	Part Number	Qty	Description	GPRS
	5407963		XRMW GRADIENT COIL MOUNTING KIT, FRU, MR750W	Yes
1	46-170684P1	.1	LOCTITE 271, 0.5CC VIAL12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER	Yes
3	5326293-5	20	Screw, M10x35 Hex Head Cap, Brass UNS C2700, per ASTM F468M	No
4	5146169	20	M10x24mm Washer	No
5	5566938	8	INWARD WEDGE, TOP	No
6	5473683-2	1	Gradient Mounting Pad	No
8	5419679	4	AXIAL END STOP, MR750W	No
9	5566941	8	INWARD WEDGE, BOTTOM	No
10	2104944-7	2	RF SHIELD RETAINER – 80.65	No

6.5.12 5431091, Revision 3: Shim Tray Locking Ring, XRMw, FRU

Item №	Part Number	Qty	Description	GPRS
	5431091		SHIM TRAY LOCKING RING, XRMW, FRU	Yes
1	5454121	1	SHIM RETAINING RING, XRMW	No
2	2162974-13	13	SCREW, THREAD CUTTING, STAINLESS	No
4	U1-SVC-520047	1	SAND PAPER, 60 GRIT, PK OF 10	No

GE Healthcare**Renewal Parts and Tooling****6.5.13 5534973, Revision 1: MR450w Safety Cover FRU Collector**

Item №	Part Number	Qty	Description	GPRS
	5534973		MR450W SAFETY COVER FRU COLLECTOR	Yes
1	5506392	1	COVER ASSEMBLY, 450W	No
2	5508411	1	HSHCS, M10-1.5 X 60MM LG, 316 SS	No
3	2370892	1	G10 FIBERGLASS FLAT WASHER	Yes

6.5.14 5534975, Revision 2: MR750w Safety Cover FRU Collector

Item №	Part Number	Qty	Description	GPRS
	5534975		MR750W SAFETY COVER FRU COLLECTOR	Yes
1	5507842	1	SAFETY COVER ASSEMBLY, MR750W	No
2	5538449	2	NYLON-INSERT HEX LOCKNUT, .250-20, BRASS	No
3	5535175	2	SET SCREW, .25-20 X 1.25, SS316	No
4	5534978	2	G10 WASHER	No
5	5415261-2	1	LOCTITE 243 10 ML	No

6.5.15 5538535, Revision 1: Collector FRU, Clamp Replacement, MR750w

Item №	Part Number	Qty	Description	GPRS
	5538535		COLLECTOR FRU, CLAMP REPLACEMENT, MR750W	Yes
1	5508113	2	CLAMP, BUSBAR, XRMW GRADIENT, MR750W	No
2	5538449	6	NYLON-INSERT HEX LOCKNUT, .250-20, BRASS	No

6.5.16 5407756, Revision 1: XRMW Coolant Manifold FRU

Item №	Part Number	Qty	Description	GPRS
	5407756		XRMW Coolant Manifold FRU	Yes
1	5401763	1	XRMW Manifold Supply	No
2	5401764	1	XRMW Manifold Return	No
3	46-252065P155	26	hose clamp, 9mm wide x 12-20mm clamping range, stainless steel, McMaster pn 5388K16	No
4	46-294167P91	1	Tape, Electrical, Rubber Splicing, .75 Inch x 30 FT. Roll, Scotch 130C	No

6.5.17 5850819, Revision 1: XRMW and VRMW Shipping Collector, FRU

Item №	Part Number	Qty	Description	GPRS
	5850819		XRMW and VRMW Shipping Collector, FRU	Yes
1	2001-M12-10	2	Large Flat Washer, ISO7093, M12, Stainless Steel 140HV, RoHS Fastener	No
2	3000-M12C-40	2	HEXAGON NUT STYLE 1, ISO 4032, M12-1.75 BRASS RoHS FASTENER	No
3	1005-M12C090-31	2	Hexagon Head Bolt, ISO 4014, M12-1.75, X90mm Long, A2-70 Stainless Steel, Passivate, RoHS Fastener	No
4	1006-M10C050-31	2	HEXAGON HEAD SCREW, ISO 4017, M10-1.5 X 50MM LONG, A2-70 STAINLESS STEEL, RoHS FASTENER	No
5	2109879-3	2	WASHER PLAIN - LARGE 10.5 MM 30 MM	Yes
6	2381267-10	2	WASHER, SS SPRING LOCK, M12	No
7	2381267-9	2	WASHER, SS SPRING LOCK, M10	No
8	5357473	2	XRMW Gradient Coil Shipping Bracket	No

GE Healthcare

Renewal Parts and Tooling

6.5.18 5620700, Revision 2: collector, Kizuna 3T busbar Y assembly replacement

Item №	Part Number	Qty	Description	GPRS
	5620700		collector, Kizuna 3T busbar Y assembly replacement	Yes
1	5720237	1	BusBar Molded Y VRMw 3TLC Assembly	No
2	3000-M10C-40	3	Hexagon Nut, ISO4032, M10-1.5, Brass, RoHS Fastener	No
3	46-294167P96	3	bushing, leadboard vibration damping, xrmw gradient	No
4	5146169	3	M10x24mm Washer	No
5	5214524-6	2	washer, m10 large od nordlock, a4 stainless	No
6	2116325-2	2	STUD SST M10 X 65	Yes
7	2274863-2	2	flange nut, m10 x 1.5, a2, per din6923	No
8	5620700	1	collector, Kizuna 3T busbar Y assembly replacement	Yes

6.5.19 5620701, Revision 2: collector, Kizuna 3T busbar XZ assembly replacement

Item №	Part Number	Qty	Description	GPRS
	5620701		collector, Kizuna 3T busbar XZ assembly replacement	Yes
1	5720236	1	BusBar Molded X-Z VRMw 3TLC Assembly	No
2	3000-M10C-40	5	Hexagon Nut, ISO4032, M10-1.5, Brass, RoHS Fastener	No
3	46-294167P96	5	bushing, leadboard vibration damping, xrmw gradient	No
4	5146169	5	M10x24mm Washer	No
5	5214524-6	4	washer, m10 large od nordlock, a4 stainless	No
6	2116325-2	4	STUD SST M10 X 65	Yes
7	2274863-2	4	flange nut, m10 x 1.5, a2, per din6923	No
8	5620701	1	collector, Kizuna 3T busbar XZ assembly replacement	Yes

6.5.20 5877116, Revision 5: Kizuna 3T Terminal Cover Kit

Item №	Part Number	Qty	Description	GPRS
	5877116	1	Kizuna 3T Terminal Cover Kit	Yes
1	5845221	1	KIZUNA Z TERMINAL COVER ASSEMBLY	No
3	5845237	2	KIZUNA XY TERMINAL COVER ASSEMBLY	No
4	5480218	3	Clamp, Cable upper Kizuna	No
5	5888021	1	KIZUNA INTERFACE BRACKET, M12 Safety Cover-Clamp Compatible	No
6	5573622	4	Standoff interface bracket Kizuna	No
7	1000-M10C035-31	4	Hexagon Socket Head Cap Screw, M10-1.5, X35mm Long, A2-70 Stainless Steel, Passivate, RoHS Fastener	No
8	1006-M8C035-31	6	Hexagon Head Screw, ISO 4017, M8-1.25, X35mm Long, A2-70 Stainless Steel, Passivate, RoHS Fastener	No

6.5.21 5483563, Revision 1: PETMR New FG Catalogs

Item №	Part Number	Qty	Description	GPRS
	5407761-2	1	FRU, XRMW Gradient with no Acoustic Barrier	Yes

6.5.22 5788902, Revision 2: Gradient termination Lug kit

Item №	Part Number	Qty	Description	GPRS
	5788902		Gradient termination Lug kit	Yes
1	5392270	26	LUG .50 STUD SIZE 1 AWG WIRE	No
2	5392407	10	LUG .50 STUD SIZE 2-0 AWG WIRE SIZE	No
3	5368384	8	LUG .25 STUD 6AWG WIRE	No

GE Healthcare**Renewal Parts and Tooling****6.5.23 5861176, Revision 1, DV Gradient Cable Lug Kit - 1 AWG (Pink Lugs)**

Item No	Part Number	Qty	Description	GPRS
	5861176		DV Gradient Cable Lug Kit - 1 AWG (Pink Lugs)	Yes
1	5799022	8	HOMAS and BETTS Terminal 54109	No
2	5799023	4	THOMAS and BETTS terminal 54155-TB	No
3	5826292	4	Copper Two-Hole Lug, Standard Barrel, Peep Hole, Max 35kV, Wire Size #1 Flex 225/24, 1/2 Inch Bolt Size, 1-3/4 Inch Hole Spacing, Tin Plated, Die Color Code Pink, Die Code 42	No
4	5368384	3	LUG .25 STUD 6AWG WIRE	No
5	5410522	1	lug .375 Stud 6 awg wire	No

6.5.24 5861177, Revision 1, DV Gradient Cable Lug Kit - 2/0 AWG (Orange Lugs)

Item No	Part Number	Qty	Description	GPRS
	5861177		DV Gradient Cable Lug Kit - 2/0 AWG (Orange Lugs)	Yes
1	5368413	8	LUG .375 STUD 3-0 AWG WIRE SIZE	No
2	5826504	4	Copper One-Hole Lug, Standard Barrel, Peep Hole, Max 35kV, Wire Size 3/0 AWG, 1/2 Inch Bolt Size, Tin Plated, Die Code 50, Die Color Code Orange	No
3	5826505	4	Copper Two-Hole Lug, Standard Barrel, Peep Hole, Max 35kV, Wire Size 3/0 AWG, 1/2 Inch Bolt Size, 1-3/4 Inch Hole Spacing, Tin Plated, Die Code 50, Die Color Code Orange	No
4	5368384	3	LUG .25 STUD 6AWG WIRE	No
5	5410522	1	lug .375 Stud 6 awg wire	No

GE Healthcare**Renewal Parts and Tooling****6.6 FRU Assemblies – VRMw gradient coil****6.6.1 5534741, Revision 5: VRMW FRU Gradient Coil Asm**

Item №	Part Number	Qty	Description	GPRS
	5534741		VRMw FRU Gradient Coil Asm	Yes
1	2134810	1	BRM Cart for Standard Room	No
2	2134810-3	1	Aluminum Cradle for BRM Cart	No
3	5534750-2	1	Bay Tested VRMw for 1.5T	No
4	5214524-6	10	Washer, m10 large od nordlock, a4 stainless	No
5	5415261	1	Loctite 243, .5ml	No
6	5357473	2	XRmW GradientCoil Shipping Bracket	No
7	46-260950P49	2	.500-13UNC X 2.00 LG., Hex HD, Stainless Steel	No
8	46-252283P81	50	Panduit Cable Tie PLT31-C	No
9	46-294167P91	1	Roll, Electrical, Rubber Splicing, .75 Inch X 30 FT. Roll, Scotch 130C	No
10	5473879	4	Shim Tray Assembly, Molded, Extractor Knob, VRMw	Yes
11	5566938	8	Inward_Wedge_Top	No
12	5566941	8	Inward_Wedge_Bottom	No
13	5543079	2	Thermal Barrier Ends VRMw	No
14	5214524-8	3	Washer, m12 large od nordlock, a4 stainless	No
15	5473683-2	1	Gradient mounting pad	No

6.6.2 5473880, Revision 2: HM Magnet, Molded Shim Tray Assembly, VRMw and XRMw

Item №	Part Number	Qty	Description	GPRS
	5473880		Molded Shim Tray Assembly, VRMW	Yes
1	5473878	1	Molded Passive Shim Tray	No
2	5262572	1	XRMW Shim Tray Retainer	No
3	5306007-3	1	SETSCREW, M8 X 12 FLAT SLOTTED FLAT POINT, SHORT GLASS FILLED NYLON	No

6.6.3 5473879, Revision 2: UA Magnet, Molded Shim Tray Assembly, VRMw and XRMw

Item №	Part Number	Qty	Description	GPRS
	5473879		Shim Tray Assembly, Molded, Extractor Knob, VRMW	Yes
1	5473878	1	Molded Passive Shim Tray	No
2	5453280	1	Knob, shim tray retainer, XRMW	No
3	5306007-3	1	SETSCREW, M8 X 12 FLAT SLOTTED FLAT POINT, SHORT GLASS FILLED NYLON	No

6.6.4 5542955, Revision 1: Collector, VRMW Coolant Manifold FRU

Item №	Part Number	Qty	Description	GPRS
	5542955		VRMW Coolant Manifold FRU	Yes
1	5500297	1	VRMW_MANIFOLD_SUPPLY	No
2	5500298	1	VRMW_MANIFOLD_RETURN	No
3	46-252065P155	26	hose clamp, 9mm wide x 12-20mm clamping range, stainless steel, McMaster pn 5388K16	No
4	46-294167P91	1	Tape, Electrical, Rubber Splicing, .75 Inch x 30 FT. Roll, Scotch 130C	No
5	5496134	3	NYLON TUBING 6MM OD 4MM ID WHITE	No
6	5542957	4	SMC Plastic Straight Union KQ2	No
7	5441374	8	Locking Clips for Plastic Tubing	No

GE Healthcare**Renewal Parts and Tooling****6.6.5 5850819, Revision 1: XRMW and VRMW Shipping Collector, FRU**

Item No	Part Number	Qty	Description	GPRS
	5850819		XRMW and VRMW Shipping Collector, FRU	Yes
1	2001-M12-10	2	Large Flat Washer, ISO7093, M12, Stainless Steel 140HV, RoHS Fastener	No
2	3000-M12C-40	2	HEXAGON NUT STYLE 1, ISO 4032, M12-1.75 BRASS RoHS FASTENER	No
3	1005-M12C090-31	2	Hexagon Head Bolt, ISO 4014, M12-1.75, X90mm Long, A2-70 Stainless Steel, Passivate, RoHS Fastener	No
4	1006-M10C050-31	2	HEXAGON HEAD SCREW, ISO 4017, M10-1.5 X 50MM LONG, A2-70 STAINLESS STEEL, RoHS FASTENER	No
5	2109879-3	2	WASHER PLAIN - LARGE 10.5 MM 30 MM	Yes
6	2381267-10	2	WASHER, SS SPRING LOCK, M12	No
7	2381267-9	2	WASHER, SS SPRING LOCK, M10	No
8	5357473	2	XRMW Gradient Coil Shipping Bracket	No

GE Healthcare**Renewal Parts and Tooling****6.7 FRU Assemblies – HRMw gradient coil****6.7.1 5768102, Revision 2: Collector, HRMw Gradient Coil FRU**

Item №	Part Number	Qty	Description	GPRS
	5768102		Collector, HRMw Gradient Coil FRU	Yes
1	5788672	1	HRMw PACKAGING FOR FRU	No
2	5727708	1	HRMw GRADIENT TESTED	No
3	5620757	2	Shim Tray Asm, HRMw	Yes
4	5837423	2	MANIFOLD SUPPORT, G3, HRMW	No
5	5822816-2	2	Male JIC 37 degree union	No
6	5837415	1	MANIFOLD BRACKET, G3, HRMW	No
7	1000-M6C012-33	2	Hexagon Socket Head Cap Screw Full threaded, ISO 4762, M6-1 X 12mm Long, Grade A4, AISI 316 Stainless Steel, RoHS Fastener	No

* This HRMWA style coil discontinued and shall be replaced with 5927118 once stock depleted

6.7.2 5863473, Revision 4: COLLECTOR, HRMW GRADIENT BRAID-WOUND HOS FRU

Item №	Part Number	Qty	Description	GPRS
	5863473		COLLECTOR, HRMW GRADIENT BRAID-WOUND HOS FRU	Yes
1	5788672	1	HRMw PACKAGING FOR FRU	No
2	5863471	1	HRMW GRADIENT BRAID-WOUND HOS TESTED	No
3	5620757	2	Shim Tray Asm, HRMw	Yes
4	5837423	2	MANIFOLD SUPPORT, G3, HRMW	No
5	5822816-2	2	Male JIC 37 degree union	No
6	5837415	1	MANIFOLD BRACKET, G3, HRMW	No
7	1000-M6C012-33	2	Hexagon Socket Head Cap Screw Full threaded, ISO 4762, M6-1 X 12mm Long, Grade A4, AISI 316 Stainless Steel, RoHS Fastener	No

* This HRMWB style coil is compatible with Signa Premier without HOS or with HOS ProDrive Driver

6.7.3 5927118, Revision 2: HRMW ATLANTIS FRU + PRO DRIVE

Item №	Part Number	Qty	Description	GPRS
	5927118		HRMW ATLANTIS FRU + PRO DRIVE	Yes
1	5863473	1	COLLECTOR, HRMW GRADIENT BRAID-WOUND HOS FRU	Yes
2	5933483	1	PRO-Drive Amp Collector for HRMw-Atlantis	No

* This HRMWB style coil is compatible with Signa Premier without HOS or with HOS NAVshim driver when 5768102 is unavailable. This kit has no stock and will be manufactured per order.

GE Healthcare**Renewal Parts and Tooling****6.7.4 5768103, Revision 1: Collector, HRMw Gradient Coil Mounting Kit, UA, FRU**

Item №	Part Number	Qty	Description	GPRS
	5768103		Collector, HRMw Gradient Coil Mounting Kit, UA, FRU	Yes
1	5721157	4	axial end stop, HRMw gradient, MR750w	No
2	5146169	20	M10x24mm Washer	No
3	5326293-5	20	Screw, M10x35 Hex Head Cap, Brass UNS C27000, per ASTM F468M	No
4	5473683-3	1	ISOLATION MAT GRADIENT	No
5	5566941	16	INWARD_WEDGE_BOTTOM	No
6	2274863-2	2	flange nut, m10 x 1.5, a2, per din6923	No
7	5214524-6	2	washer, m10 large od nordlock, a4 stainless	No
8	5214524-8	12	washer, m12 large od nordlock, a4 stainless	No
9	5415261	1	Loctite 243 .5ml	No
10	5765900-50	6	HHCS M12X30 A4-80 SMARTBOLT	No
11	5777582	1	FOAM SE 12 O CLOCK	No
12	5777583	1	FOAM SE 3 OCLOCK	No
13	5777584	1	FOAM SE 9 OCLOCK	No
14	5777585	1	FOAM SE YX	No
15	5777586	1	FOAM SE Z	No
16	5777587	1	FOAM SE 6 OCLOCK R	No
17	5777588	1	FOAM SE 6 OCLOCK L	No
18	46-252283P81	50	Panduit Cable Tie PLT3I-C	No
19	46-294167P91	1	Tape, Electrical, Rubber Splicing, .75 Inch x 30 FT. Roll, Scotch 130C	No
20	5768286	1	Ground Strap	No
21	2116325-7	1	STUD SST M10 X 35MM	No
22	1006-M10C020-40	2	Hexagon Head Screw, ISO 4017, M10-1.5, X20mm Long, Brass, RoHS Fastener	No
23	1000-M6C050-33	3	Hexagon Socket Head Cap Screw Partial threaded, ISO 4762, M6-1 X 50mm Long, Grade A4, AISI 316 Stainless Steel, RoHS Fastener	No
24	5782118	3	HRMw manifold support	No

6.7.5 5768105, Revision 2: HRMw Busbar, UA, FRU

Item №	Part Number	Qty	Description	GPRS
	5768105	1		Yes
1	5770979	1	Bus Bar YX 3TLC HRMw	No
2	5770984	1	Z Bus Bar 3TLC HRMw	No
3	2116325-10	8	STUD, M10 X 75MM, STAINLESS STEEL	No
4	3000-M10C-40	21	Hexagon Nut, ISO4032, M10-1.5, Brass, RoHS Fastener	No
5	46-294167P96	9	bushing, leadboard vibration damping, xrm b gradient	No
6	5146169	9	M10x24mm Washer	No
7	5214524-6	15	washer, m10 large od nordlock, a4 stainless	No
8	5214524-8	12	washer, m12 large od nordlock, a4 stainless	No
9	5765900-50	6	HHCS M12X30 A4-80 SMARTBOLT	No

GE Healthcare**Renewal Parts and Tooling**

6.7.6 5778991, Revision 2: HRMw Busbar, G3, FRU

Item No	Part Number	Qty	Description	GPRS
	5778991	1	HRMw Busbar, G3, FRU	Yes
1	5774647	1	XY BUSBAR ASM, HRMW, G3	No
2	5774653	1	Z BUSBAR ASM, HRMW, G3	No
3	1000-M10C050-33	8	HEXAGON SOCKET HEAD CAP SCREW, ISO 4762, M10-1.5 X 50MM LONG, A4-70 STAINLESS STEEL, RoHS FASTENER	No
4	3000-M10C-40	16	Hexagon Nut, ISO4032, M10-1.5, Brass, RoHS Fastener	No
5	46-294167P96	12	bushing, leadboard vibration damping, xrmc gradient	No
6	5146169	12	M10x24mm Washer	No
7	5214524-6	15	washer, m10 large od nordlock, a4 stainless	No
8	5214524-8	8	washer, m12 large od nordlock, a4 stainless	No
9	5765900-50	8	HHCS M12X30 A4-80 SMARTBOLT	No
10	1000-M10C060-33	4	Hexagon Socket Head Cap Screw Partial threaded, ISO 4762, M10-1.5 X 60mm Long, Grade A4, AISI 316 Stainless Steel, RoHS Fastener	No

GE Healthcare**Renewal Parts and Tooling****6.8 FRU Assemblies – IRMw GRADIENT COIL**

Item №	Part Number	Qty	Description	GPRS
1	5887011	1	IRMw GRADIENT COIL FRU	Yes
2	5919559	1	Collector, IRMw Gradient Coil	Yes
3	5919560	1	COLLECTOR, ADDITIVE MANIFOLD FRU, IRMw	Yes
4	5919561	1	IRMw Busbar, FRU	Yes

6.9 FRU Assemblies – Artist EVO GRADIENT COIL

Item №	Part Number	Qty	Description	GPRS
1	5338453-2	1	Hex Spacer Mapping Fixture	Yes
2	5338453-3	1	Hex Spacer Mapping Fixture	Yes
3	5400303	1	Cover Mount Segment Molded	Yes
4	5401509	1	COVER FRAME SUPPORT	Yes

6.10 Miscellaneous parts for gradients**6.10.1 5803053, Revision 1: Nordlock Washer Kit, FRU**

Item №	Part Number	Qty	Description	GPRS
	5803053		Nordlock Washer Kit, FRU	Yes
1	5214524-6	20	washer, m10 large od nordlock, a4 stainless	No
2	5214524-8	20	washer, m12 large od nordlock, a4 stainless	No
3	2274863-2	6	flange nut, m10 x 1.5, a2, per din6923	No
4	5765900-50	6	HHCS M12X30 A4-80 SMARTBOLT	No

6.10.2 5426726, Revision 2: IBC Temperature Sensor Replacement Kit

Item №	Part Number	Qty	Description	GPRS
	5426726		IBC TEMPERATURE SENSOR REPLACEMENT KIT	Yes
1	5161627-2	6	TEMPERATURE SENSOR, 0-70 DEGREES C, AGED TO 105 DEGREES C, GE INDUSTRIAL SENSING PN A197-A-UBC85F103W	No
2	5128728-4	1.5	TUBE, BLACK POLYOLEFIN HEAT SHRINKABLE, .125 EXPANDED ID	No
3	5128728-3	1.5	TUBE, BLACK POLYOLEFIN HEAT SHRINKABLE, .093 EXPANDED ID	No
4	46-294167P91	.09	TAPE, ELECTRICAL, RUBBER SPLICING, .75 INCH X 30 FT. ROLL, SCOTCH 130C	No
5	5262823-5	.1	SOLDER, LEAD FREE, SN3.8AGCU0.7 INDALLY 241,.062 DIA WIRE	No

6.10.3 5414896, Revision 3, Hardware FRU Kit, FRMW Gradient Asm

Item №	Part Number	Qty	Description	GPRS
	5414896		HARDWARE FRU KIT, FRMW GRADIENT ASM*	Yes
1	5214523-13	6	SCREW, M10 X 50 CUP POINT HEX SOCKET SET, A2 OR A4 STAINLESS, WITH LOCKING COMPOUND	No
2	5214524-6	6	WASHER, M10 LARGE OD NORDLOCK, A4 STAINLESS	No
3	3000-M10C-40	6	HEXAGON NUT, ISO4032, M10-1.5, BRASS, ROHS FASTENER	No
4	5344896	6	SPACER, GRADIENT TO BUSBAR, DVW	No
5	5371820	6	LUG AND RING TERMINAL INSULATOR	No
6	2381283-3	12	TIE WRAP, NYLON, 50LB MIN, 2.0 OD	No
7	2274863-2	6	FLANGE NUT, M10 X 1.5, A2, PER DIN6923	No
8	5214523-15	6	SCREW, M10 X 60 CUP POINT HEX SOCKET SET WITH LOCKING COMPOUND, A2 OR A4 STAINLESS, DIN 916	No
9	5415261-2	1	LOCTITE 243 10 ML	No

*Suitable for XRMb, XRMw terminals repair. See SN DOC0957064 for details

GE Healthcare**Renewal Parts and Tooling****6.10.4 2241521, Revision 0: GRADIENT CABLE SPLICE KIT**

Item №	Part Number	Qty	Description	GPRS
	2241521		GRADIENT CABLE SPLICE KIT*	Yes
1	2218148	8	COLD SHRINK CONNECTOR INSULATOR	Yes
2	2218149	8	COPPER COMPRESSION TWO-WAY SPLICER	Yes

* Suitable for BRM/BRM2 and TRM

6.10.5 5849225, Revision 1: M12 Gradient Asm Hardware FRU Kit

Item №	Part Number	Qty	Description	GPRS
	5849225	1	M12 Gradient Asm Hardware FRU Kit	Yes
1	5765900-50	6	HHCS M12X30 A4-80 SMARTBOLT	No
2	5214524-8	6	washer, m12 large od nordlock, a4 stainless	No
3	5849267	6	M12 Spacer, Gradient to Busbar, DVw	No
4	5371820	6	Lug and Ring Terminal Insulator	No
5	3000-M10C-40	6	Hexagon Nut, ISO4032, M10-1.5, Brass, RoHS Fastener	No
6	5214524-6	6	washer, m10 large od nordlock, a4 stainless	No
7	2381283-3	12	TIE WRAP, NYLON, 50LB MIN, 2.0 OD	No

6.10.6 5867176, Revision 1: DV Gradient Cable Lug Kit - 1 AWG (Pink Lugs)

Item №	Part Number	Qty	Description	GPRS
	5867176	1	DV Gradient Cable Lug Kit - 1 AWG (Pink Lugs)	Yes
1	5799022	8	THOMAS and BETTS Terminal 54109	No
2	5799023	4	THOMAS and BETTS terminal 54155-TB	No
3	5826292	4	Copper Two-Hole Lug, Standard Barrel, Peep Hole, Max 35kV, Wire Size #1 Flex 225/24, 1/2 Inch Bolt Size, 1-3/4 Inch Hole Spacing, Tin Plated, Die Color Code Pink, Die Code 42	No
4	5368384	3	LUG .25 STUD 6AWG WIRE	No
5	5410522	1	lug .375 Stud 6 awg wire	No

6.10.7 5867177, Revision 1: DV Gradient Cable Lug Kit - 2/0 AWG (Orange Lugs)

Item №	Part Number	Qty	Description	GPRS
	5867177	1	DV Gradient Cable Lug Kit - 2/0 AWG (Orange Lugs)	Yes
1	5368413	8	LUG .375 STUD 3-0 AWG WIRE SIZE	No
2	5826504	4	Copper One-Hole Lug, Standard Barrel, Peep Hole, Max 35kV, Wire Size 3/0 AWG, 1/2 Inch Bolt Size, Tin Plated, Die Code 50, Die Color Code Orange	No
3	5826505	4	Copper Two-Hole Lug, Standard Barrel, Peep Hole, Max 35kV, Wire Size 3/0 AWG, 1/2 Inch Bolt Size, 1-3/4 Inch Hole Spacing, Tin Plated, Die Code 50, Die Color Code Orange	No
4	5368384	3	LUG .25 STUD 6AWG WIRE	No
5	5410522	1	lug .375 Stud 6 awg wire	No

GE Healthcare

Renewal Parts and Tooling

6.11 Tool Kits – Gradient

6.11.1 2164744-9, Revision 2: Gradient Insertion Lifting Tool for BRM, TRM, XRMB and XRMw

Item №	Part Number	Qty	Description	GPRS
	2164744-9		KIT, GRADIENT INSERTION LIFTING TOOL BRM,TRM,XRMB, XRMW	No
1	5191132	1	TOOL, GRADIENT INSTALLATION JACK	No
2	2284929	1	MALE TUBE ASM - INSERTION TOOL	No
3	2284930	1	FEMALE TUBE ASM - INSERTION TOOL	No
4	2284928-2	1	TOOL, SUPPORT PLATE, BRM, TRM, XRMB INSTALLATION	No
5	2172196	1	CRATE - BRM INSERTION TOOL	No
6	2121111	4	BODY COIL WHEEL HOUSING ASSEMBLY	Yes
7	2188986	1	ROPE AND SNAP ASSEMBLY	No
8	2326753	1	GRADIENT COIL LIFT KIT	Yes
9	2280510	1	BRM INSERTION TOOL IDENT LABEL	No
10	2187590	2	MOUNTING PLATE, CRM INSERTION TOOL	No
11	46-318508P59	10	HEX HEAD CAP SCREW, M10 X 80MM, FULL THREAD, SS	No
12	5308402	2	TOOL, TUBE GUIDE ROLLER ASSEMBLY, FIELD SERVICE	No
13	5161983	1	TOOL, SERVICE END INSTALLATION PLATE, DUAL USE TRM ,XRMB	No
14	5161985	1	TOOL, PATIENT END INSTALLATION PLATE, DUAL USE TRM , XRMB	No
15	5191626	1	TOOL, STANDOFF, GRADIENT INSTALLATION TUBE,	No
16	46-318508P20	14	HEX HEAD CAP SCREW, M10 X 25, FULL THREAD, STAINLESS STEEL	No
17	2180498	16	STUD, M10 X 120MM LONG, STAINLESS	Yes
18	2109875-4	40	NUT HEXAGON M10 A2-70	Yes
19	2109878-4	24	WASHER FLAT ISO7089 M10 A2	Yes
20	46-318508P25	10	HEX HEAD CAP SCREW, M10 X 50, FULL THREAD.	No
21	46-328430P8	8	WASHER PLAIN - NORMAL 13 MM 24 MM	No
22	46-252065P118	2	HAIRPIN COTTER, 1/8 DIA WIRE, SS	No
23	2116325-2	6	STUD SST M10 X 65	Yes
24	2224994	8	3/8 - 16 STEEL NUT WING	No
25	2224996	2	3/8 X 3/4 X 3~ PLAIN STEEL SPACER	No
26	2109866-24	10	SCREW HEXAGON HEAD 10 MM 50 MM	Yes
27	2315733	8	STAINLESS STEEL SLEEVE SPACER	No
28	46-260789P89	4	SCREW HEXAGON SOCKET 9.53 MM 38.1 MM .375-16UNC X 1.50	No
29	5303993-18	4	SCREW, SOCKET HEAD CAP, SS .38-16 X 4.0	No
30	5303994-9	4	SET SCREW, SLOTTED, SS M10 X 90	No
31	5265166	1	SLEEVE, INSERTION TUBE, THREAD PROTECTOR	No
32	5266408	1	TOOL, ROTATIONAL ALIGNMENT, XRMB GRADIENT	No
33	5262238	4	SPACER, GRADIENT INSTALLATION PLATE	No
34	5338222	3	XRMW INSTALL LINK PLATE A	No
35	5338223	4	XRMW INSTALL LINK PLATE B	No
36	5338224	2	XRMW INSTALL LINK PLATE C	No
37	5339882	2	ADAPTER PLATE, GRADIENT INSERTION TOOL, DVW	No
38	46-260950P79	4	CAP SCREW, HEX HD, .75-16UNF X 1.25 LG. STAINLESS STL	No
39	5342525	1	JACK HANDLE, NON MAGNETIC	No
40	5308861	3	LABEL, MAGNETIC MATERIALS DANGER	No
41	46-318508P19	32	HEX HEAD CAP SCREW, M10 X 20, FULL THREAD.	No
42	2269326-2	1	WARNING LABEL: FERROUS MATERIAL (SM)	No
43	5389303	1	TOOL, ROTATIONAL ALIGNMENT, XRMW GRADIENT	No
44	5694574	1	tool, spacer, gradient installation alignment block	No
45	5729022	2	SPACER GRADIENT INSTALL PLATE	No
46	5729027	1	GRADIENT INSTALL TUBE, SPACER	No
47	5729471	1	GRADIENT INSTALL PLATE BRIDGE	No
48	5694574-2	1	G3 GRADIENT INSTALL, ALIGNING TOOL	No
49	1000-M12C035-33	2	Hexagon Socket Head Cap Screw Full threaded, ISO 4762, M12-1.75 X 35mm Long, Grade A4, AISI 316 Stainless Steel, RoHS Fastener	No
50	3000-M12C-40	2	HEXAGON NUT STYLE 1, ISO 4032, M12-1.75 BRASS RoHS FASTENER	No
51	5766679	9	STUD, STAINLESS STEEL M10	No
52	1000-M10C025-33	4	HEXAGON SOCKET HEAD CAP SCREW, ISO 4762, M10-1.5 X 25MM LONG, A4-70 STAINLESS STEEL, RoHS FASTENER	No

GE Healthcare**Renewal Parts and Tooling**

6.11.2 Torque wrenches kits

Item №	Part Number	Qty	Description	GPRS
	5534134		Norbar Non-Magnetic Torque Wrench Kit (Adjustable 8-50 Nm)	No
	5790177		Non-Magnetic 80Nm Fixed Torque Wrench, Kit	No
	5534134-2		Norbar Non-Magnetic Torque Wrench Kit	No
	5109890-1		titanium socket, 3_8 inch drive x 3_8 inch hex bit	Yes

6.11.3 Short Rm, Blue Coil Cart, Cradle

Item №	Part Number	Qty	Description	GPRS
	2144093-2		Short Rm, Blue Coil Cart, Cradle	No
1	2134810-2		BRM/CRM/TRM MODIFIED CART--SHORT ROOM	No
2	2134810-3		ALUMINUM CRADLE FOR BRM CART	Yes

GE Healthcare**Renewal Parts and Tooling****Chapter 7 – MISCELLANEOUS ITEMS****7. Miscellaneous Items****7.1 Non-Magnetic/Non-Ferrous Hand Tools****7.1.1 2385097, Revision 0: FE Basic Nonmagnetic Tool Kit**

Item No	Part Number	Qty	Description	GPRS
	2385097		FE BASIC NON-MAGNETIC HAND TOOL KIT	No
1	2385098	Kit	AMPCO TOOL KIT	Yes
	2357186	1	WRENCH, 8~ CRESCENT	Yes
	2357187	1	PLIERS	Yes
	2357188	1	CUTTERS, DIAGONAL	Yes
	2357186-2	1	WRENCH, 10~ CRESCENT	Yes
	2357187-2	1	NEEDLE-NOSE PLIERS	Yes
	2357187-3	1	CHANNEL LOCK PLIERS	Yes
2	2385099	Kit	IMPERIX TOOL KIT	Yes
	46-320688P1	1	PHILLIPS SCREWDRIVER #0 X 3~ LG	Yes
	46-320688P2	1	PHILLIPS SCREWDRIVER #1 X 3~ LG	Yes
	46-320688P3	1	PHILLIPS SCREWDRIVER #2 X 6~ LG	Yes
	46-320688P4	1	PHILLIPS SCREWDRIVER #3 X 6~ LG	Yes
	46-320466P1	1	3/8 SPL DRIVE RATCHET	Yes
	46-320466P2	1	4~ EXTENSION	No
	46-320467P5	1	3/8 DRIVE X 9/16	Yes
	46-320468P12	1	3/8 DRIVE X 17MM	Yes
	46-320468P15	1	3/8 DRIVE X 17MM DEEP WELL	Yes
	46-320469P12	1	5/16 ALLEN WRENCH	Yes
	46-320469P8	1	5/32 ALLEN WRENCH	No
	46-320469P9	1	3/16 ALLEN WRENCH	Yes
	46-320470P10	1	6MM ALLEN WRENCH	No
	46-320470P13	1	9MM ALLEN WRENCH	No
	46-320470P6	1	4MM ALLEN WRENCH	No
	46-320470P8	1	5MM ALLEN WRENCH	No
	46-320472P5	1	OPEN END 8MM X 10MM	No
	46-320473P10	1	STD SCREWDRIVER 1/4 DIA X 5~ LG	Yes
	46-320473P13	1	STD. SCREWDRIVER 5/16 DIA X 6~ LG	Yes
	46-320473P18	1	STD. SCREWDRIVER 1/8 DIA X 8~ LG.	No
	46-320473P4	1	STD SCREWDRIVER 1/8 DIA X 4~ LG	Yes
	46-320473P7	1	STD. SCREWDRIVER 3/16 DIA X 4~ LG	Yes
	46-320686P4	1	COMBO 7/16 X 7/16	Yes
	46-320686P5	1	COMBO 1/2 X 1/2	Yes
	46-320686P6	1	COMBO 9/16 X 9/16	Yes
	46-320686P7	1	COMBO 5/8 X 5/8	Yes
	46-320686P8	1	COMBO 11/16 X 11/16	Yes
	46-320686P9	1	COMBO 3/4 X 3/4	Yes
	46-320687P6	1	COMBO 17MM	Yes
3	46-194451P16	1	15.5X11.25X2.5 WITH OUTSIDE EAPANDABLE POCKETS BLACK CORDURA CHICAGO CASE CD-48	Yes

GE Healthcare**Renewal Parts and Tooling****7.1.2 46-320273G4, Revision 0: Nonmagnetic Tool Kit**

Item №	Part Number	Qty	Description	GPRS
	46-320273G4		INCH/METRIC NON-MAGNETIC TOOL KIT	No
2	46-216784P12	1	TOOL CASE	Yes
3	2168869	1	46-320273G4 CASE LABEL	No
4	46-320270P1	1	FOAM BOTTOM	No
5	46-320718G1	1	METRIC BASE TOOLS ENCASED	No
6	46-320719G1	1	INCH BASE TOOL ENCASED	Yes
7	46-320287P2	1	PALLET	No
8	46-320288P1	1	SINGLE	No
9	46-320466P1	1	3/8 SPL DRIVE RATCHET	Yes
10	46-320466P2	1	4~ EXTENSION	No
11	46-320466P3	1	8~ EXTENSION	No
12	46-320466P4	1	3/8 DRIVE SPEED HANDLE	Yes
13	46-320466P5	1	3/8 DRIVE ~T~ HANDLE	No
14	46-320466P6	1	3/8 DRIVE UNIVERSAL JOINT	Yes
15	46-320473P1	1	STD. SCREWDRIVER 3/32 DIA X 3~ LG	Yes
17	46-320473P7	1	STD. SCREWDRIVER 3/16 DIA X 4~ LG	Yes
18	46-320473P10	1	STD SCREWDRIVER 1/4 DIA X 5~ LG	Yes
19	46-320473P13	1	STD. SCREWDRIVER 5/16 DIA X 6~ LG	Yes
20	46-320473P3	1	STD. SCREWDRIVER 1/8 DIA X 2~ LG	Yes
21	46-320473P18	1	STD. SCREWDRIVER 1/8 DIA X 8~ LG.	No
22	46-320688P3	1	PHILLIPS SCREWDRIVER #2 X 6~ LG	Yes
23	46-320474P2	1	ADJUSTABLE WRENCH 8~	No
24	46-320474P1	1	ADJUSTABLE WRENCH 10~	No
25	46-320475P1	1	SLIPJOINT 10~	Yes
26	46-320475P2	1	NEEDLENOSE	Yes
27	46-320475P3	1	SIDE CUTTER	No
28	46-320475P4	1	STRIPPER	Yes
29	46-320476P1	1	BALLPEEN 1LB.	Yes
30	46-320477P1	1	HALF ROUND 13~LG	No
31	46-320687P3	1	COMBO 9MM	Yes
32	46-320687P4	1	COMBO 11MM	Yes
33	46-320687P5	1	COMBO 14MM	Yes
34	46-320687P6	1	COMBO 17MM	Yes
35	46-320687P7	1	COMBO 19MM	Yes
36	46-320472P1	1	OPEN END 6MM X 7MM	Yes
37	46-320472P4	1	OPEN END 8MM X 9MM	Yes
38	46-320472P5	1	OPEN END 8MM X 10MM	Yes
39	46-320472P11	1	OPEN END 12MM X 13MM	Yes
40	46-320472P13	1	OPEN END 14MM X 15MM	Yes
41	46-320472P16	1	OPEN END 16MM X 18MM	Yes
42	46-320686P2	1	COMBO 5/16 X 5/16	Yes
43	46-320686P3	1	COMBO 3/8 X 3/8	Yes
44	46-320686P4	1	COMBO 7/16 X 7/16	Yes
45	46-320686P5	1	COMBO 1/2 X 1/2	Yes
46	46-320686P6	1	COMBO 9/16 X 9/16	Yes
47	46-320686P7	1	COMBO 5/8 X 5/8	Yes
48	46-320686P8	1	COMBO 11/16 X 11/16	Yes
49	46-320686P9	1	COMBO 3/4 X 3/4	Yes
50	46-320471P1	1	OPEN END 7/32 X 1/4	Yes
51	46-320471P2	1	OPEN END 1/4 X 5/16	Yes
52	46-320687P1	1	COMBO 6MM	Yes
53	46-320687P2	1	COMBO 7MM	Yes
54	46-328349P1	4	PLASTIC ARROW CLIP (.248 BODY X .54 LG)	No
55	46-320475P5	1	PLIERS, SLIP JOINT, 15~	Yes

GE Healthcare**Renewal Parts and Tooling****7.1.3 46-320718G1, Revision 0: Nonmagnetic Metric Tool Kit**

Item No	Part Number	Qty	Description	GPRS
	46-320718G1		METRIC BASE TOOLS ENCASED	No
1	46-320738P1	1	1 COMPARTMENT BOX, 14 X 6.25 X 2.5~ FLAMBEAU PRODUCTS MIDDLEFIELD OHIO PART V601	No
2	46-320804P1	1	FOAM INSERT	No
3	46-320806P1	1	FOAM INSERT	No
4	46-320809P1	1	LABEL	No
5	46-320468P3	1	3/8 DRIVE X 8MM	Yes
6	46-320468P4	1	3/8 DRIVE X 9MM	Yes
7	46-320468P5	1	3/8 DRIVE X 10MM	Yes
8	46-320468P6	1	3/8 DRIVE X 11MM	Yes
9	46-320468P7	1	3/8 DRIVE X 12MM	Yes
10	46-320468P8	1	3/8 DRIVE X 13MM	Yes
11	46-320468P9	1	3/8 DRIVE X 14MM	Yes
12	46-320468P10	1	3/8 DRIVE X 15MM	Yes
13	46-320468P11	1	3/8 DRIVE X 16MM	Yes
14	46-320468P12	1	3/8 DRIVE X 17MM	Yes
15	46-320468P13	1	3/8 DRIVE X 18MM	Yes
16	46-320468P14	1	3/8 DRIVE X 19MM	Yes
17	46-320808G2		ALLEN TOOL KIT - METRIC	No
	46-320807P2	1	POUCH	No
	46-320470P1	1	1.5MM ALLEN WRENCH	Yes
	46-320470P2	1	2MM ALLEN WRENCH	Yes
	46-320470P3	1	2.5MM ALLEN WRENCH	Yes
	46-320470P4	1	3MM ALLEN WRENCH	Yes
	46-320470P5	1	3MM SHORT ALLEN WRENCH	Yes
	46-320470P6	1	4MM ALLEN WRENCH	No
	46-320470P7	1	4MM SHORT ALLEN WRENCH	Yes
	46-320470P8	1	5MM ALLEN WRENCH	Yes
	46-320470P9	1	5MM SHORT ALLEN WRENCH	Yes
	46-320470P10	1	6MM ALLEN WRENCH	No
	46-320470P11	1	7MM ALLEN WRENCH	Yes
	46-320470P12	1	8MM ALLEN WRENCH	Yes
	46-320470P13	1	9MM ALLEN WRENCH	No
	46-320470P14	1	10MM ALLEN WRENCH	Yes
	46-320470P15	1	11MM ALLEN WRENCH	Yes
	46-320470P16	1	12MM ALLEN WRENCH	No
	46-320470P17	1	14MM ALLEN WRENCH	No
18	46-320948P1	1	CASE LABEL	No

GE Healthcare**Renewal Parts and Tooling****7.1.4 46-320719G1, Revision 0: Nonmagnetic Inch Tool Kit**

Item No	Part Number	Qty	Description	GPRS
	46-320719G1		INCH BASE TOOL ENCASED	No
1	46-320738P1	1	1 COMPARTMENT BOX, 14 X 6.25 X 2.5~ FLAMBEAU PRODUCTS MIDDLEFIELD OHIO PART V601	No
2	46-320804P1	1	FOAM INSERT	No
3	46-320805P1	1	FOAM INSERT	No
4	46-320810P1	1	LABEL	No
5	46-320467P2	1	3/8 DRIVE X 3/8	Yes
6	46-320467P3	1	3/8 DRIVE X 7/16	Yes
7	46-320467P4	1	3/8 DRIVE X 1/2	Yes
8	46-320467P5	1	3/8 DRIVE X 9/16	Yes
9	46-320467P6	1	3/8 DRIVE X 5/8	Yes
10	46-320467P7	1	3/8 DRIVE X 11/16	Yes
11	46-320467P12	1	3/8 DRIVE X 1 1/16	No
12	46-320467P13	1	3/8 DRIVE X 1 1/4	Yes
13	46-320808G1		ALLEN TOOL KIT - INCH	No
	46-320807P1	1	POUCH	No
	46-320469P1	1	0.50 ALLEN WRENCH	Yes
	46-320469P2	1	1/16 ALLEN WRENCH	Yes
	46-320469P3	1	5/64 ALLEN WRENCH	Yes
	46-320469P4	1	3/32 ALLEN WRENCH	No
	46-320469P5	1	7/64 ALLEN WRENCH	Yes
	46-320469P6	1	1/8 ALLEN WRENCH	Yes
	46-320469P7	1	9/64 ALLEN WRENCH	Yes
	46-320469P8	1	5/32 ALLEN WRENCH	No
	46-320469P9	1	3/16 ALLEN WRENCH	Yes
	46-320469P10	1	7/32 ALLEN WRENCH	No
	46-320469P11	1	1/4 ALLEN WRENCH	Yes
	46-320469P12	1	5/16 ALLEN WRENCH	Yes
	46-320469P13	1	3/8 ALLEN WRENCH	Yes
	46-320469P14	1	7/16 ALLEN WRENCH	Yes
	46-320469P15	1	1/2 ALLEN WRENCH	Yes
14	46-320467P1	1	3/8 X 5/16	Yes
15	46-320949P1	1	LABEL CASE	No

GE Healthcare**Renewal Parts and Tooling****7.2 FRU Assemblies, Parts and Kits – Miscellaneous****7.2.1 46-271137G1, Revision 2: Cryogen Safety Equipment Kit**

Item №	Part Number	Qty	Description	GPRS
	46-271137G1		SAFETY EQUIPMENT KIT CONSISTING OF HEAD GEAR, VISION SHIELD, SAFETY GLASSES, GLOVES AND CASE WITH FOAM	Yes
1	46-271135P10	1	HEAD GEAR WITH NO VISOR	Yes
2	46-271135P11	1	VISION SHIELD TO BE USED WITH HEAD GEAR	Yes
3	46-271135P12	1	SAFETY GLASSES	Yes
4	46-258217P2	1	LEATHER FOAM INSULATED GLOVES	Yes
5	46-301453P8	1	BLOW MOLDED, 17 X 13 X 7-3/8	No
6	46-306748P1	1	FOAM FOR CASE TOP, SAFETY EQUIPMENT KIT	No
7	46-306748P2	1	FOAM FOR CASE BOTTOM, SAFETY EQUIPMENT KIT	No
8	46-306710P1	1	LABEL FOR SAFETY EQUIPMENT KIT	No

7.2.2 Miscellaneous Safety Items

Item №	Part Number	Qty	Description	GPRS
	46-194427P320		PADLOCK,BRASS WITH BRASS SHACKLE MASTER LOCK #4B	Yes
	2387081		EHS BRASS TRANSITION LOCK MASTER LOCK 4KAB 3609 BRASS PADLOCK (535384)	Yes
	46-194427P322		LOCKOUT TAG, BLACK WHITE ON RED. PKG OF 25 GLAVIN SPECIALTY NO. GE GT409-2	Yes
	2393064		TAG LOTO TRANSITION	Yes
	46-194427P313		1 INCH MULTIPLE LOCKOUT DEVICE. PRINZING MODEL T210	Yes
	46-194427P321		PLUG LOCKOUT-- LOCKS OUT ANY PLUG UP TO 3~WIDE AND 5-7/8~LONG WITH A MAXIMUM CORD DIA.OF 1-1/4~. PRINZING # PL023	Yes
	46-194427P284		DIGITAL VOLT METER. VENDOR NAME: FLUKE VENDOR PART NO. 87	No

7.2.3 46-260888G4, Revision 1: LCC Magnet Leveling Kit

Item №	Part Number	Qty	Description	GPRS
	46-260888G4*		MAGNET LEVELING KIT- LCC	Yes
1	2213945	12	SHIM PLATE,LCC MAGNET	Yes
2	2213945-2	8	SHIM PLATE,LCC MAGNET	Yes
3	2180016	24	SHIMS, 6 X 6 ALUM.,.020 THK.-CONQUEST	No
4	2180016-2	16	SHIMS, 6 X 6 ALUM.,.032 THK.-CONQUEST	No
5	2180016-3	8	SHIMS, 6 X 6 ALUM.,.040 THK.-CONQUEST	No
6	2180016-4	4	SHIMS, 6 X 6 ALUM.,.063 THK.-CONQUEST	No

7.2.4 5819317, Revision 1: MAGNET LEVELING SHIM KIT (PM and AR magnet series)

Item №	Part Number	Qty	Description	GPRS
	5819317	1	MAGNET LEVELING SHIM KIT	Yes
1	5872197-05	8	SHIM PLATE, 0.5mm Thick	No
2	5872197-15	8	SHIM PLATE, 1.5mm Thick	No
3	5872197-30	4	SHIM PLATE, 3.0mm Thick	No

GE Healthcare**Renewal Parts and Tooling****7.2.5 5926419, Revision 3: Lifting Rail FRU Kit: LCC/450/450W magnet (R, RD, CR, HM)**

Item №	Part Number	Qty	Description	GPRS
	5926419	1	Lifting Rail FRU Kit: LCC/450/450W magnet (R, RD, CR, HM)	Yes
1	5924341	16	.50-13UNC x 2.00 Long, Grade 8, Yellow Zinc, Fully Threaded, Hex Head Screw, ASTM F1941	No
2	5924341-80	30	.50-13UNC, Yellow Zinc, Lock Washer, ASTM F1941	No
3	5924341-100	30	50-13, Yellow Zinc, Flat Washer, ASTM F1941	No
4	5924341-40	16	.50-13UNC, Grade 8, Yellow Zinc, Hex Head Nut, ASTM F1941	No
5	5343359	2	Lifting Beam Angle, 1.5T Magnets	No
6	5926638	1	Bolt Torque Laminated Placard	No

7.2.6 5926420, Revision 3: Lifting Rail FRU Kit: G3/750w magnet (W, WB, UA)

Item №	Part Number	Qty	Description	GPRS
	5926420	1	Lifting Rail FRU Kit: G3/750w magnet (W, WB, UA)	Yes
1	5924341-10	20	.75-16UNF x 2.50 Long, Grade 8, Yellow Zinc, Fully Threaded, Hex Head Screw, ASTM F1941	No
2	5924341-90	40	.75-16UNF, Yellow Zinc, Lock Washer, ASTM F1941	No
3	5924341-110	40	.75-16, Yellow Zinc, Flat Washer, ASTM F1941	No
4	5924341-50	20	.75-16UNF, Grade 8, Yellow Zinc, Hex Head Nut, ASTM F194	No
5	2341993	2	BEAM, G3 LIFTING	No
6	5926638	1	Bolt Torque Laminated Placard	No

7.2.7 5926421, Revision 3: Lifting Rail FRU Kit: IPM/Starter magnet (PM, AR)

Item №	Part Number	Qty	Description	GPRS
	5926421	1	Lifting Rail FRU Kit: IPM/Starter magnet (PM, AR)	Yes
1	5924341-30	20	M12-1.75 x 65mm Long, Class 10.9, Yellow Zinc, Fully Threaded, Hex Head Screw, ISO 4042	No
2	5924341-80	40	.50-13UNC, Yellow Zinc, Lock Washer, ASTM F1941	No
3	5924341-120	40	M12, Yellow Zinc, Flat Washer, ASTM F1941	No
4	5924341-70	20	M12-1.75, Class 10, Yellow Zinc, Hex Head Nut, ISO 4042	No
5	5809386	2	PLATFORM LIFTING RAIL	No
6	5926638	1	Bolt Torque Laminated Placard	No

7.2.8 5926709, Revision 3: Hardware Collector Kit, LCC, LCCW & DVW

Item №	Part Number	Qty	Description	GPRS
	5926709*		Hardware Collector Kit, LCC, LCCW & DVW	Yes
1	5924341-40	16	.50-13UNC, Grade 8, Yellow Zinc, Hex Head Nut, ASTM F1941	No
2	5924341	16	.50-13UNC x 2.00 Long, Grade 8, Yellow Zinc, Fully Threaded, Hex Head Screw, ASTM F1941	No
3	5924341-80	30	.50-13UNC, Yellow Zinc, Lock Washer, ASTM F1941	No
4	5924341-100	30	.50-13, Yellow Zinc, Flat Washer, ASTM F1941	No
5	5926638	1	Bolt Torque Laminated Placard	No

*This kit is strictly bolts and washers. This does NOT come with lifting beams. If lifting beams are required, please order part number 5926419 that DOES include lifting beams.

GE Healthcare**Renewal Parts and Tooling****7.2.9 5926710, Revision 3: Hardware Collector Kit, 3TLC and Corona**

Item №	Part Number	Qty	Description	GPRS
	5926710*		Hardware Collector Kit, 3TLC and Corona	Yes
1	5924341-110	40	.75-16, Yellow Zinc, Flat Washer, ASTM F1941	No
2	5924341-10	20	.75-16UNF x 2.50 Long, Grade 8, Yellow Zinc, Fully Threaded, Hex Head Screw, ASTM F1941	No
3	5924341-50	20	.75-16UNF, Grade 8, Yellow Zinc, Hex Head Nut, ASTM F1941	No
4	5924341-90	40	.75-16UNF, Yellow Zinc, Lock Washer, ASTM F1941	No
5	5926638	1	Bolt Torque Laminated Placard	No

* This kit is strictly bolts and washers. This does NOT come with lifting beams. If lifting beams are required, please order part number 5926420 that DOES include lifting beams.

7.2.10 5926711, Revision 3: Hardware Collector Kit, PM and AR

Item №	Part Number	Qty	Description	GPRS
	5926711*		Hardware Collector Kit, PM and AR	Yes
1	5924341-70	40	. M12-1.75, Class 10, Yellow Zinc, Hex Head Nut, ISO 4042	No
2	5924341-30	20	M12-1.75 x 65mm Long, Class 10.9, Yellow Zinc, Fully Threaded, Hex Head Screw, ISO 4042	No
3	5924341-80	20	.50-13UNC, Yellow Zinc, Lock Washer, ASTM F1941	No
4	5924341-120	40	M12, Yellow Zinc, Flat Washer, ASTM F1941	No
5	5926638	1	Bolt Torque Laminated Placard	No

*This kit is strictly bolts and washers. This does NOT come with lifting beams. If lifting beams are required, please order part number 5926421 that DOES include lifting beams.

7.2.11 5366120, Revision 4: Stud Replacement Kit, DVW

Item №	Part Number	Qty	Description	GPRS
	5366120		STUD REPLACEMENT KIT, DVW	Yes
1	5366178	1	STUD REPLACEMENT BLOCK, DVW CYLINDER	No
2	5366177	1	STUD REPLACEMENT BLOCK, DVW FLANGE	No
3	5370258	1	SPACER ISOLATION STUD REPAIR KIT	No
4	5370258-2	1	SPACER ISOLATION STUD REPAIR KIT	No
5	5368503	1	SPACING BAR, FOR SIDE ELECTRONICS, STUD REPLACEMENT	No
6	U1-SVC-520047	1	SAND PAPER, 60 GRIT, PK OF 10	No
7	46-252065P158	1	TACK CLOTH, 18 X 36 INCHES, TRIMACO SKU #10501 OR EQUAL	No
8	5265376	1	DISPENSER, EPX PLUS II, APPLICATOR	No
9	46-252065P154	3	MIXING NOZZLE 9742, 3M SCOTCH-WELD EPX APPLICATOR, 50 ML, PN 62974299356	No
10	46-2914151P64	1	LOCTITE H3000 STRUCTURAL ADHESIVE - SPEEDBONDER, COLOR: CREAM, 50ML DUO-PAK DISPENSER WITH MIXING NOZZLE.	No
11	5268439	10	SHIM, POWER CONNECTION LEADBOARD AXIAL ADJUSTMENT, XRMB GRADIENT	No
12	5214524-6	12	WASHER, M10 LARGE OD NORDLOCK, A4 STAINLESS	No

GE Healthcare**Renewal Parts and Tooling****7.2.12 5339638, Revision 2: IO Angle Bracket Asm**

Item No	Part Number	Qty	Description	GPRS
	5339638		IO ANGLE BRACKET ASM	Yes
1	5335556	1	IO ANGLE BRACKET	No
2	5110749-2	2	3M DUAL LOCK RECLOSABLE FASTENER SJ3560	No
3	5122507	1	BULKHEAD ADAPTER, FIBER OPTIC - LC TO LC DUPLEX, PLASTIC BODY, CERAMIC SLEEVE, 0.525IN SQ. PANEL OPENING: LCOM FOA-550A-OD1	Yes
4	5110749-3	1	VELCRO	No

7.2.13. 46-318565G1, Revision 4: Vacuum Tool Kit [Emergency Vacuum Port Operator]

Item No	Part Number	Qty	Description	GPRS
	46-318565G1		SIV VACUUM TOOL KIT.	No
1	46-260852G3	1	EMERGENCY EVACUATION RELEASE TOOL	Yes

7.2.14. Miscellaneous Vacuum Maintenance Equipment

Item No	Part Number	Qty	Description	GPRS
	46-252210P1	1	VALVE OPERATOR	No
	46-265273G1	1	HELIUM LEVEL METER (GE MAGNET) CONSISTING OF METER, CABLE ASM, CASE WITH FOAM	No
	46-265387G1	1	HELIUM FLOW GAUGE ASM	Yes
	46-260267P3	1	SEAL OFF VALVE OPERATOR WITH KF16 FLANGE	Yes

GE Healthcare**Renewal Parts and Tooling****7.3 FRU Parts common – Miscellaneous**

Item №	Part Number	Qty	Description	GPRS
	5155291	1	SERVICE PLATFORM, HD	Yes
	5451642	1	MAGNET SIDE PROTECTION PLATE ASSEMBLY	Yes
	5397201	1	MAGNET SIDE KIT PROTECTION PLATE 750W	Yes
	46-281387P13	1	BELLEVILLE WASHER, 0.375 IN. NOMINAL	Yes
	2237498	1	RELIEF VALVE, 0.50 IN. NPT, 5 PSIG	Yes
	2249367	1	PRESSURE GAUGE, 0-15 PSI, 4.5 IN DIAMETER	Yes
	46-260926P1	1	SWING CHECK VALVE, 0.75IN NPT	Yes
	2219776	1	STRAIN RELIEF, COMPRESSOR LINE	Yes
	2142318	1	INSTRUMENTATION LEAD COUPLING NUT	Yes
	2214062	1	832 SLEEVE and COLD HEAD CABLE, LCC	Yes
	2259837	1	CABLE ASM,COLD HEAD CONTROL,LCC	Yes
	2242756	1	RELIEF VALVE, .375 NPT X 7 PSIG, GEN	Yes
	5114908*	1	Relief Valve, SST	Yes

* Replacement for 2242756

7.3.1 FRU Parts for R Series Magnets - Miscellaneous

Item №	Part Number	Qty	Description	GPRS
	2173354	1	RETRACTABLE SHIM LEAD ASSEMBLY	Yes
	46-252209P1	1	.25 ORIFICE E .172 PLUG (Nupro valve)	Yes
	2206930	1	PORT CONNECTOR	Yes
	46-260898P3	1	MALE CONN., .25 ODT X .125 NPT BRASS	Yes
	2102831	1	LARGE SHIM LEAD CAP/SEAL ASSEMBLY	Yes
	2206093	1	PLUG, SWAGELOCK, WITH .09 HOLE	Yes
	5130298	1	INSTRUMENTATION & ELECTRICAL LEAD ASSEMBLY (G-10 SLEEVE ASSEMBLY)	Yes
	2204623	1	INSTRUMENT BOX ASSEMBLY	Yes
	46-318060P1	1	SAV-CON CONNECTOR	Yes
	2180818	1	BAFFLE ASM,LCC	Yes
	2237254	1	DRIP TROUGH WITH SPONGE	Yes
	46-271260G4	1	WIDE OPEN DOCK ASSEMBLY	Yes
	5160991	1	DOCK BRACKET ASSEMBLY	Yes
	2226264	1	WIDE OPEN DOCK MOUNT CNTRL BOX ASM	Yes
	2156980	1	SPIKE NOISE INSULATION KIT	No
	2217622	1	BRACKET KIT, SHIM LEAD MOUNTING,LCC [Mobile]	Yes
	2349529	1	RUPTURE DISK*	Yes
	2350415	1	DISK ASM, GRAPHITE BURST (replaces 2194545-7)	Yes
	2194545-7	1	¾" BURST DISC & GASKET ASSEMBLY, 10PSIG	No
	5890495	1	LCC Warmbore Passive Shim Removal/Repair Kit Hand tools	No
	2300866	1	CABLE, DIODE TEMP READOUT, SERVICE	Yes

*2349529 is 5/16" thick, 10 psig Burst disk, used on majority of the systems

7.3.2 FRU Parts for HM/UA Series Magnets– Miscellaneous

Item №	Part Number	Qty	Description	GPRS
	5264199	1	DVW BAFFLE ASSEMBLY	Yes
	5306594	1	INSTRUMENTATION ASSEMBLY, DVW	Yes
	5366120	1	STUD REPLACEMENT KIT	Yes
	5265523	1	G10 Sleeve Asm, 2.50, w 30 Pin Latching Lemo, HDw	Yes
	2237254	1	DRIP TROUGH WITH SPONGE	Yes
	5266592	1	Top Cover Plate Asm, HDw	Yes

GE Healthcare**Renewal Parts and Tooling****7.3.3 FRU Parts for W/WB Series Magnets - Miscellaneous**

Item №	Part Number	Qty	Description	GPRS
	2173354	1	RETRACTABLE SHIM LEAD ASSEMBLY	Yes
	46-252209P1	1	.25 ORIFICE E .172 PLUG (Nupro valve)	Yes
	2206930	1	PORT CONNECTOR	Yes
	46-260898P3	1	MALE CONN., .25 ODT X .125 NPT BRASS	Yes
	5130298	1	INSTRUMENTATION & ELECTRICAL LEAD ASSEMBLY (G-10 SLEEVE ASSEMBLY)	Yes
	2237254	1	DRIP TROUGH WITH SPONGE	Yes
	5127543	1	g3 sav-con assembly	Yes
	46-271260G4	1	WIDE OPEN DOCK ASSEMBLY	Yes
	5160991	1	DOCK BRACKET ASSEMBLY	Yes
	2226264-2	1	WIDE OPEN DOCK MOUNT CNTRL BOX ASM	Yes
	46-252209P1	1	NUPRO VALVE, 0.25 IN ORIFICE (PART OF SHIM LEAD PLUMBING)	Yes
	2214062	1	832 SLEEVE and COLD HEAD CABLE, LCC	Yes

7.3.4 FRU Parts for RD Series Magnets - Miscellaneous

Item №	Part Number	Qty	Description	GPRS
	5266592	1	Top Cover Plate Asm, HDw	Yes
	5723879	1	assembly, g10 sleeve with latching lemo, lcc-w no cc magnet	Yes
	5264199	1	DVW BAFFLE ASSEMBLY	Yes
	5306594	1	INSTRUMENTATION ASSEMBLY, DVW	Yes
	5555227	1	Assembly, Plumbing, PED Certified	No

7.3.5 5265652, Revision 4: Kit, Field Spares, HDv

Item №	Part Number	Qty	Description	GPRS
	5265652		Kit, Field Spares, HDv	No
1	46-281101P9	2	1.487 ID X .103 THK., O-RING #2-128	Yes
2	46-281101P5	2	5.737 ID X .103 THK, O-RING #2-162	Yes
3	46-281101P6	2	9.487 ID X .103 THK, O-RING #2-177	Yes
4	46-281101P10	2	6.237 ID X .103 THK., #2-164	Yes
5	46-294104P1	2	NUT, VACUUM QUICK COUPLING	Yes
6	46-294105P1	2	RING, QUICK COUPLING RETAINING	Yes
7	46-260389P1	4	SILICONE, PARKER#2-116	Yes
8	46-260389P2	2	SILICONE, PARKER #2-014	Yes
9	46-260902P5	1	PLUG, .25 OD TUBE FITTING, BRASS SWAGelok CAT #B-400-P	Yes
10	2203803	2	SPLIT O-RING	Yes
11	46-318068P1	6	BRONZE WOOL USED FOR RF SHIELDING CDA 649 BRONZE 4 TO 5 IN. WIDE X 12 IN. LONG MINIMUM.	Yes
12	46-252283P58	4	CABLE TIE WRAP, PANDUIT PLT7LH-L	Yes
13	46-260342P10	5	QUICK SEAL CONNECTOR O-RING SILICONE PARKER TYPE 2-112	Yes
14	46-252065P24	1	VACUUM GREASE,DOW CORNING	Yes
15	46-260389P4	1	O-RING, SILICONE, PARKER 2-205	Yes
16	46-252065P141	1	AEROQUIP O-RING, FLAT, RUBBER	Yes
17	46-252065P157	1	Aeroquip Gasket, Catalog Part Number 22008-4, .25 in size, Flat, Rubber	No

GE Healthcare**Renewal Parts and Tooling****7.3.6 46-294744G5, Revision 10: Field spare Kit LCC**

Item №	Part Number	Qty	Description	GPRS
	46-294744G5		FIELDSPARE KIT – LCC	No
1	46-281101P9	2	1.487 ID X .103 THK., O-RING #2-128	Yes
2	46-281101P5	2	5.737 ID X .103 THK, O-RING #2-162	Yes
3	46-281101P6	2	9.487 ID X .103 THK, O-RING #2-177	Yes
4	46-281101P10	2	6.237 ID X .103 THK., #2-164	Yes
5	46-294104P1	2	NUT, VACUUM QUICK COUPLING	Yes
6	46-294105P1	2	RING, QUICK COUPLING RETAINING	Yes
7	46-260389P1	4	SILICONE, PARKER#2-116	Yes
8	46-260389P2	2	SILICONE, PARKER #2-014	Yes
9	46-260902P5	1	PLUG, .25 OD TUBE FITTING, BRASS SWAGelok CAT #B-400-P	Yes
10	2203803	2	SPLIT O-RING	Yes
11	2218148	14	COLD SHRINK CONNECTOR INSULATOR	Yes
12	2218149	2	COPPER COMPRESSION TWO-WAY SPLICER	Yes
13	46-318068P1	6	BRONZE WOOL USED FOR RF SHIELDING CDA 649 BRONZE 4 TO 5 IN. WIDE X 12 IN. LONG MINIMUM.	Yes
14	46-252283P58	4	CABLE TIE WRAP, PANDUIT PLT7LH-L	Yes
15	46-260342P10	5	QUICK SEAL CONNECTOR O-RING SILICONE PARKER TYPE 2-112	Yes
16	46-252065P24	1	VACUUM GREASE,DOW CORNING	Yes
17	46-260389P4	1	O-RING, SILICONE, PARKER 2-205	Yes
18	46-252065P141	1	AEROQUIP O-RING, FLAT, RUBBER	Yes
19	46-252065P157	1	AEROQUIP GASKET, CATALOG PART NUMBER 22008-4, .25 IN SIZE, FLAT, RUBBER	No

GE Healthcare**Renewal Parts and Tooling****7.4 Helium Vent Kits****7.4.1 46-318057G6, Revision 5: LCC Vent Kit Assembly**

Item №	Part Number	Qty	Description	GPRS
	46-318057G6		VENT KIT ASSEMBLY, LCC	Yes
1	2258644	1	VENT, VERTICAL HELIUM	No
2	46-281101P2	1	8.237 ID X .103 THK #2-172	Yes
3	2262680	1	CLAMP RING	Yes
4	46-281971P2	1	VENTGLAS, 8.00 WIDE X 80.00 LONG	Yes
5	46-281972P2	4	CLAMP, AERO-SEAL #300152H, STAINLESS STEEL, BREEZE CORPORATION, INC.	Yes
6	46-318057P1	2	SCREW, #10-16 X 1.00 LG. SELF-DRILL, PAN HEAD, PHILLIPS DRIVE, STAINLESS STEEL, TYPE 304 OR 18-8, HARDENED POINT, DRILL POINT #3, MCMASTER-CARR PART # 90415A460 OR EQUIVALENT.	No
7	46-252635P9	6	PLAIN WASHER, .375 NOM., STAIN. STL.	Yes
8	46-281046P38	6	BRASS, HEX HEAD CAP SCREW .375-16 X 1.75 LG	No

7.4.2 46-318057G2, Revision 0: Vent Kit Assembly, Mobile Sites

Item №	Part Number	Qty	Description	GPRS
	46-318057G2		VENTGLAS KIT FOR CONQUEST MOBILE	Yes
1	46-281971P2	1	VENTGLAS, 8.00 WIDE X 80.00 LONG	Yes
2	46-281972P2	4	CLAMP, AERO-SEAL #300152H, STAINLESS STEEL, BREEZE CORPORATION, INC.	Yes
3	46-318057P1	2	SCREW, #10-16 X 1.00 LG. SELF-DRILL, PAN HEAD, PHILLIPS DRIVE, STAINLESS STEEL, TYPE 304 OR 18-8, HARDENED POINT, DRILL POINT #3, MCMASTER-CARR PART # 90415A460 OR EQUIVALENT.	No

GE Healthcare**Renewal Parts and Tooling****7.5 Tool Kits – Miscellaneous****7.5.1 2230707-2, Revision 0: Service Platform, Wide Open, Mobile**

Item №	Part Number	Qty	Description	GPRS
	2230707-2		LADDER PLATFORM,WIDE OPEN,MOBILE	No

7.5.2 46-294872G3, Revision 2: Sav-Con Kit

Item №	Part Number	Qty	Description	GPRS
	46-294872G3		KIT, SAV-CON	No
1	46-318612G51	1	CASE ASSEMBLY FOR SAV-CON TOOLS	No
2	46-294310G1	1	TOOL ASSEMBLY SAV-CON INSERT REMOVE	Yes
3	46-294292G1	1	TOOL ASSEMBLY SCREW INSERT REMOVE	No
4	46-294765G1	1	LEXAN TURRET COVER	Yes
5	46-294306G1	1	FLANGE ASSEMBLY SAV-CON INSERT REMOVE TOOL	No
6	46-318561P1	1	PENETRATION ADAPTER PLATE	Yes
7	46-281934P1	1	REMOVAL-INSTRUMENTATION CONNECTOR	Yes
8	46-294310G2	1	SAV-CON INSTALLATION TOOL FOR SIII / MAX.	No
9	46-318241G1	1	LEXAN PENETRATION COVER PLATE.	Yes
10	46-260963P2	1	FLANGE BLANK	No
11	2119965	1	BELLOWS COVER ASSY, VMX	No
12	2117683	1	COVER ASM, PENETRATION, VMX	Yes

7.5.3 5158913, Revision 1: Latching Lemo Tool Kit

Item №	Part Number	Qty	Description	GPRS
	5158913		LATCHING LEMO TOOL KIT	No
1	5145459	1	TOOL, LATCHING INSTRUMENTATION LEMO REMOVAL	No
2	46-328349P1	4	PLASTIC ARROW CLIP (.248 BODY X .54 LG)	No
3	5103265-5	1	CASE, BLACK PLASTIC, PLANO MOLDING MODEL NUMBER 1535-00 OR EQUIVALENT.	No
4	5158912	1	LABEL, LATCHING LEMO TOOL KIT	No
5	5158913-2	1	INSTRUCTIONS, LATCHING LEMO TOOL KIT	No

7.5.4 5105189, Revision 1: Kit, Instrumentation Connector Reinsertion

Item №	Part Number	Qty	Description	GPRS
	5105189		KIT, INSTRUMENTATION CONNECTOR REINSERTION	No
1	46-281934P1	1	REMOVAL-INSTRUMENTATION CONNECTOR	Yes
2	5105892	1	FIXTURE, CONNECTOR REINSERTION, LCC	No
3	2303292	1	CAP ASM, VIEWER	Yes
4	46-320466P5	1	3/8 DRIVE ~T~ HANDLE	No
5	46-320466P3	1	8~ EXTENSION	No
6	46-320467P3	1	3/8 DRIVE X 7/16	Yes
7	46-320467P5	1	3/8 DRIVE X 9/16	Yes
8	46-260950P32	3	CAP SCREW, .375-16 X 1.00 LONG	No
9	46-252365P27	3	BOLT, 3/8-16 X 4.00, HEX HD, S.S.	No
10	2275388-2	1	CASE	No
11	5114892	1	LEMO CONNECTOR RESEATING	No

GE Healthcare**Renewal Parts and Tooling**

7.5.5 2262769, Revision 1: Plumbing Asm, LCC Replacement

Item No	Part Number	Qty	Description	GPRS
	2262769		PLUMBING ASM, LCC REPLACEMENT	Yes
1	2225905	1	TUBE ASM, WARM VENT, LCC	No
2	46-260898P1	2	MALE CONN., .50 ODT X .50 NPT, BRASS	Yes
3	46-252065P19	0	TEFLON TAPE .50 WIDE	Yes
4	46-252205P1	1	.50 NPT BRASS, TEE	Yes
5	2237498	1	RELIEF VALVE - .50 NPT X 5 PSIG	Yes
6	46-252557P1	1	STREET ELBOW, .50 NPT, BRASS	No
7	2239615	1	TUBE, LCC, LOOP	No
8	46-252552P4	1	.50 NPT, BRASS, PIPE PLUG	Yes
9	46-252539P3	1	ELBOW, .50FNPT - .50 TUBE, BRASS	Yes
10	46-252251P3	1	MALE RUN TEE, .50 ODT X .50 NPT, BRASS	Yes
11	46-252250P6	1	.50 STN. STL. UNION ELBOW	No

GE Healthcare

Renewal Parts and Tooling

7.5.6 5555227, Revision 12: Assembly Plumbing, PED Certified

Item №	Part Number	Qty	Description	GPRS
	5555227		Assembly, Plumbing, PED Certified	No
1	2253306-2	1	PLENUM, CAST, LCC, PED CERTIFIED	No
2	5661398	1	VENT ADAPTER, CASTING	No
3	2245511	1	FLANGE, BURST DISK	Yes
4	2349529	1	RUPTURE DISK	Yes
5	2349683-3	2	Bolt, Hex Head, SST 0.375-16 x 1.875 LG	No
6	2349683	4	HEX HEAD BOLT, SST	No
7	46-252635P9	6	PLAIN WASHER, .375 NOM., STAIN. STL.	Yes
8	46-281162P11	4	.375 NOM SPRING	Yes
9	46-252532P2	4	.375-16UNC X .625 LG STAINLESS STL	Yes
10	46-281387P15	6	BELLEVILLE WASHER,SST,.75 OD, .375 HOLE,.040 THK, ASSOC SPRING PART# K0750-C-040-S	Yes
11	2349677	6	HEX NUT, STAINLESS STEEL	No
12	46-260196P1	1	LABEL, BURST DISC WARNING - R.H.	No
13	5746243	2	BALL VALVE ASSEMBLY	No
14	5561705	1	Pressure Transducer, Wika S-10	Yes
15	46-252065P19	1	TEFLON TAPE .50 WIDE	Yes
16	46-281101P1	1	4.239 ID X .070 THK #2-046	Yes
17	46-281101P8	1	4.739 ID X .070 THK, O-RING #2-048	Yes
18	46-252065P24	0.5	VACUUM GREASE,DOW CORNING	Yes
19	46-294151P8	0.5	BOSTICK NEVER-SEEZ,REGULAR GRADE,LUBRICATING AND ANTI-SEIZE COMPOUND.	Yes
20	5732393	1	SEAMLESS COPPER TUBING	No
21	5693499	1	SEAMLESS COPPER TUBING, A	No
22	5693500	1	SEAMLESS COPPER TUBING, B	No
23	5693067-3	1	7 PSI. RELIEF VALVE	Yes
24	5692885	1	NPT MALE ELBOW	No
25	2249367-2	1	PRESSURE GAUGE, 0-15 PSI, 3.50 DIAL	Yes
26	5692887	5	PORT CONNECTOR, .50 INCH	No
27	5693125-2	1	.75 INCH MALE NPT CONNECTOR TO MALE NPT	No
28	5693122	3	TUBE STUB TO MALE .50 INCH NPT ADAPTER	Yes
29	5692886	2	TUBE CAP	Yes
30	5692883-2	1	FEMALE NPT BRANCH TEE	No
31	5693141	2	.5 INCH UNION TEE	No
32	5692949	1	.5 INCH, FEMALE NPT RUN TEE	No
33	5693143	1	.5 INCH OD X .032 INCH WALL, COPPER TUBE	No
34	5692884	1	UNION CROSS FITTING	No
35	5719674	1	MALE NPT BFRANCH TEE	No
36	5692887-2	1	PORT CONNECTOR, .75 INCH	No
37	5693067-2	1	5 PSI. RELIEF VALVE	Yes
38	5693067	1	.5 PSI. RELIEF VALVE	Yes
39	5693180	1	.375-16 SWIVEL LEVELING FOOT	No
40	5555260	1	tube, copper, stock 0.75 inch od x 0.032 inch wall	No
41	5732394	1	SEAMLESS COPPER TUBING, STRAIGHT	No
42	5732392	1	0.5INCH UNION ELBOW	No

GE Healthcare**Renewal Parts and Tooling****7.5.7 5175534, Revision 3: Plumbing Assembly, 3.0T, SST**

Item №	Part Number	Qty	Description	GPRS
	5175534		Plumbing Assembly, 3.0T, SST	No
1	46-281101P1	1	4.239 ID X .070 THK #2-046	Yes
2	2253306-2	1	PLENUM, CAST, LCC, PED CERTIFIED	No
3	2349529	1	RUPTURE DISK	Yes
4	2245511	1	FLANGE, BURST DISK	Yes
5	2253376	1	VENT ADAPTER,CAST,LCC	No
6	46-281101P8	1	4.739 ID X .070 THK, O-RING #2-048	Yes
7	2349683	4	HEX HEAD BOLT, SST	No
8	46-281162P11	4	.375 NOM SPRING	Yes
9	46-252532P2	4	.375-16UNC X .625 LG STAINLESS STL	Yes
10	46-281387P15	6	BELLEVILLE WASHER,SST,.75 OD,.375 HOLE,.040 THK, ASSOC SPRING PART# K0750-C-040-S	Yes
11	46-252635P9	6	PLAIN WASHER, .375 NOM., STAIN. STL.	Yes
12	2349677	6	HEX NUT, STAINLESS STEEL	No
13	2299843-2	1	PRESSURE TRANSDUCER, SUNX, FOR LCC	Yes
14	46-252065P19	0.05	TEFLON TAPE .50 WIDE	Yes
15	46-294151P8	0.05	BOSTICK NEVER-SEEZ,REGULAR GRADE,LUBRICATING AND ANTI-SEIZE COMPOUND.	Yes
16	46-252065P24	0.05	VACUUM GREASE,DOW CORNING	Yes
17	2204235	1	SUPPORT,VNET ADAPTER,LCC	No
18	46-252610P4	2	HEX SOC HD, .250-20UNC X .750 LG SS	Yes
19	46-260196P1	1	LABEL, BURST DISC WARNING - R.H.	No
20	46-252223P1	1	CRYOGENIC BALL VALVE, .50 NPT	Yes
21	5139403	1	Relief Valve, .50DT x 5 psig	Yes
22	5139414	1	Relief Valve, .50DT x 7psig	Yes
23	5139415	1	Pressure Gauge, .250DT, 0-15 PSI, 4.0 Dial	No
24	5139566	2	Elbow,.5NPT X .50DT	No
25	5139404	1	Hex Manifold	Yes
26	5139535	1	Tube,L-shape,Short,.50DT	No
27	5182517	1	Female Run Tee	No
28	5181099	1	Tube, U-Shape, Long	No
29	5181275	1	Tube, U-Shape, Short	No
30	5141944	1	Tube, .5 ODT	No
31	5178578	1	Check Valve, SST	Yes
32	46-252552P1	1	.25 NPT BRASS	Yes
33	5173282	2	Male Connector, SST	No
34	5173944	1	Male Branch Tee, SST	No
35	2349683-3	2	Bolt, Hex Head, SST 0.375-16 x 1.875 LG	No

GE Healthcare**Renewal Parts and Tooling****7.5.8 2293050, Revision 0: GAUSS METER W/PROBE, CABLE, CASE**

Item №	Part Number	Qty	Description	GPRS
	2293050	1	GAUSS METER W/PROBE, CABLE, CASE	No
1	2293044	1	GAUSS METER	No
2	2293045	1	Probe and cable	No
3	2293046	1	OPEN	No
4	2293047	1	CASE, FOAM, LABEL	No

7.6 MCATs – Miscellaneous**7.6.1 M1060SR, Revision 0: 2.5M Ceiling Ht Kit (M1060CH for China)**

Item №	Part Number	Qty	Description	GPRS
	M1060SR		SHIM LEAD ASM, RETRACTABLE, LOW CEILING HEIGHT	No
1	2287029	2	LEAD, FLEXIBLE CURRENT, 750A	Yes
2	2276755	1	FIXTURE,RAMP LEAD, SHORT	Yes
3	2285073	1	SHIM LEAD ASM,RETRACTABLE,LCC CEILING HEIGHT	No

7.6.2 M7000GM, Revision 2: Ceiling Height Kit-passive

Item №	Part Number	Qty	Description	GPRS
	M7000GM		Ceiling Height Kit-passive	No
1	2287029	2	LEAD, FLEXIBLE CURRENT, 750A	Yes
2	5323611	1	Ramp Lead Hold Tool, HDw	No