



PM Report Schedule C

LightSpeed VCT and Discovery CT750 HD

Customer Name: _____ Service Record #: _____

System ID/ Serial Number: _____ Customer Equipment ID: _____

Service Manual Document Number: _____ Software Version: _____

PM Frequency: _____ months

System Type (Select:)

PRELIMINARY TASKS

Last PM Inspection

Verify last PM Inspection

Last Inspection Date :

Last Schedule Completed:

List any Schedule C tasks that were completed earlier

PM Schedule Sys C

Review Error Logs

Check the tasks you completed and list any follow-up work scheduled.

Task	Completed	Comments/Follow-up
Table Hydraulic Errors	☐	
System Application Stop and Start Errors	☐	
Linux Reboot and Shutdown Errors	☐	
Record Gantry Rotations	0	

Initial Procedures

Check Temperature and Humidity

Record the temperature and humidity. Indicate if these are within spec.

Room	Temperature	Humidity	Pass/Fail	Comments/Follow-up
Scan room				
Control Room				

**Console**

Task	Completed	Comments/Follow-up
Console Monitor Cleaning & Inspection	☐	

General Gantry Inspection

Task	Pass/Fail	Comments/Follow-up
Listen for Unusual Noises		
Test Display Lights		
Test E-stops		
Test Patient Tilt Sensors		

Gantry Power Off (Covers Off)**General Tasks**

Task	Pass/Fail	Comments/Follow-up
Emergency Off Switch		
UPS Power Off Inspection		
Gantry Cover Cable Inspection		

General Console Cleaning & Inspection

Task	Pass/Fail	Comments/Follow-up
Clean Component Filters, Fans & Grills		
Clean, Inspect, Replace Desktop Components		
Visual Fan Inspection		
Seismic Anchor Check (if necessary)		

Slip Ring & Brush Block Inspection and Maintenance

Task	Pass/Fail	Comments/Follow-up
Remove Slip Ring Brush Debris		
Slip Ring Tracks Inspected		
Clean Brush Blocks		
Brush Tip Inspection		



Gantry Filter Cleaning

Task	Pass/Fail	Comments/Follow-up
Clean DAS Detector Plenum Filter		
Clean Detector Face Plate		
Clean Gantry Heater Filter		
Clean Top Cover Fans		
Inspect Rotating Path Components & Cables		

Encoder

Task	Pass/Fail	Comments/Follow-up
Inspect Rotational Encoder Wheel		
Encoder Replacement Date last replaced:		

Power On (Covers Off)

Tube Heat Exchanger & Pump

Task	Select	Comments/Follow-up
Inspect/Clean Tube Heat Exchanger		
Inspect/Clean JEDI Inverter Fan		
Tube Heat Exchanger Pump Inspection Mfg date/last replacement		

General Tasks

Task	Select	Comments/Follow-up
UPS Powered and Operational		
Plenum Fans (5) Operational		
Grease Gantry Main Bearing No. Revs since last greasing		

Power On (Covers On)

General Table Cleaning & Inspection

Task	Select	Comments/Follow-up
Inspect/Clean Table Pan		
Inspect Table anchors		
Hydraulic fluid leak Inspection		
Inspect Drive Belts		
Clean Roller		



Grease Table Rails		
Inspect Bearing Retainer		
Grease Table Leg Bearing		
Inspect Flat Tabletop Mount		
Inspect Table Covers		
Inspect Cradle Accessory Attachment (Note in comments if anything is replaced)		
Check Tape Switches		
Functional Checks Table Height Min Max Gantry Tilt S30 Min Max Gantry Tilt I30 Min Max		

General Tasks

Task	Select	Comments/Follow-up
Scan Window Inspection/Replacement		
Inspect Rotational Clearances		
Confirm Save State Completed	☐	

Finalization**System Scanning Test**

Task	QA Phantom Serial #	Location	Gantry Serial #
Scans	Scan Failures	Artifacts	Artifact Types
Series 1 (Scout/Auto Voice)			
Series 2 (Axial/Tilt)			
Series 3 (Helical/Auto Voice)			

General Tasks

Task	Select	Comments/Follow-up
X-ray On Indicators		
Scan Control Pushbuttons Test		
	Completed	
Site Log Updated	☐	
PM Paper work completed for this site	☐	

Is follow-up work required?



GE Healthcare

Proprietary to General Electric
This form is based on: **DOC0527863, Rev 1**
LightSpeed VCT and Discovery CT750 HD - PM Report Schedule C

System ID:

SR#:

Date:

Service Scheduled Date:

☐ This inspection passed and meets all GE Healthcare PM specifications

Comments:

Characters Remaining:

Tools Used:

Description	Serial Number	Bar Code	Cal Due Date	
				Add Tool Remove Tool

GE Representative:

FSEs completing this PM

FSE Name: _____

SSO #: _____

FSE SR #: _____

FSE Signature Data...

Customer Signature: ☐ Yes ☐ No



GE Healthcare

Proprietary to General Electric

This form is based on: **DOC0527863, Rev 1**

LightSpeed VCT and Discovery CT750 HD - PM Report Schedule C

System ID:

SR#:

Date:

☐ Final PM Complete

Date

Completed: _____