

Chapter 2

Table



MARK	NAME	QTY	FRU	DESCRIPTION
3	2115996-3	1.0	Yes	AXIAL HEADHOLDER ASM WITH RATING PLATE AND SELF-CATCHING LATCH
4	2115990-3	1.0	Yes	LCC HSA METALLESS SUPINE CORNAL HEADHOLDER ASM WITH RATING PLATE AND SELF-CATCHING LATCH
5	2115993-2	1.0	Yes	LCC METALLESS CRADLE EXTENDER HSA METALLESS FOOT EXTENDER WITH RATING PLATE AND SELF-CATCHING LATCH
8	2331933	1.0	Yes	METALLESS COMPATIBLE PHANTOM HOLDER
1	2269165	1.0	No	TABLE ASSEMBLY FOR H2 LCC {For further detail, see Section 1.0 , beginning on page 54 .}
n/a	2271242	—	No	TABLE ASSEMBLY FOR H2 HCC {For further detail, see Section 2.0 , beginning on page 110 .}
2	2330508	1.0	No	FOOT SWITCH ASSEMBLY {For further detail, see Section 3.0 , beginning on page 160 .}
n/a	46-296790G3	—	No	FOOT SWITCH COVER ASSEMBLY {For further detail, see Section 4.0 , beginning on page 162 .}

Table 2-1 Table w/ Footswitch (2276776-2 BOM, rev 31)

Section 1.0 — Table Assembly - LCC (2269165)





Figure 2-2 Table Assembly - LCC (2269165, sheet 1, rev 3)






Figure 2-4 Table Assembly - LCC (2269165, sheet 3, rev 0)

MARK	NAME	QTY	FRU	DESCRIPTION
4	2261682	1.0	Yes	BRKT,ACTUAT CVR MT,(REF 46-296591P1)
5	2265868	1.0	Yes	BRACKET (REF 46-296820P1)
7	2271244	1.0	Yes	RIGHT BASE COVER ASM FOR H2 TABLE {For further detail, see Section 1.2 , beginning on page 62 .}
8	2271247	1.0	Yes	LEFT BASE COVER ASM FOR H2 TABLE {For further detail, see Section 1.3 , beginning on page 64 .}
14	2261675	1.0	Yes	ENCODER COVER, (REF 46-297897P1)
15	2265857	1.0	Yes	COVER (REF 46-297374P1)
16	2261683	1.0	Yes	COVER, ACTUATOR, (REF 46-296592P1)
17	2372371	1.0	Yes	COVER TABLE TOP LEFT {For further detail, see Section 1.4 , beginning on page 66 .}
18	2372372	1.0	Yes	COVER TABLE TOP RIGHT {For further detail, see Section 1.5 , beginning on page 67 .}
19	46-296594P1	1.0	Yes	DRIVE COVER
21A	46-296679P2	1.0	Yes	ZEUS TABLE, REAR END COVER WITHOUT SWITCH OPENINGS
22	2261708	2.0	Yes	COVER GUARD, RAILS,(REF 46-296742P1)
23	46-264368G3	1.0	Yes	ENHANCED TABLE CONTROL HELIOS 2
28	46-297061P1	1.0	Yes	FISH PAPER
29	2265851	2.0	Yes	BRACKET, LR SIDE PANEL (REF 2136285)
30	2265850	2.0	Yes	BRACKET, LF SIDE PANEL (REF 2136284)
31	2261669	2.0	Yes	BRKT,U-FRONT SIDE PNL, (REF 2138103)
32	2261668	2.0	Yes	BRKT, U-REAR SIDE PNL, (REF 2138101)
33	2265870	4.0	Yes	NUT PLATE (REF 46-297171P1)
35	2261688	1.0	Yes	BRKT, LWR COV SUP (REF 46-297324P1)
36	46-297329P2	2.0	Yes	1/4 TURN RECEPTACLE - FIXED .128~ DIA. MOUNTING HOLES
38	2329950	1.0	Yes	CHANNEL,COVER SUPPORT
43	2146815-2	26.0	Yes	10-32 X 3/8 HEX SOC BUTTON HD SST MA
45	46-208562P2	1.0	Yes	SCREW SCREW OLD STL F70B5A2 ZINC PL HT TR FOR HYD EMB

Table 2-2 Table Assembly - LCC (2269165 BOM, rev 13)

MARK	NAME	QTY	FRU	DESCRIPTION
46	46-208562P3	4.0	Yes	SCREW SCREW OLD STL F70B5A2 ZINC PL HT TR FOR HYD EMB
47	46-208562P20	14.0	Yes	SCREW SCREW OLD STL F70B5A2 ZINC PL HT TR FOR HYD EMB N646P16010
49	46-170497P38	14.0	Yes	SCREW SCREW OLD HD STL F70B5A2 ZINC PL HT TR FOR HYD EMB N647P13008
52	2265882	8.0	Yes	SIDE PIVOT BUSHING (REF 46-297552P1)
55	46-297155G1	4.0	Yes	ZEUS TABLE, 9~ LONG CABLE
56	46-170012P61	2.0	Yes	0006-32 HEX .313 NUT HT. .178 STL ZINC PL ESNA 21NM-62 OR AERO 21NM-0632
59	2261690	2.0	Yes	HEX STDF,9.5MM TAP,(REF 46-297596P1)
61	46-170498P23	6.0	Yes	SCREW SCREW OLD F70B5A2 ZINC PL HT TR FOR HYD EMB
62	46-221523P10	6.0	Yes	HOLE DIA 0.437 MAT~L THK .125 MAX BUTTON
63	46-170012P63	2.0	Yes	0010-32 HEX .375 NUT HT .239 STL ZINC PL ESNA 21NM-02
64	46-220186P20	4.0	Yes	1/8~ 120 DEGREE COUNTERSUNK POP RIVET .129-.133,GRIP RANGE .126-.187 USM CORPORATION PART # AK4385
71	46-297695P1	2.0	Yes	COTTER PIN
72	46-297678P1	1.0	Yes	HINGES
73	2265878	1.0	Yes	PIN (REF 46-297676P1)
74	46-170490P50	1.0	Yes	SPRING-INCH EXTENSION 4 INCHES .65 I FREE LG 4.00 MUSIC WIRE RATE-LB/IN 1.2
78	46-197018P1	1.0	Yes	DANGER HI-VOLTAGE
80	46-220181P5	10.0	Yes	WASHER #10 (.190) SCREW SIZE, ID .194-.202 OD .334, WIDTH .062 THK .047, STEEL ZINC PLT
81	12-276	4.0	Yes	ID 0.2040 OD 0.5630 THK .0478 STL F70B5A
85	2146815-2	18.0	Yes	10-32 X 3/8 HEX SOC BUTTON HD SST MA
86	2265863	3.0	Yes	SHIM 0.80MM (REF 46-196354P1)
87	2265863-2	3.0	Yes	SHIM 0.40MM (REF 46-196354P4)
88	46-170106P1	0.0	Yes	DBL COATED FILM .007 THK 1.0W
92	1000904P358	1.0	Yes	ID 0.380 OD 1.000 THK .094 B4B33A F70D9 .0001/.0002
98	2265861-5	1.0	Yes	SPACER PLT 1.60MM (REF 46-327022P5)
100	46-208561P147	5.0	Yes	SCREW SCREW OLD F70BIAR ZINC PL HT TR FOR HYD EMB
103	2265845	5.0	Yes	PAD, LEVELING (REF 2105872)

Table 2-2 Table Assembly - LCC (2269165 BOM, rev 13) (Continued)

MARK	NAME	QTY	FRU	DESCRIPTION
105	2224296-3	1.0	Yes	VE-GE, 8MB, 2FL, GE MED, CT
108	2115992-5	1.0	Yes	METALLESS CRADLE YMS METALLESS CRADLE WITH RATING PLT AND KEEPER_FOR SELF-CATCHING LATCH
109A	2117299	1.0	Yes	METALLESS CRADLE COMPATIBLE FRONT END COVER FOR HSA TABLE, BACKWARD COMPATIBLE WITH ALL ZEUS TABLES
115	46-265874P1	1.0	Yes	50 OHM IMPEDANCE BNC FEED-THRU TERM 1ST END FEMALE-2ND END MALE
117	2117390-2	1.0	Yes	CE Compliance and WEEE Labels
118	2132620	1.0	Yes	EMC LABEL CLASS A GROUP 1
119	2103580-12	8.0	Yes	SCREW COUNTERSUNK HEA 6 MM 12 MM CAP SCREWS BLACK OXIDE, STEEL
120	46-328417P12	8.0	Yes	SCREW HEXAGON SOCKET 6 MM 12 MM ZINC PLATED, STEEL
123	46-233422P1	1.0	Yes	BNC JACK-PLUG, RIGHT ANGLE ADAPTER. 50 OHM IMPEDANCE, TFE INSULATED. (1.06~ X .56~ LENGTHS, .89~ DIA) SAME AS 215-070221 AND 215-070431
124	46-208922P5	1.0	Yes	0.5~ LONG
125	46-208984P17	3.0	Yes	BRASS STANDOFF MALE/FEMALE, 6-32 MALE THREAD, 6-32 THREADED HOLE, .25 HEX, 1.12 LONG
128	2271251	1.0	Yes	RIBBON JUMPER CABLE, ETC TO ETC I/F (CRADLE DRIVE CONTROL) FOR H2 TABLE
129	2271253	1.0	Yes	RIBBON CABLE JUMPER, ETC TO ETC I/F (FOOT SWITCHES) FOR H2 TABLE
130	2271254	1.0	Yes	RIBBON CABLE JUMPER, ETC TO ETC I/F (TILT & ELEVATION) FOR H2 TABLE
131	2271255	1.0	Yes	CABLE JUMPER, ETC TO ETC I/F (ETCDATA) FOR H2 TABLE
132A	2336447	1.0	Yes	H2 TABLE ETC/IF BD ASSY
134	46-303268P1	1.0	Yes	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH
137	P9140ET	1.0	Yes	NAME PLATE
6	2269167-2	1.0	No	SIDE PANEL KIT FOR H2 LCC CATHAY TBL {For further detail, see Section 1.1 , beginning on page 61 .}
13	2261673	2.0	No	BKT, REAR COVER, (REF 46-296832P1)
91	2132175	0.0	No	TABLE DOCUMENT COLLECTOR
93	2265862	1.0	No	BOLT, SHIPPING (REF 2184717)
94	2265861	0.0	No	SPACER PLT 0.50MM (REF 46-327022P1)
95	2265861-2	0.0	No	SPACER PLT 0.80MM (REF 46-327022P2)

Table 2-2 Table Assembly - LCC (2269165 BOM, rev 13) (Continued)

MARK	NAME	QTY	FRU	DESCRIPTION
96	2265861-3	0.0	No	SPACER PLT 1.00MM (REF 46-327022P3)
97	2265861-4	0.0	No	SPACER PLT 1.20MM (REF 46-327022P4)
101	2265840	5.0	No	ADJUSTER, LEVELING PAD (REF 2105871)
102	2265846	5.0	No	LOCK RING (REF 2106207) DIMENSIONS 1.5-12 THDS, 0.25 THK
113	2265849	1.0	No	BRACKET, ETC BRD COVER (REF 2117930)
116	2395282	1.0	No	ETL LABEL
121	2269166	1.0	No	SUBASSEMBLY FOR H2 LCC CATHAY TABLE {For further detail, see Section 1.7 , beginning on page 70 .}
122	2384711-2	1.0	No	RATING PLATE
126	2148658	1.0	No	ESD GROUND STRAP, BRAIDED FLAT WIRE

Table 2-2 Table Assembly - LCC (2269165 BOM, rev 13) (Continued)

1.1 Side Panel Kit (2269167-2)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
9	2259861-2	1.000	Yes	SIDE PANEL,RT REAR, HELIO
10	2259862-2	1.000	Yes	SIDE PANEL,LEFT REAR, HEL
11	2259863-2	1.000	Yes	SIDE PANEL,RIGHT FRONT, H
12	2259864-2	1.000	Yes	SIDE PANEL,LEFT FRONT, HE

Table 2-3 Side Panel Kit (2269167-2 BOM, rev 0)

FRU = YES

Figure 2-5 Right Base Cover Assembly (2271244, sheet 1, rev1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2271257	1.000	No	RIGHT BASE COVER FOR H2 TABLE
5	2265869-2	1.000	No	RIGHT FILLER COVER (REF 4
6	2153420	1.000	No	SPRING-INCH TORSION
7	46-297328P1	4.000	No	1/4 TURN STUD AND TWO GROMMETS
8	2329955	1.000	No	EJECTING BLADE
9	46-220186P6	2.000	No	DIA .125 HOLE .129-.133 -.390 LG
10	12-297	2.000	No	ID 0.140 OD 0.312 THK .032 STL
13	46-170024P60	1.000	No	SCREW-INCH
14	46-170012P62	1.000	No	0008-32 HEX .343 NUT HT .229

Table 2-4 Right Base Cover Assembly (2271244 BOM, rev 3)

1.3 Left Base Cover Assembly (2271247)

FRU = YES




Figure 2-6 Left Base Cover Assembly (2271247, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2271259	1.000	No	LEFT BASE COVER FOR H2 TABLE
4	2265869	1.000	No	LEFT FILLER COVER (REF 46
6	2153420	1.000	No	SPRING-INCH TORSION
7	46-297328P1	4.000	No	1/4 TURN STUD AND TWO GROMMETS
8	2329955	1.000	No	EJECTING BLADE
9	46-220186P6	2.000	No	DIA .125 HOLE .129-.133 -.390 LG
10	12-297	2.000	No	ID 0.140 OD 0.312 THK .032 STL
13	46-170024P60	1.000	No	SCREW-INCH
14	46-170012P62	1.000	No	0008-32 HEX .343 NUT HT .229

Table 2-5 Left Base Cover Assembly (2271247 BOM, rev 3)

1.4 Cover Table Top Left (2372371)

FRU = YES

MARK	NAME	QTY	FRU	DESCRIPTION
3	46-208747P1	4.0	Yes	ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)
4	46-297698P3	1.0	Yes	27~ HIGH PROFILE BEAD, RIBBON SWITCH
5	46-297329P1	2.0	Yes	1/4 TURN RECEPTACLE WITH TWO SPACERS FLOATING (.074~ ALL DIRECTIONS) .261~ DIA. MOUNTING HOLES
6	46-220186P20	4.0	Yes	1/8~ 120 DEGREE COUNTERSUNK POP RIVET .129-.133,GRIP RANGE .126-.187 USM CORPORATION PART # AK4385
7	46-297721P1	1.0	Yes	4-POS PLUG HOUSING (DRAWER CONNECT). PLUG HOUSING, DRAWER CONNECTOR. (DWG REV 7/92, JFK - OK).
8	46-158905P16	2.0	Yes	MOUNTING SCREW FOR DRAWER CONNECTOR. SAE 6-32 THD, .966~-L, STAINLESS STL MOUNTING SCREW FOR FLOATING DRAWER C ONNECTORS. SAE 6-32 THD. (DWG REV 7/92, JFK - OK).
9	46-297328P2	2.0	Yes	1/4 TURN STUD AND TWO GROMMETS STUD LENGTHS; A=.500~ B=.330~
10	46-208758P1	4.0	Yes	3.62~ X .094~ CABLE TIE 1/16~ TO 5/8~ WIRE BUNDLE RANGE 18LB TENSILE STRENGTH, NYLON MATR'L.
11	2372346	1.0	Yes	CHANNEL TAPESWITCH
14	2199652	1.0	Yes	P9235LZ LABEL
2	2372348	1.0	No	COVER, TOP LEFT
12	2363367	6.0	No	VINYL CAP 1/8 INCH
13	2363369	6.0	No	T SHAPE HEAD STEEL CLIP

Table 2-6 Cover Table Top Left (2372371 BOM, rev 1)

1.5 Cover Table Top Right (2372372)

FRU = YES

MARK	NAME	QTY	FRU	DESCRIPTION
3	46-208747P1	4.0	Yes	ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)
4	46-297698P3	1.0	Yes	27~ HIGH PROFILE BEAD, RIBBON SWITCH
5	46-297329P1	2.0	Yes	1/4 TURN RECEPTACLE WITH TWO SPACERS FLOATING (.074~ ALL DIRECTIONS) .261~ DIA. MOUNTING HOLES
6	46-220186P20	4.0	Yes	1/8~ 120 DEGREE COUNTERSUNK POP RIVET .129-.133,GRIP RANGE .126-.187 USM CORPORATION PART # AK4385
7	46-297721P1	1.0	Yes	4-POS PLUG HOUSING (DRAWER CONNECT). PLUG HOUSING, DRAWER CONNECTOR. (DWG REV 7/92, JFK - OK).
8	46-158905P16	2.0	Yes	MOUNTING SCREW FOR DRAWER CONNECTOR. SAE 6-32 THD, .966~-L, STAINLESS STL MOUNTING SCREW FOR FLOATING DRAWER C ONNECTORS. SAE 6-32 THD. (DWG REV 7/92, JFK - OK).
9	46-297328P2	2.0	Yes	1/4 TURN STUD AND TWO GROMMETS STUD LENGTHS; A=.500~ B=.330~
10	46-208758P1	4.0	Yes	3.62~ X .094~ CABLE TIE 1/16~ TO 5/8~ WIRE BUNDLE RANGE 18LB TENSILE STRENGTH, NYLON MATR'L.
11	2372346	1.0	Yes	CHANNEL TAPESWITCH
14	2199652	1.0	Yes	P9235LZ LABEL
2	2372351	1.0	No	COVER, TOP RIGHT
12	2363367	6.0	No	VINYL CAP 1/8 INCH
13	2363369	6.0	No	T SHAPE HEAD STEEL CLIP

Table 2-7 Cover Table Top Left (2372372 BOM, rev 1)

1.6 Metalless Cradle (2115992-2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
4	2268003	2.000	Yes	10-32 X .437 LONG, FLAT H
5	2222889-2	1.000	Yes	KEEPER FOR SELF-CATCHING
6	2112870-2	1.000	No	METALLESS CRADLE WITH VEL
7	46-303268P1	1.000	No	PRODUCT LOCATOR CARD FORM
3A	2285525-2	1.000	No	RATING PLATE
3B	46-302200P7	ALT	No	SERIAL #, MODEL#

Table 2-8 Metalless Cradle (2115992-2 BOM, rev 0)

This Page Intentionally Blank.

1.7 Table Subassembly (2269166)

FRU = NO




Figure 2-7 Table Subassembly (2269166, sheet 1, rev 0)




Figure 2-8 Table Subassembly (2269166, sheet 2, rev 0)




Figure 2-9 Table Subassembly (2269166, sheet 3, rev 0)




Figure 2-10 Table Subassembly (2269166, sheet 4, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2269168	1.000	Yes	CARRIAGE ASSEMBLY FOR CAT {For further detail, see Section 1.7.1 , beginning on page 76 .}
4	2269177	1.000	Yes	LATCH ASSEMBLY FOR CATHAY {For further detail, see Section 1.7.2 , beginning on page 78 .}
7	46-296370P1	2.000	Yes	ROLLER
11	2269169	1.000	Yes	CRADLE DRIVE ASM FOR CATH {For further detail, see Section 1.7.3 , beginning on page 80 .}
12	2269175	1.000	Yes	LONGITUDINAL ENCODER ASM {For further detail, see Section 1.7.4 , beginning on page 86 .}
14	46-136334P23	2.000	Yes	SPDT, 15A @125/250VAC BAS
19	2269173	1.000	Yes	ELEVATION ENCODER ASM FOR {For further detail, see Section 1.7.7 , beginning on page 96 .}
20	46-170047P7	1.000	Yes	125 PITCHES 19.635 LENG
32	2146815-2	17.000	Yes	10-32 X 3/8 HEX SOC BUTTO
34	2265877	1.000	Yes	CALIBRATION PIN (REF 46-2
48	46-276217P2	2.000	Yes	FERRITE BEADS/FILTERS
50	46-276217P3	1.000	Yes	FERRITE BEADS/FILTERS
5	2261680	1.000	No	PLATE, TABLE TOP, (REF 46
6	2265867	1.000	No	BLOCK GUIDE ROLLER (REF 4
8	2261723	2.000	No	BLOCK, BUMPER, (REF 46-27
9	46-221597P2	4.000	No	MINOR VBM-3003 NAT RUBBER
10	2265854	1.000	No	PLATE, SWITCH (REF 46-278
13	2265867-2	1.000	No	WIDE GROOVE ROLLER (REF 4
15	2261679	1.000	No	BRACE, REAR (REF 46-29681
16	2261698	1.000	No	BRACKET, SWITCH, (REF 46-
17	46-297059G1	1.000	No	SERVO ASM {For further detail, see Section 1.7.5 , beginning on page 90 .}
18	2269172	1.000	No	POWER SUPPLY ASM FOR CATH {For further detail, see Section 1.7.6 , beginning on page 92 .}

Table 2-9 Table Subassembly (2269166 BOM, rev 2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
22	46-296626P1	1.000	No	LARGE GEAR
25	46-170623P2	AR	No	BLACK, HEAVY-DUTY LUBRICA
26	46-170498P54	4.000	No	SCREW
27	46-170498P13	8.000	No	SCREW
29	12-276	5.000	No	ID 0.2040 OD 0.5630 THK .
30	46-220181P5	10.000	No	WASHER #10 (.190) SCREW S
31	46-328417P15	2.000	No	SCREW HEXAGON SOCKET 6 MM
37	2261701	1.000	No	BAR, CALIBRATION, (REF 46
38	46-170684P2	AR	No	LOCTITE 271, 10CC BOTTLE
39	46-170686P2	AR	No	LOCTITE 242, 10CC
40	2146815	5.000	No	10-32 X 1/2 HEX SOC BUTTO
41	2261724	1.000	No	CAL PLATE FOR H2 LCC CATH
43	46-328417P12	9.000	No	SCREW HEXAGON SOCKET 6 MM
44	12-275	4.000	No	WASHER, ID .2656, OD.75,
45	46-328432P3	2.000	No	WASHER LOCK - SPRING 6.1
46	2269174	1.000	No	FRAME ASSEMBLY FOR H2 LCC {For further detail, see Section 1.7.8 , beginning on page 100 .}
47	46-170015P28	1.000	No	SCREW

Table 2-9 Table Subassembly (2269166 BOM, rev 2)

1.7.1 Carriage Assembly (2269168)

FRU = YES




Figure 2-11 Carriage Assembly (2269168, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	46-218193P1	2.000	Yes	WHEEL DIA 01.21 WIDTH 0.4
5	46-208580P6	3.000	Yes	CAM FOLLOWER
6	46-208580P19	1.000	Yes	CAM FOLLOWER
10	46-170498P2	4.000	Yes	SCREW
12	46-170024P4	1.000	Yes	SCREW-INCH
2	2262161-2	1.000	No	CRADLE CARRIAGE MACH.(REF
4	2264566	2.000	No	STUD, ECCENTRIC (REF 46-2
7	2261695	1.000	No	LATCH ACTUATOR, (REF 2133
8	2264561	1.000	No	ACTU., LONG.LIMIT SWITCH(
9	46-170012P77	3.000	No	LOCK NUT, 1/4-28, .437 X
11	2264577	1.000	No	SPACER, 6.35MM THK (REF 4
14	46-170686P2	AR	No	LOCTITE 242, 10CC

Table 2-10 Carriage Assembly (2269168 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
8	46-296972G1	1.000	Yes	SOLENOID ASM}
2	2261672	1.000	No	BRKT, TABLE HOME LATCH,(R
3	2261697	1.000	No	BLK,PVT, HOME LATCH(REF 4
4	2261721	1.000	No	BAR, (REF 46-278154P1)
5	2261725	1.000	No	BRKT, LATCH ASM, (REF 46-
6	2261715	1.000	No	ARM, (REF 46-278495P1)
7	2261681	1.000	No	BRACKET, SPRING, (REF 46-
9	46-170490P54	1.000	No	SPRING-INCH EXTENSION 1 I
10	46-170497P36	3.000	No	SCREW
11	46-208732P14	1.000	No	SCREW
12	46-170012P63	1.000	No	0010-32 HEX .375 NUT HT
13	46-170020P80	1.000	No	NOMINAL DIA 0.250 X 0.750
14	46-170020P66	1.000	No	NOMINAL DIA 0.125 X 0.250
15	46-170020P85	1.000	No	NOMINAL DIA 0.125 X 0.75
16	403B335P17	1.000	No	02 POSITION TERMINAL STRI
17	403B320P94	1.000	No	2 POSITION MARKER STRIP,
18	46-170497P38	2.000	No	SCREW
19	46-170497P7	4.000	No	SCREW
20	46-170686P2	AR	No	LOCTITE 242, 10CC

Table 2-11 Latch Assembly (2269177 BOM, rev 0)

FRU = YES

Figure 2-13 Cradle Drive Assembly (2269169 sheet 1, rev 0)




Figure 2-14 Cradle Drive Assembly (2269169 sheet 2, rev 0)

MARK	NAME	QTY	FRU	DESCRIPTION
5	46-136282P30	2.0	Yes	BEARING BALL, SINGLE RO 20 MM 42 MM ID 0.7874~ OD 1.6535~ WIDTH 0.4724~ 2 SHIELDS, ABEC 1, STL SAE 52100 0.0002~-0.0008~ RAD. PLAY 09 BALLS, 1/4~ DIA, REF.
6	46-212989P2	2.0	Yes	WASHER
7	46-208562P8	2.0	Yes	SCREW 3/8X16X1~ FLAT HEAD STL F70B5A2 ZINC PL HT TR FOR HYD EMB
10	46-278575P1	1.0	Yes	SERVO MOTOR
13	46-170497P38	6.0	Yes	SCREW SCREW OLD HD STL F70B5A2 ZINC PL HT TR FOR HYD EMB N647P13008
14	403B320P86	1.0	Yes	NUMBERS 1 THU 6, LEFT TO RIGHT (DWG REV 3/91, TK)
16	46-170498P80	8.0	Yes	SCREW SCREW OLD F70B5A2 ZINC PL HT TR FOR HYD EMB
17	46-220181P3	4.0	Yes	WASHER HELICAL#6 SCREW SIZE (0.138) STEEL ZN PLT ID 0.141-0.148 OD 0.25 WIDTH 0.047 THK 0.031
18	1000904P116	5.0	Yes	ID 0.1495 OD 0.6250 THK .0598 STL B8A4B2 F70B5A
19	46-221532P18	1.0	Yes	PITCH 1/5 WIDTH 00.250 LG 009.600 NEOPRENE EXTRA LIGHT DUTY
20	46-170576P7	1.0	Yes	PITCH 1/5 TEETH 0L6 FACE 0.44 ALUM OD 00.999 ID 0.3125 LG 00.560 STOCK DRIVE 6A3-16DF02510
21	46-170576P6	1.0	Yes	PITCH 1/5 TEETH 013 FACE 0.44 ALUM OD 00.808 ID 0.250 LG 00.687 STOCK DRIVE 6A3-13DF02508
22	46-296972G1	1.0	Yes	SOLENOID ASM
26	46-170686P2	0.0	Yes	LOCTITE 242, 10CC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER
28	46-208562P7	4.0	Yes	SCREW SCREW OLD STL F70B5A2 ZINC PL HT TR FOR HYD EMB
30	2265826	1.0	Yes	SHAFT (REF 46-296048P1)
32	46-170498P79	1.0	Yes	SCREW SCREW OLD F70B5A2 ZINC PL HT TR FOR HYD EMB
34A	46-222894P1	1.0	Yes	USM CORP HDC3C-86-2BL
36	46-170498P115	6.0	Yes	SCREW SCREW OLD F65D6A BLK OX
38	46-220179P8	1.0	Yes	EXT Axial Mount 0.375 Shaft DIA 0.338 Free DIA .025 THK STL OIL DIP N900P37 BASIC
39	46-136282P33	2.0	Yes	BEARING ID 0.3750~ OD 0.8750~ WIDTH 0.2812~ 2 SHIELDS, ABEC 1, STL SAE 52100 0.0001~-0.0005~ RAD. PLAY 07 BALLS, 5/32~ DIA, REF.
40	46-136294P9	2.0	Yes	BEARING ID 2.3125 OD 2.8750 WIDTH 0.281 2 SHIELDS ABEC-1 STL SAE 52100 .0008-.0018 RAD. PLAY

Table 2-12 Cradle Drive Assembly (2269169 BOM, rev 3)

MARK	NAME	QTY	FRU	DESCRIPTION
43	46-136303P10	1.0	Yes	BEARING ID 2.00 OD 2.500 WIDTH .250 KAYDON TYPE C KA020CP0
44	46-136282P62	2.0	Yes	BEARING BALL, SINGLE RO 17 MM 35 MM ID 0.6693~ OD 1.3780~ WIDTH 0.3937~ 2 SEALS, ABEC 1, STL SAE 52100 0.0001~-0.0007~ RAD. PLAY 10 BALLS, 3/16~ DIA, REF.
45	46-220179P34	1.0	Yes	EXT AXIAL MOUNT 0.669 SHAFT DIA 0.626 FREE DIA 0.037 THK STEEL-OIL DIP
46	46-170498P84	6.0	Yes	SCREW, CAP SOCKET HEAD 8-32 X .50 F70B5A2 ZINC PL HT TR FOR HYD_EMB
48	46-296368G1	1.0	Yes	No description provided
50	46-208718P2	2.0	Yes	5/16 ID LOOP TYPE TA MFG CO. TA713D5
51	509A859P69	1.0	Yes	SPRING-INCH COMPRESSION 1.84 INCHES FREE LGTH 1.84 5 ACT COILS RATE-LB/IN 1.5 RH OR LH WOUND
55	46-297820P1	1.0	Yes	SHOULDER SCREW
59	46-311873P8	1.0	Yes	PIN, M3 X 18
61	46-221560P4	1.0	Yes	SCREW SCREW OLD PT SST 18-8
63	46-220277P6	1.0	Yes	EXT RADIAL MOUNT 0.219 SHAFT DIA 0.185 FREE DIA-.025 THK STL BLK PHOS N910P21 E-RING WALDES TRUARC #X5133-21
64	46-170012P63	1.0	Yes	0010-32 HEX .375 NUT HT .239 STL ZINC PL ESNA 21NM-02
70	46-170497P7	5.0	Yes	SCREW SCREW OLD HD STL F70B5A2 ZINC PLT HT TR FOR HYD EMB
71	46-297117G1	1.0	Yes	CRDR PWR CBL
72	46-297118G1	1.0	Yes	CRDR SIG CBL
74	46-208758P3	3.0	Yes	7.31~ X .184~ SELF-LOCKING CABLE TIE 1/16~ TO 1-3/4~ WIRE BUNDLE RANGE. 50LB TENSILE STRENGTH, NYLON MATR~L.
75	46-208758P9	1.0	Yes	11.08~X .140~ SELF-LOCKING CABLE TIE 1/16~ TO 3~ WIRE BUNDLE RANGE. 40LB TENSILE STRENGTH, NYLON MATR~L. (DWG REV 5/92, JFK - OK).
76	46-297155G1	1.0	Yes	ZEUS TABLE, 9~ LONG CABLE
77	46-220181P5	3.0	Yes	WASHER #10 (.190) SCREW SIZE, ID .194-.202 OD .334, WIDTH .062 THK .047, STEEL ZINC PLT
78	12-276	2.0	Yes	ID 0.2040 OD 0.5630 THK .0478 STL F70B5A
79	46-170623P2	0.0	Yes	BLACK, HEAVY-DUTY LUBRICATING PASTE. (1 PINT CAN) SHELF LIFE, UNDER NORMAL STORAGE CON DITIONS, HAS NO KNOWN LIMITS. (DWG REV 6/92, JFK - OK).
80	46-170498P78	4.0	Yes	SCREW SCREW OLD F70B5A2 ZINC PL HT TR FOR HYD EMB

Table 2-12 Cradle Drive Assembly (2269169 BOM, rev 3) (Continued)

MARK	NAME	QTY	FRU	DESCRIPTION
81	46-208758P10	1.0	Yes	8.00~ X .094~ SELF-LOCKING CABLE TIE 1/16~ TO 2~ WIRE BUNDLE RANGE. 18LB TENSILE STRENGTH, NYLON MATR~L. SAME AS 46-208758P7 EXCEPT RESTRICT VENDOR TO T&B. TIE MUST BE .094~ WID E. (DWG REV 5/92, JFK - OK).
2	2265831	1.0	No	END PLATE, LH (REF 46-296068P1)
3	2265832	1.0	No	END PLATE, RH (REF 46-296069P1)
9	2265818	1.0	No	IDLER SHAFT (REF 46-296038P1)
12	403B335P10	1.0	No	06 POSITION TERMINAL STRIP, 20 AMP, W/QUICK CONN TERMINALS (RIVETED). (7/16~ CENTERS, W/2 MOUNTING HOLES). TWO TABS PER POSITION. (DWG REV, TK, 4/91).
15	403B320P97	1.0	No	2 POSITION MARKER STRIP. 7 & 8 MARKED HORIZONTALY AT TOP. .437~ CENTERS, 2 MOUNTING HOLES. NUMBERS 7 & 8, LEFT TO RIGHT. (DWG R EV 6/92, JFK - OK).
23	2261718	1.0	No	BRACKET, SOLENOID, (REF 46-296107P1)
24	2265829	1.0	No	SPACER BAR (REF 46-296063P1)
25	2265830	2.0	No	PIVOT PIN (REF 46-296064P1)
27	2265833	2.0	No	SUPPORT BAR (REF 46-296070P1)
31	2265824	1.0	No	BEARING SUPPORT (REF 46-296042P1)
35	46-221451P1	0.0	No	EXXON BEACON 325 35 LB PAIL
37	2265865	1.0	No	KEY 2.38SQX15MMLG (REF 46-217252P4)
41	2265834	1.0	No	SHAFT (REF 46-296110P1)
42	2265825	1.0	No	BEARING SUPPORT (REF 46-296043P1)
49	403B335P17	1.0	No	02 POSITION TERMINAL STRIP, 20 AMP, W/QUICK CONN TERMINALS (RIVETED). (7/16~ CENTERS, W/2 MOUNTING HOLES). TWO PUSH-ON TABS PER POSITION. (DWG REV 4/91, TK)
54A	46-297397P1	1.0	No	No description provided
56	2261694	1.0	No	CLUTCH LINK BAR, (REF 46-327185P1)
60	2261717	1.0	No	BRACKET, STOP, (REF 46-296084P1)
66	2265865-2	1.0	No	KEY 4.76SQX50MMLG
73	1000904P514	2.0	No	I.D. 0.789 O.D. 1.100 THK .0785 MATL STEEL B8A3A TRUARC 5900-78 OR EQUIV.
82A	2134525	1.0	No	ZEUS TABLE CRADLE DRIVE; HIAWATHA CARBOXYLATED NITRILE
83A	2134526	1.0	No	ZEUS TABLE CRADLE DRIVE; HIAWATHA CARBOXYLATED NITRILE

Table 2-12 Cradle Drive Assembly (2269169 BOM, rev 3) (Continued)

This Page Intentionally Blank.

1.7.4 Longitudinal Encoder Asm (2269175)

FRU = YES




Figure 2-15 Longitudinal Encoder Asm (2269175, sheet 1, rev 0)




Figure 2-16 Longitudinal Encoder Asm (2269175, sheet 2, rev 0)

MARK	NAME	QTY	FRU	DESCRIPTION
2	46-296854P2	1.0	Yes	2.5 INCH, 2048 COUNT.REV INCREMENTAL W/0.027~ FLANGE ENCODER
7	46-296171P1	1.0	Yes	SPRING COIL
11	46-208562P13	5.0	Yes	SCREW SCREW OLD STL F65D6A BLK OX
12	46-208562P24	6.0	Yes	SCREW SCREW OLD STL F70B5A2 ZINC PL HT TR FOR HYD EMB N646P13006
13	46-170623P2	0.0	Yes	BLACK, HEAVY-DUTY LUBRICATING PASTE. (1 PINT CAN) SHELF LIFE, UNDER NORMAL STORAGE CON DITIONS, HAS NO KNOWN LIMITS. (DWG REV 6/92, JFK - OK).
17	46-297036G2	1.0	Yes	POT ASSEMBLY
21	46-286100P5	1.0	Yes	75 PITCHES, 7.365 INCH LENGTH 4 LBF/PIN CAPACITY
23	46-296781P1	1.0	Yes	CABLE
24	46-170497P42	1.0	Yes	SCREW SCREW OLD HD STL F65D6A BLACK OXIDE N647P13004
25	12-1723	1.0	Yes	ID 0.1590 OD 0.7500 THK .0478 STL F70B5A
26	46-296780P2	4.0	Yes	SYN-CLAMP
28	46-208735P38	1.0	Yes	SPLIT LENGTHWISE .094 DIA X 0.625 LG .022 WALL 420 S ST N533P910
29	46-170686P2	0.0	Yes	LOCTITE 242, 10CC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER
30	46-170528P23	1.0	Yes	SHAFT SIZE .250 MAT 2024-T4 ALUM ANODIZE
31	46-170528P70	1.0	Yes	CLAMP
32	46-208562P36	2.0	Yes	SCREW SCREW OLD STL F65D6A BLK OX
33	46-208758P1	1.0	Yes	3.62~ X .094~ CABLE TIE 1/16~ TO 5/8~ WIRE BUNDLE RANGE 18LB TENSILE STRENGTH, NYLON MATR'L.
3	2261719	1.0	No	BRKT, LG ENCDR ASM (REF 46-296437P1)
4	2265255	1.0	No	CABLE DRUM (REF 46-296168P1)
5	2265313	1.0	No	THUMB WHEEL (REF 46-296169P1)
6	2265260	1.0	No	SHAFT, SPRING REEL (REF 46-296893P1)
8	2265257	1.0	No	COVER, SPRING REEL (REF 46-296699P1)
9	2265256	3.0	No	SPACER (REF 46-296173P1)
10	2261716	1.0	No	MNTING PLATE ENCD,(REF 46-296436P1)
14	46-208894P1	2.0	No	SCREW SCREW OLD

Table 2-13 Longitudinal Encoder Asm (2269175 BOM, rev 1)

MARK	NAME	QTY	FRU	DESCRIPTION
15	2264577-2	1.0	No	SPACER, 4.24MM THK (REF 46-296946P5)
16	2261707	1.0	No	BRKT, POTENTIOMETER(REF 46-296435P1)
18	2329244-2	1.0	No	SPACER
19	46-296438P1	1.0	No	SPROCKET MOD
20	46-296627P4	1.0	No	SPROCKET
22	46-297539P2	1.0	No	NO. 6 HI-COLLAR SPRING, CARBON STEEL ZINC PLATED.
34	2335271	1.0	No	TEFLON SPACE
35	46-297539P2	1.0	No	NO. 6 HI-COLLAR SPRING, CARBON STEEL ZINC PLATED.

Table 2-13 Longitudinal Encoder Asm (2269175 BOM, rev 1) (Continued)

1.7.5 Servo Assembly (46-297059G1)

FRU = NO

USED ON:
2269166,
2271241




Figure 2-17 Servo Assembly (46-297059, sheet 1, rev 2)

MARK	NAME	QTY	FRU	DESCRIPTION
2	46-264370G1	1.0	Yes	CRADLE AMP FIRST MADE FOR 9800 ET
3	46-288170G1	1.0	Yes	ELEV/TILT AMP
6	46-208560P27	6.0	Yes	SCREW SCREW OLD ZINC PL
4	46-296983P1	1.0	No	CIRCUIT BOARD BRACKET
5	46-297060P1	2.0	No	WHEEL ENCODE
7	46-230523P6	2.0	No	MOLDED NYLON FLAT WASHER. .141~-I.D., .375~-O.D., .0625~ THK. (DWG REV 8/92, JFK - OK).

Table 2-14 Servo Assembly (46-297059G1 BOM, rev 9)

1.7.6 Power Supply Assembly (2269172)

FRU = NO

USED ON:
2269166




Figure 2-18 Power Supply Assembly (2269172, sheet 1, rev 0)




Figure 2-19 Power Supply Assembly (2269172, sheet 2, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
4	46-229455P1	1.000	Yes	FAN TUBEAXIAL
5	46-297093P1	3.000	Yes	SPST, AC/DC TOGGLE S
6	46-296316P2	1.000	Yes	EMC POWER SUPPLY SINGLE O
7	46-296317P3	1.000	Yes	OBC-STC & TABLE EMC POWER
10	46-221914P3	1.000	Yes	INDICATOR LIGHT W/ BUILT-
20	46-278446P1	1.000	Yes	15A, 125VAC, 2 POLE, 3 WI
28	46-170053P13	1.000	Yes	3PDT, 12V DC RELAY (EN
2	2261684	1.000	No	BRKT, PWR SUPPLY, (REF 46
3	46-170497P6	8.000	No	SCREW
8	403B335P17	1.000	No	02 POSITION TERMINAL STRI
9	46-208734P26	5.000	No	SCREW
11	2264567	1.000	No	BAR, WIRE GROUND (REF 46-
12	46-170012P63	3.000	No	0010-32 HEX .375 NUT HT
13	46-297088G1	1.000	No	POWER ASSEMBLY HARNESS
14	512A720P16	4.000	No	ID .219 OD .562 THK .11
15	1000904P329	2.000	No	ID 0.1960 OD 0.6250 THK .
16	46-220489P2	1.000	No	FAN ACCESSORIES FINGER GU
17	46-297734G1	1.000	No	RESISTOR ASM
18	403B320P94	1.000	No	2 POSITION MARKER STRIP,
19	46-170497P38	2.000	No	SCREW
21	46-297325P1	1.000	No	FISHPAPER
22	2261687	1.000	No	BRKT, PWR SUPPLY, (REF 46
23	46-208758P1	1.000	No	3.62" X .094" CABLE TIE
24	2264568	1.000	No	CVR, SERVICE OUTLET(REF 4
25	46-170012P61	3.000	No	0006-32 HEX .313 NUT HT
26	46-208759P6	1.000	No	7.81" X .184" SELF-LOCKIN
27	46-170171P22	AR	No	TAPE OLD

Table 2-15 Power Supply Assembly (2269172 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
29	46-170497P58	2.000	No	SCREW
30	1000904P33	6.000	No	ID 0.173 OD 0.375 THK .
32	46-271110P1	1.000	No	LABEL, PROTECTIVE EARTH (
33	46-208560P53	2.000	No	SCREW
35	46-170497P64	1.000	No	SCREW
36	46-220181P6	1.000	No	1/4" SCREW SIZE (0.250),
37	46-220181P5	10.000	No	WASHER #10 (.190) SCREW S
38	46-170497P7	10.000	No	SCREW
39	46-271110P6	1.000	No	"READ ALL APPLICABLE DOCU

Table 2-15 Power Supply Assembly (2269172 BOM, rev 0) (Continued)

1.7.7 Elevation Encoder Assembly (2269173)

FRU = YES




Figure 2-20 Elevation Encoder Assembly (2269173, sheet 1, rev 1)




Figure 2-21 Elevation Encoder Assembly (2269173, sheet 2, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	46-296854P2	1.000	Yes	2.5 INCH, 2048 COUNT.REV
4	46-296780P2	3.000	Yes	SYN-CLAMP
10	46-136334P57	2.000	Yes	SPDT, MINIATURE BASIC
3	2265263	1.000	No	BLOCK, (REF 46-296772P1)
5	2265262	1.000	No	SHAFT, ELEV ENCODER(REF 4
6	46-296259P4	1.000	No	ADAPTER
7	46-170048P8	1.000	No	MIN-E-PITCH, 16 T, PD 0.8
8	2265261	1.000	No	CAM, ELEV SWITCH (REF 46-
9	2261699	1.000	No	BRKT, SWT ELEV ASM,(REF 4
11	46-136282P75	2.000	No	BEARING
12	46-220277P7	1.000	No	EXT RADIAL MOUNT 0.250 SH
13	2265264	1.000	No	NUT PLATE (REF 46-296945P
14	46-170497P7	2.000	No	SCREW
15	46-170498P125	3.000	No	SCREW
16	46-170528P23	1.000	No	SHAFT SIZE .250 MAT 2024-
17	2264577-3	1.000	No	SPACER, 10.16MM THK(REF 4
18	46-170020P63	2.000	No	NOMINAL DIA 0.250 X 0.625
19	46-220181P2	1.000	No	#4 SPRING LOCK WASHER, ST
20	1000904P438	2.000	No	ID 0.12, OD 0.312, THK .0
21	2265313	1.000	No	THUMB WHEEL (REF 46-29616
22	46-208894P1	2.000	No	SCREW
23	46-297539P2	1.000	No	NO. 6 HI-COLLAR SPRING,
24	1000904P519	2.000	No	WASHER .117 ID .190 OD X
25	46-220205P1	1.000	No	ID 0.3125 OD 0.75 THK

Table 2-16 Elevation Encoder Assembly (2269173 BOM, rev 1)

This Page Intentionally Blank.




Figure 2-23 Frame Assembly (2269174, sheet 2, rev 0)




Figure 2-24 Frame Assembly (2269174, sheet 3, rev 0)




Figure 2-25 Frame Assembly (2269174, sheet 4, rev 0)




Figure 2-26 Frame Assembly (2269174, sheet 5, rev 0)

MARK	NAME	QTY	FRU	DESCRIPTION
-	46-296561P11	-	Yes	MAGNET ROD
4	2271256	1.0	Yes	SIGNAL HARNESS FOR H2 TABLE
5	2265855	1.0	Yes	PIN, ROD END (REF 46-296770P1)
6	2294322-2	1.0	Yes	Actiator, Linear, Zeus thru H16 Table, Pinned
8	2100671	2.0	Yes	CT Table gas spring
18	2259859	1.0	Yes	CARRIAGE_RAIL, (REF 46-278261P1)
19	2259860	1.0	Yes	CARRIAGE RAIL, V, (REF 46-278262P1)
23	46-297350P1	2.0	Yes	SPACER
24	46-220179P15	4.0	Yes	EXT Axial Mount 0.625 Shaft DIA 0.579 Free DIA .035 THK STL OIL DIP N900P62 BASIC
25	46-220277P15	2.0	Yes	EXT Radial Mount 0.312 Shaft DIA 0.243 Free DIA-.025 THK STL BLK Phos N910P31 E-Ring
26	46-220277P11	1.0	Yes	EXT RADIAL MOUNT 0.625 SHAFT DIA 0.480 FREE DIA-.042 THK STL BLK PHOS N910P62
28	46-328417P26	2.0	Yes	SCREW HEXAGON SOCKET 10 MM 40 MM ZINC PLATED, STEEL
29	46-328417P15	4.0	Yes	SCREW HEXAGON SOCKET 6 MM 25 MM ZINC PLATED, STEEL
30	46-170498P14	32.0	Yes	SCREW SCREW OLD F70B5A2 ZINC PL HT TR FOR HYD EMB
32	46-170012P81	1.0	Yes	005/8-18 HEX .940 NUT HT .750 STL ZINC PLATED ESNA 51NE-108
33	46-220181P6	35.0	Yes	1/4~ SCREW SIZE (0.250), STL NI PLT ID 0.255-0.263 OD 0.489 WIDTH 0.109 THK 0.062
35	46-170684P1	0.0	Yes	LOCTITE 271, 0.5CC VIAL 12 Month Shelf Life from Date of Shipment from the Manufacturer
37	2265876	2.0	Yes	PIN, GAS SPRING (REF 46-297436P1)
40	12-128	4.0	Yes	WASHER,FLAT ID 0.625 OD 1.156 THK. .0598 MATL: STEEL FINISH: PLAIN
41	46-170623P2	0.0	Yes	BLACK, HEAVY-DUTY LUBRICATING PASTE. (1 PINT CAN) SHELF LIFE, UNDER NORMAL STORAGE CON DITIONS, HAS NO KNOWN LIMITS. (DWG REV 6/92, JFK - OK).
42	2266015	1.0	Yes	TABLE POWER HARNESS
44	46-297329P2	6.0	Yes	1/4 TURN RECEPTACLE - FIXED .128~ DIA. MOUNTING HOLES
45	46-220186P20	12.0	Yes	1/8~ 120 Degree Countersunk Pop Rivet .129-.133,GRIP RANGE .126-.187 USM Corporation Part # AK4385
49	12-275	4.0	Yes	WASHER, ID .2656, OD.75, THK .050 STEEL F70B5A
50	2261663	2.0	Yes	BRKT, TBL TOP COVER, (REF 2138099)
52	46-297705P2	1.0	Yes	CONDUCTIVE TAPE

Table 2-17 Frame Assembly (2269174 BOM, rev 5)

MARK	NAME	QTY	FRU	DESCRIPTION
63	46-297687P1	1.0	Yes	8 OHM, 4W, 3-1/2~ SPEAKER. 2.9 OZ CERAMIC MAGNET, 90 DB. (0.187~ Q.C. TERMINALS) (DWG REV 3/92, JFK - OK).
64	46-297686P2	1.0	Yes	SPEAKER INSULATION FOR H2 HCC CATHAY TABLE
66	46-297155G1	4.0	Yes	ZEUS TABLE, 9~ LONG CABLE
75	46-208759P4	14.0	Yes	13.93~X .301~ SELF-LOCKING CABLE TIE 3/16~ TO 3-1/2~ WIRE BUNDLE RANGE. 120LBS TENSILE STRENGTH, NAT. NYLON. (10-32 SCREW MOUNTING HOLE SIZE) MOUNTING HOLE FOR 10-32 SCREW. (DWG REV 4/92, JFK - OK).
77	46-297648P1	1.0	Yes	CABLE CLAMP
79	46-297698P1	2.0	Yes	11~ LOW PROFILE BEAD, RIBBON SWITCH. 28 VAC OR VDC @ 2A, .12~ BEND RADIUS 8OZ ACTUATION FORCE FOR .177D ROD. (AXIAL LEADS W/ 2-PIN CONNECTORS) SEE TEXT DATA AND DRAWING FOR COMPLETE DETAILS. (DWG REV 3/92, JFK - OK) LOWER RIGHT & LOWER LEFT TAPESWITCH
81	1000904P395	4.0	Yes	ID 0.177 OD 0.312 THK .062 B4B33A F70B5A
82	46-208759P1	34.0	Yes	7.81~ X .184~ SELF-LOCKING CABLE TIE 1/16~ X 1-3/4~ WIRE BUNDLE RANGE. 50 LBS TENSILE STRENGTH, NAT. NYLON. (#10 SCREW MOUNTING HOLE SIZE) MAT~L: NYLON 6/6 PER MIL-N-20693 AND ASTM D4066 STOCK #: 08415 - NATURAL, TYPE I
83	46-208560P46	1.0	Yes	SCREW SCREW OLD ZINC PL
84	46-170497P38	4.0	Yes	SCREW SCREW OLD HD STL F70B5A2 ZINC PL HT TR FOR HYD EMB N647P13008
85	46-170012P63	24.0	Yes	0010-32 HEX .375 NUT HT .239 STL ZINC PL ESNA 21NM-02
86	46-170012P61	3.0	Yes	0006-32 HEX .313 NUT HT. .178 STL ZINC PL ESNA 21NM-62 OR AERO 21NM-0632
88	2146815-2	22.0	Yes	10-32 X 3/8 HEX SOC BUTTON HD SST MA
90	12-276	4.0	Yes	ID 0.2040 OD 0.5630 THK .0478 STL F70B5A
91	46-170497P36	5.0	Yes	SCREW SCREW OLD HD STL F70B5A2 ZINC PL HT TR FOR HYD EMB N647P13006
92	46-208758P9	2.0	Yes	11.08~X .140~ SELF-LOCKING CABLE TIE 1/16~ TO 3~ WIRE BUNDLE RANGE. 40LB TENSILE STRENGTH, NYLON MATR~L. (DWG REV 5/92, JFK - OK).
93	46-208758P5	2.0	Yes	SELF LOCKING CABLE TIE
96	46-208758P1	4.0	Yes	3.62~ X .094~ Cable Tie 1/16~ TO 5/8~ Wire Bundle Range 18LB Tensile Strength, Nylon Matr'l.
97	46-297805G1	2.0	Yes	JUMPER CONN.
98	46-297805G2	2.0	Yes	JUMPER CONN.

Table 2-17 Frame Assembly (2269174 BOM, rev 5) (Continued)

MARK	NAME	QTY	FRU	DESCRIPTION
99	46-208562P7	12.0	Yes	SCREW SCREW OLD STL F70B5A2 ZINC PL HT TR FOR HYD EMB
100	46-170012P62	4.0	Yes	0008-32 HEX .343 NUT HT .229 STL ZINC PL ESNA 21NM-82
102	46-170012P64	1.0	Yes	01/4-20 HEX .437 NUT HT .312 STL ZINC PL ESNA 41NE-040
104	2106866-2	1.0	Yes	.75 IN WIDE INSULATED CABLE W/LUGS 9.00 INCHES LONG
107A	2336449	1.0	Yes	H2 TABLE 170VDC PS BD ASSY
108	46-297335P2	2.0	Yes	6.90A, 250V 2-STAGE PWR LINE FILTER. 2 X 0.21MA LEAK CURRENT, 220V/50HZ. SOLDER LUGS/FAST-ON TERMINALS. (UL, CSA & VDE APPROVALS) 2-STAGE, GENERAL PURPOSE FILTER W/ H IGH ATTENUATION. (DWG REV 3/92, JFK - OK).
114	46-328417P12	38.0	Yes	SCREW HEXAGON SOCKET 6 MM 12 MM ZINC PLATED, STEEL
115	46-170015P3	4.0	Yes	SCREW SEMS 8-32 X .312, EX TOOTH LOCK WASHER, ZINC PLT F70B5A
119	46-170498P32	6.0	Yes	SCREW SCREW OLD F70B5A2 ZINC PL HT TR FOR HYD EMB
120	46-208562P20	6.0	Yes	SCREW SCREW OLD STL F70B5A2 ZINC PL HT TR FOR HYD EMB N646P16010
121	46-328416P23	4.0	Yes	SCREW HEXAGON HEAD 10 MM 40 MM ZINC PLATED, STEEL
122	46-328431P3	4.0	Yes	WASHER PLAIN - LARGE 10.5 MM 30 MM ZINC PLATED, STEEL
123	46-328432P4	4.0	Yes	WASHER LOCK - SPRING 10.2 MM 18.1 MM ZINC PLATED, SPRING STEEL
124	46-170686P2	0.0	Yes	LOCTITE 242, 10CC 12 Month Shelf Life from Date of Shipment from the Manufacturer
125	46-170015P27	5.0	Yes	SCREW SCREW OLD F70B5A ZINC PL
126	12-297	5.0	Yes	ID 0.140 OD 0.312 THK .032 STL F70B5A ZINC PLATE
127	46-220181P5	6.0	Yes	Washer #10 (.190) Screw Size, ID .194-.202 OD .334, WIDTH .062 THK .047, STEEL ZINC PLT
3	2269179	1.0	No	ETC/CPU ASSEMBLY FOR CATHAY-TABLE
7A	2133760-2	1.0	No	ZEUS TABLE ACTUATOR BASE MOUNT MACHINING
9	2255160-2	1.0	No	FRONT LEG FOR H2 LCC CATHAY TABLE, MACHINING
12	46-136342P135	4.0	No	BUSHING BEARING FLANGE O.D. 1.000 THK. .125
13	2265856	1.0	No	SPACER, ROD END (REF 46-296771P1)
14	2265853	1.0	No	PIN, ACTUATOR (REF 46-278310P1)
15	2255158-2	1.0	No	REAR LEG FOR H2 LCC CATHAY TABLE, MACHINING
16	2259855	1.0	No	CHANNEL,LEFT SUPPORT,(REF 2138098)
17	2259854	1.0	No	CHANNEL,RIGHT SUPPORT,(REF 2138097)

Table 2-17 Frame Assembly (2269174 BOM, rev 5) (Continued)

MARK	NAME	QTY	FRU	DESCRIPTION
20	2261670	2.0	No	PIN, TABLE TOP, FRONT (REF 2140091)
21	2261671	2.0	No	PIN, TABLE TOP, REAR, (REF 2140092)
22	2261678	1.0	No	BRACE, FRONT, (REF 46-296811P1)
34	2261685	1.0	No	BRKT, R-CHANL CON, (REF 46-297055P1)
36	2261722	1.0	No	BLOCK, MOUNTING, (REF 46-297435P1)
38	2261686	1.0	No	BRKT, L-CHANL CON, (REF 46-297056P1)
39	2261703	1.0	No	BRKT, FRT RT CHAN,(REF 46-297263P1)
43	2265871	6.0	No	BRACKET TBL BASE CVR (46-297205P1)
46	2329243	1.0	No	SHAFT
47	2329245	1.0	No	ROLLER MOUNT
48	46-170020P2	2.0	No	NOMINAL DIA 0.0938X 0.500 LG PREC. HARDENED & GROUND RC 36/40 SST #416
51	2261692	1.0	No	SPRING CLAMP, (REF 46-297677P1)
53	46-297705P14	2.0	No	METALLIC CONDUCTIVE TAPE
54	46-297705P4	1.0	No	MET-CON TAPE
55	46-170497P66	1.0	No	SCREW SCREW OLD HEAD STEEL FINISH F70B5A
61	2261689	1.0	No	CABLE TROUGH, (REF 46-297343P2)
62	2261674	1.0	No	SPEAKER GRILL, (REF 46-297684P1)
65	2265879	1.0	No	SPEAKER COVER (REF 46-297685P2)
69	2265875	3.0	No	CLAMP, GRND STRAP (REF 46-297427P1)
80	2146815	25.0	No	10-32 X 1/2 HEX SOC BUTTON HD SST MA
87	2275854	1.0	No	GROUND - TABLE TO RACEWAY GND BAR
89	46-297731P1	1.0	No	MOTOR LABEL
94	2261712	1.0	No	BRKT, LH TBL TOP CVR, (REF 2138105)
95	2261712-2	1.0	No	BRKT, RH TBL TOP CVR,(REF 2138105-2)
101	2265880	4.0	No	SHOULDER SCREW (REF 46-297858P1)
105	2265847	1.0	No	BRACKET, T/E PWR SUPPLY (REF 2114431)
106	2265848	1.0	No	COVER, T/E PWR SUPPLY (REF 2114433)
112	46-297705P10	1.0	No	METALLIC CONDUCTIVE 3M-SCOTCH BRAND

Table 2-17 Frame Assembly (2269174 BOM, rev 5) (Continued)

MARK	NAME	QTY	FRU	DESCRIPTION
113	2259856-2	1.0	No	TABLE BASE,MACHINING,(REF 2132187-2)
117	2261676	2.0	No	PIN, TBL LOWER, FRONT,(REF 2136121)
118	2261677	2.0	No	PIN, TBL LOWER, REAR, (REF 2138104)

Table 2-17 Frame Assembly (2269174 BOM, rev 5) (Continued)

Section 2.0 — Table Assembly (2271242)





Figure 2-27 Table Assembly (2271242, sheet 1, rev 1)





MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
-	46-297894G1	-	Yes	FIELD KIT FOR ACTUATOR COVER ROLLER
7	2271244	1.000	Yes	RIGHT BASE COVER ASM FOR H2 {For further detail, see Section 1.2, beginning on page 62.}
8	2271247	1.000	Yes	LEFT BASE COVER ASM FOR H2 {For further detail, see Section 1.3, beginning on page 64.}
16	46-296592P1	1.000	Yes	ACTUATOR COVER
17	2236720	1.000	Yes	COVER ASSEMBLY TABLE TOP - LEFT {For further detail, see Section 2.2, beginning on page 120.}
18	2236720-2	1.000	Yes	COVER ASSEMBLY TABLE TOP - RIGHT {For further detail, see Section 2.3, beginning on page 122.}
19	46-296594P1	1.000	Yes	DRIVE COVER
21	46-296679P2	1.000	Yes	ZEUS TABLE, REAR END COVER WITHOUT
23	46-264368G3	1.000	Yes	ENHANCED TABLE CONTROL HELIOS 2
28	46-297061P1	1.000	Yes	FISH PAPER
29	2136285	2.000	Yes	ZEUS TABLE LOWER REAR SID PNL BRACKT
30	2136284	2.000	Yes	ZEUS TABLE LOWER FRNT SID PNL BRACKT
31	2138103	2.000	Yes	ZEUS TABLE UPPER SIDE PANEL BRACKET
32	2138101	2.000	Yes	UPPER SIDE PANEL BRACKET
33	46-297171P1	4.000	Yes	NUT PLATE
35	46-297324P1	1.000	Yes	COVER SUPPORT BRACKET
36	46-297329P2	2.000	Yes	1/4 TURN RECEPTACLE - FIXED
38	46-297246P1	1.000	Yes	CHANNEL
47	46-208562P20	14.000	Yes	SCREW
52	46-297552P1	8.000	Yes	SPNL PIV TUB
64	46-220186P20	4.000	Yes	1/8" 120 DEGREE COUNTERSUNK POP
72	46-297678P1	1.000	Yes	HINGES
73	46-297676P1	1.000	Yes	PIN
74	46-170490P50	1.000	Yes	SPRING-INCH EXTENSION 4 INCHES .65 I

Table 2-18 Table Assembly (2271242 BOM, rev 5)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
81	12-276	4.000	Yes	ID 0.2040 OD 0.5630 THK .
108	2115992-2	1.000	Yes	860-941208 METALLESS CRADLE {For further detail, see Section 1.6 , beginning on page 68 .}
109	2117299	1.000	Yes	METALLESS CRADLE COMPATIBLE FRONT
113	2117930	1.000	Yes	EMC PROJECT/TABLE/CVR BRD BRKT
115	46-265874P1	1.000	Yes	50 OHM IMPEDANCE BNC FEED-THRU TERM
123	46-233422P1	1.000	Yes	BNC JACK-PLUG, RIGHT ANGLE ADAPTER
128	2271251	1.000	Yes	RIBBON JUMPER CABLE, ETC TO ETC I/F
129	2271253	1.000	Yes	RIBBON CABLE JUMPER, ETC TO ETC I/F
130	2271254	1.000	Yes	RIBBON CABLE JUMPER, ETC TO ETC I/F
131	2271255	1.000	Yes	RIBBON CABLE JUMPER, ETC TO ETC I/F
135	2287357	1.000	Yes	HANDLE, CRADLE
136	46-328417P20	2.000	Yes	SCREW HEXAGON SOCKET 6 MM
105A	2224296-3	1.000	Yes	VE-GE, 8MB, 2FL, GE MED,
105B	2224296	ALT	Yes	ARTESYN II CT COMPUTER BOARD
132A	2372485	1.000	Yes	ETC INTERFACE W/ AZ60A
132B	2229871	ALT	Yes	ETC IF BOARD
4	46-296591P1	1.000	No	COVER MOUNT
5	46-296820P1	1.000	No	BRACKET
6	2169596-2	1.000	No	SIDE PANEL KIT FOR H2 HCC {For further detail, see Section 2.1 , beginning on page 118 .}
13	46-296832P1	2.000	No	REAR COVER BRACKET
14	46-297897P1	1.000	No	ENCODER COVER
15	46-297374P1	1.000	No	COVER
22	46-296742P1	2.000	No	GUARD COVER
43	2146815-2	26.000	No	10-32 X 3/8 HEX SOC BUTTON HD SST MA
45	46-208562P2	1.000	No	SCREW
46	46-208562P3	4.000	No	SCREW

Table 2-18 Table Assembly (2271242 BOM, rev 5) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
49	46-170497P38	13.000	No	SCREW
55	46-297155G1	4.000	No	ZEUS TABLE, 9" LONG CABLE
56	46-170012P61	2.000	No	0006-32 HEX .313 NUT HT
59	46-297596P1	2.000	No	HEX STANDOFF
61	46-170498P23	6.000	No	SCREW
62	46-221523P10	6.000	No	HOLE DIA 0.437 MAT'L THK
63	46-170012P63	2.000	No	0010-32 HEX .375 NUT HT
65	46-170498P41	2.000	No	SCREW
71	46-297695P1	2.000	No	COTTER PIN
75	46-297448P1	1.000	No	HANDLE ADDITION
78	46-197018P1	1.000	No	DANGER HI-VOLTAGE
80	46-220181P5	10.000	No	WASHER #10 (.190) SCREW SIZE,
85	2146815-2	18.000	No	10-32 X 3/8 HEX SOC BUTTON HD SST MA
86	46-196354P1	3.000	No	STANDOFF/SPACER
87	46-196354P4	3.000	No	STANDOFF/SPACER
88	46-170106P1	AR	No	DBL COATED FILM .007 THK
91	2132175	0.000	No	TABLE DOCUMENT COLLECTOR
92	1000904P358	1.000	No	ID 0.380 OD 1.000 THK .
93	2184717	1.000	No	SHIPPING BOLT FOR CT TABLE
94	46-327022P1	AR	No	ZEUS TABLE, CRADLE/CARRIAGE SPACER
95	46-327022P2	AR	No	ZEUS TABLE, CRADLE/CARRIAGE SPACER
96	46-327022P3	AR	No	ZEUS TABLE, CRADLE/CARRIAGE SPACER
97	46-327022P4	AR	No	ZEUS TABLE, CRADLE/CARRIAGE SPACER
98	46-327022P5	1.000	No	ZEUS TABLE, CRADLE/CARRIAGE SPACER
100	46-208561P147	5.000	No	SCREW
101	2105871	5.000	No	GANTRY AND TABLE INSTALLATION
102	2106207	5.000	No	1 1/2 -12 THDS 2.00 O.D.
103	2105872	5.000	No	GANTRY AND TABLE INSTALLATION

Table 2-18 Table Assembly (2271242 BOM, rev 5) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
116	2265684	1.000	No	LABEL
117	2117390	1.000	No	LABEL, CE MARKING COMPLIANCE
118	2132620	1.000	No	CE MARKING LABEL CLASS A GROUP 1
119	2103580-12	8.000	No	SCREW COUNTERSUNK HEA 6 MM 12 MM
120	46-328417P12	8.000	No	SCREW HEXAGON SOCKET 6 MM
121	2271241	1.000	No	SUBASSEMBLY FOR H2 HCC CATHAY TABLE {For further detail, see Section 2.4, beginning on page 124.}
122	2285525	1.000	No	RATING PLATE
124	46-208922P5	1.000	No	SCREW
125	46-208984P17	3.000	No	BRASS STANDOFF MALE/FEMALE, 6-32
126	2148658	1.000	No	ESD GROUND STRAP, BRAIDED FLAT WIRE
134	46-303268P1	1.000	No	PRODUCT LOCATOR CARD FORM

Table 2-18 Table Assembly (2271242 BOM, rev 5) (Continued)

This Page Intentionally Blank.

2.1 Side Panel Kit (2169596-2)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
9	46-296577P2	1.000	Yes	RIGHT REAR SIDE PANEL FOR H2 HCC
10	46-296578P2	1.000	Yes	LEFT REAR SIDE PANEL FOR H2 HCC
11	46-296583P2	1.000	Yes	RIGHT FRONT SIDE PANEL FOR H2 HCC
12	46-296584P2	1.000	Yes	LEFT FRONT SIDE PANEL FOR H2 HCC

Table 2-19 Side Panel Kit (2169596-2 BOM, rev 1)

2.2 Cover Assembly Table Top - Left (2236720)

FRU = YES




Figure 2-30 Cover Assembly Table Top - Left (2236720, sheet 1, rev 0)

MARK	NAME	QTY	FRU	DESCRIPTION
3	46-208747P1	4.0	Yes	ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)
4	46-297698P3	1.0	Yes	27~ HIGH PROFILE BEAD, RIBBON SWITCH
5	46-297329P1	2.0	Yes	1/4 TURN RECEPTACLE WITH TWO SPACERS FLOATING (.074~ ALL DIRECTIONS) .261~ DIA. MOUNTING HOLES
6	46-220186P20	4.0	Yes	1/8~ 120 DEGREE COUNTERSUNK POP RIVET .129-.133,GRIP RANGE .126-.187 USM CORPORATION PART # AK4385
7	46-297721P1	1.0	Yes	4-POS PLUG HOUSING (DRAWER CONNECT). PLUG HOUSING, DRAWER CONNECTOR. (DWG REV 7/92, JFK - OK).
8	46-158905P16	2.0	Yes	MOUNTING SCREW FOR DRAWER CONNECTOR. SAE 6-32 THD, .966~-L, STAINLESS STL MOUNTING SCREW FOR FLOATING DRAWER C ONNECTORS. SAE 6-32 THD. (DWG REV 7/92, JFK - OK).
9	46-297328P2	2.0	Yes	1/4 TURN STUD AND TWO GROMMETS STUD LENGTHS; A=.500~ B=.330~
10	46-208758P1	4.0	Yes	3.62~ X .094~ CABLE TIE 1/16~ TO 5/8~ WIRE BUNDLE RANGE 18LB TENSILE STRENGTH, NYLON MATR'L.
11A	2216522	1.0	Yes	CHANNEL EXTRUSION FOR TAPE SWITCH
2	2225970	1.0	No	TOP LEFT SIDE COVER

Table 2-20 Cover Assembly Table Top - Left (2236720 BOM, rev 2)

2.3 Cover Assembly Table Top - Right (2236720-2)

FRU = YES




Figure 2-31 Cover Assembly Table Top - Right (2236720-2, sheet 1, rev 0)

MARK	NAME	QTY	FRU	DESCRIPTION
3	46-208747P1	4.0	Yes	ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)
4	46-297698P3	1.0	Yes	27~ HIGH PROFILE BEAD, RIBBON SWITCH
5	46-297329P1	2.0	Yes	1/4 TURN RECEPTACLE WITH TWO SPACERS FLOATING (.074~ ALL DIRECTIONS) .261~ DIA. MOUNTING HOLES
6	46-220186P20	4.0	Yes	1/8~ 120 DEGREE COUNTERSUNK POP RIVET .129-.133,GRIP RANGE .126-.187 USM CORPORATION PART # AK4385
7	46-297721P1	1.0	Yes	4-POS PLUG HOUSING (DRAWER CONNECT). PLUG HOUSING, DRAWER CONNECTOR. (DWG REV 7/92, JFK - OK).
8	46-158905P16	2.0	Yes	MOUNTING SCREW FOR DRAWER CONNECTOR. SAE 6-32 THD, .966~-L, STAINLESS STL MOUNTING SCREW FOR FLOATING DRAWER C ONNECTORS. SAE 6-32 THD. (DWG REV 7/92, JFK - OK).
9	46-297328P2	2.0	Yes	1/4 TURN STUD AND TWO GROMMETS STUD LENGTHS; A=.500~ B=.330~
10	46-208758P1	4.0	Yes	3.62~ X .094~ CABLE TIE 1/16~ TO 5/8~ WIRE BUNDLE RANGE 18LB TENSILE STRENGTH, NYLON MATR'L.
11A	2216522	1.0	Yes	CHANNEL EXTRUSION FOR TAPE SWITCH
2	2225968	1.0	No	TOP RIGHT SIDE COVER

Table 2-21 Cover Assembly Table Top - Right (2236720-2 BOM, rev 2)

2.4 Subassembly For H2 HCC Cathay Table (2271241)

FRU = NO




Figure 2-32 Subassembly For H2 HCC Cathay Table (2271241, sheet 1, rev 1)




Figure 2-33 Subassembly For H2 HCC Cathay Table (2271241, sheet 2, rev 1)






Figure 2-35 Subassembly For H2 HCC Cathay Table (2271241, sheet 4, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2137940	1.000	Yes	ZEUS TABLE 3 POINT CARRIAGE ASM {For further detail, see Section 2.4.1 , beginning on page 130 .}
4	46-296233G2	1.000	Yes	LATCH ASM {For further detail, see Section 2.4.2 , beginning on page 132 .}
7	46-296370P1	2.000	Yes	BAR
11	46-296235G2	1.000	Yes	ZEUS TABLE, WITH HIAWATHA {For further detail, see Section 2.4.3 , beginning on page 134 .}
12	46-296234G1	1.000	Yes	ZEUS TABLE ENCODER ASM. {For further detail, see Section 2.4.4 , beginning on page 138 .}
14	46-136334P23	2.000	Yes	SPDT, 15A @125/250VAC BASIC SWITCH.
19	46-296633G3	1.000	Yes	ELEVATION ENCODER ASM FOR H2 HCC {For further detail, see Section 2.4.6 , beginning on page 146 .}
20	46-170047P7	1.000	Yes	125 PITCHES 19.635 LENGTH
25	46-170623P2	AR	Yes	BLACK, HEAVY-DUTY LUBRICATING PASTE.
32	2146815-2	17.000	Yes	10-32 X 3/8 HEX SOC BUTTON HD SST MA
34	46-297576P1	1.000	Yes	CALIBRATION PIN
38	46-170684P2	AR	Yes	LOCTITE 271, 10CC BOTTLE
5	46-278308P1	1.000	No	TOP TABLE PL
6	46-296108P1	1.000	No	GUIDE BLOCK
8	46-278362P1	2.000	No	BUMPER BLOCK
9	46-221597P2	4.000	No	MINOR VBM-3003 NAT RUBBER 70+/-5DURO
10	46-278363P1	1.000	No	SWITCH PLATE
13	46-296108P2	1.000	No	GUIDE BLOCK
15	46-296815P1	1.000	No	REAR BRACE
16	46-296925P1	1.000	No	BRACKET, SW
17	46-297059G1	1.000	No	SERVO ASM {For further detail, see Section 1.7.5 , beginning on page 90 .}
18	46-296791G2	1.000	No	POWER ASSEMBLY {For further detail, see Section 2.4.5 , beginning on page 142 .}

Table 2-22 Subassembly For H2 HCC Cathay Table (2271241 BOM, rev 3)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
22	46-296626P1	1.000	No	LARGE GEAR
26	46-170498P54	4.000	No	SCREW
27	46-170498P13	8.000	No	SCREW
29	12-276	5.000	No	ID 0.2040 OD 0.5630 THK .0478
30	46-220181P5	9.000	No	WASHER #10 (.190) SCREW SIZE,
31	46-328417P15	2.000	No	SCREW HEXAGON SOCKET 6 MM 25 MM
37	46-297934P2	1.000	No	CALIBRATION BAR H2 HCC TABLE
39	46-170686P2	AR	No	LOCTITE 242, 10CC
40	2146815	5.000	No	10-32 X 1/2 HEX SOC BUTTON HD SST MA
41	46-297921P3	1.000	No	CALIB PLATE H2 HCC TABLE
43	46-328417P12	9.000	No	SCREW HEXAGON SOCKET 6 MM 12 MM
44	12-275	4.000	No	WASHER, ID .2656, OD.75, THK .050
45	46-328432P3	2.000	No	WASHER LOCK - SPRING 6.1 MM 11.8 MM
46	2271240	1.000	No	FRAME ASSEMBLY FOR H2 HCC CATHAY TABLE { <i>For further detail, see Section 2.4.7, beginning on page 150.</i> }
47	46-170015P28	1.000	No	SCREW
48	46-276217P2	2.000	No	FERRITE BEADS/FILTERS
50	46-276217P3	1.000	No	FERRITE BEADS/FILTERS
51	46-271110P6	1.000	No	"READ ALL APPLICABLE DOCU

Table 2-22 Subassembly For H2 HCC Cathay Table (2271241 BOM, rev 3) (Continued)

2.4.13 Point Carriage Assembly (2137940)

FRU = YES




Figure 2-36 3 Point Carriage Assembly (2137940, sheet 1, rev 2)

MARK	NAME	QTY	FRU	DESCRIPTION
3	46-218193P1	2.0	Yes	WHEEL DIA 01.21 WIDTH 0.44 BORE 0.38 CAPACITY THRUST 110 LBS RADIAL 700 LBS
4	46-296911P1	2.0	Yes	No description provided
5	46-208580P6	3.0	Yes	CAM FOLLOWER THRD 003/8-24 X 0.375 LG STEEL STUD LGTH 0.875,HEX HOLE,SEALS CROWNED
6	46-208580P19	1.0	Yes	CAM FOLLOWER THRD 001/4-28 X 0.312 LG STEEL STUD LGTH 0.750, HEX HOLE, SEALS CROWNED
7	2133779	1.0	Yes	ZEUS TABLE HOME LATCH ACTUATOR
8	2133781	1.0	Yes	LONGITUD LIMIT SWITCH ACTUATOR
9	46-170012P77	3.0	Yes	LOCK NUT, 1/4-28, .437 X .312
10	46-170498P2	4.0	Yes	SCREW SCREW OLD F70B5A2 ZINC PL HT TR FOR HYD EMB
11	46-296946P1	1.0	Yes	SPACER
12	46-170024P4	1.0	Yes	SCREW-INCH SCREW OLD SHLD DIA 0.124 X SHLD LG 0.375 THD
14	46-170686P2	0.0	Yes	LOCTITE 242, 10CC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER
2	2128887-2	1.0	No	CARRAGE FOR ZEUS TABLE-MACHINING

Table 2-23 3 Point Carriage Assembly (2137940 BOM, rev 3)

FRU = YES

Figure 2-37 Latch Assembly (46-296233, sheet 1, rev 1)

MARK	NAME	QTY	FRU	DESCRIPTION
8	46-296972G1	1.0	Yes	SOLENOID ASM
9	46-170490P54	1.0	Yes	SPRING-INCH EXTENSION 1 INCHES .25 I FREE LENGTH 1.00 MUSIC WIRE RATE LBS/IN 4.7
10	46-170497P36	3.0	Yes	SCREW SCREW OLD HD STL F70B5A2 ZINC PL HT TR FOR HYD EMB N647P13006
11	46-208732P14	1.0	Yes	SCREW SCREW OLD GE NO. N175P1612 10-32 X .750 LG. STL.
12	46-170012P63	1.0	Yes	0010-32 HEX .375 NUT HT .239 STL ZINC PL ESNA 21NM-02
17	403B320P94	1.0	Yes	2 POSITION MARKER STRIP, 1 & 2 MARKED ON BOTTOM, 0.437~ CENTERS
18	46-170497P38	2.0	Yes	SCREW SCREW OLD HD STL F70B5A2 ZINC PL HT TR FOR HYD EMB N647P13008
19	46-170497P7	4.0	Yes	SCREW SCREW OLD HD STL F70B5A2 ZINC PLT HT TR FOR HYD EMB
20	46-170686P2	0.0	Yes	LOCTITE 242, 10CC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER
2	2146883	1.0	No	ZEUSTABLE HOME LATCH BRACKET
3	46-278164P2	1.0	No	PIVOT BLOCK G2 HOME LATCH, FOR .179~ THK Z-CHANNELS
4	46-278154P1	1.0	No	BAR
5	46-278153P1	1.0	No	BRACKET
6	46-278495P1	1.0	No	ARM
7	46-278494P1	1.0	No	SPRING BRACKET
13	46-170020P80	1.0	No	NOMINAL DIA 0.250 X 0.750 LG STD HARDENED & GROUND ALLOY STL
14	46-170020P66	1.0	No	NOMINAL DIA 0.125 X 0.250 LG STD HARDENED & GROUND ALLOY STL
15	46-170020P85	1.0	No	NOMINAL DIA 0.125 X 0.75 LG HARDENED & GROUND ALLOY STL
16	403B335P17	1.0	No	02 POSITION TERMINAL STRIP, 20 AMP, W/QUICK CONN TERMINALS (RIVETED). (7/16~ CENTERS, W/2 MOUNTING HOLES). TWO PUSH-ON TABS PER POSITION. (DWG REV 4/91, TK)

Table 2-24 Latch Assembly (46-296233G2 BOM, rev 1)

FRU = YES

Figure 2-38 Zeus Table, w/Hiawatha (46-296235, sheet 1, rev 9)



MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	46-296235	0.000		CRADLE DRIVE ASSEMBLY
19	46-221532P18	1.000	2	PITCH 1/5 WIDTH 00.250
22	46-296972G1	1.000	2	SOLENOID ASM
26	46-170686P2	AR	2	LOCTITE 242, 10CC
75	46-208758P9	1.000	2	11.08"X .140" SELF-LOCKIN
78	12-276	2.000	2	ID 0.2040 OD 0.5630 THK .
79	46-170623P2	AR	2	BLACK, HEAVY-DUTY LUBRICA
81	46-208758P10	1.000	2	8.00" X .094" SELF-LOCKIN
5	46-136282P30	2.000	N	BEARING BALL, SINGLE RO 2
6	46-212989P2	2.000	N	WASHER
7	46-208562P8	2.000	N	SCREW 3/8X16X1" FLAT HEAD
10	46-278575P1	1.000	N	SERVO MOTOR
12	403B335P10	1.000	N	06 POSITION TERMINAL STRI
13	46-170497P38	6.000	N	SCREW
14	403B320P86	1.000	N	NUMBERS 1 THU 6, LEFT TO
15	403B320P97	1.000	N	2 POSITION MARKER ST
16	46-170498P80	8.000	N	SCREW
17	46-220181P3	4.000	N	WASHER HELICAL#6 SCREW SI
18	1000904P116	5.000	N	ID 0.1495 OD 0.6250 THK .
20	46-170576P7	1.000	N	PITCH 1/5 TEETH 0L6 FAC
21	46-170576P6	1.000	N	PITCH 1/5 TEETH 013 FAC
28	46-208562P7	4.000	N	SCREW
32	46-170498P79	1.000	N	SCREW
34	46-222894P1	1.000	N	USM CORP HDC3C-86-2BL
35	46-221451P1	AR	N	EXXON BEACON 325 35 LB P
36	46-170498P115	6.000	N	SCREW
38	46-220179P8	1.000	N	EXT AXIAL MOUNT 0.375 SHA
39	46-136282P33	2.000	N	BEARING
40	46-136294P9	2.000	N	BEARING

Table 2-25 Zeus Table, w/Hiawatha (46-296235G2 BOM, rev 6)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
43	46-136303P10	1.000	N	BEARING
44	46-136282P62	2.000	N	BEARING BALL, SINGLE RO 1
45	46-220179P34	1.000	N	EXT AXIAL MOUNT 0.669 SHA
46	46-170498P84	6.000	N	SCREW, CAP SOCKET HEAD 8-
49	403B335P17	1.000	N	02 POSITION TERMINAL STRI
50	46-208718P2	2.000	N	5/16 ID LOOP TYPE TA MFG
51	509A859P69	1.000	N	SPRING-INCH COMPRESSION 1
55	46-297820P1	1.000	N	SHOULDER SCREW
59	46-311873P8	1.000	N	PIN, M3 X 18
61	46-221560P4	1.000	N	SCREW
63	46-220277P6	1.000	N	EXT RADIAL MOUNT 0.219 SH
64	46-170012P63	1.000	N	0010-32 HEX .375 NUT HT
70	46-170497P7	5.000	N	SCREW
71	46-297117G1	1.000	N	CRDR PWR CBL
72	46-297118G1	1.000	N	CRDR SIG CBL
73	1000904P514	2.000	N	I.D. 0.789 O.D. 1.100
74	46-208758P3	3.000	N	7.31" X .184" SELF-LOCKIN
76	46-297155G1	1.000	N	ZEUS TABLE, 9" LONG CABLE
77	46-220181P5	2.000	N	WASHER #10 (.190) SCREW S
80	46-170498P78	4.000	N	SCREW
85	2152292	1.000	N	KIT OF PARTS FOR CRADLE D
23A	46-296107P1	1.000	N	SOLENOID BRACKET
23B	2261718	ALT	N	BRACKET, SOLENOID, (REF 4
37A	46-217252P4	1.000	N	0.0938 SQ X 0.60 LG STL
37B	2265865	ALT	N	KEY 2.38SQX15MMLG (REF 46
60A	46-296084P1	1.000	N	STOP BRACKET
60B	2261717	ALT	N	BRACKET, STOP, (REF 46-29
66A	46-140913P48	1.000	N	.1875 SQ X 2.00 LONG
66B	2265865-2	ALT	N	KEY 4.76SQX50MMLG

Table 2-25 Zeus Table, w/Hiawatha (46-296235G2 BOM, rev 6) (Continued)

2.4.4 Table Encoder Assembly (46-296234G1)

FRU = YES




Figure 2-40 Table Encoder Assembly (46-296234, sheet 1, rev 6)




Figure 2-41 Table Encoder Assembly (46-296234, sheet 2, rev 2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	46-296234	0.000		CRADLE ENCODER ASSEMBLY
2	46-296854P2	1.000	Yes	2.5 INCH, 2048 COUNT.REV
13	46-170623P2	AR	Yes	BLACK, HEAVY-DUTY LUBRICA
29	46-170686P2	AR	Yes	LOCTITE 242, 10CC
33	46-208758P1	1.000	Yes	3.62" X .094" CABLE TIE
3	46-296437P1	1.000	No	BRACKET
4	46-296168P1	1.000	No	DRUM
5	46-296169P1	1.000	No	WHEEL
6	46-296893P1	1.000	No	SHFT,SPG,REE
7	46-296171P1	1.000	No	SPRING COIL
8	46-296699P1	1.000	No	COVER
9	46-296173P1	3.000	No	SPACER
10	46-296436P1	1.000	No	ENC MTG PLT
11	46-208562P13	5.000	No	SCREW
12	46-208562P24	6.000	No	SCREW
14	46-208894P1	2.000	No	SCREW
15	46-296946P5	1.000	No	ALUMINUM TUBE, 6061-T6, .
16	46-296435P1	1.000	No	POTEN. BRACKET
17	46-297036G2	1.000	No	POT ASSEMBLY
18	46-296946P2	1.000	No	SPACER
19	46-296438P1	1.000	No	SPROCKET MOD
20	46-296627P4	1.000	No	SPROCKET
21	46-286100P5	1.000	No	75 PITCHES, 7.365 INCH LE
22	46-297539P2	1.000	No	NO. 6 HI-COLLAR SPRING,
23	46-296781P1	1.000	No	CABLE
24	46-170497P42	1.000	No	SCREW
25	12-1723	1.000	No	ID 0.1590 OD 0.7500 THK .

Table 2-26 Table Encoder Assembly (46-296234G1 BOM, rev 9)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
26	46-296780P2	4.000	No	SYN-CLAMP
28	46-208735P38	1.000	No	SPLIT LENGTHWISE .094 DIA
30	46-170528P23	1.000	No	SHAFT SIZE .250 MAT 2024-
31	46-170528P70	1.000	No	CLAMP
32	46-208562P36	2.000	No	SCREW
34	2149553	1.000	No	TEFLON SPACER .010 THK, 2
35	46-297539P2	1.000	No	NO. 6 HI-COLLAR SPRING,

Table 2-26 Table Encoder Assembly (46-296234G1 BOM, rev 9) (Continued)

USED ON:
2271241

Figure 2-42 Power Assembly (46-296791, sheet 1, rev 8)




Figure 2-43 Power Assembly (46-296791, sheet 2, rev 8)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
4	46-229455P1	1.000	Yes	FAN TUBEAXIAL
5	46-297093P1	3.000	Yes	SPST, AC/DC TOGGLE S
6	46-296316P2	1.000	Yes	EMC POWER SUPPLY SINGLE O
7	46-296317P3	1.000	Yes	OBC-STC & TABLE EMC POWER
10	46-221914P3	1.000	Yes	INDICATOR LIGHT W/ BUILT-
13	46-297088G1	1.000	Yes	POWER ASSEMBLY HARNESS
16	46-220489P2	1.000	Yes	FAN ACCESSORIES FINGER GU
23	46-208758P1	1.000	Yes	3.62" X .094" CABLE TIE
26	46-208759P6	1.000	Yes	7.81" X .184" SELF-LOCKIN
28	46-170053P13	1.000	Yes	3PDT, 12V DC RELAY (EN
2	46-296973P1	1.000	No	POWER SUPPORT BRACKET
3	46-170497P6	8.000	No	SCREW
8	403B335P17	1.000	No	02 POSITION TERMINAL STRI
9	46-208734P26	5.000	No	SCREW
11	46-297132P1	1.000	No	GROUND BAR
12	46-170012P63	3.000	No	0010-32 HEX .375 NUT HT
14	512A720P16	4.000	No	ID .219 OD .562 THK .11
15	1000904P329	2.000	No	ID 0.1960 OD 0.6250 THK .
17	46-297734G1	1.000	No	RESISTOR ASM
18	403B320P94	1.000	No	2 POSITION MARKER STRIP,
19	46-170497P38	2.000	No	SCREW
20	46-278446P1	1.000	No	15A, 125VAC, 2 POLE, 3 WI
21	46-297325P1	1.000	No	FISHPAPER
22	46-297323P1	1.000	No	POWER SUPPLY BRACKET
24	46-297340P1	1.000	No	COVER
25	46-170012P61	3.000	No	0006-32 HEX .313 NUT HT
27	46-170171P22	AR	No	TAPE OLD

Table 2-27 Power Assembly (46-296791G2 BOM, rev 2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
29	46-170497P58	2.000	No	SCREW
30	1000904P33	6.000	No	ID 0.173 OD 0.375 THK .
31	46-271110P5	1.000	No	LABEL OLD
32	46-271110P1	1.000	No	LABEL, PROTECTIVE EARTH (
33	46-208560P53	2.000	No	SCREW
35	46-170497P64	1.000	No	SCREW
36	46-220181P6	1.000	No	1/4" SCREW SIZE (0.250),
37	46-220181P5	10.000	No	WASHER #10 (.190) SCREW S
38	46-170497P7	10.000	No	SCREW

Table 2-27 Power Assembly (46-296791G2 BOM, rev 2) (Continued)

2.4.6 Elevation Encoder Assembly (46-296633G3)

FRU = YES




Figure 2-44 Elevation Encoder Assembly (46-296633G3, sheet 1, rev 8)




Figure 2-45 Elevation Encoder Assembly (46-296633G3, sheet 2, rev 8)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	46-296633	0.000		ELEVATION ENCODER ASSEMBL
2	46-296854P2	1.000	Yes	2.5 INCH, 2048 COUNT.REV
10	46-136334P57	2.000	Yes	SPDT, MINIATURE BASIC
17	46-296946P6	1.000	Yes	ALUMINUM TUBE, 6061-T6, .
3	46-296772P1	1.000	No	BLOCK
4	46-296780P2	3.000	No	SYN-CLAMP
5	46-296631P1	1.000	No	ENCODER SHAFT
6	46-296259P4	1.000	No	ADAPTER
7	46-170048P8	1.000	No	MIN-E-PITCH, 16 T, PD 0.8
9	46-296944P1	1.000	No	SWITCH BRACKET
11	46-136282P75	2.000	No	BEARING
12	46-220277P7	1.000	No	EXT RADIAL MOUNT 0.250 SH
13	46-296945P1	1.000	No	NUT PLATE
14	46-170497P7	2.000	No	SCREW
15	46-170498P125	3.000	No	SCREW
16	46-170528P23	1.000	No	SHAFT SIZE .250 MAT 2024-
18	46-170020P63	2.000	No	NOMINAL DIA 0.250 X 0.625
19	46-220181P2	1.000	No	#4 SCREW SIZE (0.112), ST
20	1000904P438	2.000	No	ID 0.12, OD 0.312, THK .0
21	46-296169P1	1.000	No	WHEEL
22	46-208894P1	2.000	No	SCREW
23	46-297539P2	1.000	No	NO. 6 HI-COLLAR SPRING,
24	1000904P519	2.000	No	WASHER .117 ID .190 OD X
25	46-220205P1	1.000	No	ID 0.3125 OD 0.75 THK
26	46-296636P3	1.000	No	ELEVATION SWITCH CAM FOR

Table 2-28 Elevation Encoder Assembly (46-296633G3 BOM, rev 2)

This Page Intentionally Blank.








Figure 2-49 Frame Assembly, Table (2271240, sheet 4, rev 1)



MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
-	46-296561P11	-	Yes	MAGNET ROD
4	2271256	1.000	Yes	SIGNAL HARNESS FOR H2 TAB
6	2294322	1.000	Yes	ACTUATOR, LINEAR
8	2100671	2.000	Yes	SUSPA # 016-03570 C (ZEUS
13	46-296771P1	1.000	Yes	ROD END SPACER
14	46-278310P1	1.000	Yes	ACTUATOR PIN
23	46-297350P1	2.000	Yes	SPACER
35	46-170684P1	AR	Yes	LOCTITE 271, 0.5CC VIAL
42	2266015	1.000	Yes	TABLE POWER HARNESS
50	2138099	2.000	Yes	ZEUS TABLE TOP COVER MTG
63	46-297687P1	1.000	Yes	8 OHM, 4W, 3-1/2" SPE
75	46-208759P4	14.000	Yes	13.93"X .301" SELF-LOCKIN
79	46-297698P1	2.000	Yes	11" LOW PROFILE BEAD, RIB
82	46-208759P1	34.000	Yes	7.81" X .184" SELF-LOCKIN
92	46-208758P9	2.000	Yes	11.08"X .140" SELF-LOCKIN
93	46-208758P5	2.000	Yes	SELF LOCKING CABLE TIE
97	46-297805G1	2.000	Yes	JUMPER CONN.
98	46-297805G2	2.000	Yes	JUMPER CONN.
104	2106866-2	1.000	Yes	.75 IN WIDE INSULATED CAB
105	2114431	1.000	Yes	T/E POWER SUPPLY
106	2114433	1.000	Yes	EMC PROJECT/COVER/TILT EL
107	2113411	1.000	Yes	170 VOLT POWER SUPPLY
108	46-297335P2	2.000	Yes	6.90A, 250V 2-STAGE PWR L
120	46-208562P20	6.000	Yes	SCREW
3	2118716	1.000	No	EMC ETC CPU
5	46-296770P1	1.000	No	ROD END PIN
7	2133760-2	1.000	No	ZEUS TABLE ACTUATOR BASE

Table 2-29 Frame Assembly (2271240 BOM, rev 6)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
9	2271296-2	1.000	No	FRONT LEG FOR H2 HCC CATH
12	46-136342P135	4.000	No	BUSHING
15	2271297-2	1.000	No	REAR LEG FOR H2 HCC CATHA
16	2138098	1.000	No	ZEUS TABLE LEFT SUPPORT C
17	2138097	1.000	No	ZEUS TABLE RIGHT SUPPORT
18	46-278261P1	1.000	No	CHANNEL
19	46-278262P1	1.000	No	V RAIL
20	2140091	2.000	No	ZEUS TABLE TOP FRONT PIVO
21	2140092	2.000	No	ZEUS TABLE TOP REAR PIVOT
22	46-296811P1	1.000	No	STAGE FIXT
24	46-220179P15	4.000	No	EXT AXIAL MOUNT 0.625 SHA
25	46-220277P15	2.000	No	EXT RADIAL MOUNT 0.312 SH
26	46-220277P11	1.000	No	EXT RADIAL MOUNT 0.625 SH
28	46-328417P26	2.000	No	SCREW HEXAGON SOCKET 10 M
29	46-328417P15	4.000	No	SCREW HEXAGON SOCKET 6 MM
30	46-170498P14	32.000	No	SCREW
32	46-170012P81	1.000	No	005/8-18 HEX .940 NUT
33	46-220181P6	35.000	No	1/4" SCREW SIZE (0.250),
34	46-297055P1	1.000	No	CONN BRACKET
36	46-297435P1	1.000	No	MTG BLOCK
37	46-297436P1	2.000	No	
38	46-297056P1	1.000	No	CONN BRACKET
39	46-297263P1	1.000	No	BRACKET
40	12-128	4.000	No	WASHER,FLAT
41	46-170623P2	AR	No	BLACK, HEAVY-DUTY LUBRICA
43	46-297205P1	6.000	No	BAS CVR BRKT
44	46-297329P2	6.000	No	1/4 TURN RECEPTACLE - FIX
45	46-220186P20	12.000	No	1/8" 120 DEGREE COUNTERSU

Table 2-29 Frame Assembly (2271240 BOM, rev 6) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
46	46-296889P1	1.000	No	BAR
47	46-296890P1	1.000	No	ROLLER MOUNT
48	46-170020P2	2.000	No	NOMINAL DIA 0.0938X 0.500
49	12-275	4.000	No	WASHER, ID .2656, OD.75,
51	46-297677P1	1.000	No	SPRING CLIP
52	46-297705P2	1.000	No	CONDUCTIVE TAPE
53	46-297705P14	2.000	No	METALLIC CONDUCTIVE TAPE
54	46-297705P4	1.000	No	MET-CON TAPE
55	46-170497P66	1.000	No	SCREW
61	46-297343P2	1.000	No	ZEUS TABLE CABLE TROUGH
62	46-297684P1	1.000	No	SPKR GRILL
64	46-297686P2	1.000	No	SPEAKER INSULATION FOR H2
65	46-297685P2	1.000	No	SPEAKER COVER FOR H2 HCC
66	46-297155G1	4.000	No	ZEUS TABLE, 9" LONG CABLE
69	46-297427P1	3.000	No	GROUNDING STRAP CLAMP
77	46-297648P1	1.000	No	CABLE CLAMP
80	2146815	25.000	No	10-32 X 1/2 HEX SOC BUTTO
81	1000904P395	4.000	No	ID 0.177 OD 0.312 THK .
83	46-208560P46	1.000	No	SCREW
84	46-170497P38	4.000	No	SCREW
85	46-170012P63	24.000	No	0010-32 HEX .375 NUT HT
86	46-170012P61	3.000	No	0006-32 HEX .313 NUT HT
87	2275854	1.000	No	GROUND - TABLE TO RACEWAY
88	2146815-2	22.000	No	10-32 X 3/8 HEX SOC BUTTO
89	46-297731P1	1.000	No	MOTOR LABEL
90	12-276	4.000	No	ID 0.2040 OD 0.5630 THK .
91	46-170497P36	5.000	No	SCREW
94	2138105	1.000	No	TOP COVER CONNECTION BRAC

Table 2-29 Frame Assembly (2271240 BOM, rev 6) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
95	2138105-2	1.000	No	TOP COVER CONNECTION BRAC
96	46-208758P1	4.000	No	3.62" X .094" CABLE TIE
99	46-208562P7	12.000	No	SCREW
100	46-170012P62	4.000	No	0008-32 HEX .343 NUT HT
101	46-297858P1	4.000	No	SHOULDER SCREW
112	46-297705P10	1.000	No	METALLIC CONDUCTIVE 3M-SC
113	2132187-2	1.000	No	EMC COST REDUCED TABLE BA
114	46-328417P12	38.000	No	SCREW HEXAGON SOCKET 6 MM
115	46-170015P3	4.000	No	SCREW SEMS 8-32 X .312, E
117	2136121	2.000	No	ZEUS TABLE LOWER LEG PIVO
118	2138104	2.000	No	REAR LEG LOWER PIVOT PIN
119	46-170498P32	6.000	No	SCREW
121	46-328416P23	4.000	No	SCREW HEXAGON HEAD 10 MM
122	46-328431P3	4.000	No	WASHER PLAIN - LARGE 10.5
123	46-328432P4	4.000	No	WASHER LOCK - SPRING 10.2
124	46-170686P2	AR	No	LOCTITE 242, 10CC
125	46-170015P27	5.000	No	SCREW
126	12-297	5.000	No	ID 0.140 OD 0.312 THK .
127	46-220181P5	6.000	No	WASHER #10 (.190) SCREW S

Table 2-29 Frame Assembly (2271240 BOM, rev 6) (Continued)

This Page Intentionally Blank.

Section 3.0 — Footswitch Assembly (2330508)

FRU = NO




Figure 2-51 Footswitch Assembly (2330508 sheet 1, rev 2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	46-296129G1	2.000	Yes	
28	2329241	4.000	Yes	FOOT PEDAL
29	2297598	1.000	Yes	FOOT SWITCH COVER
6	46-170497P8	6.000	No	SCREW
7	46-271110P1	1.000	No	LABEL, PROTECTIVE EARTH (GROUND)
8	2329242	2.000	No	PEDAL ROD
9	2329240	1.000	No	GROUND FOOT SWITCH
10	46-221828P4	2.000	No	0.063 DIA X 0.625 LG 18-8 SST
11	46-170012P63	3.000	No	0010-32 HEX .375 NUT HT .239
13	46-208560P65	11.000	No	SCREW
14	46-220181P6	11.000	No	1/4" SCREW SIZE (0.250), STL NI PLT
15	1000274P289	3.000	No	
17	46-297971P1	2.000	No	TERMINAL
19	46-328417P12	2.000	No	SCREW HEXAGON SOCKET 6 MM 12 MM
20	46-328430P3	2.000	No	WASHER PLAIN - NORMAL 6.4 MM 12 MM
21	2329239	2.000	No	GUARD BAR
22	2329246	6.000	No	SPACER GROUND BAR
23	46-170498P11	6.000	No	SCREW
24	46-170686P3	AR	No	LOCTITE 242, 50CC
25	2329247	2.000	No	SPRING PLATE
26	46-170497P20	12.000	No	SCREW
27	46-136353P135	4.000	No	SPRING-INCH COMPRESSION 1.75 INCHES
30	2297597	1.000	No	FOOT SWITCH BASE

Table 2-30 Footswitch Assembly (2330508 BOM, rev 5)

$FRU = NO$

Figure 2-52 Foot Switch Cover Assembly (46-296790, sheet 1, rev 7)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	46-296129G1	2.000	Yes	
28	46-296162P2	4.000	Yes	FOOT PEDAL FOR H2 TABLE
29	46-296161P2	1.000	Yes	FOOT SWITCH COVER FOR H2 TABLE
6	46-170497P8	6.000	No	SCREW
7	46-271110P1	1.000	No	LABEL, PROTECTIVE EARTH (GROUND)
8	46-296163P1	2.000	No	
9	46-297053P1	1.000	No	
10	46-221828P4	2.000	No	0.063 DIA X 0.625 LG 18-8 SST
11	46-170012P63	3.000	No	0010-32 HEX .375 NUT HT .239
12	46-271110P5	1.000	No	LABEL OLD
13	46-208560P65	11.000	No	SCREW
14	46-220181P6	11.000	No	1/4" SCREW SIZE (0.250), STL NI PLT
15	1000274P289	3.000	No	
17	46-297971P1	2.000	No	TERMINAL
19	46-328417P12	2.000	No	SCREW HEXAGON SOCKET 6 MM 12 MM
20	46-328430P3	2.000	No	WASHER PLAIN - NORMAL 6.4 MM 12 MM
21	2151201	2.000	No	HSA FOOT SWITCH GUARD BAR
22	2151202	6.000	No	HSA FT SW GUARD BAR SPACER
23	46-170498P11	6.000	No	SCREW
24	46-170686P3	AR	No	LOCTITE 242, 50CC
25	2153419	2.000	No	HSA FOOT SWITCH PEDAL SPRING PLATE
26	46-170497P20	12.000	No	SCREW
27	46-136353P135	4.000	No	SPRING-INCH COMPRESSION 1.75 INCHES
30	2147702-2	1.000	No	FOOT SWITCH BASE FOR H2 TABLE

Table 2-31 Foot Switch Cover Assembly (46-296790G3 BOM, rev 1)

End of Chapter.