
Chapter 5

CONSOLE

Section 1.0 : AIO console Assembly for BSD delight lite (5377708-560), Page376.

Section 2.0 : GRE CONSOLE, 16 SLICE GDAS FOR BS WITH XW8200, ETHERNET SWITCH AND JARRELL (5377708-530), Page402.

Section 3.0 : GRE CONSOLE, 4/8 SLICE GDAS FOR BS WITH XW8200, ETHERNET SWITCH AND JARRELL (5377708-520), Page425.

Section 4.0 : Sigle Piece Desktop Kit for Freedom Workspace (5166834), Page439.

Section 5.0 : Freedom Workspace Side by Side Kit (5166569), Page441.

Section 6.0 : Associated Items, Page445.

Section 7.0 : Options, Page454.

Section 1.0 — AIO console Assembly for BSD delight lite (5377708-560)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	2337414	2.0	Yes	DVD RAM MEDIA, 9.4 GB
3	2368697-2	1.0	Yes	LINUX DDR WORKSTATION 3BUTTON PS2 MOUSE
1	5378862-560	1.0	No	AIO in one console for BSD delight
4	5317250	1.0	No	AIO console rating plate for BSD lite

Table 5-1 AIO console Assembly for BSD delight lite (5377708-560 BOM, Rev 3)

1.1 All-In-One Console (5378862-560)

$FRU = NO$




Figure 5-1 All-In-One Console (5378862-560ADW, Sheet1, Rev 1)




Figure 5-2 All-In-One Console (5378862-560ADW, Sheet2, Rev1)




Figure 5-3 All-In-One Console (5378862-560ADW, Sheet3, Rev1)




Figure 5-4 All-In-One Console (5378862-560ADW, Sheet4, Rev1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	5165396-2	1.0	No	FWS GOC4 chassis,make parts.
2	5310601	1.0	No	Host Computer Assembly
3	5309829	2.0	No	Host Bracket
5	2404320-2	1.0	Yes	CONSOLE INTERCOM PRESCRIBED-TILT HUB
6	5116391-2	1.0	Yes	SCSI Tower ASSY
7	5115334	1.0	Yes	Switch, Ethernet 8-port GBe
8	N9501WS	4.0	Yes	WASHER PLAIN - LARGE 5.3 MM 10 MM WASHER M5
9	46-313575P3	4.0	Yes	HIGH COLLAR LOCK WASHER FOR M5 SCREW
11	5128181	1.0	Yes	Cover Assy - Console Front
13	2304737-2	2.0	Yes	CABLE, DVI VIDEO CABLE - OCTANE II
14	2366442-27	2.0	No	BRACKET, RACK MOUNT HOST COMPUTER
17	2324697-3	1.0	Yes	SCSI cable assembly
19	5127222-2	1.0	Yes	Cable, Ethernet, 2ft Lg - CAT5e, 4 twisted pair, 24awg solid copper, Shielded - aluminized polyester and tinned copper, gray PVC jacket, snagless PVC boot, Plugs - UL94V-0 rated polycarbonate with metallic strain relief, EIA-TIA568, UL listed, CMR rated,
25	46-136246P14	60.0	Yes	2~ WIDE, WHT, ADH BACK, VELCRO LOOP 25 YD
26	46-136246P21	60.0	Yes	2~ WIDE X 25 YD WHITE VELCRO HOOKS
27	46-208781P2	4.0	Yes	WRITE-ON LABEL FOR MARKING CABLES WHITE ADHESIVE BACKED 0.75~ X 1.90~. UNMARKED, 2,500 PC FAN-FOLD SHEETS (FULL SIZE = 1.9~ X 3.15~.)
28	5133331	1.0	No	ETL Label for OC
29	2117390-2	1.0	Yes	CE Compliance and WEEE Labels
30	2132620	1.0	Yes	EMC LABEL CLASS A GROUP 1
33	2240540	1.0	No	NETSCAPE LICENSE
41	46-312358P51	8.0	No	M5 X 0.8 X 12MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW
42	2103585-2	4.0	Yes	SCREW PAN HEAD 4 MM 10 MM SPRING WASHER, ZINC PLATED, STEEL

Table 5-2 All-In-One Console (5378862-560 BOM, Rev 2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
43	46-208758P2	5.0	Yes	5.5~ X .140~ SELF-LOCKING CABLE TIE. 1/16~ TO 1-1/4~ WIRE BUNDLE RANGE. 30LB TENSILE STRENGTH, NYLON MATR~L.
44	46-208747P1	1.0	Yes	ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)
45	5194414-2	1.0	No	Cable Labels for GRE
46	2385779	2.0	No	PANEL - COVER
47	46-208758P3	2.0	Yes	7.31~ X .184~ SELF-LOCKING CABLE TIE 1/16~ TO 1-3/4~ WIRE BUNDLE RANGE. 50LB TENSILE STRENGTH, NYLON MATR~L.

Table 5-2 All-In-One Console (5378862-560 BOM, Rev 2)

$FRU = NO$

Figure 5-5 FWS GOC4 chassis Assembly (5165396-2ADW, refer to 5165396 ADW, Sheet1, Rev3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	5164805-2	1.0	No	FWS GOC4 chassis from supplier, buy parts.
2	5160898	1.0	Yes	SubAssy, GRE Console Rear Bulkhead
3	2315117	1.0	Yes	CABLE, PS/2 MOUSE CABLE, LINUX
4	2315117-2	1.0	Yes	CABLE, PS/2 KEYBOARD CABLE, LINUX
5	5123478	1.0	Yes	Cable, USB Bulkhead, USB-Af to USB-Am x 3ft Lg
6	46-328421P3	6.0	Yes	SCREW PAN HEAD 2.5 MM 5 MM ZINC PLATED, STEEL
7	2142225	4.0	Yes	MONITOR POWER CORD
8	2142226	5.0	Yes	POWER CORD FOR EMC AREA
9	5124071	4.0	Yes	Fan, Tube Axial - 115VAC, 50-60Hz, 119mm x 119mm x 25mm x 95-115cfm, ball bearing motor - Sinwan S125AP-11-1
10	5162696	1.0	No	GOC4 Outlet Box Assy
11	2385780	1.0	No	TABLE TOP ASSEMBLY, MONITOR
13	46-208758P3	22.0	Yes	7.31~ X .184~ SELF-LOCKING CABLE TIE 1/16~ TO 1-3/4~ WIRE BUNDLE RANGE. 50LB TENSILE STRENGTH, NYLON MATR~L.
14	2103585	19.0	Yes	SCREW PAN HEAD 4 MM 8 MM SPRING WASHER, ZINC PLATED, STEEL
15	46-328431P2	8.0	Yes	WASHER PLAIN - LARGE 6.4 MM 18 MM ZINC PLATED, STEEL
16	2109881-3	8.0	Yes	WASHER LOCK - INTERNAL 6.4 MM 11 MM
17	2337083-2	8.0	No	SCREW, M6 X 20 SOC HEAD, 12.9 CLASS
20	46-221417P4	2.0	Yes	FEMALE SCREWLOCKS FOR HD CONNECTORS. 4-40 THREAD, 0.470 OVERALL LENGTH YELLOW CHROMATE, TWO PER KIT. 0.312 BOTTOM, 0.094 PANEL THK. MAX KIT INCLUDES NUTS AND WASHERS See drawing for approved vendor
21	46-208747P1	1.0	Yes	ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)
22	2150839	1.0	No	CAUTION RECEPTACLES NOT FOR GEN USE
23	46-328431P1	10.0	Yes	WASHER PLAIN - LARGE 4.3 MM 12 MM ZINC PLATED, STEEL
25	5114923-2	1.0	No	Bracket, ICOM Lower
26	5114924-2	1.0	No	Bracket, ICOM Upper
27	5122710	2.0	No	M4 Tapped Hole Speed Nut
28	46-328417P5	2.0	Yes	SCREW HEXAGON SOCKET 4 MM 16 MM ZINC PLATED, STEEL

Table 5-3 FWS GOC4 chassis Assembly (5165396-2 BOM, Rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
29	46-272301P11	4.0	Yes	FAN ACCESSORIES FINGER GUARD PAPST LZ30
30	5127629	1.0	Yes	Cable Assy - Exhaust Fan Power
31	2381089-3	16.0	Yes	SCREW, SS NYLOCK HEX NUT, M4
32	46-328430P2	16.0	Yes	WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL
33	46-312358P44	16.0	Yes	M4 X 0.7 X 35MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW
34	5128406	1.0	Yes	Assy, ICOM cable harness for GOC4
37	46-170686P3	0.0	Yes	LOCTITE 242, 50CC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER.
38	5166489	1.0	No	GOC4 chassis label

Table 5-3 FWS GOC4 chassis Assembly (5165396-2 BOM, Rev 1)

1.1.1.1 FWS GOC4 CHASSIS (5164805-2)

FRU = NO




Figure 5-6 FWS GOC4 chassis (5164805-2ADW, Sheet1, Rev2)

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	2385112-3	1.0	No	GRE,CHASSIS SUB-ASSEMBLY
2	2385763	1.0	No	REAR MTG PANEL, OUTLET BOX
3	5123480	1.0	No	Bracket, Cable Cover - Dual PS2f, Single USB-Af
5	2319043	4.0	Yes	CASTER, 2 1/2~ SWIVEL
6	2103585	24.0	Yes	SCREW PAN HEAD 4 MM 8 MM SPRING WASHER, ZINC PLATED, STEEL
9	2109877-2	16.0	Yes	NUT HEXAGON 6 MM
11	2109881-3	25.0	Yes	WASHER LOCK - INTERNAL 6.4 MM 11 MM
13	2180698-2	38.0	Yes	CLIP NUT ASSEMBLY, M5
14	2385765	4.0	No	BRACKET, SIDE COVER MTG
15	2385766	3.0	No	BRACKET, FRONT COVER MTG
16	5123484	1.0	Yes	Bracket, Exhaust Fan Mtg - Left
17	5123484-2	1.0	Yes	Bracket, Exhaust Fan Mtg - Right
18	2385773	2.0	No	PANEL, RACK MOUNT
19	2385776	10.0	No	BRACKET, RACK MOUNT
20	5123485	1.0	Yes	Rear Access Cover - VCT
21	46-328431P2	25.0	Yes	WASHER PLAIN - LARGE 6.4 MM 18 MM ZINC PLATED, STEEL
22	46-170686P3	0.0	Yes	LOCTITE 242, 50CC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER.
23	2376906	6.0	No	SOUTHCO 1/4 TURN RCVR 85-35-295-20
24	2337083-2	1.0	No	SCREW, M6 X 20 SOC HEAD, 12.9 CLASS
25	2385779	5.0	No	PANEL - COVER
26	46-312358P51	30.0	No	M5 X 0.8 X 12MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW
27	46-328430P3	4.0	Yes	WASHER PLAIN - NORMAL 6.4 MM 12 MM ZINC PLATED, STEEL
28	N9601BS	4.0	Yes	NUT HEXAGON 6 MM NUT M6 2 FE ZN

Table 5-4 FWS GOC4 chassis (5164805-2 BOM, Rev 2)

REAR ACCESS COVER - VCT (5123485)

FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	46-272301P11	4.0	Yes	FAN ACCESSORIES FINGER GUARD PAPST LZ30
3	2381089-3	16.0	Yes	SCREW, SS NYLOCK HEX NUT, M4
4	5125042	4.0	No	Gasket, Foam - Exhaust Fan
5	5124020	3.0	Yes	Foam, Sound Absorbing - 54in Wd x 1in Thick
6	5131360	1.0	No	Label, Keyboard Mouse Trackball

Table 5-5 Rear Access Cover - VCT (5123485 BOM, Rev6)

1.1.1.2 SUBASSY, GRE CONSOLE REAR BULKHEAD (5160898)

FRU = YES




Figure 5-7 SubAssy, GRE Console Rear Bulkhead (5160898ADW, Sheet1, Rev3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2265406	1.0	Yes	FIBER-OPTIC COUPLING RECEPTACLE KIT, 2.5MM BAYONET, AMP OPTIMATE
2	5160899	1.0	No	Bulkhead Panel, GRE Console Rear Cable Connections
3	2374957	1.0	No	ADAPTER - BULKHEAD BNC JACK-JACK
4	5122589	6.0	Yes	Adapter, Bulkhead RJ45 CAT5e Shielded, L-com #TDG1026KS-C5E
6	2318713-2	1.0	No	CABLE-FIBER OPTIC ST-SC X 3FT
8	5127222-3	3.0	Yes	Cable, Ethernet x 3ft Lg - CAT5e, 4 twisted pair, 24awg solid copper, Shielded - aluminized polyester and tinned copper, gray PVC jacket, snagless PVC boot, Plugs - UL94V-0 rated polycarbonate with metallic strain relief, EIA-TIA568, UL listed, CMR rated,
9	5128409	1.0	Yes	Assy, Gbe cable harness for VCT console

Table 5-6 SubAssy, GRE Console Rear Bulkhead (5160898 BOM, Rev4)

ASSY, GBE CABLE HARNESS (5128409)

FRU = YES




Figure 5-8 Assy, Gbe cable harness (5128409ADW, Sheet1, Rev3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	5127222-3	3.0	Yes	Cable, Ethernet x 3ft Lg - CAT5e, 4 twisted pair, 24awg solid copper, Shielded - aluminized polyester and tinned copper, gray PVC jacket, snagless PVC boot, Plugs - UL94V-0 rated polycarbonate with metallic strain relief, EIA-TIA568, UL listed, CMR rated,
2	5127222-5	1.0	No	Cable, Ethernet x 5ft Lg - CAT5e, 4 twisted pair, 24awg solid copper, Shielded - aluminized polyester and tinned copper, gray PVC jacket, snagless PVC boot, Plugs - UL94V-0 rated polycarbonate with metallic strain relief, EIA-TIA568, UL listed, CMR rated,
3	46-287814P2	1.0	Yes	BRAIDED EXPANDABLE SLEEVING, 1/2~

Table 5-7 Assy, Gbe cable harness (5128409 BOM, Rev4)

1.1.1.3 GOC4 Outlet Box Assy (5162696)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2174159-5	1.0	Yes	Assy, Outlet Box for GOC4
2	5162697	1.0	No	Lead Data Wire for GOC4

Table 5-8 GOC4 Outlet Box Assy (5162696 BOM, Rev 2)

Assy, Outlet Box for GOC4 (2174159-5)

FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2367803	1.0	No	WELDMENT - OUTLET BOX (GRE)
2	46-278446P1	8.0	Yes	15A, 125VAC, 2 POLE, 3 WIRE DUPLEX TYPE, GRAY IN COLOR. (REV TK, 7/91).
3	46-297335P3	1.0	Yes	16A, 250V 2-STAGE PWR LINE FILTER
4	46-243925P1	1.0	Yes	20A, 125V, 2-POLE, 3 WIRE TWIST LOCK FLANGED INLET (MALE) LOCKING TYPE (REV TK, 10/91).
6	46-328425P10	16.0	Yes	NUT HEXAGON 3 MM
7	46-328425P1	3.0	Yes	NUT HEXAGON 4 MM
8	46-312358P30	3.0	No	M3 X 0.5 X 12MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW
9	46-328430P2	3.0	Yes	WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL
11	46-208697P57	1.0	No	HEYCO PART NO.1247 SR-7W-2 BLACK
12	46-311804P2	16.0	Yes	3.2X9X.8 ZP ST DIN
13	46-311805P10	16.0	Yes	HELICAL LOCK WASHER M3 ZP DIN
14	46-208758P2	3.0	Yes	5.5~ X .140~ SELF-LOCKING CABLE TIE. 1/16~ TO 1-1/4~ WIRE BUNDLE RANGE. 30LB TENSILE STRENGTH, NYLON MATR~L.
15	2336074	1.0	No	CIRCUIT BREAKER 15 AMPS 240 VOLTS 80
16	46-328421P501	2.0	Yes	SCREW PAN HEAD 3 MM 5 MM
17	5154767	1.0	No	SWITCH BOX ASSY for GRE with Service ethernet port

Table 5-9 Assy, Outlet Box for GOC4 (2174159-5 BOM, Rev 2

1.1.2 HOST COMPUTER ASSEMBLY (5310601)

FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	5183547-50	1.0	Yes	AIO console HOST computer
2	5220229	1.0	Yes	CDIP CARD

Table 5-10 Host Computer Assembly (5310601 BOM, Rev 2)

1.1.3 CONSOLE INTERCOM PRESCRIBED-TILT HUB (2404320)

FRU = YES




Figure 5-9 Console Intercom Prescribed-tilt Hub (2404320ADW, Sheet1, Rev11)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2370567	1.0	No	POWER ENTRY MODULE 2 POLE FUSING
2	2374174	1.0	No	POWER ENTRY FUSE DRAWER
3	2374694	2.0	Yes	8A TIME DELAY FUSES
4	2399750	1.0	No	GRE INTERCOM-TILT
5	5114920-2	1.0	No	Chassis Base, ICOM
6	5114921	1.0	No	ct Intercom Box Chassis Cover
7	5115027	1.0	No	Power supply triple output 90-264V 47-63 HZ AC IN 5V 7A +12V 3A -12V 2.5A OUT
8	5115179	1.0	No	cable harness power entry to power supply
9	5115180	1.0	No	cable harness power supply to pwa
10	46-312577P14	22.0	Yes	SCREW HEAD CAP SCREW STEEL ZINC PLATED (FULL THREAD) PROPERTY CLASS 8.8 ***THIS PART IS REPLACED BY 46-328417P3***
11	46-328430P2	16.0	Yes	WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL
12	2276462	6.0	No	METRIC PLAIN FLAT WASHER 3.2MM ID X 7MM OD
13	46-312346P1	2.0	Yes	M3 X 0.5 ZINC PLATED METRIC HEX NUT WITH ATTACHED TOOTH LOCK WASHER
14	46-312577P4	2.0	No	SCREW HEAD CAP SCREW STEEL ZINC PLATED (FULL THREAD) PROPERTY CLASS 8.8
15	46-214941P2	6.0	Yes	FLAT WASHER, NYLON, M4
16	2209282-3	5.0	No	HEX STANDOFF MALE/FEMALE M4 X 25MM
17	5128595	1.0	No	M4 x 25mm male female nylon hex stand off.
18	46-328417P73	4.0	No	SCREW HEXAGON SOCKET 3 MM 6 MM
19	46-221417P1	3.0	Yes	FEMALE SCREWLOCKS, TWO PER KIT FOR HD CONNECTORS 4-40 THREAD 0.5 OVERALL LENGTH, YELLOW CHROMATE 0.312 BOTTOM, 0.062 PANEL THK. MAX KIT INCLUDES NUTS AND WASHERS

Table 5-11 Console Intercom Prescribed-tilt Hub (2404320 BOM, Rev14)

FRU = YES

Figure 5-10 SCSI Tower ASSY (5116391ADW, Sheet1, Rev1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2398863	1.0	Yes	DVD-RAM MULTI-DRIVE 4.7GB
2	2404031	1.0	No	SCSI TO IDE ADAPTER
4	2401055	1.0	No	SCSI TOWER
5	46-328421P501	1.0	Yes	SCREW PAN HEAD 3 MM 5 MM

Table 5-12 SCSI Tower ASSY (5116391-2 BOM, Rev3)

1.1.5 COVER ASSY - CONSOLE FRONT (5128181)
FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2369619	1.0	Yes	COVER, CONSOLE FRONT - GRE
2	5124020	1.0	Yes	Foam, Sound Absorbing - 54in Wd x 1in Thick
3	46-136246P18	3.21	Yes	1~ WIDE, WHT, ADH BACK, VELCRO HOOK 75 FT PER RL
4	5128180	1.0	Yes	Filter, 25 PPI Quadrafoam UL94 HF-1

Table 5-13 Cover Assy - Console Front (5128181 BOM, Rev1)

1.2 AIO 5377708-560 RATING PLATE (5317250)


FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	46-303268P1	1.0	Yes	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH

Table 5-14 AIO 5377708-560 Rating Plate (5317250 BOM, Rev 1)

Section 2.0 — GRE CONSOLE, 16 SLICE GDAS FOR BS WITH XW8200, ETHERNET SWITCH AND JARRELL (5377708-530)

FRU = NO






Figure 5-12 GRE Console, 16 Slice Gdas For Bs With Xw8200, Ethernet Switch And Jarrell (5377708-530ADW, Sheet 2, Rev 1)




Figure 5-13 GRE Console, 16 Slice Gdas For Bs With Xw8200, Ethernet Switch And Jarrell (5377708-530ADW, Sheet 3, Rev 1)




Figure 5-14 GRE Console, 16 Slice Gdas For Bs With Xw8200, Ethernet Switch And Jarrell (5377708-530ADW, Sheet 4, Rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	5117866-51	1.0	Yes	GEHC CT VCT INTEGRATED BUY LEVEL PC HP8200 WITH FX1500 GRAPHICS CARD HOST COMPUTER
3	5147442	1.0	Yes	COMPUTER - DATA ACQUISITION RECON CONTROL, 2ND GENERATION
5	2404320	1.0	Yes	CONSOLE INTERCOM PRESCRIBED-TILT HUB
6	5116391-2	1.0	Yes	SCSI TOWER ASSY
7	5115334	1.0	Yes	SWITCH, ETHERNET 8-PORT GBE
8	2368697-2	1.0	Yes	LINUX DDR WORKSTATION 3BUTTON PS2 MOUSE
11	5128181	1.0	Yes	COVER ASSY - CONSOLE FRONT
13	2304737-2	2.0	Yes	CABLE, DVI VIDEO CABLE - OCTANE II
17	2324697-2	1.0	Yes	SCSI CABLE, VHD68 TO CENT 50, LINUX
19	5127222-2	2.0	Yes	CABLE, ETHERNET, 2FT LG - CAT5E, 4 TWISTED PAIR, 24AWG SOLID COPPER, SHIELDED - ALUMINIZED POLYESTER AND TINNED COPPER, GRAY PVC JACKET, SNAGLESS PVC BOOT, PLUGS - UL94V-0 RATED POLYCARBONATE WITH METALLIC STRAIN RELIEF, EIA-TIA568, UL LISTED, CMR RATED,
25	46-136246P14	60.0	Yes	2~ WIDE, WHT, ADH BACK, VELCRO LOOP 25 YD
26	46-136246P21	60.0	Yes	2~ WIDE X 25 YD WHITE VELCRO HOOKS
27	46-208781P2	4.0	Yes	WRITE-ON LABEL FOR MARKING CABLES WHITE ADHESIVE BACKED 0.75~ X 1.90~. UNMARKED, 2,500 PC FAN-FOLD SHEETS (FULL SIZE = 1.9~ X 3.15~.)
29	2117390-2	1.0	Yes	CE COMPLIANCE AND WEEE LABELS
30	2132620	1.0	Yes	EMC LABEL CLASS A GROUP 1
36	2337414	2.0	Yes	DVD RAM MEDIA, 9.4 GB
39	46-208758P3	6.0	Yes	7.31~ X .184~ SELF-LOCKING CABLE TIE 1/16~ TO 1-3/4~ WIRE BUNDLE RANGE. 50LB TENSILE STRENGTH, NYLON MATR~L.
42	2103585-2	4.0	Yes	SCREW PAN HEAD 4 MM 10 MM SPRING WASHER, ZINC PLATED, STEEL
43	46-208747P1	1.0	Yes	ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)
44	46-208758P2	1.0	Yes	5.5~ X .140~ SELF-LOCKING CABLE TIE. 1/16~ TO 1-1/4~ WIRE BUNDLE RANGE. 30LB TENSILE STRENGTH, NYLON MATR~L.

Table 5-15 GRE Console, 16 Slice Gdas For Bs With Xw8200, Ethernet Switch And Jarrell (5377708-530 BOM, Rev 2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	5165396	1.0	No	FWS GOC4 CHASSIS,MAKE PARTS.
14	2366442-27	2.0	No	BRACKET, RACK MOUNT HOST COMPUTER
23	5184117	1.0	No	HIVALUE LIGHT 16 JARRELL RATING PLATE
28	5133331	1.0	No	ETL LABEL FOR OC
33	2240540	1.0	No	NETSCAPE LICENSE
41	46-312358P51	14.0	No	M5 X 0.8 X 12MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW

Table 5-15 GRE Console, 16 Slice Gdas For Bs With Xw8200, Ethernet Switch And Jarrell (5377708-530 BOM, Rev 2)

2.1 FWS GOC4 CHASSIS ASSEMBLY (5165396)

FRU = NO




Figure 5-15 FWS GOC4 Chassis Assembly (5165396ADW, Sheet 1, Rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION
2	5160898	1.0	Yes	SUBASSY, GRE CONSOLE REAR BULKHEAD
3	2315117	1.0	Yes	CABLE, PS/2 MOUSE CABLE, LINUX
4	2315117-2	1.0	Yes	CABLE, PS/2 KEYBOARD CABLE, LINUX
5	5123478	1.0	Yes	CABLE, USB BULKHEAD, USB-AF TO USB-AM X 3FT LG
6	46-328421P3	6.0	Yes	SCREW PAN HEAD 2.5 MM 5 MM ZINC PLATED, STEEL
7	2142225	4.0	Yes	MONITOR POWER CORD
8	2142226	5.0	Yes	POWER CORD FOR EMC AREA
9	5124071	4.0	Yes	FAN, TUBE AXIAL - 115VAC, 50-60HZ, 119MM X 119MM X 25MM X 95-115CFM, BALL BEARING MOTOR - SINWAN S125AP-11-1
13	46-208758P3	22.0	Yes	7.31~ X .184~ SELF-LOCKING CABLE TIE 1/16~ TO 1-3/4~ WIRE BUNDLE RANGE. 50LB TENSILE STRENGTH, NYLON MATR~L.
14	2103585	19.0	Yes	SCREW PAN HEAD 4 MM 8 MM SPRING WASHER, ZINC PLATED, STEEL
15	46-328431P2	8.0	Yes	WASHER PLAIN - LARGE 6.4 MM 18 MM ZINC PLATED, STEEL
16	2109881-3	8.0	Yes	WASHER LOCK - INTERNAL 6.4 MM 11 MM
20	46-221417P4	2.0	Yes	FEMALE SCREWLOCKS FOR HD CONNECTORS. 4-40 THREAD, 0.470~ OVERALL LENGTH YELLOW CHROMATE, TWO PER KIT. 0.312~ BOTTOM, 0.094~ PANEL THK. MAX (KIT INCLUDES NUTS AND WASHERS)
21	46-208747P1	1.0	Yes	ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)
23	46-328431P1	10.0	Yes	WASHER PLAIN - LARGE 4.3 MM 12 MM ZINC PLATED, STEEL
28	46-328417P5	2.0	Yes	SCREW HEXAGON SOCKET 4 MM 16 MM ZINC PLATED, STEEL
29	46-272301P11	4.0	Yes	FAN ACCESSORIES FINGER GUARD PAPST LZ30
30	5127629	1.0	Yes	CABLE ASSY - EXHAUST FAN POWER
31	2381089-3	16.0	Yes	SCREW, SS NYLOCK HEX NUT, M4
32	46-328430P2	16.0	Yes	WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL
33	46-312358P44	16.0	Yes	M4 X 0.7 X 35MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW
34	5128406	1.0	Yes	ASSY, ICOM CABLE HARNESS FOR GOC4
37	46-170686P3	0.0	Yes	LOCTITE 242, 50CC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER.

Table 5-16 FWS GOC4 Chassis Assembly (5165396 BOM, Rev 2)

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	5164805	1.0	No	FWS GOC4 CHASSIS FROM SUPPLIER, BUY PARTS.
10	5162696	1.0	No	GOC4 OUTLET BOX ASSY
11	2385780	1.0	No	TABLE TOP ASSEMBLY, MONITOR
17	2337083-2	8.0	No	SCREW, M6 X 20 SOC HEAD, 12.9 CLASS
22	2150839	1.0	No	CAUTION RECEPTACLES NOT FOR GEN USE
25	5114923-2	1.0	No	BRACKET, ICOM LOWER
26	5114924-2	1.0	No	BRACKET, ICOM UPPER
27	5122710	2.0	No	M4 TAPPED HOLE SPEED NUT
38	5166489	1.0	No	GOC4 CHASSIS LABEL

Table 5-16 FWS GOC4 Chassis Assembly (5165396 BOM, Rev 2)

2.1.1 FWS GOC4 CHASSIS (5164805)

FRU = NO

Figure 5-16 FWS GOC4 Chassis (5164805ADW, Rev 2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
5	2319043	4.0	Yes	CASTER, 2 1/2~ SWIVEL
6	2103585	24.0	Yes	SCREW PAN HEAD 4 MM 8 MM SPRING WASHER, ZINC PLATED, STEEL
9	2109877-2	16.0	Yes	NUT HEXAGON 6 MM
11	2109881-3	25.0	Yes	WASHER LOCK - INTERNAL 6.4 MM 11 MM
13	2180698-2	38.0	Yes	NUT CLIP ON 5 MM CLIP & NUT MAT~L: CARBON STEEL, 1050 PER AMS5085, HEAT TREATED CLIP FINISH: CADMIUM PLATE PER QQ-P-416 TYPE II, CLASS 1. NUT FINISH: CADMIUM PLATE PER QQ-P-416 TYPE II, CLASS 2, MOLYBDENUM BISULFIDE DRY FILM LUBRICANT PER MIL-L-46010
16	5123484	1.0	Yes	BRACKET, EXHAUST FAN MTG - LEFT
17	5123484-2	1.0	Yes	BRACKET, EXHAUST FAN MTG - RIGHT
20	5123485	1.0	Yes	REAR ACCESS COVER - VCT
21	46-328431P2	25.0	Yes	WASHER PLAIN - LARGE 6.4 MM 18 MM ZINC PLATED, STEEL
22	46-170686P3	0.0	Yes	LOCTITE 242, 50CC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER.
1	2385112	1.0	No	GRE, CHASSIS SUB-ASSEMBLY
2	2385763	1.0	No	REAR MTG PANEL, OUTLET BOX
3	5123480	1.0	No	BRACKET, CABLE COVER - DUAL PS2F, SINGLE USB-AF
14	2385765	4.0	No	BRACKET, SIDE COVER MTG
15	2385766	3.0	No	BRACKET, FRONT COVER MTG
18	2385773	2.0	No	PANEL, RACK MOUNT
19	2385776	10.0	No	BRACKET, RACK MOUNT
23	2376906	6.0	No	SOUTHCO 1/4 TURN RCVR 85-35-295-20
24	2337083-2	1.0	No	SCREW, M6 X 20 SOC HEAD, 12.9 CLASS
25	2385779	5.0	No	PANEL - COVER
26	46-312358P51	30.0	No	M5 X 0.8 X 12MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW

Table 5-17 FWS GOC4 Chassis (5164805 BOM, Rev 3)

2.1.2 SUBASSY, GRE CONSOLE REAR BULKHEAD (5160898)

FRU = NO




Figure 5-17 SubAssy, GRE Console Rear Bulkhead (5160898ADW, Rev 2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2265406	1.0	Yes	FIBER-OPTIC COUPLING RECEPTACLE KIT, 2.5MM BAYONET, AMP OPTIMATE
2	5160899	1.0	No	Bulkhead Panel, GRE Console Rear Cable Connections
3	2374957	1.0	No	ADAPTER - BULKHEAD BNC JACK-JACK
4	5122589	6.0	Yes	Adapter, Bulkhead RJ45 CAT5e Shielded, L-com #TDG1026KS-C5E
5	2377229	1.0	Yes	CABLE - COAX X 610 LG
6	2318713-2	1.0	No	CABLE-FIBER OPTIC ST-SC X 3FT
8	5127222-3	3.0	Yes	Cable, Ethernet x 3ft Lg - CAT5e, 4 twisted pair, 24awg solid copper, Shielded - aluminized polyester and tinned copper, gray PVC jacket, snagless PVC boot, Plugs - UL94V-0 rated polycarbonate with metallic strain relief, EIA-TIA568, UL listed, CMR rated,
9	5128409	1.0	Yes	Assy, Gbe cable harness for VCT console
5	2319043	4.0	Yes	CASTER, 2 1/2~ SWIVEL

Table 5-18 FWS GOC4 Chassis (5164805 BOM, Rev 3)

2.1.3 GOC4 OUTLET BOX ASSY (5162696)

FRU = NO




Figure 5-18 GOC4 Outlet Box Assy (5162696ADW refer to 2174159-3ADW, Sheet 1, Rev 3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2174159-5	1.0	Yes	ASSY, OUTLET BOX FOR GOC4
2	5162697	1.0	No	LEAD DATA WIRE FOR GOC4

Table 5-19 GOC4 Outlet Box Assy (5162696 BOM, Rev 2)

2.2 CONSOLE INTERCOM PRESCRIBED-TILT HUB (2404320)

FRU = YES

Document Information for: 2404320ADW

Type	Drawing Print
Name	2404320ADW
Revision	11
State	Release
Description	CONSOLE INTERCOM PRESCRIBED-TILT HUB
Originator	212025876_mary_j_kinsey

File List
1. 2404320ADW_s1_r11.pdf

Approval Information		
Person	Action	Date and Time
212025876_mary_j_kinsey	Approved	3/3/2005 11:15:48 AM
212025876_mary_j_kinsey	Approved	3/3/2005 11:29:43 AM

This page is generated automatically by the GEMS ePDM System
Printed documents are for Reference Only and may be out-of-date.
Check the database to ensure you have the correct revision.

State: RELEASE - Document is released and under formal Change Control. Changes are subject to the ECR/ECO Process.
See the GEMS ePDM System to determine the status of this document.

Approved Document - 2404320ADW_r11.pdf Page 1 of 2

Figure 5-19 Console Intercom Prescribed - Tilt Hub (2404320ADW, Rev 11)

MARK	NUMBER	QTY	FRU	DESCRIPTION
3	2374694	2.000	Yes	8A TIME DELAY FUSES
10	46-312577P14	22.000	Yes	SCREW HEAD CAP SCREW STEEL ZINC PLATED (FULL THREAD) PROPERTY CLASS 8.8 ***THIS PART IS REPLACED BY 46-328417P3***
11	46-328430P2	16.000	Yes	WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL
13	46-312346P1	2.000	Yes	M3 X 0.5 ZINC PLATED METRIC HEX NUT WITH ATTACHED TOOTH LOCK WASHER
15	46-214941P2	6.000	Yes	FLAT WASHER, NYLON, M4
19	46-221417P1	3.000	Yes	FEMALE SCREWLOCKS, TWO PER KIT FOR HD CONNECTORS 4-40 THREAD 0.5~ OVERALL LENGTH, YELLOW CHROMATE 0.312~ BOTTOM, 0.062~ PANEL THK. MAX (KIT INCLUDES NUTS AND WASHERS
1	2370567	1.000	N	POWER ENTRY MODULE 2 POLE FUSING
2	2374174	1.000	N	POWER ENTRY FUSE DRAWER
4	2399750	1.000	N	GRE INTERCOM-TILT
5	5114920-2	1.000	N	CHASSIS BASE, ICOM
6	5114921	1.000	N	CT INTERCOM BOX CHASSIS COVER
7	5115027	1.000	N	POWER SUPPLY TRIPLE OUTPUT 90-264V 47-63 HZ AC IN 5V 7A +12V 3A -12V 2.5A OUT
8	5115179	1.000	N	CABLE HARNESS POWER ENTRY TO POWER SUPPLY
9	5115180	1.000	N	CABLE HARNESS POWER SUPPLY TO PWA
12	2276462	6.000	N	METRIC PLAIN FLAT WASHER 3.2MM ID X 7MM OD
14	46-312577P4	2.000	N	SCREW HEAD CAP SCREW STEEL ZINC PLATED (FULL THREAD) PROPERTY CLASS 8.8
16	2209282-3	5.000	N	HEX STANDOFF MALE/FEMALE M4 X 25MM
17	5128595	1.000	N	M4 X 25MM MALE FEMALE NYLON HEX STAND OFF
18	46-328417P73	4.000	N	SCREW HEXAGON SOCKET 3 MM 6 MM

Table 5-20 Outlet Box Assembly (2404320 BOM, Rev 12)

2.3 SCSI TOWER ASSY (5116391-2)

FRU = YES

Document Information for: 5116391ADW

Type	Drawing Print
Name	5116391ADW
Revision	1
State	Release
Description	SCSI TOWER ASSY
Originator	305009506_jun_zhou

File List
1. 5116391ADW_s5116391_r1.pdf

Approval Information		
Person	Action	Date and Time
305006490_xiaoyan_zhang	Approved	6/7/2004 4:01:02 AM
305006621_yao_chen	Approved	6/6/2004 9:08:47 PM

This page is generated automatically by the GEMS ePDM System
Printed documents are for Reference Only and may be out-of-date.
Check the database to ensure you have the correct revision.

State: RELEASE - Document is released and under formal Change Control. Changes are subject to the ECRECO Process.
See the GEMS ePDM System to determine the status of this document.

Approved Document - 5116391ADW_r1.pdf Page 1 of 2

Figure 5-20 SCSI Tower Assy (5116391-2ADW refer to 5116391ADW, Rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
5	46-328421P501	1.0	Yes	SCREW PAN HEAD 3 MM 5 MM
1	2398863	1.0	N	DVD-RAM MULTI-DRIVE 4.7GB
2	2404031	1.0	No	SCSI TO IDE ADAPTER
4	2401055	1.0	No	SCSI TOWER

Table 5-21 SCSI Tower Assy (5116391-2 BOM, Rev 2)

2.4 COVER ASSY - CONSOLE FRONT (5128181)

FRU = YES




Figure 5-21 Cover Assy - Console Front (5128181ADW Sheet 1, Rev 1)



MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2369619	1.000	Yes	COVER, CONSOLE FRONT - GRE
2	5124020	1.000	Yes	Foam, Sound Absorbing - 54in Wd x 1in Thick
3	46-136246P18	3.210	Yes	1~ WIDE, WHT, ADH BACK, VELCRO HOOK 75 FT PER RL
4	5128180	1.000	Yes	Filter, 25 PPI Quadrafoam UL94 HF-1

Table 5-22 Cover Assy - Console Front (5128181 BOM, Rev 1)

Section 3.0 — GRE CONSOLE, 4/8 SLICE GDAS FOR BS WITH XW8200, ETHERNET SWITCH AND JARRELL (5377708-520)

FRU = NO




Figure 5-23 GRE Console, 4/8 Slice GDAS for BS with XW8200, Ethernet Switch and Jarrell (5377708-520ADW, Sheet 1, Rev 1)




Figure 5-24 GRE Console, 4/8 Slice GDAS for BS with XW8200, Ethernet Switch and Jarrell (5377708-520ADW, Sheet 2, Rev 1)




Figure 5-25 GRE Console, 4/8 Slice GDAS for BS with XW8200, Ethernet Switch and Jarrell (5377708-520ADW, Sheet 3, Rev 1)




Figure 5-26 GRE Console, 4/8 Slice GDAS for BS with XW8200, Ethernet Switch and Jarrell (5377708-520ADW, Sheet 4, Rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION
2	5117866-48	1.0	Yes	GEHC PET,CT 4,8,16 SLICE INTEGRATED BUY LEVEL HP8200 WITH FX1500 GRAPHICS CARD HOST COMPUTER
3	5147442	1.0	Yes	COMPUTER - DATA ACQUISITION RECON CONTROL, 2ND GENERATION
5	2404320	1.0	Yes	CONSOLE INTERCOM PRESCRIBED-TILT HUB
6	5116391-2	1.0	Yes	SCSI TOWER ASSY
7	5115334	1.0	Yes	SWITCH, ETHERNET 8-PORT GBE
8	2368697-2	1.0	Yes	LINUX DDR WORKSTATION 3BUTTON PS2 MOUSE
11	5128181	1.0	Yes	COVER ASSY - CONSOLE FRONT
13	2304737-2	2.0	Yes	CABLE, DVI VIDEO CABLE - OCTANE II
17	2324697-2	1.0	Yes	SCSI CABLE, VHD68 TO CENT 50, LINUX
19	5127222-2	2.0	Yes	CABLE, ETHERNET, 2FT LG - CAT5E, 4 TWISTED PAIR, 24AWG SOLID COPPER, SHIELDED - ALUMINIZED POLYESTER AND TINNED COPPER, GRAY PVC JACKET, SNAGLESS PVC BOOT, PLUGS - UL94V-0 RATED POLYCARBONATE WITH METALLIC STRAIN RELIEF, EIA-TIA568, UL LISTED, CMR RATED,
25	46-136246P14	60.0	Yes	2~ WIDE, WHT, ADH BACK, VELCRO LOOP 25 YD
26	46-136246P21	60.0	Yes	2~ WIDE X 25 YD WHITE VELCRO HOOKS
27	46-208781P2	4.0	Yes	WRITE-ON LABEL FOR MARKING CABLES WHITE ADHESIVE BACKED 0.75~ X 1.90~. UNMARKED, 2,500 PC FAN-FOLD SHEETS (FULL SIZE = 1.9~ X 3.15~.)
29	2117390-2	1.0	Yes	CE COMPLIANCE AND WEEE LABELS
30	2132620	1.0	Yes	EMC LABEL CLASS A GROUP 1
36	2337414	2.0	Yes	DVD RAM MEDIA, 9.4 GB
39	46-208758P3	6.0	Yes	7.31~ X .184~ SELF-LOCKING CABLE TIE 1/16~ TO 1-3/4~ WIRE BUNDLE RANGE. 50LB TENSILE STRENGTH, NYLON MATR~L.
42	2103585-2	4.0	Yes	SCREW PAN HEAD 4 MM 10 MM SPRING WASHER, ZINC PLATED, STEEL
43	46-208758P2	1.0	Yes	5.5~ X .140~ SELF-LOCKING CABLE TIE. 1/16~ TO 1-1/4~ WIRE BUNDLE RANGE. 30LB TENSILE STRENGTH, NYLON MATR~L.
44	46-208747P1	1.0	Yes	ADHESIVE MOUNT - CABLE TIE PANDUIT CORP, ABMM-A(STYLE A)

Table 5-23 GRE Console, 4/8 Slice GDAS for BS with XW8200, Ethernet Switch and Jarrell (5377708-520 BOM, Rev 4)

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	5165396	1.0	No	FWS GOC4 CHASSIS,MAKE PARTS.
14	2366442-27	2.0	No	BRACKET, RACK MOUNT HOST COMPUTER
23	5183979	1.0	No	HIVALUE 4,8 JARRELL LCGOC RATING PLATE
28	5133331	1.0	No	ETL LABEL FOR OC
33	2240540	1.0	No	NETSCAPE LICENSE
41	46-312358P51	14.0	No	M5 X 0.8 X 12MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW

Table 5-23 GRE Console, 4/8 Slice GDAS for BS with XW8200, Ethernet Switch and Jarrell (5377708-520 BOM, Rev 4)

3.1 LCGOC Chassis (5165396)

Refer to {*For further detail, see [Section 2.1 FWS GOC4 CHASSIS ASSEMBLY \(5165396\), Page408.](#)*}

3.2 Console Intercom Prescribed-Tilt Hub (2404320)

Refer to {For further detail, see [Section 2.2 CONSOLE INTERCOM PRESCRIBED-TILT HUB \(2404320\), Page418.](#)}

3.3 SCSI Tower ASSY (5116391-2)

Refer to {*For further detail, see [Section 2.3 SCSI TOWER ASSY \(5116391-2\), Page420.](#)*}

3.4 Cover Assy - Console Front (5128181)

Refer to {*For further detail, see [Section 2.4 COVER ASSY - CONSOLE FRONT \(5128181\), Page422.](#)*}

Section 4.0 — Sigle Piece Desktop Kit for Freedom Workspace (5166834)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION
2	2337534	1.0	Yes	SIDE COVERS GOC
1	2385783	1.0	No	TABLE TOP ASSEMBLY, KEYBOARD
3	2318681-2	1.0	No	LEFT ANTI-TIP LEG, CONSOLE CHASSIS
4	2318681-3	1.0	No	RIGHT ANTI-TIP LEG, CONSOLE CHASSIS
5	2109867-28	6.0	No	SCREW HEX SOCKET CAP M10X60/32 A2-70
6	2325236	2.0	No	LEVELER, CONSOLE ANTI-TIP LEGS

Table 5-24 Sigle Piece Desktop Kit for Freedom Workspace (5166834 BOM, Rev 1)

4.1 Side Covers GOC (2337534)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	2315656	1.0	Yes	RIGHT SIDE COVER, GLOBAL CONSOLE
	2315656-2	1.0	Yes	LEFT SIDE COVER, GLOBAL CONSOLE

Table 5-25 Side Covers GOC (2337534BOM, Rev 0)

Section 5.0 — Freedom Workspace Side by Side Kit (5166569)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION
4	5160577	1.0	Yes	OC LONG CABLE KIT
1	2337534-3	1.0	No	FREEDOM WORKSPACE SIDE COVER
2	5160735	1.0	No	FWS TABLE AND ARM
5	5166671	1.0	No	FWS CONSOLE CHASSIS MONITOR TABLETOP

Table 5-26 Freedom Workspace Side by Side Kit (5166569 BOM, Rev 3)

5.1 Freedom Workspace Side Cover (2337534-3)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	2315656-5	1.0	Yes	FREEDOM WORKSPACE RIGHT SIDE COVERS FOR CONSOLE CHASSIS
2	2315656-6	1.0	Yes	FREEDOM WORKSPACE LEFT SIDE COVERS FOR CONSOLE CHASSIS

Table 5-27 Freedom Workspace Side Cover (2337534-3 BOM, Rev 1)

5.2 FWS Table and ARM (5160735)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	5160735-2	1.0	Yes	FREEDOM WORKSPACE TABLE FURNITURE
2	5160735-3	1.0	Yes	FREEDOM WORKSPACE MONITOR MOUNTING ARMS

Table 5-28 FWS Table and ARM (5160735 BOM, Rev 3)

Section 6.0 — Associated Items

- For further detail, see [Section 6.1.1 LCD Monitor - 19 Inch \(B7858LC\)](#), Page446.
- For further detail, see [Section 6.2 Keyboard w/SCIM, Collectors](#), Page448.
- For further detail, see [Section 6.3 Console Slide Service Kit \(2379135\)](#), Page452.
- For further detail, see [Section 6.4 Power Switch & Cables \(2379134\)](#), Page453.
- For further detail, see [Section 6.5 Modem Kit \(B7700MG\)](#), Page453.

6.1 Monitors

6.1.1 LCD Monitor - 19 Inch (B7858LC)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	5117545	2.000	Yes	Monitor, 19 inch color LCD, silver, NEC LCD1980SXi

Table 5-29 LCD Monitor & Cable (2309592 BOM, rev 1)

6.2 Keyboard w/SCIM, Collectors



Figure 5-27 Keyboard & SCIM Assembly (Rev 11/15/00)

6.2.1 Keyboard w/SCIM Collector - English (2275758)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2275756	1.000	Yes	KEYBOARD, NORTH AMERICAN
3	2275752	1.000	Yes	SCIM COLLECTOR - ENGLISH {For further detail, see <Default ? Font>?\$paratext>, ? on page 449.}
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757	1.000	N	FUNCTION KEY TEMPLATE - E

Table 5-30 Keyboard w/SCIM Collector - English (2275758 BOM, rev 1)

6.2.1.1 SCIM Collector – English (2275752)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2262158	1.000	N	SCAN CONTROL INTERCOM MOD
2	2262158-11	1.000	N	OVERLAY, SCIM CT, TILT -
3	2262158-21	1.000	N	OVERLAY, SCIM CT, NO TILT

Table 5-31 SCIM Collector - English (2275752 BOM, rev 0)

6.2.2 French (2275758-2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2275756-2	1.000	Yes	KEYBOARD, EUROPEAN FRENCH
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757-2	1.000	N	FUNCTION KEY TEMPLATE - F
3	2275752-2	1.000	N	SCIM COLLECTOR - FRENCH {For further detail, see <Default ? Font>?\$paratext>, ? on page 450.}

Table 5-32 Keyboard w/SCIM Collector - French (2275758-2 BOM, rev 1)

6.2.2.1 SCIM Collector - French (2275752-2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2262158	1.000	N	SCAN CONTROL INTERCOM MOD
2	2262158-12	1.000	N	OVERLAY, SCIM CT, TILT -
3	2262158-22	1.000	N	OVERLAY, SCIM CT, NO TILT

Table 5-33 SCIM Collector - French (2275752-2 BOM, rev 0)

6.2.3 German (2275758-3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2275756-3	1.000	Yes	KEYBOARD, GERMAN PC PS/2
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757-3	1.000	N	FUNCTION KEY TEMPLATE - G
3	2275752-3	1.000	N	SCIM COLLECTOR - GERMAN {For further detail, see <Default ? Font>?\$paratext>, ? on page 451.}

Table 5-34 Keyboard w/SCIM Collector - German (2275758-3 BOM, rev 1)

6.2.3.1 SCIM Collector - German (2275752-3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2262158	1.000	N	SCAN CONTROL INTERCOM MOD
2	2262158-13	1.000	N	OVERLAY, SCIM CT, TILT -
3	2262158-23	1.000	N	OVERLAY, SCIM CT, NO TILT

Table 5-35 SCIM Collector - German (2275752-3 BOM, rev 0)

6.2.4 Swedish (2275758-4)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2275756-4	1.000	Yes	KEYBOARD, SWEDISH PC PS/2
3	2275752	1.000	Yes	SCIM COLLECTOR - ENGLISH {For further detail, see <Default ? Font>?\$paratext>, ? on page 449.}
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757	1.000	N	FUNCTION KEY TEMPLATE - E

Table 5-36 Keyboard w/SCIM Collector - Swedish (2275758-4 BOM, rev 1)

6.2.5 International (2275758-5)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2275756	1.000	Yes	KEYBOARD, NORTH AMERICAN
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2283830	1.000	N	FUNCTION KEY TEMPLATE COL
3	2275752-4	1.000	N	SCIM COLLECTOR - INTERNATIONAL {For further detail, see <Default ? Font>?\$paratext>, ? on page 452.}

Table 5-37 Keyboard w/SCIM Collector - International (2275758-5 BOM, rev 1)

6.2.5.1 SCIM Collector - International (2275752-4)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2262158	1.000	N	SCAN CONTROL INTERCOM MOD
2	2262158-10	1.000	N	OVERLAY, SCIM CT, TILT -
3	2262158-14	1.000	N	OVERLAY, SCIM CT, TILT -
4	2262158-20	1.000	N	OVERLAY, SCIM CT, NO TILT
5	2262158-24	1.000	N	OVERLAY, SCIM CT, NO TILT

Table 5-38 SCIM Collector - International (2275752-4 BOM, rev 0)

6.3 Console Slide Service Kit (2379135)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	2317595	1.000	N	SLIDE RAIL FOR DRAWER ASS
4	2317501	1.000	N	BRACKET, DRAWER BALL LATC
5	2317502	1.000	N	LATCH, DRAWER ASSEMBLY
6	2366442-7	4.000	N	BRACKET, RACK MOUNT

Table 5-39 Console Slide Service Kit (2379135 BOM, rev 0)

6.4 Power Switch & Cables (2379134)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2150833	1.000	N	BREAKER SWITCH
2	2108900	1.000	N	CONNECTOR 8 CONTACTS 150 VOLTS
5	2318308	1.000	N	ETHERNET CABLE, CHASSIS SVC PANEL
7	2341434-2	1.000	N	CABLE, KEYBOARD - OCTANE 2
8	2341434	1.000	N	CABLE, MOUSE - OCTANE 2

Table 5-40 Power Switch & Cables (2379134 BOM, rev 0)

6.5 Modem Kit (B7700MG)

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	2382276	1.000	N	GRE MODEM KIT

Table 5-41 GOC3 Modem Kit (B7700MG BOM, rev 0)

Section 7.0 — Options

- For further detail, see [Section 7.1 Video Upgrade For Linux Console \(B7868MC\)](#), Page454.
- For further detail, see [Section 7.2 LCD & SUSP \(Boom In Room\) LNX \(B7710LM\)](#), Page455.
- For further detail, see [Section 7.3 HS/LS Remote Color \(LNX\) Monitor \(B7530LC\)](#), Page456.

7.1 Video Upgrade For Linux Console (B7868MC)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2378278	1.000	N	VIDEO UPGRADE FOR LINUX C {For further detail, see <Default ? Font> Section 7.1.1<Default ? Font> .}

Table 5-42 Video Upgrade For Linux Console (B7868MC BOM, rev 0)

7.1.1 Video Upgrade (2378278)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	2384322	1.000	Yes	5 BNC MALE TO HD15 MALE - 75 FT PLENUM RATED
9	2302528	1.000	Yes	VIDIO CABLE FOR LCD SUSPE
13	2376772	1.000	Yes	1-IN 4-OUT RGBHV DISTRIBUTION AMP
27	2384323	10.000	Yes	75 OHM TERMINATOR
28	2384324	1.000	Yes	DVI TO HD15 CONVERTER
1	46-302200P10	1.000	N	1.72 X 1.10, MODEL#, SN, DATE, DESC
14	2376770-2	1.000	N	5 BNC MALE TO HD15 MALE (6 FT)
16	2376771	1.000	N	5 BNC MALE TO HD15 FEMALE (6 FT)
18	2200583-100	1.000	N	CT/I (SONATA) UPGRADE MANUAL
19	2211224-100	1.000	N	LIGHTSPEED QX/I OPTIONS INSTALLATION
25	2212092	1.000	N	CEILING PLATE ASSEMBLY FOR LCD

Table 5-43 Video Upgrade (2378278 BOM, rev 3)

7.2 LCD & SUSP (Boom In Room) LNX (B7710LM)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2376773	1.000	N	BOOM-IN-A-ROOM FOR LINUX {For further detail, see <Default ? Font> Section 7.2.1 <Default ? Font>.}

Table 5-44 LCD & SUSP (Boom In Room) LNX (B7710LM BOM, rev 0)

7.2.1 Boom-in-a-Room for Linux (2376773)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	2384322	1.000	Yes	5 BNC MALE TO HD15 MALE - 75 FT PLENUM RATED
3	2214293	1.000	Yes	CEILING COLUMN WITH FLANGE, 58CM
4	2214292	1.000	Yes	80-80 CM EXTENSION/SPRING ARM 10,3 KG (22.7 LBS)
5	2349792	1.000	Yes	18.1 IN LCD COLOR MONITOR
6	2211031	1.000	Yes	LCD BRACKET - SHORT
7	2212717	1.000	Yes	CABLE COVER FOR LCD SUSPENSION 3.6 METERS LONG
9	2302528	1.000	Yes	VIDIO CABLE FOR LCD SUSPE
11	2213213	1.000	Yes	GROUNDING CABLE FOR LCD SUSPENSION
12	2213218	1.000	Yes	POWER CABLE FOR LCD SUSPENSION
13	2376772	1.000	Yes	1-IN 4-OUT RGBHV DISTRIBUTION AMP
27	2384323	10.000	Yes	75 OHM TERMINATOR
28	2384324	1.000	Yes	DVI TO HD15 CONVERTER
1	46-302200P10	1.000	N	1.72 X 1.10, MODEL#, SN, DATE, DESC
8	2213216	1.000	N	GROUNDING CABLE FOR LCD-CONSOLE TO LCD
10	2213219	1.000	N	POWER CABLE FOR LCD-CONSOLE TO LCD
14	2376770-2	1.000	N	5 BNC MALE TO HD15 MALE (6 FT)
16	2376771	1.000	N	5 BNC MALE TO HD15 FEMALE (6 FT)
17	2223870	1.000	N	SONATA HARDWARE COLLECTOR KIT
18	2200583-100	1.000	N	CT/I (SONATA) UPGRADE MANUAL
19	2211224-100	1.000	N	LIGHTSPEED QX/I OPTIONS INSTALLATION

Table 5-45 Boom-in-a-Room for Linux (2376773 BOM, rev 2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
25	2212092	1.000	N	CEILING PLATE ASSEMBLY FOR LCD
26	2300059	1.000	N	HANDLE FOR LCD DISPLAY MONITOR

Table 5-45 Boom-in-a-Room for Linux (2376773 BOM, rev 2) (Continued)

7.3 HS/LS Remote Color (LNX) Monitor (B7530LC)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2377146	1.000	N	REMOTE MONITOR KIT FOR LI {For further detail, see <Default ? Font> Section 7.3.1 <Default ? Font>.}

Table 5-46 HS/LS Remote Color (LNX) Monitor (B7530LC BOM, rev 0)

7.3.1 Remote Monitor Kit (2377146)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2333243	1.000	Yes	21 IN COLOR CRT MONITOR
1A	2254212	ALT	Yes	SONY GDM-5402T9 21" COLOR

Table 5-47 Remote Monitor Kit (2377146 BOM, rev 0)