
Chapter 6

Ship Collector

Section 1.0 — BrightSpeed Elite Select SHIP COLLECTOR (5197961-2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2381045-2	1.0	Yes	HP60 FIELD BALANCE WEIGHT KIT {For further detail, see Section 1.2, beginning on page 459.}
2	2385412-2	1.0	Yes	FUSE KIT BS 90KVA {For further detail, see Section 1.3, beginning on page 460.}
4	5248959-2	1.0	Yes	BSD BASE COVER COLLECTOR, PF {For further detail, see Section 1.1, beginning on page 458.}
3	5213061-2	1.0	No	BS Lite Affiliate Buy Collector {For further detail, see Section 1.4, beginning on page 461.}

Table 6-1 BrightSpeed Elite Select Ship Collector (5197961-2 BOM, Rev 1)

1.1 BSD BASE COVER COLLECTOR (5248959-2)

FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	46-328417P4	15.0	Yes	SCREW HEXAGON SOCKET 4 MM 12 MM ZINC PLATED, STEEL
2	46-328430P2	21.0	Yes	WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL
3	46-328432P2	25.0	Yes	WASHER LOCK - SPRING 4.1 MM 7.6 MM ZINC PLATED, SPRING STEEL
4	46-328425P1	4.0	Yes	NUT HEXAGON 4 MM
5	5248961	1.0	Yes	BSD BASE FRONT ASSY, PF
6	5248965-2	1.0	Yes	BSD BASE LEFT SIDE ASSY, PF
7	5248968-2	1.0	Yes	BSD BASE RIGHT SIDE ASSY, PF
8	5248970	1.0	Yes	BSD BASE LEFT REAR ASSY, PF
9	5248973	1.0	Yes	BSD BASE RIGHT REAR ASSY, PF
16	46-328417P5	2.0	Yes	SCREW HEXAGON SOCKET 4 MM 16 MM ZINC PLATED, STEEL
10	5144160	2.0	No	SL4 BRACKET SIDE BASE COVER
11	2261325	1.0	No	BRACKET - REAR BASE COVER
12	5144162-2	1.0	No	SL4 BRACKET, BASE COVER MOUNTING - LEFT

Table 6-2 BSD BASE COVER COLLECTOR, PF (5248959-2 BOM, Rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
13	5144163-2	1.0	No	SL4 BRACKET, BASE COVER MOUNTING - RIGHT
14	5310493	2.0	No	BULKHEAD - BLANK
15	46-328420P4	4.0	No	SCREW BUTTON HEAD 4 MM 16 MM BLACK OXIDE, STEEL

Table 6-2 BSD BASE COVER COLLECTOR, PF (5248959-2 BOM, Rev 1) (Continued)

1.2 HP60 FIELD BALANCE WEIGHT KIT (2381045-2)

FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2378669	11.0	Yes	P2033DU 78WEIGHT-M
2	2378670	4.0	Yes	P2033DV 78WEIGHT-S
4	2317865-2	2.0	Yes	SPACER,180 DEG,25.4MM
5	2317994	7.0	Yes	WEIGHT,TRIM
8	2342895-2	3.0	Yes	WEIGHT,180 DEG,6MM
9	2342895-3	4.0	Yes	WEIGHT,180 DEG,1.5MM
12	46-328425P4	4.0	Yes	NUT HEXAGON 16 MM
13	46-328417P70	2.0	Yes	SCREW HEXAGON SOCKET 12 MM 100 MM
14	2376683	1.0	Yes	MCMaster CARR PN 6023A21 SOCKET 0.9375 IN DEEP IMPACT 0.500 IN DRV
18	46-328432P5	4.0	Yes	WASHER LOCK - SPRING 16.2 MM 27.4 MM ZINC PLATED, SPRING STEEL
20	2378668	1.0	Yes	P2033DT 78WEIGHT-L
21	46-328432P7	2.0	Yes	WASHER LOCK - SPRING 12.2 MM 21.1 MM
22	2236649-5	2.0	Yes	WASHER, FLAT, ISO7089,M12, ST-HB200
24	2317863	9.0	Yes	SHIM,180 DEG,3.2MM
10	2370053	1.0	No	WEIGHT,TRIAL,NEGATIVE Z
11	2370227	1.0	No	WEIGHT,TRIAL,POSITIVE Z, H-POWER
15	2377399	1.0	No	LEFT SIDE DAS COVER TOOL

Table 6-3 Balance Kit, Field, H-Power (2381045-2 BOM, Rev 2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
16	2377400	1.0	No	RIGHT SIDE DAS COVER TOOL
17	2377401	1.0	No	CENTER DAS COVER TOOL
23	5128067	2.0	No	CAP SCREW U0480AA

Table 6-3 Balance Kit, Field, H-Power (2381045-2 BOM, Rev 2) (Continued)

1.3 FUSE KIT BS 90KVA (2385412-2)

FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2351493	3.0	Yes	100A FUSE
2	2364059	2.0	Yes	GLASS FUSE
3	46-170021P50	2.0	Yes	FUSE 12 AMPS 250 VOLTS BUSSMAN MDA12 DUAL ELEMENT.
4	46-170021P15	2.0	Yes	FUSE 8 AMPS 250 VOLTS (REVIEWED TK, 7/91).
5	2336517-2	2.0	Yes	FUSE 25 AMPS 700 VOLTS 2.5M SECONDS
6	46-170021P52	3.0	Yes	FUSE 3 AMPS 250 VOLTS (REVIEWED TK, 7/91).
7	46-170021P10	2.0	Yes	006.000A 0250V 3AG FAST UL/C (REVIEWED TK, 7/91).
8	46-170021P101	1.0	Yes	FUSE 20 AMPS 700 VOLTS 2~ LONG X 9/16~ DIA FAST BLOW
9	46-170021P106	2.0	Yes	FUSE 8 AMPS 250 VOLTS BUSSMANN ONLY
10	46-170021P31	2.0	Yes	1/2A, 250V SLO-BLO FUSE. TYPE 3AG, 1.25~ X 0.25~ GLASS BODY. (REVIEWED TK, 7/91).
11	2106993-5	2.0	Yes	FUSE 20 AMPS 500 VOLTS .3 SECONDS
12	46-327160P1	2.0	Yes	12.0A, 125VAC, DUAL ELEMENT TIME DLY FNM-12, 0.406 X 1.5 LG, FIBER BODY
13	2379651	1.0	Yes	FUSE - 700 V, 200 A, FAST SEMI
15	2238207-3	1.0	Yes	FUS FF 14X51 2A 700VDC C 50KA ULR
16	99183979	1.0	Yes	SLOW BLOW FUSE HCP 6X32 10A 250V UL
17	U8009FF	3.0	Yes	FUSE 3.2A

Table 6-4 FUSE KIT BS 90KVA (2385412-2 BOM, Rev 2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
18	U8008FF	3.0	Yes	FUSE 2.0A
19	U8006FF	3.0	Yes	FUSE 1.0A

Table 6-4 FUSE KIT BS 90KVA (2385412-2 BOM, Rev 2) (Continued)

1.4 BS LITE AFFILIATE BUY COLLECTOR (5213061-2)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	5197957	1.0	No	BS COMMON SHIP COLLECTOR
2	5197959	1.0	No	BS LITE UNIQUE SHIP COLLECTOR

Table 6-5 BS Lite Affiliate Buy Collector (5213061-2BOM, Rev 1)

1.4.1 BS COMMON SHIP COLLECTOR (5197957)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	46-229452G1	1.0	Yes	CT-9800 (ACCESSORIES)
2	2105873	4.0	Yes	.531 I.D., 1.625 O.D., .25THK GANTRY & TABLE INSTALLATION ANCHOR WASHER, FOR BASEPLATE-LESS HSA SYSTEMS.
4	46-222436G1	4.0	Yes	CT-9800 - CHIN STRAP
7	46-222443P1	1.0	Yes	OPEN CELL POLYESTER FOAM
8	46-222444P1	1.0	Yes	CLOSED CELL LINEAR ETHYLENE
9	46-222445P1	1.0	Yes	CLOSED CELL LINEAR ETHYLENE
10	46-222921P1	6.0	Yes	OBSOLETE FOR MR SITES PER FMI 60288 POLYURETHANE FOAM 1.7#/FT DENSITY S.E. FINE CELL CHARCOAL GRAY
11	46-199388P1	1.0	Yes	HEAD HOLDER POLYURETHANE FOAM
13	5179183	1.0	Yes	DOLLY, PIVOTING, SMALL ROOM, RH
14	5179183-2	1.0	Yes	DOLLY, PIVOTING, SMALL ROOM, LH
15	2282509	1.0	Yes	HOIST-2-PC, ARM, HELIOS 2
16	2282510	1.0	Yes	HOIST-2-PC, LONG TUBE, HELIOS 2
18	2364629	1.0	Yes	LABEL INSTALLATION KIT
19	2336557	1.0	Yes	SEISMIC KIT FOR GLOBAL CONSOLE
20	2373436-7	1.0	Yes	CABLE, FE LAPTOP TO CONSOLE, 23000MM LONG
5	2319891-3	1.0	No	INSTALLATION SUPPORT KIT FOR H-POWER
6	46-237412P1	1.0	No	STRAP HEAD STRAP 3 STRAPS PER PACKAGE
12	5150476	2.0	No	NON-TILTING DOLLY FOR HIVALUE GANTRY COVER
17	5152094	4.0	No	PLUG,SIDE COVER SL4
21	2318713-2	1.0	No	CABLE-FIBER OPTIC ST-SC X 3FT

Table 6-6 BS Common Ship Collector (5197957 BOM, Rev 4)

1.4.2 BS Lite UNIQUE SHIP COLLECTOR (5197959)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	5193991	1.0	Yes	FLOOR TEMPLATE-HIVALUE LITE
2	2152504	1.0	Yes	P9150SQ BODY STRAP C
3	2152503	1.0	Yes	P9150SP BODY STRAP B
4	2152502	1.0	Yes	P9150SN BODY STRAP A
5	2115984-3	1.0	Yes	P9230JR PAD EXTENDER
6	P9155BZ	4.0	Yes	WASHER ANCHOR
7	2214306	6.0	Yes	ANCHORS STUD 120 MM 12 MM
8	5180876	1.0	Yes	TARGET,TTABLE-ALIGNMENT
9	2331933	1.0	Yes	METALLESS COMPATIBLE PHANTOM HOLDER
10	2106205	1.0	Yes	500 I.D. 1.00 O.D. 2.125 LONG, ANCHOR HOLE CENTERING ON BASEPLATE-LESS HSA SYSTEMS
11	2107863	1.0	Yes	SPECIAL SOCKET FOR LEVELING ADJUSTERS USED ON BASEPLATE- LESS HSA SYSTEMS
12	5214921	6.0	No	ANCHOR HSA M12x220_U0094BF
13	5194598	1.0	No	DRILL BUSHING NP LITE POSITIONING BSD LITE
14	5263871	1.0	No	Torque wrench tip for BSD Lite

Table 6-7 BS Lite Unique Ship Collector (5197959 BOM, Rev 5)

Section 2.0 — BrightSpeed Excel/Edge Select SHIP COLLECTOR (5198104-2)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	5248959-2	1.0	Yes	BSD BASE COVER COLLECTOR, PF {For further detail, see Section 1.1, beginning on page 458.}
4	2381045-2	1.0	Yes	HP60 FIELD BALANCE WEIGHT KIT {For further detail, see Section 1.2, beginning on page 459.}
5	2385412-2	1.0	Yes	FUSE KIT BS 90KVA {For further detail, see Section 1.3, beginning on page 460.}
1	5198105	1.0	No	HIVALUE LITE AND POWER 4/8 COMMON ITEMS {For further detail, see Section 2.1, beginning on page 465.}
2	5197959	1.0	No	BS LITE UNIQUE SHIP COLLECTOR {For further detail, see Section 1.4.2, beginning on page 464.}

Table 6-8 BrightSpeed Lite 4/8 Ship Collector (5198104-2 BOM, Rev 1)

2.1 BS LITE AND POWER 4/8 COMMON ITEMS (5198105)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	46-229452G1	1.0	Yes	CT-9800 (ACCESSORIES)
2	2105873	4.0	Yes	.531 I.D., 1.625 O.D., .25THK GANTRY & TABLE INSTALLATION ANCHOR WASHER, FOR BASEPLATE-LESS HSA SYSTEMS.
4	46-222436G1	4.0	Yes	CT-9800 - CHIN STRAP
7	46-222443P1	1.0	Yes	OPEN CELL POLYESTER FOAM
8	46-222444P1	1.0	Yes	CLOSED CELL LINEAR ETHYLENE
9	46-222445P1	1.0	Yes	CLOSED CELL LINEAR ETHYLENE
10	46-222921P1	6.0	Yes	OBSOLETE FOR MR SITES PER FMI 60288 POLYURETHANE FOAM 1.7#/FT DENSITY S.E. FINE CELL CHARCOAL GRAY
11	46-199388P1	1.0	Yes	HEAD HOLDER POLYURETHANE FOAM
12	5150476	2.0	Yes	NON-TILTING DOLLY FOR HIVALUE GANTRY COVER
13	5179183	1.0	Yes	DOLLY, PIVOTING, SMALL ROOM, RH

Table 6-9 BS Lite and Power 4/8 Common Items (5198105 BOM, Rev 3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
14	5179183-2	1.0	Yes	DOLLY, PIVOTING, SMALL ROOM, LH
15	2282509	1.0	Yes	HOIST-2-PC, ARM, HELIOS 2
16	2282510	1.0	Yes	HOIST-2-PC, LONG TUBE, HELIOS 2
18	2364629	1.0	Yes	LABEL INSTALLATION KIT
19	2336557	1.0	Yes	SEISMIC KIT FOR GLOBAL CONSOLE
20	2373436-7	1.0	Yes	CABLE, FE LAPTOP TO CONSOLE, 23000MM LONG
5	2319891-6	1.0	No	INSTALLATION SUPPORT KIT FOR HIVALUE 4 AND 8 SLICE SYSTEM
6	46-237412P1	1.0	No	STRAP HEAD STRAP 3 STRAPS PER PACKAGE
17	5152094	4.0	No	PLUG,SIDE COVER SL4
21	2318713-2	1.0	No	CABLE-FIBER OPTIC ST-SC X 3FT

Table 6-9 BS Lite and Power 4/8 Common Items (5198105 BOM, Rev 3)

Section 3.0 — CABLE COLLECTORS

3.1 Short Cable Collector (2281840-5)

FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2371450-2	1.000	1	RUN#55A-PDU TO RACEWAY GND
2	2371450-4	1.000	1	RUN#56A-CONSOLE TO RACEWAY
3	2373436-3	1.000	1	RUN #102A - GANTRY TO CON
4	2117848-7	1.000	1	20000 +/- 150 FIBER OPTIC
5	2371133-2	1.000	1	CONSOLE CABLES - SHORT KI
6	2371134-2	1.000	1	GANTRY CABLES - SHORT KIT

Table 6-10 Short Cable Collector (2281840-5 BOM, rev 0)

3.2 Long Cable Collector (2281840-4)

FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2371450	1.000	1	RUN#55-PDU TO RACEWAY GND
2	2371450-3	1.000	1	RUN#56-CONSOLE TO RACEWAY
3	2373436-2	1.000	1	RUN #102 - GANTRY TO CONS
4	2117848-2	1.000	1	24384 +/- 150 FIBER OPTIC
5	2371133	1.000	1	CONSOLE CABLES - LONG KIT
6	2371134	1.000	1	GANTRY CABLES - LONG KIT

Table 6-11 Long Cable Collector (2281840-4 BOM, rev 0)

