

Chapter 5

Ship Collector

RT4 MDAS Ship Collector: 2394003

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2381848	1.000	Yes	COLLECTOR-BASE COVER-SG {For further detail, see Section 1.0, beginning on page 311. }
10	2218570-2	1.000	Yes	H-POWER FUSE KIT {For further detail, see Section 3.0, beginning on page 313. }
20	46-222443P1	1.000	Yes	OPEN CELL POLYESTER FOAM
23	46-222921P1	6.000	Yes	OBSOLETE FOR MR SITES PER
24	46-199388P1	1.000	Yes	HEAD HOLDER POLYURETHANE
25	2282011	2.000	Yes	NON-TILT DOLLY REAR COVER
26	5101127	1.000	Yes	DOLLY PIVOTING , RIGHT HA
30	5101127-2	1.000	Yes	DOLLY, PIVOTING, LEFT HAN
31	2282509	1.000	Yes	HOIST-2-PC, ARM, HELIOS 2
32	2282510	1.000	Yes	HOIST-2-PC, LONG TUBE, HE
33	2277548	4.000	Yes	PLUG, SIDE COVER H2
36	2376682-2	1.000	Yes	BALANCE KIT, FIELD, HPOWER {For further detail, see Section 5.0, beginning on page 316. }
40	2364629	1.000	Yes	LABEL INSTALLATION KIT
41	2281840-5	1.000	Yes	H-POWER SHORT CABLE COLLECTOR {For further detail, see Section 6.0, beginning on page 317. }
-	2281840-4	-	Yes	H-POWER LONG CABLE COLLECTOR {For further detail, see Section 7.0, beginning on page 317. }
42	2336557	1.000	Yes	SEISMIC KIT FOR GLOBAL CONSOLE
48	2386850	2.000	Yes	TRIAL WEIGHT, SG GANTRY,
1	46-229452G1	1.000	No	CT-9800 (ACCESSORIES)
2	46-297482G1	1.000	No	IV BOARD ACC FOR ZEUS TAB

Table 5-1 LightSpeed RT STARGATE SHIP COLLECTOR Ship Collector (2394003 BOM, rev 2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
4	46-297592P1	1.000	No	BAG HOLDER
6	2105873	10.000	No	.531 I.D., 1.625 O.D., .2
7	2106573	10.000	No	HILTI ANCHOR STUD 1/2 DIA
8	2115984	1.000	No	PAD FOR HSA METALLESS FOO
9	2115994	1.000	No	METALLESS CRADLE PAD & ST
13	46-297960P1	1.000	No	FOAM INSERT
14	2319891-4	1.000	No	INSTALLATION SUPPORT KIT, {For further detail, see Section 4.0, beginning on page 314. }
18	46-222436G1	4.000	No	CT-9800 - CHIN STRAP
19	46-237412P1	1.000	No	STRAP
21	46-222444P1	1.000	No	CLOSED CELL LINEAR ETHYLE
22	46-222445P1	1.000	No	CLOSED CELL LINEAR ETHYLE
29	2389013	1.000	No	FLOOR TEMPLATE-SG
34	2217173-2	1.000	No	STOOL WITH FOOT REST AND
35	2299113	1.000	No	CT SERVICE CABINET
37	2309994	1.000	No	IV POLE ASSY
38	2329064	1.000	No	TABLE ACC ASSY
39	2355100	1.000	No	TRAY, H2 TABLE
46	46-237605P1	1.000	No	WRIST STRAP
47	2393793	2.000	No	BAR, FOOT SWITCH-SG

Table 5-1 LightSpeed RT STARGATE SHIP COLLECTOR Ship Collector (2394003 BOM, rev 2) (Continued)

RT4 GDAS Ship Collector: 5135132

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2381848-3	1.000	Yes	COLLECTOR-BASE COVER-RT {For further detail, see Section 2.0, beginning on page 312. }
10	2218570-2	1.000	Yes	H-POWER FUSE KIT {For further detail, see Section 3.0, beginning on page 313. }
20	46-222443P1	1.000	Yes	OPEN CELL POLYESTER FOAM
23	46-222921P1	6.000	Yes	OBSOLETE FOR MR SITES PER
24	46-199388P1	1.000	Yes	HEAD HOLDER POLYURETHANE
25	5180151	2.000	Yes	REAR COVER DOLLY
26	5179183	1.000	Yes	Dolly, pivoting, small room, RH
27	5179183-2	1.000	Yes	Dolly, pivoting, small room, LH
31	2282509	1.000	Yes	HOIST-2-PC, ARM, HELIOS 2
32	2282510	1.000	Yes	HOIST-2-PC, LONG TUBE, HE
33	2277548	4.000	Yes	PLUG, SIDE COVER H2
36	2376682-2	1.000	Yes	BALANCE KIT, FIELD, HPOWER {For further detail, see Section 5.0, beginning on page 316. }
40	2364629	1.000	Yes	LABEL INSTALLATION KIT
42	2336557	1.000	Yes	SEISMIC KIT FOR GLOBAL CONSOLE
48	2386850	2.000	Yes	TRIAL WEIGHT, SG GANTRY,
49	2373436-7	1.000	Yes	Cable, FE Laptop to Console, 23000mm Long
1	46-229452G1	1.000	No	CT-9800 (ACCESSORIES)
2	46-297482G1	1.000	No	IV BOARD ACC FOR ZEUS TAB
4	46-297592P1	1.000	No	BAG HOLDER
6	2105873	10.000	No	.531 I.D., 1.625 O.D., .2
7	2106573	10.000	No	HILTI ANCHOR STUD 1/2 DIA
8	2115984	1.000	No	PAD FOR HSA METALLESS FOO
9	2115994	1.000	No	METALLESS CRADLE PAD & ST
13	46-297960P1	1.000	No	FOAM INSERT

Table 5-2 STARGATE SHIP COLLECTOR, Without short cable collector (5135132 BOM, rev 9)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
14	2319891-4	1.000	No	INSTALLATION SUPPORT KIT, {For further detail, see Section 4.0, beginning on page 314. }
18	46-222436G1	4.000	No	CT-9800 - CHIN STRAP
19	46-237412P1	1.000	No	STRAP
21	46-222444P1	1.000	No	CLOSED CELL LINEAR ETHYLE
22	46-222445P1	1.000	No	CLOSED CELL LINEAR ETHYLE
29	2389013	1.000	No	FLOOR TEMPLATE-SG
34	2217173-2	1.000	No	STOOL WITH FOOT REST AND
35	5112283	1.000	No	CT SERVICE CABINET
37	2309994-2	1.000	No	IV POLE ASSY, with Multi-Language Labels
38	2329064-2	1.000	No	TABLE TRAY ACC ASSY, with Multi-Language Lables
39	2355100	1.000	No	TRAY, H2 TABLE
46	46-237605P1	1.000	No	WRIST STRAP
47	2393793	2.000	No	BAR, FOOT SWITCH-SG
51	2364311	2.000	No	TAPERED PLASTIC HANDLE - METAL INSERT

Table 5-2 STARGATE SHIP COLLECTOR, Without short cable collector (5135132 BOM, rev 9) (Continued)

Section 1.0 — Collector-Base Cover-SG (2381848)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2383481	1.000	Yes	COVER ASM,BASE,FRONT LEFT
3	2383483	1.000	Yes	COVER ASM,BASE,REAR LEFT-
4	2381839	1.000	Yes	COVER-BASE FRONT-SG
5	2383459	1.000	Yes	COVER ASSY, REAR BASE-SG
6	2383482	1.000	Yes	COVER ASM,BASE,FRONT RIGH
7	2383484	1.000	Yes	COVER ASM,BASE,REAR RIGHT
8	2381846	1.000	Yes	COVER-BASE RIGHT-SG
1	2381847	1.000	No	COVER-BASE LEFT-SG
9	2382421	2.000	No	BRACKET,REAR BASE COVER L
10	2382425	2.000	No	BRACKET,FRONT BASE COVER
11	2382426	1.000	No	BRACKET,FRONT LEFT BASE C
12	2382428	1.000	No	BRACKET,FRONT RIGHT BASE
13	2382439	2.000	No	BRACKET,REAR BASE COVER L
14	2382441	2.000	No	BRACKET,GANTRY OPTIONS MT
15	2288882	2.000	No	BULKHEAD - BLANK
16	46-328417P3	42.000	No	SCREW HEXAGON SOCKET 4 MM
17	46-328432P2	46.000	No	WASHER LOCK - SPRING 4.1
18	46-328430P2	24.000	No	WASHER PLAIN - NORMAL 4.3
19	46-328425P1	4.000	No	NUT HEXAGON 4 MM
21	2394127	2.000	No	BRACKET,BASE COVER MTG,UP
22	2393783	2.000	No	BRKT,BASE CVR LATCH,REAR

Table 5-3 Collector-Base Cover-SG (2381848 BOM, rev 1)

Section 2.0 — Collector-Base Cover-RT (2381848-3)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2383481-3	1.000	Yes	COVER ASM,BASE,FRONT LEFT-RT
3	2383483-2	1.000	Yes	COVER ASM,BASE,REAR LEFT- RT
4	2381839-3	1.000	Yes	COVER-BASE FRONT-RT
5	2383459-2	1.000	Yes	COVER ASSY, REAR BASE-RT
6	2383482-3	1.000	Yes	COVER ASM,BASE,FRONT RIGH-RT
7	2383484-2	1.000	Yes	COVER ASM,BASE,REAR RIGHT-RT
8	2381846-2	1.000	Yes	COVER-BASE RIGHT-RT
1	2381847-2	1.000	No	COVER-BASE LEFT-RT
9	2382421	2.000	No	BRACKET,REAR BASE COVER L
10	2382425-2	2.000	No	BRACKET,FRONT BASE COVER-RT
11	2382426	1.000	No	BRACKET,FRONT LEFT BASE C
12	2382428	1.000	No	BRACKET,FRONT RIGHT BASE
13	2382439	2.000	No	BRACKET,REAR BASE COVER L
14	2382441	2.000	No	BRACKET,GANTRY OPTIONS MT
15	2288882	2.000	No	BULKHEAD - BLANK
16	46-328417P3	42.000	No	SCREW HEXAGON SOCKET 4 MM
17	46-328432P2	46.000	No	WASHER LOCK - SPRING 4.1
18	46-328430P2	24.000	No	WASHER PLAIN - NORMAL 4.3
19	46-328425P1	4.000	No	NUT HEXAGON 4 MM
21	2394127	2.000	No	BRACKET,BASE COVER MTG,UP
22	2393783	2.000	No	BRKT,BASE CVR LATCH,REAR

Table 5-4 Collector-Base Cover-SG (2381848-3 BOM, rev 1)

Section 3.0 — H-Power Fuse Kit (2218570-2)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2351492	2.000	Yes	200A FUSE
2	2364059	2.000	Yes	GLASS FUSE
3	46-170021P50	2.000	Yes	FUSE 12 AMPS 250 VOLTS
4	46-170021P106	2.000	Yes	FUSE 8 AMPS 250 VOLTS BUSSMANN ONLY
5	2336517-2	2.000	Yes	FUSE 25 AMPS 700 VOLTS 2.
6	46-170021P52	3.000	Yes	FUSE 3 AMPS 250 VOLTS
7	46-170021P10	2.000	Yes	006.000A 0250V 3AG
8	46-170021P101	1.000	Yes	FUSE 20 AMPS 700 VOLTS
9	46-170021P11	2.000	Yes	7.0A, 250V SLO-BLO FUSE.
10	46-170021P31	2.000	Yes	1/2A, 250V SLO-BLO FUSE.
11	2106993-5	2.000	Yes	FUSE 20 AMPS 500 VOLTS .3
12	46-327160P1	2.000	Yes	12.0A, 125VAC, DUAL ELEME
13	2379651	1.000	Yes	FUSE - 700 V, 200 A, FAST

Table 5-5 H-Power Fuse Kit (2218570-2 BOM, rev 2)

Section 4.0 — Installation Support Kit, Stargate (2319891-4)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
28	2227089-2	1.000	Yes	DAS/DETECTOR SERVICE KIT
39	2214163	1.000	Yes	CONNECTOR, ANALOG LOOP BA
40	2214163-2	1.000	Yes	CONNECTOR, DIGITAL LOOP B
45	2284023-3	1.000	Yes	SERVICE TOOL COLLECTOR, H
46	46-136246P14	6.000	Yes	2~ WIDE, WHT, ADH BACK, VELCRO LOOP 25 YD
47	2278989	6.000	Yes	LANYARD, NYLON U-LOCK, TYPE 1 B=8~ C=4~
48	46-136246P15	6.000	Yes	3/4~ WIDE, WHT ADH BACK, VELCRO LOOP 75 FT PER RL
49	46-170123P1	1.000	Yes	WIDTH 0.625 THK .007 COLOR-BLACK MATERIAL-VINYL (REV TK, 3/92) - OK.
50	46-136246P21	6.000	Yes	2~ WIDE X 25 YD WHITE VELCRO HOOKS
52	46-136246P17	6.000	Yes	3/4~ WIDE, WHT ADH BACK, HOOK VELCRO 45 YARD
53	46-221561P2	5.000	Yes	VINYL TRIM W/ SEGMENTED METAL CORE. FOR 1/32~ EDGES, 500 FT. (COLOR - WHITE SOF-TONE) TRIM HAS SEGMENTED METAL CORE, .018~ THK COLD ROLLED STEEL, #2 TEMPER. (DWG REV 7/92, JFK - OK).
54	46-208753P3	6.000	Yes	ID 0.375 PAD WIDTH 0.75 LGTH 0.75 PVC MATERIAL. (REVIEWED TK, 7/91).
55	46-170123P2	3.000	Yes	3/4~ WIDTH, .007 BLACK VINYL
58	5111733	2.000	Yes	Label for front cover to explain that GE external laser lights are not in use
8	46-170923P1	1.000	No	2-PART, EPOXY PATCH
10	46-221561P3	1.000	No	CLAMP OR CLIP
11	46-208758P10	25.000	No	8.00" X .094" SELF-LOCKIN
12	46-208758P3	50.000	No	7.31" X .184" SELF-LOCKIN
13	46-208758P5	50.000	No	SELF LOCKING CABLE TIE
14	46-208759P5	15.000	No	6.00" X .140" SELF-LOCKIN
15	46-208759P6	20.000	No	7.81" X .184" SELF-LOCKIN
16	46-208759P4	20.000	No	13.93"X .301" SELF-LOCKIN
17	46-208747P1	25.000	No	CABLE TIE,ADH.BACKED
19	46-230118P1	2.000	No	PAPER WYP

Table 5-6 Installation Support Kit, Stargate (2319891-4 BOM, rev 4)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
20	2253867	1.000	No	SCOTCH BRIGHT SCOUR PAD M
21	46-170127P5	1.000	No	1.0" MASKING TYPE TAPE, 6
23	2207303-7	25.000	No	SAFE SKIN NITRILE GLOVE,
24	46-170686P3	1.000	No	LOCTITE 242, 50CC
26	46-183039P1	1.000	No	TEXWIPE 58019-506-91/9%
29	46-208781P3	1.000	No	WRITE-ON LABEL FOR MARKIN
37	46-194427P313	1.000	No	1 INCH MULTIPLE LOCKOUT D
38	46-194427P322	1.000	No	LOCKOUT TAG,BLACK & WHITE
41	46-303460P2	1.000	No	PAINT GE GRAY #1
42	46-303460P3	1.000	No	PAINT GE STEEL BLUE
43	46-303460P4	1.000	No	PAINT GE DARK STEEL BLUE
44	46-220178P10	2.000	No	GROMMET
56	5113046	1.000	No	0.5mm Wire for Gantry Rotation Characterization Procedure,SPRING TEMPER CARBON STEEL WIRE 0.020inches Diameter, 25.4 mm length Pck Qrt 100
57	5111615	2.000	No	external laser light disconnect label

Table 5-6 Installation Support Kit, Stargate (2319891-4 BOM, rev 4) (Continued)

Section 5.0 — Balance Kit, Field, H-Power (2376682-2)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2316653	6.000	Yes	WEIGHT,107 DEG,6MM
2	2316653-2	6.000	Yes	WEIGHT,107 DEG,1.5MM
3	2317863	2.000	Yes	SHIM,180 DEG,3.2MM
4	2317865-2	2.000	Yes	SPACER,180 DEG,25.4MM
5	2317994	20.000	Yes	WEIGHT,TRIM
6	2319142	2.000	Yes	SPACER,107 DEG,25.4MM
7	2319145	2.000	Yes	SHIM,107 DEG,3.3MM
8	2342895-2	6.000	Yes	WEIGHT,180 DEG,6MM
9	2342895-3	6.000	Yes	WEIGHT,180 DEG,1.5MM
12	46-328425P4	4.000	Yes	NUT HEXAGON 16 MM
13	46-328417P70	2.000	Yes	SCREW HEXAGON SOCKET 12 M
14	2376683	1.000	Yes	SOCKET, 15/16" DEEP IMPAC
18	46-328432P5	4.000	Yes	WASHER LOCK - SPRING 16.2
10	2370053	1.000	No	WEIGHT,TRIAL,NEGATIVE Z
11	2370227	1.000	No	WEIGHT,TRIAL,POSITIVE Z,
15	2377399	1.000	No	LEFT SIDE DAS COVER TOOL
16	2377400	1.000	No	RIGHT SIDE DAS COVER TOOL
17	2377401	1.000	No	CENTER DAS COVER TOOL

Table 5-7 Balance Kit, Field, H-Power (2376682-2 BOM, rev 1)

Section 6.0 — Short Cable Collector (2281840-5)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2371450-2	1.000	Yes	RUN#55A-PDU TO RACEWAY GND
2	2371450-4	1.000	Yes	RUN#56A-CONSOLE TO RACEWAY
3	2373436-3	1.000	Yes	RUN #102A - GANTRY TO CON
4	2117848-7	1.000	Yes	20000 +/- 150 FIBER OPTIC
5	2371133-2	1.000	Yes	CONSOLE CABLES - SHORT KI
6	2371134-2	1.000	Yes	GANTRY CABLES - SHORT KIT

Table 5-8 Short Cable Collector (2281840-5 BOM, rev 0)

Section 7.0 — Long Cable Collector (2281840-4)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2371450	1.000	Yes	RUN#55-PDU TO RACEWAY GND
2	2371450-3	1.000	Yes	RUN#56-CONSOLE TO RACEWAY
3	2373436-2	1.000	Yes	RUN #102 - GANTRY TO CONS
4	2117848-2	1.000	Yes	24384 +/- 150 FIBER OPTIC
5	2371133	1.000	Yes	CONSOLE CABLES - LONG KIT
6	2371134	1.000	Yes	GANTRY CABLES - LONG KIT

Table 5-9 Long Cable Collector (2281840-4 BOM, rev 0)

© 2005-2008, General Electric Company.
GE Medical Systems, a General Electric Company, going to market as GE Healthcare.
3000 N. Grandview Boulevard
Waukesha, Wisconsin 53188
USA

www.gehealthcare.com

