

GE Healthcare

LightSpeed 7.X

Illustrated Parts List



OPERATING DOCUMENTATION



5116416-100
Rev 12

Chapter 5

Miscellaneous Items

Section 1.0 — LS VCT/Pro32 Common Non-Mobile Shipping Collector (5120455)

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	2290278-4	1.0	Yes	BASE COVER COLLECTOR {For further detail, see Section 1.1, beginning on page 360. }
2	5125429	1.0	Yes	Fuse Kit - Originally developed for VCT {For further detail, see Section 1.2, beginning on page 361. }
3	5139186-2	1.0		VCT, Failsafe Balance Weight Kit for Pathfinder
4	5130493	1.0		VCT GANTRY AND HP OR GT TABLE COMMON SHIPPING COLLECTOR {For further detail, see Section 1.3, beginning on page 362. }

Table 5-1 LS VCT/Pro32 Gantry & HP/GT Table Common Non-Mobile Shipping Collector (5120455 BOM)

1.1 Base Cover Collector (2290278-4)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION
4	46-328425P1	5.0	Yes	NUT HEXAGON 4 MM
5	2251467-4	1.0	Yes	COVER, BASE-FRONT
6	2285829-7	1.0	Yes	COVER, LEFT BASE WITH CARDIAC OPT
7	2285830-7	1.0	Yes	COVER,RIGHT BASE WITH CARDIAC OPT
8	2275686-4	1.0	Yes	COVER ASSY, BASE REAR LHY
9	2275685-4	1.0	Yes	COVER ASSY, BASE REAR RH
1	46-328417P3	4.0		SCREW HEXAGON SOCKET 4 MM 10 MM ZINC PLATED, STEEL
2	46-328430P2	25.0		WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL
3	46-328432P2	29.0		WASHER LOCK - SPRING 4.1 MM 7.6 MM ZINC PLATED, SPRING STEEL
10	2277652	3.0		LATCH, H2 SIDE COVER COMPRESSION
11	2261321	2.0		BRACKET - FRONT BASE COVER
12	2261324	2.0		BRACKET - SIDE BASE COVER
13	2261325	1.0		BRACKET - REAR BASE COVER
14	2285469	2.0		BRACKET,BASE COVER MOUNTING
15	2288882	2.0		BULKHEAD - BLANK
16	46-328420P4	4.0		SCREW BUTTON HEAD 4 MM 16 MM BLACK OXIDE, STEEL
18	46-328417P4	17.0		SCREW HEXAGON SOCKET 4 MM 12 MM ZINC PLATED, STEEL

Table 5-2 Base Cover Collector (2290278-4 BOM, Rev 0)

1.2 Fuse Kit (5125429)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2142415-4	3.0	Yes	PROTECTION FUSE, 5A, TIME DELAY, 600VAC 250VDC, BOARDS USE FUSE CLIPS
2	2142415-5	3.0	Yes	PROTECTION FUSE, 8A, TIME DELAY, 600VAC 250VDC, BOARDS USE FUSE CLIPS
3	2142415-7	3.0	Yes	PROTECTION FUSE, 15A, TIME DELAY, 600VAC 250VDC, BOARDS USE FUSE CLIPS
4	2142415-9	2.0	Yes	PROTECTION FUSE, 25A, TIME DELAY, 600VAC 250VDC, CHASSIS MOUNTED-- CURRENT RATING HIGHER THAN FOR TYPICAL WITH FUSE CLIPS
5	2336517-2	2.0	Yes	FUSE 25 AMPS 700 VOLTS 2.5M SECONDS
6	2351492	3.0	Yes	200A FUSE
7	2364059	2.0	Yes	GLASS FUSE
8	2379651	2.0	Yes	FUSE - 700 V, 200 A, FAST SEMI
9	46-170021P106	2.0	Yes	FUSE 8 AMPS 250 VOLTS BUSSMANN ONLY
10	46-170021P75	2.0	Yes	015.000A 0250V 3AB SLOW UL BUSSMANN MDA 15
11	46-267217P20	4.0	Yes	FUSE METRIC SLO BLO 2.500 A 250V 5X20MM GLASS LITTLEFUSE 218 02.5
12	46-267217P23	2.0	Yes	FUSE METRIC SLO BLO 5.000 A 250V 5X20MM GLASS LITTLEFUSE 218 005
13	5118644	1.0	Yes	FAST ACTION 15 AMPS 250 VOLTS 0.01 SECONDS U0028FF
14	U8006FF	3.0	Yes	FUSE 1.0A
15	U8008FF	3.0	Yes	FUSE 2.0A
16	U8009FF	3.0	Yes	FUSE 3.2A
17	2390269	3.0	Yes	FUSE 2.5 AMPS 125 VOLTS

Table 5-3 Fuse Kit - Originally Developed for VCT (5125429 BOM, Rev 2)

1.3 VCT Gantry & HP/GT Table Common Shipping Collector (5130493)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION
100	2277548	4.0	Yes	PLUG, SIDE COVER H2
104	5180151	2.0	Yes	Rear cover dolly
105	5179183	1.0	Yes	Dolly, pivoting, small room, RH
106	5179183-2	1.0	Yes	Dolly, pivoting, small room, LH
107	5137423	1.0	Yes	hoist arm_VCT
108	5138605	1.0	Yes	tube, hoist arm support-VCT
115	5137422	1.0	Yes	bracket-hoist extension_VCT
116	46-199678P3	2.0	Yes	RING GRIP QUICK RELEASE PIN - .375 DIA - 3.50 RANGE
210	2364629	1.0	Yes	LABEL INSTALLATION KIT
300	2336557	1.0	Yes	SEISMIC KIT FOR GLOBAL CONSOLE
301	2373436-7	1.0	Yes	Cable, FE Laptop to Console, 23000mm Long
	5220799	1.0		Belt Frequency Measurement Class C Software
	5224640	1.0		Microphone to be used with belt frequency measurement tool.
102	2105873	10.0		.531 I.D., 1.625 O.D., .25THK GANTRY & TABLE INSTALLATION ANCHOR WASHER, FOR BASEPLATE-LESS HSA SYSTEMS.
103	2106573	10.0		HILTI ANCHOR STUD 1/2 DIA X 8.0 LG, KWIK-BOLT II WITH EXTENDED THREAD, FOR BASEPLATE-LESS HSA SYSTEMS STUD MAT~L: CARBON ST PER ASTM A510 WEDGE MAT~L: CARBON ST PER AISI 1010 NUT MAT~L: CARBON ST PER ASTM A563 GRADE A WASHER MAT~L: CARBON ST PER SAE 1005-
109	2217173-2	1.0		STOOL WITH FOOT REST AND NO ARM REST
110	5112283	1.0		CT Service Cabinet from Pentair
111	5198409	1.0		VCT Shipping Collector Kit
112	5111499	1.0		FLOOR TEMPLATE-VCT
200	46-237412P1	1.0		STRAP HEAD STRAP 3 STRAPS PER PACKAGE
201	46-222443P1	1.0		OPEN CELL POLYESTER FOAM
202	46-222444P1	1.0		CLOSED CELL LINEAR ETHYLENE
203	46-222445P1	1.0		CLOSED CELL LINEAR ETHYLENE

Table 5-4 VCT Gantry & HP/GT Table Common Shipping Collector (5130493 BOM, Rev 16)

MARK	NUMBER	QTY	FRU	DESCRIPTION
204	46-222921P1	6.0		OBSOLETE FOR MR SITES PER FMI 60288 POLYURETHANE FOAM 1.7#/FT DENSITY S.E. FINE CELL CHARCOAL GRAY
205	46-199388P1	1.0		HEAD HOLDER POLYURETHANE FOAM
206	46-222436G1	4.0		CT-9800 - CHIN STRAP
207	46-297960P1	1.0		FOAM INSERT
208	46-229452G1	1.0		CT-9800 (ACCESSORIES)
209	2309994-2	1.0		IV POLE ASSY, with Multi-Language Labels
212	46-297629P1	2.0		STRAP
213	2329064-2	1.0		TABLE TRAY ACC ASSY, with Multi-Language Lables
214	2355100	1.0		TRAY, H2 TABLE
217	46-236664P1	2.0		STRAP STRAP, WRIST 2.00 WIDE X 18.00 LONG WHITE

Table 5-4 VCT Gantry & HP/GT Table Common Shipping Collector (5130493 BOM, Rev 16)

MARK	NUMBER	QTY	FRU	DESCRIPTION
100	2277548	4.0	Yes	PLUG, SIDE COVER H2
104	2282011	2.0	Yes	NON-TILT DOLLY REAR COVER {For further detail, see Section 1.3.1, beginning on page 365. }
105	5101127	1.0	Yes	dolly pivoting , right hand
106	5101127-2	1.0	Yes	dolly, pivoting, left hand
107	5137423	1.0	Yes	hoist arm_VCT
108	5138605	1.0	Yes	tube, hoist arm support-VCT
115	5137422	1.0	Yes	bracket-hoist extension_VCT
116	46-199678P3	2.0	Yes	RING GRIP QUICK RELEASE PIN - .375 DIA - 3.50 RANGE
210	2364629	1.0	Yes	LABEL INSTALLATION KIT
300	2336557	1.0	Yes	SEISMIC KIT FOR GLOBAL CONSOLE
102	2105873	10.0		.531 I.D., 1.625 O.D., .25THK GANTRY & TABLE INSTALLATION ANCHOR WASHER, FOR BASEPLATE-LESS HSA SYSTEMS.

Table 5-5 VCT Gantry & HP/GT Table Common Shipping Collector (5130493 BOM, Rev 8)

MARK	NUMBER	QTY	FRU	DESCRIPTION
103	2106573	10.0		HILTI ANCHOR STUD 1/2 DIA X 8.0 LG, KWIK-BOLT II WITH EXTENDED THREAD, FOR BASEPLATE-LESS HSA SYSTEMS STUD MAT~L: CARBON ST PER ASTM A510 WEDGE MAT~L: CARBON ST PER AISI 1010 NUT MAT~L: CARBON ST PER ASTM A563 GRADE A WASHER MAT~L: CARBON ST PER SAE 1005-
109	2217173-2	1.0		STOOL WITH FOOT REST AND NO ARM REST
110	2299113	1.0		CT SERVICE CABINET
111	5133492	1.0		VCT Pro 32 Install Support Kit
112	5111499	1.0		FLOOR TEMPLATE-VCT
113	5133490	1.0		VCT Service Tools Collector
114	5149370	1.0		64 slice DAS Detector Service Kit
200	46-237412P1	1.0		STRAP HEAD STRAP 3 STRAPS PER PACKAGE
201	46-222443P1	1.0		OPEN CELL POLYESTER FOAM
202	46-222444P1	1.0		CLOSED CELL LINEAR ETHYLENE
203	46-222445P1	1.0		CLOSED CELL LINEAR ETHYLENE
204	46-222921P1	6.0		OBSOLETE FOR MR SITES PER FMI 60288 POLYURETHANE FOAM 1.7#/FT DENSITY S.E. FINE CELL CHARCOAL GRAY
205	46-199388P1	1.0		HEAD HOLDER POLYURETHANE FOAM
206	46-222436G1	4.0		CT-9800 - CHIN STRAP
207	46-297960P1	1.0		FOAM INSERT
208	46-229452G1	1.0		CT-9800 (ACCESSORIES) {For further detail, see Section 1.3.2, beginning on page 365. }
209	2309994-2	1.0		IV POLE ASSY, with Chinese label
211	46-266041P1	1.0		NARROW RESTRAINT
212	46-297629P1	2.0		STRAP
213	2329064-2	1.0		TABLE TRAY ACC ASSY, with Chinese label.
214	2355100	1.0		TRAY, H2 TABLE
215	46-237605P1	1.0		WRIST STRAP
216	46-266042P1	1.0		WIDE RESTRAINT
217	46-236664P1	2.0		STRAP STRAP, WRIST 2.00 WIDE X 18.00 LONG WHITE

Table 5-5 VCT Gantry & HP/GT Table Common Shipping Collector (5130493 BOM, Rev 8)

1.3.1 Non-Tilt Dolly Rear Cover (2282011)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION
3	2277193	2.0	Yes	DOLLY TO COVER ATTACHMENT PIN
4	2279717	2.0	Yes	WING NUT & WASHER SPECIAL FASTENER
1	2282012	1.0		NON-TILT FRAME WELDED
2	2277434	2.0		CASTER, HURON JR 2~
5	46-328417P12	8.0		SCREW HEXAGON SOCKET 6 MM 12 MM ZINC PLATED, STEEL
6	46-328425P2	8.0		NUT HEXAGON 6 MM
7	46-170686P2	0.0		LOCTITE 242, 10CsC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER

Table 5-6 Non-Tilt Dolly Rear Cover (2282011 BOM, Rev 1)

1.3.2 CT-9800 Accessories (46-229452G1)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION
6	46-229450P1	1.0	Yes	STRAP - SECURITY, ARM OBSOLETE FOR MR SITES PER FMI 60288 VINYL COATED POLYESTER FABRIC W/TAPE FASTENERS
2	46-278968P2	1.0		No description provided
4	46-278986P2	2.0		No description provided

Table 5-7 CT-9800 Accessories (46-229452G1 BOM, Rev 3)

Section 2.0 — Cable Collectors

2.1 Short Cable Collector (2281840-7)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	2371450-2	1.0	Yes	RUN#55A-PDU TO RACEWAY GND - SHORT
2	2371450-4	1.0	Yes	RUN#56A-CONSOLE TO RACEWAY GND-SHORT
3	2373436-3	1.0	Yes	RUN #102A - GANTRY TO CONSOLE LAN
4	5125259-2	1.0	Yes	Fiber Optic - Console to Gantry Short
5	2371133-4	1.0	Yes	Console Cable - Short VCT Kit {For further detail, see Section 2.1.1, below. }
6	2371134-4	1.0	Yes	Gantry Cables - Short VCT Kit {For further detail, see Section 2.1.2, below. }
7	2373436-6	1.0	Yes	LAN CABLE, CARDIAC INTERFACE PANEL TO CONSOLE, 18280mm LONG

Table 5-8 Short Cable Collector (2281840-7 BOM, Rev 3)

2.1.1 Console Cable - Short VCT Kit (2371133-4)

FRU = YES

MARK	NAME	QTY	FRU	DESCRIPTION
1	5121809-2	1.0		VCT Console Pwr
2	5120645-2	1.0		Cable - Console to MSUB J9 - EMC Edition 2

Table 5-9 Console Cable - Short VCT Kit (2371133-4 BOM, Rev 1)

2.1.2 Gantry Cables - Short VCT Kit (2371134-4)

FRU = YES

MARK	NAME	QTY	FRU	DESCRIPTION
1	2343529-2	1.0		CABLE, HVDC FROM PDU TO GANTRY
2	2343528-2	1.0		CABLE PDU TO GANTRY 120VAC POWER
3	2343530-2	1.0		CABLE, AXIAL DRIVE POWER PDU TO GANT
4	5120646-2	1.0		Cable - PDU to MSUB J11 - EMC Edition 2

Table 5-10 Gantry Cables - Short VCT Kit (2371134-4 BOM, Rev 2)

2.2 Long Cable Collector (2281840-6)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2371450	1.0	Yes	RUN#55-PDU TO RACEWAY GND - LONG
2	2371450-3	1.0	Yes	RUN#56-CONSOLE TO RACEWAY GND-LONG
3	2373436-2	1.0	Yes	RUN #102 - GANTRY TO CONSOLE LAN
4	5125259	1.0	Yes	Fiber Optic - Console to Gantry
5	2371133-3	1.0	Yes	Console Cables - Long VCT Cables {For further detail, see Section 2.2.1. }
6	2371134-3	1.0	Yes	Gantry Cables - Long VCT Kit {For further detail, see Section 2.2.2. }
7	2373436-4	1.0	Yes	LAN CABLE, CARDIAC INTERFACE PANEL TO CONSOLE, 26350mm LONG

Table 5-11 Long Cable Collector (2281840-6 BOM, Rev 3)

2.2.1 Console Cables - Long VCT Cables (2371133-3)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	5121809	1.0		VCT Console Pwr
2	5120645	1.0		Cable Console to MSUB J9 - EMC Edition 2

Table 5-12 Console Cables - Long VCT Cables (2371133-3 BOM, Rev 1)

2.2.2 Gantry Cables - Long VCT Kit (2371134-3)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2343529	1.0		CABLE, HVDC FROM PDU TO GANTRY
2	2343528	1.0		CABLE PDU TO GANTRY 120VAC POWER
3	2343530	1.0		CABLE, AXIAL DRIVE POWER PDU TO GANT
4	5120646	1.0		Cable - PDU to MSUB J11 - EMC Edition 2

Table 5-13 Gantry Cables - Long VCT Kit (2371134-3 BOM, Rev 1)

Section 3.0 — IVY 3150 ECG Monitor Collector

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	5112819	1.0	Yes	3150 IVY ECG MONITOR
2	5112819-2	1.0	Yes	ECG, 3-LEAD, LOW-NOISE, PATIENT CABLE, 3m LONG
3	5112819-3	1.0	Yes	ECG, LOW-NOISE, LEAD WIRES, 24in LONG, SET OF 3
4	5112819-4	1.0	Yes	ECG, ROLLING MONITOR STAND
6	2373436-4	1.0	Yes	LAN CABLE, CARDIAC INTERFACE PANEL TO CONSOLE, 26350mm LONG
7	2373436-5	1.0	Yes	LAN CABLE, CARDIAC INTERFACE PANEL TO ECG MONITOR, 1000mm LONG
8	5122589	1.0	Yes	Adapter, Bulkhead RJ45 CAT5e Shielded, L-com #TDG1026KS-C5E
5	5112819-5	1.0		ECG, RADIOTRANSLUCENT ELECTRODES, 20 PACKS OF 3

Table 5-14 IVY 3150 ECG Monitor Collector (5133408 BOM, Rev 1)


Section 4.0 — Small Room Catalog (5191743)

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	5179183	1.0		Dolly, pivoting, small room, RH {For further detail, see Section 4.1, beginning on page 370. }
2	5179183-2	1.0		Dolly, pivoting, small room, LH {For further detail, see Section 4.2, beginning on page 372. }
3	5184632	1.0		Small Room CT Gantry Installation Dolly
4	2373436-7	1.0		Cable, FE Laptop to Console, 23000mm Long


Table 5-15 Small Room Catalog (5191743 BOM, Rev 1)

4.1 Dolly, Pivoting, Small Room, RH (5179183)

FRU = NO



LH
5179183-2
(FOR SMALL ROOM)



RH
5179183
(FOR SMALL ROOM)

Figure 5-1 Dolly, Pivoting, Small Room (5179183 Sheet 1, Rev 3)

MARK	NUMBER	QTY	FRU	DESCRIPTION
9	5176413	2.0	Yes	CLAMPING LEVER, ADJUSTABLE M10 X 20MM LG MORTON MACHINE WORKS #AH-933
1	5154481	1.0		Post, adjustable for pivoting dolly used on VCT gantry front cover
2	5164802	1.0		Post, upper for pivoting dolly used on VCT gantry front cover
3	5173585	1.0		Leg, stationary for pivoting dolly used on VCT gantry front cover
4	5165043	1.0		Leg, pivoting for pivoting dolly used on VCT gantry front cover
6	5173375	1.0		Bracket, pivot, RH for pivoting dolly used on VCT gantry front cover
8	2290611	2.0		PIN-DOLLY TO COVER CONNECTING
10	5179332	1.0		LANYARD, .0469 X 18 TERMINAL AND LOOP, FASTENAL #0157360
11	2283554	1.0		PIN, H2 DOLLY PIVOT LOCKING
12	2277431	1.0		KNOB, PLASTIC KNURLED CAMPING KNOB
13	2277430	1.0		E-RING, BOWED RETAINING RING
14	2277432-2	1.0		STEEL COMPRESSION SPRING_.36 OD 1.50
15	2283321	2.0		BUSHING, FLANGED NYLON .26 ID .50 LG
16	2283320	2.0		SPEED NUT, SHAFT DIA. 8MM
17	46-220205P7	2.0		ID 0.812 OD 2.00 THK 0.156 F70B5A (FOR 3/4 IN.SCREW)
18	46-208949P8	4.0		003/4-10 1.125 HEX X .641 THK F70B5A
19	5179690	2.0		COTTERLESS HITCHPIN, .375 X 2.375 FASTENAL #0157085
20	5179187	2.0		LANYARD, .093 X 8 PIVOT POINT #NL93F-8-2-26
21	2278991	3.0		PRESS LOCK RIVET
22	5166930	2.0		Spacer, caster for pivoting dolly used on VCT gantry front cover
23	5177454	4.0		CASTER, SWIVEL, 3.5 HIGH, SHEPHERD CASTER CORP. #102241
24	5176832	1.0		BOLT, HCS .75-10 X 2 GR5 ZINC
25	5179441	2.0		WASHER, SAE FW .75 ZINC, ID 0.812, OD 1.47, THK 0.134
26	5176393	1.0		COTTERLESS HITCHPIN, .375 X 1, FASTENAL #0157080
27	46-170684P2	0.0		LOCTITE 271, 10CC BOTTLE 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER. (REV, TK, 1/92).
28	5178519	16.0		NUT, M6 X 1.0, 9.8 HEX X 5.8 THK, NYLOCK, ZINC

Table 5-16 Dolly, Pivoting, Small Room, RH (5179183 BOM, Rev 2)

4.2 Dolly, Pivoting, Small Room, LH (5179183-2)

FRU = NO





Figure 5-2 Dolly, Pivoting, Small Room (5179183 Sheet 1, Rev 3)

MARK	NUMBER	QTY	FRU	DESCRIPTION
9	5176413	2.0	Yes	CLAMPING LEVER, ADJUSTABLE M10 X 20MM LG MORTON MACHINE WORKS #AH-933
1	5154481	1.0		Post, adjustable for pivoting dolly used on VCT gantry front cover
2	5164802	1.0		Post, upper for pivoting dolly used on VCT gantry front cover
3	5173585	1.0		Leg, stationary for pivoting dolly used on VCT gantry front cover
4	5165043	1.0		Leg, pivoting for pivoting dolly used on VCT gantry front cover
7	5174002	1.0		Bracket, pivot, LH for pivoting dolly used on VCT gantry front cover
8	2290611	2.0		PIN-DOLLY TO COVER CONNECTING
10	5179332	1.0		LANYARD, .0469 X 18 TERMINAL AND LOOP, FASTENAL #0157360
17	46-220205P7	2.0		ID 0.812 OD 2.00 THK 0.156 F70B5A (FOR 3/4 IN.SCREW)
18	46-208949P8	4.0		003/4-10 1.125 HEX X .641 THK F70B5A
19	5179690	2.0		COTTERLESS HITCHPIN, .375 X 2.375 FASTENAL #0157085
20	5179187	2.0		LANYARD, .093 X 8 PIVOT POINT #NL93F-8-2-26
21	2278991	3.0		PRESS LOCK RIVET
22	5166930	2.0		Spacer, caster for pivoting dolly used on VCT gantry front cover
23	5177454	4.0		CASTER, SWIVEL, 3.5 HIGH, SHEPHERD CASTER CORP. #102241
24	5176832	1.0		BOLT, HCS .75-10 X 2 GR5 ZINC
25	5179441	2.0		WASHER, SAE FW .75 ZINC, ID 0.812, OD 1.47, THK 0.134
26	5176393	1.0		COTTERLESS HITCHPIN, .375 X 1, FASTENAL #0157080
27	46-170684P2	0.0		LOCTITE 271, 10CC BOTTLE 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER. (REV, TK, 1/92).
28	5178519	16.0		NUT, M6 X 1.0, 9.8 HEX X 5.8 THK, NYLOCK, ZINC

Table 5-17 Dolly, Pivoting, Small Room, LH (5179183-2 BOM, Rev 2)

End of Chapter.

© 2005, General Electric Company.

GE Medical Systems, a General Electric Company, going to market as GE Healthcare.

3000 N. Grandview Boulevard

Waukesha, Wisconsin 53188

USA

www.gehealthcare.com

