

GE Healthcare

Discovery MR450 & MR750
Field Replaceable Units (FRU)

Technical Publication

© 2010 General Electric Company
GE Medical Systems, a General Electric Company, going to market as GE Healthcare
3000 N. Grandview Boulevard
Waukesha, Wisconsin 53188
USA

www.gehealthcare.com



5500105
Rev 4

**COVER
BOTTOM
CROP LINE**

**COVER
BOTTOM
CROP LINE**

**COVER
LEFT
CROP LINE**

**COVER
FOLD
LINE** | **(2) CENTER STAPLES
FOR BINDING**

**COVER
RIGHT
CROP LINE**

**5/16" RADIUS
CORNER TRIM
TYPICAL**

**5/16" RADIUS
CORNER TRIM
TYPICAL**

SEND COMMENTS ABOUT THIS MANUAL TO:

GE MEDICAL SYSTEMS
MANAGER, SERVICE ENGINEERING
W-808
3200 N. GRANDVIEW BLVD.
WAUKESHA, WISCONSIN 53188
USA

This manual belongs to:

INTRODUCTION TO FIELD REPLACEABLE UNIT (FRU)

This direction is a guide to the FRU's.

Parts designated as a FRU are critical to maintain the safety and/or functionality of the system. FRU's are available on a Priority 1 basis from Service Parts Operation.

As a reference aid, the FRU's are listed under an upper-level part number where they are used. In most cases the upper-level part number is not a FRU. In some cases the parts listed under an upper-level assembly may have been nested in intermediate sub-assemblies which are not shown.

Where possible, FRU's which are part of Option Catalogs are shown under the area where they would be installed, rather than being listed in the Coils or Options Sections (SEE TABLE OF CONTENTS).


LIST OF EFFECTIVE PAGES


<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>
Cover Page ..	4	16 to 19	4	32	3
1	4	20 to 21	3	33 to 36	4
2 to 4	3	22	4	37 to 38	3
5 to 9	4	23 to 24	3	39 to 40	4
10 to 13	3	25 to 28	4	41 to 43	3
14	4	29 to 30	3	44	4
15	3	31	4		

[illegible]


TITLE	PAGE
INTRODUCTION TO FIELD REPLACEABLE UNIT(FRU)	1
LIST OF EFFECTIVE PAGES	1
TABLE OF CONTENTS	3
MAGNET ENCLOSURE	
COMMON ENCLOSURE COMPONENTS	5
ENCLOSURE COMPONENTS FOR IN-ROOM DISPLAY	6
ENCLOSURE COMPONENTS FOR 7-SEGMENT DISPLAY	7
ENCLOSURE ELECTRONIC COMPONENTS	8
ENCLOSURE & SHIM COMPONENTS	9
MR450 1.5T GRADIENT and RF BODY COILS	10
MR750 3.0T GRADIENT and RF BODY COILS	11
REAR PEDESTAL	12
LPCA	13
RF HUB	14
DOCK	16
PATIENT TABLE	17 & 18
CRADLE	19
TABLE TOP	19
HEAD COIL and PATIENT HANDLING	20
SHIELD COOLER	22
HEC CABINET	23
PGR CABINET	25 & 26
CRONUS MODULE	27
CAM MODULE	28
PEN CABINET	29
PEN CABINET PENETRATION PANEL	30
SPW	31
MAIN DISCONNECT PANEL	33
OPERATOR WORKSPACE	34
GOC COMPUTER CABINET with HP9300	35
SCIM AND KEYBOARD	36
MR750 3.0T COILS	37 & 38
UNIFIED SURFACE COIL PHANTOM SET	39
COIL STAGING FRU'S	39
SERVICE KITS	39
CABLE COLLECTORS	40
MNS 8kW RF AMP CABINET (OPTIONAL)	41
MNS COMPONENTS	42
REVISION HISTORY	44

[illegible]

COMMON ENCLOSURE COMPONENTS**MAGNET ENCLOSURE**

ENCLOSURE COMPONENTS for IN-ROOM DISPLAY (IRD)

**ENCLOSURE COMPONENTS for SEVEN SGMENT DISPLAY
(without IRD)**



MAGNET ENCLOSURE

ENCLOSURE ELECTRONIC COMPONENTS**COLD HEAD FRU's**


2154502-2	Run 621, Supply Gas Line, 65.6 ft or 20M long
2154505-2	Run 622, Return Gas Line, 65.6 ft or 20M long
2155316-2	Run 623, Compressor to Pen Panel Cable, 49.21 ft or 15M long
2218292	Run 624, Cold Head Control Cable
2214062	Run 832, LCC Sleeve and Cold Head Cable
2259837	LCC Cold Head Control Cable Asm


MAGNET FRU's

46-260852G3	Emergency Evacuation Release Tool
--------------------	-----------------------------------

CABLE FRU's

5304950	RF Cable FRU's for Body Hybrid to Body Coil (Contains Runs M1322 and M1323)
5183293	Microphone Cable, Run M3345
5137339-6	E-Stop Cable
5180907	Run 3350, MAG-SRI-J10 to MAG-BORE, Bore Temp,



ENCLOSURE AND SHIM COMPONENTS**MAGNET ENCLOSURE**

MR450 1.5T GRADIENT AND RF BODY COIL**2208999-14**

1.5T RF BODY COIL Includes
ORANGE & BLACK FRU CASE
AND SKID

2144093-12

GRADIENT COIL with HOS
Includes GRAY COIL CART
WITH CRADLE

5316997**COOLANT MANIFOLD COLLECTOR**

FRU contains:

- 5268005, Coolant Supply Manifold – BLU
- 5268152, Coolant Supply Manifold – BLK
- 46-208765P2Hose Clamps (25)
- 46-294167P91 Electrical Tape

5307840

BUSBAR FRU ASM contains:

- 5305786, Leadboard
- 46-294167P96, Bushing (8)
- 5268439, Axial Adj. Shim (12)
- 5214524-6, M10 SS Washer (20)
- 2341885, M10 Brass Nut (22)
- 5308479, Installation Gage Tool

5197020

HOS CONNECTOR
TERMINATION CAP

5191537 (1.5T)

Rear Turtle Bracket Asm

5335389

1.5T GRADIENT MOUNTING FRU KIT contains:
(For mounting front and rear of coil)

- 5170015 Radial Positioning Block (2)
- 5170098 Axi-Radial Positioning Block (4)
- 5170312 Inward Mounting Wedge (5)
- 5170035 Outward Mounting Wedge (4)
- 5326293-4 M10 x 30mm Hex Head Cap Screw (15)
- 5146169 M10 x 24mm Washer (15)
- 46-294151P11 Locktite 271

MR750 3.0T GRADIENT AND RF BODY COIL**2208999-13**

3T RF BODY COIL Includes
ORANGE FRU CASE AND SKID

2144093-12

GRADIENT COIL with HOS
Includes GRAY COIL CART
WITH CRADLE

5316997

COOLANT MANIFOLD COLLECTOR

FRU contains:

- 5268005, Coolant Supply Manifold – BLU
- 5268152, Coolant Supply Manifold – BLK
- 46-208765P2Hose Clamps (25)
- 46-294167P91 Electrical Tape

5307840

BUSBAR FRU ASM contains:

- 5305786, Leadboard
- 46-294167P96, Bushing (8)
- 5268439, Axial Adj. Shim (12)
- 5214524-6, M10 SS Washer (20)
- 2341885, M10 Brass Nut (22)
- 5308479, Installation Gage Tool

5197020

HOS CONNECTOR
TERMINATION CAP


5191537-2 (3.0T)


Rear Turtle Bracket Asm

5307838


3.0T G3GRADIENT MOUNTING FRU KIT contains:
(For mounting front and rear of coil)

- 5155280 G3 Radial Positioning Block (2)
- 5154485 G3 Axi-Radial Positioning Block (4)
- 5166786 G3 Inward Mounting Wedge (5)
- 5166458 G3 Outward Mounting Wedge (4)
- 5110638-64 M10 x 30 Brass Bolt (15)
- 5146169 M10 x 24mm Washer (15)
- 46-294151P11 Locktite 271


REAR PEDESTAL


LPCA


RF HUB




[illegible]

MAGNET DOCKING PARTS:


PATIENT TABLE**PATIENT TABLE**


PATIENT TABLE

PATIENT TABLE CRADLE**PATIENT TABLE TOP FRU's**

The Table Bumpers
FRU Kit, **5306407**,
contains:

5168057
5168057-2
5168326
5168326-2
5168193 (2)

**PATIENT TABLE**

SPLIT HEAD COIL:**Note: The following parts are not shown:**

5263175 (3.0T) QHTR Circuit Board Asm Kit

5328732 (1.5T) QHTR Circuit Board Asm Kit

5267144 Head Coil Cushion – Flat

5271434 Head Coil Cushion – Head

2355094 Replacement Spring Contacts Kit

PATIENT HANDLING FRU KIT's:

46-306574P7 0.5 x 9 x 14.5 Gray Rectangular Pad, Nylon Covered

46-306574P8 1.0 x 9 x 14.5 Gray Rectangular Pad, Nylon Covered


46-306574P9 0.5 x 9 x 36 Gray Rectangular Pad, Nylon Covered

46-306574P10 1.0 x 9 x 36 Gray Rectangular Pad, Nylon Covered

46-317758P1 Pneumatic Patient Alert System

5122500 ECG Lead Kit

[illegible]

SHIELD COOLER

2188440-2
Sumitomo CSW-71D
Shield Cooler Compressor

OR

5314887
Sumitomo F-50
Shield Cooler Compressor

SHIELD COOLER FRU's Not Shown:

2154502-2 Run 621, Supply Gas Line, 65.6 ft or 20M long

2154505-2 Run 622, Return Gas Line, 65.6 ft or 20M long

2155316-2 Run623, Compressor to Pen Panel Cable, 49.21 ft or 15M long

2218192 Run624, Cold Head Control Cable

2214062 Run 832, LCC Sleeve and Cold Head Cable

2259837 LCC Cold Head Control Cable Asm

2200835 Power Cable for Sumitomo CSW-71D Shield Cooler Compressor

2211446-3 Compressor Cable Kit consists of:

2154502-2, Run 621


2154505-2, Run 622

2259837, LCC Cold Head Control Cable

2200835, Shield Cooler Compressor Power Cable

5335002 Flow Control Adapter, used with F-50 Compressor only

HEC CABINET



HEC (Heat Exchanger Cabinet) FRU KIT's


5264701	Deionization Filter Cartridges Kit
5156315	3hp Blower Kit
5156316	Cryogen Compressor Flow Meter Kit
5156318	Mixing Valve Kit for Gradient Coil and PGR Cooling Loops
5156319	Power Box Kit
5156320	Mixing Valve Kit for Chilled Air Blower
5156321	5hp Pump Kit
5262029	Equipment Room Hoist Support Tube
5266827	Fuse Kit
5267643	Liquid Pressure Sensor Kit
5267644	Level Switch Kit
5267648	Pressure Difference Sensor Board Kit
5267649	Resistivity Sensor Kit
5268665	Liquid Temperature Sensor Kit
5269683	Coolant Removal Kit
5269685	Cryogen Compressor Filter Kit
5269689	Air Temperature Sensor Kit
5269691	Signal Box Kit

HEC COOLING SYSTEM FRU's

5267596	Liquid Cooling Valve Kit (Contains ball valves and check valves)
5174313-4	Gradient Coil Coolant, 4 Gallons – non-copper corrosive, algae resistant

[illegible]

PGR CABINET




PGR CABINET

5168052
7-16 FEMALE TO 7-16
FEMALE BULKHEAD

5263113
ICN 1 TO ICN 2
CONNECT CABLE

5257687-2
ICN 2 POWER CABLE

2379186
TERM SERVER POWER
SUPPLY/CORD

**REAR VIEW****Internal Transmit Cable FRU's**

5167234-2 3.0T Body TX cable, LMR600
5167234-3 1.5T Body TX cable, LMR600
5167234-4 1.5T Body TX cable, LMR600
5192690-2 3.0T Head TX cable, LMR600
5192690-3 1.5T Head TX cable, LMR600
5192690-4 1.5T Head TX cable, LMR600

5304065, Internal Signal Cable FRU Kit contains:

5178931	5165870	5165970	5176566	5182238	5182238-2
5182238-3	5182238-4	5184303	5257682	5257683	5257683-2
5257684	5257685	5257685-2	5257685-3	5257685-4	5257685-5
5257685-6	5257685-10	5174739	5167344	5167005	5257680
5257681	5166354	5166222-13	5166222-18	5257685-11	5257685-12
5257685-13	5176990	5176262	5176262-2	5176262-3	


5304386, Internal Fiber Optic Cable FRU Kit contains:



5182888	5182317	5182317-2	5182317-3	5182317-4	5182317-5
5182317-6	5257686				



5304635, Internal Power Cable FRU Kit contains:

5173399	5164904	5165800	5172391	5174242	5174242-2
5183977					

CRONUS MODULE



CAM MODULE (FRONT)**CAM MODULE (REAR)**

PEN CABINET**5199042
TOP IO PANEL****Top of Pen Cabinet****INTERNAL CABLE FRU KIT'S****5308001-2, Power Cable FRU Kit contains:**

5167143	5167143-2	5172838	5155030	5172849	5171737
5154994					


5308001-3, RF Cable FRU Kit contains:

5181464	5166512-3	5166512-2	5172955-2	5172955	5166749
5166749-9	5172679-10	5172679-11	5172679		


5308001-4, Power Distribution Cable FRU Kit contains:

2322434	5168641	5168043	5168317	5167882	5166461
---------	---------	---------	---------	---------	---------

PENETRATION CABINET

PEN CABINET PENTRATION PANEL**REAR VIEW****INTERNAL CABLE FRU KIT'S (Continued from previous page)****5308001-5, Signal Cable FRU Kit contains:**


5165866	5166997	5169443	5172726	5166354-2	5171626
5170521	5172868	5167365	5171583	5172959	5169381
5168461	5168597	5172812	5172595-2	5169345	5175777
5172595					


SPW (Secondary PEN Wall)

SPW

31

[illegible]

MAIN DISCONNECT PANEL (MDP)**MDP FRU's Not Shown****5161287-2** 1.5A, 380/480V Fuse**5161287-3** 60A, 380/480V Fuse**5161287-4** Green PB Lens**5161287-5** Green PB LED 120V Power Supply**2375862-2** 200A, 600V Fuse**2276547-14** 12V Battery**2276547-16** 12V Battery Charger

OPERATOR WORKSPACE**OPERATOR WORKSPACE FRU's Not Shown**

2156158-63 1.5T Software Collector for 20.1 Release

2156158-74 3.0T Software Collector for 20.x Release M4

5165713 DVD Media for Linux

GOC COMPUTER CABINET WITH HP8400;**5267342**

GOC CableCollector
(See details below)

5133768-2
FIBER OPTIC
CONVERTER
BOX

5139480-2
8-PORT SWITCH

5262750
PDU ASM

2395043 or 2395043-2
GOC AUDIO ASM

5183547-30

HOSTPC

5183547-15
73GB SAS HARD
DRIVE FOR PC

5183547-20
SINGLE 1GB DIMM

5183547-36
SINGLE 2GB DIMM

5183547-22
DVD-ROM DRIVE

5184500
USB to RS422
CONVERTER
(Connected to
Ethernet Port)

46-271382P6
1 POLE, 15 AMP
CIRCUIT BREAKER

5265358

One-Wire Collector
(See details below)

5267342 GOC Cable Collector FRU consists of;

5270774 Host Power Cable

5270803 Main LCD Monitor Power Cable

5263413 Run E3046 Main Monitor DVI-D Cable

5264533 Run E3047 SCIM-GOCAA Cable

5267826 Audio Cable Asm

5270723 GOCAA-PDU DC Power Cable

5267389 Run E3045 Modem Cable

5272839 Host Hub Ethernet Cable

5272839-2 Hub-Front Port Ethernet Cable

5264521 Run E3048 USB Extension Cable

5169105 Run E3016 Y6 Cable


5265059 Host to IRDC DVI-D Cable

5265784 Host to IRDC USB Cable

5270246 IRDC Power Cable

5265358 One-Wire Collector FRU consists of:

Humidity & Temperature iButton, iButton Holder, USB to RS422 Converter, USB and Sensor cables, Run E3019, and Dual Locking Tape

SCIM & KEYBOARD:**NOTE:**

If Keyboard is to be replaced, also order **46-317475P2**, 2 inch x 4 inch adhesive backed velcro to attach to bottom of new keyboard.

KEYBOARDS:

- 5143798-7** Standard USB Keyboard – English
- 5143798-8** Standard USB Keyboard – French
- 5143798-9** Standard USB Keyboard – German
- 5143798-10** Standard USB Keyboard – Spanish
- 5143798-11** Standard USB Keyboard – Italian
- 5143798-12** Standard USB Keyboard – Portuguese
- 5143798-13** Standard USB Keyboard – Scandinavian
- 5183547-46** Standard USB Keyboard – Danish
- 5183547-47** Standard USB Keyboard – Dutch
- 5183547-48** Standard USB Keyboard – Norwegian

SCIM OVERLAYS:



- 5148389-2** MR SCIM – Black with English Overlay
- 5148389-9** MR SCIM – Black Coordinated Global Overlay Set

MR750 COILS:






3T GP FLEX COIL	M1385AF
<u>2373366</u> GP Flex Coil Asm	
<u>2379330</u> GP Flex Coil Cable Kit	
<u>2383558</u> GP Flex Coil Hardware Kit	
3T 8 CHANNEL CARDIAC ARRAY COIL	M3334TS
<u>2414705</u> Anterior Coil	
<u>2414706</u> Posterior Coil	
<u>2414707</u> System Cable Asm	
<u>2414708</u> Phantom Positioner Kit	
3T 8 CHANNEL NEURO ARRAY COIL	M3334TP
<u>2414507</u> Staging Collector FRU	
3T 8 CHANNEL BREAST COIL	M3335LA
<u>2417367</u> Breast Coil without Cable	
<u>2417366</u> Output Cable Asm	
<u>2417234</u> Phantom Asm	
<u>2408503</u> Chamber Cover Pad	
<u>2414541</u> Biopsy Mount	
<u>2412285</u> Biopsy Arm	
<u>2412281</u> Breast Coil Headrest	
3T 8 CHANNEL T/R KNEE ARRAY COIL	M3335LB
<u>5160254-2</u> T/R Knee Array with Cable and Baseplate	
<u>5160254-4</u> T/R Knee Array Phantom	
<u>5160254-5</u> T/R Knee Array Cable	
<u>5114356-4</u> T/R Knee Array Baseplate	
<u>5114356-12</u> T/R Knee Array Phantom Positioner	
<u>5114356-13</u> T/R Knee Array Latch Kit	
3T 8 CHANNEL BRAIN SPINE COIL (Forward Production)	M3335LC
<u>5417351</u> Head Neck Unit	
<u>5417355</u> HN System Connector	
3T 8 CHANNEL BRAIN SPINE COIL (Upgrade)	M3335LD
<u>5417352</u> HN System Connector, Hypertronics Version	
3T 8 CHANNEL BRAIN ARRAY – INVIVO	M3335LL
<u>5147134-2</u> HiRes Brain Coil with Mirrors and Baseplate Asm	
<u>5147134-3</u> HiRes Brain Coil Cable	
<u>5146634-5</u> HiRes Brain Coil Inferior Mirror	
<u>5146634-6</u> HiRes Brain Coil Superior Mirror	
<u>2317112-3</u> HiRes Brain Phantom Positioner	

COILS CONTINUED ON NEXT PAGE.

MR750 COILS (Continued):**3T 8 CHANNEL CTL SPINE ARRAY COIL (HD Connector) M3335LM****2417540** Coil Asm**2417542** Cable Asm**2375136-5** Patient Cover Pad**3T 8 CHANNEL HD TORSO ARRAY COIL (HD Connector) M3335LN****2418063** Anterior Coil**2418064** Posterior Coil**2416958** Cable Asm**3T QUAD KNEE/FOOT COIL – INVIVO (HD Connector) M3335LP****5147137-2** Quad Extremity Coil Asm**5147225-3** Quad Extremity Phantom**5147225-6** Quad Extremity Cable Asm**5147225-7** Quad Extremity Phantom Positioner**2357819-10** Extremity Repair Kit**3T 3 CHANNEL SHOULDER COIL (HD Connector) M3335LR****2416877** Cable Asm**2416878** Coil Asm**3T 8 CHANNEL WRIST COIL (HD Connector) M3340CA****5160986-5** Phantom**5160986-12** Baseplate and Baseplate Riser**5268024-2** Coil with Baseplate and Cable**5268024-3** Coil Cable**3T 8 CHANNEL FOOT/ANKLE COIL (HD Connector) M3340CB****5273235-2** Coil with Cable**5273235-3** Coil Cable**5273235-5** Phantom Asm**5273235-6** Coil Baseplate Asm**3T 8 CHANNEL SHOULDER COIL (HD Connector) M3340CC****5271265-2** Coil with Cable**5271265-3** Coil Cable Asm**5271265-5** Phantom Holder Asm

1.5T UNIFIED SURFACE COIL PHANTOM SET (MR450 Only)**5342679**Large Cylindrical
Unified Phantom**5343347**

TL Unified Phantom

**5342680**Small Cylindrical
Unified Phantom**5342681**Cubical
Unified Phantom**COIL STAGING FRU'S****2424308**1.5T 12-Channel BA
Output Staging Cable
(P-Connector)**2424269**1.5T 16-Channel HNS
System Staging Cable
(P-Connector)**2424387**1.5T Phantom Positioner
for Torso Coil**SERVICE KITS****CABLE FRU KIT's****5306523**MR450/MR750 Service Cable Hardware Kit
(Includes BNC components, studs, loopbacks, gender changers)

CABLE COLLECTOR FRU'S**5192499, 3.0T Rear Pedestal Cable FRU Kit contains:**

Run M3320 Run M3321 Run M1315 Run M1335 Run M2300 Run M3351
Run M1324

5308827, 1.5T Rear Pedestal Cable FRU Kit contains:

Run M3320 Run M3321 Run M1315 Run M1335 Run M2300 Run M3351

5192500, Rear Pedestal 32-Channel Add-on Cable FRU Kit contains:

Run M3322 Run M3323 Run M3330 Run M3331 Run M3332 Run M3333
Run M3371 Run M3372 Run M3373

5191899, Scan Room Base SHORT Cable FRU Kit contains:

Run M1302 Run M1303 Run M1304 Run M1310 Run M1311 Run M1312
Run M2025 Run M3030 Run M3031 Run M3300 Run M3301 Run M3302
Run M3303 Run M3304 Run M3306 Run M3307 Run M3308 Run M3309
Run M3310 Run M3311 Run M3312 Run M3313 Run M3315 Run M3316
Run M4005 13W6 Shunt Connector for 16 & 32 Channel Only

5191899-2, Scan Room Base LONG Cable FRU Kit contains:

Run M1302 Run M1303 Run M1304 Run M1310 Run M1311 Run M1312
Run M2025 Run M3030 Run M3031 Run M3300 Run M3301 Run M3302
Run M3303 Run M3304 Run M3306 Run M3307 Run M3308 Run M3309
Run M3310 Run M3311 Run M3312 Run M3313 Run M3315 Run M3316
Run M4005 13W6 Shunt Connector for 16 & 32 Channel Only

5191278, 3.0T Scan Room Transmit SHORT Cable FRU Kit contains:

Run M1305 Run M1306 Run M1001 Run M1002 Run M1006 Run M1015

5191278-2, 3.0T Scan Room Transmit LONG Cable FRU Kit contains:

Run M1305 Run M1306 Run M1001 Run M1002 Run M1006 Run M1015

5191368, 1.5T Scan Room Transmit SHORT Cable FRU Kit contains:

Run M1305 Run M1306 Run M1001 Run M1002 Run M1341 Run M1342

5191368-2, 1.5T Scan Room Transmit LONG Cable FRU Kit contains:

Run M1305 Run M1306 Run M1001 Run M1002 Run M1341 Run M1342


5343119, Fiber Optic Cable FRU Kit contains:

Run P2003 Run P2004

5343120, In-Room Display (IRD) Cable FRU Kit contains:


Run P2026 Run M3366

5169566 Run 3367, MAG-TDM-J1 to MAG-RX1, 44-pin HD sub-D 560mm
5169566-2 Run 3368, MAG-TDM-J2 to MAG-RX2, 44-pin HD sub-D 1360mm
5180907 Run 3350, MAG-SRI-J10 to MAG-BORE, Bore Temp, 3000mm

8KW MNS CABINET

NOTE: The following RF POWER AMP internal parts are **Not Shown**:

- 2206456-2** RF Power Amp Preamp/Splitter
- 2206456-3** RF Amplifier (5)
- 2206456-4** RF Power Amp Combiner/Coupler
- 2206456-5** RF Power Amp Coupler (4)
- 2206456-6** RF Power Amp Bias Tee
- 2206456-13** Cable – Splitter To Amplifier (5)
- 2206456-20** Coupler–bias Tee
- 2206456-21** Power Bias Tee
- 2206456-22** Option 10, PMI Interface Kit
- 2206456-23** Option 11, PMI Interface Kit
- 2206456-25** Option 12, PMI Interface Kit
- 46-170021P29** 5 Amp Fuses (3)

MNS COMPONENTS

5250032
8kW BB Amp I/F BD.
Slots 17/18 U

5250034
RF DETECTOR BD.
Slots 19 U/L

2329754
50 OHM SMB
TERMINATOR

- 5263595** MNS Cable Kit – Internal
5271929 MNS Cable Kit – Interconnects
46-208990P1 50 Ohm BNC Bulkhead Adaptor

- 5300204-2** MR750 3.0T MNS Exciter (located in PEN Cabinet)
5189655-2 MNS Upconverter Module (located in Rear Ped)
2109880-2 Plain Washer, 6.4mm x 12mm

- 2354050-6** 3T 13C TR Switch (located on LPCA)

- 2109880-2** Lock Washer, 4.1mm x 7.6mm (6)

- 2109873-15** Pan Head Screw, M4 x 12mm (6)

- 2109878-2** Plain Washer, 4.3mm x 9mm (6)

- 5114799-2** 3T 13C MNS Preamplifier

- 2359525** 3T MNS QD Box Linear Coils

- 2359525-2** 3T MNS QD Box Quadrature Coils

- 2379500** 2kW MNS QD Box with One BNC

- 2379500-2** 2kW MNS QD Box with Two BNC

- 5114403-3** MNS QA Probe 3T 13C 32.125 MHz


- 5114978-4** MNS QA Probe 31P

- 5114978-3** 13C Sample Vial




2152220

MNS PHANTOM



2386800

MNS SERVICE COIL



2389665

MNS PHANTOM
POSITIONER

MNS

[illegible]

REVISION HISTORY

REV/DATE REASON FOR CHANGE

- | | | |
|---|------------|---|
| 1 | Jun 27, 08 | Initial release |
| 2 | Oct 11, 08 | Reformatted FRU manual into pocket size and added illustrations to help identify parts. |
| 3 | May 15, 09 | Added MR450 to title to reflect addition of MR450 1.5T components to manual. Update PAC part # and add PPG cable (ECO 2069724). Added cable FRU's (PQR 13250844). Updated XRMB FRU to 2144093-12 (PQR13223405). Added OW HP8400 FRU's (PCN's 2038518 & 2059074) |
| 4 | Feb 3, 10 | Added Preamp & Combiner FRU 2377703 (ECO 2062777). Updated RF Hub Backplane to 5159474-32 (ECO 2078155). Updated ODU Dock cable to 5328862 (ECO 2083025). Added new FRU's (ECO 2053514, 2096463). Added Dock Hook, Extension Spring, E-Clip, Hydraulic Reservoir Sponge, and Uplimit Cable to Patient Table FRU's (PQR 13278899 & 13280426). Added Sumitomo F-50 Compressor & flow control adapter (ECO 2059723). Corrected Cronus module board name to XGA (MRIhc46702). Corrected SPW FRU's (MRIhc46065). Added new APS2 bd, 2294300-23 (ECO 2064898). Corrected fuse part number in MDP (MRIhc46829). Added new GOCAA 2395043-2 (ECO 2084000). Added 1.5T Unified Surface Coil Phantom Set (ECO 2081133). Updated Scan Room short and long cable FRU kit part #'s (ECO 2089895). |