

HLA

SCHEDULE B (HIGH USAGE SCANNER)

SITE:

DATE:

SYSTEM ID:

| <u>Procedure</u>                                  | <u>Inspect</u> | <u>Clean</u> | <u>Split</u> |
|---------------------------------------------------|----------------|--------------|--------------|
| <b>GENERATOR</b>                                  |                |              |              |
| Replace MOV on transformer                        | _____          | _____        | _____        |
| Check generator lamps                             | _____          | _____        | _____        |
| Calibrate 80 kw internal bleeder                  | _____          | _____        | _____        |
| Change oil in transformer wells                   | _____          | _____        | _____        |
| Change oil in tube wells                          | _____          | _____        | _____        |
| Clean generator filter/check fans                 | _____          | _____        | _____        |
| Check PDU fans/ clean filters                     | _____          | _____        | _____        |
| <b>SRU</b>                                        |                |              |              |
| Check for proper fan operation of all fans in SRU | _____          | _____        | _____        |
| For older style ST-CP AP, clean AP blower filter  | _____          | _____        | _____        |
| Clean dirt/dust/debris from bottom of SRU cabinet | _____          | _____        | _____        |
| <b>SYSTEM</b>                                     |                |              |              |
| Update Site Log                                   | _____          | _____        | _____        |
| Check PM manual revision                          | _____          | _____        | _____        |

Comments:

---



---



---



---



---



---



---

Field Engineer:

Id#:

Dispatch#: