



Site Name: _____

FE Name: _____

SMR CT Planned Maintenance Checklist

Month: _____

Procedure

	Q1	Q2	Q3	Q4
Table Function (Up/Down, In/Out) Use All Input Pads	☐	☐	☐	☐
Check Scan and Laser Windows	☐	☐	☐	☐
Inspect Patient Images and Talk With Tech About Issues	☐	☐	☐	☐
Monitor System Log	☐	☐	☐	☐
Check Number of Gantry Rotations	☐	☐	☐	☐
Slice / mAs Count	☐	☐	☐	☐
Verify Emergency Stop Buttons (table/console)	☐	☐	☐	☐
Clean Gantry Ventilation Holes and Inspect Cooling Fans	☐	☐	☐	☐
Clean Tube Fan Filter	☐	☐	☐	☐
Clean SDAS/MDAS Filters	☐	☐	☐	☐
Vacuum / Clean Slip Rings	☐	☐	☐	☐
Vacuum / Clean Slip Brushes Inspect for Wear	☐	☐	☐	☐
Inspect the Axial Motor Brushes		☐		☐
Grease Main Bearing				☐
Check Axial Drive Belt for Wear	☐	☐	☐	☐
Vacuum / Clean Table Inspect Cooling Fans	☐	☐	☐	☐
Check Cradle Drive Motor for Wear and Alignment	☐	☐	☐	☐
Inspect / Clean PDU Assembly	☐	☐	☐	☐
Inspect PDU Cooling Fans	☐	☐	☐	☐
General Cleaning	☐	☐	☐	☐
Check for Smooth Rolling Trackball / Mouse CLEAN	☐	☐	☐	☐
Perform Generator Calibration	☐	☐	☐	☐
Run HHS scans	☐	☐	☐	☐
Check X-Ray On Lights	☐	☐	☐	☐
Check Operation of Final Scan Abort Push Buttons	☐	☐	☐	☐
Image Performance QA	☐	☐	☐	☐
FINAL				
Save System State to MOD/DVD	☐	☐	☐	☐
Check Alignment Lights, Store Phantoms, Install Covers, Update SMR PM/FSR Book, Update Customer Communication Book.	☐	☐	☐	☐
Shutdown System Per Site Procedure	☐	☐	☐	☐

N
o
t
e
s