

GE MRI System PM Checklist

Site:		Date:									
System ID #:		Service Report #:									
Service Engineer:											
Procedure	Serial Number for tools	Pass	Cleaned	Adjusted	Replaced	Repaired	Needs Adjustment	Needs Replacement	Needs Repair	Failed	
Safety and Regulatory											
Check Pneumatic Patient Alert System		☐	☐	☐	☐	☐	☐	☐	☐	☐	
Check Emergency Release of Cradle and Patient Transport		☐	NA	NA	NA	☐	NA	NA	☐	☐	
Check Patient Transport Casters, Armboard Set Screws and Bumper Strips		☐	NA	☐	NA	NA	☐	NA	NA	☐	
Check Patient Transport Caster Locks		☐	NA	NA	NA	☐	NA	NA	☐	☐	
Check Cradle Longitudinal Drive Clutch		☐	NA	☐	NA	NA	☐	NA	NA	☐	
Verify Magnet Emergency Rundown Unit (MRU/ERU) is attached to Magnet via Four Step Test		☐	NA	NA	NA	☐	NA	NA	☐	☐	
Image Quality											
Check ISO Center		☐	NA	☐	NA	NA	☐	NA	NA	☐	
DC Offset Calibration		☐	NA	☐	NA	NA	☐	NA	NA	☐	
Check Gradient Calibration: X		☐	NA	☐	NA	NA	☐	NA	NA	☐	
Y		☐	NA	☐	NA	NA	☐	NA	NA	☐	
Z		☐	NA	☐	NA	NA	☐	NA	NA	☐	
Check Shim/Grad Shim		☐	NA	☐	NA	NA	☐	NA	NA	☐	
Check Head and Body SNR		☐	NA	☐	NA	NA	☐	NA	NA	☐	
Other											
Check Patient Blower & Filter		☐	☐	NA	☐	NA	NA	☐	NA	☐	
Check Patient Transport Docking Alignment		☐	NA	☐	NA	NA	☐	NA	NA	☐	
Clean Lightweight Cradle Wheels		NA	☐	NA	NA	NA	NA	NA	NA	☐	
Check Host Computer Fans and Clean Air Intake Grills		☐	NA	NA	☐	NA	NA	☐	NA	☐	
Check Host Computer Power Supply Fan		☐	NA	NA	☐	NA	NA	☐	NA	☐	
Check Operator Workspace Fans And Clean Air Intake		☐	NA	NA	☐	NA	NA	☐	NA	☐	
Set System Clock		☐	NA	☐	NA	NA	☐	NA	NA	☐	
Check and Delete Error/Message Log & Check SaveInfo		☐	NA	NA	NA	NA	NA	NA	NA	☐	

GE MRI System PM Checklist

Procedure	Serial Number for tools	Pass	Cleaned	Adjusted	Replaced	Repaired	Needs Adjustment	Needs Replacement	Needs Repair	Failed
Check RF Cabinet and Clean Filters		☐	☐	NA	☐	NA	NA	☐	NA	☐
Check Fluid Level of Heat Exchanger for Gradient Coil Cooling (MOBILES ONLY)		☐	NA	☐	☐	NA	☐	☐	NA	☐
Check Water Chiller for Gradient Coil Cooling (FIXED SYSTEMS ONLY)		☐	NA	☐	☐	NA	☐	☐	NA	☐
Check Gradient Cabinet and Filters		☐	NA	NA	☐	NA	NA	☐	NA	☐
Check System Cabinet and Filters		☐	NA	NA	☐	NA	NA	☐	NA	☐
Record Cryostat Pressure		☐	NA	NA	NA	☐	☐	NA	NA	☐
Inspect and Record Sumitomo System		☐	NA	NA	NA	☐	☐	NA	NA	☐
Check Cryogen Levels and MM3 Coldhead Temperature, if applicable		☐	NA	NA	NA	☐	☐	NA	NA	☐
Check each coil for artifacts and operation		☐	☐	NA	☐	☐	NA	☐	☐	☐

Notes:

Service Engineer Signature and Date: _____