

Excite II and Newer MR PM Report



Document #: GE-0000-00-PM1

Revision: C

Affected Products

Excite II and newer MR

Reference Material

MR Planned Maintenance Data Sheets GE Excite II and newer

History

Date	Name	Rev	Change Descriptions
3/16/15	Ruben Beltran	A	Initial release.
6/6/17	Jhovany Hernandez	B	Updated to standard, fillable pdf.
9/18/17	Jhovany Hernandez	C	Added NA Field moved signature to the bottom.

Approval

Department	Name	Date
Specialist	Ruben Beltran	
Document Control	Jhovany Hernandez	9/18/17

Preface

This document is valid proof of a Preventative Maintenance performed on the system. The Preventative Maintenance must be performed in the specified intervals, and the results are to be documented herein.

Refer to the **Reference Material listed above** for specific procedures.

This section addresses planned maintenance that should be performed by service personnel. Some of it is also listed in the Operator's Manual for the user to perform on daily and weekly basis.

However, service personnel should make sure that these procedures have been followed and if not, include them on the periodic visits.

Warning: All steps must be performed in accordance with the manufacturer's service manual.

This document must be completely filled out by the Field Service Engineer, i.e.:

- All boxes must be filled out. Use N/A where appropriate.
- Enter the Serial Number and the Date in the header of each page.
- After completing the check, fill out and sign the last page.

--

--

System Information

--

--

--

--

--

Measuring Equipment Used

Enter the measuring equipment used in the table below (meters, scopes, phantoms, etc.):

--

--

--

--

--

Customer: Operator Reported Problems

--

Engineer: FSE Reported Problems

--

System Serial #:

Date:

Procedure.*

*Record values in data table.

A=Bi-monthly
B=6 months

Passed

Cleaned

Adjusted

Repaired

Replaced

Needs Adjustment

Needs Repair

Needs Replacement

Failed

N/A

Safety and Regulatory

Check O2 monitor operation, if applicable.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check peripheral pulse gating cable.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check RF power output.*	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check emergency release of cradle.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check emergency release of patient transport.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check patient transport casters.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check patient armboard setscrews.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Verify that cradle must be at home to lower table.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check cradle longitudinal drive clutch.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check all cradle and carriage wheels.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check emergency off circuits.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check power monitor functional test.	B	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check magnet run down unit (ERU, MRU).	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

Image Quality

Check X gradient calibration.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check Y gradient calibration.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check Z gradient calibration.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check shim.*	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check head SNR.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check body SNR.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check all receive path functionality.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

System Serial #:

Date:

Procedure.*	A=Bi-monthly B=6 months	Passed	Cleaned	Adjusted	Repaired	Replaced	Needs Adjustment	Needs Repair	Needs Replacement	Failed	N/A
*Record values in data table.											

Other (Available during patient scanning)

Check and clean RF cabinet filters, if applicable.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check and clean ACGD/HGD cabinet filters.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check and clean host fans and vents.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check and clean TAC cabinet filters, if applicable.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check and clean Excite/system cabinet filters.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Verify GOC/system clock accuracy.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check computer room air temperature.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

Other (Not available during patient scanning)

Check RF shielding and door seals.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check patient fans/blowers and clean filter.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check gradient coil cooling unit.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Clean patient table.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check dock pedal springs.	B	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Run AGP, APS, SCP diagnostics.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Clean GOC and HP computer fans and air intake.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check and clean GOC mouse.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check and delete error log.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Remove excess "core" files.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check gradient cable connections and support.	B	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check magnet room RF cables and connections.	B	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check patient transport hydraulic fluid level.	B	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Inspect PDU power connection.	B	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check patient alignment lights.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Check ISO Z table movement.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

System Serial #:

Date:

Procedure.*	A=Bi-monthly B=6 months	Passed	Cleaned	Adjusted	Repaired	Replaced	Needs Adjustment	Needs Repair	Needs Replacement	Failed	V/N
*Record values in data table.											

Other (Not available during patient scanning, cont'd)

Record magnet pressure.*	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Record shield cooler hours.*	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Calculate lost hours of shield cooler.*	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Record helium level.*	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Back-up system level information.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Back-up customer protocols.	A	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

Calibration Data

Inhomogeneity:	
X Gradcal:	
Y Gradcal:	
Z Gradcal:	
Head SNR:	
Body SNR:	
Shim Phantom TG:	
DQA Isocenter TG:	

Magnet Data

Magnet Pressure:	
Lhe Level:	
Compressor Pressure:	
Compressor Run Hours:	
Time Run Hours Logged:	
Lost Hours:	

Note: If system is Twinspeed record whole and zoom values.

System Serial #:

Date:

System Condition

After performing scheduled PM, an evaluation of the system is performed with the following result (check one):

No deficiencies. The image quality test resulted in no differences from required reference values.

Slight deficiencies. There is no effect on operation; however they should be corrected preventively. The image quality test resulted in no differences from required reference values.

Serious deficiencies. This must be corrected before operation of the system is permitted.

Notes (Additional work to be performed, partial PM, customer requests):

Field Service Engineer

Signature Field