

**ADMINISTRATIVE INFORMATION****FMI No: 25904****FMI Title: HiSpeed QX/i Table Rating Plate Correction**

FMI TYPE	CHARGE CODES	IMPLEMENTATION DATES
☒ Mandatory Safety ☐ Mandatory ☐ Mandatory on Request	GEMS-AM: GEMS-E: GEMS-A:	Issue Date: Close Date:

☒ Staged FMI	Customer List provided by FMI Administration. If you identify additional units within effectivity range, but they are not shown on the customer list, contact FMI Administration in your Pole.
☒	Parts and any special tools are furnished automatically. To address shortages and ordering exceptions, contact FMI Administration in your Pole.
☐	Parts Ordering: GEMS - Americas, contact your Region logistics coordinator GEMS - Europe, Order FMI kits from GPO/EPL Buc GEMS - Asia, Contact your region logistics coordinator
☐	No Parts are required to complete this FMI.
☐ Batch with	FMI # To minimize costs, combine implementation of this FMI with the FMIs listed for all customers where the effectivity is the same (Please see your customer lists).
☐ Non-Staged	No Customer List is provided. Each service location must identify affected units and place part orders as indicated in the Special Instructions.

Technical Support

GEMS - Americas	GEMS - Europe	GEMS - Asia
GE Medical Systems Insite Center (414) 524-5300 800-321-7973 Milwaukee, WI 53201, USA	GEMS - E / ESC BP34 F 78533 BUC CEDEX, FRANCE Tel: (33 1) 39 20 00 07 or ESC Fax: (33 1) 30 70 99 70	Asia Support Center Phone: 81-426-56-0033 Fax: 81-426-48-2905

Special Instructions:**Unit Tracked:****FMI Instructions Part No.: 2402446-100****Do NOT leave unsecured on site.**



Completion Certification Report

Mail To

GEMS - Americas

GE Medical Systems
Attn: FMI Admin. W-523
P.O. Box 414
Milwaukee, WI 53201

GEMS - Europe

GEMS - E / ESC
FMI Admin. 322
BP34
F 78533 Buc Cedex, France

Send to your
Country FMI or Service
Administrator

FMI No.: 25904 FMI Title: HiSpeed QX/i Table Rating Plate Correction

Customer Full Name:

Number and Street Address:

City and Country:

Country Code

Site Number

System I.D.

Unit Tracked:

Model Number	Serial Number	Compl. Date DD/MM/YY	FMI Comp. Code	Service Hours	
				Labor	Travel

FMI Completion Code:

- 1.) Modification Complete or previously modified.
- 2.) FMI not required per FMI instruction.
- 3.) Not modified. Unit in storage. FMI kit returned to stock.
- 4.) Modification refused or equipment altered. Customer signature and title required.
- 5.) Unit location unknown or unit not in area. Product Locator notified.

**Customer Refusal or Equipment Altered
- Modification Pending.**

Installer
Comments

Effectivity

HiSpeed QX/i systems with Table (Model number: 2269165-3) with the following serial numbers.
(In the list, SN: Serial Number; Date: Manufactured Date; Serial numbers are listed in the ascending order, from top to bottom and from left to right.)

Note: When you spot the serial number of your site table, note the manufactured date.

SN	Date	SN	Date	SN	Date
38240HM8	June-02	44763HM1	October-02	52225HM0	January-03
38241HM6	June-02	44764HM9	October-02	52226HM8	January-03
38242HM4	June-02	44765HM6	October-02	52227HM6	January-03
38243HM2	June-02	45975HM0	October-02	52228HM4	January-03
38244HM0	June-02	45976HM8	October-02	53054HM3	January-03
38245HM7	June-02	45977HM6	October-02	53055HM0	January-03
38246HM5	June-02	45978HM4	October-02	53056HM8	January-03
38247HM3	June-02	46696HM1	October-02	53057HM6	January-03
38248HM1	June-02	46697HM9	November-02	53177HM2	January-03
38249HM9	June-02	46698HM7	November-02	53178HM0	January-03
40007HM7	July-02	46699HM5	November-02	53179HM8	January-03
40008HM5	July-02	46700HM1	November-02	53180HM6	January-03
40009HM3	July-02	46701HM9	December-02	53829HM8	February-03
40010HM1	July-02	46702HM7	December-02	53830HM6	February-03
40011HM9	July-02	46703HM5	December-02	53831HM4	February-03
40012HM7	July-02	46704HM3	December-02	53832HM2	February-03
40013HM5	July-02	50735HM0	December-02	54510HM3	February-03
40014HM3	July-02	50736HM8	December-02	54511HM1	February-03
42660HM1	September-02	50737HM6	December-02	54512HM9	February-03
42661HM9	September-02	50738HM4	December-02	54513HM7	March-03
42662HM7	September-02	50739HM2	December-02	54514HM5	March-03
42663HM5	September-02	50740HM0	December-02	54515HM2	February-03
44754HM0	September-02	50741HM8	December-02	54516HM0	March-03
44755HM7	September-02	50742HM6	December-02	54517HM8	March-03
44756HM5	September-02	51305HM1	December-02	54518HM6	March-03
44757HM3	September-02	51306HM9	December-02	55137HM4	March-03
44758HM1	September-02	51307HM7	December-02	55138HM2	March-03
44759HM9	September-02	51308HM5	December-02	55139HM0	March-03
44760HM7	September-02	51309HM3	December-02	55140HM8	March-03
44761HM5	September-02	52221HM9	December-02	55141HM6	March-03
44762HM3	September-02	52222HM7	December-02	55334HM7	March-03
		52223HM5	December-02	55335HM4	March-03
		52224HM3	December-02	55336HM2	March-03
				55337HM0	March-03

(Continued on next page)

SN	Date	SN	Date	SN	Date
55338HM8	April-03	62542HM6	July-03	65464HM0	September-03
55339HM6	April-03	62543HM4	July-03	65465HM7	September-03
55340HM4	April-03	62544HM2	July-03	65466HM5	September-03
55341HM2	April-03	62545HM9	July-03	65467HM3	September-03
55342HM0	April-03	62546HM7	July-03	66196HM7	September-03
55343HM8	April-03	62547HM5	July-03	66197HM5	September-03
56187HM8	April-03	63460HM0	July-03	66198HM3	September-03
56188HM6	April-03	63461HM8	July-03	66199HM1	September-03
56189HM4	April-03	63462HM6	July-03	66200HM7	September-03
56190HM2	April-03	63463HM4	August-03	66201HM5	September-03
57524HM1	April-03	63464HM2	August-03	66202HM3	September-03
57525HM8	April-03	63465HM9	August-03	66203HM1	September-03
57526HM6	April-03	64221HM5	August-03	66204HM9	September-03
57527HM4	April-03	64222HM3	August-03	66205HM6	September-03
57994HM6	April-03	64222HM4	August-03	67366HM5	September-03
57995HM3	April-03	64223HM1	August-03	67367HM3	September-03
57996HM1	April-03	64224HM9	August-03	67368HM1	September-03
58532HM3	May-03	64619HM0	August-03	67369HM9	September-03
58533HM1	May-03	64620HM8	August-03	67370HM7	September-03
58534HM9	May-03	64621HM6	August-03	67371HM5	September-03
58535HM6	May-03	64623HM2	August-03	67372HM3	September-03
58554HM7	May-03	64624HM0	August-03	67373HM1	September-03
58555HM4	May-03	64625HM7	August-03	68519HM7	October-03
59546HM2	June-03	64626HM5	August-03		
59547HM0	June-03	65460HM8	August-03		
59548HM8	June-03	65461HM6	August-03		
60184HM9	June-03	65462HM4	August-03		
60185HM6	June-03	65463HM2	August-03		
60186HM4	June-03				
60851HM3	June-03				
60852HM1	June-03				
60853HM9	June-03				
60854HM7	June-03				
61272HM1	June-03				
61273HM9	June-03				
61274HM7	June-03				
61275HM4	June-03				

Purpose

The purpose of this FMI is to correct the rating plate attached on the table.

It was found that during a certain period of time tables with rating plates that the manufactured date of the table was not printed on were shipped from factory with HiSpeed QX/i systems. (FDA non-compliant issue)

In this FMI, you will attach a label which the manufactured date is printed on to a location on the rating plate where the date should be shown.

Prerequisite FMIs

None.

Related FMIs

None.

Furnished Materials

FMI 25904 Kit (Part Number 2404459) includes the following items:

Qty	Description	Part #
1	Label Sheet	2404023
1	FMI 25904 Document	2402446-100

Required Tools

None.

People Power

0.5 hours labor and 1 hour of travel

Procedure

1 Attach the Label

- 1.) Confirm that the rating plate on the table does not have the manufactured date printed.

See [Figure 1](#) and [Figure 2](#).




Figure 1 - Rating Plate on Table




Figure 2 - Rating Plate with NO Manufactured Date (Example)

- 2.) Peel a label with the correct manufactured date off the sheet (2404023) provided, and attach it to the location on the rating plate. See [Figure 3](#).



Figure 3 - Rating Plate with Manufactured Date (Example)

2 Finishing Up

Close FMI per regional procedure.

Remove FMI Document from customer site.

