

SECTION 8 – SHIPPING AND DELIVERY DATA

TABLE OF CONTENTS

<u>SECTION</u>	<u>TITLE</u>	<u>PAGE</u>
8-1	INTRODUCTION	8-3
8-2	UPGRADE CATALOGS SHIPPING INFORMATION	8-5

8-1 INTRODUCTION

Domestic transportation for the Signa Horizon upgrade will be via an air-ride moving van. Export transportation for overseas will be via air shipment in a pressurized cargo hold.

The tables in Section 8-2 UPGRADE CATALOGS SHIPPING INFORMATION (arranged by catalog numbers) provide estimates of shipping weights and dimensions for upgrade catalog major components. Refer to the Table 8-1 for page number location for a specific catalog.

Note

The weight and dimension information is estimate only, actual shipment packaging configurations for a catalog may vary.

TABLE 8-1
LIST OF CATALOG TABLES AND PAGE NUMBER LOCATION

TABLE	TITLE	PAGE
8-2	CATALOG M1000MW, TABLE FOR USE WITH ALL MONITORS	8-5
8-3	CATALOG M1000NT, LCD COLOR MONITOR	8-5
8-4	CATALOG M1033K, INDIGO TO OCTANE UPGRADE	8-5
8-5	CATALOG M1033MD, LX MNS KIT	8-5
8-6	CATALOG M1060JW, SHIELD/CRYO COOLER COMPRESSOR	8-6
8-7	CATALOG M1090RW, 512 MBYTE BAM	8-6
8-8	CATALOG M1090TL, PATIENT TABLE	8-6
8-9	CATALOG M1090WG, FAST RECEIVER KIT	8-6
8-10	CATALOG M1090WJ, RECEIVER "0" KIT	8-6
8-11	CATALOG M1090WL, LOW PROFILE CARRIAGE COVER KIT	8-6
8-12	CATALOG M3000RA, REFLEX 100 WITH 512 BAM	8-7
8-13	CATALOG M3000TA, TWINSPEED UPGRADE SYSTEM ELECTRONICS	8-7
8-14	CATALOG M3043TA, SITE COLLECTOR KIT	8-8
8-15	CATALOG M3044TA, SITE COLLECTOR KIT	8-9
8-16	CATALOG M3060TA, 1.5T MAGNET WITH TWINSPEED ENCLOSURE	8-10
8-17	CATALOG M3063TA, SITE COLLECTOR KIT	8-10
8-18	CATALOG M3080TA, 2ND ORDER RESISTIVE SHIM KIT	8-10
8-19	CATALOG M3088TA, OUTDOOR AIR COOLED (DUAL LOOP) TSCC	8-11
8-20	CATALOG M3088TE, INDOOR WATER COOLED (SINGLE LOOP) TGWC	8-11
(Continued)		

8-1 INTRODUCTION (Continued)

TABLE 8-1 (Continued)
LIST OF CATALOG TABLES AND PAGE NUMBER LOCATION

TABLE	TITLE	PAGE
8-21	CATALOG M3088TF, OUTDOOR AIR COOLED (SINGLE LOOP) TGWC	8-11
8-22	CATALOG M3088TG, INDOOR AIR COOLED (SINGLE LOOP) TGWC	8-11
8-23	CATALOG M3090AH, HORIZON TO TWINSPEED UPGRADE SYSTEM ELECTRONICS	8-12
8-24	CATALOG M3090BB, ACGD UPGRADE FOR SYSTEM WITH SGD GRADIENT CABINET & RF/PDU CABINET	8-12
8-25	CATALOG M3090BC, 1.5T ACGD UPGRADE FOR SYSTEM WITH SGD GRADIENT CABINET & RF/PENETRATION CABINET	8-12
8-26	CATALOG M3090BE, ECHOSPEED PLUS PERFORMANCE OPTION	8-13
8-27	CATALOG M3090DB, 5.5 MNS UPGRADE KIT	8-13
8-28	CATALOG M3090DC, LXSPECTRO PRE-RF/PDU TO AGCD MNS UPGRADE KIT	8-13
8-29	CATALOG M3090JA, IPG2 BOARD	8-13
8-23	CATALOG M3090AH, HORIZON TO TWINSPEED UPGRADE SYSTEM ELECTRONICS	8-12
8-30	CATALOG M3090TA, TWIN UPGRADE FOR LCC HORIZON FIXED & T/R	8-14
8-31	CATALOG M1085BG, 1.5T BODY COIL FOR SIGNA HORIZON	8-14

8-2 UPGRADE CATALOGS SHIPPING INFORMATION

TABLE 8-2
CATALOG M1000MW, TABLE FOR USE WITH ALL MONITORS

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Operator Workspace Table for Color Monitor	45 x 54 x 37 (1143 x 1372 x 940)	180 (82)	box

TABLE 8-3
CATALOG M1000NT, LCD COLOR MONITOR

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Operator Workspace Color Monitor	27 x 33 x 27 (686 x 838 x 686)	125 (57)	skid

TABLE 8-4
CATALOG M1033K, INDIGO TO OCTANE UPGRADE

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Upgrade Kit	20 x 20 x 20 (508 x 508 x 508)	37 (17)	box
System Computer	25 x 38 x 34 (635 x 965 x 864)	200 (91)	skid
Octane Computer	22 x 30 x 24 (559 x 762 x 610)	50 (23)	box on skid

TABLE 8-5
CATALOG M1033MD, LX MNS KIT

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Multi Nuclear Spectro MOD Option Key for 8X/LX	12 x 12 x 4 (305 x 305 x 102)	3 (1.5)	box
1.5T MNS Phosphorus Flex Coil Kit	18 x 24 x 6 (457 x 610 x 152)	5 (2.5)	box
MRS Phantom	16 x 16 x 16 (406 x 406 x 406)	10 (5)	box
LX MNSpectro Magnet Enclosure Kit for HI/LO Carriage Covers	24 x 24 x 24 (610 x 610 x 610)	20 (9)	box
UCERD to TNS and I/F Panel Cable	16 x 16 x 6 (406 x 406 x 152)	8 (4)	box
Heliac Cable	36 x 36 x 4 (914 x 914 x 102)	20 (9)	box
Heliac Cable	36 x 36 x 4 (914 x 914 x 102)	20 (9)	box
MNS Cabinet	24 x 42 x 71 (610 x 1067 x 1803)	430 (195)	on cabinet casters, wrapped with plastic
MNS Interconnect Cables	16 x 16 x 6 (406 x 406 x 152)	8 (4)	box

8-2 UPGRADE CATALOGS SHIPPING INFORMATION (Continued)

TABLE 8-6
CATALOG M1060JW, SHIELD/CRYO COOLER COMPRESSOR

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Shield/Cryo Cooler Compressor Cabinet	26 x 28 x 42 (660 x 711 x 1067)	240 (109)	skid with box cover

TABLE 8-7
CATALOG M1090RW, 512 MBYTE BAM

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Memory Board	16 x 16 x 5 (406 x 406 x 127)	5 (2)	box

TABLE 8-8
CATALOG M1090TL, PATIENT TABLE

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Patient Table	95 x 36 x 50 (2413 x 914 x 1270)	670 (305)	pallet

TABLE 8-9
CATALOG M1090WG, FAST RECEIVER KIT

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Fast Receiver Kit	26 x 26 x 12 (660 x 660 x 305)	35 (16)	box
Technical Publication & Software	12 x 12 x 4.5 (305 x 305 x 114)	10 (5)	box

TABLE 8-10
CATALOG M1090WJ, RECEIVER "0" KIT

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Receiver 0 Kit	26 x 26 x 4 (660 x 660 x 102)	25 (11)	box

TABLE 8-11
CATALOG M1090WL, LOW PROFILE CARRIAGE COVER KIT

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Low Profile Carriage Cover	18 x 30 x 14 (457 x 762 x 356)	50 (23)	box

8-2 UPGRADE CATALOGS SHIPPING INFORMATION (Continued)

TABLE 8-12
CATALOG M3000RA, REFLEX 100 WITH 512 BAM

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
ReFlex 100 with 512 BAM	16 x 16 x 5 (406 x 406 x 127)	5 (2)	box

TABLE 8-13
CATALOG M3000TA, TWINSPEED UPGRADE SYSTEM ELECTRONICS

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Patient Table	95 x 36 x 50 (2413 x 914 x 1270)	670 (305)	pallet
Ship Collector	18 x 30 x 24 (457 x 762 x 610)	50 (23)	box
Non-magnetic chair	36 x 24 x 38 (914 x 610 x 965)	50 (23)	box
Octane Computer	24 x 36 x 23 (610 x 914 x 584)	30 (14)	box
Head Coil Assembly	24 x 24 x 24 (610 x 610 x 610)	35 (16)	box
Phantom Set	16 x 16 x 16 (406 x 406 x 406)	20 (9)	box
Operator Workspace Equipment	24 x 36 x 23 (610 x 914 x 584)	60 (27)	box
Operator Workspace Cabinet	25 x 38 x 34 (635 x 965 x 864)	200 (91)	skid
ACGD/PDU Cabinet without Fan Module installed	24 x 37 x 75 (610 x 940 x 1905)	1810 (823)	on cabinet casters, wrapped with plastic
Fan Module for ACGD/PDU Cabinet	28 x 38 x 15 (711 x 965 x 381)	160 (73)	on pallet
SRF Cabinet	24 x 42 x 71 (610 x 1067 x 1803)	505 (230)	on cabinet casters, wrapped with plastic
System Control Cabinet	24 x 33 x 78 (610 x 838 x 1981)	494 (225)	on cabinet casters, wrapped with plastic
Twin Accessory Cabinet	23 x 39 x 50 (584 x 991 x 1270)	600 (272)	on cabinet casters, wrapped with plastic

8-2 UPGRADE CATALOGS SHIPPING INFORMATION (Continued)

TABLE 8-14
CATALOG M3043TA, SITE COLLECTOR KIT

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Technical Publications	12 x 20 x 12 (305 x 508 x 305)	40 (18)	box
Penetration Panel	22 x 40 x 25 (559 x 1016 x 635)	70 (32)	box
SPT Phantom Set	34 x 32.5 x 60 (864 x 826 x 1524)	350 (159)	on cart casters with box cover
Head & Body Transmit Cable Kit	36 x 36 x 6 (914 x 914 x 152)	20 (9)	box
Power Cables Kit	17 x 26 x 18 (432 x 660 x 457)	100 (45)	box
Signal Cables Kit	40 x 40 x 46 (1016 x 1016 x 1168)	100 (45)	4 boxes on skid
Fiber Optic Cables Kit	27.5 x 31 x 33 (699 x 787 x 838)	50 (23)	box
Gradient Cables	40 x 40 x 18 (1016 x 1016 x 457)	620 (282)	box on skid
Gradient Hose Assembly	24 x 12 x 12 (610 x 305 x 305)	30 (14)	box
Vacuum Hose Assembly	24 x 12 x 12 (610 x 305 x 305)	30 (14)	box
Cover Kit	30 x 66 x 7 (762 x 1676 x 178)	50 (23)	box
Cryo Monitor Support Bracket Kit	30 x 30 x 5 (762 x 762 x 127)	20 (9)	box
Left & Right Trim Cover Assembly	37 x 52 x 24 (940 x 1321 x 610)	125 (57)	box on skid
Blower Assembly	21 x 15 x 27 (533 x 381 x 686)	50 (23)	box
Pedstal Cover Kit	30 x 45 x 12 (762 x 1143 x 305)	50 (23)	box
Top Front & Top Rear Cover Kits	55 x 45 x 40 (1397 x 1143 x 1016)	150 (68)	2 boxes strapped to skid
Blower Hose Kit	44 x 44 x 10 (1118 x 1118 x 254)	25 (11)	box
Rear Pedestal Assembly	23 x 49 x 36 (584 x 1245 x 914)	100 (45)	box on skid
PAC/SRI Platform Assembly, PAC I/F Assembly	18 x 45 x 12 (457 x 1143 x 305)	40 (18)	box

8-2 UPGRADE CATALOGS SHIPPING INFORMATION (Continued)

TABLE 8-15
CATALOG M3044TA, SITE COLLECTOR KIT

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Technical Publications	12 x 20 x 12 (305 x 508 x 305)	40 (18)	box
Penetration Panel	22 x 40 x 25 (559 x 1016 x 635)	70 (32)	box
SPT Phantom Set	34 x 32.5 x 60 (864 x 826 x 1524)	350 (159)	on cart casters with box cover
Head & Body Transmit Cable Kit	36 x 36 x 6 (914 x 914 x 152)	20 (9)	box
Power Cables Kit	17 x 26 x 18 (432 x 660 x 457)	100 (45)	box
Signal Cables Kit	40 x 40 x 46 (1016 x 1016 x 1168)	100 (45)	4 boxes on skid
Fiber Optic Cables Kit	27.5 x 31 x 33 (699 x 787 x 838)	50 (23)	box
Gradient Cables	40 x 40 x 18 (1016 x 1016 x 457)	620 (282)	box on skid
Gradient Hose Assembly	24 x 12 x 12 (610 x 305 x 305)	30 (14)	box
Vacuum Hose Assembly	24 x 12 x 12 (610 x 305 x 305)	30 (14)	box
Cover Kit	30 x 66 x 7 (762 x 1676 x 178)	50 (23)	box
Cryo Monitor Support Bracket Kit	30 x 30 x 5 (762 x 762 x 127)	20 (9)	box
Left & Right Trim Cover Assembly	37 x 52 x 24 (940 x 1321 x 610)	125 (57)	box on skid
Blower Assembly	21 x 15 x 27 (533 x 381 x 686)	50 (23)	box
Pedstal Cover Kit	30 x 45 x 12 (762 x 1143 x 305)	50 (23)	box
Top Front & Top Rear Cover Kits	55 x 45 x 40 (1397 x 1143 x 1016)	150 (68)	2 boxes strapped to skid
Blower Hose Kit	44 x 44 x 10 (1118 x 1118 x 254)	25 (11)	box
Rear Pedestal Assembly	23 x 49 x 36 (584 x 1245 x 914)	100 (45)	box on skid
PAC/SRI Platform Assembly, PAC I/F Assembly	18 x 45 x 12 (457 x 1143 x 305)	40 (18)	box

8-2 UPGRADE CATALOGS SHIPPING INFORMATION (Continued)

TABLE 8-16
CATALOG M3060TA, 1.5T MAGNET WITH TWINSPEED ENCLOSURE

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
LCC Magnet with cryogenics, partial Quiet Technology Enclosure installed	93 x 144 x 107 (2362 x 3658 x 2718)	See Note 1	Domestic – Tarped International – crate/pallet
Magnet Accessory Equipment	48 x 48 x 28 (1219 x 1219 x 711)	400 (182)	crate
1 Approximate magnet shipping weight (includes packaging material) and weight of magnet with cryogenics, TRM Gradient & RF coils, Enclosure parts installed on magnet, lifting beams, crate, and skid (i.e. minus packaging material): 12,195 lbs (5532 kg). International shipments must add shipping crate/pallet of 2,200 lbs (998 kg).			

TABLE 8-17
CATALOG M3063TA, T/R SITE COLLECTOR KIT

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Technical Publications	12 x 20 x 12 (305 x 508 x 305)	40 (18)	box
Penetration Panel	22 x 40 x 25 (559 x 1016 x 635)	70 (32)	box
SPT Phantom Set	34 x 32.5 x 60 (864 x 826 x 1524)	350 (159)	on cart casters with box cover
Head & Body Transmit Cable Kit	36 x 36 x 6 (914 x 914 x 152)	20 (9)	box
Power Cables Kit	17 x 26 x 18 (432 x 660 x 457)	100 (45)	box
Signal Cables Kit	40 x 40 x 46 (1016 x 1016 x 1168)	100 (45)	4 boxes on skid
Fiber Optic Cables Kit	27.5 x 31 x 33 (699 x 787 x 838)	50 (23)	box
Gradient Cables	40 x 40 x 18 (1016 x 1016 x 457)	620 (282)	box on skid
Gradient Hose Assembly	24 x 12 x 12 (610 x 305 x 305)	30 (14)	box
Vacuum Hose Assembly	24 x 12 x 12 (610 x 305 x 305)	30 (14)	box
Cover Kit	30 x 66 x 7 (762 x 1676 x 178)	50 (23)	box
Cryo Monitor Support Bracket Kit	30 x 30 x 5 (762 x 762 x 127)	20 (9)	box
Left & Right Trim Cover Assembly	37 x 52 x 24 (940 x 1321 x 610)	125 (57)	box on skid
Blower Assembly	21 x 15 x 27 (533 x 381 x 686)	50 (23)	box
Pedestal Cover Kit	30 x 45 x 12 (762 x 1143 x 305)	50 (23)	box
Top Front & Top Rear Cover Kits	55 x 45 x 40 (1397 x 1143 x 1016)	150 (68)	2 boxes strapped to skid
Blower Hose Kit	44 x 44 x 10 (1118 x 1118 x 254)	25 (11)	box
Rear Pedestal Assembly	23 x 49 x 36 (584 x 1245 x 914)	100 (45)	box on skid
PAC/SRI Platform Assembly, PAC I/F Assembly	18 x 45 x 12 (457 x 1143 x 305)	40 (18)	box

TABLE 8-18
CATALOG M3080TA, 2ND ORDER RESISTIVE SHIM KIT

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
2nd Order Resistive Shim Kit	36 x 24 x 18 (914 x 610 x 457)	50 (23)	box

8-2 UPGRADE CATALOGS SHIPPING INFORMATION (Continued)

TABLE 8-19
CATALOG M3088TA, OUTDOOR AIR COOLED (DUAL LOOP) TSCC

MR COMPONENT	APPROXIMATE W x D x H	APPROXIMATE WEIGHT	METHOD OF SHIPMENT
TSCC unit	52 x 82 x 75 (1321 x 2083 x 1905)	1320 (599)	crate/skid
Remote Control Panel	20 x 10 x 20 (51 x 25 x 51)	33 (15)	box
DowFrost HD glycol	12 x 12 x 18 (305 x 305 x 457)	50 (23)	4 boxes
Ship loose items	32 x 32 x 32 (813 x 813 x 813)	50 (23)	box

TABLE 8-20
CATALOG M3088TE, INDOOR WATER COOLED (SINGLE LOOP) TGWC

MR COMPONENT	APPROXIMATE W x D x H	APPROXIMATE WEIGHT	METHOD OF SHIPMENT
TGWC Indoor Water Cooled	38 x 24 x 27 (97 x 61 x 69)	200 (91)	box on pallet
Remote Control Panel	12 x 12 x 12 (305 x 305 x 305)	20 (9)	box
DowFrost HD glycol	12 x 12 x 18 (305 x 305 x 457)	50 (23)	2 boxes
Ship loose items	12 x 12 x 12 (305 x 305 x 305)	30 (14)	box
Ship loose items	26 x 6 x 6 (660 x 152 x 152)	15 (7)	box

TABLE 8-21
CATALOG M3088TF, OUTDOOR AIR COOLED (SINGLE LOOP) TGWC

MR COMPONENT	APPROXIMATE W x D x H	APPROXIMATE WEIGHT	METHOD OF SHIPMENT
TGWC unit	55 x 40 x 82 (141 x 102 x 208)	918 (417)	crate/skid
Remote Control Panel	13 x 7 x 15 (330 x 16 x 37)	20 (9)	box
DowFrost HD glycol & ship loose items	24 x 24 x 18 (610 x 610 x 457)	50 (23)	2 boxes

TABLE 8-22
CATALOG M3088TG, INDOOR AIR COOLED (SINGLE LOOP) TGWC

MR COMPONENT	APPROXIMATE W x D x H	APPROXIMATE WEIGHT	METHOD OF SHIPMENT
TGWC Indoor Air Cooled	34 x 40 x 76 (59 x 76 x 186)	918 (417)	crate/skid
DowFrost HD glycol & ship loose items	24 x 24 x 12 (610 x 610 x 305)	50 (23)	box

8-2 UPGRADE CATALOGS SHIPPING INFORMATION (Continued)

TABLE 8-23
CATALOG M3090AH, HORIZON TO TWINSPEED UPGRADE SYSTEM ELECTRONICS

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Ship Collector	18 x 30 x 24 (457 x 762 x 610)	50 (23)	box
Non-magnetic chair	36 x 24 x 38 (914 x 610 x 965)	50 (23)	box
Octane Computer	24 x 36 x 23 (610 x 914 x 584)	30 (14)	box
System Control Cabinet	24 x 33 x 78 (610 x 838 x 1981)	494 (225)	on cabinet casters, wrapped with plastic
Operator Workspace Equipment	24 x 36 x 23 (610 x 914 x 584)	60 (27)	box
Operator Workspace Cabinet	25 x 38 x 34 (635 x 965 x 864)	200 (91)	skid
ACGD/PDU Cabinet without Fan Module installed	24 x 37 x 75 (610 x 940 x 1905)	1810 (823)	on cabinet casters, wrapped with plastic
Fan Module for ACGD/PDU Cabinet	28 x 38 15 (711 x 965 x 381)	160 (73)	on pallet
SRF Cabinet	24 x 42 x 71 (610 x 1067 x 1803)	505 (230)	on cabinet casters, wrapped with plastic

TABLE 8-24
CATALOG M3090BB, ACGD UPGRADE FOR SYSTEM WITH SGD GRADIENT CABINET & RF/PDU CABINET

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Power Cables Kit	24 x 24 x 18 (610 x 610 x 457)	50 (23)	box
ACGD/PDU Cabinet without Fan Module installed	24 x 37 x 75 (610 x 940 x 1905)	1810 (823)	on cabinet casters, wrapped with plastic
Fan Module for ACGD/PDU Cabinet	28 x 38 15 (711 x 965 x 381)	160 (73)	on pallet
SRF Cabinet	24 x 42 x 71 (610 x 1067 x 1803)	505 (230)	on cabinet casters, wrapped with plastic

TABLE 8-25
CATALOG M3090BC, 1.5T ACGD UPGRADE FOR SYSTEM WITH SGD GRADIENT CABINET & RF/PENETRATION CABINET

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Power Cables Kit	24 x 24 x 18 (610 x 610 x 457)	50 (23)	box
Dual Sensor Upgrade Kit	14 x 10 x 10 (356 x 254 x 254)	20 (9)	box
ACGD/PDU Cabinet without Fan Module installed	24 x 37 x 75 (610 x 940 x 1905)	1810 (823)	on cabinet casters, wrapped with plastic
Fan Module for ACGD/PDU Cabinet	28 x 38 15 (711 x 965 x 381)	160 (73)	on pallet
SRF Cabinet	24 x 42 x 71 (610 x 1067 x 1803)	505 (230)	on cabinet casters, wrapped with plastic
Gradient Cables	40 x 40 x 18 (1016 x 1016 x 457)	620 (282)	box on skid

8-2 UPGRADE CATALOGS SHIPPING INFORMATION (Continued)

TABLE 8-26
CATALOG M3090BE, ECHOSPEED PLUS PERFORMANCE OPTION

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Software & Install document	9 x 12 x 1 (229 x 305 x 25)	2 (1)	envelope

TABLE 8-27
CATALOG M3090DB, 5.5 MNS UPGRADE KIT

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Multi Nuclear Spectro MOD Option Key	12 x 12 x 4 (305 x 305 x 102)	3 (1.5)	box
1.5T RF/Pen Style Cabinet MNSpectro Upgrade Kit	24 x 24 x 18 (610 x 610 x 457)	30 (14)	box
1.5T MNS Phosphorus Flex Coil Kit	18 x 24 x 6 (457 x 610 x 152)	5 (2.5)	box
MRS Phantom	16 x 16 x 16 (406 x 406 x 406)	10 (5)	box
UCERD to TNS and I/F Panel Cable	16 x 16 x 6 (406 x 406 x 152)	8 (4)	box
MNS Cabinet	24 x 42 x 71 (610 x 1067 x 1803)	430 (195)	on cabinet casters, wrapped with plastic
MNS Interconnect Cables	16 x 16 x 6 (406 x 406 x 152)	8 (4)	box

TABLE 8-28
CATALOG M3090DC, LXSPECTRO PRE-RF/PDU TO AGCD MNS UPGRADE KIT

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
MNS Cabinet	24 x 42 x 71 (610 x 1067 x 1803)	430 (195)	on cabinet casters, wrapped with plastic
MNS Interconnect Cables	16 x 16 x 6 (406 x 406 x 152)	8 (4)	box

TABLE 8-29
CATALOG M3090JA, IPG2 BOARD

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
IPG2 Board	16 x 16 x 5 (406 x 406 x 127)	5 (2)	box

8-2 UPGRADE CATALOGS SHIPPING INFORMATION (Continued)

TABLE 8-30
CATALOG M3090TA, TWIN UPGRADE FOR LCC HORIZON FIXED & T/R

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Technical Publications	12 x 21 x 12 (305 x 533 x 305)	20 (9)	box
Left/right Trim Cover Assembly	37 x 52 x 22 (940 x 1321 x 102)	125 (57)	box on skid
Rear Pedestal with Split Bridge Rear Upgrade Assembly	23 x 49 x 32 (584 x 1245 x 813)	65 (30)	box
Pedstal Cover Kit	30 x 45 x 12 (762 x 1143 x 305)	50 (23)	box
Bridge	24 x 161 x 10 (610 x 4090 x 254)	125 (57)	box
Dock Control Box, Bore Light Assembly, Cable	18 x 45 x 11 (457 x 1143 x 279)	30 (14)	box
Top Front & Top Rear Cover Kits	55 x 45 x 40 (1397 x 1143 x 1016)	150 (68)	2 boxes strapped to skid
Dock Cover Kit	30 x 14 x 8 (762 x 356 x 203)	50 (23)	box
PAC/SRI Platform Assembly, PAC I/F Assembly	18 x 45 x 12 (457 x 1143 x 305)	40 (18)	box
Water Chiller power cable/circuit breaker	16 x 16 x 16 (406 x 406 x 406)	50 (23)	box
Front Cover Cart Assembly for return of old cover/endbell	36 x 106 x 57 (914 x 2693 x 1448)	100 (45)	cart
Rear Cover Cart Assembly for return of old cover/endbell	36 x 106 x 57 (914 x 2693 x 1448)	100 (45)	cart
Container for draining coil coolant	18x 18 x 24 (457 x 457 x 610)	60 (27)	box

TABLE 8-31
CATALOG M3085TA, 1.5T TRM BODY COIL

MR COMPONENT	APPROXIMATE W x D x H in. (mm)	APPROXIMATE WEIGHT lbs (kg)	METHOD OF SHIPMENT
Body Coil Assembly	35 x 94.5 x 55.88 (889 x 2400 x 1420)	3280 (1491)	Coil on cradle mated to installation cart (See Note 1)
Note: 1 Systems With LCC Magnet Currently Installed Systems with an LCC Magnet currently installed will require the Body Coil Insertion Tool to upgrade from BRM or CRM to the new TRM coil. The Body Coil Insertion Tool requires the TRM coil to be installed on the aluminum cradle and be mated onto the aluminum coil cart. The new TRM Coil will be shipped on a wooden disposable cradle, see Illustration 2-1. The TRM Coil must be transferred to an aluminum cradle utilizing a GE Service provided coil lifting bracket and a customer rigger provided forklift or crane/hoist rated for 4,000 lbs (1814 kg), refer to Section 10-3 RIGGER/CUSTOMER SUPPLIED EQUIPMENT for additional rigger/customer supplied equipment. Systems Installing LCC Magnet Upgrade Concurrent With TwinSpeed Upgrade Systems installing an LCC Magnet Upgrade at the same time as the TwinSpeed Upgrade Systems will receive (M3060TA) LCC Magnet with TwinSpeed Enclosure with the TRM Assembly (M3085TA) installed inside the magnet bore and enclosure.			