



GE Medical Systems

Technical Publications

Direction 5133321

Revision 10

Signa[®] EXCITE HD 3.0T Field Replaceable Units (FRU)

Copyright © 2004 – 2009 by General Electric Company, Inc.

All Rights Reserved



GE Medical Systems

GE Medical Systems; Telex 3797371
P.O. Box 414, Milwaukee, Wisconsin 53201 U.S.A.
(Asia, Pacific, Latin America, North America)

GE Medical Systems – Europe; Telex 261794
Shortlands, Hammersmith, London W6 8BX U.K.



**COVER
BOTTOM
CROP LINE**

**COVER
BOTTOM
CROP LINE**

**5/16" RADIUS
CORNER TRIM
TYPICAL**



**COVER
LEFT
CROP LINE**

**COVER
FOLD
LINE**

**(2) CENTER STAPLES
FOR BINDING**

**COVER
RIGHT
CROP LINE**



**5/16" RADIUS
CORNER TRIM
TYPICAL**

SEND COMMENTS ABOUT THIS MANUAL TO:

GE MEDICAL SYSTEMS
MANAGER, SERVICE ENGINEERING, W-808
3200 N. GRANDVIEW BLVD.
WAUKESHA, WISCONSIN 53188
USA

This manual belongs to:

INTRODUCTION TO FIELD REPLACEABLE UNIT (FRU)

This direction is a guide to the FRU's.

Parts designated as a FRU are critical to maintain the safety and/or functionality of the system. FRU's are available on a Priority 1 basis from Service Parts Operation.

As a reference aid, the FRU's are listed under an upper-level part number where they are used. In most cases the upper-level part number is not a FRU. In some cases the parts listed under an upper-level assembly may have been nested in intermediate sub-assemblies which are not shown.

Where possible, FRU's which are part of Option Catalogs are shown under the area where they would be installed, rather than being listed in the Coils or Options Sections (SEE TABLE OF CONTENTS).

LIST OF EFFECTIVE PAGES

<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>
Cover Page .	10	16	9	29 to 31	8
1	10	17	10	32	1
2	2	18	1	33 to 34	3
3	8	19	2	35 to 36	2
4	8	20 to 21	8	37	6
5 to 6	6	22	1	38	3
7	9	23 to 24	5	39	8
8	5	25	8	40 to 41	2
9 to 10	6	26	9	42	8
11 to 13	7	27	6	43	1
14	2	28	7	44	10
15	1				

2

NOTES

TITLE	PAGE
INTRODUCTION TO FIELD REPLACEABLE UNIT(FRU)	1
LIST OF EFFECTIVE PAGES	1
TABLE OF CONTENTS	3
PICTORIAL TABLE OF CONTENTS	4
MAGNET AND ENCLOSURE (ASSEMBLIES)	
FRONT COVER, MISC. PARTS, & TRM	5
BRM & TRM LIST (OF PARTS NOT SHOWN)	5
ENCLOSURE & VACUUM RELATED	6
FRONT & REAR ENDBELL, BRIDGE SUPPORT	6
REAR PEDESTAL PARTS	7
SPLIT BRIDGE, MAGNET DOCKING PARTS	8
LOW PROFILE CARRIAGE ASSEMBLY	9
SRI/PAC & TOUCH-UP PAINTS	10
CHILLER	11
PATIENT SUPPORT TABLE	12, 13
PENETRATION PANEL	14
OPERATOR WORKSPACE	16
GOC COMPUTER CABINET	17
WORKSPACE SCIM AND KEYBOARDS	18
MAIN DISCONNECT PANEL	19
HFD-S/PDU CABINET	20
HFD/PDU CABINET	21
RFS CABINET	23
SYSTEM CABINET	24
MGD, RRF, & ASC MODULES	25, 26, 27, 28
RF ACCESSORY CABINET	29
NB RF AMP CABINET	30, 31
8KW MNS BB RF AMP CABINET (OPTIONAL)	32
TWIN ACCESSORY CABINET (TAC)	33, 34, 35
REPLACEMENT CABLES	36, 37
SERVICE TOOLS	38
COILS & GRAFIDY KIT	39
MNS OPTIONS	41
SOFTWARE OPTION KEYS	42
BrainWave OPTIONS	43

TABLE OF CONTENTS	3
--------------------------	----------

INTRO TO FRU'S	PAGE 1
CHILLER PARTS	PAGE 11
CABLES	PAGE 36
COILS	PAGE 39
MNS OPTIONS	PAGE 41
S/W OPTIONS	PAGE 42
BrainWave	PAGE 43




2241999
PATIENT AIR SUPPLY
BLOWER CABINET (MG6)



PENETRATION
PANEL
PAGE 14

PATIENT
SUPPORT
TABLE
PAGE 12




MAGNET ENCLOSURE:
SUBASSEMBLIES PGS 5, 6, 7, 8, 9
TOUCH-UP PAINTS PAGE 10

TRM COIL,
MOUNTING
BRACKETS, CART,
AND CRADLE,
SEE PAGE 5

MAGNET ROOM


OPERATOR/EQUIPMENT ROOM



2219400
MAGNET
MONITOR
(MSM)



MAIN
DISCONNECT
PANEL
PAGE 19




ACCESSORY
CABINET
PAGES
33 THRU 35


EXCITE OPERATOR
WORKSPACE
PAGES 16, 17, & 18




SERVICE TOOLS
PAGE 38



HFD-S/PDU
PAGE 20
HFD/PDU
PAGE 21




RFS &
SYSTEM
CABINET
PAGES
23, THRU 28




RF ACCESSORY
CABINET
PAGE 30



NB RF AMP
CABINET
PAGE 31



8KW BB RF
AMP CABINET
PAGE 32

FRONT COSMETIC COVER & PARTS (BACK VIEW):

NOTE: The following LCC 3.0T FIXED W/TRM & Enclosure related parts Not Shown;
MISCELLANEOUS:

2203803	SPLIT O-RING
2242756	RELIEF VALVE, 3/8 NPT X 7 PSIG, GEN
2299843	MAGNET PRESSURE TRANSDUCER, SUNX DP2-41N
46-294231G1	MAGNET RUNDOWN UNIT
46-294231P1	BUSSMAN, F1, 110VAC .400 GDC
46-264758G1	REMOTE INTERCOM BOARD, MG3 A2 A1

TRM:


5113878*	RF SHIELD FOR TRM COIL
2144093-9	TRM WITH MOUNTING BRACKETS, CART, AND CRADLE
2208999-7	3.0T RF BODY COIL AND CASE
2218148	COLD SHRINK (FOR REAR ENDBELL AND TRM COIL)
2299498	COOLING CIRCUIT REPAIR KIT
2383070 or 5250104	DIODE SLOT BOARD
5123812	(1260) BORE TEMPERATURE SENSOR CABLE
5123813	(1261) OPERATOR DISPLAY TO SRI CABLE
5127441	GRADIENT CABLE CLAMP BLOCK KIT
2293612	MALE CONNECTOR 3/4" HOSE TO 3/4" NPT
2324086-2	5 GALLON PAIL w/ETHYLENE GLYCOL WARNING LABEL

***2389049** RF Shield Change Out Kit – MUST be ordered for the RF Shield replacement


MAGNET AND ENCLOSURE (ASSEMBLIES)

ENCLOSURE & VACUUM RELATED PARTS:

2293771
RADIAL SHIM SPACER, 1.5MM
2293771-2
RADIAL SHIM SPACER, 0.8MM
2241176
RUBBER ISOLATION PAD

**FRONT ENDBELL PARTS:**
(TRM ONLY)


2162554-4
FRONT COVER LASER
2295639
O-RING,
ACCESS
PANELS
GASKET,
ENDBELL
TO RF TUBE
2294358, 4MM
2294358-2, 6MM
FRONT
ENDBELL ASM
2296177
2297761
O-RING, WATERLINE
FEEDTHROUGH

**BRIDGE SUPPORT PARTS:**
(TRM ONLY)

2293539
BRIDGE
SUPPORT
BUSHING
2293956
NYLON HEX NUT,
M10 X 1.5MM
X 9.5MM HIGH

**REAR ENDBELL ASSEMBLY:**
(TRM ONLY)

2284170
REAR ENDBELL
VACUUM ASSEMBLY;
& 2218148
COLD SHRINK
(FOR REAR ENDBELL &
TRM COIL)
& 2218149
COMPRESSION
SPLICER
(FOR REAR ENDBELL)



2358470
COMFORT MODULE
ASSEMBLY
2206437
M4 X 12mm NYLON
SCREW FOR
SPEAKER GRILL


GRADIENT GASKETS:
(TRM ONLY)


2297540
(4MM)
2297540-2
(6MM)



NOTE: The following parts are Not Shown.


2354050	SPECTRO T/R ASM (OPTION)
2381249	MNS LPCA CABLE HARNESS (OPTION)
2297760	VACUUM GREASE, MSDS# 8366442
2297672	GLYCOL, 1 GALLON/3.8 LITERS

REAR PEDESTAL**46-264758G1 OR 6250012**REMOTE INTERCOM
BOARD, MG3 A2 A1**46-226552G2**
ENCODER LIMIT
SWITCH BD**5224109**
LONGITUDINAL
ENCODER,
MG3 A6**3.0T HYBRID SPLITTER:****2370221**
3T DUAL BODY
HYBRID**2369611-2**
IMAGE REJECT
FILTER**2377703**
3T BODY PREAMP
& COMBINER
MODULE**MAGNET AND ENCLOSURE (ASSEMBLIES)**


BRIDGE PARTS:

NOTE: The following part is **Not Shown**.

2386348	FRONT BRIDGE FIBERGLASS
5117304	FRONT BRIDGE PULLEY ASM
2293538	BRIDGE EDGE FOAM
46-282436P1	CARRIAGE TRAVEL SENSOR FLAG


MAGNET DOCKING PARTS:

8 MAGNET AND ENCLOSURE (ASSEMBLIES)


HD LOW PROFILE CARRIAGE ASM**MAGNET AND ENCLOSURE (ASSEMBLIES)**

SRI/PAC PLATFORM:

5112668
NON-LCC
PAC-2B MODULE



5122500
ECG LEAD KIT



5105270
SRI-3
MODULE

5121793
SRI-3 U/G KIT
w/CABLES


5112668
PAC-2B MODULE

TOUCH-UP PAINTS:

2183491-3 GRAY #1 0.6 OZ. NON-AEROSOL

46-303460P3 STEEL BLUE 0.6 OZ. NON-AEROSOL


46-303460P4 DARK STEEL BLUE 0.6 OZ. NON-AEROSOL




NOTE: To obtain enclosure part numbers not found in this manual, refer to Direction 2138247, *Signa Illustrated Parts List*, for enclosure covers, endbells, and other parts.

SYSTEM COOLING CABINETS:

REFER TO ELLIS & WATTS SERVICE MANUALS FOR FRU INFORMATION





REFER TO AIRSYS SERVICE MANUALS FOR FRU INFORMATION

**2343325** AIRSYS HEAT EXCHANGER (GWHX)**2341077-2** AIRSYS QUICK DISCONNECT KIT**2341077-3** AIRSYS REMOTE CONTROL PANEL**2298100-4** E & W MODEL LTL10 HEAT EXCHANGER CABINET**2298100-5** CONTROL PANEL FOR LTL10 ONLY**2297672** GLYCOL USED WITH E & W OR AIRSYS


PATIENT TABLE – INFINITY and HD

NOTE: Refer to page 13 for information on bumper FRU's for HD, Infinity, and Pre-Infinity patient tables.


**PATIENT SUPPORT TABLE**

PATIENT TABLE – BUMPERS AND PANELS

INFINITY EXISTING PART	INFINITY REPLACEMENT FRU
46-271022G1, 2393106, Panel – Left 46-271021G1, 2393106-2, Panel – Right	5179706 , New panels and side bumpers FRU kit for Infinity table
5154908-4, Panel – Left 5154908-3, Panel – Right	5179706 , New panels and side bumpers FRU kit for Infinity table
46-271025P5, 22 Side Bumper	5161172-3 , Side Bumper 22, Color P17
46-271025P6, 16.25 Side Bumper	5161172-4 , Side Bumper 22, Color P17
46-271027P3, Foot End Bumper	46-271027P3 , Foot End Bumper
46-271026P3, Head End Bumper	46-271026P3 , Head End Bumper
Pre-INFINITY EXISTING PART	Pre-INFINITY REPLACEMENT FRU
5154908-4, Panel – Left (Grey) 5154908-3, Panel – Right (Grey)	5179706 , New panels and side bumpers FRU kit for Infinity table
22 Side Bumper (Accent Grey or Blue)	5161172-3 , Side Bumper 22, Color P8
16.25 Side Bumper (Accent Grey or Blue)	5161172-4 , Side Bumper 16.25, Color P8
Foot End Bumper (Accent Grey or Blue)	46-271027P3 , Foot End Bumper
Head End Bumper (Accent Grey or Blue)	46-271026P3 , Head End Bumper



[illegible]




NOTE: The following parts are Not Shown.


2156158-50	S/W COLLECTOR FOR EXCITE HD (12.x) M4
2389540	SCAN TOOLS MOD

GOC COMPUTER CABINET WITH HP8200:**NOTES:**


- The HP8200 has visible USB ports located on the right side of front of PC cover.
- The HP8200 cannot be replaced with the HP8000.
- If substituting HP8200 for HP8000, must order 2378152 video cable for waveform LCD monitor

**GOC COMPUTER CABINET WITH HP8000:****NOTE:**

- The HP8000 can be replaced with the HP8200.

**NOTES:**

- * If WIM is not available, order part 2141978-5 and replace cable between motherboard J1 and CIF J3.
- ** If replacing WIM with GOC AA, order cable 2395735.

LATEST VERSIONS:**KEYBOARDS ONLY:**

2275756 (English)
*2275756-2 (French)
*2275756-3 (German)
*2275756-4 (Scandinavian)
*2275756-5 (Italian)
*2275756-6 (Portuguese)
*2275756-7 (Spanish)

*** NOTE:**

These types of Keyboard/SCIM Assemblies were not available at the time this book was released, see below for earlier types.

NOTE: The following parts are Not Shown.

2262158-9	CABLE, 50 PIN SCIM KEYBOARD
2289697-2	SCIM ^W /ENGLISH OVERLAY
2289697-8	MR SCIM GLOBAL OVERLAY SET


NOTE: Refer to the Directions on CD–2160623 (Signa Service Methods) for additional product information.



TABLE 1–13

LOCATION	AMPS	VOLTS	GEMS PART NUMBER	MANUFACTURER	
				BUSS	G–S
FUSE F1	1.25	480	2276547–2	SNQ–R–1–1/4	ATQR1–1/4
FUSE F2	3	120	2276547–3	FNQ–R–3	ATQR3
FUSE F3	1	12VDC	2276547–8	LP–CC–1	ATDR1
FUSE F4	0.12	480	2276547–9	SNQ–R–8	ATQR1/8
FUSE F5	125	600	2276547–19	LPJ125SP	AJT125
FUSE F6	50	600	2276547–20	LPJ–50SP	AJT50
FUSE F7	20	600	2276547–22	LPJ–20SP	AJT20
FUSE F8	5	600	2276547–21	FRS–R–5	TRS5R
FUSE F9	8	400	2276547–6	FRS–R–8	TRS8R
FUSE F10	10	120	2276547–11	FNQ–R–10	ATQR10
FUSE F11	10	120	2276547–11	FNQ–R–10	ATQR10

2276547–14	12 VDC BATTERY
2276547–16	12 VDC SUPPLY

HFD-S/PDU CABINET ASSEMBLY W/SGA/SGAPS/GP

NOTE: The following Part is Not Shown.

2303317 CAN Terminator

HFD/PDU CABINET

5116803
HFD CABINET
FAN ASM

2284074
ALUMINUM
AIR FILTER

5118444
HFD CABINET
WITH PDU

2331400-2
GP3 MODULE

2393906-4
HFA ASSEMBLY
2396050-2
HFA CONTROL BD

2334862-2
SGA PS ASSEMBLY

46-170021P87
SGA PS FUSES


2301200-3
TVSS MODULE
2374930-2
PDU CONTROL
BOARD

HFD/PDU CABINET

21

22


NOTES



NOTE: The following items are Not Shown.

5117087 RRF SERVICE DIAGNOSTIC KIT


46-296817P1 50 OHM BNC TERMINATOR

EXISTING SYSTEM CABINET UPGRADED TO EXCITE HD

NOTE: The following items are **Not Shown**.

5117087 RRF SERVICE DIAGNOSTIC KIT

46-296817P1 50 OHM BNC TERMINATOR





NOTE: The following Array Processor ReFlex board Options are Not Shown.

2325947-12	R400 MINERVA AP (MASTER)
2325947-13	R400 MINERVA AP (SLAVE)
2325947-14	MINERVA PIM
2325947-15	R800 MINERVA AP (MASTER)
2325947-16	R800 MINERVA AP (SLAVE)

NOTE: The following MGD Module parts are Not Shown.

2294300-3	AGP BOARD
2294300-4	APS BOARD WITH IT DAUGHTER BD
2294300-5	BRIDGE BOARD
2294300-6	MGD CHASSIS W/BACKPLANE & AC ENTRY (No Boards)
2294300-7	MGD FAN
2294300-8 or	
5181495-3	POWER SUPPLIES
2294300-9	INDICATOR BOARD
2294300-19	SCP2 BOARD

2305884 REMOTE RF MOD ASM (CHASSIS W/MIDPLANE) (FRONT)**2305884 REMOTE RF MOD ASM (CHASSIS W/MIDPLANE) (REAR)**

NOTE: Additional REMOTE RF MOD ASM Parts/Cables are located on page 27.

* A part that does not have a "-2" can be replaced with either the same part number or the same number ending with a "-2". If replacing the "-2" version, then only the part number with a "-2" can be used.



** 2274188 and 2274188-2 are replaced by 2274188-3. Also, all 2274188 series MUX boards can be replaced by backwards compatible MUX2 board, 5250106.

NOTE: The following RRF MODULE (page 26) Parts/Cables are Not Shown.


46-306565P1	CABINET FILTER
2294728	COMPACT PCI FRONT PANEL BLANK (8HP) (2 SLOT)
2319152	COMPACT PCI FRONT PANEL BLANK (4HP) (1 SLOT)
2308729	BNC TO SMB CABLE (MUX J8 TO CAB I/F J1)
2308728	8W8 CABLE ASSEMBLY (MUX J5)
2322578	J15/23/25/26 DATA INPUT CABLE (CAB I/F)
2322580	CABLE J18, UNBLANK (CAB I/F)
2322581	CABLE J24, MAGNET MONITOR (CAB I/F)
2322582	CABLE J30/31, SERIAL EXP (CAB I/F)
2322583	CABLE J62 (15 PIN SUB-D) (CAB I/F)
2322586	CABLE J61 (37 PIN SUB-D) (CAB I/F)
2315093	ETHERNET CABLE
2315093-2	CABLE J63 ETHERNET (CAB I/F)
2308729-2	BNC TO SMB CABLE (MUXJ13 TO CAB I/F J2)
2308729-3	BNC TO SMB CABLE (MUXJ14 TO CAB I/F J3)
2394870-3	RRF slot14R-J5 I/O J200
5110796	RRF slot14R-J13 I/O J202
5110795	RRF slot 14R-J13 I/O J202
5115417-3	RRF slot14R-J8 RF Amp J14
5115417	RRF slot12R-J31 Splitter-1
5113193	RRF slot12R-J7 RRF slot11R-J6/7
2394870-2	RRF slot5R-J5 I/O J201
5113193-2	RRF slot3R-J7 RRF slot2R-J6/7
5115417-2	RRF slot3F-J31 Splitter-2
5115630	J7 MUX (2) TO J6 MUX (1) CABLE
5115630-2	RRF slot14R-J11 RRF slot9R-J6

NOTE: The following ASC MODULE (page 28) Parts/Cables are Not Shown.

5115417-4	ASC slot14-J1 RF Amp J8
5115417-5	ASC slot14-J4 RF Amp J7
5115417-6	ASC slot11-J1 RF Amp J6
5120955	ASC POWER SUPPLY
5250032	3.0T 8KW BB2 AMP I/F ASM

2363041 or 2363041-2* ASC CHASSIS ASSEMBLY**2363041 or 2363041-2* ASC CHASSIS ASSEMBLY WITH MNS**

NOTE: Additional ASC MODULE Parts/Cables are located on page 27.



NOTE: The following RF ACCESSORY CABINET parts Not Shown:

46-306565P1 10" FRONT COVER FILTER

2365797 CABLE - I/O J127 TO ASC SLOT 11/14 J1/2/4/5

2368552 CABLE - ASC SLOT 14-J3 TO I/O J100

2369552-8 CABLE - ASC SLOT 13-J4 TO I/O J102

2369552-9 CABLE - ASC SLOT 13-J5 TO I/O J103


2369552-10 CABLE - ASC SLOT 10-J4 TO I/O J104

2369552-11 CABLE - ASC SLOT 10-J5 TO I/O J105

2379507 DRSM JUMPER PLUG

2383066 DRSM JUMPER CABLE (ASC-DRSM-I/O)

This page for G3 MAGNET (Short Bore)
35KW NARROW BAND RF AMP



NOTE: The following NB RF AMP CABINET parts are Not Shown:

2336064-8 35KW RF AMP FUSE KIT

2336064-15 35KW RF AMP TEMPERATURE SENSOR

NOTE: The following RF TRANSMIT FLEX CABLES are Not Shown:

2384850 RUN 1156 – MR8-J4 TO CONNECT TO RUN 1125

2384850-2 RUN 1160 – MG3-A3-J1 TO CONNECT TO RUN 1127


2384851 RUN 1157 – MR8-J3 TO CONNECT TO RUN 1126

2384851-2 RUN 1158 – MR1-CPC-J2 TO CONNECT TO RUN 1137

2384851-3 RUN 1159 – MR9-J4 TO CONNECT TO RUN 1146

This page for 3T94 MAGNET (Long Bore)*

25kW NB RF AMPLIFIER CABINET (for sites WITHOUT LCD Display)




2278810-4
NB RF AMP
CABINET

NOTE: The following 25kW NB RF AMP CABINET (without LCD) parts are Not Shown:

2336064-4	35KW RF AMP IPA TUBE & TUBE PULLER KIT
2336064-5	35KW RF AMP PA TUBE & HEX WRENCH KIT
5111003	55-50dB DUAL DIRECTIONAL BODY RF COUPLER FOR E2
5111004	45-40dB DUAL DIRECTIONAL HEAD RF COUPLER FOR E2

25/35kW NB RF AMPLIFIER CABINET (for sites WITH LCD Display)



2336064-14
35KW RF
DECK
(W/CASE)

5125603-2
25KW RF AMP
POWER
SUPPLY DECK


5125603
NB RF AMP
CABINET

2336064-13
RF AND PS AIR
FILTER KIT

NOTE: The following 25/35kW NB RF AMP CABINET parts are Not Shown:


2336064-8	35KW RF AMP FUSE KIT
2336064-15	35KW RF AMP TEMPERATURE SENSOR

* **Note:** If the site has a 3T94 Long Bore magnet, it MUST use a 25KW NB RF Amplifier on this page.




NOTE: The following RF POWER AMP internal parts are **Not Shown**:


2206456-2	RF POWER AMP PREAMP/SPLITTER
2206456-3	RF AMPLIFIER (5)
2206456-4	RF POWER AMP COMBINER/COUPLER
2206456-5	RF POWER AMP COUPLER (4)
2206456-6	RF POWER AMP BIAS TEE
2206456-13	CABLE – SPLITTER TO AMPLIFIER (5)
2206456-20	COUPLER-BIAS TEE
2206456-21	POWER BIAS TEE
2206456-22	PMI INTERFACE KIT
46-170021P29	5 AMP FUSES

2239500 TAC

NOTE: The following parts are Not Shown.

2282045-2	VACUUM PUMP TIP SEAL KIT
2293608	QF25 CENTERING RING & O-RING
2293608-2	QF16 CENTERING RING & O-RING
2293609	QF25 VACUUM ELBOW
2293616	QF25 ADAPTER TO ONE INCH HOSE
2293617	VACUUM REDUCING 'T'
2293618	QF16 CLAMP
2293618-2	QF25 CHAIN CLAMP
2294010-3	HOSE CLAMP, SST ONE INCH
2297760	VACUUM GREASE, MSDS# 8366442
2329620-4	VACUUM GAUGE CONTROLLER CABLE
2284301	VACUUM GAUGE TO TAC J2 CABLE


5118162 TAC



NOTE: The following parts are Not Shown.

2282045-2	VACUUM PUMP TIP SEAL KIT
5118818	CLAMP NW 20-25
2297760	VACUUM GREASE, MSDS# 8366442

5127147
SERVICE KIT FOR TAC VACUUM
SYSTEM UPGRADE
(Used to replace old style vacuum
system, with solenoid)



TAC VACUUM SYSTEM

The table below provides information for generic replacement cables. If any of the cables are ordered, make sure to order write-on labels to be used for identifying cable part number and connection designations at the end of each cable.

GENERIC REPLACEMENT CABLES

Description	Length	
	65 Feet	100 Feet
9 Pin Sub-D Male/Female	46-244714G30	46-244714G37
9 Pin Sub-D Female/Female	46-244713G30	46-244713G37
15 Pin Sub-D Male/Female	46-244712G30	46-244712G37
15 Pin Sub-D Female/Female	46-244711G30	46-244711G37
25 Pin Sub-D Male/Female	46-244710G30	46-244710G37
25 Pin Sub-D Female/Female	46-244709G30	46-244709G37
37 Pin Sub-D Male/Female	46-258877G30	46-258877G37
	50 Feet	65 Feet
BNC RG 223	46-251710G1	46-251710G18
Cable Write-on Labels (60 labels per booklet)	46-208781P3	

DESIGNATED CABLE REPLACEMENTS

CABLE PART #	CABLE RUN
46-244710G30	Run 1227 Run 1236
46-244712G30	Run 1188
46-244713G30	Run 1181
46-244714G30	Run 1180 Run 1181 Run 1233
46-258877G30	Run 1189 Run 1190

CABLE KITS

FRU Kit Description	Part #
Receive Cables for Equipment Room, 34.5 ft (10.56m)	5119953
Receive Cables for Equipment Room, 65.6 ft (20m)	5119953-2
Receive Cables for Magnet Room, 28.2 ft (8.6m)	5119954
Receive Cables for Magnet Room, 49.2 ft (15m)	5119954-2
Transportable/Relocatable MC Rcv Cable Kit (Eqmt. Rm)	5119953-3
Body Receive Cables (both rooms), 34.5 ft (10.52m)	5119956
Body Receive Cables (both rooms), 65.6 ft (20m)	5119956-2

FIBER OPTIC CABLE KIT – FIXED SITE


Description	Part #
Includes clamps, cable ties, mounts, and Runs 1045, 1083, & 1179.	5115713

TIMES MICROWAVE CABLE PARTS

Description	Part #
Straight Male N Connector for LMR 600	2352195
7/8 EIA Right Angle Plug for LMR 1200 & Center Pin	2384859

ETHERNET CABLE KITS

Description	Part #
Ethernet cables	5160775



COILS:**3T 8 CHANNEL TORSO ARRAY COIL**

2381683-2	COIL WITH CABLE
2381683-3	PHANTOM
2382683-4	PHANTOM POSITIONER

3T 8 CHANNEL BRAIN COIL

2380637-2	BRAIN COIL/MIRROR/BASEPLATE ASM
2317112-3	PHANTOM POSITIONER
2380637-5	CABLE ASM
5146634-6	MIRROR ASM

3T GP FLEX COIL

2382817	GP FLEX COIL
2389979	GP FLEX COIL CABLE KIT

3T QUAD EXTREMITY COIL

2380635-2	KNEE/FOOT ASM W/CABLE
2380635-5	KNEE/FOOT CABLE ASM
2357819-10	EXTREMITY REPAIR KIT
2357819-11	PHANTOM POSITIONER
2357819-12	PHANTOM

3T SPLIT TOP HEAD COIL

2376114	SPLIT HEAD COIL ASM
5191687	HEAD COIL LATCHING FRU KIT
2145391	MIRROR ASM
2376104	QD BOX ASM
2355094	REPLACEMENT SPRING CONTACTS KIT

3T NV ARRAY COIL

2414383	E3 NV ARRAY W/O CABLE
2414385	E3 NV SYSTEM COIL
2414386	ANTERIOR CHEST WITH CHEST CABLE
2414389	NV PHANTOM SET
2414391	E2 NV SYSTEM CABLE

3T CTL ARRAY COIL

2413107-1	CTL ARRAY W/CABLE
2413107-4	CTL ARRAY CABLE
2381682-9	CTL PHANTOM CT SECTION
2381682-10	CTL PHANTOM TL SECTION

GRAFIDY KIT:

2270695
CuSO4 SAMPLE VIAL

40

NOTES

MNS OPTIONS:

Description	Part #
3.0T 13C TR Switch	2354050-6
3.0T 13C MNS Preamp	5114799-2
MNS QA Probe 3.0T 13C	5114403-3
MNS QA Probe 3.0T 31P	5114978-4
3.0T 13C Sample Vial	5114978-3
3.0T 31P Sample Vial	5114978-4

OPTIONS:

CONFIRMA MAMMO CAD	M3033NF
5117875	COMPLETE PC WITH SOFTWARE
ReportCARD HARDWARE	
5117866-26	ReportCARD COMPUTER
2390367-4	ReportCARD SOFTWARE COLLECTOR
2390369	ReportCARD OPTON KEY

EXCITE SOFTWARE OPTION KEYS:

2156158-60	Release 12.x Software Collector EXCITE HD (M4)
2383494	Fiesta C
2390037	Asset
2390040	ConnectPro
2390075	3D FRFSE
2390078	EchoPlus
2390079	Diffusion Tensor
2390084	Probe 2000
2390085	Probe 200 Upgrade
2390086	Sage 7
2390091	2D Delayed Enhancement
2390095	iDrive ProPlus
2390096	Cardiac Tagging
2390097	SmartPrep 2000
2390098	Flouro Trigger MRA
2390103	Performed Proc Step
5110228	Single Voxel Probe-P
5110229	Single Voxel Probe-S and Probe-P
5110230	Probe 2D – CSI
5110231	Probe 3D – CSI
5110232	MNS, SAGE, Probe-P
5110319	Research Proton Spectro
5120544	Release 12.0 M4 Apps DVD
5121548	ReportCard on AW
5121550	Tractography on AW
5122363	LAVA
5122364	PFO Evaluation
5122365	PROPELLOR Brain Suite
5122366	Tractography
5122367	2D FatSat FIESTA
5122368	3D Flair
5122370	RealTime FIESTA
5122371	VIBRANT-II
5122372	8 Channel EXCITE 3
5122373	16 Channel EXCITE 3
5122375	3.0T EXCITE ScanTools 12.0

BrainWave OPTIONS:

2363727	15" LCD Display with 5 foot VGA Cable for Computer connection
2394034-5	(Cedrus) Response Pad Asm. Left & right hand units with two buttons for each hand. Includes fiber optic cable (20 ft) and optical converter box
2394034-6	(Cedrus) "Lumina" response controller box
2394034-7	(Cedrus) AC Power module for Lumina Controller Box
2394034-9	FRU usage only Kit containing all cables [Cedrus-4, -5, -6, -7, 2xxxxxx-6, 46-265505P1 (qty 2)]
2394034-10	Penetration Sub-Panel Asm
2394034-16	PC compatible Keyboard (narrow)
2394034-17	PC compatible Mouse
5110767	BrainWave Application Software CD-ROM
5110736	EXCITE II Scanner License Key MOD for the BrainWave H/W Option
5110946	PC and Platform Software Asm
5110968	Patch CD for BrainWave Stimulus PC

REVISION HISTORY

REV/DATE REASON FOR CHANGE

- 1 Nov 16, 04 Initial release
- 2 May 5, 05 Corrected RF coil, temp sensor cable, and split bridge part numbers. Added SRI-3 upgrade kit. Added CAN terminator to PHPS. Corrected Host monitor and SCSI tower part numbers. Updated HFD & HFD-S/PDU cabinet modules. Added CAN terminator & vacuum u/g kit to TAC cabinet. Updated MDP & cable replacement information. Added CTL coil w/FRU's.
- 3 Nov 16, 05 Corrected front split bridge asm part number. For OW, added DVD upgrade kit, video cable, and notes for HP800 & HP8200 substitution rules. Updated RF Mux board part number. Revised 35KW RF Amp FRU's. Updated HOS Supply part number. Added fiber optic cable kit. Updated LV Shim phantom part numbers & added DQA Phantom location bracket.
- 4 Aug 1, 06 Corrected part #'s for LPCA cover. Added FRU's to Patient Table. Added ethernet hub to GOC and ethernet cable kit. Added Cardiac Tagging s/w option key.
- 5 Jan 26, 07 Added RF Shield Change Out Kit, Speaker Grill nylon screw, LPCA Wheel Housing FRU kit (MRIhc21559), & Dock pedal FRU parts (MRIhc17629). Added new LCD Host Monitor, Term Server cable FRU kit, and RRF Module MUX2 board. Removed MCR Tool.
- 6 Mar 13, 07 Added front cover laser, 5 gallon coolant pail, and LPCA LED cable. Added note for referring to 2138247 for parts not in this manual. Added center pin to description of Times Microwave right angle connector (CQA 13099582). Removed ECMT kit (CQA 13097585). Added Head Coil Latch Kit.
- 7 Sep 17, 07 Moved Revision History to last page. Updated replacement information for HP8000 in GOC (PQR 13121799) and HFD-S & HFD cabinet part #'s (PQR 13082817). Added LPCA bracket, pillar, and screws FRU kit (PQR 13121160), Patient Table bumper and panel FRU kit (PQR 13127363), Grafidy CuSO4 vial (PQR 13098795). Updated board part #'s in MGD & ASC modules. Added NV coil FRU parts (PQR 13113173)
- 8 May 28, 08 Rear Pedestal – updated Longitudinal Encoder part number (PQR 13175398) and added Encoder board part # 46–226552G2. HFD: Updated FHA Control Bd part #. MGD: Updated SCP2 Bd part #. Added 5181495–3 MGD P.S. Tools: Added RF Transmit Termination kit. NB RF Amp: Added Astex and MKS style Amps (PQR 13142203). Coils: Added Mirror assembly to Brain Coil (ECR 203882). Option Keys: Updated EXCITE HD (M4) software part #.
- 9 Jun 15, 09 Added Flexible Coupling (PQR 13206931) to Rear Pedestal. Updated Monitor offerings on OW (PCN 290188). Also added HP8200 DIMM (ECR 2078720) and updated to 73 Gb HDD (ECR 2065991). In RRF Chassis, added 5250028–2 RF–DIF (ECO 2080729).
- 10 Nov 3, 09 Added new GOC AA 2395043–22 (ECO 2084000).