


Chapter 5

H3 Console

Note: LightSpeed 3.X CT scanners may use any one of five Operator's Consoles. These include:

- **H3 Console** - See *This* Chapter
- **GOC1 (w/ Octane 2)** - See Chapter 6
- **GOC2 (w/ Linux PC)** - See Chapter 7
- **GOC3 (a.k.a GRE Xtream)** - See the *GOC3 Illustrated Parts List*, Direction 2394950-100. It is available for order in paper, or may be found as a PDF on the *GOC3 Installation & Service Info* CD-ROM (2384451-200) that came with the GRE Console. The PDF may also be found on Rev 12 or later of the *LightSpeed 3.X Service Information* CD-ROM (2330551-200) and on GE's Common Documentation Library.
- **GOC4 (a.k.a GRE Xtream)** - See the *GOC4 Illustrated Parts List*, Direction 5134893-100. It is available for order in paper, or may be found as a PDF on the *GOC4 Installation & Service Info* CD-ROM (5134894-200) that came with the GRE Console. The PDF may also be found on Rev 14 or later of the *LightSpeed 3.X Service Information* CD-ROM (2330551-200) and on GE's Common Documentation Library.



rev. 12/27/01

Figure 5-1 Console BOM Tree

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2333243	2.000	Yes	21 IN COLOR CRT MONITOR
2A	2254212	ALT	Yes	SONY GDM-5402T9 21" COLOR {For further detail, see Section 2.1 , beginning on page 489 .}
1A	2266832-3*	1.000	No	CONSOLE ASSY - H3 {For further detail, see Section 1.0 , beginning on page 422 .}
2B	2309592	ALT		LCD MONITOR & CABLE {For further detail, see Section 2.2 , beginning on page 489 .}

* Dual Processor Console (B7858ND) uses console assembly 2266832-4.

Table 5-1 Console w/Monitor, Octane1 (B7858NC BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2304732-3**	1.000	No	CONSOLE ASSY - H3, OCTANE II For further detail, see Section 1.0 , beginning on page 422 .}

** Dual Processor Console (B7858CD) uses console assembly 2304732-4.

Table 5-2 Octane 2 Console (B7858CC BOM, rev 0)

Section 1.0 — Console Assembly (2266832-3, 2304732-3)




Figure 5-2 Console Assembly (2266832-3, sheet 1, rev 5)




Figure 5-3 Console Assembly (2266832-3, sheet 2, rev 5)





Figure 5-4 Console Assembly (2266832-3, sheet 3, rev 5)








Figure 5-7 Console Assembly (2266832-3, sheet 6, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
7	46-328421P5	4.000	Yes	SCREW PAN HEAD 2.5 MM 8 M
8	2142006	4.000	Yes	HEX SOC SCREW M2,5 X 16
9	2216467	1.000	Yes	PEGASUS IMAGE GENERATOR P
21	2206183	1.000	Yes	NEXGEN CABLE OCTANE TO MO
25	2174128	2.000	Yes	HELIOS ETHERNET SWITCH TO
27	2286340	1.000	Yes	WALL MOUNT POWER ADAPTER
28	2281056	1.000	Yes	CABLE, PMC TO HARD DRIVE
40	2214280	1.000	Yes	REPLACEMENT HARDWARE COLL
42	2277402	1.000	Yes	SCSI CABLE, OCTANE TO BU
45	2197229-5	1.000	Yes	ETHERNET SWITCH, 8 PORT
49	2197234-2	1.000	Yes	VME BASED POWER PC CIRCUI
52	2148559	1.000	Yes	10BASE-2 TO UTP CONVERTER
53	46-297332P1	1.000	Yes	50 OHM, BNC TEE ADAP
54	2154019	1.000	Yes	50 OHM BNC GROUNDED TERMI
56	2141161	1.000	Yes	RHAPSODE BACKPLANE
60	2294357	1.000	Yes	FAN ASSEMBLY FOR POWER SUPPLY {For further detail, see Section 1.9, beginning on page 482.}
72	2199633	1.000	Yes	HELIOS RIP BD TO SIDE BUL
73	2199634	1.000	Yes	HELIOS SIDE BULKHEAD TO E
79	2142226	4.000	Yes	POWER CORD FOR EMC AREA
85	2260421	1.000	Yes	VME CONSOLE POWER SUPPLY {For further detail, see Section 1.11, beginning on page 486.}
88	2141469-2	1.000	Yes	CABLE, POWER SUPPLY TO FA
89	2191444	1.000	Yes	NEXGEN HARNESS FROM POWER
95	2280022	2.000	Yes	HANDSCREW, MOTOROLA BOARD
96	2142726	2.000	Yes	DICOM MAGNETO OPTICAL DIS
* Used on Console Assembly 2266832-3 only				
** Used on Console Assembly 2266832-4 only				

Table 5-3 Console Assembly (2266832-3, -4 BOM, rev 7)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
97	2141620	1.000	Yes	GE PC-TRAC TRACKBALL
99	2227726	1.000	Yes	MOUSE PAD FOR LIGHTSPEED
100	2256482	2.000	Yes	CABLE, 13W3-DB15 RGB CRT
101	2142225	2.000	Yes	MONITOR POWER CORD
112	46-221417P3	2.000	Yes	FEMALE SCREWLOCKS FOR HD
116	2265406	1.000	Yes	FIBER-OPTIC COUPLING RECE
117	2286051-2	1.000	Yes	CABLE, INTERNAL FIBER OPT
118	2265396	1.000	Yes	CIRCUIT BOARD-SCSI DAUGHT
119	2259530	1.000	Yes	H3 PMC DIP
125	2269876	1.000	Yes	CABLE - H2 INTERNAL ABORT
127	2245794	1.000	Yes	GLOBAL MODEM KIT
128	46-221852P2	10.000	Yes	SWITCH
1	46-302200P1	1.000	No	RATING PLATE W/ SN, MN_&
2	2174061-2	1.000	No	CONSOLE USER DRIVE ASSEMBLY {For further detail, see Section 1.1, beginning on page 442.}
3	2205488-4*	1.000	No	OCTANE COMPUTER ASSEMBLY* {For further detail, see Section 1.2, beginning on page 446.}
	2205488-5**	1.000	No	OCTANE COMPUTER ASSEMBLY- W/ D3D** {For further detail, see Section 1.2, beginning on page 446.}
4	2190560	1.000	No	FRONT EMC COVER - NEXGEN
5	2297734	1.000	No	REAR EMC COVER WITH FAN A {For further detail, see Section 1.5, beginning on page 458.}
6	2287121	2.000	No	SPACER, IG CHASSIS TO PEG
10	2132370	1.000	No	CSA LABEL
12	2102550	AR	No	TRANSPARENT PROTECTIVE TA
13	2117390	1.000	No	LABEL, CE MARKING COMPLIA
14	2132620	1.000	No	CE MARKING LABEL CLASS A
* Used on Console Assembly 2266832-3 only				
** Used on Console Assembly 2266832-4 only				

Table 5-3 Console Assembly (2266832-3, -4 BOM, rev 7) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
15	2169998	1.000	No	NEXGEN/HELIOS CONSOLE SUB {For further detail, see Section 1.6, beginning on page 460.}
19	2190506	1.000	No	REMOVEABLE CARD CAGE ASSEMBLY {For further detail, see Section 1.7, beginning on page 476.}
29	46-220184P54	2.000	No	SCREW
30	1000904P475	4.000	No	ID 0.281 OD 0.625 THK .
31	46-208935P8	4.000	No	1/4-20, 7/16" HEX X 3/16
32	46-221449P6	6.000	No	SCREW SIZE 1/4 (0.250) PH
39	46-269568P1	1.000	No	INFORMIX SOFTWARE-DATA BA
41	2169590-2	1.000	No	CR INTERCOM/PRESCRIBED TILT {For further detail, see Section 1.8, beginning on page 480.}
58	2280806	1.000	No	BRACKET, ETHERNET SWITCH
63	2103585-2	27.000	No	SCREW PAN HEAD 4 MM 10 MM
65	46-328431P1	20.000	No	WASHER PLAIN - LARGE 4.3
66	46-328421P4	1.000	No	SCREW PAN HEAD 2.5 MM 6 M
67	46-328432P1	1.000	No	WASHER LOCK - SPRING 2.6
68	46-208758P9	4.000	No	11.08"X .140" SELF-LOCKIN
74	46-136246P21	AR	No	2" WIDE X 25 YD WHITE VEL
75	46-136246P14	AR	No	2" WIDE, WHT, ADH BACK, V
78	2174056-6	1.000	No	SIDE BULKHEAD ASSY-WARP III {For further detail, see Section 1.10, beginning on page 484.}
80	2270040	1.000	No	PLATE, ACCESS H2
86	2141465	1.000	No	4 AWG 62" LONG
87	2141466	1.000	No	4 AWG 62" LONG
91	46-328432P2	3.000	No	WASHER LOCK - SPRING 4.1
93	46-208759P1	1.000	No	7.81" X .184" SELF-LOCKIN
102	46-208781P2	8.000	No	WRITE-ON LABEL FOR MARKIN
* Used on Console Assembly 2266832-3 only				
** Used on Console Assembly 2266832-4 only				

Table 5-3 Console Assembly (2266832-3, -4 BOM, rev 7) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
110	2240540	1.000	No	NETSCAPE LICENSE
113	46-170686P3	AR	No	LOCTITE 242, 50CC
120	46-328421P3	10.000	No	SCREW PAN HEAD 2.5 MM 5 M
121	2209277	2.000	No	STANDOFF/SPACER FEMALE/FE
122	46-208758P1	1.000	No	3.62" X .094" CABLE TIE
123	2255847	1.000	No	VISIBROKER RUN-TIME LICEN
124	46-220183P1	4.000	No	0.500 SQ X 0.225 THK BLK
126	2284884	1.000	No	VXWORKS 5.4 LICENCE
129	2315552	1.000	No	OCTANE MOUNTING BRACKET
130	2315552-2	1.000	No	OCTANE MOUNTING BRACKET
* Used on Console Assembly 2266832-3 only				
** Used on Console Assembly 2266832-4 only				

Table 5-3 Console Assembly (2266832-3, -4 BOM, rev 7) (Continued)






Figure 5-9 Console Assembly, w/Octane2 (2304732, sheet 2, rev 4)








Figure 5-12 Console Assembly, w/Octane2 (2304732, sheet 5, rev 0)



MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
9	2216467	1.000	Yes	PEGASUS IMAGE GENERATOR P
21	2206183	1.000	Yes	NEXGEN CABLE OCTANE TO MO
25	2174128	2.000	Yes	HELIOS ETHERNET SWITCH TO
27	2286340	1.000	Yes	WALL MOUNT POWER ADAPTER
28	2281056	1.000	Yes	CABLE, PMC TO HARD DRIVE
30	1000904P475	4.000	Yes	ID 0.281 OD 0.625 THK .
40	2214280	1.000	Yes	REPLACEMENT HARDWARE COLL
42	2277402	1.000	Yes	SCSI CABLE, OCTANE TO BU
45	2197229-5	1.000	Yes	ETHERNET SWITCH, 8 PORT
49	2197234-2	1.000	Yes	VME BASED POWER PC CIRCUI
52	2148559	1.000	Yes	10BASE-2 TO UTP CONVERTER
53	46-297332P1	1.000	Yes	50 OHM, BNC TEE ADAP
54	2154019	1.000	Yes	50 OHM BNC GROUNDED TERMI
56	2141161	1.000	Yes	RHAPSODE BACKPLANE
60	2294357	1.000	Yes	FAN ASSEMBLY FOR POWER SUPPLY {For further detail, see Section 1.9 , beginning on page 482 .}
65	46-328431P1	16.000	Yes	WASHER PLAIN - LARGE 4.3
72	2199633	1.000	Yes	HELIOS RIP BD TO SIDE BUL
73	2199634	1.000	Yes	HELIOS SIDE BULKHEAD TO E
79	2142226	4.000	Yes	POWER CORD FOR EMC AREA
85	2260421	1.000	Yes	VME CONSOLE POWER SUPPLY {For further detail, see Section 1.11 , beginning on page 486 .}
88	2141469-2	1.000	Yes	CABLE, POWER SUPPLY TO FA
89	2191444	1.000	Yes	NEXGEN HARNESS FROM POWER
95	2280022	2.000	Yes	HANDSCREW, MOTOROLA BOARD
96	2142726	2.000	Yes	DICOM MAGNETO OPTICAL DIS

* Used on Console Assembly 2304732-3 only

** Used on Console Assembly 2304732-4 only

Table 5-4 Console Assembly, w/Octane2 (2304732-3, -4 BOM, rev 6)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
97	2141620	1.000	Yes	GE PC-TRAC TRACKBALL
99	2227726	1.000	Yes	MOUSE PAD FOR LIGHTSPEED
101	2142225	2.000	Yes	MONITOR POWER CORD
112	46-221417P3	2.000	Yes	FEMALE SCREWLOCKS FOR HD
116	2265406	1.000	Yes	FIBER-OPTIC COUPLING RECE
117	2286051-2	1.000	Yes	CABLE, INTERNAL FIBER OPT
118	2265396	1.000	Yes	CIRCUIT BOARD-SCSI DAUGHT
119	2259530	1.000	Yes	H3 PMC DIP
125	2269876	1.000	Yes	CABLE - H2 INTERNAL ABORT
127	2245794	1.000	Yes	GLOBAL MODEM KIT
128	46-221852P2	10.000	Yes	SWITCH
1	46-302200P1	1.000	No	RATING PLATE W/ SN, MN_ &
2	2174061-2	1.000	No	CONSOLE USER DRIVE ASSY W {For further detail, see Section 1.1 , beginning on page 442 .}
3	2295618*	1.000	No	OCTANE2 COMPUTER ASSEMBLY* {For further detail, see Section 1.3 , beginning on page 454 .}
	2295618-2**	1.000	No	OCTANE2 COMPUTER ASSEMBLY, W D3D** {For further detail, see Section 1.3 , beginning on page 454 .}
4	2190560	1.000	No	FRONT EMC COVER - NEXGEN
5	2297734	1.000	No	REAR EMC COVER WITH FAN A {For further detail, see Section 1.5 , beginning on page 458 .}
6	2287121	2.000	No	SPACER, IG CHASSIS TO PEG
7	46-328421P5	4.000	No	SCREW PAN HEAD 2.5 MM 8 M
8	2142006	4.000	No	HEX SOC SCREW M2,5 X 16
10	2132370	1.000	No	CSA LABEL
12	2102550	AR	No	TRANSPARENT PROTECTIVE TA
13	2117390	1.000	No	LABEL, CE MARKING COMPLIA
* Used on Console Assembly 2304732-3 only				
** Used on Console Assembly 2304732-4 only				

Table 5-4 Console Assembly, w/Octane2 (2304732-3, -4 BOM, rev 6) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
14	2132620	1.000	No	CE MARKING LABEL CLASS A
15	2169998	1.000	No	NEXGEN/HELIOS CONSOLE SUB {For further detail, see Section 1.6, beginning on page 460.}
19	2190506	1.000	No	REMOVEABLE CARD CAGE ASSY {For further detail, see Section 1.7, beginning on page 476.}
29	46-220184P54	2.000	No	SCREW
31	46-208935P8	4.000	No	1/4-20, 7/16" HEX X 3/16
32	46-221449P6	4.000	No	SCREW SIZE 1/4 (0.250) PH
39	46-269568P1	1.000	No	INFORMIX SOFTWARE-DATA BA
41	2169590-2	1.000	No	CR INTERCOM/PRESCRIBED TILT {For further detail, see Section 1.8, beginning on page 480.}
58	2280806	1.000	No	BRACKET, ETHERNET SWITCH
63	2103585-2	24.000	No	SCREW PAN HEAD 4 MM 10 MM
66	46-328421P4	11.000	No	SCREW PAN HEAD 2.5 MM 6 M
67	46-328432P1	1.000	No	WASHER LOCK - SPRING 2.6
68	46-208758P9	5.000	No	11.08"X .140" SELF-LOCKIN
74	46-136246P21	AR	No	2" WIDE X 25 YD WHITE VEL
75	46-136246P14	AR	No	2" WIDE, WHT, ADH BACK, V
78	2174056-6	1.000	No	SIDE BULKHEAD ASSY-WARP III For further detail, see Section 1.10, beginning on page 484.
80	2270040	1.000	No	PLATE, ACCESS H2
86	2141465	1.000	No	4 AWG 62" LONG
87	2141466	1.000	No	4 AWG 62" LONG
93	46-208759P1	1.000	No	7.81" X .184" SELF-LOCKIN
102	46-208781P2	8.000	No	WRITE-ON LABEL FOR MARKIN
110	2240540	1.000	No	NETSCAPE LICENSE
113	46-170686P3	AR	No	LOCTITE 242, 50CC
* Used on Console Assembly 2304732-3 only				
** Used on Console Assembly 2304732-4 only				

Table 5-4 Console Assembly, w/Octane2 (2304732-3, -4 BOM, rev 6) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
120	46-328421P3	10.000	No	SCREW PAN HEAD 2.5 MM 5 M
121	2209277	2.000	No	STANDOFF/SPACER FEMALE/FE
122	46-208758P1	13.000	No	3.62" X .094" CABLE TIE
123	2255847	1.000	No	VISIBROKER RUN-TIME LICEN
124	46-220183P1	4.000	No	0.500 SQ X 0.225 THK BLK
126	2284884	1.000	No	VXWORKS 5.4 LICENCE
129	2315552	1.000	No	OCTANE MOUNTING BRACKET
130	2315552-2	1.000	No	OCTANE MOUNTING BRACKET

* Used on Console Assembly 2304732-3 only

** Used on Console Assembly 2304732-4 only

Table 5-4 Console Assembly, w/Octane2 (2304732-3, -4 BOM, rev 6) (Continued)




Figure 5-15 Console User Drive Assembly (2174061-2, sheet 2, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
7	2318598	1.000	Yes	SONY MOD SMO-F551-SD
7A	2280678	ALT	Yes	MOD, STAR T5-4.5 FIRMWARE
7B	2260734-2	ALT	Yes	MAGNETO OPTICAL DISK DRIVE
8	2315456	1.000	Yes	CD ROM DRV WITH 256K BUFFER
8A	2253371	ALT	Yes	CDROM DRIVE 32X
8B	2363007	ALT	Yes	CD-ROM DRIVE, INTERNAL, 4
11	46-208908P33	6.000	Yes	SCREW STEEL, FLAT HEAD 8-
14	2141473	1.000	Yes	12VDC FAN WITH CONNECTOR
18	2191434	1.000	Yes	NEXGEN CABLE FROM SIDE BU
20	2141463	1.000	Yes	FOR 4 DRIVES, 81" LONG
22	46-269601P1	1.000	Yes	ACTIVE REGULATED SCSI TER
23	2141468	1.000	Yes	FAN POWER EXTENSION
24	46-221417P4	2.000	Yes	FEMALE SCREWLOCKS FOR HD
26	46-170686P3	AR	Yes	LOCTITE 242, 50CC
28	46-403307P1	1.000	Yes	FAN ACCESSORIES FINGER GU
31	46-311804P2	8.000	Yes	3.2X9X.8 ZP ST DIN
38	46-220191P5	1.000	Yes	PANDUIT FCM3-A-14 GRAY
40	46-208758P2	2.000	Yes	5.5" X .140" SELF-LOCKING
41	46-311805P10	12.000	Yes	HELICAL LOCK WASHER
42	2225004-4	1.000	Yes	HARD DRIVE 18GB 15KRPM SCAN DATA
42A	2225004-3	ALT	Yes	DISK DRIVE, 18GB 1500 RP
46	2174107	1.000	Yes	HELIOS HARD RIVE BOX TO HAR
1	2138605	1.000	No	USER DRIVER ENCLOSURE TOP
5	2295626	1.000	No	LOWER PANEL OF USER DRIVE
6	46-328421P501	8.000	No	SCREW PAN HEAD 3 MM 5 MM
10	2138604-2	1.000	No	USER DRIVE FRONT BEZEL
12	2202860	1.000	No	USER DRIVE TRAY NEXGEN/HE

Table 5-5 Console User Drive Assembly (2174061-2 BOM, rev 6)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
13	46-170015P14	12.000	No	PAN HD SCREW 8-32 X 3/8
15	2141457	1.000	No	USER DRIVE FAN MOUNTING P
16	46-208897P16	8.000	No	SCREW
30	12-277	12.000	No	WASHER, FLAT
35	2153015	1.000	No	CAUTION DO NOT RMVE ASSY
37	46-170209P1	AR	No	0.034 X 1.00 X 100 FT
43	2190538	1.000	No	HARD DRIVE MOUNTING BRACKET
44	46-170015P10	4.000	No	SCREW
45	46-312679P2	4.000	No	INTERNALLY RIBBED GROMMET
47	46-328425P1	2.000	No	NUT HEXAGON 4 MM
48	46-328432P2	2.000	No	WASHER LOCK - SPRING 4.1
49	46-328431P1	6.000	No	WASHER PLAIN - LARGE 4.3
50	2213190	1.000	No	HARD DRIVE BOX GROUND STRAP
51	2103585-2	4.000	No	SCREW PAN HEAD 4 MM 10 MM
52	46-208560P53	2.000	No	SCREW

Table 5-5 Console User Drive Assembly (2174061-2 BOM, rev 6) (Continued)

1.2 Octane Computer Assembly

FRU = NO




Figure 5-16 Octane Computer Assembly (2205488, sheet 4, rev 7)

FRU = NO

USED ON:
2266832-3

MARK	NUMBER	QTY	FRU	DESCRIPTION
2	2169940-36	1.000	Yes	DUAL HEAD SE GRAPHICS OPTION/SPARE
3	2169940-37	1.000	Yes	4MB SE TEXTURE MEMORY OPTION/SPARE
5	2169940-7	1.000	Yes	OPTION PCI CARD CAGE
6A	2199806-5	2.000	Yes	2 X 128 MB OCTANE KIT/PAIR
6B	2280645	ALT	Yes	256MB MEMORY KIT, KINGSTON
7	2169940-8	1.000	Yes	OPTION PCI ETHERNET CARD
8	2169940-12	1.000	Yes	OPTION 9GB DISK DRIVE
8A	2169940-52	ALT	Yes	9 GB "PISCES" DISK DRIVE
8B	2169940-46	ALT	Yes	PCA NEW 9GB DISK DRIVE
9	2282000	1.000	Yes	PCI 4 PORT SERIAL CARD
10	2169940-9	1.000	Yes	OPTION PCI SCSI CARD
1	2169940-43	1.000	No	R12K BASE CT/MR OCTANE COMPUTER

Table 5-6 Octane Computer Assembly (2205488-4 BOM, rev 2)




Figure 5-17 Octane Computer Assembly w/D3D (2205488, sheet 5, rev 7)


FRU = NO

USED ON:
2266832-4

MARK	NUMBER	QTY	FRU	DESCRIPTION
2	2169940-36	1.000	Yes	DUAL HEAD SE GRAPHICS OPTION
3	2169940-37	1.000	Yes	4MB SE TEXTURE MEMORY OPTION
5	2169940-7	1.000	Yes	OPTION PCI CARD CAGE
6A	2199806-8	2.000	Yes	MEMORY, 2X256MB, LOW PROF
6B	2280645-2	ALT	Yes	512MB MEMORY KIT, KINGSTON
7	2169940-8	1.000	Yes	OPTION PCI ETHERNET CARD
9	2282000	1.000	Yes	PCI 4 PORT SERIAL CARD
10	2169940-9	1.000	Yes	OPTION PCI SCSI CARD
1	2169940-50	1.000	No	GEMS BASE DUAL R12K OCTANE


Table 5-7 Octane Computer Assembly w/D3D (2205488-5 BOM, rev 0)

Octane Computer



Revision 11/15/00

Figure 5-18 Octane Computer, (sheet 1, rev 11/14/00)



Revision 10/27/99

Figure 5-19 Octane Computer, (sheet 2, rev 10/27/99)

MARK	NUMBER	QTY	FRU	DESCRIPTION
3	2169940- 3	1	Yes	GRAPHICS CARD 1 AND CARD 2 / SGI# 030- 0938- 003
3A	2169940- 36	ALT	Yes	GRAPHICS CARD 1 AND CARD 2 / SGI PART# HU- DHSE OR 030- 1241- 002
4	2169940- 4	1	Yes	TEXTURE MEMORY (4MB) / SGI# 030- 0961- 002
4A	2169940- 37	ALT	Yes	TEXTURE MEMORY (4MB) / SGI PART# HU- SETEX OR 030- 1263- 002
7	2169940- 7	1	Yes	PCI MODULE / SGI# PCI- CARD CAGE
10A	2169940-52	1	Yes	9 GB "PISCES" DISK DRIVE
10B	2169940-46	ALT	Yes	PCA NEW 9GB DISK DRIVE
13	2169940- 45	1	Yes	SYSTEM MODULE WITH OUT PROCESSOR OR DIMMS / SGI# 030- 0887- 003
14	2169940- 14	1	Yes	POWER SUPPLY / SGI# 060- 0038- 001
15	2169940- 15	1	Yes	XIO MODULE / SGI# 013- 1600- 005
21	2169940- 21	1	Yes	COMPRESSION CAP / SGI# 050- 0358- 002
24	2169940- 24	1	Yes	FRONT PLANE ASSEMBLY / SGI# 013- 1597- 001
25	2169940- 25	1	Yes	SYSTEM FAN / SGI# 013- 2000- 001
26	2169940- 26	1	Yes	SYSTEM IDENTIFICATION MODULE / SGI# 013- 2501- 001
27	2169940- 27	1	Yes	DRIVE BAYS
28	2169940- 28	1	Yes	MOUSE / SGI# 063- 0009- 001
32	2169940- 32	1	Yes	LIGHT MODULE / SGI# 013- 1918- 001
40	2169940- 8	1	Yes	PCI ETHERNET CARD

Table 5-8 OCTANE Computer Parts

MARK	NUMBER	QTY	FRU	DESCRIPTION
13	2169940- 45	1	Yes	SYSTEM MODULE WITHOUT PROCESSOR OR DIMMS - MOTHERBOARD
18A	2169940- 44	1	Yes	270MHZ R12K SINGLE PROCESSOR / SGI PART # 030-1432-001
18B	2169940- 48	1	Yes	270MHZ R12K DUAL PROCESSOR / SGI PART # 030-1433-001 (DIRECT-3D OPTION)
42	2199806- 8	1 OR 3	Yes	256MB SDRAM DIMMS (PAIR) (STANDARD & OPTION - DIRECT-3D)

Table 5-9 OCTANE Computer (270 MHz R12K and Dual Processor)

MARK	NUMBER	QTY	FRU	DESCRIPTION
18C	2169940-49	1	Yes	300MHZ R12K DUAL PROCESSOR / SGI PART # 030-1356-001 (DIRECT 3D OPTION)

Table 5-10 OCTANE Computer (300 MHz D3D Option)

1.3 Octane2 Computer Assembly

FRU = NO




Figure 5-20 Octane2 Computer Assembly (2295618, sheet 1, rev 0)

FRU = NO

USED ON:
2304732-3

MARK	NUMBER	QTY	FRU	DESCRIPTION
5	2282000	1.000	Yes	PCI 4 PORT SERIAL CARD
7	2199806-8	2.000	Yes	MEMORY, 2X256MB, LOW PROF
7A	2280645-2	ALT	Yes	512MB MEMORY KIT, KINGSTO
8	2304737-2	2.000	Yes	CABLE, DVI VIDEO CABLE -
1	2169940-76	1.000	No	CT OCTANE 2 BASE BUY LEVE
1A	2169940-61	ALT	No	CT LIGHTSPEED OCTANE II B
2	2169940-63	1.000	No	DUAL CHANNEL DSPLY DAUGHT
3	2169940-7	1.000	No	OPTION PCI CARD CAGE
4	2169940-8	1.000	No	OPTION PCI ETHERNET CARD
6	2169940-9	1.000	No	OPTION PCI SCSI CARD

Table 5-11 Octane2 Computer Assembly (2295618 BOM, rev 1)


FRU = NO

USED ON:
2304732-4

MARK	NUMBER	QTY	FRU	DESCRIPTION
5	2282000	1.000	Yes	PCI 4 PORT SERIAL CARD
7A	2199806-8	2.000	Yes	MEMORY, 2X256MB, LOW PROF
7B	2280645-2	ALT	Yes	512MB MEMORY KIT, KINGSTO
8	2304737-2	2.000	Yes	CABLE, DVI VIDEO CABLE -
1	2169940-77	1.000	No	CT OCTANE 2 D3D BUY LEVE
1	2169940-62	ALT	No	CT D3D OPTION OCTANE II B
2	2169940-63	1.000	No	DUAL CHANNEL DSPLY DAUGHT
3	2169940-7	1.000	No	OPTION PCI CARD CAGE
4	2169940-8	1.000	No	OPTION PCI ETHERNET CARD
6	2169940-9	1.000	No	OPTION PCI SCSI CARD

Table 5-12 Octane2 Computer Assembly (2295618-2 BOM, rev 1)

1.4 Octane2 Replacement Parts



Revision 10/09/01

Figure 5-21 Octane2 Replacement Parts (sheet 1, rev 10/09/01)

MARK	NUMBER	QTY	FRU	DESCRIPTION
1A	2300169	1	Yes	OCTANE II FRU, SINGLE PROC UNIT
1B	2300169-2	ALT	Yes	OCTANE II FRU, DUAL PROC UNIT
7	2300173	1	Yes	CT PCI CARD CAGE ASSEMBLY
10	2169940-78	1	Yes	18GB 10K RPM DISK DRIVE (X2)

Table 5-13 Octane2 Replacement Parts

1.5 Rear Cover W/Fan (2297734)

FRU = NO

USED ON:
2266832-3,
2266832-4,
2304732-3,
2304732-4




Figure 5-22 Rear Cover w/Fan (2297734, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2297827	1.000	Yes	80X80X15 MM FAN ASSEMBLY
1	2190561	1.000	No	REAR EMC COVER - NEXGEN
3	46-208897P16	8.000	No	SCREW
4	2325173	2.000	No	USER DRIVE FAN GASKET

Table 5-14 Rear EMC Cover w/Fan Assembly (2297734 BOM, rev 1)




Figure 5-24 Nexgen/Helios Console Sub-Assy (2169998, sheet 2, rev 7)




Figure 5-25 Nexgen/Helios Console Sub-Assy (2169998, sheet 3, rev 6)




Figure 5-26 Nexgen/Helios Console Sub-Assy (2169998, sheet 4, rev 7)

MARK	NAME	QTY	FRU	DESCRIPTION
4	2169593	1.0	Yes	CR RHAPSODE MONITOR TABLE TOP
5	46-170015P38	12.0	Yes	SCREW SCREW OLD TOOTH F70B5A ZINC PL
9	12-275	12.0	Yes	WASHER, ID .2656, OD.75, THK .050 STEEL F70B5A
10	46-208695P73	1.0	Yes	HEYCO PART NO. 2198 SB-1250-15-BLACK
12	46-220205P11	12.0	Yes	#10 .219 ID, .500 O.D. .049 THK
13	46-170012P40	2.0	Yes	1/4-20 HEX 7/16 NUT THK 3/16 O~ALL HT 15/64 WASH OD 1/2 WASH MAT THK .025 ZINC PL KEPS
14	46-328431P1	5.0	Yes	WASHER PLAIN - LARGE 4.3 MM 12 MM ZINC PLATED, STEEL
27	2191429	1.0	Yes	NEXGEN CABLE FROM COMPUTER TO SERVICE KEY
28	2169593-2	1.0	Yes	CR RHAPSODE KEYBOARD TABLE TOP
30	46-170686P3	0.0	Yes	LOCTITE 242, 50CC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER.
31	46-221417P5	1.0	Yes	FEMALE SCREWLOCKS FOR HD CONNECTORS. 4-40 THREAD, 0.438~ OVERALL LENGTH YELLOW CHROMATE, TWO PER KIT. 0.250~ BOTTOM, 0.062~ PANEL THK. MAX (KIT INCLUDES NUTS AND WASHERS)
32	2103585	5.0	Yes	SCREW PAN HEAD 4 MM 8 MM SPRING WASHER, ZINC PLATED, STEEL
33	46-170015P6	12.0	Yes	SCREW SCREW OLD F70B5A ZINC PL
36	46-170498P48	4.0	Yes	SCREW SCREW OLD F70B5A2 ZINC PL HT TR FOR HYD EMB
37	46-220205P2	4.0	Yes	ID 0.375 OD 0.875 THK 0.083 F70B5A (FOR 5/16 IN.SCREW)
38	2150525	2.0	Yes	CT RHAPSODE DOLLY
40	46-208695P75	1.0	Yes	HEYCO NO.2210 SB-1375-16 BLACK
41	46-208758P2	6.0	Yes	5.5~ X .140~ SELF-LOCKING CABLE TIE. 1/16~ TO 1-1/4~ WIRE BUNDLE RANGE. 30LB TENSILE STRENGTH, NYLON MATR~L.
43	46-170108P1	0.0	Yes	DBL COAT NEOPRENE,CLOSED CELL 0.031 X 0.500 BLACK 50 LB/CU FT (REV TK, 3/92). OK
45	2141473	1.0	Yes	12VDC FAN WITH CONNECTOR
46	46-403307P1	2.0	Yes	FAN ACCESSORIES FINGER GUARD (USED WITH ETRI 98, 125, 129, 141 AND 144 SERIES FANS). USED WITH ETRI 98, 125, 129, 141 AND 144 SERIES FANS (4-1/2~). (DWG REV 2/92, JFK - OK).

Table 5-15 Nexgen/Helios Console Sub-Assy (2169998 BOM, rev 20)

MARK	NAME	QTY	FRU	DESCRIPTION
48	2196381	1.0	Yes	FAN AND CONNECTOR
49	2196974	2.0	Yes	FAN ACCESSORIES FINGER GUARD
50	2208553	1.0	Yes	FAN POWER EXTENSION CABLE
51	46-208897P16	16.0	Yes	SCREW SCREW OLD
52	46-328425P1	4.0	Yes	NUT HEXAGON 4 MM
53	46-328432P2	4.0	Yes	WASHER LOCK - SPRING 4.1 MM 7.6 MM ZINC PLATED, SPRING STEEL
54	12-281	2.0	Yes	ID 0.2813 OD 0.625 THK .050 B11H89C BRS OLD #61
57	46-170684P3	0.0	Yes	LOCTITE 271, 50ML BOTTLE. 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER. (REV, TK, 1/92).
60	2318465	1.0	Yes	OCTANE MNTG PLATE
61	2103585-2	3.0	Yes	SCREW PAN HEAD 4 MM 10 MM SPRING WASHER, ZINC PLATED, STEEL
62	46-328431P1	3.0	Yes	WASHER PLAIN - LARGE 4.3 MM 12 MM ZINC PLATED, STEEL
1	2169588	1.0	No	NEXGEN/HELIOS CABINET ASSY {For further detail, see Section 1.6.1 , beginning on page 466 .}
2	2170651	1.0	No	NEXGEN/HELIOS CONSOLE CHASSIS
3	2174159	1.0	No	OUTLET BOX ASSEMBLY {For further detail, see Section 1.6.2 , beginning on page 470 .}
39	46-208695P76	1.0	No	HEYCO PART NO.2211 SB-1375-16 WHITE
42	2154226	1.0	No	RAPSODE CONSOLE SHIPPING PALLET
44	2208538	1.0	No	FAN MOUNTING BRACKET-LOWER LEFT SIDE
47	2208540	1.0	No	FAN MOUNTING BRACKET-UPPER RIGHT
55	2241388	1.0	No	SURFACE RACEWAY 3 FEET
56	2307744	1.0	No	KEYBOARD SHELF SPPT BKT
58	2315554	7.0	No	DOUBLE THREADED SCREW SIZE 9/10
59	2307744	2.0	No	KEYBOARD SHELF SPPT BKT

Table 5-15 Nexgen/Helios Console Sub-Assy (2169998 BOM, rev 20) (Continued)

1.6.1 Cabinet Assembly (2169588)

FRU = NO

USED ON:
2169998




Figure 5-27 Cabinet Assembly (2169588, sheet 1, rev 2)




Figure 5-28 Cabinet Assembly (2169588, sheet 2, rev 1)

MARK	NAME	QTY	FRU	DESCRIPTION
1	2138395	1.0	Yes	RHAPSODE CONSOLE RIGHT SIDE PANEL
2	2138395-2	1.0	Yes	RHAPSODE CONSOLE LEFT SIDE PANEL
3	2138394	2.0	Yes	RHAPSODE CONSOLE LOWER BASE CVR
4	2173114	1.0	Yes	NEXGEN/HELIOS CONSOLE FRONT PANEL
6	2138390	2.0	Yes	RHAPSODE CONSOLE STABILIZER FOOT
7	2138388	2.0	Yes	RHAPSODE CONSOLE MONITOR TABLE BRKT
9	2200318	1.0	Yes	KEYBRD TBLTOP BRACE WITH LATCH
10	2144080	2.0	Yes	RHAPSODE CONSOLE SEISMIC/SHIP BKT
11	2173113	2.0	Yes	CR-CT/I MONITOR TABLE TOP BRACE
5	2283157	1.0	No	REAR ACCESS PANEL, CONSOLE
8	2144033	1.0	No	RHAPSODE CONSOLE FRAME

Table 5-16 Cabinet Assembly (2169588 BOM, rev 2)

This Page Intentionally Blank.

1.6.2 Outlet Box Assembly (2174159)

FRU = NO

USED ON:
2169998




Figure 5-29 Outlet Box Assembly (2174159, sheet 1, rev 1)

MARK	NAME	QTY	FRU	DESCRIPTION
4	46-278446P1	8.0	Yes	15A, 125VAC, 2 POLE, 3 WIRE DUPLEX TYPE, GRAY IN COLOR. (REV TK, 7/91).
5	46-297335P3	2.0	Yes	16A, 250V 2-STAGE PWR LINE FILTER
6	46-327213P1	1.0	Yes	AMPS 30 VOLTS 120/208 5 WIRES FLANGED INLET, MALE, U.L. LISTED
24	46-328425P10	24.0	Yes	NUT HEXAGON 3 MM
25	46-328425P1	3.0	Yes	NUT HEXAGON 4 MM
28	46-328430P2	3.0	Yes	WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL
30	46-208695P61	1.0	Yes	HEYCO PART NO.2126 SB-875-11 BLACK
31	46-271110P1	1.0	Yes	LABEL, PROTECTIVE EARTH (GROUND) IEC STANDARD GRAPHIC SYMBOL
37	46-311804P2	20.0	Yes	3.2X9X.8 ZP ST DIN
38	46-311805P10	24.0	Yes	HELICAL LOCK WASHER M3 ZP DIN
39	46-208758P2	3.0	Yes	5.5~ X .140~ SELF-LOCKING CABLE TIE. 1/16~ TO 1-1/4~ WIRE BUNDLE RANGE. 30LB TENSILE STRENGTH, NYLON MATR~L.
1	2172419	1.0	No	OUTLET BOX COVER
2	2173324	1.0	No	DIVIDER - OUTLET BOX
9	2200241-3	1.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
10	2200241-2	1.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
12	2200241-5	1.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
15	2200241-8	1.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
16	2200241-9	1.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
17	2200241-10	1.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
18	2200241-11	3.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
19	2200241-12	2.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
20	2200241-13	6.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
21	2200241-14	1.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
22	2200241-15	3.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
23	2200241-16	1.0	No	WIRE DATA FOR NEXGEN OUTLET BOX
26	46-312358P30	3.0	No	M3 X 0.5 X 12MM METRIC FREEDRIV ZINC PLATED CROSS RECESSED PAN HEAD MACHINE SCREW

Table 5-17 Outlet Box Assembly (2174159 BOM, rev 1)

MARK	NAME	QTY	FRU	DESCRIPTION
32	2190750	1.0	No	CONSOLE POWER - NEXGEN HELIOS {For further detail, see “Breaker Box Assembly (2190750),” on page 474.}
33	46-313024P5	1.0	No	No description provided
34	46-208697P57	2.0	No	HEYCO PART NO.1247 SR-7W-2 BLACK
35	2150839	1.0	No	CAUTION RECEPTACLES NOT FOR GEN USE
36	46-313024P2	0.0	No	No description provided

Table 5-17 Outlet Box Assembly (2174159 BOM, rev 1) (Continued)

This Page Intentionally Blank.

Breaker Box Assembly (2190750)

FRU = NO

USED ON:
2174159




Figure 5-30 Breaker Box Assembly (2190750, Sheet 1, rev 2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2150833	1	2	BREAKER SWITCH
5	2200242	2	2	POWER CABLE FROM OUTLET B
11	46-170012P39	1	2	10-32, 3/8" HEX KEPS NUT,
1	2200158	1	N	BREAKER BOX COVER PLATE
2	2200159	1	N	FRONT CIRCUIT BREAKER BOX
4	2150830	1	N	POWER BREAKER BEZEL
6	2190749	1	N	CIRCIUT BREAKER MOUNTING
7	46-208697P57	2	N	HEYCO PART NO.1247 SR-7W
8	46-328430P2	6	N	WASHER PLAIN - NORMAL 4.3
9	2103585	4	N	SCREW PAN HEAD 4 MM 8 MM
10	46-328425P1	1	N	NUT HEXAGON 4 MM
13	46-313024P2	AR	N	
14	2108900	1	N	CONNECTOR 8 CONTACTS 150

Table 5-18 Breaker Box Assembly (2190750 BOM, rev 3)

MARK	NAME	QTY	FRU	DESCRIPTION
2	2190504	1.0	Yes	REM CARD CAGE FAN ASSY-NEXGEN HELIOS {For further detail, see "NexGen Fan Assembly (2190504)," on page 478.}
4	46-221513P7	4.0	Yes	CLAMP OR CLIP 8.0~ X 0.312~ NYLON CARD GUIDE. 3 TABS SNAP INTO A .187~ MGT HOLE. (CARD THICKNESS-0.076~) GUIDES SNAP MOUNT INTO A 0.187~ HOLE . (DWG REV 9/91, JFK).
5	46-221513P3	12.0	Yes	CLAMP OR CLIP 6.0~ X 0.312~ NYLON CARD GUIDE. 2 TABS SNAP INTO A .187~ MGT HOLE. (CARD THICKNESS-0.077~) GUIDES SNAP MOUNT INTO A 0.187~ HOLE . (DWG REV 9/91, JFK).
6	46-296128P2	4.0	Yes	CARD GUIDE_9U
7	46-208758P2	1.0	Yes	5.5~ X .140~ SELF-LOCKING CABLE TIE. 1/16~ TO 1-1/4~ WIRE BUNDLE RANGE. 30LB TENSILE STRENGTH, NYLON MATR~L.
9	2103585	2.0	Yes	SCREW PAN HEAD 4 MM 8 MM SPRING WASHER, ZINC PLATED, STEEL
11	46-220178P10	0.0	Yes	GROMMET 12.75~ LONG GROMMET STRIP, OFF WHITE FOR 0.090~ THICK APPLICATIONS
12	2141472	2.0	Yes	DRIVE POWER EXTENSION
13	46-328421P5	6.0	Yes	SCREW PAN HEAD 2.5 MM 8 MM ZINC PLATED, STEEL
17	46-328432P1	6.0	Yes	WASHER LOCK - SPRING 2.6 MM 5.1 MM ZINC PLATED, SPRING STEEL
18	46-221513P7	2.0	Yes	CLAMP OR CLIP 8.0~ X 0.312~ NYLON CARD GUIDE. 3 TABS SNAP INTO A .187~ MGT HOLE. (CARD THICKNESS-0.076~) GUIDES SNAP MOUNT INTO A 0.187~ HOLE . (DWG REV 9/91, JFK).
19	46-221513P3	2.0	Yes	CLAMP OR CLIP 6.0~ X 0.312~ NYLON CARD GUIDE. 2 TABS SNAP INTO A .187~ MGT HOLE. (CARD THICKNESS-0.077~) GUIDES SNAP MOUNT INTO A 0.187~ HOLE . (DWG REV 9/91, JFK).
20	46-170108P1	0.0	Yes	DBL COAT NEOPRENE,CLOSED CELL 0.031 X 0.500 BLACK 50 LB/CU FT (REV TK, 3/92). OK
1	2190507	1.0	No	REMOVEABLE CARD CAGE - HELIOS
3	2280020	1.0	No	PANEL, FILLER - PEGASUS
10	2141475	1.0	No	RHAP CONSOLE VME FILLER PANEL
16	2210165	1.0	No	FILLER PANEL - HELIOS CARD CAGE ASSY

Table 5-19 Removable Card Cage Assembly (2190506 BOM, rev 5)

1.7.1 NexGen Fan Assembly (2190504)

FRU = YES




Figure 5-32 NexGen Fan Assembly (2190504, Sheet 1, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2141473	3	Yes	12VDC FAN WITH CONNECTOR
7	2141467	1	Yes	FAN DISTRIBUTION
11	46-403307P1	3	Yes	FAN ACCESSORIES FINGER GU
1	2190512	1	No	FAN MOUNTING BRACKET-REM
3	2140987	1	No	CONSOLE,3 CONNECTOR FAN B
4	2140988	1	No	CONSOLE,SINGLE CONN. BRKT
5	46-208897P16	24	No	SCREW
6	2103585	4	No	SCREW PAN HEAD 4 MM 8 MM
9	46-208758P2	2	No	5.5" X .140" SELF-LOCKING

Table 5-20 NexGen Fan Assembly (2190504 BOM, rev 1)

1.8 CR Intercom/Prescribed Tilt Circuit (2169590-2)

FRU = NO

USED ON:
2266832-3,
2266832-4,
2304732-3,
2304732-4




Figure 5-33 CR Intercom/Prescribed Tilt Circuit (2169590-2, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2167014	1.000	Yes	RHAPSODE INTERCOM
3	46-208560P48	6.000	Yes	SCREW
4	1000904P279	6.000	Yes	ID 0.1910 OD 0.3125 THK
6	46-221417P4	1.000	Yes	FEMALE SCREWLOCKS FOR HD
8	2269601	1.000	Yes	PRESCRIBED TILT BOARD
9	46-170015P41	6.000	Yes	SCREW
2	2268979	1.000	No	BRACKET - CIRCUIT BOARD M
5	46-287187P4	2.000	No	SCREW, SLT HD FLT, 2.5 x 6.0 MM
7	46-170686P3	AR	No	LOCTITE 242, 50CC

Table 5-21 CR Intercom/Prescribed Tilt Circuit (2169590-2 BOM, rev 1)

1.9 Box Assy - Fan and Cables (2294357)

FRU = YES




Figure 5-34 Box Assy - Fan and Cables (2294357, sheet 1, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
5	46-296212P1	1.000	Yes	NEXGEN HARNESS FROM HARD
1	2294356	1.000	No	HARD DRIVE ENCLOSURE HELI
2	2196972	2.000	No	FAN AND CONNECTOR
3	2196974	4.000	No	FAN ACCESSORIES FINGER GU
4	46-208897P16	16.000	No	HELIOS HARD DR BOX TO HAR
6	511A590P316	2.000	No	SCREW

Table 5-22 Box Assy - Fan and Cables (2294357 BOM, rev 0)

1.10 Side Bulkhead Assy (2174056-5)

FRU = NO

USED ON:
2266832-3,
2266832-4,
2304732-3,
2304732-4




Figure 5-35 Bulkhead Assembly, Side (2174056-5, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2269875-2	1.000	Yes	CABLE-H2 S/BLK TO RIP BD,
3	2269874-2	1.000	Yes	CABLE-H2 OCTOPUS
8	2191450	1.000	Yes	NEXGEN HARNESS FROM POWER
11	2142224	1.000	Yes	VME POWER SUPPLY TO RECEP
12	2256483-2	1.000	Yes	SCSI CABLE, S/BLKHD TO PP
1	2174374-4	1.000	No	BRACKET, SIDE BULKHEAD
6	46-316540P3	2.000	No	SCREW PAN HEAD 3 MM 10 MM
7	46-221417P4	2.000	No	FEMALE SCREWLOCKS FOR HD
9	46-170686P3	AR	No	LOCTITE 242, 50CC
13	2109168-3	1.000	No	FERRITE CLIP, 13.4MM

Table 5-23 Side Bulkhead Assembly - Helios 2 (2174056-5 BOM, rev 1)

1.11 VME Console Power Supply Assembly (2260421)

FRU = YES




Figure 5-36 VME Console Power Supply Assembly (2260421, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2255666	1.000	No	POWER SUPPLY TRAY NEXGEN/
2	2255986	1.000	No	POWER SUPPLY, MODULAR,600
3	2255660	2.000	No	5V POWER SUPPLY ADAPTER
4	46-220181P5	2.000	No	WASHER #10 (.190) SCREW
5	12-276	2.000	No	ID 0.2040 OD 0.5630 THK
6	46-208941P9	2.000	No	00010-32 .375 HEX X .124
7	12-277	2.000	No	WASHER, FLAT
8	2103585-2	2.000	No	SCREW PAN HEAD 4 MM 10 MM
9	46-170015P8	4.000	No	SCREW

Table 5-24 VME Console Power Supply Assembly (2260421 BOM, rev 0)

1.12 Global Modem Kit (2245794)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2216697	1.000	Yes	CABLE, RJ11 MODULAR, 6-CO
5	2249450	1.000	Yes	ADAPTER ASSY GLOBAL MODEM
1	2245773	1.000	No	MULTI-TECH GLOBAL MODEM
2	2253088	1.000	No	LABEL,CUSTOMER PROPERTY F
4	2255874	1.000	No	CABLE, RJ11 MODULAR,6-CON

Table 5-25 Global Modem Kit (2245794 BOM, rev 3)

Section 2.0 — Console Display Monitors

2.1 CRT Monitors

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1A	2189192	2	Yes	PHILIPS 20 INCH MONITOR FOR LIGHTSPEED
1B	2244724	2	Yes	SONY GDM-500PS 21" COLOR MONITOR (RGB INPUTS)
1C	2254212	2	Yes	SONY GDM-5402T9 21" COLOR MONITOR (DB15 INPUT)

Table 5-26 Console Display Monitors

2.2 LCD Monitor & Cable (2309592)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2349792	1.000	Yes	18.1 IN LCD COLOR MONITOR
2	2295019	1.000	Yes	13W3-BD15 RGB LCD STRAIGHT CONN CBL

Table 5-27 LCD Monitor & Cable (2309592 BOM, rev 1)

Section 3.0 — DASMs

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1A	2249894	1	Yes	Analog DASM
1B	2249895	1	Yes	Digital DASM

Table 5-28 DASMs

Section 4.0 — Keyboard w/SCIM, Collectors



Figure 5-37 Keyboard & SCIM Assembly (Rev 11/15/00 - English shown)

4.1 English (2275758)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2275756	1.000	Yes	KEYBOARD, NORTH AMERICAN
3	2275752	1.000	Yes	SCIM COLLECTOR - ENGLISH {For further detail, see SCIM Collector – English (2275752) , below.}
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757	1.000	No	FUNCTION KEY TEMPLATE - E

Table 5-29 Keyboard w/SCIM Collector - English (2275758 BOM, rev 1)

SCIM Collector – English (2275752)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2262158	1.000	No	SCAN CONTROL INTERCOM MOD
2	2262158-11	1.000	No	OVERLAY, SCIM CT, TILT -
3	2262158-21	1.000	No	OVERLAY, SCIM CT, NO TILT

Table 5-30 SCIM Collector - English (2275752 BOM, rev 0)

4.2 French (2275758-2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2275756-2	1.000	Yes	KEYBOARD, EUROPEAN FRENCH
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757-2	1.000	No	FUNCTION KEY TEMPLATE - F
3	2275752-2	1.000	No	SCIM COLLECTOR - FRENCH {For further detail, see SCIM Collector - French (2275752-2) , below.}

Table 5-31 Keyboard w/SCIM Collector - French (2275758-2 BOM, rev 1)

SCIM Collector - French (2275752-2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2262158	1.000	No	SCAN CONTROL INTERCOM MOD
2	2262158-12	1.000	No	OVERLAY, SCIM CT, TILT -
3	2262158-22	1.000	No	OVERLAY, SCIM CT, NO TILT

Table 5-32 SCIM Collector - French (2275752-2 BOM, rev 0)

4.3 German (2275758-3)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2275756-3	1.000	Yes	KEYBOARD, GERMAN PC PS/2
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757-3	1.000	No	FUNCTION KEY TEMPLATE - G
3	2275752-3	1.000	No	SCIM COLLECTOR - GERMAN {For further detail, see SCIM Collector - German (2275752-3) , below.}

Table 5-33 Keyboard w/SCIM Collector - German (2275758-3 BOM, rev 1)

SCIM Collector - German (2275752-3)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2262158	1.000	No	SCAN CONTROL INTERCOM MOD
2	2262158-13	1.000	No	OVERLAY, SCIM CT, TILT -
3	2262158-23	1.000	No	OVERLAY, SCIM CT, NO TILT

Table 5-34 SCIM Collector - German (2275752-3 BOM, rev 0)

4.4 Swedish (2275758-4)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2275756-4	1.000	Yes	KEYBOARD, SWEDISH PC PS/2
3	2275752	1.000	Yes	SCIM COLLECTOR - ENGLISH {For further detail, see " SCIM Collector – English (2275752) ," on page 491.}
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757	1.000	No	FUNCTION KEY TEMPLATE - E

Table 5-35 Keyboard w/SCIM Collector - Swedish (2275758-4 BOM, rev 1)

4.5 International (2275758-5)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2275756	1.000	Yes	KEYBOARD, NORTH AMERICAN
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2283830	1.000	No	FUNCTION KEY TEMPLATE COL
3	2275752-4	1.000	No	SCIM COLLECTOR - INTERNATIONAL {For further detail, see SCIM Collector - International (2275752-4) , below.}

Table 5-36 Keyboard w/SCIM Collector - Swedish (2275758-5 BOM, rev 1)

SCIM Collector - International (2275752-4)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2262158	1.000	No	SCAN CONTROL INTERCOM MOD
2	2262158-10	1.000	No	OVERLAY, SCIM CT, TILT -
3	2262158-14	1.000	No	OVERLAY, SCIM CT, TILT -
4	2262158-20	1.000	No	OVERLAY, SCIM CT, NO TILT
5	2262158-24	1.000	No	OVERLAY, SCIM CT, NO TILT

Table 5-37 SCIM Collector - International (2275752-4 BOM, rev 0)