



GE Medical Systems

gehealthcare.com

Technical Publication

Direction 2394950-100

Revision 04

**GE Medical Systems
GOC3 Illustrated Parts List**

Copyright © 2003, 2004, 2005 by General Electric Company.
All rights reserved.



LEGAL NOTES

TRADEMARKS

All products and their name brands are trademarks of their respective holders.

COPYRIGHTS

All Material, Copyright © 2003, 2004, 2005 by General Electric Company. All rights reserved.

IMPORTANT PRECAUTIONS

LANGUAGE

WARNING

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

AVERTISSEMENT

- CE MANUEL DE MAINTENANCE N'EST DISPONIBLE QU'EN ANGLAIS.
- SI LE TECHNICIEN DU CLIENT A BESOIN DE CE MANUEL DANS UNE AUTRE LANGUE QUE L'ANGLAIS, C'EST AU CLIENT QU'IL INCOMBE DE LE FAIRE TRADUIRE.
- NE PAS TENTER D'INTERVENTION SUR LES ÉQUIPEMENTS TANT QUE LE MANUEL SERVICE N'A PAS ÉTÉ CONSULTÉ ET COMPRIS.
- LE NON-RESPECT DE CET AVERTISSEMENT PEUT ENTRAÎNER CHEZ LE TECHNICIEN, L'OPÉRATEUR OU LE PATIENT DES BLESSURES DUES À DES DANGERS ÉLECTRIQUES, MÉCANIQUES OU AUTRES.

WARNUNG

- DIESES KUNDENDIENST-HANDBUCH EXISTIERT NUR IN ENGLISCHER SPRACHE.
- FALLS EIN FREMDER KUNDENDIENST EINE ANDERE SPRACHE BENÖTIGT, IST ES AUFGABE DES KUNDEN FÜR EINE ENTSPRECHENDE ÜBERSETZUNG ZU SORGEN.
- VERSUCHEN SIE NICHT, DAS GERÄT ZU REPARIEREN, BEVOR DIESES KUNDENDIENST-HANDBUCH ZU RATE GEZOGEN UND VERSTANDEN WURDE.
- WIRD DIESE WARNUNG NICHT BEACHTET, SO KANN ES ZU VERLETZUNGEN DES KUNDENDIENSTTECHNIKERS, DES BEDIENERS ODER DES PATIENTEN DURCH ELEKTRISCHE SCHLÄGE, MECHANISCHE ODER SONSTIGE GEFAHREN KOMMEN.

LANGUAGE (CONTINUED)

AVISO

- ESTE MANUAL DE SERVICIO SÓLO EXISTE EN INGLÉS.
- SI ALGÚN PROVEEDOR DE SERVICIOS AJENO A GEMS SOLICITA UN IDIOMA QUE NO SEA EL INGLÉS, ES RESPONSABILIDAD DEL CLIENTE OFRECER UN SERVICIO DE TRADUCCIÓN.
- NO SE DEBERÁ DAR SERVICIO TÉCNICO AL EQUIPO, SIN HABER CONSULTADO Y COMPRENDIDO ESTE MANUAL DE SERVICIO.
- LA NO OBSERVANCIA DEL PRESENTE AVISO PUEDE DAR LUGAR A QUE EL PROVEEDOR DE SERVICIOS, EL OPERADOR O EL PACIENTE SUFRAN LESIONES PROVOCADAS POR CAUSAS ELÉCTRICAS, MECÁNICAS O DE OTRA NATURALEZA.

ATENÇÃO

- ESTE MANUAL DE ASSISTÊNCIA TÉCNICA SÓ SE ENCONTRA DISPONÍVEL EM INGLÊS.
- SE QUALQUER OUTRO SERVIÇO DE ASSISTÊNCIA TÉCNICA, QUE NÃO A GEMS, SOLICITAR ESTES MANUAIS NOUTRO IDIOMA, É DA RESPONSABILIDADE DO CLIENTE FORNECER OS SERVIÇOS DE TRADUÇÃO.
- NÃO TENHA TENTAR REPARAR O EQUIPAMENTO SEM TER CONSULTADO E COMPREENDIDO ESTE MANUAL DE ASSISTÊNCIA TÉCNICA.
- O NÃO CUMPRIMENTO DESTE AVISO PODE POR EM PERIGO A SEGURANÇA DO TÉCNICO, OPERADOR OU PACIENTE DEVIDO A CHOQUES ELÉTRICOS, MECÂNICOS OU OUTROS.

AVVERTENZA

- IL PRESENTE MANUALE DI MANUTENZIONE È DISPONIBILE SOLTANTO IN INGLESE.
- SE UN ADDETTO ALLA MANUTENZIONE ESTERNO ALLA GEMS RICHIEDE IL MANUALE IN UNA LINGUA DIVERSA, IL CLIENTE È TENUTO A PROVVEDERE DIRETTAMENTE ALLA TRADUZIONE.
- SI PROCEDA ALLA MANUTENZIONE DELL'APPARECCHIATURA SOLO DOPO AVER CONSULTATO IL PRESENTE MANUALE ED AVERNE COMPRESO IL CONTENUTO.
- NON TENERE CONTO DELLA PRESENTE AVVERTENZA POTREBBE FAR COMPIERE OPERAZIONI DA CUI DERIVINO LESIONI ALL'ADDETTO ALLA MANUTENZIONE, ALL'UTILIZZATORE ED AL PAZIENTE PER FOLGORAZIONE ELETTRICA, PER URTI MECCANICI OD ALTRI RISCHI.

LANGUAGE (CONTINUED)

警告

このサービスマニュアルには英語版しかありません。
GEMS以外でサービスを担当される業者が英語以外の言語を要求される場合、翻訳作業はその業者の責任で行うものとさせていただきます。
このサービスマニュアルを熟読し理解せずに、装置のサービスを行わないで下さい。
この警告に従わない場合、サービスを担当される方、操作員あるいは患者さんが、感電や機械的又はその他の危険により負傷する可能性があります。

注意:

本维修手册仅存有英文本。
非 GEMS 公司的维修员要求非英文本的维修手册时，
客户需自行负责翻译。
未详细阅读和完全了解本手册之前，不得进行维修。
忽略本注意事项会对维修员，操作员或病人造成触电，机械伤害或其他伤害。

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://us44hdd21/sctq/InstallFulfill/InstalFulfillment.htm>
- Contact your local service coordinator for more information on this process.

Rev. Jan. 5, 2005

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

LITHIUM BATTERY CAUTIONARY STATEMENTS



CAUTION
Risk of
Explosion

Danger of explosion if battery is incorrectly replaced. Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.

OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GEMS PQR Process to report all omissions, errors, and defects in this publication.

Revision History

Revision	Date	Reason for change
0	09/24/03	Initial Release for GRE Xstream Console Upgrade.
1	10/08/03	Added P/Ns for 4- & 8-slice DARCs
2	01/06/04	Added drawings for 2378862 and 2362591.
3	10/22/04	EMC Ed2 Updates. New LCD Monitor.
4	4/18/05	PQR 13028473.

Table of Contents


GOC3 (GRE Xstream Console)	13
GOC3 Console BOM Tree	13
1 GOC3 (2378862)	14
GRE Console (2378862, sheet 1, rev 04)	14
GOC3 (2378862, sheet 2, rev 02)	15
GOC3 (2378862, sheet 3, rev 04)	16
GOC3 16-Slice MDAS (2378862 BOM, rev 3)	17
2 Console Chassis Assembly (2362591)	20
Console Chassis Assembly (2362591 sheet 1, rev 1)	20
Console Chassis Assembly (2362591 sheet 2, rev 1)	21
Console Chassis Assembly (2362591 sheet 3, rev 1)	22
Console Chassis Assembly (2362591 BOM, rev 1)	23
2.1 Assembly - Outlet Box (2174159-3)	26
Assembly - Outlet Box (2174159-3, sheet 1, rev 1)	26
Outlet Box Assembly (2174159-3 BOM, rev 1)	27
3 SCSI Tower, DVD-RAM, MOD (2368838)	28
SCSI Tower, DVD-RAM, MOD (2368838, sheet 1, rev 0)	28
SCSI Tower, DVD-RAM, MOD (2368838 BOM, rev 0)	29
4 Associated Items	31
4.1 Monitors	32
Monitor Selection Flowchart	32
4.1.1 – Color CRT Monitors (B7850MM)	33
Console Display Monitors (B7850MM BOM, rev 1)	33
4.1.2 – LCD Monitor & Cable (2309592)	33
LCD Monitor & Cable (2309592 BOM, rev 1)	33
4.1.3 – LCD Monitor - 19 Inch (B7858LC)	33
LCD Monitor & Cable (2309592 BOM, rev 1)	33
4.2 Keyboard w/SCIM, Collectors	34
Keyboard & SCIM Assembly (Rev 11/15/00)	34

4.2.1 – Keyboard w/SCIM Collector - English (2275758)	35
Keyboard w/SCIM Collector - English (2275758 BOM, rev 1)	35
4.2.1.1 – SCIM Collector – English (2275752)	35
SCIM Collector - English (2275752 BOM, rev 0)	35
4.2.2 – French (2275758-2)	36
Keyboard w/SCIM Collector - French (2275758-2 BOM, rev 1)	36
4.2.2.1 – SCIM Collector - French (2275752-2)	36
SCIM Collector - French (2275752-2 BOM, rev 0)	36
4.2.3 – German (2275758-3)	37
Keyboard w/SCIM Collector - German (2275758-3 BOM, rev 1)	37
4.2.3.1 – SCIM Collector - German (2275752-3)	37
SCIM Collector - German (2275752-3 BOM, rev 0)	37
4.2.4 – Swedish (2275758-4)	37
Keyboard w/SCIM Collector - Swedish (2275758-4 BOM, rev 1)	37
4.2.5 – International (2275758-5)	38
Keyboard w/SCIM Collector - International (2275758-5 BOM, rev 1)	38
4.2.5.1 – SCIM Collector - International (2275752-4)	38
SCIM Collector - International (2275752-4 BOM, rev 0)	38
4.3 GOC3, Hardware Kit (2379133)	39
GOC3, Hardware Kit (2379133 BOM, rev 0)	39
4.4 GOC3 Cable Kit (2384112)	41
GOC3 Cable Kit (2384112 BOM, rev 0)	41
4.5 GOC3, SCSI Cable Kit (2379129)	41
GOC3, SCSI Cable Kit (2379129 BOM, rev 0)	41
4.6 GOC3, Ethernet (2379130)	41
GOC3, Ethernet (2379130 BOM, rev 0)	41
4.7 GOC3, Video/Power Cables (2379131)	42
GOC3, Video/Power Cables (2379131 BOM, rev 0)	42
4.8 Console Slide Service Kit (2379135)	42
Console Slide Service Kit (2379135 BOM, rev 0)	42
4.9 Power Switch & Cables (2379134)	42
Power Switch & Cables (2379134 BOM, rev 0)	42

4.10	GOC3, Electrical Component (2379132)	43
	GOC3, Electrical Component (2379132 BOM, rev 0)	43
4.11	GOC3 Modem Kit (B7700MG)	43
	GOC3 Modem Kit (B7700MG BOM, rev 0)	43
5	Options	44
5.1	Video Upgrade For Linux Console (B7868MC)	44
	Video Upgrade For Linux Console (B7868MC BOM, rev 0)	44
	5.1.1 – Video Upgrade (2378278)	44
	Video Upgrade (2378278 BOM, rev 0)	44
5.2	LCD & SUSP (Boom In Room) LNX (B7710LM)	45
	LCD & SUSP (Boom In Room) LNX (B7710LM BOM, rev 0)	45
	5.2.1 – Boom-in-a-Room for Linux (2376773)	45
	Boom-in-a-Room for Linux (2376773 BOM, rev 0)	45
5.3	HS/LS Remote Color (LNX) Monitor (B7530LC)	46
	HS/LS Remote Color (LNX) Monitor (B7530LC BOM, rev 0)	46
	5.3.1 – Remote Monitor Kit (2377146)	46
	Remote Monitor Kit (2377146 BOM, rev 0)	46

GOC3 (GRE Xstream Console)

Figure 1: GOC3 Console BOM Tree



1 GOC3 (2378862)

FRU = NO

Figure 2: GRE Console (2378862, sheet 1, rev 04)




Figure 3: GOC3 (2378862, sheet 2, rev 02)




Figure 4: GOC3 (2378862, sheet 3, rev 04)




Table 1: GOC3 16-Slice MDAS (2378862 BOM, rev 3)

MARK	NUMBER	QTY	FRU	DESCRIPTION
-	2249895	1	Yes	Digital DASM
-	2379129	-	Yes	GRE CONSOLE, SCSI CABLE KIT
-	2379130	-	Yes	GRE CONSOLE, ETHERNET /
-	2379131	-	Yes	GRE CONSOLE, VIDEO/POWER CABLES
-	2379132	-	Yes	GRE CONSOLE, ELECTRICAL COMPONENT
-	2379133	-	Yes	GRE CONSOLE, HARDWARE KIT
-	2379134	-	Yes	GRE CONSOLE, PWR SWITCH & CABLES
-	2379135	-	Yes	GRE CONSOLE SLIDE SERVICE KIT
-	2382274	-	Yes	150' OPTION CABLE FOR REMOTE LCD
-	B7530LC	-	Yes	HS/LS REMOTE COLOR (LNX) MONITOR
-	B7710LM	-	Yes	LCD & SUSP (BOOM IN ROOM) LNX
-	B7868MC	-	Yes	VIDEO UPGRADE FOR LINUX CONSOLE
2	2368697-3	1.000	Yes	Hewlett-Packard XW8000 PC (CT GRE Only)
3	2362870-10 2362870-12* 2362870-13**	1.000	Yes	DARC - GRE (16 SLICE MDAS GDIP) DARC - GRE (4/8 SLICE MDAS GDIP)* DARC - GRE (4 SLICE SDAS GDIP)**
4	2362872	1.000***	Yes	IMAGE GENERATION COMPUTER - GRE
5	2362873	1.000	Yes	SCAN DATA DISK ARRAY - GR
6	2368838	1.000	Yes	SCSI TOWER, DVD-RAM, MOD {For further detail, see Section 3. }
8	2368697-2	2.000	Yes	LINUX DDR WORKST'N 3BUTTO
9	2399753	1.000	Yes	TRACKBALL, USB
11	2369619	1.000	Yes	COVER, CONSOLE FRONT - GR
12	2337534	1.000	Yes	SIDE COVERS GOC
13	2304737-2	2.000	Yes	CABLE, DVI VIDEO CABLE -

* Used on GOC3 4/8 MDAS (2378862-4, 2378862-5 or 2378862-6)

** Used on GOC3 4-Slice SDAS (2378862-7, 2378862-8 or 2378862-9)

*** 2378862-2, 2378862-5 and 2378862-8 use 2 IGs; 2378862-3, 2378862-6 and 2378862-9 use 3 IGs

Table 1: GOC3 16-Slice MDAS (2378862 BOM, rev 3) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION
15	2377227	1.000	Yes	UTP TRANSCEIVER - BNC TO
16	2323942	1.000	Yes	CABLE, HD50 TO HD50
17	2324697-2	1.000	Yes	SCSI CABLE, VHD68 TO CENT
18	2374965	1.000	Yes	CABLE - ETHERNET CAT5E X
19	2374965-2	2.000	Yes	CABLE - ETHERNET CAT5E X
20	2376295	1.000	Yes	CABLE, SCSI U320 HD68M/M
J20	2376296	1.0	Yes	Cable Data DB25 M/F SHLD
21	2377380	2.000	Yes	CABLE - AUDIO 3.5MM JACK
22	2377597	1.000	Yes	CABLE - DB9F TO DB9F X 30
23	2383458	1.000	Yes	CABLE - SERIAL RHARD (GRE
1	2362591	1.000	No	CONSOLE CHASSIS ASSY - GR {For further detail, see Section 2. }
14	2366442-27	2.000	No	BRACKET, RACK MOUNT HOST
24	2335376	1.000	No	CABLE - PDU TO CONSOLE FO
25	46-136246P14	1.000 ft	No	2" WIDE, WHT, ADH BACK, V
26	46-136246P21	1.000 ft	No	2" WIDE X 25 YD WHITE VEL
27	46-208781P2	4.000	No	WRITE-ON LABEL FOR MARKIN
28	2383464	1.000	No	LABEL ETL MARK - ITS 8352
29	2117390	1.000	No	LABEL, CE MARKING COMPLIA
30	2132620	1.000	No	CE MARKING LABEL CLASS A
31	46-279769P1	1.000	No	IEC TYPE B SYMBOL
32	46-269337P6	1.000	No	LABEL OLD
33	2240540	1.000	No	NETSCAPE LICENSE
34	46-269568P1	1.000	No	INFORMIX SOFTWARE-DATA BA
35	2214164	1.000	No	5.25" 2.3GB 512 BYTES/SEC

* Used on GOC3 4/8 MDAS (2378862-4, 2378862-5 or 2378862-6)

** Used on GOC3 4-Slice SDAS (2378862-7, 2378862-8 or 2378862-9)

*** 2378862-2, 2378862-5 and 2378862-8 use 2 IGs; 2378862-3, 2378862-6 and 2378862-9 use 3 IGs

Table 1: GOC3 16-Slice MDAS (2378862 BOM, rev 3) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION
36	2337414	2.000	No	DVD RAM MEDIA, 9.4 GB
37	2227726	1.000	No	MOUSE PAD FOR LIGHTSPEED
39	46-208758P3	6.000	No	7.31" X .184" SELF-LOCKIN
41	46-312358P51	14.000	No	M5 X 0.8 X 12MM METRIC FR
42	2103585-2	4.000	No	SCREW PAN HEAD 4 MM 10 MM
44	46-221856P13	AR	No	TAPE, ELECTRICAL, PVC, CO
45	46-221856P14	AR	No	TAPE, ELECTRICAL, PVC, CO
46	5123478	1.000	No	ADAPTER, USB TO USB
				ASSOCIATED OPTIONS {For further detail, see Section 5 .}

* Used on GOC3 4/8 MDAS (2378862-4, 2378862-5 or 2378862-6)

** Used on GOC3 4-Slice SDAS (2378862-7, 2378862-8 or 2378862-9)

*** 2378862-2, 2378862-5 and 2378862-8 use 2 IGs; 2378862-3, 2378862-6 and 2378862-9 use 3 IGs

FRU = YES

Figure 5: Console Chassis Assembly (2362591 sheet 1, rev 1)





Figure 6: Console Chassis Assembly (2362591 sheet 2, rev 1)



REVISIONS

REV	DATE	BY	CHK	DESCRIPTION
1	220981			ADDED K, L & M POWER CABLES. PERTINENT INFO IN EACH VIEW. ADDED CABLE LENGTHS FOR CONNECTION LOCATIONS. ADDED VIEW OF POWER SWITCH AT D15 FOR CLARITY.

FIG. 1

FIG. 2

FIG. 3

FIG. 4

FIG. 5

FIG. 6

FIG. 7

FIG. 8

FIG. 9

FIG. 10

FIG. 11

FIG. 12

FIG. 13

FIG. 14

FIG. 15

FIG. 16

FIG. 17

FIG. 18

FIG. 19

FIG. 20

FIG. 21

FIG. 22

FIG. 23

FIG. 24

FIG. 25

FIG. 26

FIG. 27

FIG. 28

FIG. 29

FIG. 30

FIG. 31

FIG. 32

FIG. 33

FIG. 34

FIG. 35

FIG. 36

FIG. 37

FIG. 38

FIG. 39

FIG. 40

FIG. 41

FIG. 42

FIG. 43

FIG. 44

FIG. 45

FIG. 46

FIG. 47

FIG. 48

FIG. 49

FIG. 50

FIG. 51

FIG. 52

FIG. 53

FIG. 54

FIG. 55

FIG. 56

FIG. 57

FIG. 58

FIG. 59

FIG. 60

FIG. 61

FIG. 62

FIG. 63

FIG. 64

FIG. 65

FIG. 66

FIG. 67

FIG. 68

FIG. 69

FIG. 70

FIG. 71

FIG. 72

FIG. 73

FIG. 74

FIG. 75

FIG. 76

FIG. 77

FIG. 78

FIG. 79

FIG. 80

FIG. 81

FIG. 82

FIG. 83

FIG. 84

FIG. 85

FIG. 86

FIG. 87

FIG. 88

FIG. 89

FIG. 90

FIG. 91

FIG. 92

FIG. 93

FIG. 94

FIG. 95

FIG. 96

FIG. 97

FIG. 98

FIG. 99

FIG. 100

FIG. 101

FIG. 102

FIG. 103

FIG. 104

FIG. 105

FIG. 106

FIG. 107

FIG. 108

FIG. 109

FIG. 110

FIG. 111

FIG. 112

FIG. 113

FIG. 114

FIG. 115

FIG. 116

FIG. 117

FIG. 118

FIG. 119

FIG. 120

FIG. 121

FIG. 122

FIG. 123

FIG. 124

FIG. 125

FIG. 126

FIG. 127

FIG. 128

FIG. 129

FIG. 130

FIG. 131

FIG. 132

FIG. 133

FIG. 134

FIG. 135

FIG. 136

FIG. 137

FIG. 138

FIG. 139

FIG. 140

FIG. 141

FIG. 142

FIG. 143

FIG. 144

FIG. 145

FIG. 146

FIG. 147

FIG. 148

FIG. 149

FIG. 150

FIG. 151

FIG. 152

FIG. 153

FIG. 154

FIG. 155

FIG. 156

FIG. 157

FIG. 158

FIG. 159

FIG. 160

FIG. 161

FIG. 162

FIG. 163

FIG. 164

FIG. 165

FIG. 166

FIG. 167

FIG. 168

FIG. 169

FIG. 170

FIG. 171

FIG. 172

FIG. 173

FIG. 174

FIG. 175

FIG. 176

FIG. 177

FIG. 178

FIG. 179

FIG. 180

FIG. 181

FIG. 182

FIG. 183

FIG. 184

FIG. 185

FIG. 186

FIG. 187

FIG. 188

FIG. 189

FIG. 190

FIG. 191

FIG. 192

FIG. 193

FIG. 194

FIG. 195

FIG. 196

FIG. 197

FIG. 198

FIG. 199

FIG. 200

FIG. 201

FIG. 202

FIG. 203

FIG. 204

FIG. 205

FIG. 206

FIG. 207

FIG. 208

FIG. 209

FIG. 210

FIG. 211

FIG. 212

FIG. 213

FIG. 214

FIG. 215

FIG. 216

FIG. 217

FIG. 218

FIG. 219

FIG. 220

FIG. 221

FIG. 222

FIG. 223

FIG. 224

FIG. 225

FIG. 226

FIG. 227

FIG. 228

FIG. 229

FIG. 230

FIG. 231

FIG. 232

FIG. 233

FIG. 234

FIG. 235

FIG. 236

FIG. 237

FIG. 238

FIG. 239

FIG. 240

FIG. 241

FIG. 242

FIG. 243

FIG. 244

FIG. 245

FIG. 246

FIG. 247

FIG. 248

FIG. 249

Table 2: Console Chassis Assembly (2362591 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION
4	2367805	1.000	Yes	CABLE BULKHEAD ASSY (GRE)
9	2142225	4.000	Yes	MONITOR POWER CORD
10	2142226	5.000	Yes	POWER CORD FOR EMC AREA
22	2180698-2	40.000	Yes	NUT CLIP ON 5 MM
35	2376909	2.000	Yes	FAN ASSY - DUAL 120MM FAN
36	2174159-3	1.000	Yes	ASSY - OUTLET BOX (GRE)
38	2315657	1.000	Yes	TABLE TOP ASSEMBLY, MONIT
39	2330566	1.000	Yes	TABLE TOP ASSEMBLY, KEYBO
1	2366442	1.000	No	CONSOLE CHASSIS SUB-ASSY,
2	2366442-11	1.000	No	REAR MTG PANEL, OUTLET BO
3	2366442-14	1.000	No	COVER - CABLE
5	2315117	1.000	No	CABLE, PS/2 MOUSE CABLE,
6	2315117-2	1.000	No	CABLE, PS/2 KEYBOARD CABL
7	5123478	1.000	No	CABLE, USB TO USB
8	46-328421P3	6.000	No	SCREW PAN HEAD 2.5 MM 5 M
13	2318681	2.000	No	ANTI-TIP LEG, CONSOLE CHA
14	2319043	4.000	No	CASTER, 2 1/2" SWIVEL
15	2103585	38.000	No	SCREW PAN HEAD 4 MM 8 MM
16	2325236	2.000	No	LEVELER, CONSOLE ANTI-TIP
17	2109867-25	4.000	No	SCREW HEXAGON SOCKET
18	2109877-2	16.000	No	NUT HEXAGON 6 MM
19	2109878-4	4.000	No	WASHER FLAT ISO7089 M10 A
20	2109881-3	25.000	No	WASHER LOCK - INTERNAL 6.
21	2109881-4	4.000	No	WASHER LOCK - INTERNAL 10
23	2366442-12	4.000	No	BRACKET, SIDE COVER MTG
24	2366442-13	3.000	No	BRACKET, FRONT COVER MTG
25	2366442-29	1.000	No	BRACKET - FAN LEFT

Table 2: Console Chassis Assembly (2362591 BOM, rev 1) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION
26	2366442-30	1.000	No	BRACKET - FAN RIGHT
27	2366442-6	2.000	No	PANEL, RACK MOUNT
28	2366442-7	10.000	No	BRACKET, RACK MOUNT
29	2366442-8	1.000	No	REAR ACCESS COVER, CONSOL
31	46-328431P2	25.000	No	WASHER PLAIN - LARGE 6.4
32	46-170686P3	AR	No	LOCTITE 242, 50CC
33	2376906	6.000	No	SOUTHCO 1/4 TURN RCVR 85-
34	2337083-2	9.000	No	SCREW, M6 X 20 SOC HEAD,
37	2366442-31	3.000	No	PANEL - COVER
40	46-312358P51	26.000	No	M5 X 0.8 X 12MM METRIC FR
41	46-208758P3	19.000	No	7.31" X .184" SELF-LOCKIN

This Page Intentionally Blank.

FRU = YES




Table 3: Outlet Box Assembly (2174159-3 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION
14	46-208758P2	3.000	Yes	5.5" X .140" SELF-LOCKING
1	2367803	1.000	No	WELDMENT - OUTLET BOX (GR
2	46-278446P1	7.000	No	15A, 125VAC, 2 POLE, 3 WI
3	46-297335P3	1.000	No	16A, 250V 2-STAGE PWR LIN
4	46-243925P1	1.000	No	20A, 125V, 2-POLE, 3 WIRE
6	46-328425P10	16.000	No	NUT HEXAGON 3 MM
7	46-328425P1	3.000	No	NUT HEXAGON 4 MM
8	46-312358P30	3.000	No	M3 X 0.5 X 12MM METRIC FR
9	46-328430P2	3.000	No	WASHER PLAIN - NORMAL 4.3
11	46-208697P57	1.000	No	HEYCO PART NO.1247 SR-7W
12	46-311804P2	16.000	No	3.2X9X.8 ZP ST DIN
13	46-311805P10	16.000	No	HELICAL LOCK WASHER
15	2336074	1.000	No	CIRCUIT BREAKER 15 AMPS 2
16	46-328421P501	2.000	No	SCREW PAN HEAD 3 MM 5 MM
17	2367804	1.000	No	SWITCH BOX ASSY (GRE)

3 SCSI Tower, DVD-RAM, MOD (2368838)

FRU = YES

Figure 9: SCSI Tower, DVD-RAM, MOD (2368838, sheet 1, rev 0)




Table 4: SCSI Tower, DVD-RAM, MOD (2368838 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2319046	1.000	No	DVDRAM PLAYER
2	2318598	1.000	No	SONY MOD SMO-F551-SD
4	2315118	1.000	No	BAY EXTERNAL SCSI TOWER
5	46-328421P501	8.000	No	SCREW PAN HEAD 3 MM 5 MM
8	2350404	1.000	No	TERMINATOR, CENT50, ACTIV


This Page Intentionally Blank.

4 Associated Items

- [Color CRT Monitors \(B7850MM\)](#)
- [LCD Monitor & Cable \(2309592\)](#)
- [Keyboard w/SCIM, Collectors](#)
- [GOC3, Hardware Kit \(2379133\)](#)
- [GOC3 Cable Kit \(2384112\)](#)
- [GOC3, SCSI Cable Kit \(2379129\)](#)
- [GOC3, Ethernet \(2379130\)](#)
- [GOC3, Video/Power Cables \(2379131\)](#)
- [Console Slide Service Kit \(2379135\)](#)
- [Power Switch & Cables \(2379134\)](#)
- [GOC3, Electrical Component \(2379132\)](#)
- [GOC3 Modem Kit \(B7700MG\)](#)

4.1 Monitors

Figure 10: Monitor Selection Flowchart



4.1.1 Color CRT Monitors (B7850MM)

Table 5: Console Display Monitors (B7850MM BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2333243	2.000	No	21 IN COLOR CRT MONITOR
1A	2254212	ALT	No	SONY GDM-5402T9 21" COLOR MONITOR

4.1.2 LCD Monitor & Cable (2309592)

FRU = NO

Table 6: LCD Monitor & Cable (2309592 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2349792	1.000	Yes	18.1 INCH LCD COLOR MONITOR
2	2295019	1.000	Yes	13W3-BD15 RGB LCD STRAIGHT CONN CBL

4.1.3 LCD Monitor - 19 Inch (B7858LC)

Table 10-1 LCD Monitor & Cable (2309592 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	5117545	2.000	Yes	Monitor, 19 inch color LCD, silver, NEC LCD1980SXi

4.2 Keyboard w/SCIM, Collectors

Figure 11: Keyboard & SCIM Assembly (Rev 11/15/00)



4.2.1 Keyboard w/SCIM Collector - English (2275758)

FRU = YES

Table 7: Keyboard w/SCIM Collector - English (2275758 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2275756	1.000	Yes	KEYBOARD, NORTH AMERICAN
3	2275752	1.000	Yes	SCIM COLLECTOR - ENGLISH {For further detail, see Section 4.2.1.1 , below.}
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757	1.000	No	FUNCTION KEY TEMPLATE - E

4.2.1.1 SCIM Collector – English (2275752)

FRU = YES

Table 8: SCIM Collector - English (2275752 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2262158	1.000	No	SCAN CONTROL INTERCOM MOD
2	2262158-11	1.000	No	OVERLAY, SCIM CT, TILT -
3	2262158-21	1.000	No	OVERLAY, SCIM CT, NO TILT

4.2.2 French (2275758-2)

Table 9: Keyboard w/SCIM Collector - French (2275758-2 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2275756-2	1.000	Yes	KEYBOARD, EUROPEAN FRENCH
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757-2	1.000	No	FUNCTION KEY TEMPLATE - F
3	2275752-2	1.000	No	SCIM COLLECTOR - FRENCH {For further detail, see Section 4.2.2.1 , below.}

4.2.2.1 SCIM Collector - French (2275752-2)

Table 10: SCIM Collector - French (2275752-2 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2262158	1.000	No	SCAN CONTROL INTERCOM MOD
2	2262158-12	1.000	No	OVERLAY, SCIM CT, TILT -
3	2262158-22	1.000	No	OVERLAY, SCIM CT, NO TILT

4.2.3 German (2275758-3)

Table 11: Keyboard w/SCIM Collector - German (2275758-3 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2275756-3	1.000	Yes	KEYBOARD, GERMAN PC PS/2
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757-3	1.000	No	FUNCTION KEY TEMPLATE - G
3	2275752-3	1.000	No	SCIM COLLECTOR - GERMAN {For further detail, see Section 4.2.3.1 , below.}

4.2.3.1 SCIM Collector - German (2275752-3)

Table 12: SCIM Collector - German (2275752-3 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2262158	1.000	No	SCAN CONTROL INTERCOM MOD
2	2262158-13	1.000	No	OVERLAY, SCIM CT, TILT -
3	2262158-23	1.000	No	OVERLAY, SCIM CT, NO TILT

4.2.4 Swedish (2275758-4)

Table 13: Keyboard w/SCIM Collector - Swedish (2275758-4 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2275756-4	1.000	Yes	KEYBOARD, SWEDISH PC PS/2
3	2275752	1.000	Yes	SCIM COLLECTOR - ENGLISH {For further detail, see Section 4.2.1.1 , beginning on page 35 .}
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2275757	1.000	No	FUNCTION KEY TEMPLATE - E

4.2.5 International (2275758-5)

Table 14: Keyboard w/SCIM Collector - International (2275758-5 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2275756	1.000	Yes	KEYBOARD, NORTH AMERICAN
4	2323942	1.000	Yes	CABLE, HD50 TO HD50
2	2283830	1.000	No	FUNCTION KEY TEMPLATE COL
3	2275752-4	1.000	No	SCIM COLLECTOR - INTERNATIONAL {For further detail, see Section 4.2.5.1 , below.}

4.2.5.1 SCIM Collector - International (2275752-4)

Table 15: SCIM Collector - International (2275752-4 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2262158	1.000	No	SCAN CONTROL INTERCOM MOD
2	2262158-10	1.000	No	OVERLAY, SCIM CT, TILT -
3	2262158-14	1.000	No	OVERLAY, SCIM CT, TILT -
4	2262158-20	1.000	No	OVERLAY, SCIM CT, NO TILT
5	2262158-24	1.000	No	OVERLAY, SCIM CT, NO TILT

4.3 GOC3, Hardware Kit (2379133)

Table 16: GOC3, Hardware Kit (2379133 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	2103585	8.000	No	SCREW PAN HEAD 4 MM 8 MM
3	2142006	8.000	No	HEX SOC SCREW M2,5 X 16
4	2260836	8.000	No	SCREW PAN HEAD 3 MM 6 MM
6	2324209	8.000	No	SCREW SOCKET HEAD M4 X 60
7	2325236	2.000	No	LEVELER, CONSOLE ANTI-TIP
8	2331393	8.000	No	SCREW, POLYTOP SCREW, M4
10	2103585-2	8.000	No	SCREW PAN HEAD 4 MM 10 MM
11	2109867-25	8.000	No	SCREW HEXAGON SOCKET
12	2109873-11	8.000	No	SCREW PAN HEAD 4 MM X 5 M
13	2109877-2	8.000	No	NUT HEXAGON 6 MM
14	2109878-4	4.000	No	WASHER FLAT ISO7089 M10 A
15	2109881-3	4.000	No	WASHER LOCK - INTERNAL 6.
16	2109881-4	4.000	No	WASHER LOCK - INTERNAL 10
18	46-170015P14	8.000	No	PAN HD SCREW 8-32 X 3/8
19	46-170015P23	4.000	No	SCREW, 1/4"-20 X 1/2"
20	46-208935P8	4.000	No	1/4-20, 7/16" HEX X 3/16
26	46-287187P4	8.000	No	SCREW, CSNK 2.5MM X 6LG
27	46-311804P2	8.000	No	3.2X9X.8 ZP ST DIN
28	46-311805P10	8.000	No	HELICAL LOCK WASHER
29	46-312346P2	8.000	No	NUT, M4 X 0.7 HEX WITH LO
30	46-328421P11	8.000	No	SCREW PAN HEAD 4 MM 5 MM
31	46-328421P3	8.000	No	SCREW PAN HEAD 2.5 MM 5 M
32	46-328421P4	8.000	No	SCREW PAN HEAD 2.5 MM 6 M
33	46-328421P501	8.000	No	SCREW PAN HEAD 3 MM 5 MM
35	46-328431P1	8.000	No	WASHER PLAIN - LARGE 4.3
36	46-328431P2	8.000	No	WASHER PLAIN - LARGE 6.4

Table 16: GOC3, Hardware Kit (2379133 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION
37	46-328432P1	8.000	No	WASHER LOCK - SPRING 2.6
38	46-312679P2	4.000	No	INTERNALLY RIBBED GROMMET
39	46-208695P61	2.000	No	HEYCO PART NO.2126SB-87
40	46-170686P3	1.000	No	LOCTITE 242, 50CC
41	46-208897P16	8.000	No	SCREW
42	46-312358P51	8.000	No	M5 X 0.8 X 12MM METRIC FR
43	2337083-2	8.000	No	SCREW, M6 X 20 SOC HEAD,
44	2180698-2	8.000	No	NUT CLIP ON 5 MM
45	46-328425P10	8.000	No	NUT HEXAGON 3 MM
46	46-328425P1	4.000	No	NUT HEXAGON 4 MM
47	46-328430P2	4.000	No	WASHER PLAIN - NORMAL 4.3
48	46-312358P30	4.000	No	M3 X 0.5 X 12MM METRIC FR
49	46-208758P2	4.000	No	5.5" X .140" SELF-LOCKING

4.4 GOC3 Cable Kit (2384112)

Table 17: GOC3 Cable Kit (2384112 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	2323945	1.000	No	CABLE - MODEM TO HOST COM
2	2377597	1.000	No	CABLE - DB9F TO DB9F X 30
3	2383458	1.000	No	CABLE - SERIAL RHARD (GRE
4	2377229	1.000	No	CABLE - COAX X 610 LG
5	2368504	1.000	No	CONVERTER, DB25 TO USB

4.5 GOC3, SCSI Cable Kit (2379129)

Table 18: GOC3, SCSI Cable Kit (2379129 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2324697-2	1.000	Yes	SCSI CABLE, VHD68 TO CENT 50, LINUX
2	2323942	1.000	Yes	CABLE, HD50 TO HD50
3	2376295	1.000	Yes	CABLE, SCSI U320 HD68M/M 305MM

4.6 GOC3, Ethernet (2379130)

Table 19: GOC3, Ethernet (2379130 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2374965	3.000	Yes	CABLE - ETHERNET CAT5E X 305MM LG
2	2374965-2	2.000	Yes	CABLE - ETHERNET CAT5E X 610MM LG
3	2318713-2	1.000	No	CABLE - FIBER OPTIC ST-SC X 3FT
	2379130			COMMUNICATION CABLE KIT

4.7 GOC3, Video/Power Cables (2379131)

Table 20: GOC3, Video/Power Cables (2379131 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2142225	2.000	No	MONITOR POWER CORD
2	2142226	2.000	No	POWER CORD FOR EMC AREA
3	2304737-2	2.000	No	CABLE, DVI VIDEO CABLE -

4.8 Console Slide Service Kit (2379135)

Table 21: Console Slide Service Kit (2379135 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2317595	1.000	No	SLIDE RAIL FOR DRAWER ASS
4	2317501	1.000	No	BRACKET, DRAWER BALL LATC
5	2317502	1.000	No	LATCH, DRAWER ASSEMBLY
6	2366442-7	4.000	No	BRACKET, RACK MOUNT

4.9 Power Switch & Cables (2379134)

Table 22: Power Switch & Cables (2379134 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2150833	1.000	No	BREAKER SWITCH
2	2108900	1.000	No	CONNECTOR 8 CONTACTS 150 VOLTS
5	2318308	1.000	No	ETHERNET CABLE, CHASSIS SVC PANEL
7	2341434-2	1.000	No	CABLE, KEYBOARD - OCTANE 2
8	2341434	1.000	No	CABLE, MOUSE - OCTANE 2

4.10 GOC3, Electrical Component (2379132)

Table 23: GOC3, Electrical Component (2379132 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
4	46-297335P3	1.000	Yes	16A, 250V 2-STAGE PWR LIN
5	2350404	1.000	Yes	TERMINATOR, CENT50, ACTIV
3	2154019	2.000	No	50 OHM BNC GROUNDED TERMI
6	46-243925P1	1.000	No	20A, 125V, 2-POLE, 3 WIRE
7	46-208697P57	1.000	No	HEYCO PART NO.1247 SR-7W
8	2336074	1.000	No	CIRCUIT BREAKER 15 AMPS 2

4.11 GOC3 Modem Kit (B7700MG)

Table 24: GOC3 Modem Kit (B7700MG BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	2382276	1.000	No	GRE MODEM KIT

5 Options

- [Video Upgrade For Linux Console \(B7868MC\)](#)
- [LCD & SUSP \(Boom In Room\) LNX \(B7710LM\)](#)
- [HS/LS Remote Color \(LNX\) Monitor \(B7530LC\)](#)

5.1 Video Upgrade For Linux Console (B7868MC)

Table 25: Video Upgrade For Linux Console (B7868MC BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2378278	1.000	No	VIDEO UPGRADE FOR LINUX C {For further detail, see Section 5.1.1. }

5.1.1 Video Upgrade (2378278)

Table 26: Video Upgrade (2378278 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
9	2302528	1.000	Yes	VIDIO CABLE FOR LCD SUSPE
1	46-302200P10	1.000	No	1.72 X 1.10, MODEL#, SN,
2	2376770	2.000	No	5 BNC MALE TO HD15 MALE (
13	2376772	1.000	No	1-IN 4-OUT RGBHV DISTRIBUTU
14	2376770-2	1.000	No	5 BNC MALE TO HD15 MALE (
16	2376771	1.000	No	5 BNC MALE TO HD15 FEMALE
18	2200583-100	1.000	No	CT/I (SONATA) UPGRADE MAN
19	2211224-100	1.000	No	LIGHTSPEED QX/I OPTIONS I
25	2212092	1.000	No	CEILING PLATE ASSEMBLY FO

5.2 LCD & SUSP (Boom In Room) LNX (B7710LM)

Table 27: LCD & SUSP (Boom In Room) LNX (B7710LM BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2376773	1.000	No	BOOM-IN-A-ROOM FOR LINUX {For further detail, see Section 5.2.1. }

5.2.1 Boom-in-a-Room for Linux (2376773)

Table 28: Boom-in-a-Room for Linux (2376773 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2214293	1.000	Yes	CEILING COLUMN WITH FLANG
4	2214292	1.000	Yes	80-80 CM EXTENSION/SPRING
5	2349792	1.000	Yes	18.1 IN LCD COLOR MONITOR
6	2211031	1.000	Yes	LCD BRACKET - SHORT
7	2212717	1.000	Yes	CABLE COVER FOR LCD SUSPE
9	2302528	1.000	Yes	VIDIO CABLE FOR LCD SUSPE
11	2213213	1.000	Yes	GROUNDING CABLE FOR LCD S
12	2213218	1.000	Yes	POWER CABLE FOR LCD SUSPE
1	46-302200P10	1.000	No	1.72 X 1.10, MODEL#, SN,
2	2376770	1.000	No	5 BNC MALE TO HD15 MALE (
8	2213216	1.000	No	GROUNDING CABLE FOR LCD-C
10	2213219	1.000	No	POWER CABLE FOR LCD-CONSO
13	2376772	1.000	No	1-IN 4-OUT RGBHV DISTRIBU
14	2376770-2	1.000	No	5 BNC MALE TO HD15 MALE (
16	2376771	1.000	No	5 BNC MALE TO HD15 FEMALE
17	2223870	1.000	No	SONATA HARDWARE COLLECTOR
18	2200583-100	1.000	No	CT/I (SONATA) UPGRADE MAN
19	2211224-100	1.000	No	LIGHTSPEED QX/I OPTIONS I
25	2212092	1.000	No	CEILING PLATE ASSEMBLY FO
26	2300059	1.000	No	HANDLE FOR LCD DISPLAY MO

5.3 HS/LS Remote Color (LNX) Monitor (B7530LC)

Table 29: HS/LS Remote Color (LNX) Monitor (B7530LC BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2377146	1.000	No	REMOTE MONITOR KIT FOR LI {For further detail, see Section 5.3.1. }

5.3.1 Remote Monitor Kit (2377146)

Table 30: Remote Monitor Kit (2377146 BOM, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2333243	1.000	Yes	21 IN COLOR CRT MONITOR
1A	2254212	ALT	Yes	SONY GDM-5402T9 21" COLOR
2	2376772	1.000	No	1-IN 4-OUT RGBHV DISTRIBUTOR
3	2376770	1.000	No	5 BNC MALE TO HD15 MALE (
4	2376770-2	1.000	No	5 BNC MALE TO HD15 MALE (
5	2376771	1.000	No	5 BNC MALE TO HD15 FEMALE



GE MEDICAL SYSTEMS

**GE MEDICAL SYSTEMS-AMERICAS: FAX 262.312.7434
3000 N. GRANDVIEW BLVD., WAUKESHA, WI 53188, U.S.A.**

**GE MEDICAL SYSTEMS-EUROPE: FAX 33.1.40.93.33.33
PARIS, FRANCE**

**GE MEDICAL SYSTEMS-ASIA: FAX 65.291.7006
SINGAPORE**