

System Configuration Data Sheets

1 Requirement

Record valuable system information in the data sheets that follow. Consult with your customer or network administrator to obtain the information. Understanding how the customer plans to use their LightSpeed 4.X scanner, and their network and filming expectation reduces the time required to reconfigure the system.

- [Table 1: Manual Film Composer Options, on page 2](#)
- [Table 2: System Network Configuration, on page 3](#)
- [Table 3: Networking Application \(Image transfer\) Configuration, on page 4](#)
- [Table 4: DASM Laser Camera Configuration, on page 6](#)
- [Table 5: DICOM Print Camera Configuration, on page 6](#)
- [Table 6: DICOM Print Camera Advanced Configuration, on page 7](#)

2 Manual Film Composer Options

Table 1: Manual Film Composer Options

MANUAL FILM COMPOSER OPTIONS	
Slide Format (if available):	
Greyscale:	
Auto Printing:	
Auto Clear Page:	
Icon Labels:	
Expose Order:	
No. of Copies:	

3 System Network Configuration

Table 2: System Network Configuration

SYSTEM NETWORK CONFIGURATION			
	FIELD NAME:	SETENV NAME:	FIELD VALUE:
System Settings:	Service ID	SERVICE_ID	
	Hospital Name	HOSPITAL_NAME	
	Exam Number *	* ask customer or checklog	
	DAS Type	DASTYPE	
	PDU Type	PDUTYPE	
Network Settings:	Gateway Host Name	GATEWAY_HOSTNAME	
	Gateway IP	GATEWAY_IP	
	Gateway Net Mask	GATEWAY_NETMASK	
	Gateway Broadcast Mask	GATEWAY_BROADCAST	
Option	Suite Name	SUITEID	
	Network Printer IP Address		
Option	HIS Server IP Address		
Option	HIS Server AE Title		
Option	HIS server AE Port		
Option	CT Server AE Title		
Option	Connect Pro IP Address		

4 Network Application (Image Transfer) Configuration

Record the network application (image transfer) configuration.

Table 3: Networking Application (Image transfer) Configuration

NETWORKING APPLICATION (IMAGE TRANSFER) CONFIGURATION				
AE TITLE OR HOST NAME	NETWORK ADDRESS	NETWORK PROTOCOL	PORT NUMBER	COMMENTS

5 Host Ethernet Address

_____._____._____._____._____._____

6 Camera Application Configuration

Record the camera application configuration for the DASM or DICOM print camera.

Table 4: DASM Laser Camera Configuration

DASM LASER CAMERA CONFIGURATION	
Camera Type:	
DASM Type:	
Film Smooth/Sharp Setting:	
Options:	
Valid Film Formats:	
Default Film Formats:	

Table 5: DICOM Print Camera Configuration

DICOM PRINT CAMERA CONFIGURATION	
Camera Type:	
Host Name:	
IP Address:	
AE Title:	
TCP/IP Listen Port:	
Comments (Optional):	
Valid Film Formats:	
Default Film Formats:	
Destination:	
Orientation:	
Medium Type:	
Magnification Type:	

Table 6: DICOM Print Camera Advanced Configuration

DICOM PRINT CAMERA ADVANCE CONFIGURATION

Smoothing Type:	
Configuration:	
Minimum Density:	
Maximum Density:	
Empty Density:	
Border Density:	
Association Timeout:	
Session Timeout:	
N-Set Timeout:	
N-Action Timeout:	
N-Create Timeout:	
N-Delete Timeout:	
N-Get Timeout:	