

Chapter 2

Table

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	46-303268P1	1.0	Yes	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH
2	5139485	1.0	No	WEEE LABEL WITH BIN ONLY
3	5132810	1.0	No	ETL LABEL COMPONENT CT BEIJING
4	5194739	1.0	No	NP LITE TABLE TOP ASSY POSITIONING BSD LITE

Table 2-1 Np Lite Table Assy (5182488 BOM, Rev 10)

Section 1.0 — NP LITE TABLE TOP ASSY POSITIONING BSD LITE (5194739)


FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	5194735	1.0	No	NP LITE TABLE WITHOUT DORRY AND ACC ASSY POSITIONING BSD LITE
2	5194741	1.0	No	NP LITE ACCESSORY COLLECTOR POSITIONING BSD LITE
3	5192574	1.0	No	NP LITE SERVICE TOOL COLLECTOR POSITIONING BSD LITE
4	2219924	1.0	No	P9870FB VE TABLE DORRY ASSY

Table 2-2 Np Lite Table Top Assy Positioning Bsd Lite (5194739 BOM, Rev 2)

1.1 NP LITE TABLE WITHOUT DORRY AND ACC ASSY POSITIONING BSD LITE (5194735)

FRU = NO








Figure 2-3 NP Lite Table (5194735ADW refer to 2200193ADW Sheet 3, Rev 17)








Figure 2-6 NP Lite Table (5194735ADW refer to 2200193ADW Sheet 6, Rev 15)



Page 31




Figure 2-8 NP Lite Table (5194735ADW refer to 2200193ADW Sheet 8, Rev 15)

製造基準

組立製作図

5194735は2200193(NP Table)の一部の部品構成を変更したものである。

このADWでは2200193との変更点について示している。
このADWに指示がない風船番号の部品については、
2200193ADW参照のこと。

5194735 has similar construction with
2200193(NP table).

This ADW shows different point from 2200193.
For the mark number which is not shown in it,
refer to 2200193ADW.

△		P4 Erroneous recording correction (Rmemo-161-07-X078)	2042997	07-05-17	花島	中村	
△		P5 Scannable range label position change(Rmemo-4161-07-X045)		07-04-20	花島	中村	
△		P2.4 Update several items due to feed back from GEHW Mfg		07-02-21	花島	中村	
△		Drawing release for M3		06-12-19	花島	中村	
△		Initial Release		06-11-13	花島	中村	
REV.	PAGE	DESCRIPTION		DATE	REV.BY	CHK	
△		DR. 花島	機構設計	NP LITE TABLE WITHOUT DORRY AND ACC ASSY			
△		CH. 中村					
△		APP. 中村	2006-11-13				
△	全	06-12-19	花島	5194735 ADW			
				Page 1	of 5		



YOKOGAWA

GE Yokogawa Medical Systems

Approved Document - 5194735ADW_16.pdf Page 2 of 6

State: RELEASE - Document is released and under formal Change Control. Changes are subject to the ECRECO Process.
 See the GEHC Myworkshop System to determine the status of this document.

Figure 2-9 NP Lite Table (5194735ADW Sheet 1, Rev 5)



Section 1.0 — - NP LITE TABLE TOP ASSY POSITIONING BSD LITE (5194739)




Figure 2-11 NP Lite Table (5194735ADW Sheet 3, Rev 2)



Figure 2-12 NP Lite Table (5194735ADW Sheet 4, Rev 5)



Figure 2-13 NP Lite Table (5194735ADW Sheet 5, Rev 4)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	2113651	1.0	Yes	P9235CG ROLLER
19	2304402	1.0	Yes	P9237PJ CRADLE LOCK ASSY
63	2243960	1.0	Yes	P9235RA TA CYL ASSY
97	5194770	2.0	Yes	FOOT SW PEDAL ASSY Positioning BSD LITE
128	2141669	1.0	Yes	P9235EL BRKT ROLLER ASSY
184	5191040	1.0	Yes	CRADLE NP LITE ASSY
213	2210366	1.0	Yes	P9235NZ FRONT COVER BOT
214	2159531	1.0	Yes	P9235GQ RAER BOT COVER 2
218	5192576	1.0	Yes	TOP SIDE COVER R ASSY Positioning BSD Lite
221	5192577	1.0	Yes	TOP SIDE COVER L ASSY Positioning BSD Lite
222	2198339-2	1.0	Yes	P9235LS TOP COVER FRONT
223	2173530	1.0	Yes	P9235HN TOP COVER REAR
224	2198336-2	1.0	Yes	P9235LP IMS COVER REAR
225	2198332-2	1.0	Yes	P9235LM IMS COVER RIGHT
226	2198335-2	1.0	Yes	P9235LN IMS COVER LEFT
247	2210368	1.0	Yes	P9235PA R COVER STAY6
249	2181778	2.0	Yes	P9235JN R COVER STAY 7
251	2181780	1.0	Yes	P9235JQ COVER STAY 8
281	2192909	1.0	Yes	P9235KP COVER STAY 10
282	2192910	1.0	Yes	P9235KQ COVER STAY 11
306	2236089-2	1.0	Yes	P9235QL SPONGE L Positioning BSD Lite
307	2236088-2	1.0	Yes	P9235QK SPONGE R Positioning BSD Lite
319	2232363	4.0	Yes	P9235QC SHIM 0.5
320	2232364	12.0	Yes	P9235QD SHIM 1
336	2201349	1.0	Yes	P9235MA PUMP ASSY

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
337	2201350	1.0	Yes	P9235MP E VALVE ASSY
343	2201798	1.0	Yes	U0034AP BELT
345	2201814	1.0	Yes	U0036AP PULLEY
363	2201839	2.0	Yes	U0426AH BEARING ATTACHMENT BEARING UNIT
369	2201850	4.0	Yes	U0430AH CAM FOLLOWER WITH SHAFT 5MM 5 MM 13 MM
370	2201851	2.0	Yes	U0431AH CAM FOLLOWER WITH SHAFT 10MM 10 MM 22 MM
371	2201854	1.0	Yes	U0041AP BELT
401	N9020YX	6.0	Yes	SCREW HEXAGON SOCKET 10 MM 20 MM CAP SCREW 10X20 FE
404	N9040YX	4.0	Yes	SCREW HEXAGON SOCKET 10 MM 40 MM CAP SCREW M10X40 FE
406	N9101BS	2.0	Yes	NUT HEXAGON 10 MM NUT M10 2 FE ZN
412	N9304JS	14.0	Yes	SCREW
413	N9306ES	54.0	Yes	SCREW COUNTERSUNK HEA 3 MM 6 MM
414	N9306HR	10.0	Yes	SCREW TP 3X6 FE ZN
418	N9305JS	3.0	Yes	SCREW NABE 3X5 Fe Zn
432	N9401BS	4.0	Yes	NUT HEXAGON 4 MM NUT M4 FE ZN
433	N9401WS	40.0	Yes	WASHER PLAIN _ LARGE 4.3 MM 9 MM WASHER 4.3x9x0.8 FE ZN
434	N9405JS	6.0	Yes	SCREW NABE 4X5 Fe Zn
436	N9406SJ	6.0	Yes	SCREW SET SCREW 4 MM 6 MM
438	N9408HR	30.0	Yes	SCREW TP
439	N9408JS	13.0	Yes	SCREW NABE 4X8 Fe Zn
440	N9408LK	42.0	Yes	SCREW PAN HEAD 4 MM 8 MM
441	N9408ZX	27.0	Yes	SCREW HEXAGON SOCKET 4 MM 8 MM CAP SCREW 4X8 FE
442	N9410ES	2.0	Yes	SCREW SCREW SARA 4X10 Fe Zn * SCREW s[4X10 FE ZN
443	N9412ZX	16.0	Yes	SCREW HEXAGON SOCKET 4 MM 12 MM SCREW
446	N9501BS	2.0	Yes	NUT HEXAGON 5 MM NUT M5 FE ZN
447	N9501WS	76.0	Yes	WASHER PLAIN - LARGE 5.3 MM 10 MM WASHER M5

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
448	N9416ZX	2.0	Yes	SCREW HEXAGON SOCKET 4 MM 16 MM CAP 4X16
451	N9420JS	4.0	Yes	SCREW PAN HEAD 4 MM 20 MM SCREW osX 4X20 FE ZN
463	N9508SJ	8.0	Yes	SCREW SET SCREW 5 MM 8 MM SET SCREW M5X0.8X8htXvsg
464	N9510HR	3.0	Yes	SCREW TP
466	N9510ZX	56.0	Yes	SCREW HEXAGON SOCKET 5 MM 10 MM CAP SCREW 5X10 FE
468	U0030AA	5.0	Yes	SCREW Y_fhj_qDjhXr/X
469	N9516ZX	12.0	Yes	CAP SCREW 5X16 Fe
470	N9520ZX	27.0	Yes	SCREW HEXAGON SOCKET 5 MM 20 MM CAP SCREW 5X20 FE
471	N9525JS	4.0	Yes	SCREW SCREW NABE 5X25 Fe Zn * SCREW osX 5X25 FE ZN
473	N9601BS	4.0	Yes	NUT HEXAGON 6 MM NUT M6 2 FE ZN
474	N9601WS	6.0	Yes	WASHER PLAIN - LARGE 6.4 MM 12.5 MM WASHER 6.4X12.5X1.6FEZN
476	N9612YX	29.0	Yes	SCREW HEXAGON SOCKET 6 MM 12 MM CAP SCREW 6X12 FE
478	N9518TM	2.0	Yes	WASHER OUT OF CATALOG WASHER jj/sX 4.5X13X1 BS
480	N9510JS	7.0	Yes	SCREW NABE 5X10 Fe Zn
493	N9620YX	52.0	Yes	SCREW HEXAGON SOCKET 6 MM 20 MM CAP SCREW 6X20 FE
494	N9625YX	4.0	Yes	SCREW CAP 6X25
498	N9700ET	1.0	Yes	RETAINING RING E-RING 8 MM E-RING M7 SUS
501	N9801WS	16.0	Yes	WASHER PLAIN - LARGE 8.4 MM 17 MM
502	N9816YX	16.0	Yes	SCREW HEXAGON SOCKET 8 MM 16 MM CAP SCREW 8X16 FE
525	P9105SH	2.0	Yes	70W-1815 WASHER
526	P9105SR	1.0	Yes	SMR1-18 ROD END
527	P9123PC	6.0	Yes	CAP
528	U0129AE	2.0	Yes	CUSHION SJ-5008 _Xwa/
529	P9145CC	5.0	Yes	SHAFT ARM
530	2205823-2	1.0	Yes	TABLE CYLINDER SHAFT Positioning NP
531	P9145CE	1.0	Yes	CYLINDER BLOCK

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
532	P9145CF	2.0	Yes	RING
542	P9107SS	0.02	Yes	BLIND TAPE 2M
553	P9155BP	2.0	Yes	ADJUSTER BLOCK
554	P9155BQ	2.0	Yes	ADJ BOLT
555	P9155CH	1.0	Yes	SHAFT ARM
556	P9155CL	2.0	Yes	ROD POT
559	5175348	1.0	Yes	Elevation POT For NP lite
581	U0096AG	4.0	Yes	STANDOFF/SPACER MALE/FEMALE 4 MM 8 M
582	U0111AE	2.0	Yes	ADJUSTER
583	U0206AH	8.0	Yes	BEARING BALL, SINGLE RO 30 MM 55 MM 6006NR
584	U0212AC	42.0	Yes	CLAMP TMS_20
586	U0232AA	1.0	Yes	NUT HEXAGON 5 MM
587	U0320AK	1.0	Yes	VALVE
588	U0402AH	4.0	Yes	CAM FOLLOWER WITH SHAFT 6 MM 16 MM CF 6 VBUU
590	U0406AA	2.0	Yes	SCREW CAP SCREW
591	U0423AE	4.0	Yes	CATCH uhXr_n C-100A-2E
592	5224875	1.0	Yes	ELBOW U0493AK
597	U0462AK	2.0	Yes	JOINT CIN-3/8PT1/8
598	U0202AC	64.0	Yes	CLAMP CV_150
599	U0208AC	3.0	Yes	CLAMP TC-142 TAIRATSUPU KOTEI
600	U0203AC	3.0	Yes	CLAMP CV_200
601	U0509AC	0.807	Yes	PROTECTOR CE_032
602	U0507AC	0.228	Yes	PROTECTOR EDGING CE-016
603	U0506AC	28.4	Yes	PROTECTOR EDGING CE-012
604	U0534AC	0.06	Yes	PROTECTOR SPP-11 TUBE
612	U0611AK	1.0	Yes	U2-1-3/8

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
616	U0030AC	1.0	Yes	CLAMP OR CLIP CLAMP
617	U0031AC	4.0	Yes	CLAMP
701	5173372	1.0	Yes	NPCN board for NP Lite Table
706	5178099	1.0	Yes	MOTOR_NPLITE Positioning BSD Lite
712	46-297687P1	1.0	Yes	8 OHM, 4W, 3-1/2~ SPEAKER. 2.9 OZ CERAMIC MAGNET, 90 DB. (0.187~ Q.C. TERMINALS) (DWG REV 3/92, JFK - OK).
714	5176286	2.0	Yes	TAPE SWITCH NP LITE Positioning BSD Lite
817	46-297971P1	4.0	Yes	TERMINAL GROUND LUG FOR #4 TO 2/0 CABLE. ALUMINUM BASE, .25~ BOLT HOLE. SCREW COLORED GREEN FOR GROUND IDENT IFI- CATION. (DWG REV 12/92, JFK - OK) .
900	5175778	1.0	Yes	From Gantry to GTCB2 J9
901	5175734	1.0	Yes	From Home Latch ASSY to NPCN J12
902	5172410	1.0	Yes	NPSUB PWA for NP lite Table
903	5175748	1.0	Yes	From NPSUB board J3 to Valve and pump
904	5175846	1.0	Yes	From GTCB2-J4 to Sub-board J2
905	5175274	1.0	Yes	From GTCB2-J13 to NPCN J14
906	5175214	1.0	Yes	From GTCB2-J11 to NPCN board J2
907	5175602	1.0	Yes	From GTCB2-J12 to NPCN board J3
908	5175096	1.0	Yes	NPCN Board J10 to ABS ENCODER
909	5121382	1.0	Yes	FOOT SW R CABLE ASSY Positioning GT
910	5121381	1.0	Yes	FOOT SW L CABLE ASSY Positioning GT
911	5121380	4.0	Yes	FOOT SW END ASSY Positioning GT
912	5175968	1.0	Yes	From GTCB2 J10 to Speaker ASSY
913	5175990	1.0	Yes	From NPCN board J6 to GT Touch Sensor
914	5175934	1.0	Yes	From Step Motor driver CN2 to Motor
915	5173958	1.0	Yes	ENCODER for NP Lite table
917	5175897	1.0	Yes	From GANTRY 120VAC to Breaker For NP lite

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
918	5175539	1.0	Yes	From Filter to NPSUB board J1
919	5177613	1.0	Yes	From channel to channel support
920	5177634	1.0	Yes	From front leg to NP lite Base
921	5177667	1.0	Yes	From channel support to back leg
922	5177663	1.0	Yes	From back leg to NP lite base
923	5177373	1.0	Yes	From Pump to back leg
924	5127867	1.0	Yes	TABLE GND TO RACEWAY GND BAR CABLE Positioning GT
	2192909	1.0	Yes	P9235KP COVER STAY 10
	2326244	1.0	Yes	P9239FD TGAGE ASSY
1	2316894-2	1.0	No	CARRIAGE NP LITE Positioning BSD Lite
3	2113688-2	1.0	No	UPPER ARM R FOR NPLITE Positioning BSD Lite
4	2113691-2	1.0	No	UPPER ARM L FOR NPLITE Positioning BSD Lite
5	2113730	1.0	No	P9235BB RIGHT SUPPORT
6	2113731	1.0	No	P9235BC LEFT SUPPORT
8	2113733	1.0	No	P9235BY CRADLE RAIL
9	2113734	1.0	No	P9235BD LINK BAR F
10	2113774	1.0	No	P9235BE LINK BAR R
11	2113775	1.0	No	P9235CF BEARING NUT
12	2113776	1.0	No	P9235CE HANDLE NUT
15	2113779	1.0	No	P9235CJ SLIDE PLATE
16	2113780	1.0	No	P9235CK SLIDE SUP
17	2113781	1.0	No	P9235BZ RAIL PLATE
20	2336123	1.0	No	P9237QE STOPPER
21	2304379	1.0	No	P9237NP SHIM 1.0
22	2304378	2.0	No	P9237NN SHIM 0.5
25	2340403	4.0	No	P9237QJ SHIM 0.1

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
31	2316896	2.0	No	P9237MH BAR1
34	2316897	1.0	No	P9237MJ BAR3
35	2113916	1.0	No	P9235CP STOPPER
37	2113947	1.0	No	P9235CL DRIVE BASE
38	5183262	1.0	No	MOTOR SUP NPLITE Positioning BSD Lite
39	2113982	1.0	No	P9235BN REINFORCE1
40	2113983	1.0	No	P9235BP REINFORCE2
41	2113984-2	1.0	No	REINFORCE3 NP LITE TABLE Positioning BSD Lite
42	2113985	1.0	No	P9235BR REINFORCE4
43	2113986	1.0	No	P9235CD BEARING SUP
44	2257919	2.0	No	P9237KP RAIL COVER R W/TAPE ASSY
45	2114162	6.0	No	P9235DA SHIM1
46	2114163	6.0	No	P9235DB SHIM2
47	2114164	6.0	No	P9235DC SHIM4
48	2210365	2.0	No	P9235NY BAR CENTER
61	2114165	3.0	No	P9235DD SHIM8
65	2136523	1.0	No	P9225CS SHAFT
66	2141091-2	1.0	No	NP LITE TABLE BOT COVER SPRT Positioning BSD Lite
67	2141092-2	1.0	No	NP LITE TABLE COVER BRKT Positioning BSD Lite
68	2141093-2	1.0	No	NP LITE TABLE R BASE Positioning BSD Lite
69	2141094-2	1.0	No	LEFT BASE Positioning BSD Lite
70	2141095-2	1.0	No	RIGHT BASE Positioning BSD Lite
73	2141098-2	1.0	No	RAIL COVER F Positioning BSD Lite
76	2304368	1.0	No	P9237MN TRAY-F
77	2332998	1.0	No	P9237MD TRAY R W/O IMS
81	2304366	1.0	No	P9237ML HOLE PLATE

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
92	2141106	4.0	No	P9235LG COVER SUP R
93	2141109	1.0	No	P9235DH CONDENSER BRKT
94	2141110	1.0	No	P9235DJ POMP BRKT
95	2141181-2	1.0	No	NP LITE TABLE OUTER FRONT ARM Positioning BSD Lite
96	2141182	1.0	No	P9235DG INNER REAR ARM
126	2141667	1.0	No	P9235EF SUPPORT R
127	2141668	1.0	No	P9235EG SUPPORT L
136	2199989	3.0	No	P9235JA SHIM
158	2202925-2	1.0	No	SHAFT PULLY Positioning BSD Lite
159	2141869	1.0	No	P9235FL BEARING SUP
167	2142079-2	1.0	No	NP LITE TABLE F BASE Positioning BSD Lite
183	2143454	1.0	No	U0215AH BEARING BALL,SINGLE RO 15MM
185	2158471	1.0	No	P9235GE SUP COVER L
186	2158472	1.0	No	P9235GF SUP COVER R
187	2158473-2	1.0	No	NP LITE TABLE SPEAKER SPACER Positioning BSD Lite
189	2158475-2	1.0	No	NP LITE TABLE SPEAKER NET Positioning BSD Lite
190	2207566	1.0	No	P9235NS PUMP BASE
191	2158478	1.0	No	P9235CA HANDLE NUT L
192	2158479	2.0	No	P9235FM BEARING SPACER
212	2304377	1.0	No	P9237NM UNDER COVER
215	2304382	2.0	No	P9237NS COVER BRKT
216	2164759	1.0	No	P9235GU SIDE BASE L
217	2164760	1.0	No	P9235GV SIDE BASE R
241	2181770	1.0	No	P9235JE R COVER STAY 1
242	2181771	1.0	No	P9235JF L COVER STAY 1
243	2181772	2.0	No	P9235JG COVER STAY 2

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
245	2181774	2.0	No	P9235JJ COVER STAY 4
252	2181781	1.0	No	P9235JR COVER STAY 9
254	2181997	2.0	No	P9235KT IMS COVER STAY2
257	2182375	1.0	No	P9235JV COVER SUPPORT I1
271	2182376	1.0	No	P9235JW COVER SUPPORT I2
272	2182377	1.0	No	P9235JX COVER SUPPORT I3
273	2182378	1.0	No	P9235JY COVER SUPPORT I4
274	2207422	1.0	No	P9235NJ STOPPER_R
278	2185517	1.0	No	P9235HW FRAME SUPPORT2
279	2185518	2.0	No	P9235HX FRAME SUPPORT1
285	2194327-2	1.0	No	FRONT COVER BOT FOR NPLITE Positioning BSD Lite
286	2194328-3	1.0	No	P9235HZ SHEET
287	2194329	4.0	No	P9235BG COVER RAIL BRKT
288	2182815	2.0	No	P9235DE SHAFT ARM BOT
301	2194330-2	1.0	No	COVER STAY Positioning BSD Lite
302	2194333	1.0	No	P9235KM IMS COVER STAY10
303	2194334	1.0	No	P9235KN IMS COVER STAY11
304	2194335	2.0	No	P9235HS RAIL COVER CAP
305	2194336	1.0	No	P9235HV RAIL COVER F R
308	2194339	1.0	No	P9235JZ COVER SUP I5
311	2195672	1.0	No	P9235KC PLATE POT
312	2198680	1.0	No	P9235KG WIRE CLAMP
315	2199987	1.0	No	P9235LK COVER STAY 13
316	2200025	1.0	No	P9235MC F BOT COVER STAY L
317	2200026	1.0	No	P9235MD F BOT COVER STAY R
318	2194650	2.0	No	P9235GZ COVER RAIL BRKT F

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
331	2200027	2.0	No	P9235ME CARBON BRKT1
332	2200028	2.0	No	P9235MF CARBON BRKT2
333	2200530	2.0	No	P9235HA SIDE PL
334	2200531	1.0	No	P9235HD DRIVER COVER
335	2200532	1.0	No	P9235LY FRAME SUPPORT3
339	2201656	1.0	No	U0479AK HOSE ADAPTER
340	5229873	1.0	No	TEE U0495AK
342	2201668	1.0	No	U0623AK HOSE
348	2207420	1.0	No	P9235NG REAR_PUMP_CO
362	2201819	1.0	No	U0029AH BEARING BALL,SINGLE RO 9MM
372	2201855	1.0	No	P9235MW IMS BOT COVER F
376	2202293	1.0	No	P9235MY LEG COVER STAY1
377	2202302	1.0	No	P9235MZ LEG COVER
379	2203765	14.0	No	U0452AA BOLT
402	N9025YX	4.0	No	SCREW CAP SCREW 10X25FE * CAP SCREW 10X25FE
405	N9100WU	4.0	No	WASHER FLAT WASHER10.5X21X2 SUS q[?/_ - 10.5X21X2 SUS
407	N9120WS	10.0	No	WASHER FLAT WASHER 13X24X2.5 Fe q[?/_ - 13X24X2.5 FE ZN
419	N9301WS	16.0	No	WASHER PLAIN - LARGE 3.2 MM 7 MM
431	N9314JS	16.0	No	SCREW PAN HEAD 3 MM 14 MM SCREW osX 3X14 FE ZN
435	N9406ES	6.0	No	SCREW SARA 4X6 Fe Zn
444	N9410ZX	2.0	No	SCREW CAP SCREW 4X10 Fe CAP SCREW 4X10 FE
461	N9505JS	1.0	No	SCREW SCREW NABE 5X5 Fe Zn * SCREW osX 5X5 FE ZN
462	N9506ZX	4.0	No	SCREW CAP 5X6 * CAP 5X6
477	N9614YX	8.0	No	SCREW CAP SCREW 6X14 Fe * CAP SCREW 6X14 FE
482	N9514ZX	5.0	No	SCREW HEXAGON SOCKET 5 MM 14 MM
483	N9525EM	3.0	No	SCREW SCREW SARA 5X25 SUS316 * SCREW s[5X25 SUS316

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)


MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
491	N9616YX	16.0	No	SCREW HEXAGON SOCKET 6 MM 16 MM CAP SCREW 6X16 FE
492	N9620VX	12.0	No	SCREW SARA CAP SCREW 6X20 Fe s[CAP SCREW 6X20 FE
495	N9625ZX	10.0	No	SCREW CAP SCREW 12X25 Fe * CAP SCREW 12X25 FE
497	N9640YX	6.0	No	SCREW HEXAGON SOCKET 6 MM 40 MM CAP SCREW 6X40 FE
499	N9800SP	2.0	No	WASHER SPRING WASHER 8.2X15.4X2 .~DrjhX?/_ - 8.2X15.4X2
500	N9801BS	2.0	No	NUT HEXAGON 8 MM FE ZN
504	N9830YX	4.0	No	SCREW HEXAGON SOCKET 8 MM 30 MM CAP SCREW 8X30 FE
505	N9820YX	12.0	No	SCREW HEXAGON SOCKET 8 MM 20 MM CAP SCREW 8X20 FE
506	2299437	1.0	No	U0472AA_FLANGE_SOCKET
507	002-6050	1.0	No	BIG WASHER FOR M5 GB96
521	P9103FH	4.0	No	SPACER
522	P9103FJ	4.0	No	SPACER
523	P9103FK	7.0	No	SPACER
534	P9145DC	1.0	No	KEY
538	2205824	1.0	No	P9235ND CYLINDER BLOCK
557	P9185BV	1.0	No	BASE_ANGLE BASE ANGLE
567	U0022BN	8.0	No	NAME PLATE NAME PLATE-EARTH MARK
568	2225033	1.0	No	P9235QB BASE ANGLE
585	5192249	4.0	No	NUT C440 U0279AA
589	5224936	1.0	No	ELBOW U0494AK
596	U0461BA	2.0	No	TUBE_D16/16 INCH
605	U0241AA	2.0	No	NUT HEXAGON 4 MM NUT
611	U0472AK	1.0	No	ELBOW
800	5194600	1.0	No	POWER DRIVER ASSY Positioning BSD Lite
801	5194601	1.0	No	NP LITE BOTTOM COVER ASSY Positioning BSD Lite
810	5177186	1.0	No	ENCODER MOUNTING Positioning BSD Lite

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
811	5183495	1.0	No	PULLEY1 NPLITE Positioning BSD Lite
812	5183264	1.0	No	PULLEY2 NPLITE Positioning BSD Lite
813	2113733-2	1.0	No	CRADLE RAIL Positioning BSD Lite
816	5112992	1.0	No	GROUND Positioning GT
818	5183591	2.0	No	GROUND SUP Positioning BSD Lite
819	5193303	1.0	No	ENCODER COVER Positioning BSD Lite
825	U0098AG	4.0	No	STANDOFF/SPACER STANDOFF STUD FF-13
826	U0036BN	3.0	No	NAME PLATE
827	2379698	0.008	No	U0064BA_EDGING_TAPE

Table 2-3 NP Lite Table (5194735 BOM, Rev 8)

FRU = YES



MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
14	N9510ZX	1.0	Yes	SCREW HEXAGON SOCKET 5MM 10MM CAP SCREW 5X10 FE
16	N9306ZX	2.0	Yes	SCREW HEXAGON SOCKET 3MM 6MM CAP SCREW 3X6 FE
17	N9408ZX	4.0	Yes	SCREW HEXAGON SOCKET 4MM 8MM CAP SCREW 4X8 FE
19	N9516ZX	4.0	Yes	CAP SCREW 5X16 FE
20	U0509AC	0.104	Yes	PROTECTOR CE_032
1	2334963	1.0	No	P9237QC SOLENOID ASSY
2	2304367	1.0	No	P9237MM SOLENOID PLATE
3	2304373	1.0	No	P9237MT ARM
4	2304374	1.0	No	P9237MU STOPPER BLOCK
5	2316895	1.0	No	P9237MG STOPPER LEG
6	2304375	1.0	No	P9237MV BASE
7	2304376	2.0	No	P9237NL GUIDE BLOCK
8	2304381	2.0	No	P9237NR NUT PLATE
9	2304404	1.0	No	P9237PY SOLENOID SPONGE
10	2304405	1.0	No	P9237PZ STOPPER RUBBER
11	2320464	1.0	No	P9237PF SOLENOID RUBBER
12	2327140	1.0	No	P9237PC CANCEL LEVER-A
13	2336051	1.0	No	U0037AQ SPRING
15	N9310ZX	5.0	No	SCREW HEXAGON SOCKET 3MM 10MM CAP SCREW 3X10 FE
18	N9316PS	1.0	No	PIN

Table 2-4 P9237PJ CRADLE LOCK ASSY (2304402 BOM, Rev 2)

1.1.2 TA CYL ASSY (2243960)

FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
	U0462AK	1.0	Yes	JOINT CIN-3/8PT1/8
	U0611AK	1.0	Yes	U2-1-3/8
	2245857	1.0	No	U0709AK CYLINDER
	2245852	1.0	No	U0486AK ELBOW

Table 2-5 TA CYL ASSY (2243960 BOM, Rev 5)

FRU = YES

Figure 2-15 Foot SW Pedal Assy Positioning BSD Lite (5194770ADW refer to 2141647ADW, Rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	2174553	2.0	Yes	P9235JU PEDAL
4	5120586-2	2.0	Yes	PEDAL COVER POSITIONING GT
6	P9155EQ	8.0	Yes	STOP RUBBER
13	5120587-2	2.0	Yes	PEDAL COVER LID POSITIONING GT
21	N9401WS	16.0	Yes	WASHER PLAIN _ LARGE 4.3 MM 9 MM WASHER 4.3X9X0.8 FE ZN
22	N9408JS	14.0	Yes	SCREW NABE 4X8 FE ZN
23	U0023AA	8.0	Yes	SCREW
24	N9408ZX	8.0	Yes	SCREW HEXAGON SOCKET 4 MM 8 MM CAP SCREW 4X8 FE
25	N9501WS	4.0	Yes	WASHER PLAIN - LARGE 5.3 MM 10 MM WASHER M5
26	N9510JS	4.0	Yes	SCREW NABE 5X10 FE ZN
27	N9520JS	4.0	Yes	SCREW PAN HEAD 5 MM 20 MM SCREW OSX 5X20 FE ZN
28	N9608YX	2.0	Yes	SCREW CAP SCREW 6X8 FE
30	U0233AA	2.0	Yes	NUT HEXAGON 6 MM M3-17-9
1	2174552	2.0	No	P9235JT SW PEDAL
3	2174555	1.0	No	P9235JS COVER SUPPORT
5	2201351	1.0	No	P9235EW FOOT SW BASE
7	P9185ES	4.0	No	SW_BRKT2
8	P9225EM	2.0	No	PLATE_WASHER
9	P9225EN	4.0	No	SHAFT
10	P9225EP	4.0	No	STOPPER
11	P9225ER	2.0	No	SW_BRKT1
12	P9225EW	2.0	No	SW_BRKT3
29	U0004AQ	4.0	No	SPRING COMPRESSION 47 MM SC-085-120-0470 4017
31	U0445AA	2.0	No	CAP SCREW

Table 2-6 Foot SW Pedal Assy Positioning BSD Lite (5194770 BOM, Rev 1)

FRU = YES



Figure 2-16 BRKT Roller Assy (2141669ADW, Rev1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	2199992	1.0	Yes	P9235LB BRKT ROLLER
4	2201853	2.0	Yes	U0226AH BEARING BALL,DOUBLE RO 20MM
10	P9145CZ	4.0	Yes	WASHER ROLLER
12	N9612YX	4.0	Yes	SCREW HEXAGON SOCKET 6 MM 12 MM CAP SCREW 6X12 FE
1	2199990	2.0	No	P9235JB PIN
2	2199991	4.0	No	P9235LA SPACER INNER

Table 2-7 BRKT Roller Assy (2141669 BOM, Rev 4)

1.1.5 TOP SIDE COVER R ASSY POSITIONING BSD LITE (5192576)

FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2198337-2	1.0	Yes	P9235LQ TOP COVER RIGHT
2	5176286	1.0	Yes	TAPE SWITCH NP LITE POSITIONING BSD LITE
3	2199652	1.0	Yes	P9235LZ LABEL
4	2199976	1.0	Yes	P9235LW SR MEAN LABEL
6	U0008BD	5.0	Yes	ADHESIVE B.-D-30(ZJH)
5	5181948	1.0	No	SWITCH PLUG POSITIONING BSD LITE

Table 2-8 TOP SIDE COVER R ASSY POSITIONING BSD LITE (5192576 BOM, Rev 3)

1.1.6 TOP SIDE COVER L ASSY POSITIONING BSD LITE (5192577)

FRU = YES

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2198338-2	1.0	Yes	P9235LR TOP COVER LEFT
2	5176286	1.0	Yes	TAPE SWITCH NP LITE POSITIONING BSD LITE
3	2199652	1.0	Yes	P9235LZ LABEL
4	2199976	1.0	Yes	P9235LW SR MEAN LABEL
6	U0008BD	5.0	Yes	ADHESIVE B.-D-30(ZJH)
5	5181948	1.0	No	SWITCH PLUG POSITIONING BSD LITE

Table 2-9 TOP SIDE COVER L ASSY POSITIONING BSD LITE (5192577 BOM, Rev 3)

1.1.7 RAIL COVER R W/TAPE ASSY (2257919)

FRU = NO




Figure 2-17 Rail Cover R W/Tape Assy (2257919ADW, Rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2258937	2.0	Yes	P9237KW RAIL COVER R
2	2257923	4.0	No	P9237SE MAGIC TAPE FEMALE

Table 2-10 Rail Cover R W/Tape Assy (2257919 BOM, Rev 0)

1.1.8 POWER DRIVER ASSY POSITIONING BSD LITE (5194600)

FRU = NO




Figure 2-18 Power Driver Assy Positioning BSD Lite (5194600ADW, Rev 3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	5182383-2	1.0	Yes	CRADLE MOTOR DRIVER ROHS POSITIONING BSD LITE
4	46-297335P3	1.0	Yes	16A, 250V 2-STAGE PWR LINE FILTER
5	5151535	1.0	Yes	CIRCUIT BREAKER 15 AMPS 250V U0061EK
7	5159097	1.0	Yes	24V SWITCHING POWER SUPPLY U0298EE
8	5142177-3	1.0	Yes	GTCB2 BOARD ASSY POSITIONING GT
20	U0212AC	4.0	Yes	CLAMP TMS_20
21	U0096AG	3.0	Yes	STANDOFF/SPACER MALE/FEMALE 4 MM 8 M
22	N9408JS	16.0	Yes	SCREW NABE 4X8 FE ZN
23	N9306JS	4.0	Yes	SCREW
24	N9408HR	4.0	Yes	SCREW TP
25	N9510HR	9.0	Yes	SCREW TP
26	N9405JS	3.0	Yes	SCREW NABE 4X5 FE ZN
27	N9306ES	4.0	Yes	SCREW COUNTERSUNK HEA 3 MM 6 MM
28	N9508JS	1.0	Yes	SCREW NABE 5X8 FE ZN
30	5175131	1.0	Yes	FROM POWER SUPPLIER TO STEP MOTOR DRIVER
31	5175527	1.0	Yes	FROM PS TO GTCB2-J31
32	5175663	1.0	Yes	FROM AC FILTER TO PS
33	5174484	1.0	Yes	NP LITE TABLE CABLE, FROM GTCB2 J5 TO STEP MOTOR DRIVER CN1.
34	5183805	1.0	Yes	FROM GTCB2 J14 TO HART J1 POSITIONING LITETABLE
1	5178523	1.0	No	POWER DRIVER SUPPORT POSITIONING BSD LITE
2	5178248	1.0	No	POWER DRIVER BASE PLATE POSITIONING BSD LITE
6	5176418	1.0	No	PROTECTOR BOX POSITIONING BSD LITE
9	5183171	1.0	No	HART BOARD PROTECTOR POSITIONING BSD LITE
10	5179856	1.0	No	HART BOARD ASSY
11	5192250	1.0	No	NAME PLATE U0074BN
12	U0022BN	1.0	No	NAME PLATE-EARTH MARK

Table 2-11 Power Driver Assy Positioning BSD Lite (5194600 BOM, Rev 7)

1.1.9 NP LITE BOTTOM COVER ASSY POSITIONING BSD LITE (5194601)

FRU = NO




Figure 2-19 NP Lite Bottom Cover Assy Positioning BSD Lite (5194601ADW, Rev 3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	2142068	1.0	Yes	P9235GA FLAPPER
4	2173525-2	1.0	Yes	P9235HH BOTTOM COVER R
5	5193093-2	1.0	Yes	FOOT SW COVER INJ NPLITE POSITIONING BSD Lite
6	U0908AH	2.0	Yes	BUSHING 80F-0608
7	2401506	6.0	Yes	P0793HN HANDLE CAP
8	5180382	1.0	Yes	BOTTOM SKIRT COVER POSITIONING BSD LITE
9	5180501	1.0	Yes	BTM FRONT BLIND COVER POSITIONING BSD LITE
20	N9510ZX	4.0	Yes	SCREW HEXAGON SOCKET 5 MM 10 MM CAP SCREW 5X10 FE
21	N9510LK	6.0	Yes	SCREW PAN HEAD 5 MM 10 MM
22	N9408JS	4.0	Yes	SCREW NABE 4X8 FE ZN
1	2142069-2	1.0	No	P9235GB FLAP BRKT R POSITIONING BSD LITE
2	2142070-2	1.0	No	P9235GC FLAP BRKT LH POSITIONING BSD LITE
10	5180384	2.0	No	BTM COVER SET PIN POSITIONING BSD LITE

Table 2-12 NP Lite Bottom Cover Assy Positioning BSD Lite (5194601 BOM, Rev 4)

1.2 NP LITE ACCESSORY COLLECTOR POSITIONING BSD LITE (5194741)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2115996-3	1.0	Yes	AXIAL HEADHOLDER ASM WITH RATING PLATE AND SELF-CATCHING LATCH
2	2201803-2	1.0	Yes	P9235RJ EXTENDER ASSY
	2115985-3	1.0	Yes	P9230JU CRADLE PAD
3	2115985-5	1.0	Yes	LCC CRADLE PAD POSITIONING BSD LITE

Table 2-13 NP Lite Accessory Collector Positioning BSD Lite (5194741 BOM, Rev 2)

1.3 NP LITE SERVICE TOOL COLLECTOR POSITIONING BSD LITE (5192574)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	5192575	1.0	Yes	FEK SCREW KIT POSITIONING BSD LITE
3	5177354	1.0	Yes	CHAR CALIBRATIONI JIG1 POSITIONING BSD LITE
4	5177619	1.0	Yes	CHAR CALIBRATION JIG2 POSITIONING BSD LITE
5	2199974	2.0	Yes	P9235LU C ARM MARK
6	2200775	2.0	Yes	P9235MG C ARM MARK2
7	P9123PC	2.0	Yes	CAP
8	P9145LN	2.0	Yes	HANDLE CAP
9	2401506	2.0	Yes	P0793HN HANDLE CAP
10	U0202AC	20.0	Yes	CLAMP CV_150
11	U0203AC	20.0	Yes	CLAMP CV_200
12	5191747	4.0	Yes	DUMMY CONNECTOR FOR TOUCH SWITCH
13	5192391	1.0	Yes	CABLE, OUTLET TO TABLE PUMP FOR LITE TABLE SERVICE TOOL
1	2212862	1.0	No	U0043BF TOOL BOX
14	2199989	10.0	No	P9235JA SHIM
15	5220232	1.0	No	MSDS COLLECTOR POSITIONING BSD LITE

Table 2-14 NP LITE SERVICE TOOL COLLECTOR POSITIONING BSD LITE (5192574 BOM, Rev 4)

1.3.1 FEK SCREW KIT POSITIONING BSD LITE (5192575)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	N9306ES	5.0	Yes	SCREW COUNTERSUNK HEA 3 MM 6 MM
2	N9408HR	5.0	Yes	SCREW TP
3	N9408JS	5.0	Yes	SCREW NABE 4X8 FE ZN
4	N9408LK	5.0	Yes	SCREW PAN HEAD 4 MM 8 MM
5	N9412ZX	5.0	Yes	SCREW HEXAGON SOCKET 4 MM 12 MM SCREW
6	N9416ZX	5.0	Yes	SCREW HEXAGON SOCKET 4 MM 16 MM CAP 4X16
7	N9510ZX	5.0	Yes	SCREW HEXAGON SOCKET 5 MM 10 MM CAP SCREW 5X10 FE
8	N9612YX	5.0	Yes	SCREW HEXAGON SOCKET 6 MM 12 MM CAP SCREW 6X12 FE
9	N9620YX	5.0	Yes	SCREW HEXAGON SOCKET 6 MM 20 MM CAP SCREW 6X20 FE
12	N9401WS	5.0	Yes	WASHER PLAIN _ LARGE 4.3 MM 9 MM WASHER 4.3X9X0.8 FE ZN
13	N9501WS	5.0	Yes	WASHER PLAIN - LARGE 5.3 MM 10 MM WASHER M5
14	N9518TM	5.0	Yes	WASHER OUT OF CATALOG WASHER J]/SX 4.5X13X1 BS
15	N9601WS	5.0	Yes	WASHER PLAIN - LARGE 6.4 MM 12.5 MM WASHER 6.4X12.5X1.6FEZN
16	U0023AA	5.0	Yes	SCREW
17	U0030AA	5.0	Yes	SCREW Y_FHJ_QD]HXR/X
18	U0232AA	5.0	Yes	NUT HEXAGON 5 MM
19	U0233AA	5.0	Yes	NUT HEXAGON 6 MM M3-17-9
10	N9820YX	5.0	No	SCREW HEXAGON SOCKET 8 MM 20 MM CAP SCREW 8X20 FE
11	N9830YX	5.0	No	SCREW HEXAGON SOCKET 8 MM 30 MM CAP SCREW 8X30 FE

Table 2-15 FEK SCREW KIT POSITIONING BSD LITE (5192575 BOM, Rev 1)