

MI&CT ESERVICE NOTE
SN07030710a

Class C information is Property of GE.

This short-term informational note does not apply after April 1st, 2007. If setting time after April 1st, use “setdate” command to change time!

Date: 03/21/2007

Subject: Configuring Daylight Saving Time for IRIX O/S (Spring 2007).

Effectivity: This Service Note only applies to sites in the United States and Canada. And only those sites that observe Daylight Saving Time.

The following systems with the software revisions listed below are affected: LightSpeed CT Products (Plus, Ultra, LS16). Effectivity for this applies to all LightSpeed 2X – 4X systems with IRIX based GOC consoles and CT systems using **Sgi / IRIX – O/S**.

- 405I.3_H4.0M5.1
- 308I.2_H3.1M5
- 101I.2_HSOXIM5

Note: This is not applicable to CT/PET Systems, and can be done remotely.

Details: Daylight Saving Time (DST) in the US and Canada arrives three weeks earlier and finishes one week later, beginning in 2007. This Service Note disables DST so that the time does not change on the previous DST dates, (e.g. April 1st). **This procedure must be run anytime after running 'reconfig'.**

Tool Requirements: No special tools required.

Manpower Requirements: One (1) Service Engineer for 20 minutes.

Resolution: Note: The initial verification steps can be done with Applications up, but prior to running “setdate”, Applications will need to be shut down!

1. Check your software version.

Type: **check config**

Press **[Enter]**

Note: If the result does not contain one of the software versions listed in the Effectivity section, this service note does not apply.

2. Check your existing Time Zone parameters.

3. Click *Unix Shell* under the *Tool Chest* in the upper right of the monitor and put the cursor in the Shell window. If done remotely, just enter the following in your telnet window.

Type: **grep TZ /etc/TIMEZONE**

Press **[Enter]**

If the result does not contain one of the following outputs, **STOP**, this procedure does not apply. There are several reasons why you may not have an output matching the following table. For example, if the area in which your system is located does not recognize Daylight Saving Time,

such as parts of Arizona, or if your system has already been changed to Universal Time (e.g. GMT-6). In either case, do not continue with this procedure.

TZ=AST4ADT3	TZ=EST5EDT	TZ=CST6CDT	TZ=MST7MDT
TZ=PST8PDT	TZ=AKST9AKDT	TZ=HAST10HADT	

4. Please note the results from step 3, as they will be needed in Step 7 below.

5. Check the Time and Zone info with the “**date**” command:
Type: **date** Press [Enter]

(Output will be similar to below)

Sun Mar 11 01:55:02 CST 2007

Note: Standard Time is indicated by: CST or AST, EST, MST, PST, etc... If still set to Standard Time and it's before April 1st 2007 the following steps will turn off old DST rules, so 'setdate' is not needed till Oct. 28th 2007, or executing this service note again. (Eliminates customer RFS call, April 1st or 2nd.)

6. The following steps can be done with Applications up, but a reboot must occur before “**setdate**” changes take effect!
 - a. If Applications are up select the *Service Desktop* Icon, *Utilities*, and then *Application Shutdown*.
 - b. Click *Unix Shell* under the *Tool Chest* in the upper right of the monitor and put the cursor in the Shell window. Enter the following in the telnet window.
 - c. Type: **su -** Press [Enter]
 - d. Type the appropriate password. Press [Enter]
 - e. Type: **configTime -p** Press [Enter]
Starting timezone configuration
Please enter correct timezone >
7. Choose the correct Time Zone line for your US/Canada DST 2007 time zones, and turn off all DST rules to only use Standard Time.

<u>Time Zone</u>	<u>Fix / Solution</u>
For AST4ADT3 (Atlantic Time Zone)	Type: AST4 Press [Enter]
For EST5EDT (Eastern Time Zone)	Type: EST5 Press [Enter]
For CST6CDT (Central Time Zone)	Type: CST6 Press [Enter]
For MST7MDT (Mountain Time Zone)	Type: MST7 Press [Enter]
For PST8PDT (Pacific Time Zone)	Type: PST8 Press [Enter]
For AKST9AKDT (Alaska Time Zone)	Type: AKST9 Press [Enter]
For HAST10HADT (Hawaii Time Zone)	Type: HAST10 Press [Enter]

8. Check results:
Type: **grep TZ /etc/TIMEZONE** Press [Enter]

(One of the following should appear in the output).

TZ=AST4	TZ=AKST9	TZ=EST5
TZ=CST6	TZ=MST7	TZ=PST8
TZ=HAST10		

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved.
Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare
Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE
Healthcare is prohibited.
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

9. Check Time and Zone info with the “**date**” command:

Type: **date** Press **[Enter]**

Sun Mar 11 01:55:02 CST 2007

Note: Standard Time is indicated by: CST or AST, EST, MST, PST, etc...

- **If the date is correct at this point, you are done.** Another setdate command is not needed until the fall of 2007 when DST changes again, or the time drifts and requires resetting.

The changes made to disable DST will become effective at the next reboot!

- **If the time and date are not correct in step 9, then continue with ‘reboot’ in step 10, and correct the time with “setdate”!**

10. Reboot system with the following command:

Type: **rebootsys** Press **[Enter]**

Expected results from “rebootsys”:

```
/usr/g/scripts/rebootsys[7]: source: not found
Shutdown started.      Fri Mar  2 01:37:49 CST 2007
UX:wall: INFO: Broadcast Message from root (???) on msecqxi3 Fri Mar  2
01:37:49 2007
THE SYSTEM IS BEING SHUT DOWN! Log off now.
Shutdown started.      Fri Mar  2 01:37:49 CST 2007
Broadcast Message from root (ttyq1) on msecqxi3 Fri Mar  2 01:37:50 2007
The system will be shut down in 10 seconds.
Please log off now.
Broadcast Message from root (ttyq1) on msecqxi3 Fri Mar  2 01:38:00 2007
THE SYSTEM IS BEING SHUT DOWN! Log off now.
```

11. Check displayed time after the reboot with the “**date**” command:

Type: **date** Press **[Enter]**

12. You will have to perform the ‘**setdate**’ command, if the time displayed is incorrect.

- If Applications are up and site allows you to shutdown, select the *Service Desktop Icon*, *Utilities*, and then *Application Shutdown*.
- Click *Unix Shell* under the *Tool Chest* in the upper right of the monitor and put the cursor in the Shell window. If done remotely, enter the following in the telnet window.
- Type: **su -** Press **[Enter]**
- Type the appropriate password. Press **[Enter]**
- Type: **setdate** Press **[Enter]**

Press **[Enter]** if the displayed field is correct, or type in the appropriate date/time parameters.
Also, if time change is entered, another reboot should be performed.

13. Reboot the system with the following command:

Type: **rebootsys** Press **[Enter]**

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved.
Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare
Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE
Healthcare is prohibited.
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

14. To re-check, and verify all Time Zone settings:

Type: **grep TZ /usr/g/config/INFO /etc/TIMEZONE; date** Press [Enter]

Note: When running 'reconfig' changes are made to the file named: /etc/TIMEZONE, therefore you must execute this Service Note, every time 'reconfig' is performed. If service GUI – "reconfig" is run, delay the system reboot by clicking **NO** on the exit box, and return to beginning of this service note.

Additionally, whenever 'reconfig' is performed, if site is Insite enabled, the 'ConfigLink' command must also be done to add the OLC hosts. To run the 'ConfigLink' command do the following:

Type: **cd ~insite; ConfigLink;** Press [Enter]

Kirby Louque
OnLine Engineer

MI&CT Engineering Manager Signature: Mark Nyman