



GE Healthcare

Technical Publications

Direction 2409677

Revision 5

**Signa® 3.0T with EXCITE Technology
Field Replaceable Units (FRU)**

Copyright© 2004 – 2008 by General Electric Company, Inc.

All Rights Reserved



GE Healthcare

GE Medical Systems; Telex 3797371
P.O. Box 414, Milwaukee, Wisconsin 53201 U.S.A.
(Asia, Pacific, Latin America, North America)

GE Medical Systems – Europe; Telex 261794
Shortlands, Hammersmith, London W6 8BX U.K.



**COVER
BOTTOM
CROP LINE**

**COVER
BOTTOM
CROP LINE**

**5/16" RADIUS
CORNER TRIM
TYPICAL**



**COVER
LEFT
CROP LINE**

**COVER
FOLD
LINE**

**(2) CENTER STAPLES
FOR BINDING**

**COVER
RIGHT
CROP LINE**



**5/16" RADIUS
CORNER TRIM
TYPICAL**

REVISION HISTORY

<u>REV/DATE</u>	<u>REASON FOR CHANGE</u>
-----------------	--------------------------

- | | | |
|---|------------|--|
| 0 | May 25, 04 | Initial release |
| 1 | Jan 06, 05 | Added cable to LPCA. Added new monitor and GOC AA to OW. Added HFD-S cabinet. Added latest RF Module parts. Corrected part number in RF Accessory cabinet and added new ASC Chassis parts. Added 2nd 25kW NB RF amp for sites without LCD display. Added 25/35kW NB RF amp. Added BrainWave options and Fiesta C option. |
| 2 | May 23, 05 | Added information for systems upgraded to EXCITE HD: Enclosure Modules, LPCA, Penetration Panel, OW Modules, HFD/PDU, RFS Cabinet, and cables. |
| 3 | Dec 30, 05 | Updated Body Coil Blower and SRI parts. Added 6 & 8 Channel Resistive Shim filters to Pen Panel. Updated OW monitor info. Updated RRF module and UPM Rf Detector Board part numbers. Added NB RF Amp Cabinet temp sensor FRU. |
| 4 | Jan 26, 07 | Added Speaker Grill nylon screw, LPCA Wheel Housing FRU kit (MRIhc21559), & Dock pedal FRU parts (MRIhc17629). Updated Host monitor FRU. Added RRF Module MUX2 board and Term Server cable FRU kit. Added 10cm Silicon Phantom Service Tool. |
| 5 | May 28, 08 | Rear Pedestal: Updated Longitudinal Encoder part # (PQR 13175398) and added Encoder board part # 46-226552G2. Chiller: Added Lytron pump coupler part # (PQR 13136000). MGD: Updated SCP2 Bd part #. Added 5181495-3 MGD P.S. NB RF Amp: Removed 25/35kw MKS RF Amp (PQR 13142203). |

SEND COMMENTS ABOUT THIS MANUAL TO:

GE HEALTHCARE
MANAGER, SERVICE ENGINEERING
W-808
3200 N. GRANDVIEW BLVD.
WAUKESHA, WISCONSIN 53188
USA

This manual belongs to:

INTRODUCTION TO FIELD REPLACEABLE UNIT (FRU)

This direction is a guide to the FRU's for 3.0T with EXCITE Technology upgrade. A complete parts list of the 3T Eclipse system can be found in Direction 2236792.

Parts designated as a FRU are critical to maintain the safety and/or functionality of the system. FRU's are available on a Priority 1 basis from Service Parts Operation.

As a reference aid, the FRU's are listed under an upper-level part number where they are used. In most cases the upper-level part number is not a FRU. In some cases the parts listed under an upper-level assembly may have been nested in intermediate sub-assemblies which are not shown.

Where possible, FRU's which are part of Option Catalogs are shown under the area where they would be installed, rather than being listed in the Coils or Options Sections (SEE TABLE OF CONTENTS).


LIST OF EFFECTIVE PAGES

<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>
Cover Page	5	10 to 11	2	26	4
1	5	12 to 13	3	27 to 28	2
2	2	14	5	29	3
3 to 4	5	15	2	30	2
5	2	16	3	31	5
6	4	17 to 22	2	32 to 34	2
7	5	23 to 24	4	35	4
8 to 9	4	25	5	36	2

This image shows a single page of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page, leaving small margins at the top and bottom. There is no handwriting or other markings on the page.

TITLE	PAGE
INTRODUCTION TO FIELD REPLACEABLE UNIT(FRU)	1
LIST OF EFFECTIVE PAGES	1
TABLE OF CONTENTS	3
PICTORIAL TABLE OF CONTENTS	4
MAGNET AND ENCLOSURE (ASSEMBLIES)	5, 6, 7, 8, & 9
MISCELLANEOUS LIST (OF PARTS NOT SHOWN)	5
BRM & TRM LIST (OF PARTS NOT SHOWN)	5
FRONT COSMETIC COVER & PARTS (BACK VIEW)	5
BLOWER CABINET	6
SRI MODULE	6
COMFORT MODULE	6
REAR PEDESTAL PARTS	7
MAGNET DOCKING PARTS	7
LOW PROFILE CARRIAGE ASSEMBLY	8 & 9
PATIENT SUPPORT TABLE	10
PENETRATION PANEL	12 & 13
CHILLER	14
OPERATOR WORKSPACE – EXCITE	16
OPERATOR WORKSPACE COMPUTER CABINET – LINUX	17
WORKSPACE SCIM AND KEYBOARDS	18
MAIN DISCONNECT PANEL	19
ACGD/PDU CABINET	20
HFD–S/PDU CABINET	21
HFD/PDU CABINET	22
RFS CABINET	23
SYSTEM CABINET	24
RFS & SYSTEM CABINET MODULES	25, 26, & 27
RF ACCESSORY CABINET	28 & 29
NB RF AMP CABINET	30
8KW MNS BB RF AMP CABINET (OPTIONAL)	32
ACCESSORIES & CABLES	34
COILS	35
EXCITE OPTIONS	36

INTRO TO FRU'S	PAGE 1
CHILLER PARTS	PAGE 14
ACCESSORIES	PAGE 34
CABLES	PAGE 34
COILS	PAGE 35
OPTION KEYS	PAGE 36


MAGNET ENCLOSURE:**SUBASSEMBLIES PGS 5, 6, 7, 8 & 9**

PATIENT AIR SUPPLY
BLOWER CABINET (MG6)
PAGE 6



PENETRATION
PANEL
PAGES 12 & 13

PATIENT
SUPPORT
TABLE
PAGE 10

**MAGNET ROOM****OPERATOR/EQUIPMENT ROOM**


2219400
MAGNET
MONITOR
(MSM)




MAIN DISCONNECT
PANEL
PAGE 19




EXCITE II OPERATOR
WORKSPACE WITH LINUX
PAGES 16, 17, & 18



ACGD/PDU
PAGE 20
HFD-S/PDU
PAGE 21
HFD/PDU
PAGE 22




RFS CABINET
PAGE 23
SYSTEM
CABINET
PAGE 24



RF ACCESSORY
CABINET
PAGES
28 & 29





NB RF AMP
CABINET
PAGE 30




BB RF AMP
CABINET
PAGE 32

4 PICTORIAL TABLE OF CONTENTS

FRONT COVER & PARTS (BACK VIEW):**FRONT END BELL (BACK VIEW):****NOTE: The following parts Not Shown:**

5113881*	3T CRM RF Shield and FRU Container
5109260	3T Eclipse Body Coil FRU (in case)
2319599	3T/94 RT Shim Kit
2370499	3T Eclipse CRM B0 Coil FRU


***2389049** RF Shield Change Out Kit – MUST be ordered for the RF Shield replacement

BLOWER CABINET


2121199-2
COMFORT MODULE
BLOWER



5106787
BODY COIL
BLOWER



SRI MODULE

5121793
SRI-3 MODULE
5123812
RUN 1260
5123813
RUN 1261
5115910
SRI-2 MODULE

COMFORT MODULE


2220079-2
COMFORT MODULE


2206437
M4 X 12mm NYLON
SCREW FOR
SPEAKER GRILL

3.0T HYBRID SPLITTER:**LONGITUDINAL MOTOR DRIVE PARTS:**

NOTE: The following REAR PEDESTAL parts Not Shown:


2362995	8 CHANNEL SWITCH (EXCITE II)
5109645-2	16 CHANNEL SWITCH (EXCITE HD)
5250136	REMOTE INTERCOM BOARD
46-226552G2	MAGNET INTERFACE BOARD


MAGNET DOCKING PARTS:**MAGNET AND ENCLOSURE (ASSEMBLIES)**

LOW PROFILE CARRIAGE ASSEMBLY – EXCITE II


NOTE: The following Part is Not Shown.


5114423 CABLE, RUN 1035, E2 HEAD RECEIVE

LOW PROFILE CARRIAGE ASSEMBLY – EXCITE HD**MAGNET AND ENCLOSURE (ASSEMBLIES)**




[illegible]

PENETRATION PANEL – EXCITE II

PENETRATION PANEL – EXCITE HD

PENETRATION PANEL

GRADIENT WATER CHILLER


2222564-2
60 Hz GRADIENT CHILLER

OR

2222564-21
50 Hz GRADIENT CHILLER

NOTE: The following WATER CHILLER parts are Not Shown:

2222564-6	LYTRON PUMP MOTOR
2222564-8	LYTRON 0.25A/250V FUSE
2222564-18	LYTRON 60Hz FAN MOTOR
2222564-19	LYTRON 50Hz FAN MOTOR
2222564-31	LYTRON NYLON MOTOR COUPLING
2138791	BODY COIL COOLANT, 1 GALLON





NOTE: The following parts are Not Shown.

2156158-50	S/W COLLECTOR FOR EXCITE HD (12.x) M4
2156158-43	S/W COLLECTOR FOR ECLIPSE 2 (11.x)
2385096	DVI-A TO 5-BNC MONITOR CABLE
2389540	SCAN TOOLS MOD


GOC COMPUTER CABINET WITH HP8200:

NOTE: The HP8200 PC has visible USB ports located on the right side of front of PC cover.

**GOC COMPUTER CABINET WITH HP8000:****NOTES:**

* If WIM is not available, order part 2141978-5 and replace cable between motherboard J1 and CIF J3.

** If replacing WIM with GOC AA, order cable 2395735.

LATEST VERSIONS:**KEYBOARDS ONLY:**

2275756 (English)
***2275756-2** (French)
***2275756-3** (German)
***2275756-4** (Scandinavian)
***2275756-5** (Italian)
***2275756-6** (Portuguese)
***2275756-7** (Spanish)

*** NOTE:**

These types of Keyboard/SCIM Assemblies were not available at the time this book was released, see below for earlier types.

NOTE: The following parts are Not Shown.

2262158-9	CABLE, 50 PIN SCIM KEYBOARD
2289697-2	SCIM ^W /ENGLISH OVERLAY
2289697-8	MR SCIM GLOBAL OVERLAY SET


NOTE: Refer to the Directions on CD–2160623 (Signa Service Methods) for additional product information.



TABLE 1–1

LOCATION	AMPS	VOLTS	GEMS PART NUMBER	MANUFACTURER	
				BUSS	G–S
FUSE F1	1.25	480	2276547–2	SNQ–R–1–1/4	ATQR1–1/4
FUSE F2	3	120	2276547–3	FNQ–R–3	ATQR3
FUSE F3	1	12VDC	2276547–8	LP–CC–1	ATDR1
FUSE F4	0.12	480	2276547–9	SNQ–R–8	ATQR1/8
FUSE F5	125	600	2276547–19	LPJ125SP	AJT125
FUSE F6	50	600	2276547–20	LPJ–50SP	AJT50
FUSE F7	20	600	2276547–22	LPJ–20SP	AJT20
FUSE F8	5	600	2276547–21	FRS–R–5	TRS5R
FUSE F9	8	400	2276547–6	FRS–R–8	TRS8R
FUSE F10	10	120	2276547–11	FNQ–R–10	ATQR10
FUSE F11	10	120	2276547–11	FNQ–R–10	ATQR10

2276547–14	12 VDC BATTERY
2276547–16	12 VDC SUPPLY

ACGD/PDU CABINET ASSEMBLY W/SGA/SGAPS/GP


NOTE: The following Part is Not Shown.

2261124-4 ACGD GRADIENT DATA PDL DOCUMENT

2374930-2 3T ACGD PDU CONTROL BOARD

HFD-S/PDU CABINET

NOTE: These parts are ONLY to be ordered if site has the HFD-S option.



NOTE: The following Part is Not Shown.

2303317 CAN Terminator

HFD-S/PDU CABINET

HFD/PDU CABINET

5116803
HFD CABINET
FAN ASM

2284074
ALUMINUM
AIR FILTER

5118444
HFD CABINET
WITH PDU


2331400-2
GP3 MODULE


2393906-2
HFA ASSEMBLY
2396050
HFA CONTROL BD

2334862-2
SGA PS ASSEMBLY

46-170021P87
SGA PS FUSES

2301200-3
TVSS MODULE
2374930-2
PDU CONTROL
BOARD




RFS CABINET (EXCITE HD ONLY)

NOTE: The following items are Not Shown.

5117087 RRF SERVICE DIAGNOSTIC KIT


46-296817P1 50 OHM BNC TERMINATOR

SYSTEM CABINET (EXCITE or EXCITE HD Upgrade)



NOTE: The following System Cabinet Options are Not Shown.

2325947-12	MINERVA AP (MASTER)
2325947-13	MINERVA AP (SLAVE)
2325947-14	MINERVA PIM
5117087	RRF SERVICE DIAGNOSTIC KIT
46-296817P1	50 OHM BNC TERMINATOR



NOTE: The following Array Processor ReFlex board Options are Not Shown.

2325947-12 R400 MINERVA AP (MASTER)

2325947-13 R400 MINERVA AP (SLAVE)

2325947-14 MINERVA PIM

2325947-15 R800 MINERVA AP (MASTER)

2325947-16 R800 MINERVA AP (SLAVE)

NOTE: The following MGD Module parts are Not Shown.

2294300-3 AGP BOARD

2294300-4 APS BOARD WITH IT DAUGHTER BD

2294300-5 BRIDGE BOARD

2294300-6 MGD CHASSIS W/BACKPLANE & AC ENTRY (No Boards)



2294300-7 MGD FAN

2294300-8 or

5181495-3 POWER SUPPLIES

2294300-9 INDICATOR BOARD

2294300-19 SCP2 BOARD

2305884 REMOTE RF MOD ASM (CHASSIS W/MIDPLANE) (FRONT)**2305884 REMOTE RF MOD ASM (CHASSIS W/MIDPLANE) (REAR)**


NOTE: Additional REMOTE RF MOD ASM Parts/Cables are located on page 27.

* A part that does not have a "-2" or "-3" can be replaced with either the same part number or the same number ending with a "-2" or "-3". If replacing the "-2" or "-3" version, then only the part number with a "-2" or "-3" can be used.

** 2274188 and 2274188-2 are replaced by 2274188-3. Also, all 2274188 series MUX boards can be replaced by backwards compatible MUX2 board, 5250106.

NOTE: The following RRF MODULE (page 26) Parts/Cables are Not Shown.

46-306565P1	CABINET FILTER
2294728	COMPACT PCI FRONT PANEL BLANK (8HP) (2 SLOT)
2319152	COMPACT PCI FRONT PANEL BLANK (4HP) (1 SLOT)
2308729	BNC TO SMB CABLE (MUX J8 TO CAB I/F J1)
2308728	8W8 CABLE ASSEMBLY (MUX J5)
2322578	J15/23/25/26 DATA INPUT CABLE (CAB I/F)
2322580	CABLE J18, UNBLANK (CAB I/F)
2322581	CABLE J24, MAGNET MONITOR (CAB I/F)
2322582	CABLE J30/31, SERIAL EXP (CAB I/F)
2322583	CABLE J62 (15 PIN SUB-D) (CAB I/F)
2322586	CABLE J61 (37 PIN SUB-D) (CAB I/F)
2315093	ETHERNET CABLE
2315093-2	CABLE J63 ETHERNET (CAB I/F)
2308729-2	BNC TO SMB CABLE (MUXJ13 TO CAB I/F J2)
2308729-3	BNC TO SMB CABLE (MUXJ14 TO CAB I/F J3)
2394879	RRF slot14R-J5 I/O J200
5110796	RRF slot14R-J13 I/O J202
5110795	RRF slot 14R-J13 I/O J202
5115417-3	RRF slot14R-J8 RF Amp J14
5115417	RRF slot12R-J31 Splitter-1
5113193	RRF slot12R-J7 RRF slot11R-J6/7
2394870-2	RRF slot5R-J5 I/O J201
5113193-2	RRF slot3R-J7 RRF slot2R-J6/7
5115417-2	RRF slot3F-J31 Splitter-2
5115630	J7 MUX (2) TO J6 MUX (1) CABLE
5115630-2	RRF slot14R-J11 RRF slot9R-J6



NOTE: The following RF ACCESSORY CABINET parts Not Shown:

46-306565P1 10" FRONT COVER FILTER

2365797 CABLE – I/O J127 TO ASC SLOT 11/14 J1/2/4/5

2368552 CABLE – ASC SLOT 14–J3 TO I/O J100

2369552-8 CABLE – ASC SLOT 13–J4 TO I/O J102


2369552-9 CABLE – ASC SLOT 13–J5 TO I/O J103

2369552-10 CABLE – ASC SLOT 10–J4 TO I/O J104

2369552-11 CABLE – ASC SLOT 10–J5 TO I/O J105

2379507 DRSM JUMPER PLUG

2383066 DRSM JUMPER CABLE (ASC–DRSM–I/O)


FOR 2363624-2 RF ACCESSORY CABINET WITH MNS:

* A part that does not have a "-2" can be replaced **ONLY** with the same no-dash part number. If replacing the "-2" version, then only the part number with a "-2" can be used.

NOTE: The following ASC MODULE Parts/Cables are Not Shown.

5115417-4	ASC slot14-J1 RF Amp J8
5115417-5	ASC slot14-J4 RF Amp J7
5115417-6	ASC slot11-J1 RF Amp J6
5120955	ASC POWER SUPPLY
5250032	3.0T 8KW BB2 AMP I/F ASM


25kW NB RF AMPLIFIER CABINET (for sites **WITHOUT** LCD Display)



NOTE: The following 25kW NB RF AMP CABINET (without LCD) parts are **Not Shown**:

<u>2336064-4</u>	<u>35KW RF AMP IPA TUBE & TUBE PULLER KIT</u>
<u>2336064-5</u>	<u>35KW RF AMP PA TUBE & HEX WRENCH KIT</u>
<u>5111003</u>	<u>55-50dB DUAL DIRECTIONAL BODY RF COUPLER FOR E2</u>
<u>5111004</u>	<u>45-40dB DUAL DIRECTIONAL HEAD RF COUPLER FOR E2</u>

[illegible]



NOTE: The following RF POWER AMP internal parts are **Not Shown**:

2206456-2	RF POWER AMP PREAMP/SPLITTER
2206456-3	RF AMPLIFIER (5)
2206456-4	RF POWER AMP COMBINER/COUPLER
2206456-5	RF POWER AMP COUPLER (4)
2206456-6	RF POWER AMP BIAS TEE
2206456-13	CABLE – SPLITTER TO AMPLIFIER (5)
2206456-20	COUPLER-BIAS TEE
2206456-21	POWER BIAS TEE
2206456-22	PMI INTERFACE KIT
46-170021P29	5 AMP FUSES (3)

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

ACCESSORIES

The following E-Catalog accessories are available for purchase with M-Catalog Options.

COLOR POSTSCRIPT PRINTER FOR A/W E7014LA

2112346 POSTSCRIPT COLOR NETWORK PRINTER

PATIENT ECG LEADWIRE SET E8811EC

2219452 HIGH IMPEDENCE ECG LEAD KIT

46-306867G1 PATIENT ECG CABLE ASM (E8811EC)

GENERIC REPLACEMENT CABLES

The table below provides information for generic replacement cables. If any of the cables are ordered, make sure to order write-on labels to be used for identifying cable part number and connection designations at the end of each cable.

Description	Length	
	65 Feet	100 Feet
9 Pin Sub-D Male/Female	46-244714G30	46-244714G37
9 Pin Sub-D Female/Female	46-244713G30	46-244713G37
15 Pin Sub-D Male/Female	46-244712G30	46-244712G37
15 Pin Sub-D Female/Female	46-244711G30	46-244711G37
25 Pin Sub-D Male/Female	46-244710G30	46-244710G37
25 Pin Sub-D Female/Female	46-244709G30	46-244709G37
37 Pin Sub-D Male/Female	46-258877G30	46-258877G37
	50 Feet	65 Feet
BNC RG 223	46-251710G1	46-251710G18
Cable Write-on Labels (60 labels per booklet)	46-208781P3	

COILS**3.0T 4 CHANNEL TORSO ARRAY COIL**

2366942	3T TORSO ARRAY ASSEMBLY
5110897 or -2	3T/94 EXCITE 4 CHANNEL TORSO COIL CABLE

3.0T 8 CHANNEL BRAIN COIL

2380637-2	BRAIN COIL/MIRROR/BASEPLATE ASM
2380637-3	PHANTOM POSITIONER
2380637-5	CABLE ASM

3.0T GP FLEX COIL

2382817	GP FLEX COIL
2389979	GP FLEX COIL CABLE KIT
46-2217351	PIN DIODE
2383558	HARDWARE KIT

3.0T QUAD EXTREMITY COIL

2380635-2	KNEE/FOOT ASM W/CABLE
2380635-5	KNEE/FOOT CABLE ASM
2357819-10	EXTREMITY REPAIR KIT
2357819-11	PHANTOM POSITIONER
2357819-12	PHANTOM

3.0T ECLIPSE HEAD COIL

46-282118G8	3T/94 HEAD COIL
2376104-2	HEAD QUICK DISCONNECT CABLE FOR ECLIPSE 2

3.0T NV ARRAY COIL

2383478-2	NV ARRAY W/COIL
2383477-3	NV COIL PHANTOM
2383477-4	NV COIL PHANTOM POSITIONER
2383478-5	NV CABLE ASM

3.0T 4 CHANNEL CTL COIL

2380634-2	3T/94 CTL ARRAY
2377326-3	USAI 8 CHANNEL COIL PHANTOM
5110896 or -2	3T/94 EXCITE 4 CHANNEL CTL ARRAY CABLE

SERVICE TOOLS

2258366	10cm SILICON PHANTOM
----------------	----------------------

OPTION KEYS:

2383494	Fiesta C
2390037	Asset
2390040	ConnectPro
2390075	3D FRFSE
2390078	EchoPlus
2390079	Diffusion Tensor
2390086	Sage 7
2390095	iDrive ProPlus
2390097	SmartPrep 2000
2390098	Flouro Trigger MRA
2390103	Performed Proc Step
5110228	Single Voxel Probe-P
5110229	Single Voxel Probe-S and Probe-P
5110230	Probe 2D – CSI
5110231	Probe 3D – CSI
5110232	MNS, SAGE, Probe-P
5110319	Research Proton Spectro
5125929	HFD-S Option Key

BrainWave OPTIONS:

2363727	15" LCD Display with 5 foot VGA Cable for Computer connection
2394034-5	(Cedrus) Response Pad Asm. Left & right hand units with two buttons for each hand. Includes fiber optic cable (20 ft) and optical converter box
2394034-6	(Cedrus) "Lumina" response controller box
2394034-7	(Cedrus) AC Power module for Lumina Controller Box
2394034-9	FRU usage only Kit containing all cables [Cedrus-4, -5, -6, -7, 2xxxxxx-6, 46-265505P1 (qty 2)]
2394034-10	Penetration Sub-Panel Asm
2394034-16	PC compatible Keyboard (narrow)
2394034-17	PC compatible Mouse
5110767	BrainWave Application Software CD-ROM
5110736	EXCITE II Scanner License Key MOD for the BrainWave H/W Option
5110946	PC and Platform Software Asm
5110968	Patch CD for BrainWave Stimulus PC