

GE Healthcare

Signa HDx 1.5T Series
Field Replaceable Units (FRU)

Technical Publication

© 2006 – 2011 General Electric Company
GE Medical Systems, a General Electric Company, going to market as GE Healthcare
3000 N. Grandview Boulevard
Waukesha, Wisconsin 53188
USA

www.gehealthcare.com



5166012
Rev 14

**COVER
BOTTOM
CROP LINE**

**COVER
BOTTOM
CROP LINE**

**COVER
FOLD
LINE** | **(2) CENTER STAPLES
FOR BINDING**

**COVER
RIGHT
CROP LINE**

**5/16" RADIUS
CORNER TRIM
TYPICAL**

**5/16" RADIUS
CORNER TRIM
TYPICAL**

**COVER
LEFT
CROP LINE**

SEND COMMENTS ABOUT THIS MANUAL TO:

GE MEDICAL SYSTEMS
MANAGER, SERVICE ENGINEERING
W-808
3200 N. GRANDVIEW BLVD.
WAUKESHA, WISCONSIN 53188
USA

This manual belongs to:

INTRODUCTION TO FIELD REPLACEABLE UNIT (FRU)

This direction is a guide to the FRU's for HDxt, HDx, HDi, and Vibrant.

Parts designated as a FRU are critical to maintain the safety and/or functionality of the system. FRU's are available on a Priority 1 basis from Service Parts Operation.

As a reference aid, the FRU's are listed under an upper-level part number where they are used. In most cases the upper-level part number is not a FRU. In some cases the parts listed under an upper-level assembly may have been nested in intermediate sub-assemblies which are not shown.

Where possible, FRU's which are part of Option Catalogs are shown under the area where they would be installed, rather than being listed in the Coils or Options Sections (SEE TABLE OF CONTENTS).

LIST OF EFFECTIVE PAGES


<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>
Cover Page .	14	23	11	31	11
1 to 3	14	24	14	32 to 39	13
4 to 8	11	25	12	40 to 42	11
9	14	26	14	43 to 48	7
10 to 21	11	27 to 28	11	49 to 58	11
22	14	29 to 30	14	59 to 64	14

TITLE	PAGE
INTRODUCTION, LIST OF EFFECTIVE PAGES	1
TABLE OF CONTENTS	2
PICTORIAL TABLE OF CONTENTS	3
MAGNET AND ENCLOSURE (ASSEMBLIES)	
BRM & TRM LIST (OF PARTS NOT SHOWN)	5
FRONT COVER, REAR END BELL	6
ENCLOSURE & VACUUM RELATED PARTS	7
FRONT & REAR ENDBELL, BRIDGE SUPPORT	7
REAR PEDESTAL PARTS	8, 9
SPLIT BRIDGE & DOCKING PARTS	10
LPCA	11
SRI/PAC & TOUCH-UP PAINTS	12
MONITORS, MRU, & TRANSFORMER	13
HDx MAGNET ENCLOSURE	14
HDx BRIDGE PARTS	16
HDx LPCA	17
CHILLER	18
PATIENT SUPPORT TABLE	20, 21
PENETRATION PANEL	22
OPERATOR WORKSPACE	24
GOC COMPUTER CABINET, SCIM AND KEYBOARD	25
MAIN DISCONNECT PANEL	28
HFD-S/PDU CABINET	29
HFD/PDU CABINET	30
RFS CABINET	32, 33, 34
SYSTEM CABINET	35
CAM, MGD, RRF, & ASC MODULES	36, 37, 38, 39, 40
HDx 32 CHANNEL CABINET (OPTION)	41
SRFD CABINET	43
SRFD2 CABINET	44
TWIN ACCESSORY CABINET (TAC)	45, 46, 47
REPLACEMENT CABLES	49, 50
SERVICE TOOLS	51
UNIFIED SURFACE COIL PHANTOM SET	52
TOOLS, REPAIR KITS, AND HEAD COILS	53
COILS	54 THRU 59
OPTIONS	60
SOFTWARE OPTION KEYS	61
BrainWave OPTIONS	62
REVISION HISTORY	63

INTRO TO FRU'S	PAGE 1
MONITORS	PAGE 13
CHILLER PARTS	PAGE 18
CABLES	PAGE 49
COILS	PAGE 55
OPTIONS	PAGE 60
S/W OPTIONS	PAGE 61
BrainWave	PAGE 62




PATIENT AIR SUPPLY
BLOWER CABINET (MG6)
PAGE 9



PENETRATION
PANEL
PAGE 22

PATIENT
SUPPORT
TABLE
PAGE 20




WideOpen MAGNET ENCLOSURE:
PGS 5 thru 12
HDx MAGNET ENCLOSURE:
PGS 14 thru 17

TRM COIL,
MOUNTING
BRACKETS, CART,
AND CRADLE,
SEE PAGE 5

MAGNET ROOM

OPERATOR/EQUIPMENT ROOM




SERVICE TOOLS
PAGE 51




OPERATOR
WORKSPACE
PAGES 24 & 25




MAIN
DISCONNECT
PANEL
PAGE 28




32 CHANNEL
CABINET
PAGE 41




HFD-S/PDU
PAGE 29
HFD/PDU
PAGE 30



RFS CABINET
PAGE 32 to 34
SYSTEM CABINET
PAGE 35



SRF CABINET
PAGE 43
SRFD2 CABINET
PAGE 44



ACCESSORY
CABINET
PAGES 45
THRU 47

PICTORIAL TABLE OF CONTENTS

[illegible]

NOTE: The following WideOpen & TwinSpeed common parts Not Shown:

5123812	BORE TEMPERATURE SENSOR, RUN 1260
5123813	OPERATOR DISPLAY CABLE TO SRI-3, RUN 1261
5133997	IIQ KIT, contains:
	BRIDGE SUPPORT FOR LCC MAGNET ENCLOSURES
	BRIDGE SUPPORT FOR S-SERIES MAGNET ENCLOSURES
	RADIAL SUPPORTS FOR ALL MAGNETS
5184777	RF SHIELD GROUNDING KIT
2324086-2	5 GALLON PAIL w/ETHYLENE GLYCOL WARNING LABEL

NOTE: The BRM parts Not Shown:

2379255	FOAM, GRAY, 116" LG. (Inserted between BRM and bore)
2297672	GLYCOL, 1 GALLON/3.8 LITERS
2124499-5	GRAD COIL LCC ONE-PIECE SHIELD
2144093-3	BRM GRADIENT COIL W/CART & CRADLE (Cx & S-series)
2144093-5	BRM-D GRADIENT COIL W/CART & CRADLE (LCC)
2241521	SPLICE KIT GRADIENT COIL SPLICE KIT
5110336	BRM GRADIENT CABLE LEAD BOARD KIT
5127441-2	GRADIENT CABLE CLAMP BLOCK KIT
2208999-3	SHORT RF COIL & CASE (Order for RF COIL Replacement-LCC)
2208999-6	LONG RF COIL & CASE (Order for RF COIL Replacement-Cx & S-series)
2208999-10	HDx new color SHORT RF COIL & Yellow/Black CASE (Order for RF COIL Replacement-LCC)
5113879*	RF SHIELD FOR BRM COIL



NOTE: The following TwinSpeed 1.5T FIXED W/TRM & Enclosure related parts Not Shown; MISCELLANEOUS:


2297760	VACUUM GREASE, MSDS# 8366442
2203803	SPLIT O-RING
2242756	RELIEF VALVE, 3/8 NPT X 7 PSIG, GEN
2299843	MAGNET PRESSURE TRANSDUCER, SUNX DP2-41N
46-294231G1	MAGNET RUNDOWN UNIT
46-294231P1	BUSSMAN, F1, 110VAC .400 GDC


TRM:




5113878*	RF SHIELD & FRU CONTAINER FOR TRM COIL
2144093-7	TRM WITH MOUNTING BRACKETS, CART, AND CRADLE
2144093-9	TRM WITH MOUNTING BRACKETS, CART, AND CRADLE
2218148	COLD SHRINK (FOR REAR ENDBELL AND TRM COIL)
2208999-8	RF BODY COIL AND CASE
2208999-11	HDx new color Twin Body SHORT RF COIL & Silver/Black CASE (Order for RF COIL Replacement-LCC)
2297672	GLYCOL, 1 GALLON/3.8 LITERS
2299498	COOLING CIRCUIT REPAIR KIT
2303520	DIODE SLOT BOARD
5127441	GRADIENT CABLE CLAMP BLOCK KIT
2293612	MALE CONNECTOR 3/4" HOSE TO 3/4" NPT


***2389049** RF Shield Change Out Kit – MUST be ordered for the RF Shield replacement

WideOpen FRONT COSMETIC COVER & PARTS (BACK VIEW):




REAR PEDESTAL**Note:** See next page for Hybrid Splitter information**8 MAGNET AND ENCLOSURE (ASSEMBLIES)**


REAR PEDESTAL**1.5T HYBRID SPLITTER – TELEDYNE:****1.5T HYBRID SPLITTER – NETCOM:****BLOWER CABINET****MAGNET AND ENCLOSURE (ASSEMBLIES)**

WideOpen BRIDGE PARTS:

NOTE: The following part is Not Shown.


2386348	FRONT BRIDGE FIBERGLASS
5117304	FRONT BRIDGE PULLEY ASM
2293538	BRIDGE EDGE FOAM
46-282436P1	CARRIAGE TRAVEL SENSOR FLAG

MAGNET DOCKING PARTS:


HD LOW PROFILE CARRIAGE ASM**MAGNET AND ENCLOSURE (ASSEMBLIES)**

SRI/PAC PLATFORM:

5112668
NON-LCC
PAC-2B MODULE



5122500
ECG LEAD KIT



5105270
SRI-3
MODULE

5121793
SRI-3 U/G KIT
w/CABLES


5112668-2
PAC-2B MODULE

TOUCH-UP PAINTS:


2183491-3 GRAY #1 0.6 OZ. NON-AEROSOL

46-303460P3 STEEL BLUE 0.6 OZ. NON-AEROSOL

46-303460P4 DARK STEEL BLUE 0.6 OZ. NON-AEROSOL




NOTE: To obtain enclosure part numbers not found in this manual, refer to Direction 2138247, *Signa Illustrated Parts List*, for enclosure covers, endbells, and other parts.

CRYOGEN PANEL STEP-DOWN TRANSFORMER FRU'S:**MAGNET RUNDOWN UNIT (MRU) FRU'S:**


5196918	MAGNET RUNDOWN UNIT (Primary)
5196919	MAGNET RUNDOWN UNIT (Remote)
5363631	SPARE FUSE KIT includes: 1.25A, 250V Fuse; 2.5A, 250V Fuse; 400ma, 120V Fuse; 200ma, 250V fuse
5199490	BATTERY PACK ASM
5180324	ADAPTER PLATE
5196921	REMOTE TO PEN PANEL CABLE
5196920	PRIMARY TO PEN PANEL CABLE
2128699	MAGNET TO MRU CABLE (Run 606)

OXYGEN MONITOR (O2) FRU'S:

5390676	OXYGEN MONITOR FRU – Contains all of the parts included with the O2 monitor, box, sensor and cables
5399042	OXYGEN MONITOR CALIBRATION KIT – Contains all of the parts needed to calibrate either the nw or old O2 monitor
5390678-3	REMOTE 5-YEAR SENSOR – Sensor for the ISA-60M O2 monitor


MAGNET MONITOR:


2219400
MAGNET MONITOR (MSM)

HDx MAGNET ENCLOSURE

BRM SPECIFIC HDx ENCLOSURE PARTS	
5142289	Gradient Cable Lead Board BRM
5142288	Gradient Cable Lead Board Cover BRM
5160993	Front End Bell BRM Assembly
5160995	Rear End Bell BRM Assembly
5161267	Feed Thru Cover Plate, Rear End Bell BRM
5143732	Rear Cover Assembly, 1.5T BRM
5142278	Rear Turtle Arm BRM
5142278-2	Rear Turtle Arm with Notch BRM
5161263	Turtle Bracket Gasket BRM
2222549	BRM Lumitex Light Panel, 726mm
TRM SPECIFIC HDx ENCLOSURE PARTS	
5161262	Turtle Bracket Gasket TRM
5150552	Rear Turtle Arm, TRM Single Boss
5150552-2	Rear Turtle Arm, TRM Double Boss
5150552-3	Rear Turtle Arm, TRM Single Boss with Relief Cut
5160997	Front End Bell TRM Assembly
5160999	Filler Cover Left, Front End Bell TRM
5160999-2	Filler Cover Right, Front End Bell TRM
5161000	Rear End Bell TRM Assembly
5161003	Filler Cover Left, Rear End Bell TRM
5161003-2	Filler Cover Right, Rear End Bell TRM
5143732-2	Rear Cover Assembly, TRM
5146696	TRM Comfort Module Assembly
2222549-4	Twin Lumitex Light Panel, 1085mm
Vibrant SPECIFIC HDx ENCLOSURE PARTS	
46-271407G4	Door Pedal Assembly with Bearings
5151727-2	Front Right Skirt Cover Assembly Vibrant
5151729-2	Front Left Skirt Cover Assembly Vibrant
5151731-2	Rear Right Skirt Cover Assembly Vibrant
5151732-2	Rear Left Skirt Cover Assembly Vibrant
5309702	Vibrant Front Right Arc Cover Assembly with Graphic
5309703	Vibrant Front Left Arc Cover Assembly with Graphic
5309704	Vibrant Rear Right Arc Cover Assembly with Graphic
5309705	Vibrant Rear Left Arc Cover Assembly with Graphic

BRM & TRM COMMON HDx ENCLOSURE PARTS	
5164786	HDx Enclosure Hardware Spares Kit
5151721	Front Right Arc Cover Assembly
5151722	Front Left Arc Cover Assembly
5151723	Rear Right Arc Cover Assembly
5151725	Rear Left Arc Cover Assembly
5151718-2	Right Side Top Cover, White
5151771-2	Right Side Bottom Cover, White
5161378-2	Right Latch Hook
5151720-2	Left Side Bottom Cover, White
5161378	Left Latch Hook
5151734-2	Left Side Top Cover, White
5151729	Front Left Skirt Cover
5151727	Front Right Skirt Cover
5151731	Rear Right Skirt Cover
5151732	Rear Left Skirt Cover
5143554	Left Rear Pedestal Cover
5143555	Right Rear Pedestal Cover
2215716-51	Center Display Panel With Overlay
2215716-10	Center Display Panel
2215716-60	Left Control Panel with all Overlays
2215716-11	Left Control Panel
2215716-70	Right Control Panel with all Overlays
2215716-12	Right Control Panel
5160276	Bracket I and Q Lines
5150097	Turtle Bracket Front Arm
5160991	Dock Bracket Assembly
5150410	Front Trim Ring Assembly
5150411	Rear Trim Ring Assembly
5143642	Front Cover Assembly
5141500	Microphone Block and Cable Assembly
5142985	Vacuum Port Cover Assembly
5149066	Remote PAC Base
5149068	Remote PAC Inner Mounting Tray
5149065	Remote PAC Accessory Hook
2344509-2	Bore Light End Cap
5160206	Side Cover Mounting Bracket Assembly RH
5160206-2	Side Cover Mounting Bracket Assembly LH
5160992	Rear Bracket Assembly, Left
5160992-2	Rear Bracket Assembly, Right
2206437	M4 x 12mm Nylon Screw for Speaker Grill


HDx BRIDGE PARTS:

HDx LOW PROFILE CARRIAGE ASM**MAGNET AND ENCLOSURE (ASSEMBLIES)**

SYSTEM COOLING CABINETS:

REFER TO ELLIS & WATTS SERVICE
MANUALS FOR FRU INFORMATION

REFER TO AIRSYS
SERVICE MANUALS
FOR FRU INFORMATION



2343325 AIRSYS HEAT EXCHANGER (GWHX)

2341077-2 AIRSYS QUICK DISCONNECT KIT

2341077-3 AIRSYS REMOTE CONTROL PANEL

2230235 E & W MODEL LTL-4 MOBILE HEAT EXCHANGER

2230235-3 E & W PUMP/MOTOR SET

2230235-12 E & W LEVEL SWITCH

2230235-24 QUICK-DISCONNECT FACILITY

2230235-25 QUICK-DISCONNECT COIL

2298100-4 E & W MODEL LTL10 HEAT EXCHANGER CABINET

2298100-5 CONTROL PANEL FOR LTL10 ONLY

2297672 GLYCOL USED WITH E & W OR AIRSYS

NOTE: The following parts are **Not Shown**.

2222564-22 LYTRON CHILLER, 50Hz OR 60Hz

2222564-23 LYTRON CIRCUIT BREAKER

2222564-24 LYTRON PUMP

2222564-25 LYTRON MOTOR

2222564-26 LYTRON KEY PAD

2222564-27 LYTRON FLOW SWITCH

2222564-28 LYTRON HOSE FITTINGS


2222564-29 LYTRON PUMP with NYLON COUPLING


2222564-31 LYTRON NYLON MOTOR COUPLING

2260427 COOLANT (1 GALLON) (LYTRON ONLY)

PATIENT TABLE – INFINITY and HDx

NOTE: Refer to page 21 for information on bumper FRU's for HDx, Infinity, and Pre-Infinity patient tables.




PATIENT TABLE – BUMPERS AND PANELS


HDx EXISTING PART	HDx REPLACEMENT FRU
5145946, Panel – Left 5145947, Panel – Right	5179590 , New panels and side bumpers FRU kit for HDx table
5145939, 22 Side Bumper	5161172 , Side Bumper 22, Color P17
5145939–2, 16.25 Side Bumper	5161172–2 , Side Bumper 16.25, P17
5145941, Foot End Bumper	5145941 , Foot End Bumper
5145935, Head End Bumper	5145935 , Head End Bumper

INFINITY EXISTING PART	INFINITY REPLACEMENT FRU
46–271022G1, 2393106, Panel – Left 46–271021G1, 2393106–2, Panel – Right	5179706 , New panels and side bumpers FRU kit for Infinity table
5154908–4, Panel – Left 5154908–3, Panel – Right	5179706 , New panels and side bumpers FRU kit for Infinity table
46–271025P5, 22 Side Bumper	5161172–3 , Side Bumper 22, Color P17
46–271025P6, 16.25 Side Bumper	5161172–4 , Side Bumper 22, Color P17
46–271027P3, Foot End Bumper	46–271027P3 , Foot End Bumper
46–271026P3, Head End Bumper	46–271026P3 , Head End Bumper

Pre-INFINITY EXISTING PART	Pre-INFINITY REPLACEMENT FRU
5154908–4, Panel – Left (Grey) 5154908–3, Panel – Right (Grey)	5179706 , New panels and side bumpers FRU kit for Infinity table
22 Side Bumper (Accent Grey or Blue)	5161172–3 , Side Bumper 22, Color P8
16.25 Side Bumper (Accent Grey or Blue)	5161172–4 , Side Bumper 16.25, Color P8
Foot End Bumper (Accent Grey or Blue)	46–271027P3 , Foot End Bumper
Head End Bumper (Accent Grey or Blue)	46–271026P3 , Head End Bumper




[illegible]

OPERATOR WORKSPACE – HDxt

* When ordering a monitor replacement, be aware that models may be updated without notice.

NOTE: The following parts are Not Shown.


2156158–73	S/W COLLECTOR FOR 15.0 M4
2156158–82	S/W COLLECTOR FOR 15.0 M4A
2156158–84	S/W COLLECTOR FOR 15.0 M4B
2156158–86	S/W COLLECTOR FOR 16.0 M3
2156158–89	S/W COLLECTOR FOR 16.0v2
2156158–91	S/W COLLECTOR FOR HD23.0
2389540	SCAN TOOLS MOD

GOC COMPUTER CABINET WITH HP Z400:

* The Z400 Host computer is only compatible with the 15M4B software or higher.

NOTE: The following parts are Not Shown.

- 5373029** Cable FRU kit for GOC contains:
Run E3016, Run E3046, Run E3047, Run E3048, Host to IRDC DVI-D cable, Host to IRDC USB cable, IRDC Power cable
- 5376172** Dust Filter FRU Kit for Z400
- 5364704** One Wire Collector for Z400

SCIM & KEYBOARD:**NOTE:**


If Keyboard is to be replaced, also order **46-317475P2**, 2 inch x 4 inch adhesive backed velcro to attach to bottom of new keyboard.

KEYBOARDS ONLY:

- 5143798-7** (English)
5143798-8 (French)
5143798-9 (German)
5143798-10 (Spanish)
5143798-11 (Italian)
5143798-12 (Portuguese)
5143798-13 (Scandinavian)

NOTE: The following parts are Not Shown.

- 2262158-9** CABLE, 50 PIN SCIM KEYBOARD
- 5148389-2** SCIM ^W/ENGLISH OVERLAY (Black)
- 5148389-9** MR SCIM GLOBAL OVERLAY SET (Black)

GOC COMPUTER CABINET WITH HP9300 OR HP8400:

* **5143798-14** is only compatible with 14M5 software or later

** Sites with SCSI Tower **MUST** order 5183547-41 as a replacement.
Sites without a SCSI Tower can order 5183547-41 or 5183547-60 as a replacement.


NOTE: The following parts are Not Shown.

5151612 USB CABLE FOR OPERATOR WORKSPACE

5141564 AUTOVOICE CABLE

[illegible]

NOTE: Refer to the Directions on applicable MR Service Methods DVD for additional product information.



NOTE: The following items are Not Shown.

5338648 Upgrade kit for MDP Battery and Charger


5375947 Power Sonic Battery FRU

5316807 Power Supply FRU

LOCATION	AMPS	VOLTS	GEMS PART NUMBER	MANUFACTURER	
				BUSS	G-S
FUSE F1	1.25	480	2276547-2	SNQ-R-1-1/4	ATQR1-1/4
FUSE F2	3	120	2276547-3	FNQ-R-3	ATQR3
FUSE F3	1	12VDC	2276547-8	LP-CC-1	ATDR1
FUSE F4	0.12	480	2276547-9	SNQ-R-8	ATQR1/8
FUSE F5	125	600	2276547-19	LPJ125SP	AJT125
FUSE F6	50	600	2276547-20	LPJ-50SP	AJT50
FUSE F7	20	600	2276547-22	LPJ-20SP	AJT20
FUSE F8	5	600	2276547-21	FRS-R-5	TRS5R
FUSE F9	8	400	2276547-6	FRS-R-8	TRS8R
FUSE F10	10	120	2276547-11	FNQ-R-10	ATQR10
FUSE F11	10	120	2276547-11	FNQ-R-10	ATQR10


2276547-14 12 VDC BATTERY


2276547-16 12 VDC SUPPLY

HFD-S/PDU CABINET ASSEMBLY W/SGA/SGAPS/GP

NOTE: The following Part is Not Shown.

2303317 CAN Terminator

HFD/PDU CABINET

HDx RFS CABINET with CAM Module

- * If a cabinet contains a 5127452 ICN, a replacement 5127452 must be ordered. If a 5127452-3 is in the cabinet, a 5127542-3 must be ordered.
- ** The SUN 4170 ICN is only compatible with 15M4B software or higher
- ***The ORACLE/SUN 4170M2 ICN is only compatible with 16.0v02 software or higher.

CAM MODULE SEE PAGE 36

5139592
16 PORT
ETHERNET
SWITCH

OR 5162073
24 PORT ETHERNET SWITCH
(32 CHANNEL OPTION)

5127452*
IMAGE
COMPUTE
NODE

OR 5337894-3**
SUN 4170 ICN
(16 Channel)

OR 5911000-4***
ORACLE/SUN
4170M2 ICN
(8 & 16 Channel)

OR 5127452-3*
SUN 4100 ICN
5127452-4
SUN HDD

or 5337894-4**
SUN 4170 ICN
(32 Channel)
5337894-5
4170 HDD

or 5911000-5***
ORACLE/SUN
4170M2 ICN
(32 Channel)
5337894-5
4170 HDD

**REMOTE RF MODULE
SEE PAGE 38**

5110612
DRIVER
MODULE

2351573
RF AMP FRU

5120780
TERM SERVER
CABLE FRU KIT
(4 CABLES)


2379186
TERMINAL
SERVER

5115759
SINGLE PORT
TERM SERVER
(OPTION)

5149828
6 OUTLET
POWER STRIP

5143103
AC DISTRIBUTION
BOX

5125176
UTNS
SPLITTER BOX




NOTE: The following items are Not Shown.

5272613 NMR SERVICE CABLE KIT

46-296817P150 OHM BNC TERMINATOR

5306523 DVMR SERVICE CABLE KIT (Supports IRF2b, RRF-DIF3 & later
(Order if DVMR Loopback adaptor is needed; Supports HD 16.0v02 or later))

HDx RFS CABINET with MGD Module

- * If a cabinet contains a 5127452 ICN, a replacement 5127452 must be ordered. If a 5127452-3 is in the cabinet, a 5127542-3 must be ordered.
- ** The SUN 4170 ICN is only compatible with 15M4B software or higher
- ***The ORACLE/SUN 4170M2 ICN is only compatible with 16.0v02 software or higher.

MGD MODULE
 SEE PAGE 37

5139592
 16 PORT
 ETHERNET
 SWITCH

OR 5162073
 24 PORT ETHERNET SWITCH
 (32 CHANNEL OPTION)

5127452*
 IMAGE
 COMPUTE
 NODE
OR
5337894-3**
 SUN 4170 ICN
 (16 Channel)

5911000-4***
 ORACLE/SUN
 4170M2 ICN
 (8 & 16 Channel)

OR
5127452-3*
 SUN 4100 ICN
5127452-4
 SUN HDD

or
5337894-4**
 SUN 4170 ICN
 (32 Channel)
5337894-5
 4170 HDD

OR
5911000-5***
 ORACLE/SUN
 4170M2 ICN
 (32 Channel)
5337894-5
 4170 HDD

REMOTE RF MODULE
 SEE PAGE 38

5110612
 DRIVER MODULE

2351573
 RF AMP FRU



5120780
 TERM SERVER
 CABLE FRU KIT
 (4 CABLES)

2379186
 TERMINAL
 SERVER

5115759
 SINGLE PORT
 TERM SERVER
 (OPTION)


5149828
 6 OUTLET
 POWER STRIP

2363041-2
 ASC CHASSIS
 SEE PAGE 40


5125176
 UTNS SPLITTER BOX

5143103
 AC DISTRIBUTION
 BOX

NOTE: The following items are Not Shown.
5272613 NMR SERVICE CABLE KIT

46-296817P150 OHM BNC TERMINATOR

5306523 DVMR SERVICE CABLE KIT (Supports IRF2b, RRF-DIF3 & later
 (Order if DVMR Loopback adaptor is needed; Supports HD 16.0v02 or later)

HD RFS CABINET UPGRADED TO HDx**NOTE: The following items are Not Shown.****5272613** NMR SERVICE CABLE KIT**46-296817P150** OHM BNC TERMINATOR**5306523** DVMR SERVICE CABLE KIT (Supports IRF2b, RRF-DIF3 & later

(Order if DVMR Loopback adaptor is needed; Supports HD 16.0v02 or later)

EXISTING SYSTEM CABINET UPGRADED TO HDx


NOTE: The following items are Not Shown.

5272613 NMR SERVICE CABLE KIT

46-296817P150 OHM BNC TERMINATOR

5306523 DVMR SERVICE CABLE KIT (Supports IRF2b, RRF-DIF3 & later
(Order if DVMR Loopback adaptor is needed; Supports HD 16.0v02 or later))


SYSTEM CABINET



* For 32 Channel systems, use **2358122** IRF2 Bd.

** Replace the following boards only with the same P/N board:
2358122-2, 5250066-3, 2274196, and 2274196-2.

NOTE: In the MGD, CAM, and RRF Chassis, the RRF-DIF3 (6250184), IRF2b (5250066-3), and IRF-I/O2 (2274196-2) Boards only work with each other in the system. Do not intermingle other P/N boards with these. Also, they are only compatible with 16.0v02 software and higher.





* For 32 Channel systems, use **2358122** IRF2 Bd.

** Replace the following boards only with the same P/N board:
2358122-2, 5250066-3, 2274196, and 2274196-2.

NOTE: In the MGD, CAM, and RRF Chassis, the RRF-DIF3 (6250184), IRF2b (5250066-3), and IRF-I/O2 (2274196-2) Boards only work with each other in the system. Do not intermingle other P/N boards with these. Also, they are only compatible with 16.0v02 software and higher.

NOTE: The following MGD Module parts are **Not Shown**.

2294300-5	BRIDGE BOARD
2294300-6	MGD CHASSIS W/BACKPLANE & AC ENTRY (No Boards)
2294300-7	MGD FAN
2294300-8 or	
5181495-3	POWER SUPPLIES
2294300-9	INDICATOR BOARD

2305884 REMOTE RF MOD ASM (CHASSIS W/MIDPLANE) (FRONT)**2305884 REMOTE RF MOD ASM
(CHASSIS W/MIDPLANE) (REAR)**

NOTE: In the MGD, CAM, and RRF Chassis, the RRF-DIF3 (6250184), IRF2b (5250066-3), and IRF-I/O2 (2274196-2) Boards only work with each other in the system. Do not intermingle other P/N boards with these. Also, they are only compatible with 16.0v02 software and higher.

NOTE: Additional REMOTE RF MOD ASM Parts/Cables are located on page 39.

* A part that does not have a "-2" can be replaced with either the same part number or the same number ending with a "-2". If replacing the "-2" version, then only the part number with a "-2" can be used.

** 2274188 and 2274188-2 are replaced by 2274188-3. Also, all 2274188 series MUX boards can be replaced by backwards compatible MUX2 board, 5250106.

*** For 32 Channel systems, use 5250144 Exciter 2.1, in main RFS Cabinet only (not in 32 Channel cabinet)

**** Replace the following boards only with the same P/N board:
2274184, 5250028, 5250028-2, and 6250184.



NOTE: The following RRF MODULE (page 38) Parts/Cables are Not Shown.

46-306565P1	CABINET FILTER
2294728	COMPACT PCI FRONT PANEL BLANK (8HP) (2 SLOT)
2319152	COMPACT PCI FRONT PANEL BLANK (4HP) (1 SLOT)
2308729	BNC TO SMB CABLE (MUX J8 TO CAB I/F J1)
2308728	8W8 CABLE ASSEMBLY (MUX J5)
2322578	J15/23/25/26 DATA INPUT CABLE (CAB I/F)
2322580	CABLE J18, UNBLANK (CAB I/F)
2322581	CABLE J24, MAGNET MONITOR (CAB I/F)
2322582	CABLE J30/31, SERIAL EXP (CAB I/F)
2322583	CABLE J62 (15 PIN SUB-D) (CAB I/F)
2322586	CABLE J61 (37 PIN SUB-D) (CAB I/F)
2308729-2	BNC TO SMB CABLE (MUXJ13 TO CAB I/F J2)
2308729-3	BNC TO SMB CABLE (MUXJ14 TO CAB I/F J3)
2394870-3	RRF slot14R-J5 I/O J200
5110796	RRF slot14R-J13 I/O J202
5110795	RRF slot 14R-J13 I/O J202
5115417-3	RRF slot14R-J8 RF Amp J14
5115417	RRF slot12R-J31 Splitter-1
5113193	RRF slot12R-J7 RRF slot11R-J6/7
2394870-2	RRF slot5R-J5 I/O J201
5113193-2	RRF slot3R-J7 RRF slot2R-J6/7
5115417-2	RRF slot3F-J31 Splitter-2
5115630	J7 MUX (2) TO J6 MUX (1) CABLE
5115630-2	RRF slot14R-J11 RRF slot9R-J6
5116213-2	RF Extension cable to Run 1238
5121269	RF Extension cable to Run 1259
5145157	Fiber Optic cable, MGD Slot 11R-J20 to RF Slot 10F-J21
5145157-2	Fiber Optic cable, MGD Slot 11F-J7 to RF Slot 1F-J21
5150527	Y Power cable, 10 feet
5160775	HDx, FRU ICN Ethernet Cable Collector


5429284	Fiber Optic Cable FRU Kit for RRF-DIF3, IRF2b & IRF-I/O2
	contains:
	5257686-2, MGD Slot 11F J10 to RF Slot 10F J21
	5416142-2, MGD Slot 11F J8 to RF Slot 10F CLK
	5257686-3, MGD Slot 11F J11 to RF Slot 1F J21

NOTE: The following ASC MODULE (page 40) Parts/Cables are Not Shown.


5115417-4	ASC slot14-J1 RF Amp J8
5115417-5	ASC slot14-J4 RF Amp J7
5115417-6	ASC slot11-J1 RF Amp J6
5120955	ASC POWER SUPPLY
2391148	1.5T BB2 AMP I/F ASM


2363041 or 2363041-2* ASC CHASSIS ASSEMBLY**2363041 or 2363041-2* ASC CHASSIS ASSEMBLY WITH MNS****NOTES:**


* Additional ASC MODULE Parts/Cables are located on page 39.

HDx 32 CHANNEL CABINET (Option)

[illegible]


EXISTING SRF CABINET UPGRADED TO HDx**SRF CABINET**


EXISTING SRFD II CABINET UPGRADED TO HDx

2239500 TAC

NOTE: The following parts are Not Shown.

2282045-2	VACUUM PUMP TIP SEAL KIT
2293608	QF25 CENTERING RING & O-RING
2293608-2	QF16 CENTERING RING & O-RING
2293609	QF25 VACUUM ELBOW
2293616	QF25 ADAPTER TO ONE INCH HOSE
2293617	VACUUM REDUCING 'T'
2293618	QF16 CLAMP
2293618-2	QF25 CHAIN CLAMP
2294010-3	HOSE CLAMP, SST ONE INCH
2297760	VACUUM GREASE, MSDS# 8366442
2329620-4	VACUUM GAUGE CONTROLLER CABLE
2284301	VACUUM GAUGE TO TAC J2 CABLE


5118162 TAC



NOTE: The following parts are Not Shown.

2282045-2	VACUUM PUMP TIP SEAL KIT
5118818	CLAMP NW 20-25
2297760	VACUUM GREASE, MSDS# 8366442

5127147
SERVICE KIT FOR TAC VACUUM
SYSTEM UPGRADE
(Used to replace old style vacuum
system, with solenoid)



TAC VACUUM SYSTEM

[illegible]

The table below provides information for generic replacement cables. If any of the cables are ordered, make sure to order write-on labels to be used for identifying cable part number and connection designations at the end of each cable.

GENERIC REPLACEMENT CABLES

Description	Length	
	65 Feet	100 Feet
9 Pin Sub-D Male/Female	46-244714G30	46-244714G37
9 Pin Sub-D Female/Female	46-244713G30	46-244713G37
15 Pin Sub-D Male/Female	46-244712G30	46-244712G37
15 Pin Sub-D Female/Female	46-244711G30	46-244711G37
25 Pin Sub-D Male/Female	46-244710G30	46-244710G37
25 Pin Sub-D Female/Female	46-244709G30	46-244709G37
37 Pin Sub-D Male/Female	46-258877G30	46-258877G37
	50 Feet	65 Feet
BNC RG 223	46-251710G1	46-251710G18
Cable Write-on Labels (60 labels per booklet)	46-208781P3	

DESIGNATED CABLE REPLACEMENTS

CABLE PART #	CABLE RUN
46-244710G30	Run 1227 Run 1236
46-244712G30	Run 1188
46-244713G30	Run 1181
46-244714G30	Run 1180 Run 1181 Run 1233
46-258877G30	Run 1189 Run 1190

CABLE KITS

FRU Kit Description	Part #
Receive Cables for Equipment Room, 34.5 ft (10.56m)	5119953
Receive Cables for Equipment Room, 65.6 ft (20m)	5119953-2
Receive Cables for Magnet Room, 28.2 ft (8.6m)	5119954
Receive Cables for Magnet Room, 49.2 ft (15m)	5119954-2
Transportable/Relocatable MC Rcv Cable Kit (Eqmt. Rm)	5119953-3
Body Receive Cables (both rooms), 34.5 ft (10.52m)	5119956
Body Receive Cables (both rooms), 65.6 ft (20m)	5119956-2
Term Server Cables (4 cables)	5120780

FIBER OPTIC CABLE KIT – FIXED SITE


Description	Part #
Includes clamps, cable ties, mounts, and Runs 1045 & 1179. (Fixed Site)	5159412
Includes clamps, cable ties, mounts, and Runs 1045 & 1179. (Mobile & T/R sites)	5159413


TIMES MICROWAVE CABLE PARTS

Description	Part #
Straight Male N Connector for LMR 600	2352195
7/8 EIA Right Angle Plug & Center Pin	2319986
7/8 EIA Right Angle Plug for LMR 1200 & Center Pin	2384859


VRE CABLE KITS

Description	Part #
Infiniband cables – Infiniband to MGD and ICN	5160774
Ethernet cables – ICN to Ethernet, Runs 1371 & 1345	5160775



UNIFIED SURFACE COIL PHANTOM SET


5342679
Large Cylindrical
Unified Phantom



5342680
Small Cylindrical
Unified Phantom




5342681
Cubical
Unified Phantom



5343347
TL Unified Phantom

MCR3 TOOL:

For use with systems that have Hypertronics connectors only in LPCA.

**RF TRANSMIT TERMINATION KIT:**

5131924 1.5T AND 3.0T RF TRANSMIT CABLE SERVICE KIT

GRAFIDY KIT:

2270695
CuSO4 SAMPLE VIAL

HDx TEMPERATURE AND SENSOR KIT:

5147136-2 ENVIRONMENTAL REPAIR KIT

HEAD COILS:

1.5T SPLIT HEAD COIL ASM

2341973 1.5T SPLIT HEAD COIL ASM

5191687 HEAD COIL LATCHING FRU KIT

2145391 MIRROR ASM

5140886 HDx HEAD COIL QHTR ASSEMBLY

HDx SINGLE CHANNEL SURFACE COIL ADAPTOR G3335LJ
(Part of M3335JA or M3335JB)

5140888 SURFACE COIL SINGLE CHANNEL ASSEMBLY

HDx DUAL ARRAY SURFACE COIL ADAPTOR M3335MG

5140889 SURFACE COIL DUAL ARRAY ASSEMBLY

1.5T GENERAL PURPOSE COIL CABLE W/COIL ID M1085NJ

2383623 SINGLE CHANNEL RECEIVE ONLY CABLE W/COIL ID

1.5T LINEAR EXTREMITY COIL CABLE W/COIL ID M1085PC

2383622 LINEAR EXTREMITY COIL CABLE W/COIL ID

1.5T LEGACY HEAD COIL CABLE W/COIL ID M1085PJ

2384268 LEGACY HEAD COIL CABLE W/COIL ID

COILS CONTINUED ON NEXT PAGE.

EXTREMITY COILS:

1.5T 3 CHANNEL SURFACE COIL	M3335MN
2416757	COIL ASM
2417165	CABLE ASM
2375136-3	PHANTOM POSITIONER

1.5T 8 CHANNEL CTL COIL	M3335MA
2416758	COIL ASM
2225545-9	ADAPTOR BLOCK
2377426-8	CERVICAL BRIDGE
2417163	CABLE ASM
2377426-3	PHANTOM KIT

1.5T 8 CHANNEL BODY ARRAY COIL	M3335MC
2416759	COIL ASM (POSTERIOR)
2418093	COIL ASM (ANTERIOR)
2417162	CABLE ASM

1.5T 8 CHANNEL CARDIAC COIL (HD CONNECTOR)	M3335MD
2416710	SYSTEM CABLE
2416711	POSTERIOR CABLE
2416897	POSTERIOR COIL
2416896	ANTERIOR COIL

1.5T 12 CHANNEL BODY ARRAY COIL – SINGLE	M3335LE
2417327	COIL ASM
2417328	SYSTEM CABLE
2377425-3	PHANTOM KIT

1.5T 12 CHANNEL BODY ARRAY COIL – DOUBLE	M3335LF
2417331	COIL ASM
2417329	ANTERIOR CABLE
2417330	POSTERIOR CABLE
2377425-3	PHANTOM KIT

1.5T BRAIN/SPINE COIL	M3335LG (FP) & M3335LH (U/G)
2417362	HEAD NECK UNIT (A & P)
2418074	HNU ADAPTOR BLOCK
2417360	NECK CHEST UNIT
2417361	TL UNIT
2417356	OUTPUT CABLE – HYPERTRONICS (UPGRADE)
2417357	OUTPUT CABLE – BENDIX (UPGRADE)
2417359	OUTPUT CABLE (FORWARD PRODUCTION)
2417358	TL UNIT CABLE
2418080	MIRROR

1.5T BRAIN/SPINE ACCESSORIES	
2418076	TL PHANTOM KIT
2418078	HN PHANTOM KIT
2418081	PAD KIT
5344844	PHANTOM POSITIONER, HNS-UNIFIED

COILS CONTINUED ON NEXT PAGE.

1.5T 16 CHANNEL SINGLE CONNECTOR PV COIL	M3335LM
2414066	RIGHT LEG ASM
2414067	LEFT LEG ASM
2414068	BASE PLATE
2417920	PRE-AMP BOX w/o SYSTEM CABLE
2417785	CABLE KIT
2414071	PHANTOM KIT
1.5T MNS PHOS FLEX COIL	M1040PH
2219090	MNS P31 COIL ASM.
2241998	QUICK DISCONNECT BOX W/CABLE TRAP
46-317287G7	COAXIAL CABLE, 32"
1.5T EXTREMITY SURFACE COIL	M1085AC
46-265442G1	COIL, CYLINDRICAL
46-317278P1	CABLE ASM.
1.5T BODY FLEX COIL	M1085BE
46-328215P1	DUAL-COIL CONNECTOR AND SLEEVE
2151291	ANTERIOR COIL W/CABLE
2151291-2	POSTERIOR COIL W/CABLE
1.5T GENERAL PURPOSE FLEX COIL	M1085GR
2128554	COIL ASM.
2128557	CABLE KIT
1.5T BILATERAL BREAST COIL	M1085BR
46-320785P3	MEDRAD BREAST COIL W/O CABLE
46-320765P2	CABLE TO HEAD PORT
46-320756P2	PA CABLE
1.5T OPEN BREAST COIL	M1087BT
2246360	BREAST COIL
2246360-3	CABLE ASM.
1.5T SHOULDER COIL	M1085AN
46-301058P9	MEDRAD SURFACE COIL
1.5T SHOULDER ARRAY COIL	M1087SD
2225252-7	SMALL COIL ASM.
2225252-6	LARGE COIL ASM.
1.5T SHOULDER PHASED ARRAY COIL	M1087SH
2100937-17	MEDRAD ANTERIOR COIL W/CABLE
2100937-18	MEDRAD ANTERIOR CABLE ASM.
2100937-19	MEDRAD POSTERIOR COIL W/CABLE
2100937-20	MEDRAD POSTERIOR CABLE ASM.
1.5T HIGH RES WRIST ARRAY	M1087HW
2293670-2	HIGH RES WRIST ARRAY W/BASEPLATE ASM.
2293670-3	WRIST PHANTOM
2293670-4	PHANTOM HOLDER
1.5T EXTREMITY COIL W/COIL ID	M1085NA
2383613-2	MAI QUAD KNEE COIL W/CABLE
2383613-3	MAI QUAD KNEE CABLE

COILS CONTINUED ON NEXT PAGE.

1.5T QUAD EXTREMITY KNEE/FOOT COIL	M1085E
<u>2225478-7</u>	VELCRO STRAP/LATCH
<u>2383613-2</u>	COIL ASM.
<u>2383613-3</u>	CABLE ASM.
<u>2225478-10</u>	PHANTOM
1.5T PA KNEE COIL	M1087AA
<u>46-320406P1</u>	COIL ASM.
1.5T PA TR KNEE/FOOT IMAGING OPTION	M1087E
<u>2293674-14</u>	COIL W/CABLE
<u>2293674-13</u>	CABLE
<u>2293674-3</u>	PHANTOM
<u>2293674-5</u>	PHANTOM POSITIONER
1.5T CARDIAC PHASED ARRAY IMAGING OPTION	M1087CA
<u>2174499</u>	CARDIAC ARRAY
<u>2189499</u>	COIL ASM. – ANTERIOR
<u>2189499-2</u>	COIL ASM. – POSTERIOR
<u>2215748</u>	CABLE ASM.
1.5T CARDIAC PERIPHERAL VASCULAR COIL	M1087PR
<u>2293676-2</u>	CVNV PV ARRAY COIL
<u>2293676-3</u>	CABLE ASM.
<u>2293676-4</u>	PHANTOM
<u>2192406-7</u>	LOADING PHANTOM
<u>2293676-6</u>	PHANTOM POSITIONER
<u>2285178-4</u>	PHANTOM STRAP
1.5T PERIPHERAL VASCULAR ARRAY COIL	M1087PW
<u>2192406-2</u>	CABLE ASM
<u>2192406-12</u>	COIL ASM.
<u>2192406-5</u>	PIN DIODES
<u>2192406-6</u>	SNR PHANTOM
<u>2192406-7</u>	LOADING PHANTOM
<u>2192406-9</u>	PHANTOM POSITIONER
1.5T NEUROVASCULAR ARRAY COIL	M1087N
<u>2225476-2</u>	COIL ASM.
<u>2225476-3</u>	CABLE ASM.
<u>2225476-6</u>	COIL LATCHES
<u>2225476-7</u>	PHANTOM POSITIONER
1.5T NEUROVASCULAR ARRAY COIL W/ NEW TOP	M1087NB
<u>2293668-2</u>	COIL W/CABLE
<u>2225476-3</u>	CABLE ASM.
<u>2225476-6</u>	COIL LATCHES
<u>2293668-7</u>	PHANTOM POSITIONER
<u>2293668-8</u>	NV ARRAY CABLE – HD COMPATIBLE
1.5T NEUROVASCULAR HEAD/NECK COIL	M1085HN
<u>2118138-5</u>	COIL ASM.
<u>2383614</u>	QUAD NV CABLE
1.5T ANTERIOR NECK COIL	M1085AK
<u>46-287396P10</u>	MEDIUM SURFACE COIL

COILS CONTINUED ON NEXT PAGE.

REVISION 11

DIRECTION 5166012

1.5T POSTERIOR CERVICAL SPINE COIL	M1085AM
46-287807P12	MEDRAD NECK COIL
1.5T QUADRATURE T/L SPINE COIL	M1085AP
46-307387P10	MEDRAD COIL W/POSITIONER
1.5T QUADRATURE CERVICAL SPINE COIL	M1085AR
46-320422P1	COIL W/CABLE
1.5T CTL COIL	M1087CT
2104799	CTL MULTICOIL ASM.
1.5T CTL SPINE ARRAY COIL	M1087SP
2225545-6	COIL W/CABLE & BRIDGE
2225545-2	CABLE ASM.
2225545-8	ANTERIOR NECK COIL ASM.
2225545-9	ADAPTER BLOCK
1.5T TORSO ARRAY COIL	M1087TB
2215261	COIL
2215748-2	CABLE ASM.
1.5T PHASED ARRAY PELVIC COIL	M1087PA
2144100-2	PELVIC COIL
46-307821G1	CABLE ASM.
2144806	PHASE ADJUST
1.5T TORSO PELVIS COIL	M1087TP
2298043-2	COIL ASM.
2298043-3	CABLE ASM.
2298043-4	PHANTOM
2298043-5	PHANTOM POSITIONER
1.5T HI RES BRAIN COIL	M3087JA
2317112-2	8CHANNEL HR BRAIN CL/BASEPLATE
2317112-3	8CHANNEL HR BRAIN PHANTOM POSITIONER
2317112-8	8CHANNEL HR HEAD CABLE ASSEMBLY
2317112-12	8CHANNEL HR BRAIN LATCH KIT
5146634-6	MIRROR ASM
1.5T 8 CHANNEL NV ARRAY COIL	M3087JB
2317115-2	8CHANNEL NV ARRAY
2317115-3	8CHANNEL PHANTOM KIT
2317115-10	8CHANNEL RIGHT/LEFT LATCH KIT
2317115-11	8CHANNEL SUPERIOR LATCH KIT
5344806	PHANTOM POSITIONER, NV-UNIFIED
1.5T 8 CHANNEL CTL COIL	M3087JC
2377426-2	8CHANNEL COIL W/CABLE
2377426-3	8CHANNEL PHANTOM
2377426-7	8CHANNEL CABLE
2377426-10	8 CHANNEL PREAMP
2418267	ANTERIOR BRIDGE
2415553	CTL COIL TABLE LOCATOR FRU KIT
1.5T 8 CHANNEL CARDIAC COIL	M3087JD
2304255-2	8 CHANNEL CARDIAC ARRAY COIL W/CABLE
2304255-3	8 CHANNEL CARDIAC ARRAY PHANTOM KIT
2304255-7	8 CHANNEL CARDIAC ARRAY CABLE
2377426-10	8 CHANNEL PREAMP

COILS CONTINUED ON NEXT PAGE.

1.5T 8 CHANNEL FOOT ANKLE COIL	
5308342-2	1.5T 8 CH FOOT ANKLE COIL BY INVIVO COIL
5308342-3	1.5T 8 CH FOOT ANKLE COIL CABLE by INVIVO COIL
5273235-5	1.5T 8 CH FOOT ANKLE COIL PHANTOM by INVIVO COIL

NOTES:

[illegible]

HD23 COILS:

1.5T FLEX SUITE PREMIUM-HCON	M7000SG
5430000-2	GEM Flex Coil 16-L Array, 1.5T Receive Only
5430000-3	GEM Flex Coil 16-M Array, 1.5T Receive Only
5430000-4	GEM Flex Coil 16-S Array, 1.5T Receive Only
5430000-8	GEM Flex Knee Stabilization Device for Curved & Flat Tables
5430000-9	GEM Flex General Purpose Strap and Interface Cover Kit
5430000-10	GEM Flex 16-L Coil Covering Replacement Kit
5430000-11	GEM Flex 16-M Coil Covering Replacement Kit
5430000-12	GEM Flex 16-S Coil Covering Replacement Kit
5430000-13	GEM Flex Stabilization Device for Curved & Flat Tables
5430000-14	GEM Flex System Cable Take Up Pad
5430000-17	GEM Flex Interface 16 Channel Fixed, 1.5T HD-Connector
5430000-21	GEM Flex 16 Channel System Cable, HD Connector
1.5T FLEX SUITE STANDARD-HCON	M7000SH
5430000-2	GEM Flex Coil 16-L Array, 1.5T Receive Only
5430000-3	GEM Flex Coil 16-M Array, 1.5T Receive Only
5430000-8	GEM Flex Knee Stabilization Device for Curved & Flat Tables
5430000-9	GEM Flex General Purpose Strap and Interface Cover Kit
5430000-10	GEM Flex 16-L Coil Covering Replacement Kit
5430000-13	GEM Flex Stabilization Device for Curved & Flat Tables
5430000-14	GEM Flex System Cable Take Up Pad
5430000-17	GEM Flex Interface 16 Channel Fixed, 1.5T HD-Connector
5430000-21	GEM Flex 16 Channel System Cable, HD Connector
1.5T FLEX SM w/INTER-HCON	M7000SJ
5430000-4	GEM Flex Coil 16-S Array, 1.5T Receive Only
5430000-9	GEM Flex General Purpose Strap and Interface Cover Kit
5430000-12	GEM Flex 16-S Coil Covering Replacement Kit
5430000-17	GEM Flex Interface 16 Channel Fixed, 1.5T HD-Connector
5430000-21	GEM Flex 16 Channel System Cable, HD Connector
1.5T 16 CHANNEL CARDIAC COIL	M3340CZ
5418908-2	1.5T 16 Channel Cardiac Coil – Anterior
5418908-3	1.5T 16 Channel Cardiac Coil – Posterior
5418908-4	1.5T 16 Channel Cardiac Coil Connector Block Assembly
5418908-5	Straps Set includes both Anterior and Posterior Velcro Straps
5418908-6	16 Channel Cardiac Coil Pads – Head Pad, Foot Pad, Anterior Pad, Posterior Top Pad, Posterior Bottom Pad and Straps
5418908-7	16 Channel Cardiac Coil Phantom Positioner

MNS OPTIONS

Description	Part #
1.5T 31P TR Switch	2354050-2
1.5T 13C TR Switch	2354050-3
1.5T 13C MNS Preamp	5114799
MNS QA Probe 1.5T 13C	5114403
MNS QA Probe 1.5T 31P	5114978-2
1.5T 13C Sample Vial	5114978-3
1.5T 31P Sample Vial	5114978-4

OPTIONS:

CADSTREAM 4.1 w/SURELOC **M3335KL**
5164701 CONFIRMA 4.1 SERVER PC & S/W w/MONITOR & MOUSE
5160894 CADSTREAM S/W & DOCS CD-ROM6
5154402 CADSTREAM 3.X TO 4.1 UPGRADE

CONFIRMA MAMMO CAD **M3033NF**
5117875 COMPLETE PC WITH SOFTWARE

EXCITE HD TWIN RESISTIVE SHIM POWER SUPPLY **M3000PP**
5118323 1.5T HIGHER ORDER SHIM POWER SUPPLY KIT

CLARIVIEW IMAGE FILTERING **M1090PN**
2203037 AW1.2 CLARIVIEW FLOPPY ASM.
2203038 AW2.0 CLARIVIEW FLOPPY ASM.
2203039 SECURITY KEY FLOPPY ASM.

BAR CODE READER KIT **M1090PM**
2180690-2 READER ONLY
2180690-3 Y CABLE ONLY

ReportCARD HARDWARE
5117866-26 ReportCARD COMPUTER
5160890 ReportCARD SOFTWARE COLLECTOR
2390369 ReportCARD OPTON KEY

FUNCTOOL
5160880 FUNCTOOL 4.0 APPLICATIONS CD

EXCITE SOFTWARE OPTION KEYS:

2156158-64	Release 15.x Software Collector EXCITE HDxt (M3)
2156158-65	Release 14.x Software Collector EXCITE HDx (M5)
2383494	Fiesta C
2383495	3D Delayed Enhancement
2383496	Navigator 3D
2383497	3D FatSat Fiesta
2383498	Tricks
2390037	Asset
2390040	ConnectPro
2390075	3D FRFSE
2390078	EchoPlus
2390079	Diffusion Tensor
2390084	Probe 2000
2390085	Probe 200 Upgrade
2390086	Sage 7
2390091	2D Delayed Enhancement
2390095	iDrive ProPlus
2390096	Cardiac Tagging
2390097	SmartPrep 2000
2390098	Flouro Trigger MRA
2390103	Performed Proc Step
5110228	Single Voxel Probe-P
5110229	Single Voxel Probe-S and Probe-P
5110230	Probe 2D – CSI
5110231	Probe 3D – CSI
5110232	MNS, SAGE, Probe-P
5110319	Research Proton Spectro
5120544	Release 12.0 M4 Apps DVD
5160890	ReportCard on AW
5121550	Tractography on AW
5122363	LAVA
5122364	PFO Evaluation
5122365	PROPELLOR Brain Suite
5122366	Tractography
5122367	2D FatSat FIESTA
5122368	3D Flair
5122370	RealTime FIESTA
5122371	VIBRANT-II
5122372	8 Channel HDx
5122373	16 Channel HDx
5122374	1.5T HDx ScanTools 14.0

SOFTWARE OPTIONS

BrainWave OPTIONS:

2363727	15" LCD Display with 5 foot VGA Cable for Computer connection
2394034-5	(Cedrus) Response Pad Asm. Left & right hand units with two buttons for each hand. Includes fiber optic cable (20 ft) and optical converter box
2394034-6	(Cedrus) "Lumina" response controller box
2394034-7	(Cedrus) AC Power module for Lumina Controller Box
2394034-9	FRU usage only Kit containing all cables [Cedrus-4, -5, -6, -7, 2xxxxxx-6, 46-265505P1 (qty 2)]
2394034-10	Penetration Sub-Panel Asm
2394034-16	PC compatible Keyboard (narrow)
2394034-17	PC compatible Mouse
5160892	BrainWave 2.0 Application Software CD-ROM
5173371-2	PC and Platform Software Asm (For HDx Only)
5110968	Patch CD for BrainWave Stimulus PC

REVISION HISTORY

REV/DATE REASON FOR CHANGE

- 1 Feb 28, 06 Initial release
- 2 Jun 26, 06 Updated gradient cable block kit part #. Added new HDx RF coils with new colors. Updated to ethernet switching hub on OW. Updated ethernet switch part # in 32 Channel cabinet. Added HD NV Array coil & Cardiac Tagging option key and CTL Table Locator FRU kit. Corrected other part numbers as needed.
- 3 Jan 26, 07 Added RF Shield Change Out Kit, Speaker Grill nylon screw, LPCA Wheel Housing FRU kit (MRIhc21559), & Dock pedal FRU parts (MRIhc17629). Corrected Linux 3-button mouse part number. For RFS cabinet, added Term Server cable FRU kit and optional 24 port Ethernet Switch. Added new ASC Module UPM board and RRF Module MUX2 board. Removed Infiniband Switch in 32 Channel cabinet. Removed MCR2 Tool.
- 4 Mar 13, 07 Added front cover laser, 5 gallon coolant pail, and LPCA LED cable. Added note for referring to 2138247 for parts not in this manual. Updated Ethernet cable info. Added center pin to description of Times Microwave right angle connector (CQA 13099582). Removed ECMT kit (CQA 13097585). Added Head Coil Latch Kit. Updated BrainWave s/w asm part # to 5173371.
- 5 Jun 7, 07 Added RF Shield Grounding Kit. Removed erroneous LED cable in LPCA. Updated APS & APG Module parts numbers in MGD to 2294300-16. Corrected part number for Single Coil Surface Coil Adaptor. (CQA 13098155).
- 6 Sep 17, 07 Moved Revision History to last page. Added new Netcom Hybrid Splitter and FRU's, revised existing Teledyne splitter FRU's. Updated HFA-S Asm part # in HFD-S cabinet (ECO 20428110). Added RFS Cabinet with CAM Module and CAM Module board details. Revised MGD & ASC Module part #'s. Added BRM gray insulation foam (PQR 13047306), LPCA bracket, pillar, and screws FRU kit (PQR 13121160) & updated LPCA cover part # (PQR 13083953), added Patient Table bumper and panel FRU kit (PQR 13127363), Grafidy CuSO4 vial (PQR 13098795), Brain Spine accessories, and 8-Channel CTL anterior bridge (PQR 13115149).
- 7 May 28, 08 Changed name of manual to *Signa HDx 1.5T Series FRU's* to reflect the addition of the HDxt & Vibrant products to HDx and HDi program. Added Ground Strap kit part #. Rear Pedestal: Updated Longitudinal Encoder part number (PQR 13175398) and added Encoder board part # 46-226552G2. Added Vibrant specific magnet enclosure cover part #'s. Chiller: Added E&W mobile Heat Exchanger parts (PQR13133952). Added Lytron pump coupler part # (PQR 13136000). HFD-S: Updated FHA Control Bd part #. RFS Cab: Added Sun ICN part # and replacement note. Updated SCP2 Bd part #. Added 5181495-3 MGD P.S. Tools: Added RF Transmit Termination kit. Coils: Added Mirror assembly to Brain Coil (ECR 203882). Options: Updated BrainWave software part #.
- 8 Sep 23, 08 Added HDxt parts to OW, coils, and software sections. Updated Remote Intercom Board part number. Added EMI Filter Asm to Pen Panel (ECO 2060894). Updated PHPS FRU's on Pen Panel (ECO 2040374). Added Linux PC module for 14M5 s/w. Added Sun hard drive FRU. Updated part # for Environmental Repair Kit.
- 9 Jun 15, 09 Added Flexible Coupling (PQR 13206931) and Depopulated 16-CH Switch to Rear Pedestal (ECO 2065359). Updated Filter FRU's on PEN Panel (ECO 2084029). Updated Monitor offerings on OW (ECO 2047560). Also added HP9300 DIMM (ECR 2078642), HP8400 DIMM (ECR 2038518) & DVD-RW (ECR 2059074). In RRF Chassis, added 5250028-2 RF-DIF (ECO 2080729).

CONTINUED ON NEXT PAGE.

REVISION HISTORY (Continued)

REV/DATE REASON FOR CHANGE

- 10 Nov 3, 09 Updated OW LCD monitor to 5343103 (ECO 2083909) and added new GOC AA 2395043–2 (ECO 2084000). Added 2294300–22 SCP3 Bd. to CAM and MGD modules (ECO 2064203, MRIhc42663). Added 2294300–23 APS2 bd to CAM (ECO 2064898). Added Unified Surface Coil Phantom Set and added Phantom positioner to Coils section (ECO 2081133).
- 11 Jun 28, 10 Revised TRM gradient coil to 2144093–9 (ECO 2102049). Added Rear Bridge Plate Cover (ECO 20635630. Added O2 Monitor and MRU (ECO 2076814) and Step–Down Transformer (ECO 2091712). Updated Host Monitor to 5245466. Added HP Z400 computer (ECO 2056587). Added HP8400 altermate replacement (ECO 2101300). Added MDP battery and charger FRU's (ECO 2097255). Added new 16 & 32 channel ICN's (ECO 2101939).
- 12 Jan 20, 11 Removed BrainWave option key, 5110736 (ECO 2109423). Added 2156158–86 S/W collector for 16.0 M3 (DOC0853690). Added S/W collectors for 15.x. Added notes for Z400 computer and 4170 ICN are only compatible for use with 15M4B s/w or higher (MRIhc54664).
- 13 Aug 2, 11 IChanges made to reflect HD 16.0v2 update, ICN end–of–life, and IRF2b & RRF–DIF3 hardware upgrade configurations.
- 14 Dec 5 11 IChanges made to reflect HD23.x update (SPR HCSDM00101291). Revised Blower Box FRUs (PCN 207370, ECO 2123843). Updated PEN panel FRUs (ECR 2113479). Updated HFA Assembly FRU (ECR 2115417)