



Service Request Number _____

Site Name: _____

FE Name: _____

N*250k PM Cycle

SwMR CT Planned Maintenance Checklist

Month: _____

Procedure

	250k	500k	750k	1000k
GANTRY - Table Function (Up/Down, In/Out) Use All Input Pads	☐	☐	☐	☐
GANTRY - Check Scan and Laser Windows	☐	☐	☐	☐
SYSTEM - QA Phantom Scan	☐	☐	☐	☐
SYSTEM - Monitor System Log	☐	☐	☐	☐
GANTRY - Check Number of Gantry Rotations _____	☐	☐	☐	☐
SYSTEM - Slice/mA Count _____	☐	☐	☐	☐
SYSTEM - Verify Emergency Stop Buttons (table/console) _____	☐	☐	☐	☐
GANTRY - Clean Ventilation Holes and Check Fans	☐	☐	☐	☐
GANTRY - Clean Tube Fan Filter	☐	☐	☐	☐
GANTRY - SDAS DAS/Duct Cleaning (only for assemblies with filter)	☐	☐	☐	☐
SYSTEM - Shutdown System at Console and Power off at A1 Panel	☐	☐	☐	☐
GANTRY - Handling and Removal of Slip Ring Brush Debris	☐	☐	☐	☐
GANTRY - Slip Ring Vacuuming / Brushes Clean Inspect	☐	☐	☐	☐
GANTRY - Inspect the Axial Motor Brushes		☐		☐
GANTRY - Grease Main Bearing				☐
GANTRY - Check Drive Belt for Wear				☐
TABLE - Upper Assembly PM	☐	☐	☐	☐
TABLE - Lower Assembly PM	☐	☐	☐	☐
TABLE - Check for Wear or Damage on Table Extension	☐	☐	☐	☐
PDU/GPDU/CRPDU - Inspect / Clean	☐	☐	☐	☐
PDU/GPDU/CRPDU - Inspect Fan Operation	☐	☐	☐	☐
PDU - Inspect 4 Green Lamps (L1, L2, L3, DCRGS status)	☐	☐	☐	☐
PDU - Inspect 6 Green Lamps on Transformer Assembly	☐	☐	☐	☐
GPDU - Inspect 7 Green Indicator Lamps	☐	☐	☐	☐
CRPDU - PM Checks	☐	☐	☐	☐
PDU/GPDU/CRPDU - Check Tightness of Wire Lugs	☐	☐	☐	☐
SYSTEM - General Cleaning	☐	☐	☐	☐
CONSOLE - Check for Smooth Rolling Trackball / Mouse CLEAN	☐	☐	☐	☐
SYSTEM - Power on and Boot System to Applications	☐	☐	☐	☐
HIGH VOLTAGE - Auto mA Calibration			☐	
SYSTEM - Check X-Ray On Lights	☐	☐	☐	☐
SYSTEM - Check Operation of Final Scan Abort Push Buttons	☐	☐	☐	☐
SYSTEM - Image Performance QA	☐	☐	☐	☐
SYSTEM - Update saveSTATE	☐	☐	☐	☐
FINAL				
Perform Site Restoration - DQA Phantom, Check Alignment Lights, Remove Test Scans, and Check Cabinet Doors / Covers, Update vBLOG and Customer Communication Book.	☐	☐	☐	☐
SYSTEM - Shutdown System Per Site Procedure	☐	☐	☐	☐

N
o
t
e
s