

MAGNETOM

Maintenance Protocol System

Preventive Maintenance Protocol

Customer:	
Address:	
Fax / E-Mail:	
Fax / E-Mail:	

Department:		Room:	
Material-No.:		Serial-No.:	
Contract-No.:		Expire date:	
Order-No.:		System - ID:	

Avanto

Espreo

The instructions M6-020.831.05.03.XX are required for this protocol

Evaluating the Condition of the System

The system has no deficiencies. *	☐
The system has slight deficiencies that have no effect on continued operation of the system. The deficiencies should be corrected preventively. *	☐
The system has serious deficiencies. For safety reasons, continued operation of the system is permitted only after successfully correcting the deficiencies.	☐

*) For imaging systems: The result of the image quality check shows no deviation from the reference values.

The evaluation was performed after completing all work steps.

Signature: _____

Date:

Name:

If required by country-specific regulations:

The customer or a representative has taken note of the result of the evaluation of the system condition.

Signature: _____

Date:

Name:

Remarks:

Remarks Regarding the Protocol

The chapter numbers refer to the chapters in the instruction, which is referenced on the cover page.

All pages have to have the serial number of the system and the date of maintenance in the page header.

The assignment n.a. (not applicable) indicates that the checkpoint or measured value is not used for this system.

On page 2 the completeness and the results of the maintenance work is confirmed.

Explanation of Abbreviations in the Protocol

Abbrev.	Explanation	Abbrev.	Explanation
SI	Safety Inspection	PMF	Preventive Maintenance, Operating Value Check, Function Check
SIE	Electrical Safety Inspection	Q	System Quality, Image Quality
SIM	Mechanical Safety Inspection	QIQ	Image Quality
PM	Preventive Maintenance	QSQ	System Quality Check
PMP	Periodic Preventive Maintenance	SW	Software Maintenance
PMA	Preventive Maintenance Adjustments	CSE	Customer Service Engineer

Measuring Equipment and Measuring Instruments Used

Measuring instruments and measuring devices (phantoms, MR coils, etc.) may not be entered in the table if they have already been entered in the mobile device.

Measuring equipment / -instruments	Serial No.	Date used

2 Maintenance Instructions

2.2 Cooling

PM Cooling system general checks

☐ ☐ ☐

2.3 ACC Cabinet

PMP Checking the fan assembly

☐ ☐ ☐

2.4 Gradient filters

PM Checking the gradient filter fans

☐ ☐ ☐

2.5 Vacuum pump (comfort kit option)

PMP Replacing the vacuum pump filter (every 2 years)

☐ ☐ ☐

Startup Date:

☐ ☐ . ☐ ☐ . ☐ ☐

Date of Last Replacement:

☐ DD ☐ ☐ ☐ ☐ YYYY

2.6 SEP (option)

PMP Cleaning the strainer of the primary water circuit

☐ ☐ ☐

PMP Cleaning the strainer of the secondary water circuit

☐ ☐ ☐

2.7 Magnet System

2.7.1 Cold Head

PM Cold head is monitored (cold head replacement on demand)

☐ ☐ ☐

PM Replacing the cold head (w/o cold head monitoring replacement every 24 months)

☐ ☐ ☐

Startup Date:

☐ ☐ . ☐ ☐ . ☐ ☐

Date of Last Replacement:

☐ DD ☐ ☐ ☐ ☐ YYYY

2.7.2 Adsorber

2.7.2.1 Replacing the adsorber (Sumitomo CSW 71; every two years)

PM Replacing the adsorber (Sumitomo CSW 71)

☐ ☐ ☐

Startup Date:

☐ ☐ . ☐ ☐ . ☐ ☐

Date of Last Replacement:

☐ DD ☐ ☐ ☐ ☐ YYYY

OK not
OK OK n.a.**2.7.2.2 Replacing the adsorber (F70; every three years)****PM** Replacing the F70 adsorber (every three years)☐ ☐ ☐

Startup Date:

☐ ☐ . ☐ ☐ . ☐ ☐Date of Last Replacement: ☐ DD ☐ MMM ☐ YYYY**2.8 Computer MRC (MRAWP) - MRSC (MRWP)****PM** Cleaning the computer☐ ☐ ☐**2.9 TFT monitor****2.9.2 Check****PMA** Checking the TFT monitor☐ ☐ ☐**2.10 Patient Table****2.10.1 Patient Table - Drive Unit****PM** Checking the patient table drive unit☐ ☐ ☐**2.11 Patient Trolley****2.11.1 Patient Trolley (option)****PM** Replacing the patient trolley battery (option)☐ ☐ ☐**2.12 Mobile - Maintenance****2.12.1 RD298 Shock-log****2.12.1.6 Downloading the RD298 Shock-log Data****PM** Downloading the Shock-log data☐ ☐ ☐**2.12.1.7 Shock-log RD298 Battery Replacement****PM** Replacing the Shock-log battery☐ ☐ ☐**2.12.2 298™ Shock-log****2.12.2.6 Downloading the 298™ Shock-log Data****PM** Downloading the shock-log™ Data☐ ☐ ☐**2.12.2.7 Shock-log 298™ Battery Replacement****PM** Replacing the Shock-log Battery☐ ☐ ☐**2.13 Software****2.13.1 Cleaning up Directories****2.13.1.1 Cleaning up directories****SW** Deleting the MR save log files☐ ☐ ☐

OK not
OK OK n.a.**2.13.2 Site-specific data****SW** Saving site-specific data☐ ☐ ☐**2.14 End of Maintenance Instructions**

Work steps of this section were performed by:

Signature: _____

Date:

Name:

All documents may only be used for rendering services on Siemens Healthcare Products. Any document in electronic form may be printed once. Copy and distribution of electronic documents and hardcopies is prohibited. Offenders will be liable for damages. All other rights are reserved.

healthcare.siemens.com/services

Siemens Healthineers Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Telephone: +49 9131 84-0
siemens.com/healthineers

Print No.: M6-020.832.05.03.02 | Replaces: M6-020.832.05.02.02
Doc. Gen. Date: 08.17 | Language: English
© Siemens Healthcare GmbH, 2015

siemens.com/healthineers