

Chapter 2

Table



KEY:

rev. 04/23/03

= No FRUs at this level

= No drawing available

= No structure below this level

(#)

= Mark number (from parent BOM)

Page ## +

= Starting page number of info

Figure 2-1 Table BOM Tree

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
8	2331933	1.000	1	METALLESS COMPATIBLE PHANTOM HOLDER
1	2269165	1.000	N	TABLE ASSEMBLY FOR H2 LCC {For further detail, see Section 1.0 , beginning on page 40 .}
2	2330508	1.000	N	FOOT SWITCH ASSEMBLY {For further detail, see Section 2.0 , beginning on page 100 .}
3	2115996-3	1.000	N	AXIAL HEADHOLDER ASM WITH
4	2115990-3	1.000	N	LCC HSA METALLESS SUPINE
5	2115993-2	1.000	N	LCC METALLESS CRADLE EXTE

Table 2-1 Table w/ Footswitch (2276776-2 BOM, rev 2)

$$FRU = NO$$





Figure 2-3 Table Assembly - LCC (2269165, sheet 2, rev 3)




Figure 2-4 Table Assembly - LCC (2269165, sheet 3, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
-	46-297894G1	-	1	FIELD KIT FOR ACTUATOR COVER ROLLER
7	2271244	1.000	1	RIGHT BASE COVER ASM FOR {For further detail, see Section 1.2 , beginning on page 50 .}
8	2271247	1.000	1	LEFT BASE COVER ASM FOR H {For further detail, see Section 1.3 , beginning on page 52 .}
14	2261675	1.000	1	ENCODER COVER, (REF 46-29
16	2261683	1.000	1	COVER, ACTUATOR, (REF 46-
17	2372371	1.000	1	COVER TABLE TOP LEFT {For further detail, see Section 1.4 , beginning on page 54 .}
18	2372372	1.000	1	COVER TABLE TOP RIGHT {For further detail, see Section 1.5 , beginning on page 56 .}
19	46-296594P1	1.000	1	DRIVE COVER
21A	46-296679P2	1.000	1	ZEUS TABLE, REAR END COVE
21B	2332150	ALT	1	TOP REAR COVER
23	46-264368G3	1.000	1	ENHANCED TABLE CONTROL HE
28	46-297061P1	1.000	1	FISH PAPER
29	2265851	2.000	1	BRACKET, LR SIDE PANEL (R
30	2265850	2.000	1	BRACKET, LF SIDE PANEL (R
31	2261669	2.000	1	BRKT,U-FRONT SIDE PNL, (R
32	2261668	2.000	1	BRKT, U-REAR SIDE PNL, (R
36	46-297329P2	2.000	1	1/4 TURN RECEPTACLE - FIX
38	2329950	1.000	1	CHANNEL,COVER SUPPORT
47	46-208562P20	14.000	2	SCREW
52	2265882	8.000	1	SIDE PIVOT BUSHING (REF 4
72	46-297678P1	1.000	1	HINGES
74	46-170490P50	1.000	1	SPRING-INCH EXTENSION 4 I
81	12-276	4.000	1	ID 0.2040 OD 0.5630 THK .
85	2146815-2	18.000	1	10-32 X 3/8 HEX SOC BUTTO

Table 2-2 Table Assembly - LCC (2269165 BOM, rev 9)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
103	2265845	5.000	1	PAD, LEVELING (REF 210587
105	2224296-3	1.000	1	VE-GE, 8MB, 2FL, GE MED,
108	2115992-2	1.000	1	860-941208 METALLESS CRAD {For further detail, see Section 1.6 , beginning on page 58 .}
115	46-265874P1	1.000	1	50 OHM IMPEDANCE BNC FEED
109A	2117299	1.000	1	METALLESS CRADLE COMPATIB
109B	2332149	ALT	1	FRONT END COVER
128	2271251	1.000	1	RIBBON JUMPER CABLE, ETC
129	2271253	1.000	1	RIBBON CABLE JUMPER, ETC
130	2271254	1.000	1	RIBBON CABLE JUMPER, ETC
131	2271255	1.000	1	CABLE JUMPER, ETC TO ETC
132A	2372485	1.000	1	ETC INTERFACE W/ AZ60A
132B	2229871	ALT	1	ETC IF BOARD
135	2287357	1.000	1	HANDLE, CRADLE
4	2261682	1.000	N	BRKT,ACTUAT CVR MT,(REF 4
5	2265868	1.000	N	BRACKET (REF 46-296820P1)
6	2269167-2	1.000	N	SIDE PANEL KIT FOR H2 LCC {For further detail, see Section 1.1 , beginning on page 48 .}
13	2261673	2.000	N	BKT, REAR COVER, (REF 46-
15	2265857	1.000	N	COVER (REF 46-297374P1)
22	2261708	2.000	N	COVER GUARD, RAILS,(REF 4
33	2265870	4.000	N	NUT PLATE (REF 46-297171P
35	2261688	1.000	N	BRKT, LWR COV SUP (REF 46
43	2146815-2	26.000	N	10-32 X 3/8 HEX SOC BUTTO
45	46-208562P2	1.000	N	SCREW
46	46-208562P3	4.000	N	SCREW
49	46-170497P38	14.000	N	SCREW
55	46-297155G1	4.000	N	ZEUS TABLE, 9" LONG CABLE

Table 2-2 Table Assembly - LCC (2269165 BOM, rev 9) (Continued)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
56	46-170012P61	2.000	N	0006-32 HEX .313 NUT HT
59	2261690	2.000	N	HEX STDF,9.5MM TAP,(REF 4
61	46-170498P23	6.000	N	SCREW
62	46-221523P10	6.000	N	HOLE DIA 0.437 MAT'L THK
63	46-170012P63	2.000	N	0010-32 HEX .375 NUT HT
64	46-220186P20	4.000	N	1/8" 120 DEGREE COUNTERSU
65	46-170498P41	2.000	N	SCREW
71	46-297695P1	2.000	N	COTTER PIN
73	2265878	1.000	N	PIN (REF 46-297676P1)
75	2265858	1.000	N	HANDLE ADDITION (REF 46-2
78	46-197018P1	1.000	N	DANGER HI-VOLTAGE
80	46-220181P5	10.000	N	WASHER #10 (.190) SCREW S
86	2265863	3.000	N	SHIM 0.80MM (REF 46-19635
87	2265863-2	3.000	N	SHIM 0.40MM (REF 46-19635
88	46-170106P1	AR	N	DBL COATED FILM .007 THK
91	2132175	0.000	N	TABLE DOCUMENT COLLECTOR
92	1000904P358	1.000	N	ID 0.380 OD 1.000 THK .
93	2265862	1.000	N	BOLT, SHIPPING (REF 21847
94	2265861	AR	N	SPACER PLT 0.50MM (REF 46
95	2265861-2	AR	N	SPACER PLT 0.80MM (REF 46
96	2265861-3	AR	N	SPACER PLT 1.00MM (REF 46
97	2265861-4	AR	N	SPACER PLT 1.20MM (REF 46
98	2265861-5	1.000	N	SPACER PLT 1.60MM (REF 46
100	46-208561P147	5.000	N	SCREW
101	2265840	5.000	N	ADJUSTER, LEVELING PAD (R
102	2265846	5.000	N	LOCK RING (REF 2106207) D
113	2265849	1.000	N	BRACKET, ETC BRD COVER (R
116	2265684	1.000	N	LABEL

Table 2-2 Table Assembly - LCC (2269165 BOM, rev 9) (Continued)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
117	2117390	1.000	N	LABEL, CE MARKING COMPLIA
118	2132620	1.000	N	CE MARKING LABEL CLASS A
119	2103580-12	8.000	N	SCREW COUNTERSUNK HEA 6 M
120	46-328417P12	8.000	N	SCREW HEXAGON SOCKET 6 MM
121	2269166	1.000	N	SUBASSEMBLY FOR H2 LCC CA {For further detail, see Section 1.7 , beginning on page 60 .}
122	2296332	1.000	N	RATING PLATE
123	46-233422P1	1.000	N	BNC JACK-PLUG, RIGHT ANGL
124	46-208922P5	1.000	N	SCREW
125	46-208984P17	3.000	N	BRASS STANDOFF MALE/FEMAL
126	2148658	1.000	N	ESD GROUND STRAP, BRAIDED
134	46-303268P1	1.000	N	PRODUCT LOCATOR CARD FORM
136	46-328417P20	2.000	N	SCREW HEXAGON SOCKET 6 MM
137	P9140ET	1.000	N	NAME PLATE

Table 2-2 Table Assembly - LCC (2269165 BOM, rev 9) (Continued)

This Page Intentionally Blank.

1.1 Side Panel Kit (2269167-2)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
9	2259861-2	1.000	1	SIDE PANEL,RT REAR, HELIO
10	2259862-2	1.000	1	SIDE PANEL,LEFT REAR, HEL
11	2259863-2	1.000	1	SIDE PANEL,RIGHT FRONT, H
12	2259864-2	1.000	1	SIDE PANEL,LEFT FRONT, HE

Table 2-3 Side Panel Kit (2269167-2 BOM, rev 0)

This Page Intentionally Blank.

1.2 Right Base Cover Assembly (2271244)

FRU = 1




Figure 2-5 Right Base Cover Assembly (2271244, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	2271257	1.000	N	RIGHT BASE COVER FOR H2 TABLE
5	2265869-2	1.000	N	RIGHT FILLER COVER (REF 4
6	2153420	1.000	N	SPRING-INCH TORSION
7	46-297328P1	4.000	N	1/4 TURN STUD AND TWO GROMMETS
8	2329955	1.000	N	EJECTING BLADE
9	46-220186P6	2.000	N	DIA .125 HOLE .129-.133 -.390 LG
10	12-297	2.000	N	ID 0.140 OD 0.312 THK .032 STL
13	46-170024P60	1.000	N	SCREW-INCH
14	46-170012P62	1.000	N	0008-32 HEX .343 NUT HT .229

Table 2-4 Right Base Cover Assembly (2271244 BOM, rev 3)

$$FRU = 1$$

Figure 2-6 Left Base Cover Assembly (2271247, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	2271259	1.000	N	LEFT BASE COVER FOR H2 TABLE
4	2265869	1.000	N	LEFT FILLER COVER (REF 46
6	2153420	1.000	N	SPRING-INCH TORSION
7	46-297328P1	4.000	N	1/4 TURN STUD AND TWO GROMMETS
8	2329955	1.000	N	EJECTING BLADE
9	46-220186P6	2.000	N	DIA .125 HOLE .129-.133 -.390 LG
10	12-297	2.000	N	ID 0.140 OD 0.312 THK .032 STL
13	46-170024P60	1.000	N	SCREW-INCH
14	46-170012P62	1.000	N	0008-32 HEX .343 NUT HT .229

Table 2-5 Left Base Cover Assembly (2271247 BOM, rev 3)

1.4 Cover Table Top Left (2372371)

FRU = 1




Figure 2-7 Cover Table Top Left (2372371 sheet 1, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	46-208747P1	4.000	2	CABLE TIE,ADH.BACKED
4	46-297698P3	1.000	2	27" HIGH PROFILE BEAD, RI
5	46-297329P1	2.000	2	1/4 TURN RECEPTACLE WITH
6	46-220186P20	4.000	2	1/8" 120 DEGREE COUNTERSU
7	46-297721P1	1.000	2	4-POS PLUG HOUSING (DRAWE
8	46-158905P16	2.000	2	MOUNTING SCREW FOR DRAWER
9	46-297328P2	2.000	2	1/4 TURN STUD AND TWO GRO
10	46-208758P1	4.000	2	3.62" X .094" CABLE TIE
11	2372346	1.000	1	CHANNEL TAPESWITCH
2	2372348	1.000	N	COVER, TOP LEFT
12	2363367	6.000	N	VINYL CAP 1/8 INCH
13	2363369	6.000	N	T SHAPE HEAD STEEL CLIP
14	2199652	1.000	N	P9235LZ LABEL

Table 2-6 Cover Table Top Left (2372371 BOM, rev 1)

1.5 Cover Table Top Right (2372372)

FRU = 1




Figure 2-8 Cover Table Top Right (2372372 sheet 1, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	46-208747P1	4.000	2	CABLE TIE,ADH.BACKED
4	46-297698P3	1.000	2	27" HIGH PROFILE BEAD, RI
5	46-297329P1	2.000	2	1/4 TURN RECEPTACLE WITH
6	46-220186P20	4.000	2	1/8" 120 DEGREE COUNTERSU
7	46-297721P1	1.000	2	4-POS PLUG HOUSING (DRAWE
8	46-158905P16	2.000	2	MOUNTING SCREW FOR DRAWER
9	46-297328P2	2.000	2	1/4 TURN STUD AND TWO GRO
10	46-208758P1	4.000	2	3.62" X .094" CABLE TIE
11	2372346	1.000	1	CHANNEL TAPESWITCH
2	2372351	1.000	N	COVER, TOP RIGHT
12	2363367	6.000	N	VINYL CAP 1/8 INCH
13	2363369	6.000	N	T SHAPE HEAD STEEL CLIP
14	2199652	1.000	N	P9235LZ LABEL

Table 2-7 Cover Table Top Left (2372372 BOM, rev 1)

1.6 Metalless Cradle (2115992-2)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
4	2268003	2.000	1	10-32 X .437 LONG, FLAT H
5	2222889-2	1.000	1	KEEPER FOR SELF-CATCHING
6	2112870-2	1.000	N	METALLESS CRADLE WITH VEL
7	46-303268P1	1.000	N	PRODUCT LOCATOR CARD FORM
3A	2285525-2	1.000	N	RATING PLATE
3B	46-302200P7	ALT	N	SERIAL #, MODEL#

Table 2-8 Metalless Cradle (2115992-2 BOM, rev 0)

This Page Intentionally Blank.

1.7 Table Subassembly (2269166)

FRU = NO




Figure 2-9 Table Subassembly (2269166, sheet 1, rev 0)




Figure 2-10 Table Subassembly (2269166, sheet 2, rev 0)




Figure 2-11 Table Subassembly (2269166, sheet 3, rev 0)



MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	2269168	1.000	1	CARRIAGE ASSEMBLY FOR CAT {For further detail, see Section 1.7.1 , beginning on page 66 .}
4	2269177	1.000	1	LATCH ASSEMBLY FOR CATHAY {For further detail, see Section 1.7.2 , beginning on page 68 .}
7	46-296370P1	2.000	1	ROLLER
11	2269169	1.000	1	CRADLE DRIVE ASM FOR CATH {For further detail, see Section 1.7.3 , beginning on page 70 .}
12	2269175	1.000	1	LONGITUDINAL ENCODER ASM {For further detail, see Section 1.7.4 , beginning on page 76 .}
14	46-136334P23	2.000	1	SPDT, 15A @125/250VAC BAS
19	2269173	1.000	1	ELEVATION ENCODER ASM FOR {For further detail, see Section 1.7.7 , beginning on page 86 .}
20	46-170047P7	1.000	1	125 PITCHES 19.635 LENG
32	2146815-2	17.000	1	10-32 X 3/8 HEX SOC BUTTO
34	2265877	1.000	1	CALIBRATION PIN (REF 46-2
48	46-276217P2	2.000	1	FERRITE BEADS/FILTERS
50	46-276217P3	1.000	1	FERRITE BEADS/FILTERS
5	2261680	1.000	N	PLATE, TABLE TOP, (REF 46
6	2265867	1.000	N	BLOCK GUIDE ROLLER (REF 4
8	2261723	2.000	N	BLOCK, BUMPER, (REF 46-27
9	46-221597P2	4.000	N	MINOR VBM-3003 NAT RUBBER
10	2265854	1.000	N	PLATE, SWITCH (REF 46-278
13	2265867-2	1.000	N	WIDE GROOVE ROLLER (REF 4
15	2261679	1.000	N	BRACE, REAR (REF 46-29681
16	2261698	1.000	N	BRACKET, SWITCH, (REF 46-
17	46-297059G1	1.000	N	SERVO ASM {For further detail, see Section 1.7.5 , beginning on page 80 .}
18	2269172	1.000	N	POWER SUPPLY ASM FOR CATH {For further detail, see Section 1.7.6 , beginning on page 82 .}

Table 2-9 Table Subassembly (2269166 BOM, rev 2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
22	46-296626P1	1.000	N	LARGE GEAR
25	46-170623P2	AR	N	BLACK, HEAVY-DUTY LUBRICA
26	46-170498P54	4.000	N	SCREW
27	46-170498P13	8.000	N	SCREW
29	12-276	5.000	N	ID 0.2040 OD 0.5630 THK .
30	46-220181P5	10.000	N	WASHER #10 (.190) SCREW S
31	46-328417P15	2.000	N	SCREW HEXAGON SOCKET 6 MM
37	2261701	1.000	N	BAR, CALIBRATION, (REF 46
38	46-170684P2	AR	N	LOCTITE 271, 10CC BOTTLE
39	46-170686P2	AR	N	LOCTITE 242, 10CC
40	2146815	5.000	N	10-32 X 1/2 HEX SOC BUTTO
41	2261724	1.000	N	CAL PLATE FOR H2 LCC CATH
43	46-328417P12	9.000	N	SCREW HEXAGON SOCKET 6 MM
44	12-275	4.000	N	WASHER, ID .2656, OD.75,
45	46-328432P3	2.000	N	WASHER LOCK - SPRING 6.1
46	2269174	1.000	N	FRAME ASSEMBLY FOR H2 LCC {For further detail, see Section 1.7.8 , beginning on page 90 .}
47	46-170015P28	1.000	N	SCREW

Table 2-9 Table Subassembly (2269166 BOM, rev 2) (Continued)

1.7.1 Carriage Assembly (2269168)

FRU = 1




Figure 2-13 Carriage Assembly (2269168, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	46-218193P1	2.000	1	WHEEL DIA 01.21 WIDTH 0.4
5	46-208580P6	3.000	1	CAM FOLLOWER
6	46-208580P19	1.000	1	CAM FOLLOWER
10	46-170498P2	4.000	1	SCREW
12	46-170024P4	1.000	1	SCREW-INCH
2	2262161-2	1.000	N	CRADLE CARRIAGE MACH.(REF
4	2264566	2.000	N	STUD, ECCENTRIC (REF 46-2
7	2261695	1.000	N	LATCH ACTUATOR, (REF 2133
8	2264561	1.000	N	ACTU., LONG.LIMIT SWITCH(
9	46-170012P77	3.000	N	LOCK NUT, 1/4-28, .437 X
11	2264577	1.000	N	SPACER, 6.35MM THK (REF 4
14	46-170686P2	AR	N	LOCTITE 242, 10CC

Table 2-10 Carriage Assembly (2269168 BOM, rev 1)

1.7.2 Latch Assembly (2269177)

FRU = 1




Figure 2-14 Latch Assembly (2269177, sheet 1, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
8	46-296972G1	1.000	1	SOLENOID ASM}
2	2261672	1.000	N	BRKT, TABLE HOME LATCH,(R
3	2261697	1.000	N	BLK,PVT, HOME LATCH(REF 4
4	2261721	1.000	N	BAR, (REF 46-278154P1)
5	2261725	1.000	N	BRKT, LATCH ASM, (REF 46-
6	2261715	1.000	N	ARM, (REF 46-278495P1)
7	2261681	1.000	N	BRACKET, SPRING, (REF 46-
9	46-170490P54	1.000	N	SPRING-INCH EXTENSION 1 I
10	46-170497P36	3.000	N	SCREW
11	46-208732P14	1.000	N	SCREW
12	46-170012P63	1.000	N	0010-32 HEX .375 NUT HT
13	46-170020P80	1.000	N	NOMINAL DIA 0.250 X 0.750
14	46-170020P66	1.000	N	NOMINAL DIA 0.125 X 0.250
15	46-170020P85	1.000	N	NOMINAL DIA 0.125 X 0.75
16	403B335P17	1.000	N	02 POSITION TERMINAL STRI
17	403B320P94	1.000	N	2 POSITION MARKER STRIP,
18	46-170497P38	2.000	N	SCREW
19	46-170497P7	4.000	N	SCREW
20	46-170686P2	AR	N	LOCTITE 242, 10CC

Table 2-11 Latch Assembly (2269177 BOM, rev 0)

$$FRU = 1$$

Figure 2-15 Cradle Drive Assembly (2269169 sheet 1, rev 0)



MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
19	46-221532P18	1.000	1	PITCH 1/5 WIDTH 00.250
22	46-296972G1	1.000	1	SOLENOID ASM
59	46-311873P8	1.000	2	PIN, M3 X 18
2	2265831	1.000	N	END PLATE, LH (REF 46-296
3	2265832	1.000	N	END PLATE, RH (REF 46-296
5	46-136282P30	2.000	N	BEARING BALL, SINGLE RO 2
6	46-212989P2	2.000	N	WASHER
7	46-208562P8	2.000	N	SCREW 3/8X16X1" FLAT HEAD
9	2265818	1.000	N	IDLER SHAFT (REF 46-29603
10	46-278575P1	1.000	N	SERVO MOTOR
12	403B335P10	1.000	N	06 POSITION TERMINAL STRI
13	46-170497P38	6.000	N	SCREW
14	403B320P86	1.000	N	NUMBERS 1 THU 6, LEFT TO
15	403B320P97	1.000	N	2 POSITION MARKER ST
16	46-170498P80	8.000	N	SCREW
17	46-220181P3	4.000	N	WASHER HELICAL#6 SCREW SI
18	1000904P116	5.000	N	ID 0.1495 OD 0.6250 THK .
20	46-170576P7	1.000	N	PITCH 1/5 TEETH 0L6 FAC
21	46-170576P6	1.000	N	PITCH 1/5 TEETH 013 FAC
23	2261718	1.000	N	BRACKET, SOLENOID, (REF 4
24	2265829	1.000	N	SPACER BAR (REF 46-296063
25	2265830	2.000	N	PIVOT PIN (REF 46-296064P
26	46-170686P2	AR	N	LOCTITE 242, 10CC
27	2265833	2.000	N	SUPPORT BAR (REF 46-29607
28	46-208562P7	4.000	N	SCREW
30	2265826	1.000	N	SHAFT (REF 46-296048P1)
31	2265824	1.000	N	BEARING SUPPORT (REF 46-2

Table 2-12 Cradle Drive Assembly (2269169 BOM, rev 2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
32	46-170498P79	1.000	N	SCREW
35	46-221451P1	AR	N	EXXON BEACON 325 35 LB P
36	46-170498P115	6.000	N	SCREW
37	2265865	1.000	N	KEY 2.38SQX15MMLG (REF 46
38	46-220179P8	1.000	N	EXT AXIAL MOUNT 0.375 SHA
39	46-136282P33	2.000	N	BEARING
40	46-136294P9	2.000	N	BEARING
41	2265834	1.000	N	SHAFT (REF 46-296110P1)
42	2265825	1.000	N	BEARING SUPPORT (REF 46-2
43	46-136303P10	1.000	N	BEARING
44	46-136282P62	2.000	N	BEARING BALL, SINGLE RO 1
45	46-220179P34	1.000	N	EXT AXIAL MOUNT 0.669 SHA
46	46-170498P84	6.000	N	SCREW, CAP SOCKET HEAD 8-
48	46-296368G1	1.000	N	
49	403B335P17	1.000	N	02 POSITION TERMINAL STRI
50	46-208718P2	2.000	N	5/16 ID LOOP TYPE TA MFG
51	509A859P69	1.000	N	SPRING-INCH COMPRESSION 1
55	46-297820P1	1.000	N	SHOULDER SCREW
56	2261694	1.000	N	CLUTCH LINK BAR, (REF 46-
60	2261717	1.000	N	BRACKET, STOP, (REF 46-29
61	46-221560P4	1.000	N	SCREW
63	46-220277P6	1.000	N	EXT RADIAL MOUNT 0.219 SH
64	46-170012P63	1.000	N	0010-32 HEX .375 NUT HT
66	2265865-2	1.000	N	KEY 4.76SQX50MMLG
70	46-170497P7	5.000	N	SCREW
71	46-297117G1	1.000	N	CRDR PWR CBL
72	46-297118G1	1.000	N	CRDR SIG CBL
73	1000904P514	2.000	N	I.D. 0.789 O.D. 1.100

Table 2-12 Cradle Drive Assembly (2269169 BOM, rev 2) (Continued)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
74	46-208758P3	3.000	N	7.31" X .184" SELF-LOCKIN
75	46-208758P9	1.000	N	11.08"X .140" SELF-LOCKIN
76	46-297155G1	1.000	N	ZEUS TABLE, 9" LONG CABLE
77	46-220181P5	3.000	N	WASHER #10 (.190) SCREW S
78	12-276	2.000	N	ID 0.2040 OD 0.5630 THK .
79	46-170623P2	AR	N	BLACK, HEAVY-DUTY LUBRICA
80	46-170498P78	4.000	N	SCREW
81	46-208758P10	1.000	N	8.00" X .094" SELF-LOCKIN
34A	46-222894P1	1.000	N	USM CORP HDC3C-86-2BL
34B	2361751	ALT	N	HARMONIC DRIVER
54A	46-297397P1	1.000	N	
54B	2265852-2	ALT	N	GEAR, LEVER, MACHINING (4
82A	2134525	1.000	N	ZEUS TABLE CRADLE DRIVE;
82B	2134525-2	ALT	N	ZEUS TABLE CRADLE DRIVE;
82C	2265827-2	ALT	N	IDLER ROLLER, NEW (REF 21
82D	2265827-3	ALT	N	IDLER ROLLER RECYCLD (REF
83A	2134526	1.000	N	ZEUS TABLE CRADLE DRIVE;
83B	2134526-2	ALT	N	ZEUS TABLE CRADLE DRIVE;
83C	2265820-2	ALT	N	DRIVE ROLLER,NEW (REF 213
83D	2265820-3	ALT	N	DRIVE ROLLER RECYCLD (REF

Table 2-12 Cradle Drive Assembly (2269169 BOM, rev 2) (Continued)

This Page Intentionally Blank.

1.7.4 Longitudinal Encoder Asm (2269175)

FRU = 1




Figure 2-17 Longitudinal Encoder Asm (2269175, sheet 1, rev 0)




Figure 2-18 Longitudinal Encoder Asm (2269175, sheet 2, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
17	46-297036G2	1.000	1	POT ASSEMBLY
21	46-286100P5	1.000	2	75 PITCHES, 7.365 INCH LE
23	46-296781P1	1.000	2	CABLE
32	46-208562P36	2.000	2	SCREW
2	46-296854P2	1.000	N	2.5 INCH, 2048 COUNT.REV
3	2261719	1.000	N	BRKT, LG ENCDR ASM (REF 4
4	2265255	1.000	N	CABLE DRUM (REF 46-296168
5	2265313	1.000	N	THUMB WHEEL (REF 46-29616
6	2265260	1.000	N	SHAFT, SPRING REEL (REF 4
7	46-296171P1	1.000	N	SPRING COIL
8	2265257	1.000	N	COVER, SPRING REEL (REF 4
9	2265256	3.000	N	SPACER (REF 46-296173P1)
10	2261716	1.000	N	MNTING PLATE ENCD,(REF 46
11	46-208562P13	5.000	N	SCREW
12	46-208562P24	6.000	N	SCREW
13	46-170623P2	AR	N	BLACK, HEAVY-DUTY LUBRICA
14	46-208894P1	2.000	N	SCREW
15	2264577-2	1.000	N	SPACER, 4.24MM THK (REF 4
16	2261707	1.000	N	BRKT, POTENTIOMETER(REF 4
18	2329244-2	1.000	N	SPACER
19	46-296438P1	1.000	N	SPROCKET MOD
20	46-296627P4	1.000	N	SPROCKET
22	46-297539P2	1.000	N	NO. 6 HI-COLLAR SPRING,
24	46-170497P42	1.000	N	SCREW
25	12-1723	1.000	N	ID 0.1590 OD 0.7500 THK .
26	46-296780P2	4.000	N	SYN-CLAMP
28	46-208735P38	1.000	N	SPLIT LENGTHWISE .094 DIA

Table 2-13 Longitudinal Encoder Asm (2269175 BOM, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
29	46-170686P2	AR	N	LOCTITE 242, 10CC
30	46-170528P23	1.000	N	SHAFT SIZE .250 MAT 2024-
31	46-170528P70	1.000	N	CLAMP
33	46-208758P1	1.000	N	3.62" X .094" CABLE TIE
34	2335271	1.000	N	TEFLON SPACE
35	46-297539P2	1.000	N	NO. 6 HI-COLLAR SPRING,

Table 2-13 Longitudinal Encoder Asm (2269175 BOM, rev 1) (Continued)

1.7.5 Servo Assembly (46-297059G1)

FRU = NO




Figure 2-19 Servo Assembly (46-297059, sheet 1, rev 2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	46-264370G1	1.000	2	CRADLE AMP
3	46-288170G1	1.000	2	ELEV/TILT AMP
8	46-297059T	0.000		TEST SPECIFICATION FOR TH
7	46-230523P6	2.000	N	MOLDED NYLON FLAT WA
6	46-208560P27	6.000	N	SCREW
5	46-297060P1	2.000	N	WHEEL ENCODE
4	46-296983P1	1.000	N	CIRCUIT BOARD BRACKET
1	46-297059	0.000		

Table 2-14 Servo Assembly (46-297059G1 BOM, rev 3)

1.7.6 Power Supply Assembly (2269172)

FRU = NO




Figure 2-20 Power Supply Assembly (2269172, sheet 1, rev 0)




Figure 2-21 Power Supply Assembly (2269172, sheet 2, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
4	46-229455P1	1.000	1	FAN TUBEAXIAL
5	46-297093P1	3.000	1	SPST, AC/DC TOGGLE S
6	46-296316P2	1.000	1	EMC POWER SUPPLY SINGLE O
7	46-296317P3	1.000	1	OBC-STC & TABLE EMC POWER
10	46-221914P3	1.000	1	INDICATOR LIGHT W/ BUILT-
20	46-278446P1	1.000	1	15A, 125VAC, 2 POLE, 3 WI
28	46-170053P13	1.000	1	3PDT, 12V DC RELAY (EN
2	2261684	1.000	N	BRKT, PWR SUPPLY, (REF 46
3	46-170497P6	8.000	N	SCREW
8	403B335P17	1.000	N	02 POSITION TERMINAL STRI
9	46-208734P26	5.000	N	SCREW
11	2264567	1.000	N	BAR, WIRE GROUND (REF 46-
12	46-170012P63	3.000	N	0010-32 HEX .375 NUT HT
13	46-297088G1	1.000	N	POWER ASSEMBLY HARNESS
14	512A720P16	4.000	N	ID .219 OD .562 THK .11
15	1000904P329	2.000	N	ID 0.1960 OD 0.6250 THK .
16	46-220489P2	1.000	N	FAN ACCESSORIES FINGER GU
17	46-297734G1	1.000	N	RESISTOR ASM
18	403B320P94	1.000	N	2 POSITION MARKER STRIP,
19	46-170497P38	2.000	N	SCREW
21	46-297325P1	1.000	N	FISHPAPER
22	2261687	1.000	N	BRKT, PWR SUPPLY, (REF 46
23	46-208758P1	1.000	N	3.62" X .094" CABLE TIE
24	2264568	1.000	N	CVR, SERVICE OUTLET(REF 4
25	46-170012P61	3.000	N	0006-32 HEX .313 NUT HT
26	46-208759P6	1.000	N	7.81" X .184" SELF-LOCKIN
27	46-170171P22	AR	N	TAPE OLD

Table 2-15 Power Supply Assembly (2269172 BOM, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
29	46-170497P58	2.000	N	SCREW
30	1000904P33	6.000	N	ID 0.173 OD 0.375 THK .
32	46-271110P1	1.000	N	LABEL, PROTECTIVE EARTH (
33	46-208560P53	2.000	N	SCREW
35	46-170497P64	1.000	N	SCREW
36	46-220181P6	1.000	N	1/4" SCREW SIZE (0.250),
37	46-220181P5	10.000	N	WASHER #10 (.190) SCREW S
38	46-170497P7	10.000	N	SCREW
39	46-271110P6	1.000	N	"READ ALL APPLICABLE DOCU

Table 2-15 Power Supply Assembly (2269172 BOM, rev 0) (Continued)

1.7.7 Elevation Encoder Assembly (2269173)

FRU = 1




Figure 2-22 Elevation Encoder Assembly (2269173, sheet 1, rev 1)




Figure 2-23 Elevation Encoder Assembly (2269173, sheet 2, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	46-296854P2	1.000	1	2.5 INCH, 2048 COUNT.REV
4	46-296780P2	3.000	1	SYN-CLAMP
10	46-136334P57	2.000	1	SPDT, MINIATURE BASIC
3	2265263	1.000	N	BLOCK, (REF 46-296772P1)
5	2265262	1.000	N	SHAFT, ELEV ENCODER(REF 4
6	46-296259P4	1.000	N	ADAPTER
7	46-170048P8	1.000	N	MIN-E-PITCH, 16 T, PD 0.8
8	2265261	1.000	N	CAM, ELEV SWITCH (REF 46-
9	2261699	1.000	N	BRKT, SWT ELEV ASM,(REF 4
11	46-136282P75	2.000	N	BEARING
12	46-220277P7	1.000	N	EXT RADIAL MOUNT 0.250 SH
13	2265264	1.000	N	NUT PLATE (REF 46-296945P
14	46-170497P7	2.000	N	SCREW
15	46-170498P125	3.000	N	SCREW
16	46-170528P23	1.000	N	SHAFT SIZE .250 MAT 2024-
17	2264577-3	1.000	N	SPACER, 10.16MM THK(REF 4
18	46-170020P63	2.000	N	NOMINAL DIA 0.250 X 0.625
19	46-220181P2	1.000	N	#4 SPRING LOCK WASHER, ST
20	1000904P438	2.000	N	ID 0.12, OD 0.312, THK .0
21	2265313	1.000	N	THUMB WHEEL (REF 46-29616
22	46-208894P1	2.000	N	SCREW
23	46-297539P2	1.000	N	NO. 6 HI-COLLAR SPRING,
24	1000904P519	2.000	N	WASHER .117 ID .190 OD X
25	46-220205P1	1.000	N	ID 0.3125 OD 0.75 THK

Table 2-16 Elevation Encoder Assembly (2269173 BOM, rev 1)

This Page Intentionally Blank.


$$FRU = NO$$


Figure 2-24 Frame Assembly (2269174, sheet 1, rev 0)





Figure 2-25 Frame Assembly (2269174, sheet 2, rev 0)






Figure 2-27 Frame Assembly (2269174, sheet 4, rev 0)




Figure 2-28 Frame Assembly (2269174, sheet 5, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
-	46-296561P11	-	1	MAGNET ROD
4	2271256	1.000	1	SIGNAL HARNESS FOR H2 TAB
6	2294322	1.000	1	ACTUATOR, LINEAR
8	2100671	2.000	1	SUSPA # 016-03570 C (ZEUS
18	2259859	1.000	1	CARRIAGE_RAIL, (REF 46-27
19	2259860	1.000	1	CARRIAGE RAIL, V, (REF 46
23	46-297350P1	2.000	1	SPACER
25	46-220277P15	2.000	1	EXT RADIAL MOUNT 0.312 SH
26	46-220277P11	1.000	1	EXT RADIAL MOUNT 0.625 SH
35	46-170684P1	AR	1	LOCTITE 271, 0.5CC VIAL
37	2265876	2.000	1	PIN, GAS SPRING (REF 46-2
42	2266015	1.000	1	TABLE POWER HARNESS
50	2261663	2.000	1	BRKT, TBL TOP COVER, (REF
63	46-297687P1	1.000	1	8 OHM, 4W, 3-1/2" SPE
79	46-297698P1	2.000	1	11" LOW PROFILE BEAD, RIB
93	46-208758P5	2.000	1	SELF LOCKING CABLE TIE
97	46-297805G1	2.000	1	JUMPER CONN.
98	46-297805G2	2.000	1	JUMPER CONN.
104	2106866-2	1.000	1	.75 IN WIDE INSULATED CAB
107	2113411	1.000	1	170 VOLT POWER SUPPLY
108	46-297335P2	2.000	1	6.90A, 250V 2-STAGE PWR L
120	46-208562P20	6.000	2	SCREW
3	2269179	1.000	N	ETC/CPU ASSEMBLY FOR CATH
5	2265855	1.000	N	PIN, ROD END (REF 46-2967
7	2133760-2	1.000	N	ZEUS TABLE ACTUATOR BASE
9	2255160-2	1.000	N	FRONT LEG FOR H2 LCC CATH
12	46-136342P135	4.000	N	BUSHING

Table 2-17 Frame Assembly (2269174 BOM, rev 4)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
13	2265856	1.000	N	SPACER, ROD END (REF 46-2
14	2265853	1.000	N	PIN, ACTUATOR (REF 46-278
15	2255158-2	1.000	N	REAR LEG FOR H2 LCC CATHA
16	2259855	1.000	N	CHANNEL,LEFT SUPPORT,(REF
17	2259854	1.000	N	CHANNEL,RIGHT SUPPORT,(RE
20	2261670	2.000	N	PIN, TABLE TOP, FRONT (RE
21	2261671	2.000	N	PIN, TABLE TOP, REAR, (RE
22	2261678	1.000	N	BRACE, FRONT, (REF 46-296
24	46-220179P15	4.000	N	EXT AXIAL MOUNT 0.625 SHA
28	46-328417P26	2.000	N	SCREW HEXAGON SOCKET 10 M
29	46-328417P15	4.000	N	SCREW HEXAGON SOCKET 6 MM
30	46-170498P14	32.000	N	SCREW
32	46-170012P81	1.000	N	005/8-18 HEX .940 NUT
33	46-220181P6	35.000	N	1/4" SCREW SIZE (0.250),
34	2261685	1.000	N	BRKT, R-CHANL CON, (REF 4
36	2261722	1.000	N	BLOCK, MOUNTING, (REF 46-
38	2261686	1.000	N	BRKT, L-CHANL CON, (REF 4
39	2261703	1.000	N	BRKT, FRT RT CHAN,(REF 46
40	12-128	4.000	N	WASHER,FLAT
41	46-170623P2	AR	N	BLACK, HEAVY-DUTY LUBRICA
43	2265871	6.000	N	BRACKET TBL BASE CVR (46-
44	46-297329P2	6.000	N	1/4 TURN RECEPTACLE - FIX
45	46-220186P20	12.000	N	1/8" 120 DEGREE COUNTERSU
46	2329243	1.000	N	SHAFT
47	2329245	1.000	N	ROLLER MOUNT
48	46-170020P2	2.000	N	NOMINAL DIA 0.0938X 0.500
49	12-275	4.000	N	WASHER, ID .2656, OD.75,
51	2261692	1.000	N	SPRING CLAMP, (REF 46-297

Table 2-17 Frame Assembly (2269174 BOM, rev 4) (Continued)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
52	46-297705P2	1.000	N	CONDUCTIVE TAPE
53	46-297705P14	2.000	N	METALLIC CONDUCTIVE TAPE
54	46-297705P4	1.000	N	MET-CON TAPE
55	46-170497P66	1.000	N	SCREW
61	2261689	1.000	N	CABLE TROUGH, (REF 46-297
62	2261674	1.000	N	SPEAKER GRILL, (REF 46-29
64	46-297686P2	1.000	N	SPEAKER INSULATION FOR H2
65	2265879	1.000	N	SPEAKER COVER (REF 46-297
66	46-297155G1	4.000	N	ZEUS TABLE, 9" LONG CABLE
69	2265875	3.000	N	CLAMP, GRND STRAP (REF 46
75	46-208759P4	14.000	N	13.93"X .301" SELF-LOCKIN
77	46-297648P1	1.000	N	CABLE CLAMP
80	2146815	25.000	N	10-32 X 1/2 HEX SOC BUTTO
81	1000904P395	4.000	N	ID 0.177 OD 0.312 THK .
82	46-208759P1	34.000	N	7.81" X .184" SELF-LOCKIN
83	46-208560P46	1.000	N	SCREW
84	46-170497P38	4.000	N	SCREW
85	46-170012P63	24.000	N	0010-32 HEX .375 NUT HT
86	46-170012P61	3.000	N	0006-32 HEX .313 NUT HT
87	2275854	1.000	N	GROUND - TABLE TO RACEWAY
88	2146815-2	22.000	N	10-32 X 3/8 HEX SOC BUTTO
89	46-297731P1	1.000	N	MOTOR LABEL
90	12-276	4.000	N	ID 0.2040 OD 0.5630 THK .
91	46-170497P36	5.000	N	SCREW
92	46-208758P9	2.000	N	11.08"X .140" SELF-LOCKIN
94	2261712	1.000	N	BRKT, LH TBL TOP CVR, (RE
95	2261712-2	1.000	N	BRKT, RH TBL TOP CVR,(REF
96	46-208758P1	4.000	N	3.62" X .094" CABLE TIE

Table 2-17 Frame Assembly (2269174 BOM, rev 4) (Continued)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
99	46-208562P7	12.000	N	SCREW
100	46-170012P62	4.000	N	0008-32 HEX .343 NUT HT
101	2265880	4.000	N	SHOULDER SCREW (REF 46-29
102	46-170012P64	1.000	N	01/4-20 HEX .437 NUT HT
105	2265847	1.000	N	BRACKET, T/E PWR SUPPLY (R
106	2265848	1.000	N	COVER, T/E PWR SUPPLY (REF
112	46-297705P10	1.000	N	METALLIC CONDUCTIVE 3M-SC
113	2259856-2	1.000	N	TABLE BASE,MACHINING,(REF
114	46-328417P12	38.000	N	SCREW HEXAGON SOCKET 6 MM
115	46-170015P3	4.000	N	SCREW SEMS 8-32 X .312, E
117	2261676	2.000	N	PIN, TBL LOWER, FRONT,(RE
118	2261677	2.000	N	PIN, TBL LOWER, REAR, (RE
119	46-170498P32	6.000	N	SCREW
121	46-328416P23	4.000	N	SCREW HEXAGON HEAD 10 MM
122	46-328431P3	4.000	N	WASHER PLAIN - LARGE 10.5
123	46-328432P4	4.000	N	WASHER LOCK - SPRING 10.2
124	46-170686P2	AR	N	LOCTITE 242, 10CC
125	46-170015P27	5.000	N	SCREW
126	12-297	5.000	N	ID 0.140 OD 0.312 THK .
127	46-220181P5	6.000	N	WASHER #10 (.190) SCREW S

Table 2-17 Frame Assembly (2269174 BOM, rev 4) (Continued)

This Page Intentionally Blank.

Section 2.0 — Foot Switch Assembly (2330508)

FRU = NO




Figure 2-29 Foot Switch Assembly (2330508, sheet 1, rev 2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	46-296129G1	2.000	1	
28	2329241	4.000	1	FOOT PEDAL
29	2297598	1.000	1	FOOT SWITCH COVER
6	46-170497P8	6.000	N	SCREW
7	46-271110P1	1.000	N	LABEL, PROTECTIVE EARTH (GROUND)
10	46-221828P4	2.000	N	0.063 DIA X 0.625 LG 18-8 SST
11	46-170012P63	3.000	N	0010-32 HEX .375 NUT HT .239
13	46-208560P65	11.000	N	SCREW
14	46-220181P6	11.000	N	1/4" SCREW SIZE (0.250), STL NI PLT
17	46-297971P1	2.000	N	TERMINAL
19	46-328417P12	2.000	N	SCREW HEXAGON SOCKET 6 MM 12 MM
20	46-328430P3	2.000	N	WASHER PLAIN - NORMAL 6.4 MM 12 MM
21	2329239	2.000	N	GUARD BAR
22	2329246	6.000	N	SPACER GROUND BAR
23	46-170498P11	6.000	N	SCREW
24	46-170686P3	AR	N	LOCTITE 242, 50CC
25	2329247	2.000	N	SPRING PLATE
26	46-170497P20	12.000	N	SCREW
27	46-136353P135	4.000	N	SPRING-INCH COMPRESSION 1.75 INCHES
30	2297597	1.000	N	FOOT SWITCH BASE
8	2329242	2.000	N	PEDAL ROD
9	2329240	1.000	N	GROUND FOOT SWITCH
15	1000274P289	3.000	N	

Table 2-18 Footswitch Assembly (2330508 BOM, rev 5)

End of Chapter.