

Chapter 6

Global Console - Linux (GOC2)

Note: LightSpeed 4.X CT scanners may use any one of four Operator Consoles. These include:

- **GOC1 (w/ Octane)** - See Chapter 5
- **GOC2 (w/ Linux PC)** - See *This* Chapter
- **GOC3 or GOC4** - See the *GOC3 IPL*, 2394950-100, or the *GOC4 IPL*, 5134893-100, respectively.




Figure 6-1 Global Console - Linux (GOC2) BOM Tree




Figure 6-2 Global Console Assembly (2366438, sheet 1, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2350046	1.000	1	CONSOLE ASSY, GOC, LINUX {For further detail, see Section 1.0 , beginning on page 420 .}
2	2337534	1.000	1	SIDE COVERS GOC {For further detail, see Section 2.0 , beginning on page 471 .}
3	46-302200P1	1.000	N	RATING PLATE W/ SN, MN_&

Table 6-1 Global Console Collector (2366438 BOM, rev 0)

Section 1.0 — Console Assembly (2350046)

FRU = 1




Figure 6-3 Console Assembly (2350046, sheet 1, rev 1)




Figure 6-4 Console Assembly (2350046, sheet 2, rev 0)






Figure 6-6 Console Assembly (2350046, sheet 4, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	2368697	1.000	1	LINUX DDR WORKSTATION - G
3	2315655	1.000	1	FRONT COVER, CONSOLE CHAS
5	2341924	1.000	1	RECON BOX ASSEMBLY, H16 {For further detail, see Section 1.2 , beginning on page 432 .}
7	2368838	1.000	1	SCSI TOWER, DVD-RAM, MOD {For further detail, see Section 1.5 , beginning on page 468 .}
8	2197229-5	1.000	1	ETHERNET SWITCH, 8 PORT
9	2154019	1.000	1	50 OHM BNC GROUNDED TERMI
10	2148559	1.000	1	10BASE-2 TO UTP CONVERTER
11	2245794	1.000	1	GLOBAL MODEM KIT {For further detail, see Section 1.6 , beginning on page 470 .}
27	2336906	1.000	1	TRACKBALL FOR LINUX CONSO
28	2368697-2	2.000	1	LINUX DDR WORKST'N 3BUTTO
28A	2330676-2	ALT	1	COMMON WORKST'N 3-BUTTON
1	2330564	1.000	N	CONSOLE TABLE TOP COLLECTOR {For further detail, see Section 1.1 , beginning on page 426 .}
6	2333194	1.000	N	ASSY-CONSOLE CHASSIS {For further detail, see Section 1.3 , beginning on page 454 .}
6A	2317687	ALT	N	CONSOLE CHASSIS ASSEMBLY {For further detail, see Section 1.4 , beginning on page 464 .}
12	2317488	1.000	N	OCTOPUS CABLE, LINUX COMP
13	2142225	3.000	N	MONITOR POWER CORD
14	2142226	3.000	N	POWER CORD FOR EMC AREA
15	46-328421P4	1.000	N	SCREW PAN HEAD 2.5 MM 6 M
16	46-328432P1	1.000	N	WASHER LOCK - SPRING 2.6
17	2318311	2.000	N	ETHERNET CABLE, RJ45 TO R
18	2323942	1.000	N	CABLE, HD50 TO HD50
19	2304737-2	2.000	N	CABLE, DVI VIDEO CABLE -
21	2323945	1.000	N	CABLE - MODEM TO HOST COM

Table 6-2 Console Assembly (2350046 BOM, rev 3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
24	46-136246P14	1.000 ft	N	2" WIDE, WHT, ADH BACK, V
25	46-136246P21	1.000 ft	N	2" WIDE X 25 YD WHITE VEL
26	46-208758P9	3.000	N	11.08"X .140" SELF-LOCKIN
29	2227726	1.000	N	MOUSE PAD FOR LIGHTSPEED
30	2318712	1.000	N	VME ETHERNET CABLE, LINUX
33	46-297332P1	1.000	N	50 OHM, BNC TEE ADAP
34	46-208781P2	7.000	N	WRITE-ON LABEL FOR MARKIN
35	2324697-2	1.000	N	SCSI CABLE, VHD68 TO CENT
36	2335376	1.000	N	CABLE - PDU TO CONSOLE FO
39	2132370	1.000	N	CSA LABEL
40	2117390	1.000	N	LABEL, CE MARKING COMPLIA
42	2284884	1.000	N	VXWORKS 5.4 LICENCE
43	46-269568P1	1.000	N	INFORMIX SOFTWARE-DATA BA
44	2240540	1.000	N	NETSCAPE LICENSE
45	2214164	1.000	N	5.25" 2.3GB 512 BYTES/SEC
46	2337083-2	8.000	N	SCREW, M6 X 20 SOC HEAD,
47	46-328430P3	8.000	N	WASHER PLAIN - NORMAL 6.4
48	46-328432P3	8.000	N	WASHER LOCK - SPRING 6.1
49	2337414	1.000	N	DVD RAM MEDIA, 9.4 GB
50	2132620	1.000	N	CE MARKING LABEL CLASS A
51	2103585-2	4.000	N	SCREW PAN HEAD 4 MM 10 MM
52	46-328431P1	4.000	N	WASHER PLAIN - LARGE 4.3
53	46-279769P1	1.000	N	IEC TYPE B SYMBOL
54	46-269337P6	1.000	N	LABEL OLD
55	46-208758P3	16.000	N	7.31" X .184" SELF-LOCKIN
56	46-208758P2	3.000	N	5.5" X .140" SELF-LOCKING

Table 6-2 Console Assembly (2350046 BOM, rev 3) (Continued)

1.1 Console Table Top Collector (2330564)

FRU = NO

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2315657	1.000	1	TABLE TOP ASSEMBLY, CONSOLE CHASSIS {For further detail, see Section 1.1.1 , beginning on page 428 .}
2	2330566	1.000	1	KEYBOARD TABLE TOP {For further detail, see Section 1.1.2 , beginning on page 430 .}

Table 6-3 Console Table Top Collector (2330564 BOM, rev 0)

This Page Intentionally Blank.

$$FRU = 1$$


MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2317496	1.000	N	MONITOR TABLE TOP
2	2317588	2.000	N	ADJUSTMENT BRKT, MONITOR
3	2317503	2.000	N	SUPPORT BRACKET, MONITOR
4	46-328422P501	4.000	N	SCREW COUNTERSUNK HEA 5 M
5	46-170015P23	6.000	N	SCREW
6	46-170686P3	AR	N	LOCTITE 242, 50CC

Table 6-4 Table Top Assembly, Monitor (2315657 BOM, rev 0)

$$FRU = 1$$

Figure 6-8 Table Top Assembly, Keyboard (2330566, sheet 1, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2317494	1.000	N	KEYBOARD TABLE TOP
3	2317512	2.000	N	SUPPORT BRACKET, KEYBD TA
6	2317498	2.000	N	KEYBOARD ADJUSTMENT BRACK
7	2317593	2.000	N	LATCH HINGE, LEFT HANDED
8	2317594	2.000	N	LATCH HINGE, RIGHT HANDED
9	2317497	2.000	N	ALIGNMENT BRACKET
10	2317500	2.000	N	TUBE, KEYBD ALIGNMENT DEV
11	2317595	1.000	N	SLIDE RAIL FOR DRAWER ASS
12	2317499	1.000	N	DRAWER
13	2317513	2.000	N	DRAWER BRACKET
14	2317501	1.000	N	BRACKET, DRAWER BALL LATC
15	2317502	1.000	N	LATCH, DRAWER ASSEMBLY
16	46-170015P23	12.000	N	SCREW
17	46-170015P14	8.000	N	PAN HD SCREW 8-32 X 3/8
18	46-328421P11	8.000	N	SCREW PAN HEAD 4 MM 5 MM
19	46-312346P2	8.000	N	NUT, M4 X 0.7 HEX WITH LO
21	46-328421P3	4.000	N	SCREW PAN HEAD 2.5 MM 5 M
22	46-170686P3	AR	N	LOCTITE 242, 50CC
23	2103585	8.000	N	SCREW PAN HEAD 4 MM 8 MM
24	2337083	13.000	N	SCREW, M6 X 35 SOC HEAD
25	46-328430P3	13.000	N	WASHER PLAIN - NORMAL 6.4
26	46-328432P3	13.000	N	WASHER LOCK - SPRING 6.1

Table 6-5 Table Top Assembly, Keyboard (2330566 BOM, rev 0)

$$FRU = 1$$

Figure 6-9 Recon Box Assembly (2341924, sheet 1, rev 3)




Figure 6-10 Recon Box Assembly (2341924, sheet 2, rev 0)




Figure 6-11 Recon Box Assembly (2341924, sheet 3, rev 3)



MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2333179	1.000	1	CHASSIS,VME RECON BOX {For further detail, see Section 1.2.1 , beginning on page 438 .}
1A	2307794	ALT	1	VME CHASSIS ASSEMBLY {For further detail, see Section 1.2.2 , beginning on page 448 .}
2	2300482	1.000	1	PEGASUS/NANCY ASSEMBLY
4	2225004-5	2.000	1	HARD DRIVE 36GB 15KRPM SC
5	2167014	1.000	1	RHAPSODE INTERCOM
6	2269601	1.000	1	PRESCRIBED TILT BOARD
7	2339284	1.000	1	PPC MVME2400 ASSY, HAWK3
9	2339285	1.000	1	PMC SCSI CONTROLLER CARD
10	2326523	1.000	1	16 SLICE DIP
35	2349539	1.000	1	SCSI CABLE, SCAN DATA, H1
36	2341109	1.000	1	CABLE, SCSI, VHD68 TO HD6
8	2209277	2.000	N	STANDOFF/SPACER FEMALE/FE
11	46-316540P2	12.000	N	SCREW PAN HEAD 3 MM 8 MM
12	46-221417P5	1.000	N	FEMALE SCREWLOCKS FOR HD
13	46-170686P3	AR	N	LOCTITE 242, 50CC
14	46-287187P4	2.000	N	
16	2280022	2.000	N	HANDSCREW, MOTOROLA BOARD
17	2142006	4.000	N	HEX SOC SCREW M2,5 X 16
28	2119591-2	2.000	N	25 PIN SUB-D EMI FILTER G
29	46-328421P4	10.000	N	SCREW PAN HEAD 2.5 MM 6 MM
31	46-170015P10	8.000	N	SEMS,6-32 X .25
32	46-136246P14	0.250 ft	N	2" WIDE, WHT, ADH BACK, VELCRO
33	46-136246P21	0.250 ft	N	2" WIDE X 25 YD WHITE VELCRO
34	46-208758P2	4.000	N	5.5" X .140" SELF-LOCKING CABLE TIE
37	2322199-2	1.000	N	HD POWER CABLE, VME RECON

Table 6-6 Recon Box Assembly (2341924 BOM, rev 4)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
38	2215884	2.000	N	CONNECTOR HARDWARE 9.02 M
39	46-328430P1	2.000	N	WASHER PLAIN - NORMAL 2.7

Table 6-6 Recon Box Assembly (2341924 BOM, rev 4) (Continued)

$$FRU = 1$$






MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2333180	1.000	N	VME RECON BOX
2	2324186	1.000	N	FAN ASSEMBLY, VME CHASSIS {For further detail, see Section 1.2.1.1 , beginning on page 444 .}
3	46-221513P3	4.000	N	CLAMP OR CLIP
4	46-221513P7	2.000	N	CLAMP OR CLIP
5	2260421-2	1.000	N	VME CONSOLE POWER SUPPLY {For further detail, see Section 1.2.1.2 , beginning on page 446 .}
6	2333181	1.000	N	BRACKET
7	2321164	2.000	N	FAN GUARD, 172 MM FANS
8	46-220191P5	1.000	N	PANDUIT FCM3-A-14 GRAY
9	2324209	8.000	N	SCREW SOCKET HEAD M4 X 60
10	2103585	12.000	N	SCREW PAN HEAD 4 MM 8 MM
11	2109873-11	8.000	N	SCREW PAN HEAD 4 MM X 5 M
12	2319051	1.000	N	5V POWER CABLE (+)
13	2319050	1.000	N	5V POWER CABLE, (-)
14	2318710	1.000	N	POWER CABLE, BACKPLANE TO
15	2307792	1.000	N	POWER CABLE, VME FANS
16	2337509	1.000	N	CABLE, ETHERNET
17	2322199	1.000	N	HD POWER CABLE, VME RECON
18	2315116	1.000	N	SCSI CABLE, 68 PIN RIBBON
19	2318714	1.000	N	CABLE, INTERNAL ABORT
20	2265406	1.000	N	FIBER-OPTIC COUPLING RECE
21	2142224	1.000	N	VME POWER SUPPLY TO RECEP
23	2318713	1.000	N	VME FIBER OPTIC CABLE
24	2318711	1.000	N	VME OCTOPUS CABLE
25	2319648	1.000	N	FIBER OPTIC CABLE, SC TO
26	2119591-3	1.000	N	37 PIN SUB-D EMI FILTER G
27	2119591-4	1.000	N	9 PIN SUB-D EMI FILTER GA

Table 6-7 Chassis, VME Recon Box (2333179 BOM, rev 4)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
28	46-328421P3	2.000	N	SCREW PAN HEAD 2.5 MM 5 M
29	2319650	1.000	N	FIBER OPTIC BULKHEAD RECE
30	46-221417P4	5.000	N	FEMALE SCREWLOCKS FOR HD
31	N9406JS	2.000	N	SCREW
32	46-208758P2	8.000	N	5.5" X .140" SELF-LOCKING
33	U0507AC	1.000	N	PROTECTOR
36	2141161	1.000	N	RHAPSODE BACKPLANE
37	46-208935P8	4.000	N	1/4-20, 7/16" HEX X 3/16
38	46-220184P54	2.000	N	SCREW
39	1000904P475	8.000	N	ID 0.281 OD 0.625 THK .
40	46-221449P6	4.000	N	SCREW SIZE 1/4 (0.250) PH
41	2103585-2	10.000	N	SCREW PAN HEAD 4 MM 10 MM
42	46-328431P1	10.000	N	WASHER PLAIN - LARGE 4.3
43	46-221852P2	5.000	N	SWITCH
44	2327361	1.000	N	SINGLE SIDED SILICON RUBB
45	2323985	1.000	N	BRACKET, VME CHASSIS-BACK
46	2323987	1.000	N	BRACKET, VME CHASSIS, FRO
47	2332913	2.000	N	FOAM, POLYETHLENE, SELF A
48	46-312679P2	4.000	N	INTERNALLY RIBBED GROMMET
49	12-276	4.000	N	ID 0.2040 OD 0.5630 THK .

Table 6-7 Chassis, VME Recon Box (2333179 BOM, rev 4) (Continued)

This Page Intentionally Blank.

1.2.1.1 Fan Assembly, VME Chassis (2324186)

FRU = 1




Figure 6-16 Fan Assembly, VME Chassis (2324186, sheet 1, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2307791	2.000	N	FAN,172MMX150X51MM, 120 C
2	2324197	1.000	N	DB9 CONNECTOR, SOLDER CUP
3	2324196	1.000	N	DB9 BACKSHELL KIT WITH SC

Table 6-8 Fan Assembly, VME Chassis (2324186 BOM, Rev 0)

1.2.1.2 VME Console Power Supply (2260421-2)

FRU = 1




Figure 6-17 VME Console Power Supply (2260421, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	2255986	1.000	N	POWER SUPPLY, MODULAR,600
3	2255660	2.000	N	5V POWER SUPPLY ADAPTER
4	46-220181P5	2.000	N	WASHER #10 (.190) SCREW S
5	12-276	2.000	N	ID 0.2040 OD 0.5630 THK .
6	46-208941P9	2.000	N	00010-32 .375 HEX X .124

Table 6-9 VME Console Power Supply (2260421-2 BOM, rev 1)

1.2.2 VME Chassis Assembly (2307794)

FRU = 1




Figure 6-18 VME Chassis Assembly (2307794, sheet 1, rev 2)




Figure 6-19 VME Chassis Assembly (2307794, sheet 2, rev 2)




Figure 6-20 VME Chassis Assembly (2307794, sheet 3, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	2324186	1.000	1	FAN ASSEMBLY, VME CHASSIS {For further detail, see Section 1.2.1.1 , beginning on page 444 .}
5	2260421-2	1.000	1	VME CONSOLE POWER SUPPLY {For further detail, see Section 1.2.1.2 , beginning on page 446 .}
1	2304975	1.000	N	CHASSIS, VME RECON BOX
3	46-221513P3	4.000	N	CLAMP OR CLIP
4	46-221513P7	2.000	N	CLAMP OR CLIP
6	2320003	1.000	N	POWER SUPPLY TIE DOWN BRA
7	2321164	2.000	N	FAN GUARD, 172 MM FANS
8	46-220191P5	1.000	N	PANDUIT FCM3-A-14 GRAY
9	2324209	8.000	N	SCREW SOCKET HEAD M4 X 60
10	2103585	16.000	N	SCREW PAN HEAD 4 MM 8 MM
11	2109873-11	8.000	N	SCREW PAN HEAD 4 MM X 5 M
12	2319051	1.000	N	5V POWER CABLE (+)
13	2319050	1.000	N	5V POWER CABLE, (-)
14	2318710	1.000	N	VME POWER CABLE, 110 VOLT
15	2307792	1.000	N	POWER CABLE, VME FANS
16	2337509	1.000	N	CABLE, ETHERNET
17	2322199	1.000	N	HD POWER CABLE, VME RECON
18	2315116	1.000	N	SCSI CABLE, 68 PIN RIBBON
19	2318714	1.000	N	CABLE, INTERNAL ABORT
20	2265406	1.000	N	FIBER-OPTIC COUPLING RECE
21	2142224	1.000	N	VME POWER SUPPLY TO RECEP
23	2318713	1.000	N	VME FIBER OPTIC CABLE
24	2318711	1.000	N	VME OCTOPUS CABLE
25	2319648	1.000	N	FIBER OPTIC CABLE, SC TO
26	2119591-3	1.000	N	37 PIN SUB-D EMI FILTER G
27	2119591-4	1.000	N	9 PIN SUB-D EMI FILTER GA

Table 6-10 VME Chassis Assembly (2307794 BOM, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
28	46-328421P3	2.000	N	SCREW PAN HEAD 2.5 MM 5 M
29	2319650	1.000	N	FIBER OPTIC BULKHEAD RECE
30	46-221417P4	3.000	N	FEMALE SCREWLOCKS FOR HD
31	46-328422P12	2.000	N	SCREW COUNTERSUNK HEAD 4 MM
32	46-208758P2	8.000	N	5.5" X .140" SELF-LOCKING
33	46-220178P10	AR	N	GROMMET
34	46-170071P1	AR	N	MULTI-PURPOSE ADHESIVE US
35	46-170686P3	AR	N	LOCTITE 242, 50CC
36	2141161	1.000	N	RHAPSODE BACKPLANE
37	46-208935P8	4.000	N	1/4-20, 7/16" HEX X 3/16
38	46-220184P54	2.000	N	SCREW
39	1000904P475	8.000	N	ID 0.281 OD 0.625 THK .
40	46-221449P6	4.000	N	SCREW SIZE 1/4 (0.250) PH
41	2103585-2	10.000	N	SCREW PAN HEAD 4 MM 10 MM
42	46-328431P1	10.000	N	WASHER PLAIN - LARGE 4.3
43	46-221852P2	5.000	N	SWITCH
44	2327361	1.000	N	SINGLE SIDED SILICON RUBB
45	2323985	1.000	N	BRACKET, VME CHASSIS-BACK
46	2323987	1.000	N	BRACKET, VME CHASSIS, FRO
47	2332913	2.000	N	FOAM, POLYETHLENE, SELF A

Table 6-10 VME Chassis Assembly (2307794 BOM, rev 0) (Continued)

This Page Intentionally Blank.

$$FRU = NO$$

Figure 6-21 Console Chassis Assembly (2333194 sheet 1, rev 1)




Figure 6-22 Console Chassis Assembly (2333194 sheet 2, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2335086	1.000	N	CONSOLE CHASSIS FROM SUPP
2	2174159-2	1.000	N	OUTLET BOX ASSEMBLY GLOBA {For further detail, see Section 1.3.1 , beginning on page 458 .}
3	2323937	4.000	N	STRAP, HOLD DOWN
4	2323986	2.000	N	BRACKET, HOST COMPUTER
5	2325234	1.000	N	PULL STRAP, HOST COMPUTER
6	2331988	2.000	N	STRAP, HOLDING STRAP
7	2332913	2.000	N	FOAM, POLYETHYLENE, SELF A
8	2323943	1.000	N	FAN ASSEMBLY, GLOBAL CONS {For further detail, see Section 1.3.2 , beginning on page 462 .}
9	46-208897P16	12.000	N	SCREW
10	46-328431P1	9.000	N	WASHER PLAIN - LARGE 4.3
11	2103585	9.000	N	SCREW PAN HEAD 4 MM 8 MM
12	2342018	1.000	N	LABEL, KEYBOARD/MOUSE I/O

Table 6-11 Chassis Assembly (2333194 BOM, rev 2)

This Page Intentionally Blank.

1.3.1 Outlet Box Assembly (2174159-2)

FRU = 1




Figure 6-23 Outlet Box Assembly (2174159-2, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2334322	1.000	N	OUTLET BOX COVER, GOC
4	46-278446P1	8.000	N	15A, 125VAC, 2 POLE, 3 WI
5	46-297335P3	1.000	N	16A, 250V 2-STAGE PWR LIN
6	46-243925P1	1.000	N	20A, 125V, 2-POLE, 3 WIRE
9	2335486	1.000	N	LEAD DATA FOR OUTLET BOX
24	46-328425P10	18.000	N	NUT HEXAGON 3 MM
25	46-328425P1	3.000	N	NUT HEXAGON 4 MM
26	46-312358P30	3.000	N	M3 X 0.5 X 12MM METRIC FR
28	46-328430P2	3.000	N	WASHER PLAIN - NORMAL 4.3
30	46-208695P61	1.000	N	HEYCO PART NO.2126 SB-87
31	46-271110P1	1.000	N	LABEL, PROTECTIVE EARTH (
32	2190750-2	1.000	N	POWER SWITCH BOX-GLOB CON {For further detail, see "Power Switch Box (2190750-2)," on page 460.}
33	46-313024P5	1.650	N	
34	46-208697P57	1.000	N	HEYCO PART NO.1247 SR-7W
35	2150839	1.000	N	CAUTION RECEPTACLES NOT F
37	46-311804P2	18.000	N	3.2X9X.8 ZP ST DIN
38	46-311805P10	18.000	N	HELICAL LOCK WASHER
39	46-208758P2	3.000	N	5.5" X .140" SELF-LOCKING
40	2336074	1.000	N	CIRCUIT BREAKER 15 AMPS 2
41	46-328421P501	2.000	N	SCREW PAN HEAD 3 MM 5 MM

Table 6-12 Outlet Box Assembly (2174159-2 BOM, rev 0)

Power Switch Box (2190750-2)

FRU = NO




Figure 6-24 Power Switch Box (2190750-2, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2317186	1.000	N	COVER-BACK COVER FOR SERV
2	2317185	1.000	N	BRACKET-FRONT SERVICE PAN
3	2336081	1.000	N	ROCKER SWITCH, 2PST
4	2150830	1.000	N	POWER BREAKER BEZEL
5	2200242-2	1.000	N	PWR CBL OUTLET BOX TO SER
7	46-208697P57	2.000	N	HEYCO PART NO.1247 SR-7W
8	46-328430P2	4.000	N	WASHER PLAIN - NORMAL 4.3
9	2103585	6.000	N	SCREW PAN HEAD 4 MM 8 MM
10	46-328425P1	2.000	N	NUT HEXAGON 4 MM
11	2108900	1.000	N	CONNECTOR 8 CONTACTS 150
15	46-328421P3	4.000	N	SCREW PAN HEAD 2.5 MM 5 M
16	2315117	1.000	N	PS/2 BULKHEAD CABLE
17	2318308	1.000	N	ETHERNET CABLE, CHASSIS S
18	46-170012P34	2.000	N	4-40, 1/4" HEX, 3/32 THIC
19	2368504	1.000	N	CONVERTER, DB25 TO USB
20	2315117-2	1.000	N	CABLE, PS/2 KEYBOARD CABL

Table 6-13 Power Switch Box (2190750-2 BOM, rev 2)

1.3.2 Fan Assembly, Global Console Chassis (2323943)

FRU = 1




Figure 6-25 Fan Assembly, Console Chassis (2323943, sheet 1, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2326685	1.000	N	POWER CABLE FOR 110V FAN
2	2326686	1.000	N	FAN ASSEMBLY,110V,120X120
3	2293755	3.000	N	TWIN COOLING FAN GRILL
4	46-208897P16	12.000	N	SCREW

Table 6-14 Fan Assembly, Global Console Chassis (2323943 BOM, rev 0)

1.4 Console Chassis Assembly (2317687)

FRU = NO




Figure 6-26 Console Chassis Assembly (2317687, sheet 1, rev 1)




Figure 6-27 Console Chassis Assembly (2317687, sheet 2, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	2174159-2	1.000	1	OUTLET BOX ASSEMBLY GLOBAL CONSOLE {For further detail, see Section 1.3.1 , beginning on page 458 .}
3	2318681	2.000	1	ANTI-TIP LEG, CONSOLE CHA
4	2319043	4.000	1	CASTER, 2 1/2" SWIVEL
11	2323943	1.000	1	FAN ASSEMBLY, GLOBAL CONSOLE {For further detail, see Section 1.3.2 , beginning on page 462 .}
1	2310119	1.000	N	GLOBAL CONSOLE CHASSIS
5	2317787	1.000	N	BRACKET, ETHERNET SWITCH
6	2323937	4.000	N	COMPUTER HOLD DOWN STRAP
7	2332913	2.000	N	FOAM, POLYETHYLENE, SELF A
8	2323986	2.000	N	BRACKET, HOST COMPUTER
9	2325234	1.000	N	PULL STRAP, HOST COMPUTER
10	2331988	2.000	N	STRAP, HOLDING STRAP
12	2109877-2	18.000	N	NUT HEXAGON 6 MM
13	46-208897P16	12.000	N	SCREW
14	2283156	2.000	N	RUBBER BOOT, REAR ACCESS
15	2331392	2.000	N	BRACKET FOR RUBBER BOOT
16	2331393	8.000	N	SCREW, POLYTOP SCREW, M4
17	2109867-25	4.000	N	SCREW HEXAGON SOCKET
18	2109878-4	4.000	N	WASHER PLAIN - NORMAL 10.
19	2109881-4	4.000	N	WASHER LOCK - INTERNAL 10
20	2325236	2.000	N	LEVELER, CONSOLE ANTI-TIP
21	46-328431P1	9.000	N	WASHER PLAIN - LARGE 4.3
22	2103585	9.000	N	SCREW PAN HEAD 4 MM 8 MM
24	46-328431P2	19.000	N	WASHER PLAIN - LARGE 6.4
25	2109881-3	18.000	N	WASHER LOCK - INTERNAL 6.
26	46-170686P3	AR	N	LOCTITE 242, 50CC
27	2342018	1.000	N	LABEL, KEYBOARD/MOUSE I/O

Table 6-15 Console Chassis Assembly (2317687 BOM, rev 1)

This Page Intentionally Blank.

$$FRU = 1$$

Figure 6-28 SCSI Tower Assembly (2368838, sheet 1, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2319046	1.000	N	DVDRAM PLAYER
2	2318598	1.000	N	SONY MOD SMO-F551-SD
4	2315118	1.000	N	2 BAY EXTERNAL SCSI TOWER
5	46-328421P501	8.000	N	SCREW PAN HEAD 3 MM 5 MM
8	2350404	1.000	N	TERMINATOR, CENT50, ACTIV

Table 6-16 SCSI Tower, DVD-RAM, MOD (2368838 BOM, rev 0)

1.6 Global Modem Kit (2245794)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2245773	1.000	N	MULTI-TECH GLOBAL MODEM
2	2253088	1.000	N	LABEL,CUSTOMER PROPERTY F
3	2216697	1.000	2	CABLE (RJ11 25 FEET)
4	2255874	1.000	N	CABLE, RJ11 MODULAR, 6-CON
5	2249450	1.000	2	GLOBAL MODEM ADAPTER ASSY

Table 6-17 Global Modem Kit (2245794 BOM, rev 3)

Section 2.0 — Side Covers (2337534)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2315656	1.000	N	RIGHT SIDE COVER, GLOBAL
2	2315656-2	1.000	N	LEFT SIDE COVER, GLOBAL C

Table 6-18 Side Covers (2337534 BOM, rev 0)

Section 3.0 — Miscellaneous (Associated Items)

NUMBER	FRU CODE	DESCRIPTION/NAME
2244724	1	SONY GDM-500PS 21" COLOR MONITOR
2254212	1	SONY GDM-5402T9 21" COLOR MONITOR
2349792	1	18.1 IN LCD COLOR MONITOR BEIGE
2295019	1	13W3-DB15 RGB LCD STRAIGHT CONN CBL
2249894	1	ANALOG DASM II VDB
2249895	1	DIGITAL DASM II CAM
2324697-2	1	SCSI CABLE, VHD68 TO CENT 50, LINUX
2336557	1	SEISMIC KIT FOR GLOBAL CONSOLE {For further detail, see Section 3.1 , beginning on page 473 .}
2275758	1	KEYBOARD W/SCIM COLLECTOR-ENGLISH {For further detail, see Section 5.1 , beginning on page 481 .}
2275752	1	SCIM COLLECTOR - ENGLISH {For further detail, see " SCIM Collector – English (2275752) ," on page 481 .}

Table 6-19 Miscellaneous Items (2366438 AILs, rev 2003-04-25)

3.1 Seismic Kit (2336557)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
2	2323984	1.000	N	BRACKET, SCSI TOWER MOUNT
3	2323937-2	1.000	N	STRAP, HOLD DOWN, 1" WIDE
4	46-136246P14	2.000 ft	N	2" WIDE, WHT, ADH BACK, V
5	46-136246P21	2.000 ft	N	2" WIDE X 25 YD WHITE VEL
6	46-208759P1	6.000	N	7.81" X .184" SELF-LOCKIN
7	46-328431P1	4.000	N	WASHER PLAIN - LARGE 4.3
8	46-328421P21	4.000	N	SCREW PAN HEAD 4 MM 40 MM
9	46-312346P2	4.000	N	NUT, M4 X 0.7 HEX WITH LO

Table 6-20 Seismic Kit (2336557 BOM, rev 1)

Section 4.0 — Service Kits (Associated Items)

NUMBER	FRU CODE	DESCRIPTION/NAME
2333829	1	GLOBAL CONSOLE, SCSI CABLE KIT {For further detail, see Section 4.1 , beginning on page 475 .}
2333831	1	GLOBAL CONSOLE, ETHERNET / {For further detail, see Section 4.2 , beginning on page 475 .}
2333832	1	GLOBAL CONSOLE, VIDEO/POWER CABLES {For further detail, see Section 4.3 , beginning on page 476 .}
2333833	1	GLOBAL CONSOLE, ELECTRICAL COMPONENT {For further detail, see Section 4.4 , beginning on page 476 .}
2333834	1	GLOBAL CONSOLE, BRACKETS & STRAPS {For further detail, see Section 4.5 , beginning on page 476 .}
2333835	1	GLOBAL CONSOLE, MISC PM HARDWARE {For further detail, see Section 4.6 , beginning on page 477 .}
2333836	1	GLOBAL CONSOLE, HARDWARE KIT {For further detail, see Section 4.7 , beginning on page 477 .}
2333837	1	GLOBAL CONSOLE, PWR SWITCH & CABLES {For further detail, see Section 4.8 , beginning on page 479 .}
2336912	1	CONSOLE SLIDE SERVICE KIT {For further detail, see Section 4.9 , beginning on page 479 .}

Table 6-21 Service Kits (2350046 AILs, rev. 2003-04-25)

4.1 SCSI Cable Kit (2333829)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2324697-2	1.000	N	SCSI CABLE, VHD68 TO CENT
3	2351053	1.000	N	CABLE, SCSI HD68 TO CENT5
4	2191429	1.000	N	NEXGEN CABLE FROM COMPUTE
5	2323942	1.000	N	CABLE, HD50 TO HD50
6	2317488	1.000	N	OCTOPUS CABLE, LINUX COMP
7	2315116	1.000	N	SCSI CABLE, 68 PIN RIBBON
8	2341433	1.000	N	CABLE, OCTOPUS - OCTANE 2
9	2341109	1.000	N	CABLE, SCSI, VHD68 TO HD6
10	2349539	1.000	N	SCSI CABLE, SCAN DATA, H1

Table 6-22 SCSI Cable Kit (2333829 BOM, rev 2)

4.2 Ethernet / Communication Cable Kit (2333831)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2318712	1.000	N	VME ETHERNET CABLE, LINUX
2	2318308	1.000	N	ETHERNET CABLE, CHASSIS S
3	2323945	1.000	N	CABLE - MODEM TO HOST COM
4	2318311	1.000	N	ETHERNET CABLE, RJ45 TO R
5	2318713	1.000	N	VME FIBER OPTIC CABLE
7	2319648	1.000	N	FIBER OPTIC CABLE, SC TO
8	2341436	1.000	N	CABLE, VME ETHERNET - OCT

Table 6-23 Communication Cable Kit (2333831 BOM, rev 0)

4.3 Video / Power Cables (2333832)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2142225	2.000	N	MONITOR POWER CORD
2	2142226	2.000	N	POWER CORD FOR EMC AREA
3	2304737-2	2.000	N	CABLE, DVI VIDEO CABLE -

Table 6-24 Video/Power Cables (2333832 BOM, rev 0)

4.4 Electrical Components Kit (2333833)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
4	46-297335P3	1.000	1	16A, 250V 2-STAGE PWR LIN
5	2350404	1.000	1	TERMINATOR, CENT50, ACTIV
1	46-269601P6	2.000	N	HIGH DENSITY 50-PIN SLICK
2	46-297332P1	2.000	N	50 OHM, BNC TEE ADAP
3	2154019	2.000	N	50 OHM BNC GROUNDED TERMI

Table 6-25 Electrical Components (2333833 BOM, rev 1)

4.5 Brackets and Straps (2333834)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2323984	1.000	N	BRACKET, SCSI TOWER MOUNT
2	2323937-2	1.000	N	STRAP, HOLD DOWN, 1" WIDE
3	2317787	1.000	N	BRACKET, ETHERNET SWITCH
4	2323937	4.000	N	STRAP, HOLD DOWN
5	2323986	2.000	N	BRACKET, HOST COMPUTER
6	2325234	1.000	N	PULL STRAP, HOST COMPUTER
7	2331988	2.000	N	STRAP, HOLDING STRAP

Table 6-26 Brackets & Straps (2333834 BOM, rev 0)

4.6 Miscellaneous PM Hardware (2333835)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2227726	1.000	N	MOUSE PAD FOR LIGHTSPEED
2	2325236	2.000	N	LEVELER, CONSOLE ANTI-TIP
3	2283156	2.000	N	RUBBER BOOT, REAR ACCESS
4	46-221513P3	4.000	N	CLAMP OR CLIP
5	46-221513P7	2.000	N	CLAMP OR CLIP
6	2119591-2	1.000	N	25 PIN SUB-D EMI FILTER G
7	2332913	1.000	N	FOAM, POLYETHYLENE, SELF A
8	2119670-2	4.000	N	EMC CONDUCTIVE FOAM TAPE
9	2336907	1.000	N	CONSOLE, HARD RUBBER EDGE

Table 6-27 Miscellaneous PM Hardware (2333835 BOM, rev 0)

4.7 Console Hardware Kit (2333836)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2103585	8.000	N	SCREW PAN HEAD 4 MM 8 MM
2	2104490	8.000	N	NYLON WASHER KEYSTONE #33
3	2142006	8.000	N	HEX SOC SCREW M2,5 X 16
4	2260836	8.000	N	SCREW PAN HEAD 3 MM 6 MM
5	2280022	8.000	N	HANDSCREW, MOTOROLA BOARD
6	2324209	8.000	N	SCREW SOCKET HEAD M4 X 60
7	2325236	2.000	N	LEVELER, CONSOLE ANTI-TIP
8	2331393	8.000	N	SCREW, POLYTOP SCREW, M4
9	1000904P475	8.000	N	ID 0.281 OD 0.625 THK .
10	2103585-2	8.000	N	SCREW PAN HEAD 4 MM 10 MM
11	2109867-25	8.000	N	SCREW HEXAGON SOCKET

Table 6-28 Hardware Kit (2333836 BOM, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
12	2109873-11	8.000	N	SCREW PAN HEAD 4 MM X 5 M
13	2109877-2	8.000	N	NUT HEXAGON 6 MM
14	2109878-4	4.000	N	WASHER FLAT ISO7089 M10 A
15	2109881-3	4.000	N	WASHER LOCK - INTERNAL 6.
16	2109881-4	4.000	N	WASHER LOCK - INTERNAL 10
17	2119591-2	2.000	N	25 PIN SUB-D EMI FILTER G
18	46-170015P14	8.000	N	PAN HD SCREW 8-32 X 3/8
19	46-170015P23	4.000	N	SCREW, 1/4"-20 X 1/2"
20	46-208935P8	4.000	N	1/4-20, 7/16" HEX X 3/16
21	46-220178P10	2.000	N	GROMMET
22	46-220184P54	4.000	N	SCREW
23	46-221417P4	8.000	N	FEMALE SCREWLOCKS FOR HD
24	46-221449P6	4.000	N	SCREW SIZE 1/4 (0.250) PH
25	46-221852P2	8.000	N	SWITCH
26	46-287187P4	8.000	N	SCREW, CSNK 2.5MM X 6LG
27	46-311804P2	8.000	N	3.2X9X.8 ZP ST DIN
28	46-311805P10	8.000	N	HELICAL LOCK WASHER
29	46-312346P2	8.000	N	NUT, M4 X 0.7 HEX WITH LO
30	46-328421P11	8.000	N	SCREW PAN HEAD 4 MM 5 MM
31	46-328421P3	8.000	N	SCREW PAN HEAD 2.5 MM 5 M
32	46-328421P4	8.000	N	SCREW PAN HEAD 2.5 MM 6 M
33	46-328421P501	8.000	N	SCREW PAN HEAD 3 MM 5 MM
34	46-328422P501	8.000	N	SCREW COUNTERSUNK HEA 5 M
35	46-328431P1	8.000	N	WASHER PLAIN - LARGE 4.3
36	46-328431P2	8.000	N	WASHER PLAIN - LARGE 6.4
37	46-328432P1	8.000	N	WASHER LOCK - SPRING 2.6
38	46-312679P2	4.000	N	INTERNALLY RIBBED GROMMET
39	46-208695P61	2.000	N	HEYCO PART NO.2126 SB-87
40	46-170686P3	1.000	N	LOCTITE 242, 50CC

Table 6-28 Hardware Kit (2333836 BOM, rev 0) (Continued)

4.8 Power Switch & Cables (2333837)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2150833	1.000	N	BREAKER SWITCH
2	2108900	1.000	N	CONNECTOR 8 CONTACTS 150
3	2315117	1.000	N	CABLE, PS/2 MOUSE CABLE,
4	2315117-2	1.000	N	CABLE, PS/2 KEYBOARD CABL
5	2318308	1.000	N	ETHERNET CABLE, CHASSIS S
6	2191429	1.000	N	NEXGEN CABLE FROM COMPUTE
7	2341435	1.000	N	CABLE, KEYBOARD - OCTANE
8	2341434	1.000	N	CABLE, MOUSE - OCTANE 2

Table 6-29 Power Switch & Cables (2333837 BOM, rev 0)

4.9 Slide Service Kit (2336912)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2336909	1.000	N	COMMUNICATION TRAY SLIDE
2	2336911	1.000	N	RECON / COMPUTER TRAY SLI
3	2317595	1.000	N	SLIDE RAIL FOR DRAWER ASS
4	2317501	1.000	N	BRACKET, DRAWER BALL LATC
5	2317502	1.000	N	LATCH, DRAWER ASSEMBLY

Table 6-30 Slide Service Kit (2336912 BOM, rev 0)

Section 5.0 — Keyboard w/SCIM, Collectors



Figure 6-29 Keyboard & SCIM Assembly (Rev 11/15/00)

5.1 English (2275758)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2275756	1.000	1	KEYBOARD, NORTH AMERICAN
3	2275752	1.000	1	SCIM COLLECTOR - ENGLISH {For further detail, see "SCIM Collector – English (2275752)," on page 481.}
4	2323942	1.000	1	CABLE, HD50 TO HD50
2	2275757	1.000	N	FUNCTION KEY TEMPLATE - E

Table 6-31 Keyboard w/SCIM Collector - English (2275758 BOM, rev 1)

SCIM Collector – English (2275752)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2262158	1.000	N	SCAN CONTROL INTERCOM MOD
2	2262158-11	1.000	N	OVERLAY, SCIM CT, TILT -
3	2262158-21	1.000	N	OVERLAY, SCIM CT, NO TILT

Table 6-32 SCIM Collector - English (2275752 BOM, rev 0)

5.2 French (2275758-2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2275756-2	1.000	1	KEYBOARD, EUROPEAN FRENCH
4	2323942	1.000	1	CABLE, HD50 TO HD50
2	2275757-2	1.000	N	FUNCTION KEY TEMPLATE - F
3	2275752-2	1.000	N	SCIM COLLECTOR - FRENCH {For further detail, see "SCIM Collector - French (2275752-2)," on page 482.}

Table 6-33 Keyboard w/SCIM Collector - French (2275758-2 BOM, rev 1)

SCIM Collector - French (2275752-2)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2262158	1.000	N	SCAN CONTROL INTERCOM MOD
2	2262158-12	1.000	N	OVERLAY, SCIM CT, TILT -
3	2262158-22	1.000	N	OVERLAY, SCIM CT, NO TILT

Table 6-34 SCIM Collector - French (2275752-2 BOM, rev 0)

5.3 German (2275758-3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2275756-3	1.000	1	KEYBOARD, GERMAN PC PS/2
4	2323942	1.000	1	CABLE, HD50 TO HD50
2	2275757-3	1.000	N	FUNCTION KEY TEMPLATE - G
3	2275752-3	1.000	N	SCIM COLLECTOR - GERMAN {For further detail, see "SCIM Collector - German (2275752-3)," on page 483.}

Table 6-35 Keyboard w/SCIM Collector - German (2275758-3 BOM, rev 1)

SCIM Collector - German (2275752-3)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2262158	1.000	N	SCAN CONTROL INTERCOM MOD
2	2262158-13	1.000	N	OVERLAY, SCIM CT, TILT -
3	2262158-23	1.000	N	OVERLAY, SCIM CT, NO TILT

Table 6-36 SCIM Collector - German (2275752-3 BOM, rev 0)

5.4 Swedish (2275758-4)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2275756-4	1.000	1	KEYBOARD, SWEDISH PC PS/2
3	2275752	1.000	1	SCIM COLLECTOR - ENGLISH {For further detail, see “SCIM Collector – English (2275752),” on page 481.}
4	2323942	1.000	1	CABLE, HD50 TO HD50
2	2275757	1.000	N	FUNCTION KEY TEMPLATE - E

Table 6-37 Keyboard w/SCIM Collector - Swedish (2275758-4 BOM, rev 1)

5.5 International (2275758-5)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2275756	1.000	1	KEYBOARD, NORTH AMERICAN
4	2323942	1.000	1	CABLE, HD50 TO HD50
2	2283830	1.000	N	FUNCTION KEY TEMPLATE COL
3	2275752-4	1.000	N	SCIM COLLECTOR - INTERNATIONAL {For further detail, see "SCIM Collector - International (2275752-4)," on page 485.}

Table 6-38 Keyboard w/SCIM Collector - Swedish (2275758-5 BOM, rev 1)

SCIM Collector - International (2275752-4)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2262158	1.000	N	SCAN CONTROL INTERCOM MOD
2	2262158-10	1.000	N	OVERLAY, SCIM CT, TILT -
3	2262158-14	1.000	N	OVERLAY, SCIM CT, TILT -
4	2262158-20	1.000	N	OVERLAY, SCIM CT, NO TILT
5	2262158-24	1.000	N	OVERLAY, SCIM CT, NO TILT

Table 6-39 SCIM Collector - International (2275752-4 BOM, rev 0)

End of Chapter.

© 2002-2007, General Electric Company.
GE Medical Systems, a General Electric Company, going to market as GE Healthcare.
3000 N. Grandview Boulevard
Waukesha, Wisconsin 53188
USA

www.gehealthcare.com

