
Chapter 6

Miscellaneous Items

Section 1.0 — Pro16 Ship Collector (Non-Mobile) (2373160)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2366670	1.000	1	BASE COVER COLLECTOR H4 {For further detail, see Section 1.1 , beginning on page 324.}
2	2218570-2	1.000	1	H-POWER FUSE KIT {For further detail, see Section 1.2 , beginning on page 325.}
3	2376682	1.000	1	BALANCE KIT, FIELD, HPOWER {For further detail, see Section 1.3 , beginning on page 326.}
4	2385413	1.000	N	P2030BD SHIP C HP-100/60 {For further detail, see Section 1.4 , beginning on page 327.}

Table 6-1 Pro16 Ship Collector (Non-Mobile) (2373160 BOM, rev 4)

1.1 Base Cover Collector H4 (2366670)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	46-328417P3	4.000	N	SCREW HEXAGON SOCKET 4 MM
2	46-328430P2	25.000	N	WASHER PLAIN - NORMAL 4.3
3	46-328432P2	29.000	N	WASHER LOCK - SPRING 4.1
4	46-328425P1	5.000	N	NUT HEXAGON 4 MM
5	2333997	1.000	1	COVER, FRONT BASE - H4
6	2333996	1.000	1	COVER, BASE LEFT - H4
7	2333995	1.000	1	COVER, BASE RIGHT - H4
8	2275686	1.000	1	COVER ASSY, BASE REAR - L
9	2275685	1.000	1	COVER ASSY, BASE REAR - R
10	2277652	3.000	N	LATCH, H2 SIDE COVER COMP
11	2261321	2.000	N	BRACKET - FRONT BASE COVE
12	2261324	2.000	N	BRACKET - SIDE BASE COVER

Table 6-2 Base Cover Collector H4 (2366670 BOM, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
13	2261325	1.000	N	BRACKET - REAR BASE COVER
14	2285469	2.000	N	BRACKET,BASE COVER MOUNTI
15	2288882	2.000	N	BULKHEAD - BLANK
16	46-328420P4	4.000	N	SCREW BUTTON HEAD 4 MM 16
17	2366670IDW	1.000	N	BASE COVER COLLECTOR H4
18	46-328417P4	17.000	N	SCREW HEXAGON SOCKET 4 MM

Table 6-2 Base Cover Collector H4 (2366670 BOM, rev 0) (Continued)

1.2 H-Power Fuse Kit (2218570-2)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2351492	2.000	1	200A FUSE
2	2364059	2.000	1	GLASS FUSE
3	46-170021P50	2.000	1	FUSE 12 AMPS 250 VOLTS
4	46-170021P15	2.000	1	FUSE 8 AMPS 250 VOLTS
5	2336517-2	2.000	1	FUSE 25 AMPS 700 VOLTS 2.
6	46-170021P52	3.000	1	FUSE 3 AMPS 250 VOLTS
7	46-170021P10	2.000	1	006.000A 0250V 3AG
8	46-170021P101	1.000	1	FUSE 20 AMPS 700 VOLTS
9	46-170021P11	2.000	1	7.0A, 250V SLO-BLO FUSE.
10	46-170021P31	2.000	1	1/2A, 250V SLO-BLO FUSE.
11	2106993-5	2.000	1	FUSE 20 AMPS 500 VOLTS .3
12	46-327160P1	2.000	1	12.0A, 125VAC, DUAL ELEME
13	2379651	1.000	1	FUSE - 700 V, 200 A, FAST

Table 6-3 H-Power Fuse Kit (2218570-2 BOM, rev 1)

1.3 Balance Kit, Field, H-Power (2376682)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2316653	6.000	N	WEIGHT,107 DEG,6MM
2	2316653-2	6.000	N	WEIGHT,107 DEG,1.5MM
3	2317863	2.000	N	SHIM,180 DEG,3.2MM
4	2317865-2	2.000	N	SPACER,180 DEG,25.4MM
5	2317994	20.000	N	WEIGHT,TRIM
6	2319142	2.000	N	SPACER,107 DEG,25.4MM
7	2319145	2.000	N	SHIM,107 DEG,3.3MM
8	2342895-2	6.000	N	WEIGHT,180 DEG,6MM
9	2342895-3	6.000	N	WEIGHT,180 DEG,1.5MM
10	2370053	1.000	N	WEIGHT,TRIAL,NEGATIVE Z
11	2370227	1.000	N	WEIGHT,TRIAL,POSITIVE Z,
12	46-328425P4	4.000	N	NUT HEXAGON 16 MM
13	46-328417P70	2.000	N	SCREW HEXAGON SOCKET 12 M
14	2376683	1.000	N	SOCKET, 15/16" DEEP IMPAC
15	2377399	1.000	N	LEFT SIDE DAS COVER TOOL
16	2377400	1.000	N	RIGHT SIDE DAS COVER TOOL
17	2377401	1.000	N	CENTER DAS COVER TOOL
18	46-328432P5	4.000	N	WASHER LOCK - SPRING 16.2

Table 6-4 Balance Kit, Field, H-Power (2376682 BOM, rev 0)

1.4 P2030BD Ship C HP-100/60 Common (2385413)

FRU = N

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	46-229452G1	1.000	N	CT-9800 (ACCESSORIES)
2	46-297482G1	1.000	N	IV BOARD ACC FOR ZEUS TAB
4	46-297592P1	1.000	N	BAG HOLDER
6	2105873	10.000	N	.531 I.D., 1.625 O.D., .2
7	2106573	10.000	N	HILTI ANCHOR STUD 1/2 DIA
8	2115984	1.000	N	PAD FOR HSA METALLESS FOO
9	2115994	1.000	N	METALLESS CRADLE PAD & ST
12	46-222436G1	1.000	N	CT-9800 - CHIN STRAP
13	46-297960P1	1.000	N	FOAM INSERT
14	2319891-3	1.000	N	INSTALLATION SUPPORT KIT
18	46-222436G1	3.000	N	CT-9800 - CHIN STRAP
19	46-237412P1	1.000	N	STRAP
20	46-222443P1	1.000	N	OPEN CELL POLYESTER FOAM
21	46-222444P1	1.000	N	CLOSED CELL LINEAR ETHYLE
22	46-222445P1	1.000	N	CLOSED CELL LINEAR ETHYLE
23	46-222921P1	6.000	N	OBSOLETE FOR MR SITES PER
24	46-199388P1	1.000	N	HEAD HOLDER POLYURETHANE
25	2282011	2.000	1	NON-TILT DOLLY REAR COVER
26	5101127	1.000	1	DOLLY PIVOTING , RIGHT HA
26A	2289712	ALT	1	DOLLY, PIVOTING - RH
29	2363851	1.000	N	TEMPLATE, INSTALLATION HP
30	5101127-2	1.000	1	DOLLY, PIVOTING, LEFT HAN
30A	2289712-2	ALT	1	DOLLY, PIVOTING -_LH
31	2282509	1.000	1	HOIST-2-PC, ARM, HELIOS 2
32	2282510	1.000	1	HOIST-2-PC, LONG TUBE, HE
33	2277548	4.000	1	PLUG, SIDE COVER H2

Table 6-5 P2030BD Ship C HP-100/60 Common (2385413 BOM, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
34	2217173-2	1.000	N	STOOL WITH FOOT REST AND
35	2299113	1.000	N	CT SERVICE CABINET
37	2309994	1.000	N	IV POLE ASSY
38	2329064	1.000	N	TABLE ACC ASSY
39	2355100	1.000	N	TRAY, H2 TABLE
40	2364629	1.000	1	LABEL INSTALLATION KIT
41	2281840-5	1.000	1	H-POWER SHORT CABLE COLLECTOR {For further detail, see Section 3.1 , beginning on page 334 .}
-	2281840-4	-	1	H-Power LONG CABLE COLLECTOR {For further detail, see Section 3.2 , beginning on page 334 .}
42	2336557	1.000	1	SEISMIC KIT FOR GLOBAL CO

Table 6-5 P2030BD Ship C HP-100/60 Common (2385413 BOM, rev 1) (Continued)

Section 2.0 — LightSpeed RT Ship Collector (2394003)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	46-229452G1	1.000	N	CT-9800 (ACCESSORIES)
2	46-297482G1	1.000	N	IV BOARD ACC FOR ZEUS TAB
3	2381848	1.000	1	COLLECTOR-BASE COVER-SG {For further detail, see Section 2.1 , beginning on page 330 .}
4	46-297592P1	1.000	N	BAG HOLDER
6	2105873	10.000	N	.531 I.D., 1.625 O.D., .2
7	2106573	10.000	N	HILTI ANCHOR STUD 1/2 DIA
8	2115984	1.000	N	PAD FOR HSA METALLESS FOO
9	2115994	1.000	N	METALLESS CRADLE PAD & ST
10	2218570-2	1.000	1	H-POWER FUSE KIT {For further detail, see Section 1.2 , beginning on page 325 .}
13	46-297960P1	1.000	N	FOAM INSERT
14	2319891-4	1.000	N	INSTALLATION SUPPORT KIT, {For further detail, see Section 2.2 , beginning on page 331 .}
18	46-222436G1	4.000	N	CT-9800 - CHIN STRAP
19	46-237412P1	1.000	N	STRAP
20	46-222443P1	1.000	2	OPEN CELL POLYESTER FOAM
21	46-222444P1	1.000	N	CLOSED CELL LINEAR ETHYLE
22	46-222445P1	1.000	N	CLOSED CELL LINEAR ETHYLE
23	46-222921P1	6.000	2	OBSOLETE FOR MR SITES PER
24	46-199388P1	1.000	2	HEAD HOLDER POLYURETHANE
25	2282011	2.000	1	NON-TILT DOLLY REAR COVER
26	5101127	1.000	1	DOLLY PIVOTING , RIGHT HA
26A	2289712	ALT	1	DOLLY, PIVOTING - RH
29	2389013	1.000	N	FLOOR TEMPLATE-SG
30	5101127-2	1.000	1	DOLLY, PIVOTING, LEFT HAN
30A	2289712-2	ALT	1	DOLLY, PIVOTING -_LH

Table 6-6 LightSpeed RT Ship Collector (2394003 BOM, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
31	2282509	1.000	1	HOIST-2-PC, ARM, HELIOS 2
32	2282510	1.000	1	HOIST-2-PC, LONG TUBE, HE
33	2277548	4.000	1	PLUG, SIDE COVER H2
34	2217173-2	1.000	N	STOOL WITH FOOT REST AND
35	2299113	1.000	N	CT SERVICE CABINET
36	2376682	1.000	1	BALANCE KIT, FIELD, HPOWE
37	2309994	1.000	N	IV POLE ASSY
38	2329064	1.000	N	TABLE ACC ASSY
39	2355100	1.000	N	TRAY, H2 TABLE
40	2364629	1.000	1	LABEL INSTALLATION KIT
41	2281840-5	1.000	1	H-POWER SHORT CABLE COLLECTOR {For further detail, see Section 3.1 , beginning on page 334 .}
-	2281840-4	-	1	H-POWER LONG CABLE COLLECTOR {For further detail, see Section 3.2 , beginning on page 334 .}
42	2336557	1.000	1	SEISMIC KIT FOR GLOBAL CO
46	46-237605P1	1.000	N	WRIST STRAP
47	2393793	2.000	N	BAR, FOOT SWITCH-SG
48	2386850	2.000	1	TRIAL WEIGHT, SG GANTRY,

Table 6-6 LightSpeed RT Ship Collector (2394003 BOM, rev 0) (Continued)

2.1 Collector-Base Cover-SG (2381848)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2381847	1.000	N	COVER-BASE LEFT-SG
2	2383481	1.000	1	COVER ASM,BASE,FRONT LEFT
3	2383483	1.000	1	COVER ASM,BASE,REAR LEFT-
4	2381839	1.000	1	COVER-BASE FRONT-SG

Table 6-7 Collector-Base Cover-SG (2381848 BOM, rev 1)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
5	2383459	1.000	1	COVER ASSY, REAR BASE-SG
6	2383482	1.000	1	COVER ASM,BASE,FRONT RIGH
7	2383484	1.000	1	COVER ASM,BASE,REAR RIGHT
8	2381846	1.000	1	COVER-BASE RIGHT-SG
9	2382421	2.000	N	BRACKET,REAR BASE COVER L
10	2382425	2.000	N	BRACKET,FRONT BASE COVER
11	2382426	1.000	N	BRACKET,FRONT LEFT BASE C
12	2382428	1.000	N	BRACKET,FRONT RIGHT BASE
13	2382439	2.000	N	BRACKET,REAR BASE COVER L
14	2382441	2.000	N	BRACKET,GANTRY OPTIONS MT
15	2288882	2.000	N	BULKHEAD - BLANK
16	46-328417P3	42.000	N	SCREW HEXAGON SOCKET 4 MM
17	46-328432P2	46.000	N	WASHER LOCK - SPRING 4.1
18	46-328430P2	24.000	N	WASHER PLAIN - NORMAL 4.3
19	46-328425P1	4.000	N	NUT HEXAGON 4 MM
20	2381848IDW	1.000		COLLECTOR-BASE COVER-SG
21	2394127	2.000	N	BRACKET,BASE COVER MTG,UP
22	2393783	2.000	N	BRKT,BASE CVR LATCH,REAR

Table 6-7 Collector-Base Cover-SG (2381848 BOM, rev 1) (Continued)

2.2 Installation Support Kit, Stargate (2319891-4)

FRU = N

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
3	46-316372P1	1.000	2	FORMULA 409 CLEANER.
5	46-208592P1	1.000	N	WINDEX GLASS CLEANER.
7	2253865	1.000	N	DOW BATH CLEANER PUMP SPR

Table 6-8 Installation Support Kit, Stargate (2319891-4 BOM, rev 0)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
8	46-170923P1	1.000	N	2-PART, EPOXY PATCH
10	46-221561P3	1.000	N	CLAMP OR CLIP
11	46-208758P10	25.000	N	8.00" X .094" SELF-LOCKIN
12	46-208758P3	50.000	N	7.31" X .184" SELF-LOCKIN
13	46-208758P5	50.000	N	SELF LOCKING CABLE TIE
14	46-208759P5	15.000	N	6.00" X .140" SELF-LOCKIN
15	46-208759P6	20.000	N	7.81" X .184" SELF-LOCKIN
16	46-208759P4	20.000	N	13.93"X .301" SELF-LOCKIN
17	46-208747P1	25.000	N	CABLE TIE,ADH.BACKED
19	46-230118P1	2.000	N	PAPER WYP
20	2253867	1.000	N	SCOTCH BRIGHT SCOUR PAD M
21	46-170127P5	1.000	N	1.0" MASKING TYPE TAPE, 6
23	2207303-7	25.000	N	SAFE SKIN NITRILE GLOVE,
24	46-170686P3	1.000	N	LOCTITE 242, 50CC
26	46-183039P1	1.000	N	TEXWIPE 58019-506-91/9%
27	46-183000P164	1.000	N	802 ALCOHOL 92%
28	2227089-2	1.000	1	DAS/DETECTOR SERVICE KIT
29	46-208781P3	1.000	N	WRITE-ON LABEL FOR MARKIN
36	46-194427P314	1.000	N	#3 PADLOCK AND 2 KEYS,9/3
37	46-194427P313	1.000	N	1 INCH MULTIPLE LOCKOUT D
38	46-194427P322	1.000	N	LOCKOUT TAG,BLACK & WHITE
39	2214163	1.000	1	CONNECTOR, ANALOG LOOP BA
40	2214163-2	1.000	1	CONNECTOR, DIGITAL LOOP B
41	46-303460P2	1.000	N	PAINT GE GRAY #1
42	46-303460P3	1.000	N	PAINT GE STEEL BLUE
43	46-303460P4	1.000	N	PAINT GE DARK STEEL BLUE
44	46-220178P10	2.000	N	GROMMET
45	2284023-3	1.000	1	SERVICE TOOL COLLECTOR, H

Table 6-8 Installation Support Kit, Stargate (2319891-4 BOM, rev 0) (Continued)

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
46	46-136246P14	6.000	N	2" WIDE, WHT, ADH BACK, V
47	2278989	6.000	N	LANYARD, NYLON U-LOCK, TY
48	46-136246P15	6.000	N	3/4" WIDE, WHT ADH BACK,
49	46-170123P1	1.000	N	WIDTH 0.625 THK .007 C
50	46-136246P21	6.000	N	2" WIDE X 25 YD WHITE VEL
52	46-136246P17	6.000	N	3/4" WIDE, WHT ADH BACK,
53	46-221561P2	5.000	N	VINYL TRIM W/ SEGMENTED M
54	46-208753P3	6.000	N	ID 0.375 PAD WIDTH 0.75
55	46-170123P2	3.000	N	3/4" WIDTH, .007 BLACK VI

Table 6-8 Installation Support Kit, Stargate (2319891-4 BOM, rev 0) (Continued)

Section 3.0 — Cable Collectors

3.1 Short Cable Collector (2281840-5)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2371450-2	1.000	1	RUN#55A-PDU TO RACEWAY GND
2	2371450-4	1.000	1	RUN#56A-CONSOLE TO RACEWAY
3	2373436-3	1.000	1	RUN #102A - GANTRY TO CON
4	2117848-7	1.000	1	20000 +/- 150 FIBER OPTIC
5	2371133-2	1.000	1	CONSOLE CABLES - SHORT KI
6	2371134-2	1.000	1	GANTRY CABLES - SHORT KIT

Table 6-9 Short Cable Collector (2281840-5 BOM, rev 0)

3.2 Long Cable Collector (2281840-4)

FRU = 1

MARK	NUMBER	QTY	FRU CODE	DESCRIPTION/NAME
1	2371450	1.000	1	RUN#55-PDU TO RACEWAY GND
2	2371450-3	1.000	1	RUN#56-CONSOLE TO RACEWAY
3	2373436-2	1.000	1	RUN #102 - GANTRY TO CONS
4	2117848-2	1.000	1	24384 +/- 150 FIBER OPTIC
5	2371133	1.000	1	CONSOLE CABLES - LONG KIT
6	2371134	1.000	1	GANTRY CABLES - LONG KIT

Table 6-10 Long Cable Collector (2281840-4 BOM, rev 0)



Imagination at work

www.gehealthcare.com