

GE Healthcare

LightSpeed 5.X Pro16

Illustrated Parts List

OPERATING DOCUMENTATION



2349187-100
Rev 18

Chapter 3

PDU

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2326492	1.0	No	NGPDU (H4 POWER DISTRIBUT {For further detail, see Section 1.0 , beginning on page 78 .}

Table 3-1 H Power PDU (B7816PD BOM, rev A)

Section 1.0 — NGPDU (2326492)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2348074	1.0	No	NGPDU ASSEMBLY1 {For further detail, see Section 1.1 , beginning on page 80 .}
2	P9336BZ	1.0	No	EMBLEM
3	2111577	1.0	No	P9810AM RATING PLATE FILM
4	2357432	1.0	No	RATING PLATE FOR NGPDU
5	2362890	1.0	No	LABEL
6	2117390	1.0	No	LABEL, CE MARKING COMPLIANCE
7	5139485	1.0	No	WEEE LABEL WITH BIN ONLY

Table 3-2 NGPDU (H4 Power Distribution Unit) (2326492 BOM, rev 3)

1.1 NGPDU Assembly (2348074)

FRU = NO




Figure 3-1 NGPDU Assembly 1 (2348074, sheet 1, rev 0)





Figure 3-2 NGPDU Assembly 1 (2348074, sheet 2, rev 1)



MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2354578-2	1.0	Yes	FRONT COVER
4	2354579	2.0	Yes	SIDE COVER
5	2354580	1.0	Yes	REAR COVER
7	2354581	1.0	Yes	TOP COVER
10	2354559-2	1.0	Yes	PANEL PLATE
48	N9408HR	25.0	Yes	SCREW SCREW TP
49	N9420JS	8.0	Yes	SCREW PAN HEAD 4 MM 20 MM SCREW osX 4X20 FE ZN
50	N9510HR	39.0	Yes	SCREW TP
51	N9501WS	41.0	Yes	WASHER PLAIN - LARGE 5.3 MM 10 MM WASHER M5
53	N9501BS	10.0	Yes	NUT HEXAGON 5 MM NUT M5 FE ZN
54	N9508JS	6.0	Yes	SCREW NABE 5X8 FE ZN
60	2103585	38.0	Yes	SCREW PAN HEAD 4 MM 8 MM SPRING WASHER, ZINC PLATED, STEEL
63	N9816YX	13.0	Yes	SCREW HEXAGON SOCKET 8 MM 16 MM CAP SCREW 8X16 FE
65	N9801WS	25.0	Yes	WASHER PLAIN - LARGE 8.4 MM 17 MM
68	N9514JS	44.0	Yes	SCREW PAN HEAD 5 MM 14 MM SCREW osX 5X14 FE ZN
70	N9414JS	8.0	Yes	SCREW PAN HEAD 4 MM 14 MM SCREW osX 4X14 FE ZN
73	N9306HR	6.0	Yes	SCREW TP 3X6 FE ZN
75	U0203AC	15.0	Yes	CLAMP CV_200
76	U0208AC	15.0	Yes	CLAMP TC-142 TAIRATSUPU KOTEI
77	U0212AC	10.0	Yes	CLAMP TMS_20
78	U0204AC	50.0	Yes	CLAMP CV_300
79	U0205AC	15.0	Yes	CLAMP CV_350
80	U0508AC	1.0	Yes	PROTECTOR EDGING CE-024
81	U0502AC	2.0	Yes	PROTECTOR BUSHING KG-024
85	N9825YX	15.0	Yes	SCREW CAP SCREW 8X25 Fe
88	46-220181P5	4.0	Yes	WASHER #10 (.190) SCREW SIZE, ID .194-.202 OD .334, WIDTH .062 THK .047, STEEL ZINC PLT

Table 3-3 NGPDU Assembly (2348074 BOM, Rev 16)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
90	2334820-2	1.0	Yes	PDU CONTROL BOARD
92	2405168	1.0	Yes	CB ASSY
93	2351484	1.0	Yes	15A BREAKER
94	2351485	1.0	Yes	40A BREAKER
95	2351486	1.0	Yes	10A BREAKER
96	2351487	1.0	Yes	40A BREAKER
97	2295544	1.0	Yes	BRK-16A
98	2295545	1.0	Yes	BRK-20A
99	2295543	1.0	Yes	BRK-32A
100	2295546	1.0	Yes	BRK-4A
101	2351488	1.0	Yes	100A CONTACTOR
102	2351489	1.0	Yes	50A CONTACTOR
103	2351490	1.0	Yes	18-32A CONTACTOR
104	2351491	1.0	Yes	14-25A CONTACTOR
105	2351492	3.0	Yes	200A FUSE
107	2113764-3	6.0	Yes	CAPACITOR, 6UF, 370VAC PM,FILTER CAP C1,C2,C3,C7,C8,C9
108	2351494	2.0	Yes	4600UF/450V CAPACITOR
109	2356924	1.0	Yes	GROUND BLOCK
110	2356926	3.0	Yes	TERMINAL CONNECTOR
111	2295531	2.0	Yes	TERMINAL CONNECTOR
112	2298763	1.0	Yes	TERMINAL CONNECTOR
113	2356927	19.0	Yes	TERMINAL CONNECTOR
114	2298762	6.0	Yes	TERMINAL CONNECTOR
115	2356928	10.0	Yes	TERMINAL CONNECTOR
119	2351495	1.0	Yes	DIODE BRIDGE
120	2361774	1.0	Yes	+24V PS
121	2391041	1.0	Yes	SS RESISTER ASSY
122	2388230	1.0	Yes	U0053SP SWITCH

Table 3-3 NGPDU Assembly (2348074 BOM, Rev 16) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
123	2286488	1.0	Yes	LED ASSY
125	2356929	12.0	Yes	END BRAKET
204	2356907-2	1.0	Yes	SIGNAL CABLE ASSY
214	2356917-2	1.0	Yes	WIRE 18#ASSY
226	2396924	1.0	Yes	REAL SWITCH
1	2354582	1.0	No	NGPDU FRAME & TRANSFORMER ASSY
13	2354566	1.0	No	TOP PROTECT PANEL
14	2354564	1.0	No	FRONT PROTECT PANEL
15	2354567	1.0	No	REAR PROTECT PANEL
16	2355078	3.0	No	SUPPORT SCREW
17	2354563	2.0	No	FIXTURE
20	2354571	1.0	No	ALUMINIUM PANEL
21	2295889	1.0	No	INSULATED PANEL-1
22	2295890	6.0	No	COPPER BLOCK
25	2354572	1.0	No	INSULATED PANEL-2
26	2354569	2.0	No	CONDUCTION BAR-1
27	2354570	1.0	No	CONDUCTION BAR-2
30	2354557	1.0	No	BUS BAR
33	2356930	1.0	No	RAIL1
34	2356931	1.0	No	RAIL2
35	2356932	1.0	No	RAIL3
37	2354574	7.0	No	SCREW
38	2356179	1.0	No	PS SUPPORT PANEL
39	2354576	6.0	No	CLAMP
47	N9516JL	8.0	No	SCREW
64	N9830YX	12.0	No	SCREW HEXAGON SOCKET 8 MM 30 MM CAP SCREW 8X30 FE
66	N9530JS	4.0	No	SCREW PAN HEAD 5 MM 30 MM SCREW 5X30
71	N9800SP	1.0	No	WASHER SPRING WASHER 8.2X15.4X2 .~DrjhX?/- 8.2X15.4X2

Table 3-3 NGPDU Assembly (2348074 BOM, Rev 16) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
72	N9516JS	4.0	No	SCREW NABE 5X16 FE ZN *
74	2356925	8.0	No	CABLE CLAMP
82	N9820NS	3.0	No	SCREW HEXAGON HEAD 8 MM 20 MM YfhtXyn 8X20 FE ZN
83	N9801BS	3.0	No	NUT HEXAGON 8 MM FE ZN
84	2249348	4.0	No	U0503AC BUSHING
86	2390699	1.0	No	END COVER
87	2401535	4.0	No	SCREW COUNTERSUNK(HEA 5 MM 10 MM)
89	N9528TM	4.0	No	WASHER OUT OF CATALOG WASHER jj/sX 5.5X14X1 BS
124	2356923	1.0	No	FERRITE CORE
126	2357125	2.0	No	CABLE CONNECTOR
127	2357127	1.0	No	LABEL
128	2357128	1.0	No	LABEL
129	2357130	1.0	No	LABEL
130	2335556	1.0	No	WIRE 1#ASSY
131	2335557	1.0	No	WIRE 1#ASSY
132	2335558	1.0	No	WIRE 1#ASSY
133	2335559	1.0	No	WIRE 10#ASSY
134	2335560	1.0	No	WIRE 1#ASSY
135	2335561	1.0	No	WIRE 1#ASSY
136	2335562	1.0	No	WIRE 1#ASSY
137	2335563	1.0	No	WIRE 6#ASSY
138	2335564	1.0	No	WIRE 6#ASSY
139	2335565	1.0	No	WIRE 6#ASSY
140	2335566	1.0	No	WIRE 6#ASSY
141	2335567	1.0	No	WIRE 6#ASSY
142	2335568	1.0	No	WIRE 6#ASSY
143	2335569	1.0	No	WIRE 6#ASSY
144	2335570	1.0	No	WIRE 6#ASSY

Table 3-3 NGPDU Assembly (2348074 BOM, Rev 16) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
145	2335571	1.0	No	WIRE 6#ASSY
146	2335573	1.0	No	WIRE 6#ASSY
147	2335572	1.0	No	WIRE 10#ASSY
150	2335574	1.0	No	WIRE 8#ASSY
151	2335575	1.0	No	WIRE 8#ASSY
152	2335576	1.0	No	WIRE 8#ASSY
153	2335577	1.0	No	WIRE 8#ASSY
154	2335578	1.0	No	WIRE 8#ASSY
155	2335579	1.0	No	WIRE 8#ASSY
156	2335580	1.0	No	WIRE 8#ASSY
157	2335581	1.0	No	WIRE 8#ASSY
158	2335582	1.0	No	WIRE 8#ASSY
159	2335583	1.0	No	WIRE 10#ASSY
160	2335584	1.0	No	WIRE 8#ASSY
161	2335585	1.0	No	WIRE 8#ASSY
162	2335586	1.0	No	WIRE 8#ASSY
163	2335587	1.0	No	WIRE 8#ASSY
164	2335588	1.0	No	WIRE 10#ASSY
165	2335589	1.0	No	WIRE 10#ASSY
166	2335590	1.0	No	WIRE 10#ASSY
167	2335591	1.0	No	WIRE 10#ASSY
168	2335592	1.0	No	WIRE 8#ASSY
169	2335593	1.0	No	WIRE 10#ASSY
170	2335594	1.0	No	WIRE 10#ASSY
171	2335595	1.0	No	WIRE 10#ASSY
172	2335596	1.0	No	WIRE 10#ASSY
173	2335597	1.0	No	WIRE 10#ASSY
174	2335598	1.0	No	WIRE 10#ASSY

Table 3-3 NGPDU Assembly (2348074 BOM, Rev 16) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
175	2335599	1.0	No	WIRE 10#ASSY
176	2335601	1.0	No	WIRE 10#ASSY
180	2335602	1.0	No	WIRE 14#ASSY
181	2335603	1.0	No	WIRE 14#ASSY
182	2335604	1.0	No	WIRE 14#ASSY
183	2335605	1.0	No	WIRE 16#ASSY
184	2335606	1.0	No	WIRE 16#ASSY
185	2335607	1.0	No	WIRE 16#ASSY
187	2335609	1.0	No	WIRE 16#ASSY
188	2335610	1.0	No	WIRE 16#ASSY
189	2335612	1.0	No	WIRE 16#ASSY
190	2335613	1.0	No	WIRE 16#ASSY
191	2335614	1.0	No	WIRE 16#ASSY
192	2335615	1.0	No	WIRE 16#ASSY
193	2335616	1.0	No	WIRE 16#ASSY
194	2335617	1.0	No	WIRE 16#ASSY
195	2335618	1.0	No	WIRE 16#ASSY
196	2335619	1.0	No	WIRE 16#ASSY
197	2335620	1.0	No	WIRE 16#ASSY
200	2335622	1.0	No	WIRE 18#ASSY
201	2335623	1.0	No	WIRE 18#ASSY
202	2335624	1.0	No	WIRE 18#ASSY
205	2356908	1.0	No	WIRE 18#ASSY
206	2356909	1.0	No	WIRE 18#ASSY
207	2356910	1.0	No	WIRE 18#ASSY
208	2356911	1.0	No	WIRE 18#ASSY
209	2356912	1.0	No	WIRE 18#ASSY
210	2356913	1.0	No	WIRE 18#ASSY

Table 3-3 NGPDU Assembly (2348074 BOM, Rev 16) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
211	2356914	1.0	No	WIRE 18#ASSY
212	2356915	1.0	No	WIRE 18#ASSY
213	2356916	1.0	No	WIRE 18#ASSY
215	2356918	1.0	No	WIRE 18#ASSY
216	2356919	1.0	No	WIRE 18#ASSY
217	2356920	1.0	No	WIRE 18#ASSY
218	2356921	1.0	No	WIRE 10#ASSY
219	2356922	1.0	No	WIRE 10#ASSY
220	2357126	1.0	No	WIRE 10#ASSY
221	2335625	1.0	No	WIRE 10#ASSY
222	2357131	1.0	No	WIRE 10#ASSY
223	2357132	1.0	No	WIRE 10#ASSY
224	2357129	1.0	No	WIRE 14#ASSY
227	2396802	4.0	No	SPONGE
228	2395645	6.0	No	A TERMINAL NUMBER LABEL

Table 3-3 NGPDU Assembly (2348074 BOM, Rev 16) (Continued)

Section 2.0 — Major NGPDU Components

ITEM	PART NO.	NAME / DESCRIPTION	COMMENTS
1	2326492	NGPDU assemble 1	2326492(H-Power) non stockable
2	2326492-3	NGPDU assemble 3	2326492-3(H-Power 60) non stockable
3	2354578	Front Cover	
4	2354579	Side Cover	
5	2354580	Rear Cover	
6	2354581	Top Cover	
7	2334820	PDU Control Board	
9	2351483	110A CB ASSY	CB1
10	2351484	15A Breaker	CB2
11	2351485	40A Breaker	CB3
12	2351486	10A Breaker	CB4
13	2351487	40A Breaker	CB5
14	2295544	16A Breaker	CB6
15	2295545	20A Breaker	CB7
16	2295543	32A Breaker	CB5(2326492-2 Qxi), CB8
17	2295546	4A Breaker	CB9
18	2351488	100A Contactor	Kxg
19	2351489	50A Contactor	Ktg
20	2351490	18-32A Contactor	Ksv
21	2351491	14-25A Contactor	Kss
22	2351492	200A Fuse	F1-3(2326492 H-Power)
23	2351493	100A Fuse	F1-3(2326492-3 H-Power 60)
24	2113764-3	6uf / 370Vac Capacitor	C1-6
25	2351494	4600uf / 450V Capacitor	C7, C8
26	2351495	Diode Bridge	BR1
27	2201229-2	Power Supply	PS

Table 3-4 Major Components

ITEM	PART NO.	NAME / DESCRIPTION	COMMENTS
28	2391041	Resistor Assy	R
29	2388230	Emergency Switch	SW
30	2286488	LED Assy	LED
31	2396924	Gantry Enable Switch	SW including LED

Table 3-4 Major Components (Continued)

Table of Contents - Chapter 4 (Gantry)

Section 1.0	Vertical Assembly (2357738-8, -10)	100
1.1	S/R 120V Power and Intercom Board Assembly (2357521, -2)	126
1.1.1	Transmitter/Receiver Power Supply (2354605)	129
1.2	DAS-Detector Plate Assembly (2358133-4)	130
1.2.1	GDAS EMC Collector Asm (2358139-3)	132
1.2.2	Watson16 Det Assembly with Diode (2354053)	134
1.2.3	DAS Duct Fan Assembly (2333651)	138
1.3	Drive Assy Hpower Axial (2217930-3)	140
1.3.1	Axial Motor with Sprocket (2235342-3)	144
1.3.2	Brake Assy, H2 AX Drive (2277225)	146
1.3.3	AX Drive Motor Controller (2233402)	148
1.4	ORP Assembly (2358141, -2)	152
1.5	850 MB H-Power Slipring (2349589)	154
1.5.1	850 MB Platter HSDCD Slip (2320265)	156
1.5.1.1	850 MB Receiver-Antenna (2333609)	158
1.5.1.2	Brushblock Assembly For H-power (2349598)	160
	Power Brush Tip Collector (2349603)	161
	Signal Brush Tip Collector (2349604)	161
1.5.1.3	Slipring Assembly For H-power (2349587)	162
1.6	Gantry Tilt Feedback ASM, EMC ED2 (2356863-2)	164
1.7	TGPU ED2 Assembly (5111425-2)	166
1.8	Service Outlet/EStop Assembly (2349257)	168
1.9	Power Pan Assembly 550V (2387385)	170
1.10	Power Pan Assembly - 120V (2334726, -2)	174
1.11	GDAS-16 Main Harness, With ASTEC PS's (2394982-5)	178
Section 2.0	Vertical Assembly (2357738, -3, -4)	180
2.1	Detector Plate Assembly (2358133, -2, -3)	212
2.1.1	DAS Collector Assembly (2358139, -2)	216
2.1.1.1	BP/Chassis LT Assembly - H4 (2371625)	220

	<i>Chassis LT Assembly - H4 (2371628)</i>	222
	<i>2.1.1.2 – BP/Chassis CTR Assembly - H4 (2371626)</i>	224
	<i>Chassis CTR Assembly - H4 (2371630)</i>	226
	<i>2.1.1.3 – BP/Chassis RT Assembly - H4 (2371627)</i>	228
	<i>Chassis RT Assembly - H4 (2371631)</i>	230
2.2	H-Power Axial Dynamic Brake Module (2357524)	232
	<i>2.2.1 – Step Up Transformer (2354743)</i>	235
2.3	Encoder Assembly (2369747)	236
2.4	Power Supply ASM, 5V Logic (2334458)	238
2.5	Tilt, Feedback Assembly (2356863)	240
2.6	MSUB/TGP Assembly (2349209)	242
	<i>2.6.1 – MSUB/TGP Power Supply Assembly (2349205)</i>	244
2.7	Harness, GDAS DET/ORP/P.S./DHCB/ (2394982)	246
Section 3.0	Horizontal Assembly (2357737-8, -10)	248
3.1	H-power Gantry High Level (2357736)	262
	<i>3.1.1 – Hydraulic Tilt Motor Assembly (2334727)</i>	266
3.2	Helios 4 collimator with notched filler bar and rating plate (5140001-4)	270
	<i>3.2.1 – Helios 4 collimator with notched filler bar (5140000-4)</i>	272
	<i>3.2.1.1 – H4 Collimator Base Assembly (2273645-3)</i>	276
3.3	Jedi-HP Collector For CT (2219132)	278
	<i>3.3.1 – JH4_Inverter (2281950)</i>	280
	<i>3.3.2 – Jedi H4 MPHV Tank (2266521)</i>	290
	<i>3.3.3 – Inverter/auxiliary Parts (2368331)</i>	300
3.4	Hercules 2 Fan Heat Exchanger (2334365)	302
	<i>3.4.1 – Hercules H-Power Heat Exchanger Assembly (5125371)</i>	304
3.5	Home Flag, PWB MTG, H4 Gantry (2227383-2)	312
Section 4.0	Horizontal Assembly (2357737, -3, -4)	314
4.1	Helios 4 Collimator (2214768-4)	326
	<i>4.1.1 – Helios 4 Collimator (2120100-4)</i>	328
4.2	Power Supply Assembly, 24V DC (2334457)	332
4.3	Fuse Box Assy - H4 (2331963-2)	334

4.3.1 – Harness, FSEBX/ORP/DAS/INV (2361658).....	338
4.4 Power Supply Assembly, 12V DC Parallel (2334459).....	340
Section 5.0 Cover Assembly, Front Gantry (2361872)	342
5.1 Cover Sub-assembly, Front Gantry (2361871).....	344
5.2 P9213YX MIC Assembly (2111955)	348
Section 6.0 Cover Assembly, Rear Gantry (2361874).....	350
6.1 Cover Sub-Assembly, Rear Gant (2361873).....	352
Section 7.0 Assembly, H Power Top Cover (2357525)	356
Section 8.0 Cover Asm, Side Gantry Right (2276690, -5).....	358
Section 9.0 Cover Asm, Side Gantry Left (2276691-7)	360
Section 10.0 Cover Asm, Side Gantry Left (2276691-5)	362

End of TOC.