

Chapter 6

Ship Collector (5323140)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	46-229452G1	1.00	Yes	CT-9800 (ACCESSORIES)
2	46-236664P1	1.00	Yes	STRAP STRAP, WRIST 2.00 WIDE X 18.00 LONG WHITE
3	2381848-4	1.000	Yes	COLLECTOR-BASE COVER-LARGER OPENNING-RT
4	46-297592P1	1.000	Yes	BAG HOLDER
6	2105873	10.000	Yes	.531 I.D., 1.625 O.D., .25THK GANTRY & TABLE INSTALLATION ANCHOR WASHER, FOR BASEPLATE-LESS HSA SYSTEMS.
7	2106573	10.000	Yes	HILTI ANCHOR STUD 1/2 DIA X 8.0 LG, KWIK-BOLT II WITH EXTENDED THREAD, FOR BASEPLATE-LESS HSA SYSTEMS STUD MAT~L: CARBON ST PER ASTM A510 WEDGE MAT~L: CARBON ST PER AISI 1010 NUT MAT~L: CARBON ST PER ASTM A563 GRADE A WASHER MAT~L: CARBON ST PER SAE 1005-1033 FINISH: ALL PARTS ZINC ELECTROPLATED CLEAR CHROMATE PER ASTM B633 (1998), TYPE 2, CLASS 5.
10	2218570-2	1.000	Yes	H-POWER FUSE KIT {For further detail, see Section 1.0, beginning on page 411.}
13	46-297960P1	1.000	Yes	FOAM INSERT
18	46-222436G1	4.000	Yes	CT-9800 - CHIN STRAP
20	46-222443P1	1.000	Yes	OPEN CELL POLYESTER FOAM
21	46-222444P1	1.000	Yes	CLOSED CELL LINEAR ETHYLE
22	46-222445P1	1.000	Yes	CLOSED CELL LINEAR ETHYLE
23	46-222921P1	6.000	Yes	OBSOLETE FOR MR SITES PER
24	46-199388P1	1.000	Yes	HEAD HOLDER POLYURETHANE
25	5180151	2.000	Yes	REAR Cover DOLLY
26	5179183	1.000	Yes	Dolly, pivoting, small room, RH
27	5179183-2	1.000	Yes	Dolly, pivoting, small room, LH
31	2282509-2	1.000	Yes	substitute of 2282509, HOIST-2-PC ARM

Table 6-1 LightSpeed RT¹⁶ Vision SHIP COLLECTOR (5323140 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
32	2282510-2	1.000	Yes	substitute of 2282510Hoist-2-PC, LONG TUBE
33	2277548	4.000	Yes	PLUG, SIDE COVER H2
34	2217173-2	1.000	Yes	STOOL WITH FOOT REST AND
36	2376682-2	1.000	Yes	BALANCE KIT FOR PATHFINDER, FIELD, HPOWER {For further detail, see Section 3.0, beginning on page 414.}
37	2309994-2	1.000	Yes	IV POLE ASSY WITH CHINESE LABEL
40	2364629	1.000	Yes	LABEL INSTALLATION KIT
42	2336557	1.000	Yes	SEISMIC KIT FOR GLOBAL CO
48	2386850	2.000	Yes	TRIAL WEIGHT, SG GANTRY
49	5118556	1.000	Yes	Jumper flex tool
50	2373436-7	1.000	Yes	Cable, FE Laptop to Console, 23000mm Long
52	5314211	5.000	Yes	HILTI ANCHOR STUD 0.625inch DIA X 8Inch LG KWIK-BOLT III WITH EXTENDED THREAD, FOR BASEPLATE-LESS HSA SYSTEMS
53	5315833	1.000	Yes	drill bushing
54	5314487	5.000	Yes	washer
14	2319891-5	1.000	No	INSTALLATION SUPPORT KIT FOR RT16 {For further detail, see Section 2.0, beginning on page 412.}
19	46-237412P1	1.000	No	STRAP
29	5164728	1.000	No	RT FLOOR TEMPLATE
35	5112283	1.000	No	CT Service Cabinet from Pentair
38	2329064-2	1.000	No	TABLE ACC ASSY WITH CHINESE LABEL
39	2355100	1.000	No	TRAY, H2 TABLE
46	46-237605P1	1.000	No	WRIST STRAP
47	2393793	2.000	No	BAR, FOOT SWITCH-SG
51	2364311	2.000	No	TAPERED PLASTIC HANDLE - METAL INSERT

Table 6-1 LightSpeed RT¹⁶ Vision SHIP COLLECTOR (5323140 BOM, rev 1) (Continued)

Section 1.0 — H-Power Fuse Kit (2218570-2)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2351492	2.000	Yes	200A FUSE
2	2364059	2.000	Yes	GLASS FUSE
3	46-170021P50	2.000	Yes	FUSE 12 AMPS 250 VOLTS
4	46-170021P106	2.000	Yes	FUSE 8 AMPS 250 VOLTS
5	2336517-2	2.000	Yes	FUSE 25 AMPS 700 VOLTS 2.
6	46-170021P52	3.000	Yes	FUSE 3 AMPS 250 VOLTS
7	46-170021P10	2.000	Yes	006.000A 0250V 3AG
8	46-170021P101	1.000	Yes	FUSE 20 AMPS 700 VOLTS
9	46-170021P11	2.000	Yes	7.0A, 250V SLO-BLO FUSE.
10	46-170021P31	2.000	Yes	1/2A, 250V SLO-BLO FUSE.
11	2106993-5	2.000	Yes	FUSE 20 AMPS 500 VOLTS .3
12	46-327160P1	2.000	Yes	12.0A, 125VAC, DUAL ELEME
13	2379651	1.000	Yes	FUSE - 700 V, 200 A, FAST

Table 6-2 H-Power Fuse Kit (2218570-2 BOM, rev 3)

Section 2.0 — Installation Support Kit - LightSpeed RT¹⁶ (2319891-5)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
28	2344539	1.000	Yes	16 SLICE DAS/DETECTOR SERVICE KIT
39	2214163	1.000	Yes	CONNECTOR, ANALOG LOOP BA
45	2284023-3	1.000	Yes	SERVICE TOOL COLLECTOR, H
59	2106205	1.000	Yes	500 I.D. 1.00 O.D. 2.125 LONG, ANCHOR HOLE CENTERING ON BASEPLATE- LESS HSA SYSTEMS
60	2107863	1.000	Yes	SPECIAL SOCKET FOR LEVELING ADJUSTERS USED ON BASEPLATE-LESS HSA SYSTEMS.
8	46-170923P1	1.000	No	2-PART, EPOXY PATCH
10	46-221561P3	1.000	No	CLAMP OR CLIP
11	46-208758P10	25.000	No	8.00" X .094" SELF-LOCKIN
12	46-208758P3	50.000	No	7.31" X .184" SELF-LOCKIN
13	46-208758P5	50.000	No	SELF LOCKING CABLE TIE
14	46-208759P5	15.000	No	6.00" X .140" SELF-LOCKIN
15	46-208759P6	20.000	No	7.81" X .184" SELF-LOCKIN
16	46-208759P4	20.000	No	13.93"X .301" SELF-LOCKIN
17	46-208747P1	25.000	No	CABLE TIE,ADH.BACKED
19	46-230118P1	2.000	No	PAPER WYP
20	2253867	1.000	No	SCOTCH BRIGHT SCOUR PAD M
21	46-170127P5	1.000	No	1.0" MASKING TYPE TAPE, 6
24	46-170686P2	1.000	No	LOCTITE 242, 10CC 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER
26	46-183039P1	1.000	No	TEXWIPE 58019-506-91/9%
29	46-208781P3	1.000	No	WRITE-ON LABEL FOR MARKIN
41	46-303460P2	1.000	No	PAINT GE GRAY #1
42	46-303460P3	1.000	No	PAINT GE STEEL BLUE
43	46-303460P4	1.000	No	PAINT GE DARK STEEL BLUE

Table 6-3 Installation Support Kit for RT¹⁶ (2319891-5 BOM, rev 3)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
44	46-220178P10	2.000	No	GROMMET
46	46-136246P14	6.000	No	2" WIDE, WHT, ADH BACK, V
47	2278989	6.000	No	LANYARD, NYLON U-LOCK, TY
48	46-136246P15	6.000	No	3/4" WIDE, WHT ADH BACK,
49	46-170123P1	1.000	No	WIDTH 0.625 THK .007 C
50	46-136246P21	6.000	No	2" WIDE X 25 YD WHITE VEL
52	46-136246P17	6.000	No	3/4" WIDE, WHT ADH BACK,
53	46-221561P2	5.000	No	VINYL TRIM W/ SEGMENTED M
54	46-208753P3	6.000	No	ID 0.375 PAD WIDTH 0.75
55	46-170123P2	3.000	No	3/4" WIDTH, .007 BLACK VI
56	5113046	1.000	No	0.5mm Wire for Gantry Rotation Characterization Procedure, SPRING TEMPER CARBON STEEL WIRE 0.020inches Diameter, 25.4 mm length Pck Qrt 100
57	5111615	2.000	No	external laser light disconnect label
58	5111733	2.000	No	Label for front cover to explain that GE external laser lights are not in use

Table 6-3 Installation Support Kit for RT¹⁶ (2319891-5 BOM, rev 3) (Continued)

Section 3.0 — BALANCE KIT, FIELD, H-POWER (2376682-2)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2316653	6.000	Yes	WEIGHT,107 DEG,6MM
2	2316653-2	6.000	Yes	WEIGHT,107 DEG,1.5MM
3	2317863	11.000	Yes	SHIM,180 DEG,3.2MM
4	2317865-2	2.000	Yes	SPACER,180 DEG,25.4MM
5	2317994	20.000	Yes	WEIGHT,TRIM
6	2319142	2.000	Yes	SPACER,107 DEG,25.4MM
7	2319145	2.000	Yes	SHIM,107 DEG,3.3MM
8	2342895-2	6.000	Yes	WEIGHT,180 DEG,6MM
9	2342895-3	6.000	Yes	WEIGHT,180 DEG,1.5MM
12	46-328425P4	4.000	Yes	NUT HEXAGON 16 MM
13	46-328417P70	2.000	Yes	SCREW HEXAGON SOCKET 12 M
14	2376683	1.000	Yes	SOCKET, 15/16" DEEP IMPAC
18	46-328432P5	4.000	Yes	WASHER LOCK - SPRING 16.2
10	2370053	1.000	No	WEIGHT,TRIAL,NEGATIVE Z
11	2370227	1.000	No	WEIGHT,TRIAL,POSITIVE Z,
15	2377399	1.000	No	LEFT SIDE DAS COVER TOOL
16	2377400	1.000	No	RIGHT SIDE DAS COVER TOOL
17	2377401	1.000	No	CENTER DAS COVER TOOL

Table 6-4 Balance Kit, Field, H-Power (2376682-2 BOM, rev 2)

Section 4.0 — CABLE COLLECTORS

4.1 Short Cable Collector (2281840-5)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2371450-2	1.000	Yes	RUN#55A-PDU TO RACEWAY GND
2	2371450-4	1.000	Yes	RUN#56A-CONSOLE TO RACEWAY
3	2373436-3	1.000	Yes	RUN #102A - GANTRY TO CON
4	2117848-7	1.000	Yes	20000 +/- 150 FIBER OPTIC
5	2371133-2	1.000	Yes	CONSOLE CABLES - SHORT KIT
6	2371134-2	1.000	Yes	GANTRY CABLES - SHORT KIT

Table 6-5 Short Cable Collector (2281840-5 BOM, rev 0)

4.2 Long Cable Collector (2281840-4)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2371450	1.000	Yes	RUN#55-PDU TO RACEWAY GND
2	2371450-3	1.000	Yes	RUN#56-CONSOLE TO RACEWAY
3	2373436-2	1.000	Yes	RUN #102 - GANTRY TO CONS
4	2117848-2	1.000	Yes	24384 +/- 150 FIBER OPTIC
5	2371133	1.000	Yes	CONSOLE CABLES - LONG KIT
6	2371134	1.000	Yes	GANTRY CABLES - LONG KIT

Table 6-6 Long Cable Collector (2281840-4 BOM, rev 0)

