
Chapter 5

Miscellaneous Itemsf

Section 1.0 — HP60 Ship Collector (2385414, -2, -3, -4 / 5135131, -2)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2290278-4	1.0	Yes	BASE COVER COLLECTOR {For further detail, see Section 1.1, beginning on page 479. }
3	2381045-2	1.0	Yes	HP60 FIELD BALANCE WEIGHT KIT {For further detail, see Section 1.3, beginning on page 480. }
4	2385412	1.0	Yes	P2030BE FUSE KIT HP60 90KVA {For further detail, see Section 1.2, beginning on page 479. }
5	2281840-5	1.0	Yes	H-POWER SHORT CABLE COLLECTOR {For further detail, see Section 2.1, beginning on page 484. }
1	2385413	1.0	No	P2030BD SHIP C HP-100-60 COMMON {For further detail, see Section 1.4, beginning on page 481. }

Table 5-1 HP60-16 SHIP COLLECTOR (2385414 BOM, rev 3)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2290278-4	1.0	Yes	BASE COVER COLLECTOR {For further detail, see Section 1.1, beginning on page 479. }
3	2381045-2	1.0	Yes	HP60 FIELD BALANCE WEIGHT KIT {For further detail, see Section 1.3, beginning on page 480. }
4	2385412	1.0	Yes	P2030BE FUSE KIT HP60 90KVA {For further detail, see Section 1.2, beginning on page 479. }
5	2227089-2	1.0	Yes	DAS/DETECTOR SERVICE KIT {For further detail, see Section 1.5, beginning on page 483. }
6	2281840-5	1.0	Yes	H-POWER SHORT CABLE COLLECTOR {For further detail, see Section 2.1, beginning on page 484. }
1	2385413	1.0	No	P2030BD SHIP C HP-100-60 COMMON {For further detail, see Section 1.4, beginning on page 481. }

Table 5-2 HP60 4 and 8 slice Ship Collector (2385414-2 BOM, rev 4)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2290278-4	1.0	Yes	BASE COVER COLLECTOR {For further detail, see Section 1.1, beginning on page 479. }
3	2381045-2	1.0	Yes	HP60 FIELD BALANCE WEIGHT KIT {For further detail, see Section 1.3, beginning on page 480. }
4	2385412	1.0	Yes	P2030BE FUSE KIT HP60 90KVA {For further detail, see Section 1.2, beginning on page 479. }
	5146898	1.0	Yes	AMC CONFIGURATION TOOL
1	2385413	1.0	No	P2030BD SHIP C HP-100-60 COMMON {For further detail, see Section 1.4, beginning on page 481. }

Table 5-3 HP60-16 Ship Collector - PF1.0 (2385414-3 BOM, rev 3)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2290278-4	1.0	Yes	BASE COVER COLLECTOR {For further detail, see Section 1.1, beginning on page 479. }
3	2381045-2	1.0	Yes	HP60 FIELD BALANCE WEIGHT KIT {For further detail, see Section 1.3, beginning on page 480. }
4	2385412	1.0	Yes	P2030BE FUSE KIT HP60 90KVA {For further detail, see Section 1.2, beginning on page 479. }
5	2227089-2	1.0	Yes	DAS/DETECTOR SERVICE KIT {For further detail, see Section 1.5, beginning on page 483. }
	5146898	1.0	Yes	AMC CONFIGURATION TOOL
1	2385413	1.0	No	P2030BD SHIP C HP-100-60 COMMON {For further detail, see Section 1.4, beginning on page 481. }

Table 5-4 HP60 4, 8 Ship Collector - PF1.0 (2385414-4 BOM, rev 3)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2290278-4	1.0	Yes	BASE COVER COLLECTOR {For further detail, see Section 1.1, beginning on page 479. }
3	2381045-2	1.0	Yes	HP60 FIELD BALANCE WEIGHT KIT {For further detail, see Section 1.3, beginning on page 480. }
4	2385412	1.0	Yes	P2030BE FUSE KIT HP60 90KVA {For further detail, see Section 1.2, beginning on page 479. }
1	2385413	1.0	No	P2030BD SHIP C HP-100-60 COMMON {For further detail, see Section 1.4, beginning on page 481. }

Table 5-5 HP60-16 SHIP COLLECTOR, Without short cable collector - PF1.5 (5135131 BOM, rev 2)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2290278-4	1.0	Yes	BASE COVER COLLECTOR {For further detail, see Section 1.1, beginning on page 479. }
3	2381045-2	1.0	Yes	HP60 FIELD BALANCE WEIGHT KIT {For further detail, see Section 1.3, beginning on page 480. }
4	2385412	1.0	Yes	P2030BE FUSE KIT HP60 90KVA {For further detail, see Section 1.2, beginning on page 479. }
5	2227089-2	1.0	Yes	DAS/DETECTOR SERVICE KIT {For further detail, see Section 1.5, beginning on page 483. }
1	2385413	1.0	No	P2030BD SHIP C HP-100-60 COMMON {For further detail, see Section 1.4, beginning on page 481. }

Table 5-6 HP60 4 and 8 SLICE SHIP COLLECTOR, Without short cable collector - PF1.5 (5135131-2 BOM, rev 2)

1.1 Base Cover Collector (2290278-4)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	46-328417P3	4.0	Yes	SCREW HEXAGON SOCKET 4 MM 10 MM ZINC PLATED, STEEL
2	46-328430P2	25.0	Yes	WASHER PLAIN - NORMAL 4.3 MM 9 MM ZINC PLATED, STEEL
3	46-328432P2	29.0	Yes	WASHER LOCK - SPRING 4.1 MM 7.6 MM ZINC PLATED, SPRING STEEL
4	46-328425P1	5.0	Yes	NUT HEXAGON 4 MM
5	2251467-4	1.0	Yes	COVER, BASE-FRONT
6	2285829-7	1.0	Yes	COVER, LEFT BASE WITH CARDIAC OPT
7	2285830-7	1.0	Yes	COVER, RIGHT BASE WITH CARDIAC OPT
8	2275686-4	1.0	Yes	COVER ASSY, BASE REAR LHY
9	2275685-4	1.0	Yes	COVER ASSY, BASE REAR RH
10	2277652	3.0	Yes	LATCH, H2 SIDE COVER COMPRESSION
11	2261321	2.0	No	BRACKET - FRONT BASE COVER
12	2261324	2.0	No	BRACKET - SIDE BASE COVER
13	2261325	1.0	No	BRACKET - REAR BASE COVER
14	2285469	2.0	No	BRACKET, BASE COVER MOUNTING
15	2288882	2.0	No	BULKHEAD - BLANK
16	46-328420P4	4.0	No	SCREW BUTTON HEAD 4 MM 16 MM BLACK OXIDE, STEEL
18	46-328417P4	17.0	Yes	SCREW HEXAGON SOCKET 4 MM 12 MM ZINC PLATED, STEEL

Table 5-7 BASE COVER COLLECTOR (2290278-4 BOM, rev 0)

1.2 Fuse Kit HP60 90kVA (2385412)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2351493	3.0	Yes	100A FUSE
2	2364059	2.0	Yes	GLASS FUSE

Table 5-8 Fuse Kit HP60 90kVA (2385412 BOM, rev 2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	46-170021P50	2.0	Yes	FUSE 12 AMPS 250 VOLTS BUSSMAN MDA12 DUAL ELEMENT.
4	46-170021P15	2.0	Yes	FUSE 8 AMPS 250 VOLTS (REVIEWED TK, 7/91).
5	2336517-2	2.0	Yes	FUSE 25 AMPS 700 VOLTS 2.5M SECONDS
6	46-170021P52	3.0	Yes	FUSE 3 AMPS 250 VOLTS (REVIEWED TK, 7/91).
7	46-170021P10	2.0	Yes	006.000A 0250V 3AG FAST UL/C (REVIEWED TK, 7/91).
8	46-170021P101	1.0	Yes	FUSE 20 AMPS 700 VOLTS 2~ LONG X 9/16~ DIA FAST BLOW
9	46-170021P106	2.0	Yes	FUSE 8 AMPS 250 VOLTS BUSSMANN ONLY
10	46-170021P31	2.0	Yes	1/2A, 250V SLO-BLO FUSE. TYPE 3AG, 1.25~ X 0.25~ GLASS BODY. (REVIEWED TK, 7/91).
11	2106993-5	2.0	Yes	FUSE 20 AMPS 500 VOLTS .3 SECONDS
12	46-327160P1	2.0	Yes	12.0A, 125VAC, DUAL ELEMENT TIME DLY FNM-12, 0.406 X 1.5 LG, FIBER BODY
13	2379651	1.0	Yes	FUSE - 700 V, 200 A, FAST SEMI
15	2238207-3	1.0	Yes	FUS FF 14X51 2A 700VDC C 50KA ULR
16	99183979	1.0	Yes	SLOW BLOW FUSE HCP 6X32 10A 250V UL
	46-170021P11	2.0	Yes	7.0A, 250V SLO-BLO FUSE. * TYPE 3AG, 1.25~ X 0.25~ GLASS BODY. (REVIEWED TK, 7/91).

Table 5-8 Fuse Kit HP60 90kVA (2385412 BOM, rev 2) (Continued)

1.3 HP60 Field Balance Weight Kit (2381045-2)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2378669	11.0	Yes	P2033DU 78WEIGHT-M
2	2378670	4.0	Yes	P2033DV 78WEIGHT-S
4	2317865-2	2.0	Yes	SPACER,180 DEG,25.4MM
5	2317994	7.0	Yes	WEIGHT,TRIM
8	2342895-2	3.0	Yes	WEIGHT,180 DEG,6MM
9	2342895-3	4.0	Yes	WEIGHT,180 DEG,1.5MM
12	46-328425P4	4.0	Yes	NUT HEXAGON 16 MM

Table 5-9 HP60 FIELD BALANCE WEIGHT KIT (2381045-2 BOM, rev 2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
13	46-328417P70	2.0	Yes	SCREW HEXAGON SOCKET 12 MM 100 MM
14	2376683	1.0	Yes	MCMaster CARR PN 6023A21 SOCKET 0.9375 IN DEEP IMPACT 0.500 IN DRV
18	46-328432P5	4.0	Yes	WASHER LOCK - SPRING 16.2 MM 27.4 MM ZINC PLATED, SPRING STEEL
20	2378668	1.0	Yes	P2033DT 78WEIGHT-L
21	46-328432P7	2.0	Yes	WASHER LOCK - SPRING 12.2 MM 21.1 MM
22	2236649-5	2.0	Yes	WASHER, FLAT, ISO7089,M12, ST-HB200
24	2317863	9.0	Yes	SHIM,180 DEG,3.2MM
10	2370053	1.0	No	WEIGHT,TRIAL,NEGATIVE Z
11	2370227	1.0	No	WEIGHT,TRIAL,POSITIVE Z, H-POWER
15	2377399	1.0	No	LEFT SIDE DAS COVER TOOL
16	2377400	1.0	No	RIGHT SIDE DAS COVER TOOL
17	2377401	1.0	No	CENTER DAS COVER TOOL
23	5128067	2.0	No	CAP SCREW U0480AA

Table 5-9 HP60 FIELD BALANCE WEIGHT KIT (2381045-2 BOM, rev 2) (Continued)

1.4 P2030BD Ship C HP-100/60 Common (2385413)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	46-229452G1	1.0	Yes	CT-9800 (ACCESSORIES)
2	46-297482G1	1.0	Yes	IV BOARD ACC FOR ZEUS TABLE
4	46-297592P1	1.0	Yes	BAG HOLDER
6	2105873	10.0	Yes	.531 I.D., 1.625 O.D., .25THK GANTRY & TABLE INSTALLATION ANCHOR WASHER, FOR BASEPLATE-LESS HSA SYSTEMS.
7	2106573	10.0	Yes	HILTI ANCHOR STUD 1/2 DIA X 8.0 LG, KWIK-BOLT II WITH EXTENDED THREAD, FOR BASEPLATE-LESS HSA SYSTEMS STUD MAT~L: CARBON ST PER ASTM A510 WEDGE MAT~L: CARBON ST PER AISI 1010 NUT MAT~L: CARBON ST PER ASTM A563 GRADE A WASHER MAT~L: CARBON ST PER SAE 1005-

Table 5-10 P2030BD Ship C HP-100/60 Common (2385413 BOM, rev 7)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
8	2115984	1.0	Yes	PAD FOR HSA METALLESS FOOT EXTENDER
12	46-222436G1	1.0	Yes	CT-9800 - CHIN STRAP
13	46-297960P1	1.0	Yes	FOAM INSERT
18	46-222436G1	3.0	Yes	CT-9800 - CHIN STRAP
20	46-222443P1	1.0	Yes	OPEN CELL POLYESTER FOAM
21	46-222444P1	1.0	Yes	CLOSED CELL LINEAR ETHYLENE
22	46-222445P1	1.0	Yes	CLOSED CELL LINEAR ETHYLENE
23	46-222921P1	6.0	Yes	OBSOLETE FOR MR SITES PER FMI 60288 POLYURETHANE FOAM 1.7#/FT DENSITY S.E. FINE CELL CHARCOAL GRAY
24	46-199388P1	1.0	Yes	HEAD HOLDER POLYURETHANE FOAM
25	5180151	2.0	Yes	Rear cover dolly
26	5179183	1.0	Yes	Dolly, pivoting, small room, RH
29	2363851	1.0	Yes	TEMPLATE, INSTALLATION HPOWER
30	5179183-2	1.0	Yes	Dolly, pivoting, small room, LH
31	2282509	1.0	Yes	HOIST-2-PC, ARM, HELIOS 2
32	2282510	1.0	Yes	HOIST-2-PC, LONG TUBE, HELIOS 2
33	2277548	4.0	Yes	PLUG, SIDE COVER H2
34	2217173-2	1.0	Yes	STOOL WITH FOOT REST AND NO ARM REST
40	2364629	1.0	Yes	LABEL INSTALLATION KIT
42	2336557	1.0	Yes	SEISMIC KIT FOR GLOBAL CONSOLE
43	2373436-7	1.0	Yes	Cable, FE Laptop to Console, 23000mm Long
9	2115994	1.0	No	METALLESS CRADLE PAD & STRAPS
14	2319891-3	1.0	No	INSTALLATION SUPPORT KIT FOR H-POWER
19	46-237412P1	1.0	No	STRAP HEAD STRAP 3 STRAPS PER PACKAGE
35	5112283	1.0	No	CT Service Cabinet from Pentair
37	2309994-2	1.0	No	IV POLE ASSY, with Multi-Language Labels
38	2329064-2	1.0	No	TABLE TRAY ACC ASSY, with Multi-Language Lables
39	2355100	1.0	No	TRAY, H2 TABLE

Table 5-10 P2030BD Ship C HP-100/60 Common (2385413 BOM, rev 7) (Continued)

1.5 DAS/Detector Service Kit (2227089-2)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	46-183039P1	24.0	Yes	TEXWIPE 58019-506-91/9% ALCOHOL PADS
2	2169110-2	1.0	Yes	ONE PKG OF 25 PC ELASTOMER SHIN-ETSU
3	2207303-6	1.0	Yes	SAFE SKIN NITRILE GLOVE, LARGE
4	2207303-7	1.0	Yes	SAFE SKIN NITRILE GLOVE, EXTRA LARGE
6	2335064	1.0	Yes	DUSTER-AIR SPRAY SYSTEM
8	2203205	20.0	Yes	ANTISTATIC BAG, 7 X 20~
9	2254888-2	1.0	Yes	GB TWIN ELAST. (25PCS) IN CONTAINER
5	2185881	8.0	No	ESD BOOT-NP-DET-TOOL

Table 5-11 DAS/Detector Service Kit (2227089-2 BOM, rev 0)

Section 2.0 — Cable Collectors

2.1 Short Cable Collector (2281840-5)

FRU = 1

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2371450-2	1.0	Yes	RUN#55A-PDU TO RACEWAY GND - SHORT
2	2371450-4	1.0	Yes	RUN#56A-CONSOLE TO RACEWAY GND-SHORT
3	2373436-3	1.0	Yes	RUN #102A - GANTRY TO CONSOLE LAN
4	2117848-7	1.0	Yes	20000 +/- 150 FIBER OPTIC CABLE
5	2371133-2	1.0	Yes	CONSOLE CABLES - SHORT KIT
6	2371134-2	1.0	Yes	GANTRY CABLES - SHORT KIT

Table 5-12 Short Cable Collector (2281840-5 BOM, rev 0)

2.2 Long Cable Collector (2281840-4)

FRU = 1

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2371450	1.0	Yes	RUN#55-PDU TO RACEWAY GND
2	2371450-3	1.0	Yes	RUN#56-CONSOLE TO RACEWAY
3	2373436-2	1.0	Yes	RUN #102 - GANTRY TO CONS
4	2117848-2	1.0	Yes	24384 +/- 150 FIBER OPTIC
5	2371133	1.0	Yes	CONSOLE CABLES - LONG KIT
6	2371134	1.0	Yes	GANTRY CABLES - LONG KIT

Table 5-13 Long Cable Collector (2281840-4 BOM, rev 0)

Section 3.0 — Counter Weight Kit

3.1 HP60 Gantry 78Weight Upgrade Kit (5179280)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2378668	2.0	-	P2033DT 78WEIGHT-L
2	2378669	5.0	-	P2033DU 78WEIGHT-M
3	2378670	2.0	-	P2033DV 78WEIGHT-S
12	46-328432P7	6.0	-	WASHER LOCK - SPRING 12.2 MM 21.1 MM
13	2236649-5	6.0	-	WASHER, FLAT, ISO7089,M12, ST-HB200
14	46-328425P13	2.0	-	NUT HEXAGON 12 MM
20	46-170684P3	0.0	-	LOCTITE 271, 50ML BOTTLE. 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER. (REV, TK, 1/92).
10	5128067	4.0	-	CAP SCREW U0480AA
11	5126559	2.0	-	STUD BOLT U0479AA

Table 5-14 Long Cable Collector (5179280 BOM, rev 1)



GE MEDICAL SYSTEMS

**GE MEDICAL SYSTEMS-AMERICAS: FAX 262.312.7434
3000 N. GRANDVIEW BLVD., WAUKESHA, WI 53188 U.S.A.**

**GE MEDICAL SYSTEMS-EUROPE: FAX 33.1.40.93.33.33
PARIS, FRANCE**

**GE MEDICAL SYSTEMS-ASIA: FAX 65.291.7006
SINGAPORE**