

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

**MI&CT ESERVICE NOTE
SN07101212b**

Class C information is Property of GE.

This short-term informational note is invalid 6 months after the issue date.

Date: October 17, 2007

Subject: Configuring Daylight Saving Time (2007) for IRIX O/S.

Effectivity: This Service Note only applies to sites in the United States and Canada.

The following LightSpeed CT Products (Qxi, Plus, Ultra, LS16) systems with the software revisions listed below are affected:. Effectivity for this applies to all LightSpeed 1X – 4X systems with **Sgi** computer based GOC consoles and CT systems running **IRIX – O/S**.

- | | | |
|-------------------|----------------|--------------------|
| • 405I.3_H4.0M5.1 | 405I.2_H4.0M5 | 308I.2_H3.1M5 |
| • 2.4.2_H2.4M5 | 101I.2_HSQXIM5 | 10.5_2.8.2I_H1.3M4 |

Note: This is not applicable to CTi, nor CT/PET Systems.

Details: Daylight Saving Time (DST) in the US and Canada has changed to meet the New Federal Code. Starting in 2007, the US Daylight Saving Time is observed from the second Sunday in March to the first Sunday in November, adding about a month to Daylight Saving Time. This Service Note adjusts DST parameters so that the time changes on the correct DST 2007 dates. **This procedure must be run every time 'reconfig' is run.**

Tool Requirements: No special tools required.

Manpower Requirements: One (1) Service Engineer for 30 minutes.

Resolution: Note: The initial verification steps can be done with Applications up.

1. Check your software version. Click *Unix Shell* under the *Tool Chest* in the upper right of the monitor and put the cursor in the Shell window.

Type: **check_config**

Press **[Enter]**

Note: If the result does not contain one of the software versions listed in the Effectivity section above, this service note does not apply.

2. Check your existing Time Zone parameters.

Type: **grep TZ /etc/TIMEZONE**

Press **[Enter]**

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

Note: The result should contain one of the following outputs in Tables 1-3

- o If your system has the original Time Zone values, output should match **Table 1**.
- o If your system already has **DST rules disabled**, output should match **Table 2**.
- o If your system has been changed to **Universal Time**, as listed in **Table 3**, you must run '**reconfig**' first and choose an Americas – Canada / US supported Time Zone, and reboot before executing this Service Note to adjust DST rules. If DST rules apply, one cannot use **Universal Time** parameters.

Table 1 below lists the original Time Zone values.

TZ=AST4ADT3	TZ=EST5EDT	TZ=CST6CDT	TZ=MST7MDT
TZ=PST8PDT	TZ=AKST9AKDT	TZ=HAST10HADT	

Table 2 below lists Time Zone values, DST rules disabled.

TZ=AST4	TZ=EST5	TZ=CST6	TZ=MST7
TZ=PST8	TZ=AKST9	TZ=HAST10	

- **Note your result from Tables 1 or 2 listed above, as needed in step 6 below.**

Table 3 below lists Universal Time Zone values, requiring 'reconfig'!

TZ=local+4:00	TZ=local+5:00	TZ=local+6:00	TZ=local+7:00
TZ=local+8:00	TZ=local+9:00	TZ=local+10:00	

- **If Universal Time is set, run 'reconfig' and choose Americas, and your Canada / US - Time Zone, then return to Step 2 of this Service Note.**

3. Check the Time and Zone info:

Type: **date** Press **[Enter]**

Output will be similar to this: Sun Oct 21 07:55:02 CST 2007

Note: Standard Time is indicated by: CST or EST, MST, PST, etc...
The following steps will adjust DST rules, so '**setdate**' is not needed. (Eliminates customer RFS 'Request For Service' call Oct. 28th, or Oct. 29th and again on Nov. 4th, or 5th)

4. If Applications are up, select the *Service Desktop* Icon, *Utilities*, and then *Application Shutdown*.

- a. Click *Unix Shell* under the *Tool Chest* in the upper right of the monitor and put the cursor in the Shell window. Enter the following in the telnet window.

Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

- b. Type: **su** - Press **[Enter]**
- c. Type the appropriate password. Press **[Enter]**
- d. Type: **configTime -p** Press **[Enter]**

Starting timezone configuration
Please enter correct timezone >

5. Choose the Time Zone line for your Canada/US - DST 2007 adjustments from chart listed below or skip to Step 6 if DST rules are not applicable to your area.

<u>Time Zone</u>	<u>Adjustment / Solution</u>	
For AST4ADT or AST4 (Atlantic)	Type: AST4ADT,M3.2.0,M11.1.0	Press [Enter]
For EST5EDT or EST5 (Eastern)	Type: EST5EDT,M3.2.0,M11.1.0	Press [Enter]
For CST6CDT or CST6 (Central)	Type: CST6CDT,M3.2.0,M11.1.0	Press [Enter]
For MST7MDT or MST7 (Mountain)	Type: MST7MDT,M3.2.0,M11.1.0	Press [Enter]
For PST8PDT or PST8 (Pacific)	Type: PST8PDT,M3.2.0,M11.1.0	Press [Enter]
For AKST9AKDT or AKST9 (Alaska)	Type: AKST9AKDT,M3.2.0,M11.1.0	Press [Enter]

Note: Reading the characters one would say for CST6CDT: Type all in Upper Case letters: 'CST', number 'SIX', 'CDT', put a **COMMA**, 'M', number 'THREE', the **period**, number 'TWO', the **period**, number 'ZERO', another **COMMA**, then, 'M', number 'ELEVEN', the **period**, number 'ONE', the **period**, number 'ZERO'.

6. **Some areas do not use Daylight Saving Time**, For the U.S. and its territories, Daylight Saving Time is NOT observed in Hawaii, American Samoa, Guam, Puerto Rico, the Virgin Islands, and most of Arizona. Choose your Time Zone line, and **disable all DST rules** to only use Standard Time for that area.

<u>Time Zone</u>	<u>DISABLE all DST rules</u>
For AST4ADT (Atlantic Time Zone)	Type: AST4 Press [Enter]
For EST5EDT (Eastern Time Zone)	Type: EST5 Press [Enter]
For CST6CDT (Central Time Zone)	Type: CST6 Press [Enter]
For MST7MDT (Mountain Time Zone)	Type: MST7 Press [Enter]
For PST8PDT (Pacific Time Zone)	Type: PST8 Press [Enter]
For AKST9AKDT (Alaska Time Zone)	Type: AKST9 Press [Enter]
For HAST10HADT (Hawaii Time Zone)	Type: HAST10 Press [Enter]

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

7. Check Adjustment or Disable results:

Type: **grep TZ /etc/TIMEZONE** Press **[Enter]**

One of the following should be seen in the output result:

Table 4 below lists the adjusted DST 2007 Time Zone values.

TZ=AST4ADT,M3.2.0,M11.1.0	TZ=EST5EDT,M3.2.0,M11.1.0	TZ=CST6CDT,M3.2.0,M11.1.0
TZ=MST7MDT,M3.2.0,M11.1.0	TZ=PST8PDT,M3.2.0,M11.1.0	TZ=AKST9AKDT,M3.2.0,M11.1.0

- Refer to Step 2, Table 2, page 2 result, if DST rules where disabled.

8. If console contains a **Sgi-O2** sbc, ie: LightSpeed 1x, you must copy /etc/TIMEZONE to the sbc, use the following command line:

Type: **rcp -p /etc/TIMEZONE root@sbc:/etc** Press **[Enter]**

9. Check Time and Zone info with the '**date**' command:

Type: **date** Press **[Enter]**

Output will be similar to this: Mon Oct 22 07:55:22 CDT 2007

Note: Standard Time is indicated by: CST or EST, MST, PST, etc...

- If Daylight Saving Time is indicated ie: CDT at this point, AND date/time is correct, you are done! The changes made to correct DST will become effective at the next reboot!
- If the time and date are not correct in Step 9, then continue with 'rebootsys' in step 10, and correct the time with 'setdate'!

10. Put the cursor in the Shell window. If done remotely, enter the following in the telnet window.

Type: **rebootsys** Press **[Enter]**

11. Check displayed time after the reboot with the "**date**" command:

Type: **date** Press **[Enter]**

Output will be similar to this: Mon Oct 22 07:55:22 CDT 2007

12. You will have to perform the '**setdate**' command, if the time displayed is incorrect.
- a. If Applications are up, select the Service Desktop Icon, *Utilities*, and then *Application Shutdown*.

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.

Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

- b. Click *Unix Shell* under the *Tool Chest* in the upper right of the monitor and put the cursor in the Shell window. If done remotely, enter the following in the telnet window.
- c. Type: **su -** Press **[Enter]**
- d. Type the appropriate password. Press **[Enter]**
- e. Type: **setdate** Press **[Enter]**

Press **[Enter]** if the displayed field is correct, or type in the appropriate date/time parameters.

13. To re-check, and verify all Time Zone settings:

Type: **grep TZ /usr/g/config/INFO /etc/TIMEZONE; date** Press **[Enter]**

Output will be similar to below:

```
/usr/g/config/INFO:setenv MENUTZ "Americas"  
/usr/g/config/INFO:setenv SUNTZ "US/Central"  
/usr/g/config/INFO:setenv TZ "CST6CDT US"  
/usr/g/config/INFO:setenv TZTEXT "United States/Central"  
/etc/TIMEZONE:TZ=CST6CDT,M3.2.0,M11.1.0  
Wed Sep 26 20:59:41 CDT 2007
```

Note: When running '**reconfig**' changes are made to the file named: /etc/TIMEZONE, therefore you must execute this Service Note, every time '**reconfig**' is performed. If Service GUI - '**reconfig**' is run, delay the system reboot by clicking **NO** on the exit box, and return to beginning of this service note.

Additionally, whenever '**reconfig**' is performed, if site is Insite enabled, the '**ConfigLink**' command must also be done to add the OLC hosts. To run the '**ConfigLink**' command, do the following:

Type: **cd ~insite; ConfigLink** Press **[Enter]**

Affected Service Publication: Qxi/CTi

PQR/SPR #: PQR 13098836

Kirby Louque
OnLine Engineer

MI&CT Engineering Manager Reviewer: _____ Mark Nyman _____