

SIEMENS

MAGNETOM

MR

Diagrams

System K2200

Harmony *syngo*MR
Symphony *syngo*MR
Sonata *syngo*MR
Upgrades to *syngo*MR:
from Impact/Expert
from Vision

© Siemens 1996
- For internal use only - All documents may only be used by authorized personnel for rendering services on Siemens Healthcare Products. Any document in electronic form may be printed once. Copy and distribution of electronic documents and hardcopies is prohibited. Offenders will be liable for damages. All other rights are reserved.

Copyright

"© Siemens, 2010" refers to the copyright of a Siemens entity such as Siemens Aktiengesellschaft - Germany, Siemens Shenzhen Magnetic Resonance Ltd. - China, Siemens Shanghai Medical Equipment Ltd. - China, Siemens Medical Solutions USA Inc. - USA and/or Siemens Healthcare Diagnostics Inc. - USA.

Document Version

Siemens reserves the right to change its products and services at any time.

In addition, manuals are subject to change without notice. The hard-copy documents correspond to the version at the time of system delivery and/or printout. Versions to hardcopy documentation are not automatically distributed. Please contact your local Siemens office to order current version or refer to our website <http://www.healthcare.siemens.com>.

Disclaimer

Siemens provides this documentation "as is" without the assumption of any liability under any theory of law.

The service of equipment described herein is to be performed by qualified personnel who are employed by Siemens or one of its affiliates or who are otherwise authorized by Siemens or one of its affiliates to provide such services.

Assemblers and other persons who are not employed by or otherwise directly affiliated with or authorized by Siemens or one of its affiliates are not entitled to use this documentation without prior written authority.

Table of Contents

The Table of Contents is found in the Bookmarks tab card which is shown/hidden with key **F5**.

Links

The Blue boxes found throughout the diagrams are hyperlinks.

Useful Acrobat Reader Shortcuts


There are a few keyboard shortcuts which will make using this PDF file much easier.


Alt + left arrow key go back to previous page/view


Alt + "+" zoom in


Alt + "-" zoom out


Ctrl + 0 fit page in window (useful if you zoomed)



























Goto
NEW

**Goto
OLD**




**to
Page 2**





CCS Plumbing Diagram

to CCS Wiring Diagram



Legend


- 1 Heat exchanger
- 2 blowers M1 - M3
- 3 3-way valve
- 4 Drain valve
- 5 Temperature sensor
- 6 3/4" connector
- 7 Strainer
- 8 Shut-off valve
- 9 Flow glass

R1 Thermal sensor



CCS Wiring Diagram

to CCS Plumbing Diagram



X800



X2




T1



X3

X4




R1











Y1






Block
Diagram

GPA K2209

Board	Settings	
D1 K2209	X16	2-3
	S1	both DOWN
	S2	both DOWN
	S3	7
	S4	2
D2 K2209	S1	ON 
	S2	ON 
D2 K2209	S1	ON 
	S2	ON 
D2 K2209	S1	ON 
	S2	ON 
D3 K2209	X3	3A - 3B 6A - 6B 9A - 9B
D4 K2209	S1	both DOWN
	S2	both DOWN
	S3	ON 
	S4	ON 
D5 K2209	S1	ON 
	S2	ON 
D6 K2209	X601	1 - 2

Board	Settings	
D10 K2209 Driver Supply	X4	1 - 2

Patient Table

Board	Settings	
D9 K2209	S1	ON 

Gradient Coil Temp Monitor

Block
DiagramFront View Description
GPA K2209

D5 - Turbomodulator (only for Upgrade)			
REF			Description
V2		red	Overtemperature Powerstage X2, Y2 or Z2
		red	Driver Error X2, Y2 or Z2
V3		red	Spare
		red	Overvoltage Powerstage X2
V4		red	Overvoltage Powerstage Y2
		red	Overvoltage Powerstage Z2
V5		red	not used
		red	not used
V6		red	not used
		red	not used
V7		red	not used
		red	not used


D6 - CAN/Service Board			
REF			Description
S401			Select gain, filter
S402			RESCON ON
S403			RESCON OFF
S404			MANUAL PROCESSOR RESET
V403		green	LOADWARE ACTIV
		green	FIRMWARE ACTIV
V407		green	RESCON ON
		green	RESCON ON
V408		yellow	ext. Channel on
		yellow	low pass filter on
V412		yellow	gain 1
		yellow	gain 2
		yellow	gain 5
		yellow	gain 10
		yellow	gain 10
V413		red	EP1 FAIL
		red	COILERROR
V607		red	POWERFAIL
X301			AUX IN
X401			CAN
X402			CAN
X500			N_LOOP OUT
X501			LOOP OUT
X502			SCOPE

D1 - DAC Board		
REF		Description
R153		OFFSET X
R177		OFFSET Y
R200		OFFSET Z
S1		DAC X OFF
		DAC X ON
		DAC Y OFF
		DAC Y ON
S2		DAC Z OFF
		DAC Z ON
		no function
		no function
X3		input: waveform - data and clock
X4		Current Nominal Value X-Axis
X5		Current Nominal Value Y-Axis
X6		Current Nominal Value Z-Axis
X81		Current Nominal Value X-Axis
X82		Current Nominal Value Y-Axis
X83		Current Nominal Value Z-Axis
X84		no function
X85		DAC out X-Axis
X86		DAC out Y-Axis
X87		DAC out Z-Axis
X88		no function
X9		ANA_GND

D2 - Regulator		
REF		Description
R191		Current Sensor Offset
V24		red CAV ERROR
V27		yellow Regulator Loop Error
X3		CAV (Current Actual Value) - (<40A/V)
X5		GND
X6		GND
X8		ANA_GND
X9		ANA_GND
X401		CX/Y/Z)1-Switch 1
X402		CX/Y/Z)1-Switch 2
X403		CX/Y/Z)1-Switch 3
X404		CX/Y/Z)1-Switch 4
X405		CX/Y/Z)2-Switch 1
X406		CX/Y/Z)2-Switch 2
X407		CX/Y/Z)2-Switch 3
X408		CX/Y/Z)2-Switch 4
X409		CX/Y/Z)3-Switch 1
X410		CX/Y/Z)3-Switch 2
X411		CX/Y/Z)3-Switch 3
X412		CX/Y/Z)3-Switch 4
X701		Regulator error signal
X702		not connected
X703		Regulator output
X704		not connected
X705		CAV (Current Actual Value) - (<40A/V)
X706		not connected
X707		Current Sensor 1
X708		not connected
X709		Current Sensor 2
X710		not connected
X711		not connected
X712		input: Control-test

D3 - Synchronisation		
REF		Description
R66		Phase or free wheeling frequency
S1		RESET
V30		green Vcc o.k.
V31		red Error N_15V
V32		red Error P_15V
V33		red Spare error 1
V34		red Spare error 2
V35		red Spare error 3
V36		red Spare error 4
V37		red ERROR SYNC
V38		red ERROR TRIA
X10		TRIA
X11		not connected
X12		CLK 50kHz
X13		not connected
X14		SYNC_CLK 50kHz
X15		not connected
X16		P_15V
X17		not connected
X18		N_15V
X19		not connected
X20		VCC
X21		not connected
X22		ANA_GND
X23		ANA_GND

D4 - Modulator		
REF		Description
S1		Modulator X - off
		Modulator X - on
		SERVICE MODE - on
S2		SERVICE MODE - off
		Modulator Y - off
		Modulator Y - on
S2		Modulator Z - off
		Modulator Z - on
V1		yellow Service Mode
		at least one Modulator disabled
V2		red Overtemperature Powerstage X1, Y1 or Z1
		red Driver Error X1, Y1 or Z1
V3		red Spare
		red Overvoltage Powerstage X1
V4		red Overvoltage Powerstage Y1
		red Overvoltage Powerstage Z1
V5		red Overcurrent X - Axis
		red Duty Cycle Error Powerstage X1
V6		red Overcurrent Y - Axis
		red Duty Cycle Error Powerstage Y1
V7		red Overcurrent Z - Axis
		red Duty Cycle Error Powerstage Z1
V8		yellow Current X positiv
		yellow Current X negativ
V9		yellow Current Y positiv
		yellow Current Y negativ
V10		yellow Current Z positiv
		yellow Current Z negativ



NUMARIS 4 Cascade Quantum

Board		Settings
D11	X16	2-3
	S1	both DOWN
	S2	both DOWN
	S3	3
	S6	6
D12	S1	ON
	S2	ON
	S3	ON
	S4	ON
	S5	ON
	S6	ON
	S7	ON
	S8	ON
	S9	ON
	X16	1 - 2
	X17	1 - 2
	X18	1 - 2
	X19	1 - 2
	X20	1 - 2
	X21	1 - 2
	X22	1 - 2
	X23	1 - 2
	X24	1 - 2
	X25	18 - 19 21 - 22 24 - 25

Board	Settings	
D14	X9	1 - 2
	X10	1 - 2
	X11	6 - 7 18 - 19
	X12	1 - 2
	X13	1 - 2
D16	X3	1 - 2
	X601	1 - 2
D17	X9	11 - 12 19 - 20 21 - 22
	S1	both DOWN
	S2	both DOWN

NUMARIS 4 Cascade Light

Board		Settings
D11	X16	2-3
	S1	both DOWN
	S2	both DOWN
	S3	3
	S6	6
D12	S1	ON
	S2	ON
	S3	ON
	S4	ON
	S5	ON
	S6	ON
	S7	ON
	S8	ON
	S9	ON
	X16	1 - 2
	X17	1 - 2
	X18	1 - 2
	X19	1 - 2
	X20	1 - 2
	X21	1 - 2
	X22	1 - 2
	X23	1 - 2
	X24	1 - 2
	X25	18 - 19 21 - 22 24 - 25

Board		Settings
D14	X9	1 - 2
	X10	1 - 2
	X11	6 - 7 18 - 19
	X12	1 - 2
	X13	1 - 2
D16	X3	1 - 2
	X601	1 - 2
D17	X9	11 - 12 15 - 16 19 - 20 21 - 22
	S1	both DOWN
	S2	both DOWN

ME03
1000V - 300ABlock
Diagram

E5 K2217 : GSSU -L ight			
Connected with	Printed Board	Jumper	
D18 K2217 X11 , X12	D111 K2217	X16	2 - 3
		S1	both down
		S2	both down
		S3	3
		S4	6
D18 K2217 X21 , X22	D124 K2217	S1	
		S2	
		S3	
		S4	
		S5	
		S6	
		S7	
		S8	
		S9	
		X16	1 - 2
		X17	1 - 2
		X18	1 - 2
		X19	1 - 2
		X20	1 - 2
		X21	1 - 2
		X24	1 - 2
		X22	without Jumper
		X23	without Jumper
		X25	18 - 19
		X25	21 - 22
X25	24 - 25		
D18 K2217 X31 , X32	D142 K2217	X9	without Jumper
		X10	2 - 3
		X11	6 - 7
		X11	18 - 19
		X12	1 - 2
		X13	1 - 2
D18 K2217 X41 , X42	D16 K2217	X3	2 - 3
		X601	1 - 2
D18 K2217 X51 , X52	not used		
D18 K2217 X61 , X62	D17 K2217	X9	11 - 12
		X9	15 - 16
		X9	17 - 18
		X9	19 - 20
		X9	21 - 22
		S1	both down
		S2	both down

Jumper Settings GPA K2217 for NUMARIS 3.5

Block
Diagram

Sonata / 2000V

E5 K2217 : GSSU			
Connected with	Printed board	Jumper	
D18 K2217 X11 , X12	D11 K2217	X16	2 - 3
		S1	both down
		S2	both down
		S3	3
		S4	6
D18 K2217 X21 , X22	D12 K2217	S1	
		S2	
		S3	
		S4	
		S5	
		S6	
		S7	
		S8	
		S9	
D18 K2217 X31 , X32	D14 K2217	X16	1 - 2
		X17	1 - 2
		X18	1 - 2
		X19	1 - 2
		X20	1 - 2
		X21	1 - 2
		X24	1 - 2
		X22	2 - 3
		X23	2 - 3
		X25	18 - 19
D18 K2217 X41 , X42	D16 K2217	X9	2 - 3
		X10	2 - 3
		X11	6 - 7
		X12	18 - 19
D18 K2217 X51 , X52	D13 K2217	X100	1 - 2
		X101	1 - 2
		X102	1 - 2
		X103	1 - 2
		X104	1 - 2
		X105	1 - 2
		X106	2 - 3
		X107	1 - 2
		X108	1 - 2
		X109	1 - 2
D18 K2217 X61 , X62	D17 K2217	X9	11 - 12
		X9	19 - 20
		X9	21 - 22
		S1	both down
		S2	both down

039T / 2000V

E5 K2217 : GSSU			
Connected with	Printed board	Jumper	
D18 K2217 X11 , X12	D11 K2217	X16	2 - 3
		S1	both down
		S2	both down
		S3	3
		S4	6
D18 K2217 X21 , X22	D12 K2217	S1	
		S2	
		S3	
		S4	
		S5	
		S6	
		S7	
		S8	
		S9	
D18 K2217 X31 , X32	D14 K2217	X16	1 - 2
		X17	1 - 2
		X18	1 - 2
		X19	1 - 2
		X20	1 - 2
		X21	1 - 2
		X24	1 - 2
		X22	2 - 3
		X23	without Jumper
		X25	18 - 19
D18 K2217 X41 , X42	D16 K2217	X9	2 - 3
		X10	2 - 3
		X11	6 - 7
		X12	18 - 19
D18 K2217 X51 , X52	D13 K2217	X100	1 - 2
		X101	1 - 2
		X102	1 - 2
		X103	1 - 2
		X104	1 - 2
		X105	1 - 2
		X106	2 - 3
		X107	1 - 2
		X108	1 - 2
		X109	1 - 2
D18 K2217 X61 , X62	D17 K2217	X9	11 - 12
		X9	19 - 20
		X9	21 - 22
		S1	both down
		S2	both down

039T / 1000V

E5 K2217 : GSSU			
Connected with	Printed board	Jumper	
D18 K2217 X11 , X12	D11 K2217	X16	2 - 3
		S1	both down
		S2	both down
		S3	3
		S4	6
D18 K2217 X21 , X22	D12 K2217	S1	
		S2	
		S3	
		S4	
		S5	
		S6	
		S7	
		S8	
		S9	
D18 K2217 X31 , X32	D14 K2217	X16	1 - 2
		X17	1 - 2
		X18	1 - 2
		X19	1 - 2
		X20	1 - 2
		X21	1 - 2
		X24	1 - 2
		X22	2 - 3
		X23	without Jumper
		X25	18 - 19
D18 K2217 X41 , X42	D16 K2217	X9	without Jumper
		X10	2 - 3
		X11	6 - 7
		X12	18 - 19
D18 K2217 X51 , X52	not used	X100	1 - 2
		X101	1 - 2
		X102	1 - 2
		X103	1 - 2
		X104	1 - 2
		X105	1 - 2
		X106	2 - 3
		X107	1 - 2
		X108	1 - 2
		X109	1 - 2
D18 K2217 X61 , X62	D17 K2217	X9	11 - 12
		X9	15 - 16
		X9	19 - 20
		X9	21 - 22
		S1	both down
		S2	both down

ME03 / 2000V / 300A

E5 K2217 : GSSU			
Connected with	Printed board	Jumper	
D18 K2217 X11 , X12	D11 K2217	X16	2 - 3
		S1	both down
		S2	both down
		S3	3
		S4	6
D18 K2217 X21 , X22	D12 K2217	S1	
		S2	
		S3	
		S4	
		S5	
		S6	
		S7	
		S8	
		S9	
D18 K2217 X31 , X32	D14 K2217	X16	1 - 2
		X17	1 - 2
		X18	1 - 2
		X19	1 - 2
		X20	1 - 2
		X21	1 - 2
		X24	1 - 2
		X22	without Jumper
		X23	without Jumper
		X25	18 - 19
D18 K2217 X41 , X42	D16 K2217	X9	2 - 3
		X10	2 - 3
		X11	6 - 7
		X12	18 - 19
D18 K2217 X51 , X52	not used	X100	1 - 2
		X101	1 - 2
		X102	1 - 2
		X103	1 - 2
		X104	1 - 2
		X105	1 - 2
		X106	2 - 3
		X107	1 - 2
		X108	1 - 2
		X109	1 - 2
D18 K2217 X61 , X62	D17 K2217	X9	11 - 12
		X9	13 - 14
		X9	15 - 16
		X9	17 - 18
		S1	both down
		S2	both down

Allegra/ 2000V

E5 K2217 : GSSU			
Connected with	Printed board	Jumper	
D18 K2217 X11 , X12	D11 K2217	X16	2 - 3
		S1	both down
		S2	both down
		S3	3
		S4	6
D18 K2217 X21 , X22	D12 K2217	S1	
		S2	
		S3	
		S4	
		S5	
		S6	
		S7	
		S8	
		S9	
D18 K2217 X31 , X32	D14 K2217	X16	1 - 2
		X17	1 - 2
		X18	1 - 2
		X19	1 - 2
		X20	1 - 2
		X21	1 - 2
		X24	1 - 2
		X22	without Jumper
		X23	2 - 3
		X25	18 - 19
D18 K2217 X41 , X42	D16 K2217	X9	2 - 3
		X10	2 - 3
		X11	6 - 7
		X12	18 - 19
D18 K2217 X51 , X52	D13 K2217	X100	1 - 2
		X101	1 - 2
		X102	1 - 2
		X103	1 - 2
		X104	1 - 2
		X105	1 - 2
		X106	2 - 3
		X107	1 - 2
		X108	1 - 2
		X109	1 - 2
D18 K2217 X61 , X62	D17 K2217	X9	11 - 12
		X9	19 - 20
		X9	21 - 22
		S1	both down
		S2	both down

GPA K2217


Block Diagram

only on NUM3.5 systems
with Quantum gradients

Front View Description


GPA K2217

[illegible]



Sonata/Trio/Allegra

R = flow reducer
R* = 1st production
run cabinets
R = Schroff cabinets



For internal use only. Copyright Siemens, 2010

RCA Cabinet Wiring Leybold Compressor

SIEMENS AG
Healthcare Sector



RCA - Cabinet Wiring Diagram

Project
024/047

Print No.
M1-010.844.01.06

Date
06.2011
Page
27

**Goto
APD**



RCA Cabinet Wiring APD Compressor

SIEMENS AG
Healthcare Sector

RCA - Cabinet Wiring Diagram


Project
024/047

Print No.
M1-010.844.01.06


Date
06.2011

Page
28

**Goto
Leybold**



Top Panel Layout (Schroff)



Secondary Water Cooler Block Diagram

power ---
coax ---
signal ---
fiber optic -.-

SIEMENS AG
Healthcare Sector


RCA - Cabinet Wiring Diagram


Project
024/047


Print No.
M1-010.844.01.06


Date
06.2011


Page
30





















Breast Array



Head Array




Spine Array



Neck Array Body Array


X5




A1 is connected to the RX3 and RX3_Gnd signals of X3. The PIN3 and PIN2_Gnd1 signals are superimposed on A1 as well.


A2 is connected to the RX4 and RX4_Gnd signals of X3 and also includes the +10V supply.

Extremity Coil




X2







Board Layout




**Wiring
Diagram**



X100 Mains input
X2 Patient Fan
X3 Vertical Motor
X4 transformer secondary
X5 connection to PTAB
X51 Freq Converter
X52 Freq Converter






<W008>

Micro Out
Intercom

Micro In
HOST Audio Port (pink)




Length	Part Number	Notes
3m	75 46 125 K2201	

<W009>

Speaker Out
HOST Sound Card (green)


CV in
Intercom



Length	Part Number	Notes
3m	75 46 133 K2201	


X50-X57
CCA roof

X01-X08
Filter panel



CCA roof

TAS




Length	Part Number	Notes
1.5m	47 59 713 K2201	RG214 cable
5m	57 51 743 K2201	
10m	47 66 056 K2201	LCF 7/8" cable
13m	47 60 695 K2201	LCF 7/8"


<W020>

CCA roof

SAMI





Length	Part Number	Notes
1,5m	71 21 689 K2201	
5m	75 44 630 K2201	
10m	75 44 484 K2201	
13m	75 44 450 K2201	



<W026a>

Port 7
Serial Distrib.

X5
Alarm Box



D-sub, 25 pin, receptacle


D-sub, 9 pin, receptacle

Length	Part Number	Notes
15m	57 53 350 K2200	

<W027a>

Video Out
Host Video Card


Monitor




Length	Part Number	Notes
3m	57 53 343 K2220	5-conductor 75 Ohm coax cable


(AEM-Nr.: 105 755 REV:01)


<W029a>



Length	Part Number	Notes
4m	57 53 327 K2220	







A

B

C

D

E

F

G

H

A

B

C

D

E

F


G

H

<W033a>

MMC
RX4

Imager
Slot 7 (D35 K5004)




Length	Part Number	Notes
25m	57 53 319 K2200	


<W033b>


MMC
RX4

Imager
Slot 1 (D35 K5004)



Length	Part Number	Notes
25m	71 22 794 K2200	








<W044a>

X31
CCA roof

2X
Switch (MRC)




D-sub, 9 pin, plug



RJ 45, 8 pin


Length	Part Number	Notes
23m	57 53 301 K2200	ETHERNET Class 5 + performance type

<W047>



Length	Part Number	Notes
23m	47 55 133 K2203	

<W048a>





Length	Part Number	Notes
4m	57 53 293 K2220	

<W057>

X23
CCA roof

S4
RF Door contact




D-sub, 9 pin,
receptacle

Length	Part Number	Notes
20m	47 62 030 K2201	

X11
Intercom

<W058>

D15-4 / X1
Filter Panel




Length	Part Number	Notes
25m	47 62 048 K2201	

<W061>

X103
CCA roof

UPS




Length	Part Number	Notes
9m	57 50 273 K2220	


<W062>

X104
CCA roof

X101
RCA




Length	Part Number	Notes
15m	47 59 564 K2200	




<W064>

X107
CCA roof

Printer



Length	Part Number	Notes
23m	47 62 881 K2203	




Length	Part Number	Notes
23m	47 64 077 K2201	

<W066a>

X109
CCA roof

Line Distributor for MRC




Length	Part Number	Notes
23m	57 53 277 K2200	

<W068>

X105
CCA roof

Z4, Z5
Filter Panel



Main power to
Patient Table

Vacuum Pump


AMP - Metrimate, 6 pin,
Receptacle

Length	Part Number	Notes
2m	47 59 598 K2200	
13m	47 59 606 K2200	

<W069>


X110
CCA roof

MRIR (Imager)



Length	Part Number	Notes
23m	57 53 269 K2200	

<W072a>




Length	Part Number	Notes
1.5m	71 21 002 K2200	

<W076a>

MRC
Power Distibutor

Host



IEC 320 C13
plug


IEC 320 C14
socket

Length	Part Number	Notes
3m	57 53 244 K2200	

<W077a>

X111
Mains Box

X100
GPA




Length	Part Number	Notes
2m	57 53 426 K2200	

X01-X08
Filter Panel

7m	47	59	895	K2201
10m	47	66	296	K2201
15m	57	51	636	K2201

X04
N-type

X08
N-type




X04
QLA-type

X08
QLA-type


<W118>

Filter Panel

LCCS



Length	Part Number	Notes
7m	47 59 754 K2201	
10m	47 66 098 K2201	
15m	57 51 784 K2201	



SIEMENS AG
Healthcare Sector

Cable W119

Project
024/047

Print No.
M1-010.844.01.06

Date
06.2011

Page
80

<W119>

TAS

BCCS

X4
7/16-type

SSB1_4

X1
7/16-type

Length	Part Number	Notes
--------	-------------	-------

7m	47 59 762 K2201	
----	-----------------	--

10m	47 66 106 K2201	
-----	-----------------	--

15m	57 51 792 K2201	
-----	-----------------	--

1

2

3

4

Δ


5

6

Δ


7

8



QLA-type

Length	Part Number	Notes
7m	47 59 770 K2201	
10m	47 66 114 K2201	
15m	57 51 800 K2201	



A

B

C

D

E

F

G

H

A

B

C


D

E


F


G


H



Length	Part Number	Notes
0.5m	71 21 697 K2201	







A

B

C

D

E

F

G

H

A

B

C

D

E

F


G

H

<W143>

X35
LCCS

X7
TAS




Length	Part Number	Notes
7m	47 61 529 K2204	
10m	47 66 189 K2204	
15m	57 51 875 K2201	

<W145>


Filter Panel

X23

Filter Panel




Length	Part Number	Notes
	47 55 463 K2204	



<W149>

X1
Camera


X56A
Filter Panel



Sühner - Type 11
BNT-50-3-1c

Sühner - Type 11
BNT-50-3-1c

Length	Part Number	Notes
		10m Sühner Triax cable G03333




Length	Part Number	Notes
0,1m	47 55 455 K2204	

<W151>

X101
Filter Panel

Filters




Length	Part Number	Notes
0,1m	47 55 604 K2204	


<W152>

Filters

X101A
Filter Panel




Length	Part Number	Notes
0,1m	47 55 612 K2204	




Length	Part Number	Notes
7m	47 60 547 K2200	
12m	47 66 197 K2200	
15m	57 51 842 K2201	

For internal use only. Copyright Siemens, 2010



For internal use only. Copyright Siemens, 2010

For internal use only. Copyright Siemens, 2010




<W200>
618-430

X3
CCA roof


X12
Filter Panel

710-920 WIRING DETAIL

PIN NO.	COLOUR
1	BK
2	RD
3	WH
4	GN
5	OR
6	BU
7	BN
8	YL
9	VI
10	GY



Thomas&Betts,
Triad 20 circular, 14 pin



Thomas&Betts,
Triad 20 circular, 14 pin


Length	Part Number	Notes
2m	101 04 382 K2200	
13m	72 78 323 K2200	
15m		

<W201>
618-431

X12_A
Filter Panel


X1
Gradient Coil

Shim Filter




Thomas&Betts,
Triad 20 circular, 14 pin plug

AMP 206150-1, receptacle



710-647 WIRING DETAIL

PIN	COLOUR
14	BLACK
1	RED
10	WHITE
26	GREEN
24	ORANGE
3	BLUE
12	BROWN
28	YELLOW
35	VIOLET
37	GREY




711-197 WIRING DETAIL

PIN	COLOUR
1	BLACK
2	RED
3	WHITE
4	GREEN
5	ORANGE
6	BLUE
7	BROWN
8	YELLOW
9	VIOLET
10	GREY

Length	Part Number	Notes
10m	86 68 654 K2261	

D15-3 X1
Filter Panel



D-sub, 25 pin
plug

Length	Part Number
2m	101 17 063 K2200
13m	101 04 387 K2200
15m	

<W205>

Side Connector D


D15-3 X2
Filter Panel

Magnet
A

B

C

ERDU
Button 1




Length	Part Number	Notes
10m	57 62 492 K2200	

<W206>
618-448

X1
CCA Roof

X1
Alarm Box




Length	Part Number	Notes
20m	70 91 916 K2200	

<W210>
618-975

X102
RCA Roof

X101
Filter Panel




Length	Part Number	Notes
10m	101 04 295 K2200	

<W211>

X101A
Filter Panel

X100
Cold Head




Length	Part Number	Notes
15m	101 04 566 K2257	old SMT pn: 618-975

<W552>


X2
Gradient Coil


X20
PTAB



Length	Part Number	Notes
	81 10 053 K2257	

										SIEMENS AG Healthcare Sector		Cable List		Project 024/047		Print No. M1-010.844.01.06		Date 06.2011		Page 106	
A	W001-	1	11 40 842	K2200	SIG	Grad_X,Y,Z	W110	7	47 59 739	K2201	COAX	MRSIG_1									
	W003	10	11 41 436	K2200				10	47 66 072	K2201											
	K2217	13	11 40 859	K2200	1	59 39 520	W111	7	47 59 747	K2201	COAX	MRSIG_2									
			for Sonata only						10	47 66 080	K2201										
		10	59 39 314	K2217		for Sonata only	W112	7	47 62 139	K2201	COAX	MRSIG_3									
		13	57 27 420	K2217		for Sonata only		10	47 66 288	K2201											
	W008	3	75 46 125	K2201	SIG	Intercom to Host (Mic)	W113	7	47 59 895	K2201	COAX	MRSIG_4									
	W009	3	75 46 133	K2201	SIG	Host to Intercom (CV_In)		10	47 66 296	K2201											
	W010	1.5	47 59 689	K2201	COAX	MRSIG 1	W118	7	47 59 754	K2201	COAX	TTX									
		10	47 66 023	K2201				10	47 66 098	K2201											
B		13	47 59 796	K2201			W119	7	47 59 762	K2201	COAX	SSB									
	W011	1.5	47 59 697	K2201	COAX	MRSIG 2		10	47 66 106	K2201											
		10	47 66 031	K2201			W120	7	47 62 410	K2200	COAX	TX_LC									
		13	47 59 804	K2201				10	47 66 213	K2200											
	W012	1.5	47 60 893	K2203	COAX	MRSIG 3	W121	7	47 59 770	K2201	COAX	CV									
		10	47 66 262	K2203				10	47 66 114	K2201											
		13	47 59 877	K2203			W123	7	47 59 788	K2201	COAX	Dummy									
	W013	1.5	47 62 147	K2203	COAX	MRSIG 4		10	47 66 247	K2201											
		10	47 66 270	K2203			W124	0.5	71 21 697	K2201	COAX	SSB2									
		13	47 60 885	K2203			W125	7	75 44 476	K2201	COAX	TX_LC									
C	W014	1.5	57 52 832	K2201	COAX	MRSIG 5		10	71 22 349	K2201											
		10	57 52 873	K2201				15	75 45 648	K2201											
	W015	1.5	57 52 295	K2201			W140	7	47 60 802	K2201	SIG	intercom to PTab									
		10	57 52 840	K2201	COAX	MRSIG 6		1.5	47 66 155	K2201											
		10	57 52 881	K2201			W142	7	47 62 055	K2204	PWR	P_N_15 to RFAS									
		13	57 52 303	K2201				10	47 66 163	K2204											
	W016	1.5	57 52 857	K2201	COAX	MRSIG 7	W143	7	47 61 529	K2204	SIG	TAS control									
		10	57 52 899	K2201				10	47 66 189	K2204											
	W017	1.5	57 52 311	K2201			W145	47 55 463	K2204	PWR											
		10	57 52 865	K2201	COAX	MRSIG 8	W146	7	47 59 812	K2201	SIG	Coil Control									
D	W018	1.5	57 52 907	K2201				10	47 66 171	K2201											
		13	57 52 329	K2201			W150	0.1	47 55 455	K2204	PWR	X21 to Z1-Z3									
	W019	1.5	47 59 705	K2201	COAX	TTX	W151	0.1	47 55 604	K2204	PWR	X101 to Z8-Z10									
		10	47 66 049	K2201			W152	0.1	47 55 612	K2204	PWR	Z8-Z10 to X101A									
		13	47 60 687	K2201			W160	7	47 60 547	K2200	PWR	AC to PTab									
	W020	1.5	47 59 713	K2201	COAX	SSB1_3		12	47 66 197	K2200											
		10	47 66 056	K2201			W182	7	47 60 596	K2200	GND	Filter Panel to Magnet									
		13	47 59 695	K2201			W183	0.5	47 60 604	K2200	GND	Filter Panel to Cabin									
	W021	5.0	71 21 689	K2201	COAX	SSB2_3	W200	2	101 04 382	K2200	SIG	Shim									
		13	75 44 630	K2201				13	72 78 323	K2200											
E	W022a	5	75 44 484	K2201			W201	10	86 68 654	K2261	SIG	Shim									
	W023a	15	75 44 450	K2201			W204	2	101 17 063	K2200	SIG	Magnet monitoring									
	W026a	15	47 60 711	K2201	COAX	CV		13	101 04 387	K2200											
	W027a	15	47 66 221	K2201			W205	10	57 62 492	K2200	SIG	Magnet monitoring									
	W029a	4	47 60 745	K2201			W206	20	70 91 916	K2200	SIG	CCA to Alarm Box									
	W030	12	57 53 384	K2200	SIG	LAN Host to switch	W210	10	101 04 295	K2200	PWR	Line to RCA Controller									
		15	57 53 350	K2200	SIG	Host to Alarm Box	W211	15	101 04 566	K2257	PWR	Cold Head									
		15	57 53 343	K2200	SIG	D-Sub-BNC video cable															
	W031	12	57 53 327	K2200	SIG	Host to Mouse	W541	0.2													
		15	47 60 711	K2201	FOC	PTAb, RF control	W542	0.2													
F	W032	22	47 66 221	K2201			W552	2	81 10 053	K2257	SIG	NTC temp sensor to Ptab									
		22	47 60 745	K2201																	
	W033a	25	47 60 729	K2201	FOC	Grad Coil Temp Error															
	W033b	25	47 66 239	K2201																	
	W039a	2	47 60 752	K2201																	
	W040	1.5	47 60 737	K2201	FOC	CAN - MSUP to RCA															
		10	57 53 319	K2201	FOC	Rawdata RX4 to MRIR															
		13	71 22 794	K2201	FOC	Rawdata RX4 to MRIR															
		2	57 53 459	K2200	SIG	BNC-BNC video cable															
		10	47 60 760	K2201	SIG	LC plug 1-2															
G	W041	1.5	47 66 122	K2201																	
		13	47 60 810	K2201																	
		10	47 60 778	K2201	SIG	LC plug 3-4															
	W042	1.5	47 66 130	K2201																	
		13	47 62 063	K2201																	
		10	47 60 786	K2201	PWR	P_N_15V															
		13	47 66 148	K2201	PWR																
	W044a	23	47 59 820	K2201	PWR																
	W045	23	57 53 301	K2200	SIG	LAN - Host to switch (MRC)															
	W046	23	54 70 505	K2075	SIG	rawdata cable to SMI-5															
H	W047	23	54 70 513	K2075	SIG																
	W048a	4	47 55 133	K2203	SIG	Physio beeper to intercom															
	W051a	23	57 53 293	K2200	SIG	Host to Keyboard															
	W051b	23	57 53 285	K2200	SIG	LAN Switch (MRC) to MRIR															
	W057	20	71 22 786	K2200	SIG	LAN Switch (MRC) to MRIR															
	W058	25	47 62 030	K2201	SIG	door contact															
	W061	9	47 62 048	K2201	SIG	Intercom to filter panel															
	W062	15	57 50 273	K2200	PWR	UPS															
	W063	15	47 59 564	K2200	PWR	to RCA Controller															
	W064	23	47 59 572	K2200	PWR	3phase to RCA															





1		2		3		4 ▽		5		6 ▽		7		8					
Changes to previous version: All drawings converted to vector and text to allow search queries. Links have been added for better navigation. Block diagrams for Cooling system added Block diagrams for RFAS added										SIEMENS AG Healthcare Sector		Changes to previous version Project 024/047		Print No. M1-010.844.01.06		Date 06.2011		Page 109	