

GE Healthcare

LightSpeed 5.X Pro16

Illustrated Parts List

OPERATING DOCUMENTATION



2349187-100
Rev 18

Chapter 5

Miscellaneous Items

This chapter contains the following:

- Section 1.0 — Pro16 Ship Collector.page 368
 - 1.1 Base Cover Collector H4 (2366670)page 369
 - 1.2 H-Power Fuse Kit (2218570-2)page 370
 - 1.3 Balance Kit, Field, H-Power (2376682)page 371
 - 1.4 P2030BD Ship C HP-100/60 Common (2385413)page 372
 - 1.4.1 Short Cable Collector (2281840-5)page 374
 - 1.4.2 Long Cable Collector (2281840-4)page 374
- Section 2.0 — Small Room Catalog (5191743)page 375
 - 2.1 Dolly, Pivoting, Small Room, RH (5179183)page 376
 - 2.2 Dolly, Pivoting, Small Room, LH (5179183-2)page 378

Section 1.0 — Pro16 Ship Collector

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2366670	1.000	Yes	BASE COVER COLLECTOR H4 {For further detail, see Section 1.1 , beginning on page 369 .}
2	2218570-2	1.000	Yes	H-POWER FUSE KIT {For further detail, see Section 1.2 , beginning on page 370 .}
3	2376682	1.000	Yes	BALANCE KIT, FIELD, HPOWER {For further detail, see Section 1.3 , beginning on page 371 .}
4	2385413	1.000	No	P2030BD SHIP C HP-100/60 {For further detail, see Section 1.4 , beginning on page 372 .}

Table 5-1 Pro16 Ship Collector (Non-Mobile) (2373160 BOM, rev 4)

1.1 Base Cover Collector H4 (2366670)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
5	2333997	1.000	Yes	COVER, FRONT BASE - H4
6	2333996	1.000	Yes	COVER, BASE LEFT - H4
7	2333995	1.000	Yes	COVER, BASE RIGHT - H4
8	2275686	1.000	Yes	COVER ASSY, BASE REAR - L
9	2275685	1.000	Yes	COVER ASSY, BASE REAR - R
1	46-328417P3	4.000	No	SCREW HEXAGON SOCKET 4 MM
2	46-328430P2	25.000	No	WASHER PLAIN - NORMAL 4.3
3	46-328432P2	29.000	No	WASHER LOCK - SPRING 4.1
4	46-328425P1	5.000	No	NUT HEXAGON 4 MM
10	2277652	3.000	No	LATCH, H2 SIDE COVER COMP
11	2261321	2.000	No	BRACKET - FRONT BASE COVE
12	2261324	2.000	No	BRACKET - SIDE BASE COVER
13	2261325	1.000	No	BRACKET - REAR BASE COVER
14	2285469	2.000	No	BRACKET, BASE COVER MOUNTI
15	2288882	2.000	No	BULKHEAD - BLANK
16	46-328420P4	4.000	No	SCREW BUTTON HEAD 4 MM 16
17	2366670IDW	1.000	No	BASE COVER COLLECTOR H4
18	46-328417P4	17.000	No	SCREW HEXAGON SOCKET 4 MM

Table 5-2 Base Cover Collector H4 (2366670 BOM, rev 0)

1.2 H-Power Fuse Kit (2218570-2)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2351492	2.000	Yes	200A FUSE
2	2364059	2.000	Yes	GLASS FUSE
3	46-170021P50	2.000	Yes	FUSE 12 AMPS 250 VOLTS
4	46-170021P15	2.000	Yes	FUSE 8 AMPS 250 VOLTS
5	2336517-2	2.000	Yes	FUSE 25 AMPS 700 VOLTS 2.
6	46-170021P52	3.000	Yes	FUSE 3 AMPS 250 VOLTS
7	46-170021P10	2.000	Yes	006.000A 0250V 3AG
8	46-170021P101	1.000	Yes	FUSE 20 AMPS 700 VOLTS
9	46-170021P11	2.000	Yes	7.0A, 250V SLO-BLO FUSE.
10	46-170021P31	2.000	Yes	1/2A, 250V SLO-BLO FUSE.
11	2106993-5	2.000	Yes	FUSE 20 AMPS 500 VOLTS .3
12	46-327160P1	2.000	Yes	12.0A, 125VAC, DUAL ELEME
13	2379651	1.000	Yes	FUSE - 700 V, 200 A, FAST

Table 5-3 H-Power Fuse Kit (2218570-2 BOM, rev 1)

1.3 Balance Kit, Field, H-Power (2376682)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2316653	6.000	No	WEIGHT,107 DEG,6MM
2	2316653-2	6.000	No	WEIGHT,107 DEG,1.5MM
3	2317863	2.000	No	SHIM,180 DEG,3.2MM
4	2317865-2	2.000	No	SPACER,180 DEG,25.4MM
5	2317994	20.000	No	WEIGHT,TRIM
6	2319142	2.000	No	SPACER,107 DEG,25.4MM
7	2319145	2.000	No	SHIM,107 DEG,3.3MM
8	2342895-2	6.000	No	WEIGHT,180 DEG,6MM
9	2342895-3	6.000	No	WEIGHT,180 DEG,1.5MM
10	2370053	1.000	No	WEIGHT,TRIAL,NEGATIVE Z
11	2370227	1.000	No	WEIGHT,TRIAL,POSITIVE Z,
12	46-328425P4	4.000	No	NUT HEXAGON 16 MM
13	46-328417P70	2.000	No	SCREW HEXAGON SOCKET 12 M
14	2376683	1.000	No	SOCKET, 15/16" DEEP IMPAC
15	2377399	1.000	No	LEFT SIDE DAS COVER TOOL
16	2377400	1.000	No	RIGHT SIDE DAS COVER TOOL
17	2377401	1.000	No	CENTER DAS COVER TOOL
18	46-328432P5	4.000	No	WASHER LOCK - SPRING 16.2

Table 5-4 Balance Kit, Field, H-Power (2376682 BOM, rev 0)

1.4 P2030BD Ship C HP-100/60 Common (2385413)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
-	2281840-4	-	Yes	H-Power LONG CABLE COLLECTOR {For further detail, see Section 1.4.2 , beginning on page 374 .}
25	2282011	2.000	Yes	NON-TILT DOLLY REAR COVER
26	5101127	1.000	Yes	DOLLY PIVOTING , RIGHT HA
30	5101127-2	1.000	Yes	DOLLY, PIVOTING, LEFT HAN
31	2282509	1.000	Yes	HOIST-2-PC, ARM, HELIOS 2
32	2282510	1.000	Yes	HOIST-2-PC, LONG TUBE, HE
33	2277548	4.000	Yes	PLUG, SIDE COVER H2
40	2364629	1.000	Yes	LABEL INSTALLATION KIT
41	2281840-5	1.000	Yes	H-POWER SHORT CABLE COLLECTOR {For further detail, see Section 1.4.1 , beginning on page 374 .}
42	2336557	1.000	Yes	SEISMIC KIT FOR GLOBAL CO
26A	2289712	ALT	Yes	DOLLY, PIVOTING - RH
30A	2289712-2	ALT	Yes	DOLLY, PIVOTING - _LH
1	46-229452G1	1.000	No	CT-9800 (ACCESSORIES)
2	46-297482G1	1.000	No	IV BOARD ACC FOR ZEUS TAB
4	46-297592P1	1.000	No	BAG HOLDER
6	2105873	10.000	No	.531 I.D., 1.625 O.D., .2
7	2106573	10.000	No	HILTI ANCHOR STUD 1/2 DIA
8	2115984	1.000	No	PAD FOR HSA METALLESS FOO
9	2115994	1.000	No	METALLESS CRADLE PAD & ST
12	46-222436G1	1.000	No	CT-9800 - CHIN STRAP
13	46-297960P1	1.000	No	FOAM INSERT
14	2319891-3	1.000	No	INSTALLATION SUPPORT KIT
18	46-222436G1	3.000	No	CT-9800 - CHIN STRAP
19	46-237412P1	1.000	No	STRAP

Table 5-5 P2030BD Ship C HP-100/60 Common (2385413 BOM, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
20	46-222443P1	1.000	No	OPEN CELL POLYESTER FOAM
21	46-222444P1	1.000	No	CLOSED CELL LINEAR ETHYLE
22	46-222445P1	1.000	No	CLOSED CELL LINEAR ETHYLE
23	46-222921P1	6.000	No	OBSOLETE FOR MR SITES PER
24	46-199388P1	1.000	No	HEAD HOLDER POLYURETHANE
29	2363851	1.000	No	TEMPLATE, INSTALLATION HP
34	2217173-2	1.000	No	STOOL WITH FOOT REST AND
35	2299113	1.000	No	CT SERVICE CABINET
37	2309994	1.000	No	IV POLE ASSY
38	2329064	1.000	No	TABLE ACC ASSY
39	2355100	1.000	No	TRAY, H2 TABLE

Table 5-5 P2030BD Ship C HP-100/60 Common (2385413 BOM, rev 1) (Continued)

1.4.1 Short Cable Collector (2281840-5)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2371450-2	1.000	Yes	RUN#55A-PDU TO RACEWAY GND
2	2371450-4	1.000	Yes	RUN#56A-CONSOLE TO RACEWAY
3	2373436-3	1.000	Yes	RUN #102A - GANTRY TO CON
4	2117848-7	1.000	Yes	20000 +/- 150 FIBER OPTIC
5	2371133-2	1.000	Yes	CONSOLE CABLES - SHORT KI
6	2371134-2	1.000	Yes	GANTRY CABLES - SHORT KIT

Table 5-6 Short Cable Collector (2281840-5 BOM, rev 0)

1.4.2 Long Cable Collector (2281840-4)

FRU = YES

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
1	2371450	1.000	Yes	RUN#55-PDU TO RACEWAY GND
2	2371450-3	1.000	Yes	RUN#56-CONSOLE TO RACEWAY
3	2373436-2	1.000	Yes	RUN #102 - GANTRY TO CONS
4	2117848-2	1.000	Yes	24384 +/- 150 FIBER OPTIC
5	2371133	1.000	Yes	CONSOLE CABLES - LONG KIT
6	2371134	1.000	Yes	GANTRY CABLES - LONG KIT

Table 5-7 Long Cable Collector (2281840-4 BOM, rev 0)

Section 2.0 — Small Room Catalog (5191743)

MARK	NUMBER	QTY	FRU	DESCRIPTION
1	5179183	1.0		Dolly, pivoting, small room, RH {For further detail, see Section 2.1 , beginning on page 376 .}
2	5179183-2	1.0		Dolly, pivoting, small room, LH {For further detail, see Section 2.2 , beginning on page 378 .}
3	5184632	1.0		Small Room CT Gantry Installation Dolly
4	2373436-7	1.0		Cable, FE Laptop to Console, 23000mm Long

Table 5-8 Small Room Catalog (5191743 BOM rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION
9	5176413	2.0	Yes	CLAMPING LEVER, ADJUSTABLE M10 X 20MM LG MORTON MACHINE WORKS #AH-933
1	5154481	1.0		Post, adjustable for pivoting dolly used on VCT gantry front cover
2	5164802	1.0		Post, upper for pivoting dolly used on VCT gantry front cover
3	5173585	1.0		Leg, stationary for pivoting dolly used on VCT gantry front cover
4	5165043	1.0		Leg, pivoting for pivoting dolly used on VCT gantry front cover
6	5173375	1.0		Bracket, pivot, RH for pivoting dolly used on VCT gantry front cover
8	2290611	2.0		PIN-DOLLY TO COVER CONNECTING
10	5179332	1.0		LANYARD, .0469 X 18 TERMINAL AND LOOP, FASTENAL #0157360
11	2283554	1.0		PIN, H2 DOLLY PIVOT LOCKING
12	2277431	1.0		KNOB, PLASTIC KNURLED CAMPING KNOB
13	2277430	1.0		E-RING, BOWED RETAINING RING
14	2277432-2	1.0		STEEL COMPRESSION SPRING_.36 OD 1.50
15	2283321	2.0		BUSHING, FLANGED NYLON .26 ID .50 LG
16	2283320	2.0		SPEED NUT, SHAFT DIA. 8MM
17	46-220205P7	2.0		ID 0.812 OD 2.00 THK 0.156 F70B5A (FOR 3/4 IN.SCREW)
18	46-208949P8	4.0		003/4-10 1.125 HEX X .641 THK F70B5A
19	5179690	2.0		COTTERLESS HITCHPIN, .375 X 2.375 FASTENAL #0157085
20	5179187	2.0		LANYARD, .093 X 8 PIVOT POINT #NL93F-8-2-26
21	2278991	3.0		PRESS LOCK RIVET
22	5166930	2.0		Spacer, caster for pivoting dolly used on VCT gantry front cover
23	5177454	4.0		CASTER, SWIVEL, 3.5 HIGH, SHEPHERD CASTER CORP. #102241
24	5176832	1.0		BOLT, HCS .75-10 X 2 GR5 ZINC
25	5179441	2.0		WASHER, SAE FW .75 ZINC, ID 0.812, OD 1.47, THK 0.134
26	5176393	1.0		COTTERLESS HITCHPIN, .375 X 1, FASTENAL #0157080
27	46-170684P2	0.0		LOCTITE 271, 10CC BOTTLE 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER. (REV, TK, 1/92).
28	5178519	16.0		NUT, M6 X 1.0, 9.8 HEX X 5.8 THK, NYLOCK, ZINC

Table 5-9 Dolly, Pivoting, Small Room, RH (5179183 BOM, rev 2)

2.2 Dolly, Pivoting, Small Room, LH (5179183-2)

FRU = NO




Figure 5-2 Dolly, Pivoting, Small Room (5179183 sheet 1, rev 3)

MARK	NUMBER	QTY	FRU	DESCRIPTION
9	5176413	2.0	Yes	CLAMPING LEVER, ADJUSTABLE M10 X 20MM LG MORTON MACHINE WORKS #AH-933
1	5154481	1.0		Post, adjustable for pivoting dolly used on VCT gantry front cover
2	5164802	1.0		Post, upper for pivoting dolly used on VCT gantry front cover
3	5173585	1.0		Leg, stationary for pivoting dolly used on VCT gantry front cover
4	5165043	1.0		Leg, pivoting for pivoting dolly used on VCT gantry front cover
7	5174002	1.0		Bracket, pivot, LH for pivoting dolly used on VCT gantry front cover
8	2290611	2.0		PIN-DOLLY TO COVER CONNECTING
10	5179332	1.0		LANYARD, .0469 X 18 TERMINAL AND LOOP, FASTENAL #0157360
17	46-220205P7	2.0		ID 0.812 OD 2.00 THK 0.156 F70B5A (FOR 3/4 IN.SCREW)
18	46-208949P8	4.0		003/4-10 1.125 HEX X .641 THK F70B5A
19	5179690	2.0		COTTERLESS HITCHPIN, .375 X 2.375 FASTENAL #0157085
20	5179187	2.0		LANYARD, .093 X 8 PIVOT POINT #NL93F-8-2-26
21	2278991	3.0		PRESS LOCK RIVET
22	5166930	2.0		Spacer, caster for pivoting dolly used on VCT gantry front cover
23	5177454	4.0		CASTER, SWIVEL, 3.5 HIGH, SHEPHERD CASTER CORP. #102241
24	5176832	1.0		BOLT, HCS .75-10 X 2 GR5 ZINC
25	5179441	2.0		WASHER, SAE FW .75 ZINC, ID 0.812, OD 1.47, THK 0.134
26	5176393	1.0		COTTERLESS HITCHPIN, .375 X 1, FASTENAL #0157080
27	46-170684P2	0.0		LOCTITE 271, 10CC BOTTLE 12 MONTH SHELF LIFE FROM DATE OF SHIPMENT FROM THE MANUFACTURER. (REV, TK, 1/92).
28	5178519	16.0		NUT, M6 X 1.0, 9.8 HEX X 5.8 THK, NYLOCK, ZINC

Table 5-10 Dolly, pivoting, small room, LH (5179183-2 BOM, rev 2)

End of Chapter.

© 2003, General Electric Company.

GE Medical Systems, a General Electric Company, going to market as GE Healthcare.

3000 N. Grandview Boulevard

Waukesha, Wisconsin 53188

USA

www.gehealthcare.com

