

GE Healthcare

LightSpeed 5.X Pro16

Illustrated Parts List

OPERATING DOCUMENTATION



2349187-100
Rev 18

Chapter 2

Table

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2115996-3	1.000	Yes	AXIAL HEADHOLDER ASM WITH
4	2115990-3	1.000	Yes	LCC HSA METALLESS SUPINE
5	2115993-2	1.000	Yes	LCC METALLESS CRADLE EXTE
8	2331933	1.000	Yes	METALLESS COMPATIBLE PHAN
1	2335179	1.000	No	TABLE ASSY - H POWER {For further detail, see Section 1.0 , beginning on page 30 .}
2	2330508	1.000	No	FOOT SWITCH ASSEMBLY {For further detail, see Section 2.0 , beginning on page 74 .}

Table 2-1 Cat Coll/H Power Table with Footswitch (2373155 BOM, rev 0)

$$FRU = NO$$

Figure 2-1 Table Assembly - H Power (2335179, sheet 1, rev 0)




Figure 2-2 Table Assembly - H Power (2335179, sheet 2, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
7	2271244	1.000	Yes	RIGHT BASE COVER ASM FOR
8	2271247	1.000	Yes	LEFT BASE COVER ASM FOR H
16	2261683	1.000	Yes	COVER, ACTUATOR, (REF 46-
17	2372371	1.000	Yes	COVER TABLE TOP LEFT {For further detail, see Section 1.2 , beginning on page 38 .}
18	2372372	1.000	Yes	COVER TABLE TOP RIGHT {For further detail, see Section 1.3 , beginning on page 40 .}
19	46-296594P1	1.000	Yes	DRIVE COVER
21	46-296679P2	1.000	Yes	ZEUS TABLE, REAR END COVE
29	2265851	2.000	Yes	BRACKET, LR SIDE PANEL (R
30	2265850	2.000	Yes	BRACKET, LF SIDE PANEL (R
31	2261669	2.000	Yes	BRKT,U-FRONT SIDE PNL, (R
32	2261668	2.000	Yes	BRKT, U-REAR SIDE PNL, (R
36	46-297329P2	2.000	Yes	1/4 TURN RECEPTACLE - FIX
38	2329950	1.000	Yes	CHANNEL,COVER SUPPORT
43	2146815-2	26.000	Yes	10-32 X 3/8 HEX SOC BUTTO
47	46-208562P20	14.000	Yes	SCREW
52	2265882	8.000	Yes	SIDE PIVOT BUSHING (REF 4
72	46-297678P1	1.000	Yes	HINGES
74	46-170490P50	1.000	Yes	SPRING-INCH EXTENSION 4 I
81	12-276	4.000	Yes	ID 0.2040 OD 0.5630 THK .
103	2265845	5.000	Yes	PAD, LEVELING (REF 210587
108	2115992-2	1.000	Yes	860-941208 METALLESS CRADLE {For further detail, see Section 1.4 , beginning on page 42 .}
109	2117299	1.000	Yes	METALLESS CRADLE COMPATIB
135	2287357	1.000	Yes	HANDLE, CRADLE
137	2343484	1.000	Yes	COVER, CRADLE ENCODER
138	2348433	1.000	Yes	BRACKET, COVER MOUNTING-R

Table 2-2 Table Assembly - H Power (2335179 BOM, rev 4)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
139	2346979	1.000	Yes	COVER - TCB EMC
4	2261682	1.000	No	BRKT,ACTUAT CVR MT,(REF 4
5	2265868	1.000	No	BRACKET (REF 46-296820P1)
6	2269167-2	1.000	No	SIDE PANEL KIT FOR H2 LCC {For further detail, see Section 1.1 , beginning on page 36 .}
13	2261673	2.000	No	BKT, REAR COVER, (REF 46-
22	2261708	2.000	No	COVER GUARD, RAILS,(REF 4
33	2265870	4.000	No	NUT PLATE (REF 46-297171P
35	2261688	1.000	No	BRKT, LWR COV SUP (REF 46
45	46-208562P2	1.000	No	SCREW
46	46-208562P3	4.000	No	SCREW
55	46-297155G1	4.000	No	ZEUS TABLE, 9" LONG CABLE
59	2261690	2.000	No	HEX STDF,9.5MM TAP,(REF 4
61	46-170498P23	6.000	No	SCREW
62	46-221523P10	6.000	No	HOLE DIA 0.437 MAT'L THK
63	46-170012P63	2.000	No	0010-32 HEX .375 NUT HT
64	46-220186P20	4.000	No	1/8" 120 DEGREE COUNTERSU
65	46-170498P41	4.000	No	SCREW
71	46-297695P1	2.000	No	COTTER PIN
73	2265878	1.000	No	PIN (REF 46-297676P1)
75	2265858	1.000	No	HANDLE ADDITION (REF 46-2
80	46-220181P5	8.000	No	WASHER #10 (.190) SCREW S
85	2146815-2	18.000	No	10-32 X 3/8 HEX SOC BUTTO
86	2265863	3.000	No	SHIM 0.80MM (REF 46-19635
87	2265863-2	3.000	No	SHIM 0.40MM (REF 46-19635
88	46-170106P1	AR	No	DBL COATED FILM .007 THK
92	1000904P358	1.000	No	ID 0.380 OD 1.000 THK .
93	2265862	1.000	No	BOLT, SHIPPING (REF 21847

Table 2-2 Table Assembly - H Power (2335179 BOM, rev 4) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
94	2265861	AR	No	SPACER PLT 0.50MM (REF 46
95	2265861-2	AR	No	SPACER PLT 0.80MM (REF 46
96	2265861-3	AR	No	SPACER PLT 1.00MM (REF 46
97	2265861-4	AR	No	SPACER PLT 1.20MM (REF 46
98	2265861-5	1.000	No	SPACER PLT 1.60MM (REF 46
100	46-208561P147	5.000	No	SCREW
101	2265840	5.000	No	ADJUSTER, LEVELING PAD (R
102	2265846	5.000	No	LOCK RING (REF 2106207) D
116	2394485	1.000	No	HPOWER TABLE ETL LABEL
117	2117390	1.000	No	LABEL, CE MARKING COMPLIA
118	2132620	1.000	No	CE MARKING LABEL CLASS A
119	2103580-12	8.000	No	SCREW COUNTERSUNK HEA 6 M
120	46-328417P12	8.000	No	SCREW HEXAGON SOCKET 6 MM
121	2335178	1.000	No	TABLE SUB ASSY - H POWER {For further detail, see Section 1.5 , beginning on page 44 .}
122	2384711	1.000	No	RATING PLATE
126	2148658	1.000	No	ESD GROUND STRAP, BRAIDED
134	46-303268P1	1.000	No	PRODUCT LOCATOR CARD FORM
136	46-328417P20	2.000	No	SCREW HEXAGON SOCKET 6 MM
140	P9140ET	1.000	No	NAME PLATE
141	46-302200P7	1.000	No	SERIAL #, MODEL#
142	46-170171P22	AR	No	TAPE OLD
143	2327686	2.000	No	WING NUT M6
144	46-328430P3	2.000	No	WASHER PLAIN - NORMAL 6.4
145	46-328432P3	2.000	No	WASHER LOCK - SPRING 6.1

Table 2-2 Table Assembly - H Power (2335179 BOM, rev 4) (Continued)

1.1 Side Panel Kit (2269167-2)

FRU = NO

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
9	2259861-2	1.000	Yes	SIDE PANEL,RT REAR, HELIO
10	2259862-2	1.000	Yes	SIDE PANEL,LEFT REAR, HEL
11	2259863-2	1.000	Yes	SIDE PANEL,RIGHT FRONT, H
12	2259864-2	1.000	Yes	SIDE PANEL,LEFT FRONT, HE

Table 2-3 Side Panel Kit for H2 LCC Cathay Table (2269167-2 BOM, rev 0)

1.2 Cover Table Top - Left (2372371)

FRU = YES




Figure 2-3 Cover Table Top - Left (2372371, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	46-208747P1	4.000	Yes	CABLE TIE,ADH.BACKED
4	46-297698P3	1.000	Yes	27" HIGH PROFILE BEAD, RI
5	46-297329P1	2.000	Yes	1/4 TURN RECEPTACLE WITH
6	46-220186P20	4.000	Yes	1/8" 120 DEGREE COUNTERSU
7	46-297721P1	1.000	Yes	4-POS PLUG HOUSING (DRAWE
8	46-158905P16	2.000	Yes	MOUNTING SCREW FOR DRAWER
9	46-297328P2	2.000	Yes	1/4 TURN STUD AND TWO GRO
10	46-208758P1	4.000	Yes	3.62" X .094" CABLE TIE
11	2372346	1.000	Yes	CHANNEL TAPESWITCH
2	2372348	1.000	No	COVER, TOP LEFT
12	2363367	6.000	No	VINYL CAP 1/8 INCH
13	2363369	6.000	No	T SHAPE HEAD STEEL CLIP
14	2199652	1.000	No	P9235LZ LABEL

Table 2-4 Cover Table Top - Left (2372371 BOM, rev 01)

1.3 Cover Table Top - Right (2372372)

FRU = YES




Figure 2-4 Cover Table Top - Right (2372372, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	46-208747P1	4.000	Yes	CABLE TIE,ADH.BACKED
4	46-297698P3	1.000	Yes	27" HIGH PROFILE BEAD, RI
5	46-297329P1	2.000	Yes	1/4 TURN RECEPTACLE WITH
6	46-220186P20	4.000	Yes	1/8" 120 DEGREE COUNTERSU
7	46-297721P1	1.000	Yes	4-POS PLUG HOUSING (DRAWE
8	46-158905P16	2.000	Yes	MOUNTING SCREW FOR DRAWER
9	46-297328P2	2.000	Yes	1/4 TURN STUD AND TWO GRO
10	46-208758P1	4.000	Yes	3.62" X .094" CABLE TIE
11	2372346	1.000	Yes	CHANNEL TAPESWITCH
2	2372351	1.000	No	COVER, TOP RIGHT
12	2363367	6.000	No	VINYL CAP 1/8 INCH
13	2363369	6.000	No	T SHAPE HEAD STEEL CLIP
14	2199652	1.000	No	P9235LZ LABEL

Table 2-5 Cover Table Top - Right (2372372 BOM, rev 1)

FRU = YES



MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
4	2268003	2.000	Yes	10-32 X .437 LONG, FLAT H
5	2222889-2	1.000	Yes	KEEPER FOR SELF-CATCHING
6	2112870-2	1.000	No	METALLESS CRADLE WITH VEL
7	46-303268P1	1.000	No	PRODUCT LOCATOR CARD FORM
3A	2285525-2	1.000	No	RATING PLATE
3B	46-302200P7	ALT	No	SERIAL #, MODEL#

Table 2-6 Metalless Cradle (2115992-2 BOM, rev 0)

1.5 Table Subassembly - H Power (2335178)

FRU = NO




Figure 2-6 Table Subassembly (2335178, sheet 1, rev 2)








Figure 2-9 Table Subassembly (2335178, sheet 4, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	2269168-2	1.000	Yes	CARRIAGE ASSEMBLY FOR HPO {For further detail, see Section 1.5.1 , beginning on page 52 .}
4	2269177	1.000	Yes	LATCH ASSEMBLY FOR CATHAY {For further detail, see Section 1.5.2 , beginning on page 54 .}
7	46-296370P1	2.000	Yes	ROLLER
11	2339120	1.000	Yes	CRADLE DRIVE ASSY - H PWR {For further detail, see Section 1.5.3 , beginning on page 56 .}
12	2342727-2	1.000	Yes	ENCODER,MULTI-TURN W/ TAK
14	46-136334P23	2.000	Yes	SPDT, 15A @125/250VAC BAS
19	2337217	1.000	Yes	ENCODER ASSY, ELEVATION - {For further detail, see Section 1.5.4 , beginning on page 62 .}
20	46-170047P7	1.000	Yes	125 PITCHES 19.635 LENG
29	12-276	5.000	Yes	ID 0.2040 OD 0.5630 THK .
32	2146815-2	20.000	Yes	10-32 X 3/8 HEX SOC BUTTO
34	2265877	1.000	Yes	CALIBRATION PIN (REF 46-2
37	2261701	1.000	Yes	BAR, CALIBRATION, (REF 46
41	2261724	1.000	Yes	CAL PLATE FOR H2 LCC CATH
51	2342729-3	2.000	Yes	H-POWER, SERVO AMP DRIVE
54	2258712-4	1.000	Yes	POWER SUPPLY AC-DC SINGLE
57	2258704	1.000	Yes	POWER SUPPLY AC-DC SINGLE
61	2349543	1.000	Yes	MTCB
62	2348054	2.000	Yes	SPACER, HEX-TCB COVER 10-
79	46-297155G1	2.000	Yes	ZEUS TABLE, 9" LONG CABLE
80	46-297155G4	1.000	Yes	GROUND CABLE_24" LONG
90	2380051	1.000	Yes	FILTER COVER
5	2261680	1.000	No	PLATE, TABLE TOP, (REF 46
6	2265867	1.000	No	BLOCK GUIDE ROLLER (REF 4
8	2261723	2.000	No	BLOCK, BUMPER, (REF 46-27

Table 2-7 Table Subassembly (2335178 BOM, rev 5)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
9	46-221597P2	4.000	No	MINOR VBM-3003 NAT RUBBER
10	2265854	1.000	No	PLATE, SWITCH (REF 46-278
13	2265867-2	1.000	No	WIDE GROOVE ROLLER (REF 4
15	2261679	1.000	No	BRACE, REAR (REF 46-29681
16	2261698	1.000	No	BRACKET, SWITCH, (REF 46-
22	46-296626P1	1.000	No	LARGE GEAR
25	46-170623P2	AR	No	BLACK, HEAVY-DUTY LUBRICA
26	46-170498P54	4.000	No	SCREW
27	46-170498P13	8.000	No	SCREW
30	46-220181P5	23.000	No	WASHER #10 (.190) SCREW S
31	46-328417P15	4.000	No	SCREW HEXAGON SOCKET 6 MM
38	46-170684P2	AR	No	LOCTITE 271, 10CC BOTTLE
39	46-170686P2	AR	No	LOCTITE 242, 10CC
40	2146815	3.000	No	10-32 X 1/2 HEX SOC BUTTO
43	46-328417P12	16.000	No	SCREW HEXAGON SOCKET 6 MM
44	12-275	4.000	No	WASHER, ID .2656, OD.75,
45	46-328432P3	2.000	No	WASHER LOCK - SPRING 6.1
46	2335177	1.000	No	TABLE FRAME ASSY - H POWE {For further detail, see Section 1.5.5 , beginning on page 64 .}
47	46-170015P28	1.000	No	SCREW
52	2264567	1.000	No	BAR, WIRE GROUND (REF 46-
55	2346974	1.000	No	BRACKET - 24V PS
56	2346976	2.000	No	BRACKET - POWER SUPPLY CL
58	2346975	1.000	No	BRACKET - 5V PS
60	2346977	1.000	No	PANEL - TCB MOUNTING
63	46-170171P22	AR	No	TAPE OLD
64	46-328417P13	4.000	No	SCREW HEXAGON SOCKET 6 MM
65	46-328430P3	4.000	No	WASHER PLAIN - NORMAL 6.4

Table 2-7 Table Subassembly (2335178 BOM, rev 5) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
66	46-328417P1	3.000	No	SCREW HEXAGON SOCKET 4 MM
67	46-328417P73	3.000	No	SCREW HEXAGON SOCKET 3 MM
68	46-170497P7	12.000	No	SCREW
69	46-170497P38	8.000	No	SCREW
70	2339118	1.000	No	BRACKET, CRADLE ENCODER M
71	46-170497P4	5.000	No	SCREW
74	46-328425P1	1.000	No	NUT HEXAGON 4 MM
75	46-328432P2	1.000	No	WASHER LOCK - SPRING 4.1
76	46-208758P3	4.000	No	7.31" X .184" SELF-LOCKIN
77	1000274P3	2.000	No	SPACER, 1/2" OD 9/32" ID
78	46-328431P2	2.000	No	WASHER PLAIN - LARGE 6.4
81	12-297	8.000	No	ID 0.140 OD 0.312 THK .
82	46-220181P3	8.000	No	WASHER HELICAL#6 SCREW SI
83	46-297335P3	1.000	No	16A, 250V 2-STAGE PWR LIN
84	46-328430P2	4.000	No	WASHER PLAIN - NORMAL 4.3
85	46-328432P2	9.000	No	WASHER LOCK - SPRING 4.1
86	46-328417P4	4.000	No	SCREW HEXAGON SOCKET 4 MM
87	46-328430P3	1.000	No	WASHER PLAIN - NORMAL 6.4
88	46-328432P3	1.000	No	WASHER LOCK - SPRING 6.1
89	46-328417P11	1.000	No	SCREW HEXAGON SOCKET 6 MM

Table 2-7 Table Subassembly (2335178 BOM, rev 5) (Continued)

1.5.1 Carriage Assembly (2269168-2)

FRU = NO




Figure 2-10 Carriage Assembly (2269168-2, sheet 2, rev1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	46-218193P1	2.000	Yes	WHEEL DIA 01.21 WIDTH 0.4
5	46-208580P6	3.000	Yes	CAM FOLLOWER
6	46-208580P19	1.000	Yes	CAM FOLLOWER
10	46-170498P2	4.000	Yes	SCREW
15	2356591	1.000	Yes	BRACKET, CRADLE ENCODER D
2	2262161-2	1.000	No	CRADLE CARRIAGE MACH.(REF
4	2264566	2.000	No	STUD, ECCENTRIC (REF 46-2
7	2261695	1.000	No	LATCH ACTUATOR, (REF 2133
8	2264561	1.000	No	ACTU., LONG.LIMIT SWITCH(
9	46-170012P77	3.000	No	LOCK NUT, 1/4-28, .437 X
14	46-170686P2	AR	No	LOCTITE 242, 10CC
16	46-170497P55	1.000	No	SCREW

Table 2-8 Carriage Assembly (2269168-2 BOM, rev 1)

1.5.2 Latch Assembly (2269177)
FRU = YES




Figure 2-11 Latch Assembly (2269177, sheet 1, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
8	46-296972G1	1.000	Yes	SOLENOID ASM
2	2261672	1.000	No	BRKT, TABLE HOME LATCH,(R
3	2261697	1.000	No	BLK,PVT, HOME LATCH(REF 4
4	2261721	1.000	No	BAR, (REF 46-278154P1)
5	2261725	1.000	No	BRKT, LATCH ASM, (REF 46-
6	2261715	1.000	No	ARM, (REF 46-278495P1)
7	2261681	1.000	No	BRACKET, SPRING, (REF 46-
9	46-170490P54	1.000	No	SPRING-INCH EXTENSION 1 I
10	46-170497P36	3.000	No	SCREW
11	46-208732P14	1.000	No	SCREW
12	46-170012P63	1.000	No	0010-32 HEX .375 NUT HT
13	46-170020P80	1.000	No	NOMINAL DIA 0.250 X 0.750
14	46-170020P66	1.000	No	NOMINAL DIA 0.125 X 0.250
15	46-170020P85	1.000	No	NOMINAL DIA 0.125 X 0.75
16	403B335P17	1.000	No	02 POSITION TERMINAL STRI
17	403B320P94	1.000	No	2 POSITION MARKER STRIP,
18	46-170497P38	2.000	No	SCREW
19	46-170497P7	4.000	No	SCREW
20	46-170686P2	AR	No	LOCTITE 242, 10CC

Table 2-9 Latch Assembly (2269177 BOM, rev 0)

FRU = YES

Figure 2-12 Cradle Drive Assembly (2339120, sheet 1, rev 0)





MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
17	46-313575P1	4.000	Yes	HIGH COLLAR LOCK WASHER F
19	46-221532P18	1.000	Yes	PITCH 1/5 WIDTH 00.250
22	46-296972G1	1.000	Yes	SOLENOID ASM
59	46-311873P8	1.000	Yes	PIN, M3 X 18
84	46-328417P504	4.000	Yes	SCREW HEXAGON SOCKET 3 MM
2	2265831	1.000	No	END PLATE, LH (REF 46-296
3	2339121	1.000	No	PLATE, END - RH
5	46-136282P30	2.000	No	BEARING BALL, SINGLE RO 2
6	46-212989P2	2.000	No	WASHER
7	46-208562P8	2.000	No	SCREW 3/8X16X1" FLAT HEAD
9	2265818	1.000	No	IDLER SHAFT (REF 46-29603
10	2342731	1.000	No	H-POWER, AC CRADLE MOTOR
12	403B335P10	1.000	No	06 POSITION TERMINAL STRI
13	46-170497P38	6.000	No	SCREW
14	403B320P86	1.000	No	NUMBERS 1 THU 6, LEFT TO
15	403B320P97	1.000	No	2 POSITION MARKER ST
16	46-170498P80	8.000	No	SCREW
18	1000904P116	1.000	No	ID 0.1495 OD 0.6250 THK .
20	46-170576P7	1.000	No	PITCH 1/5 TEETH 0L6 FAC
21	46-170576P6	1.000	No	PITCH 1/5 TEETH 013 FAC
23	2261718	1.000	No	BRACKET, SOLENOID, (REF 4
24	2265829	1.000	No	SPACER BAR (REF 46-296063
25	2265830	2.000	No	PIVOT PIN (REF 46-296064P
26	46-170686P2	AR	No	LOCTITE 242, 10CC
27	2265833	2.000	No	SUPPORT BAR (REF 46-29607
28	46-208562P7	4.000	No	SCREW
30	2265826	1.000	No	SHAFT (REF 46-296048P1)

Table 2-10 Cradle Drive Assembly (2339120 BOM, rev 3)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
31	2265824	1.000	No	BEARING SUPPORT (REF 46-2
32	46-170498P79	1.000	No	SCREW
35	46-221451P1	AR	No	EXXON BEACON 325 35 LB P
36	46-170498P115	6.000	No	SCREW
37	2265865	1.000	No	KEY 2.38SQX15MMLG (REF 46
38	46-220179P8	1.000	No	EXT AXIAL MOUNT 0.375 SHA
39	46-136282P33	2.000	No	BEARING
40	46-136294P9	2.000	No	BEARING
41	2265834	1.000	No	SHAFT (REF 46-296110P1)
42	2265825	1.000	No	BEARING SUPPORT (REF 46-2
43	46-136303P10	1.000	No	BEARING
44	46-136282P62	2.000	No	BEARING BALL, SINGLE RO 1
45	46-220179P34	1.000	No	EXT AXIAL MOUNT 0.669 SHA
46	46-170498P84	6.000	No	SCREW, CAP SOCKET HEAD 8-
48	46-296368G1	1.000	No	
49	403B335P17	1.000	No	02 POSITION TERMINAL STRI
50	46-208718P2	2.000	No	5/16 ID LOOP TYPE TA MFG
51	509A859P69	1.000	No	SPRING-INCH COMPRESSION 1
55	46-297820P1	1.000	No	SHOULDER SCREW
56	2261694	1.000	No	CLUTCH LINK BAR, (REF 46-
60	2261717	1.000	No	BRACKET, STOP, (REF 46-29
61	46-221560P4	1.000	No	SCREW
63	46-220277P6	1.000	No	EXT RADIAL MOUNT 0.219 SH
64	46-170012P63	1.000	No	0010-32 HEX .375 NUT HT
66	2265865-2	1.000	No	KEY 4.76SQX50MMLG
70	46-170497P7	5.000	No	SCREW
72	46-297118G1	1.000	No	CRDR SIG CBL
73	1000904P514	2.000	No	I.D. 0.789 O.D. 1.100

Table 2-10 Cradle Drive Assembly (2339120 BOM, rev 3) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
74	46-208758P3	5.000	No	7.31" X .184" SELF-LOCKIN
76	46-297155G1	1.000	No	ZEUS TABLE, 9" LONG CABLE
77	46-220181P5	3.000	No	WASHER #10 (.190) SCREW S
78	12-276	2.000	No	ID 0.2040 OD 0.5630 THK .
79	46-170623P2	AR	No	BLACK, HEAVY-DUTY LUBRICA
85	46-328430P9	4.000	No	S N N 1 WASHER PLAIN - NORMAL 3.2
86	46-328425P10	4.000	No	S N N 2 NUT HEXAGON 3 MM
34A	46-222894P1	1.000	No	USM CORP HDC3C-86-2BL
34B	2361751	ALT	No	HARMONIC DRIVER
54A	46-297397P1	1.000	No	
54B	2265852-2	ALT	No	GEAR, LEVER, MACHINING (4
82A	2134525	1.000	No	ZEUS TABLE CRADLE DRIVE;
82B	2134525-2	ALT	No	ZEUS TABLE CRADLE DRIVE;
82C	2265827-2	ALT	No	IDLER ROLLER, NEW (REF 21
82D	2265827-3	ALT	No	IDLER ROLLER RECYCLD (REF
83A	2134526	1.000	No	ZEUS TABLE CRADLE DRIVE;
83B	2134526-2	ALT	No	ZEUS TABLE CRADLE DRIVE;
83C	2265820-2	ALT	No	DRIVE ROLLER,NEW (REF 213
83D	2265820-3	ALT	No	DRIVE ROLLER RECYCLD (REF

Table 2-10 Cradle Drive Assembly (2339120 BOM, rev 3) (Continued)

1.5.4 Encoder Assembly, Elevation - H Power (2337217)

FRU = YES




Figure 2-15 Encoder Assembly, Elevation - H Power (2337217, sheet 1, rev 1)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
2	2342726-2	1.000	Yes	ABSOLUTE ENCODER, SINGLE
4	46-296780P2	3.000	Yes	SYN-CLAMP
6	46-296259P9	1.000	Yes	ROCOM COUPLING: BORES 06M
10	46-136334P57	2.000	Yes	SPDT, MINIATURE BASIC
3	2337218	1.000	No	BLOCK, ELEVATION ENCODER
5	2346925	1.000	No	SHAFT, ELEVATION ENCODER
7	46-170048P8	1.000	No	MIN-E-PITCH, 16 T, PD 0.8
8	2265261	1.000	No	CAM, ELEV SWITCH (REF 46-
9	2261699	1.000	No	BRKT, SWT ELEV ASM,(REF 4
11	46-136282P75	2.000	No	BEARING
12	46-220277P7	1.000	No	EXT RADIAL MOUNT 0.250 SH
13	2265264	1.000	No	NUT PLATE (REF 46-296945P
14	46-170497P7	2.000	No	SCREW
15	46-170498P125	3.000	No	SCREW
16	46-170528P23	1.000	No	SHAFT SIZE .250 MAT 2024-
18	46-170020P63	2.000	No	NOMINAL DIA 0.250 X 0.625
19	46-220181P2	1.000	No	#4 SPRING LOCK WASHER, ST
20	1000904P438	2.000	No	ID 0.12, OD 0.312, THK .0
23	46-297539P2	1.000	No	NO. 6 HI-COLLAR SPRING,
25	46-220205P1	1.000	No	ID 0.3125 OD 0.75 THK

Table 2-11 Encoder Assembly, Elevation - H Power (2337217 BOM, rev 2)

1.5.5 Table Frame Assembly (2335177)

FRU = NO




Figure 2-16 Table Frame Assembly (2335177, sheet 1, rev 1)




Figure 2-17 Table Frame Assembly (2335177, sheet 2, rev 0)




Figure 2-18 Table Frame Assembly (2335177, sheet 3, rev 1)






Figure 2-20 Table Frame Assembly (2335177, sheet 5, rev 0)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
4	2347434-2	1.000	Yes	SIGNAL HARNESS ASSY - H P
6	2343479	1.000	Yes	ACTUATOR, H-POWER TABLE E
8	2100671	2.000	Yes	SUSPA # 016-03570 C (ZEUS
18	2259859	1.000	Yes	CARRIAGE_RAIL, (REF 46-27
19	2259860	1.000	Yes	CARRIAGE RAIL, V, (REF 46
23	46-297350P1	2.000	Yes	SPACER
25	46-220277P15	2.000	Yes	EXT RADIAL MOUNT 0.312 SH
26	46-220277P11	1.000	Yes	EXT RADIAL MOUNT 0.625 SH
35	46-170684P1	AR	Yes	LOCTITE 271, 0.5CC VIAL
37	2265876	2.000	Yes	PIN, GAS SPRING (REF 46-2
42	2373415	1.000	Yes	HARNESS TABLE AC PWR
44	46-297329P2	6.000	Yes	1/4 TURN RECEPTACLE - FIX
63	46-297687P1	1.000	Yes	8 OHM, 4W, 3-1/2" SPE
79	46-297698P1	2.000	Yes	11" LOW PROFILE BEAD, RIB
87	2367487	1.000	Yes	GROUND - TABLE TO RACEWAY
88	2146815-2	21.000	Yes	10-32 X 3/8 HEX SOC BUTTO
90	12-276	6.000	Yes	ID 0.2040 OD 0.5630 THK .
97	46-297805G1	2.000	Yes	JUMPER CONN.
98	46-297805G2	2.000	Yes	JUMPER CONN.
104	2106866-2	1.000	Yes	.75 IN WIDE INSULATED CAB
120	46-208562P20	6.000	Yes	SCREW
128	2342733	1.000	Yes	H-POWER, AC ELEVATION MOT
129	46-313575P3	4.000	Yes	HIGH COLLAR LOCK WASHER F
133	2373417	1.000	Yes	HARNESS BULKHEAD TO FILTE
5	2265855	1.000	No	PIN, ROD END (REF 46-2967
9	2255160-2	1.000	No	FRONT LEG FOR H2 LCC CATH
12	46-136342P135	4.000	No	BUSHING

Table 2-12 Table Frame Assembly (2335177 BOM, rev 3)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
13	2265856	1.000	No	SPACER, ROD END (REF 46-2
14	2265853	1.000	No	PIN, ACTUATOR (REF 46-278
15	2255158-2	1.000	No	REAR LEG FOR H2 LCC CATHA
16	2259855	1.000	No	CHANNEL,LEFT SUPPORT,(REF
17	2338713	1.000	No	CHANNEL, RIGHT SUPPORT
20	2261670	2.000	No	PIN, TABLE TOP, FRONT (RE
21	2261671	2.000	No	PIN, TABLE TOP, REAR, (RE
22	2261678	1.000	No	BRACE, FRONT, (REF 46-296
24	46-220179P15	4.000	No	EXT AXIAL MOUNT 0.625 SHA
28	46-328417P26	2.000	No	SCREW HEXAGON SOCKET 10 M
29	46-328417P15	4.000	No	SCREW HEXAGON SOCKET 6 MM
30	46-170498P14	32.000	No	SCREW
32	46-170012P81	1.000	No	005/8-18 HEX .940 NUT
33	46-220181P6	33.000	No	1/4" SCREW SIZE (0.250),
34	2261685	1.000	No	BRKT, R-CHANL CON, (REF 4
36	2261722	1.000	No	BLOCK, MOUNTING, (REF 46-
40	12-128	4.000	No	WASHER,FLAT
41	46-170623P2	AR	No	BLACK, HEAVY-DUTY LUBRICA
43	2265871	6.000	No	BRACKET TBL BASE CVR (46-
45	46-220186P20	12.000	No	1/8" 120 DEGREE COUNTERSU
46	2329243	1.000	No	SHAFT
47	2329245	1.000	No	ROLLER MOUNT
48	46-170020P2	2.000	No	NOMINAL DIA 0.0938X 0.500
49	12-275	2.000	No	WASHER, ID .2656, OD.75,
50	2261663	2.000	No	BRKT, TBL TOP COVER, (REF
51	2261692	1.000	No	SPRING CLAMP, (REF 46-297
52	46-297705P2	2.000	No	CONDUCTIVE TAPE
54	46-297705P4	1.000	No	MET-CON TAPE

Table 2-12 Table Frame Assembly (2335177 BOM, rev 3) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
61	2261689	1.000	No	CABLE TROUGH, (REF 46-297
62	2261674	1.000	No	SPEAKER GRILL, (REF 46-29
64	46-297686P2	1.000	No	SPEAKER INSULATION FOR H2
65	2265879	1.000	No	SPEAKER COVER (REF 46-297
66	46-297155G1	4.000	No	ZEUS TABLE, 9" LONG CABLE
69	2265875	2.000	No	CLAMP, GRND STRAP (REF 46
75	46-208759P4	6.000	No	13.93"X .301" SELF-LOCKIN
77	46-297648P1	1.000	No	CABLE CLAMP
80	2146815	20.000	No	10-32 X 1/2 HEX SOC BUTTO
81	1000904P395	4.000	No	ID 0.177 OD 0.312 THK .
82	46-208759P1	25.000	No	7.81" X .184" SELF-LOCKIN
83	46-208560P46	1.000	No	SCREW
84	46-170497P38	4.000	No	SCREW
85	46-170012P63	16.000	No	0010-32 HEX .375 NUT HT
86	46-170012P61	2.000	No	0006-32 HEX .313 NUT HT
89	46-297731P1	1.000	No	MOTOR LABEL
92	46-208758P9	2.000	No	11.08"X .140" SELF-LOCKIN
94	2261712	1.000	No	BRKT, LH TBL TOP CVR, (RE
95	2261712-2	1.000	No	BRKT, RH TBL TOP CVR,(REF
96	46-208758P1	2.000	No	3.62" X .094" CABLE TIE
99	46-208562P7	12.000	No	SCREW
100	46-170012P62	4.000	No	0008-32 HEX .343 NUT HT
101	2265880	4.000	No	SHOULDER SCREW (REF 46-29
113	2343440-2	1.000	No	BASE, TABLE - MACHINING H
114	46-328417P12	25.000	No	SCREW HEXAGON SOCKET 6 MM
117	2261676	2.000	No	PIN, TBL LOWER, FRONT,(RE
118	2261677	2.000	No	PIN, TBL LOWER, REAR, (RE
119	46-170498P32	6.000	No	SCREW

Table 2-12 Table Frame Assembly (2335177 BOM, rev 3) (Continued)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
121	46-328416P23	4.000	No	SCREW HEXAGON HEAD 10 MM
122	46-328431P3	4.000	No	WASHER PLAIN - LARGE 10.5
123	46-328432P4	4.000	No	WASHER LOCK - SPRING 10.2
124	46-170686P2	AR	No	LOCTITE 242, 10CC
127	46-220181P5	8.000	No	WASHER #10 (.190) SCREW S
130	46-328417P505	4.000	No	SCREW HEXAGON SOCKET 5 MM
131	46-328424P36	2.000	No	SCREW SET SCREW 5 MM 8 MM
132	46-297705P16	6.000	No	METALLIC CONDUCTIVE TAPE,
7A	2133760-2	1.000	No	ZEUS TABLE ACTUATOR BASE
7B	2261710-2	ALT	No	MACHINING,BASE MOUNT,(REF

Table 2-12 Table Frame Assembly (2335177 BOM, rev 3) (Continued)

Section 2.0 — Foot Switch Assembly (2330508)

FRU = NO

This Page Intentionally Blank.




Figure 2-21 Foot Switch Assembly (2330508, sheet 1, rev 2)

MARK	NUMBER	QTY	FRU	DESCRIPTION/NAME
3	46-296129G1	2.000	Yes	
28	2329241	4.000	Yes	FOOT PEDAL
29	2297598	1.000	Yes	FOOT SWITCH COVER
6	46-170497P8	6.000	No	SCREW
7	46-271110P1	1.000	No	LABEL, PROTECTIVE EARTH (
8	2329242	2.000	No	PEDAL ROD
9	2329240	1.000	No	GROUND FOOT SWITCH
10	46-221828P4	2.000	No	0.063 DIA X 0.625 LG18-
11	46-170012P63	3.000	No	0010-32HEX .375NUT HT
13	46-208560P65	11.000	No	SCREW
14	46-220181P6	11.000	No	1/4" SCREW SIZE (0.250),
15	1000274P289	3.000	No	SPACER
17	46-297971P1	2.000	No	TERMINAL
19	46-328417P12	2.000	No	SCREW HEXAGON SOCKET 6 MM
20	46-328430P3	2.000	No	WASHER PLAIN - NORMAL 6.4
21	2329239	2.000	No	GUARD BAR
22	2329246	6.000	No	SPACER GROUND BAR
23	46-170498P11	6.000	No	SCREW
24	46-170686P3	AR	No	LOCTITE 242, 50CC
25	2329247	2.000	No	SPRING PLATE
26	46-170497P20	12.000	No	SCREW
27	46-136353P135	4.000	No	SPRING-INCH COMPRESSION 1
30	2297597	1.000	No	FOOT SWITCH BASE

Table 2-13 Foot Switch Assembly (2330508 BOM, rev 5)

