



Service Request Number _____

Isocenter / Gradcal Data Sheet

Site Name: _____

Magnet Serial Number

FE Name: _____

Month: _____

Forward/Reverse Quadrature

☐ Check forward reverse quadrature (Check signal and swap the blue cables in the rear pedestal)Initial connections R1% R2% Switched connections R1% R2%

Iso Center Checks

Exam/Series/Image	Plane	Error Value
Old isoVectorZ		New isoVectorZ

Gradcal Checks

Initial Measurement	Measurement Spec 99 to 101	Final CV Value
X Gradient		cfxfull
Y Gradient		cfyfull
Z Gradient		cfzfull

New Measurement	Measurement Spec 99 to 101	Final CV Value
X Gradient		cfxfull
Y Gradient		cfyfull
Z Gradient		cfzfull

System Gain Data

Body 1100 X Curent Recon Scale Factor = New Recon Scale Factor Acceptance Specifications
Average Pixel Intensity Pixel intensity 1000-1200

Exam/Series/Image	Initial Recon scale	R1/R2/TG	System Frequency	Pixel Average	Pixel Std. Deveation	New Recon Scale

Head 1100 X Curent Recon Scale Factor = New Recon Scale Factor Acceptance Specifications
Average Pixel Intensity Pixel intensity 1000-1200

Exam/Series/Image	Initial Recon scale	R1/R2/TG	System Frequency	Pixel Average	Pixel Std. Deveation	New Recon Scale