

MRI DE-INSTALL FINAL REPORT

1290 Elm st. Suite B Ontario, Ca. 91761
Main: 951.734.3737, Fax: 951-735-3373



OEM ID	_____	SMR ID	_____
Site Name	_____	City, State	_____
Site Contact	_____	Phone Number	_____
SMR Rep	_____	Completion Date	_____

Inventory (Guided Install Hardware Tab)

Magnet	_____	Serial Number	_____	Helium	_____	Pres.	_____	Date	_____
Package Name	_____	Magnet Cover Type (Square or tight fit) _____							
Gradients	_____	Serial Number	_____	GE Part #		_____			
Body Coil Type	_____	Magnet Bore Measurement (cm)	_____	Upgrade Date		_____			
System	_____	Serial Number	_____	GE Part #		_____			
CERD	_____	Serial Number	_____	GE Part #		_____			
RF	_____	Serial Number	_____	GE Part #		_____			
PDU	_____	Serial Number	_____	GE Part #		_____			
Compressor	_____	Serial Number	_____	hrs.	_____	Pres.	_____		
Coldhead	_____	Serial Number	_____	Repair Date		_____			
Body Coil Chiller	_____	Serial Number	_____	GE Part #		_____			
Other	_____	Serial Number	_____	GE Part #		_____			
Other	_____	Serial Number	_____	GE Part #		_____			

Host Computer

Computer Type	_____	Serial Number	_____	Software Revision (using service browser)		_____			
Array Processor (using hw_config_status)	_____	Number of Recon Processors	_____	Recon Memory (in Mbytes)	_____	AMD Processor (using hinv)	_____	Main Memory	_____
Display Monitor	_____	Archive Media	_____	☐ HP DAT		☐ DVD/ROM/RAM			
Host IP Address	_____	Default Router	_____	Magnet Monitor IP		_____			

Extras and Options (Verify what comes with the system)

Main Disconnect	_____	DC Lighting	_____	Stereo System	_____
Injector system	_____	Chiller unit	_____	Other	_____

MRI DE-INSTALL FINAL REPORT

1290 Elm st. Suite B Ontario, Ca. 91761
Main: 951.734.3737, Fax: 951-735-3373



OEM ID _____
Site Name _____
Site Contact _____
SMR Rep _____

Van ID _____
City, State _____
Phone Number _____
Inspection Date _____

Options (Verify with system list)

- | | | | | | | |
|---|--|---|--|---|---|--|
| ☐ Echo Planar Imaging | ☐ Smart Prep | ☐ Multi-Nuclear Spectro | ☐ T2 Breathhold | ☐ ProbePRESS | ☐ BrainWave RT | ☐ Blood Flow and Volume Measurements |
| ☐ Fast Gradient Echo | ☐ SSFSE | ☐ Functool2 | ☐ SSFSE MRCP | ☐ CREATIV | ☐ BrainWave PA | ☐ Time Course |
| ☐ Cine | ☐ Three Plane Localizer | ☐ Voxtool | ☐ T1 Breathhold | ☐ Navigator | ☐ BrainWave HW | ☐ SWIFT |
| ☐ Fast Spin Echo & FLAIR | ☐ Modality Worklist | ☐ Interactive Vascular Imag | ☐ Line Scan Diffusion Imag | ☐ Diffusion Tensor | ☐ Propeller T2 | ☐ T2MAP |
| ☐ Time of Flight | ☐ E3DTOF | ☐ Clairview | ☐ Image Filter | ☐ 3DFRFSE | ☐ 3T Restricted Apps | ☐ 2D MERGE |
| ☐ Phase Contrast Vascular | ☐ FSE_XL | ☐ Spectro Analysis GE | ☐ ACGD Plus | ☐ Asset Plus | ☐ 2D Fat Sat Fiesta | ☐ LAVA-XV |
| ☐ Research | ☐ Bloodsupp | ☐ Preformed Procedure Step | ☐ Probe 3D Prostate | ☐ Fat Water 2001 | ☐ MR-Echo | ☐ SMART-P |
| ☐ Proactive Service | ☐ Fastcine | ☐ High Speed | ☐ Fluoro-triggered MRA | ☐ Rapid Frame | ☐ 8 Channel System | ☐ BRAVO |
| ☐ Research PSD | ☐ Tagging | ☐ fgret | ☐ MRCP3 | ☐ TRICKS | ☐ 16 Channel System | ☐ Breast Single Voxel Spectroscopy |
| ☐ Spectroscopy/PROBE | ☐ sgdperf | ☐ fgret Real Time | ☐ PSIR | ☐ FIESTA-c | ☐ Lava | ☐ COSMIC |
| ☐ Videoconferencing | ☐ DW-EPIZ | ☐ Spiral Hi-Res | ☐ Dynamic R1 | ☐ BREAST2 | ☐ FiberTrak | ☐ 32 Channel System |
| ☐ SGD_EchoSpeed | ☐ iDrive Pro | ☐ Spiral Real Time | ☐ Fiesta 2D | ☐ Propeller DWI | ☐ FLAIR3D | ☐ Integrated Patient Comfort Module |
| ☐ DW EPI | ☐ iDrive | ☐ iDrive Pro Plus | ☐ Fiesta 3D | ☐ 3D Delayed Enhancements | ☐ Multi-Phase (variable delays) | ☐ Brainwave Fusion |
| ☐ FLAIR EPI | ☐ SmartPrep 2000 Upgrade | ☐ PROBE 3d Brain | ☐ Delayed Enhancement | ☐ 3D Fat Sat FIESTA | ☐ Dual Array | |
| ☐ SPECIAL | ☐ Probe 2000 Upgrade | ☐ Ultrashort TR | ☐ Asset | ☐ MRAppStudio | ☐ Timeout Application | |

System Performance

	OK	Cleaned	Adjusted	Needs Correction
Does System Start-up Without Errors?.....	☐	☐	☐	☐
Does Table Function (Up/Dn/In/Out/Alignment).....	☐	☐	☐	☐
Check SNR and RF output- Head.....	☐	☐	☐	☐
APS -R1 _____ R2 _____ TG _____ AX _____				
Check SNR and RF output- Body.....	☐	☐	☐	☐
APS -R1 _____ R2 _____ TG _____ AX _____				
Check SNR (CTLTOP).....	☐	☐	☐	☐
APS -R1 _____ R2 _____ TG _____ AX _____				
RC1% RC2% RC3% RC4%				
RC5% RC6% RC7% RC8%				
Check SNR (CTLBOT).....	☐	☐	☐	☐
APS -R1 _____ R2 _____ TG _____ AX _____				
RC1% RC2% RC3% RC4%				
RC5% RC6% RC7% RC8%				
Check Shim (45 DSV / 22 DSV)(make sure to run auto L-coil and DC off-sets).....	☐	☐	☐	☐

MRI DE-INSTALL FINAL REPORT

1290 Elm st. Suite B Ontario, Ca. 91761
Main: 951.734.3737, Fax: 951-735-3373



OEM ID _____
Site Name _____
Site Contact _____
SMR Rep _____

Van ID _____
City, State _____
Phone Number _____
Inspection Date _____

Tested **Verify Coils With System List**

☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____

Tested **Verify Coils With System List**

☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____
☐		S/N	_____

S/C Shim Currents

AX1		T1-1		T2-1	
AX2		T1-2		T2-2	
AX3		T1-3		T2-3	
AX4		T1-4		T2-4	
AX5		T1-5		T2-5	
AX6		T1-6		T2-6	

Ramp Info

Magnet S/N _____

Parking Current Center Frequency
Ramp config Field Strength

- ☐ GE user manuals, log books and data sheets on site
☐ MOD options and current Save Info for system on site

Data sheets completed

- ☐ Ramp down data sheet ☐ Cryogen fill logs ☐ De-install final report

Supercon Shim

Shim Values Before X _____ Y _____ Z _____ Z2 _____

Zero Gradient Shim - (Test mode) 45cm DSV

Inhomogeneity	_____	Z1(1,0)	_____	X(1,1)	_____	Z2X(3,1)	_____
Spine	_____	Z2(2,0)	_____	Y(1,-1)	_____	Z2Y(3,-1)	_____
Shoulder	_____	Z3(3,0)	_____	ZX(2,1)	_____	Zxy2(3,2)	_____
	_____	Z4(4,0)	_____	ZY(2,-1)	_____	ZXY(3,-2)	_____
Average B field	_____	Z5(5,0)	_____	X2-Y2(2,2)	_____	X3(3,3)	_____
Std B field	_____	Z6(6,0)	_____	XY(2,-2)	_____	Y3(3,-3)	_____

With Gradient Shim - 22cm DSV

Shim Values After X _____ Y _____ Z _____ Z2 _____

Inhomogeneity	_____	Z1(1,0)	_____	X(1,1)	_____	Z2X(3,1)	_____
Spine	_____	Z2(2,0)	_____	Y(1,-1)	_____	Z2Y(3,-1)	_____
Shoulder	_____	Z3(3,0)	_____	ZX(2,1)	_____	Zxy2(3,2)	_____
	_____	Z4(4,0)	_____	ZY(2,-1)	_____	ZXY(3,-2)	_____
Average B field	_____	Z5(5,0)	_____	X2-Y2(2,2)	_____	X3(3,3)	_____
Std B field	_____	Z6(6,0)	_____	XY(2,-2)	_____	Y3(3,-3)	_____

MRI DE-INSTALL FINAL REPORT

1290 Elm st. Suite B Ontario, Ca. 91761
Main: 951.734.3737, Fax: 951-735-3373



OEM ID _____
Site Name _____
Site Contact _____
SMR Rep _____

Van ID _____
City, State _____
Phone Number _____
Inspection Date _____

Items list with stickers

Pallet/item 1	_____	Pallet/item 11	_____
Pallet/item 2	_____	Pallet/item 12	_____
Pallet/item 3	_____	Pallet/item 13	_____
Pallet/item 4	_____	Pallet/item 14	_____
Pallet/item 5	_____	Pallet/item 15	_____
Pallet/item 6	_____	Pallet/item 16	_____
Pallet/item 7	_____	Pallet/item 17	_____
Pallet/item 8	_____	Pallet/item 18	_____
Pallet/item 9	_____	Pallet/item 19	_____
Pallet/item 10	_____	Pallet/item 20	_____

Final Check list

- ☐ Everything packaged and wrapped to go with the system
- ☐ Make sure everything gets on the truck and secured so nothing will be damaged
- ☐ Make sure magnet shim lead is up and record the final helium and pressure reading
- ☐ Make sure customer is happy and everything is cleaned up when we leave

Helium _____

Pres. _____

Date _____

Time _____

Notes or Follow up Items (Verify that everything is completed)

Notes	_____
-------	-------

Signatures of Completion

Engineer	_____	Customer	_____	Manager	_____
Signature	_____	Signature	_____	Signature	_____
Date/Time	_____	Date/Time	_____	Date/Time	_____