

MRI INSTALL FINAL REPORT

1290 Elm st. Suite B Ontario, Ca. 91761
Main: 951.734.3737, Fax: 951-735-3373



OEM ID	_____	SMR ID	_____
Site Name	_____	City, State	_____
Site Contact	_____	Phone Number	_____
SMR Rep	_____	Completion Date	_____

Inventory (Guided Install Hardware Tab)

Magnet	_____	Serial Number	_____	Helium	_____	Pres.	_____	Date	_____
Package Name	_____	Magnet Cover Type (Square or tight fit) _____							
Gradients	_____	Serial Number	_____	GE Part #		_____			
Body Coil Type	_____	Magnet Bore Measurement (cm)	_____	Upgrade Date		_____			
System	_____	Serial Number	_____	GE Part #		_____			
CERD	_____	Serial Number	_____	GE Part #		_____			
RF	_____	Serial Number	_____	GE Part #		_____			
PDU	_____	Serial Number	_____	GE Part #		_____			
Compressor	_____	Serial Number	_____	hrs.	_____	Pres.	_____		
Coldhead	_____	Serial Number	_____	Repair Date		_____			
Body Coil Chiller	_____	Serial Number	_____	GE Part #		_____			
Other	_____	Serial Number	_____	GE Part #		_____			
Other	_____	Serial Number	_____	GE Part #		_____			

Host Computer

Computer Type	_____	Serial Number	_____	Software Revision (using service browser)		_____			
Array Processor (using hw_config_status)	_____	Number of Recon Processors	_____	Recon Memory (in Mbytes)	_____	AMD Processor (using hinv)	_____	Main Memory	_____
Display Monitor	_____	Archive Media	_____	☐ HP DAT		☐ DVD/ROM/RAM			
Host IP Address	_____	Default Router	_____	Magnet Monitor IP		_____			

Extras and Options (Verify what was supplied with the system)

Main Disconnect	_____	DC Lighting	_____	Stereo System	_____
Injector system	_____	Chiller unit	_____	Other	_____

MRI INSTALL FINAL REPORT

1290 Elm st. Suite B Ontario, Ca. 91761
Main: 951.734.3737, Fax: 951-735-3373



OEM ID _____
Site Name _____
Site Contact _____
SMR Rep _____

Van ID _____
City, State _____
Phone Number _____
Inspection Date _____

Options (Verify with system list)

- | | | | | | | |
|---|--|---|--|---|---|--|
| ☐ Echo Planar Imaging | ☐ Smart Prep | ☐ Multi-Nuclear Spectro | ☐ T2 Breathhold | ☐ ProbePRESS | ☐ BrainWave RT | ☐ Blood Flow and Volume Measurements |
| ☐ Fast Gradient Echo | ☐ SSFSE | ☐ Functool2 | ☐ SSFSE MRCP | ☐ CREATIV | ☐ BrainWave PA | ☐ Time Course |
| ☐ Cine | ☐ Three Plane Localizer | ☐ Voxtool | ☐ T1 Breathhold | ☐ Navigator | ☐ BrainWave HW | ☐ SWIFT |
| ☐ Fast Spin Echo & FLAIR | ☐ Modality Worklist | ☐ Interactive Vascular Imag | ☐ Line Scan Diffusion Imag | ☐ Diffusion Tensor | ☐ Propeller T2 | ☐ T2MAP |
| ☐ Time of Flight | ☐ E3DToF | ☐ Clairview | ☐ Image Filter | ☐ 3DFRFSE | ☐ 3T Restricted Apps | ☐ 2D MERGE |
| ☐ Phase Contrast Vascular | ☐ FSE_XL | ☐ Spectro Analysis GE | ☐ ACGD Plus | ☐ Asset Plus | ☐ 2D Fat Sat Fiesta | ☐ LAVA-XV |
| ☐ Research | ☐ Bloodsupp | ☐ Preformed Procedure Step | ☐ Probe 3D Prostate | ☐ Fat Water 2001 | ☐ MR-Echo | ☐ SMART-P |
| ☐ Proactive Service | ☐ Fastcine | ☐ High Speed | ☐ Fluoro-triggered MRA | ☐ Rapid Frame | ☐ 8 Channel System | ☐ BRAVO |
| ☐ Research PSD | ☐ Tagging | ☐ fgret | ☐ MRCP3 | ☐ TRICKS | ☐ 16 Channel System | ☐ Breast Single Voxel Spectroscopy |
| ☐ Spectroscopy/PROBE | ☐ sgdperf | ☐ fgret Real Time | ☐ PSIR | ☐ FIESTA-c | ☐ Lava | ☐ COSMIC |
| ☐ Videoconferencing | ☐ DW-EPIZ | ☐ Spiral Hi-Res | ☐ Dynamic R1 | ☐ BREAST2 | ☐ FiberTrak | ☐ 32 Channel System |
| ☐ SGD_EchoSpeed | ☐ iDrive Pro | ☐ Spiral Real Time | ☐ Fiesta 2D | ☐ Propeller DWI | ☐ FLAIR3D | ☐ Integrated Patient Comfort Module |
| ☐ DW EPI | ☐ iDrive | ☐ iDrive Pro Plus | ☐ Fiesta 3D | ☐ 3D Delayed Enhancements | ☐ Multi-Phase (variable delays) | ☐ Brainwave Fusion |
| ☐ FLAIR EPI | ☐ SmartPrep 2000 Upgrade | ☐ PROBE 3d Brain | ☐ Delayed Enhancement | ☐ 3D Fat Sat FIESTA | ☐ Dual Array | |
| ☐ SPECIAL | ☐ Probe 2000 Upgrade | ☐ Ultrashort TR | ☐ Asset | ☐ MRAppStudio | ☐ Timeout Application | |

System Performance

	OK	Cleaned	Adjusted	Needs Correction
Does System Start-up Without Errors?.....	☐	☐	☐	☐
Does Table Function (Up/Dn/In/Out/Alignment).....	☐	☐	☐	☐
Check SNR and RF output- Head.....	☐	☐	☐	☐
APS-R1 _____ R2 _____ TG _____ AX _____				
Check SNR and RF output- Body.....	☐	☐	☐	☐
APS-R1 _____ R2 _____ TG _____ AX _____				
Check SNR (CTLTOP).....	☐	☐	☐	☐
APS-R1 _____ R2 _____ TG _____ AX _____				
RC1% RC2% RC3% RC4%				
RC5% RC6% RC7% RC8%				
Check SNR (CLTBOT).....	☐	☐	☐	☐
APS-R1 _____ R2 _____ TG _____ AX _____				
RC1% RC2% RC3% RC4%				
RC5% RC6% RC7% RC8%				
Check Shim (45 DSV / 22 DSV)(make sure to run auto L-coil and DC off-sets).....	☐	☐	☐	☐
APS-R1 _____ R2 _____ TG _____ AX _____				

MRI INSTALL FINAL REPORT

1290 Elm st. Suite B Ontario, Ca. 91761
Main: 951.734.3737, Fax: 951-735-3373



OEM ID _____
Site Name _____
Site Contact _____
SMR Rep _____

Van ID _____
City, State _____
Phone Number _____
Inspection Date _____

Tested **Verify Coils With System List**

☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____

Tested **Verify Coils With System List**

☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____
☐		S/N _____

S/C Shim Currents

AX1	T1-1	T2-1
AX2	T1-2	T2-2
AX3	T1-3	T2-3
AX4	T1-4	T2-4
AX5	T1-5	T2-5
AX6	T1-6	T2-6

Ramp Info

Magnet S/N _____

Parking Current Center Frequency
Ramp config Field Strength

- ☐ GE user manuals, log books and data sheets on site
☐ MOD options and current Save Info for system on site

Data sheets completed

- | | | |
|---|--|--|
| ☐ Ramp up data sheet | ☐ Grafidy data sheet | ☐ Install checklist |
| ☐ Supercon shim data sheet | ☐ iso center gradcal data sheet | ☐ Install final report form |
| ☐ Gradient shim data sheet | ☐ Coil SNR data sheet | |

Shim Values Before X _____ Y _____ Z _____ Z2 _____

Zero Gradient Shim - (Test mode) 45cm DSV

Inhomogeneity _____	Z1(1,0) _____	X(1,1) _____	Z2X(3,1) _____
Spine _____	Z2(2,0) _____	Y(1,-1) _____	Z2Y(3,-1) _____
Shoulder _____	Z3(3,0) _____	ZX(2,1) _____	Zxy2(3,2) _____
	Z4(4,0) _____	ZY(2,-1) _____	ZXY(3,-2) _____
Average B field _____	Z5(5,0) _____	X2-Y2(2,2) _____	X3(3,3) _____
Std B field _____	Z6(6,0) _____	XY(2,-2) _____	Y3(3,-3) _____

With Gradient Shim - 22cm DSV

Shim Values After X _____ Y _____ Z _____ Z2 _____

Inhomogeneity _____	Z1(1,0) _____	X(1,1) _____	Z2X(3,1) _____
Spine _____	Z2(2,0) _____	Y(1,-1) _____	Z2Y(3,-1) _____
Shoulder _____	Z3(3,0) _____	ZX(2,1) _____	Zxy2(3,2) _____
	Z4(4,0) _____	ZY(2,-1) _____	ZXY(3,-2) _____
Average B field _____	Z5(5,0) _____	X2-Y2(2,2) _____	X3(3,3) _____
Std B field _____	Z6(6,0) _____	XY(2,-2) _____	Y3(3,-3) _____

MRI INSTALL FINAL REPORT

1290 Elm st. Suite B Ontario, Ca. 91761
Main: 951.734.3737, Fax: 951-735-3373



OEM ID _____
Site Name _____
Site Contact _____
SMR Rep _____

Van ID _____
City, State _____
Phone Number _____
Inspection Date _____

Final Check list

- ☐ Verify system operation, system boots up without errors or problems, everything is configured properly
- ☐ Verify all table functionality; table undock/dock, up/down, drive in/out, and dog house re-tracks
- ☐ All extra options and equipment installed and verified; Injector, camera, stereo, ect...
- ☐ All coils tested for SNR and verified with the purchase agreement
- ☐ All table pads, coil pads, positioning pads, wedges, patient straps and phantoms with system and put away
- ☐ Test network configuration, send images to their PACS and or Camera, Verify modality worklist (if acceptable)
- ☐ Communication log book and all calibration data sheets left in a book at the site with a copy for the office
- ☐ GE User manuals, coil books, MOD options with system and current save info done
- ☐ Final test scan ran for head, body and phased array with everything put away
- ☐ Final walk through to check that everything is done and clean with all covers in place
- ☐ Watch a patient scan to confirm image quality of the system before turn over
- ☐ Turn system over to the customer and have them sign off that everything is completed

Notes or Follow up Items (Verify that everything is completed)

N
o
t
e
s

Parts Used on Install

P
A
R
T
S

Signatures of Completion

Engineer _____
Signature _____
Date/Time _____

Customer _____
Signature _____
Date/Time _____

Manager _____
Signature _____
Date/Time _____