

Ventri PM Report

Document #: GE-07-00-001

Revision: A

Keywords: pm, Ventri, log

Affected Products

Ventri

Reference Material

5543208-1EN Rev. 1 Ventri Service Manual

Manufacturer's service documentation can be obtained at:

<http://apps.gehealthcare.com/servlet/ClientServlet?REQ=RNEW&MODALITY=NM>

History

Date	Name	Rev	Change descriptions
12/3/14	Document Control	A	Initial Release

Approval

Department	Name	Date
Document Control	Jehovany Hernandez	12/3/14

Table of Contents

Preface 3

System Information 4

Measuring Equipment Used 4

Checklist 5

System Condition 7

Preface

This document is valid proof of a Preventative Maintenance performed on the system. The Preventative Maintenance must be performed in the specified intervals, and the results are to be documented herein.

Warning: All steps must be performed in accordance with the manufacturer's service manual.

The section numbers indicated in the checklist correspond to instructions contained within the service manual. (see *Reference Material*).

This document must be completely filled out by the Field Service Engineer, i.e.:

- All boxes must be filled out. Use **N/A** where appropriate.
- Enter the Serial Number and the Date in the header of each page.
- After completing the check, fill out and sign the last page.

System Serial #: _____ *Date:* _____

System Information

Customer:	
Model:	
System S/N:	
System ID:	

Measuring Equipment Used

Enter the measuring equipment used in the table below (meters, scopes, phantoms, etc.).

[illegible]

Checklist

Safety Mechanisms

Good	Bad	N/A	Action	Frequency
☐	☐	☐	Pressure Sensitive Device Functionality Test	1 year
☐	☐	☐	Collimator Identification Check	1 year
☐	☐	☐	Emergency Button Shut-off Test	1 year
☐	☐	☐	Release Mechanism	1 year
☐	☐	☐	Patient Table Check	1 year

Mechanical Maintenance

Good	Bad	N/A	Action	Frequency
☐	☐	☐	Motion Noises	1 year
☐	☐	☐	Rotation Gear for leakage	1 year
☐	☐	☐	Timing Belts & Pulleys	1 year
☐	☐	☐	Cleaning the Rubber Rail Strip	2 years
☐	☐	☐	Check Rotor Covers For Cracks	1 year
☐	☐	☐	Check Table Covers for cracks	1 year

Hardware

Good	Bad	N/A	Action	1 year
☐	☐	☐	RCU Switches	1 year
☐	☐	☐	Hardware Sensors and Analog Port Status	1 year
☐	☐	☐	Switches: Check the ON/OFF panel	1 year
☐	☐	☐	Head Rest, Arm Rest and Leg Supports	1 year
☐	☐	☐	Mattress is not torn and is clean	1 year
☐	☐	☐	Examine monitor screen for scratches	1 year

Image QC

Good	Bad	N/A	Action	Frequency
☐	☐	☐	Periodic QC	1 year
☐	☐	☐	CROT Test	1 year
☐	☐	☐	COR Calibration, if DROT fails	1 year

Lubrication

Good	Bad	N/A	Action	Frequency
☐	☐	☐	Stator Lubrication	2 years
☐	☐	☐	Rotor Lubrication	1 year
☐	☐	☐	Automatic Table Lubrication	2 years

Cleanup

Good	Bad	N/A	Action	Frequency
☐	☐	☐	General Cleanup - Back of Gantry and Table Sides	1 year

Backup

Good	Bad	N/A	Action	Frequency
☐	☐	☐	Perform a full backup	1 year

System Condition

After performing scheduled PM, an evaluation of the system is performed with the following result (check one):

- ☐ No deficiencies. The image quality test resulted in no differences from required reference values.
- ☐ Slight deficiencies. There is no effect on operation; however they should be corrected preventively. The image quality test resulted in no differences from required reference values.
- ☐ Serious deficiencies. This must be corrected before operation of the system is permitted.

FSE

Signature: _____

Date: _____

Name: _____

Notes (Additional work to be performed, partial PM, customer requests):

Customer

Signature: _____

Date: _____

Name: _____