

INTRODUCTION



SYSTEM LEVEL



MAINS DISTRIBUTION



Magnet



PATIENT SUPPORT



GRADIENT AMPLIFIER



DATA ACQUISITION SYSTEM



PHYSIOLOGY



OPERATOR CONSOLE



COMPUTER



OPTIONS



CONTENTS:

Blockdiagram MR System Z1-2.1

Simplified Diagram Power Supply NTDAC / Filter Unit etc. Z1-2.2

Unit Location Z3-2.1

Total View System Filter Box Z..... Z3-2.2

Signal Cable Plan Magnet Z5-2.1

Signal Cable Plan Patient Support-MT Z5-2.2

Standard/Omni Signal Cable Plan Gradient Chain Z5-2.3

Power Signal Cable Plan Gradient Chain..... Z5-2.4

Master Signal Cable Plan Gradient Chain Z5-2.5

Signal Cable Plan RF Chain Z5-2.6

Signal Cable Plan Console Z5-2.7

Standard/Omni Mains Cable Plan Z5-2.8

Power/Master Mains Cable Plan Z5-2.9


Standard/Omni Earth Cable Plan Z5-2.10

Power/Master Earth Cable Plan Z5-2.11

Signal Cable Plan Synergy Platform..... Z5-2.12



Yellow Information Stickers..... Z9-2.1

Red Warning Stickers Z9-2.2



Intera 0.5 T			Intera 1.0 T			Intera 1.5 T		
Standard	Omni	Power	Omni	Power	Master			

BLOCKDIAGRAM MR SYSTEM



Q-ITEM	A.	REFRIGERATOR SYSTEM
768	AA	CRYO COMPRESSOR
768	AB	COLD HEAD
768	AD	CHILLER
762	AH	HEAT EXCHANGER
768	AS	COLD HEAD SHIELDING
768	AT	COMPRESSOR TRANSFORMER
Q-ITEM	B.	BACK-END EQUIPMENT
775	BA	AUDIO MODULE
775	BB	MAINS DISTRIBUTION
775	BC	AUDIO EQUIPMENT
775	BE	BACK-END MISCELLANEOUS BOX
775	BEA	AUDIO-SWITCH 2
775	BEB	RS232 FIBER 2
775	BED	BACK-PANEL
775	BEE	POWER SUPPLY
775	BEG	FANS
775	BEI	IIM UNIT
775	BEO	OPTICAL DISK DRIVE
786	BF	HARD COPY UNIT
775	BH	HOST COMPUTER GYROSCAN
775	BK	KEYBOARD WITH MOUSE
775	BL	PRINTER
775	BM	OPERATOR CONSOLE DISPLAY
785	BN	NETWORK
785	BNA	NETWORK SWITCH
785	BND	PRINT SERVER
785	BNE	PRINT SERVER POWER SUPPLY
781	BP	PATIENT OBSERVATION MONITOR
775	BQ	VIDEO CLAMP
775	BR	VIDEO RECORDER
775	BS	SLAVE MONITOR
781	BW	FUNCTIONAL BRAIN IMAGING BOX


Q-ITEM	D.	GRADIENT AMPLIFIER CABINET
762	DI	GRADIENT INTERFACE
762	DM1	AMPLIFIER-X
762	DM2	AMPLIFIER-Y
762	DM3	AMPLIFIER-Z
762	DP1	DC POWER SUPPLY 15kW COPLEY 234
762	DP2	28VDC SUPPLY FANS COPLEY 234
762	DP1	MASTER DC POWER SUPPLY 28kW COPLEY 274
762	DP2	SLAVE DC POWER SUPPLY 28kW COPLEY 274
762	DS	POWER DISTRIBUTION UNIT COPLEY 274
Q-ITEM	F.	NT DAC
769	FC	DATA ACQUISITION SYSTEM BDAS
769	FI	MAINS INLET UNIT
773	FL	RECONSTRUCTOR
773	FLA	RECONSTRUCTOR DISPLAY
773	FLB	RECONSTRUCTOR KEYBOARD
769	FM	DAC CABLE INTERFACE UNIT
764	FR	RF POWER AMPLIFIER
764	FT	RF POWER AMPLIFIER MULTI NUCLEI
764	FV	MULTI NUCLEI FILTER
Q-ITEM	H.	MAINS DISTRIBUTION UNIT
783	HA	EXTENDED MAINS DISTRIBUTION UNIT
Q-ITEM	L.	PATIENT SUPPORT
776	LA	VERTICAL ACTUATOR
776	LB	POWER SUPPLY
776	LC	PATIENT SUPPORT CONTROL MODULE
776	LE	ENCLOSURES
776	LH	HORIZONTAL DRIVE UNIT

Q-ITEM	M.	MAGNET ASSEMBLY
768	MA	MAGNET
768	MB	MAGNET INTERFACE BOX
768	MC	PICU FRONT (1) AND REAR (2)
768	ME	ENCLOSURES
781	MF	PHYSIOLOGY DISPLAY (OPTION)
762	MG	GRADIENT COIL
766	MGB	HYBRID BOX
766	MGQ	QUAD BODY COIL
768	MH	LIGHT VISORS
768	MK	BRIDGE SECTION
768	ML	BORE LIGHTS
769	MM	PFEI
769	MO	RF LOAD
768	MR	REMOTE MAGNET MONITORING UNIT
768	MY	RINGFRAME FRONT (1) AND REAR (2)
781	MVM	MICROPHONE UNIT
768	MVP	FINGER PROTECTION
781	MW1	CEILING SPEAKER FRONT
781	MW2	CEILING SPEAKER REAR
768	MY	EMERGENCY RUN/DOWN SWITCHES


Q-ITEM	P.	PATIENT SYSTEMS
781	PA	PATIENT OBSERVATION CAMERA
781	PC	COMFORT ZONE
781	PP	PHYSIOLOGY UNIT
781	PV	PATIENT VENTILATION
Q-ITEM	S.	MISCELLANEOUS
781	SA	ACCESSORY SET
776	SB	TROLLEY
781	SE	PATIENT MONITORING
787	SF	PHANTOM SET
768	SL	MEASURING ACCESSORIES
768	SM	NON MAGNETIC TOOLS
780	SO	OTHER
780	SR	RECEIVE COILS
780	SR	HEAD COIL
780	ST	TRANSMIT/RECEIVE COILS
782	SV	ADD. VIEWING STATION
778	SW	SOFTWARE
768	SX	TOOLS
Q-ITEM	Z.	SYSTEM FILTER BOX
779	ZC	COMBI FILTER
779	ZE	ENCLOSURES
779	ZF	GRADIENT FILTER
779	ZK	POWER SUPPLY
779	ZL	300V CONVERTER

Intera 0.5 T	Intera 1.0 T	Intera 1.5 T
Standard	Omni	Power
	Omni	Power
		Master


UNIT LOCATION




VIEW A---A




VIEW B---B




VIEW C---C




VIEW D---D




VIEW E---E



VIEW F---F




VIEW G---G




VIEW H---H


SIDE VIEW FRAME FILTERBOX




TOP VIEW FILTERBOX WITHOUT COVER




Detail A:



Detail B:




Detail C:



Intera 0.5 T	Intera 1.0 T	Intera 1.5 T
Standard	Unity	Power
	Unity	Power
	Unity	Master


TOTAL VIEW
SYSTEM FILTER BOX Z.



REMARKS:
(X) = DEPENDED ON LEVEL
ONLY LAST EIGHT NUMBERS SHOWN
4522 (12NC OF CABLE)

Intera 0.5 T		Intera 1.0 T		Intera 1.5 T	
Standard	Omni	Power	Omni	Power	Master

SIGNAL CABLE PLAN
MAGNET




(X) = DEPENDED ON LEVEL

Intera 05 T		Intera 10 T		Intera 15 T	
Standard	Omni	Power	Omni	Power	Master


SIGNAL CABLE PLAN
PATIENT SUPPORT MT

(X) = DEPENDED ON LEVEL



Inter 05 T	Inter 10 T		Inter 15 T		
Standard	Omni		Omni		


SIGNAL CABLE PLAN
GRADIENT CHAIN




(X) = DEPENDED ON LEVEL

Intera 0.5 T	Intera 1.0 T	Intera 1.5 T
	Power	Power

SIGNAL CABLE PLAN
GRADIENT CHAIN






MR/Intera-DRW-71

99.0

Z5-2.6


SIGNAL CABLE PLAN RF CHAIN


Intera 0.5 T	Intera 1.0 T	Intera 1.5 T
Standard	Omni	Power
	Omni	Power
		Master





Intera 05 T		Intera 10 T		Intera 15 T	
Standard	Omni	Power	Omni	Power	Master

SIGNAL CABLE PLAN
CONSOLE









REMARKS:

- (X) = DEPENDENT ON LEVEL
4522 ... (12NC OF CABLE)
- * = OPTIONAL UNIT
** = LOCAL FOR LOCAL
= EARTH WIRES FROM MAINS CABLES


Intera 0.5 T	Intera 1.0 T	Intera 1.5 T
	Power	Power Master

EARTH CABLE PLAN

MR/Intera-DRW-71

00.0

Z5-2.11



(X) = DEPENDED ON LEVEL

Intera 0.5 T			Intera 1.0 T			Intera 1.5 T		
Standard	Omni	Power	Omni	Power	Master			

SIGNAL CABLE PLAN SYNERGY PLATFORM

