

## Sensation

**CT**

### Troubleshooting Guide

WCS (water/air)

Troubleshooting Guide WCS (water/air)

This document is valid for:  
SOMATOM Sensation 10/16/Cardiac  
SOMATOM Sensation 40/64/Cardiac 64/Open

© Siemens 2002

The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights, including rights created by patent grant or registration of a utility model or design, are reserved.

Print No.: CT02-023.840.06.05.02  
Replaces: CT02-023.840.06.04.02

English  
Doc. Gen. Date: 06.05

**Document revision level**

The document corresponds to the version/revision level effective at the time of system delivery. Revisions to hardcopy documentation are not automatically distributed.

Please contact your local Siemens office to order current revision levels.

**Disclaimer**

The installation and service of equipment described herein is to be performed by qualified personnel who are employed by Siemens or one of its affiliates or who are otherwise authorized by Siemens or one of its affiliates to provide such services.

Assemblers and other persons who are not employed by or otherwise directly affiliated with or authorized by Siemens or one of its affiliates are directed to contact one of the local offices of Siemens or one of its affiliates before attempting installation or service procedures.

|          |                                                          |           |
|----------|----------------------------------------------------------|-----------|
| <b>1</b> | <b>General</b>                                           | <b>4</b>  |
|          | Safety information. . . . .                              | 4         |
|          | Tools . . . . .                                          | 5         |
| <b>2</b> | <b>WCS:water/air type</b>                                | <b>6</b>  |
|          | Overview. . . . .                                        | 6         |
|          | Temperature value . . . . .                              | 7         |
|          | Nominal temperature value. . . . .                       | 7         |
|          | Current temperature values . . . . .                     | 8         |
|          | If current temperature is out of tolerance . . . . .     | 10        |
|          | Electrical circuit . . . . .                             | 11        |
|          | N1 - WCS controller . . . . .                            | 12        |
|          | Other components" parts in electrical circuit. . . . .   | 16        |
|          | Water circuits . . . . .                                 | 17        |
|          | Draining and filling the cooling system . . . . .        | 17        |
|          | Reasons for a drop in pressure. . . . .                  | 19        |
|          | Testing for leaks with WCS switched off . . . . .        | 19        |
|          | Testing for leaks in individual WCS components . . . . . | 20        |
|          | Additional hints . . . . .                               | 23        |
| <b>3</b> | <b>Changes to previous version</b>                       | <b>24</b> |

## Safety information

Read the following instructions as well as the most recent version of

- **General Safety Notes** (TD00-000.860.01.xx.xx)
- **Product-specific safety Information** (CT02-023.860.01.xx.xx)

before trouble shooting the WCS.

The following instructions and the above mentioned **Safety documentation** contain important safety information as well as information about tools to be used.

|             |
|-------------|
| <b>NOTE</b> |
|-------------|

---

**Only use the <Read> button in the functions <Table load/modify>. All other functions are not necessary for service, and are normally for development staff only.**

**=> Improper handling can result in system malfunctions.**

---

|             |
|-------------|
| <b>NOTE</b> |
|-------------|

---

**Operating the WCS water pump without water (e.g. during filling/refilling the gantry) will destroy the pump.**

**=> Do not operate the water pump without water!**

---

|             |
|-------------|
| <b>NOTE</b> |
|-------------|

---

**Do not operate the system in customer mode at the increased water pressure as used for the test. Always return the water pressure to the nominal operating pressure of 2.0 bar  $\pm$  0.2 bar, with the heat exchanger switched off.**

---

## Tools

- Standard tool kit
- 2 water hoses with 1/2" connections for filling and draining the system
- Manometer (included with the delivery of WCS)
- Absorbent cloth
- Blotting paper

## Overview

The Water Cooling System has to ensure the proper temperature inside the gantry. Therefore, errors caused by a defective WCS may generate:

- temperature error messages for components within the gantry
- temperature error messages for the MCU
- image quality problems

This section provides additional troubleshooting information, in case information in either the detailed message or functional descriptions is insufficient.




Fig. 1: Cooling circuit overview

## Temperature value

### Nominal temperature value




Fig. 2: Temperature Sensor in Gantry

Pos. 1 B 305

Pos. 2 B 302

### Internal gantry temperature

- ⇒ B302 = the actual temperature value, used for regulation.  
B305 = the safety temperature switch. In case of overheating ( $> 43^{\circ}\text{C}$ ), the PDC switches off the power to PSS.
- ⇒ Nominal temperature:  
 $28^{\circ}\text{C} \pm 2^{\circ}\text{C}$
- ⇒ Thresholds  
 $< 35^{\circ}\text{C}$  (B302)  
 $< 43^{\circ}\text{C}$  (B305)

### Temperature inside the WCS




Fig. 3: Temperature Sensor in WCS

Pos. 1 R1  
Pos. 2 R3



Fig. 4: Temperature sensor in the air inlet

Pos. 1 clip  
Pos. 2 R2

- ⇒ Temperature Sensor in WCS:  
R1 = the outlet water temperature value, used for regulation.  
R2 = the inlet air temperature value.  
R3 = the refrigerant temperature value, used to protect against freezing.
- ⇒ Nominal temperature:  
R1 = 17.5°C to 19.5°C.  
R2 = -25.0°C to +40°C.

### Current temperature values

The current temperature values of the various sensors can be evaluated via the syngo service software.

#### How to obtain the current temperature values

1. From the user interface, select **Options -> Local Service**.
2. In the authentication page, enter the **Service Password** and Click **OK**.
3. In the Home menu, click the **Control** button.



4. Click **Table load/modify**:

⇒ The following page appears.




Fig. 5: Table load / modify page

## 5. To access the control system

- Select **T00** for TABLE NAME
- Select **CONTROL SYSTEM** for TABLE SOURCE

⇒ Click **GO**, and the following page appears.




Fig. 6: Example of Parameter Set 1

6. Select **Dec**, and click **Read**:

**NOTE**

Only use the <Read> button in the functions <Table load/modify>. All other functions are not necessary for service, and are normally for development staff only.

⇒ Improper handling can result in system malfunctions.

⇒ The corresponding table will be activated.

Gantry temp. (B302) = <gantry temperature>

Water temperature in secondary cooling circuit (R1) = <WCS water temp.>

More important information is available, please refer to (Tab. 2 / p. 15).

**If current temperature is out of tolerance****Gantry temperature**

Problems may exist with the following:

- Secondary cooling circuit
  - ⇒ Refer to ([Secondary cooling circuit / p. 10](#))
- Cooling fan in gantry
  - ⇒ Check or replace the fan for the heat exchanger as described in replacement instructions  
CT02-023-841.01.XX.XX
- Temperature sensor (B302)
  - ⇒ Check or replace the temperature sensor.

**Secondary cooling circuit**

Problems may exist with the following:

- WCS
  - ⇒ Refer to ([Electrical circuit / p. 11](#))
- Water leaks
  - ⇒ Refer([Water circuits / p. 17](#))
- Air flow
  - ⇒ Check the air flow circuits (e.g., radial fan, condenser, inlet/outlet of air, air hoses, air temperature sensor (R2)).
  - ⇒ Check or replace the radial fan or the temperature sensor (R2) as described in replacement instructions  
CT02-023-841.06.XX.XX
- Temperature sensor (R1)
  - ⇒ Check or replace the temperature sensor (R1) as described in replacement instructions.  
CT02-023-841.06.XX.XX

## Electrical circuit




Fig. 7: WCS Water / Air Electrical Circuit Overview 1




Fig. 8: WCS Water / Air Electrical Circuit Overview 2

## N1 - WCS controller




Fig. 9: N1 - WCS Controller (schematic layout)

For further investigations within the WCS, it is necessary to know the current status of the N1 WCS controller. Obtain the actual status from the LED on N1 ([LED"s on the N1 controller / p. 12](#)) or the system control table via the service software ([The system control table / p. 13](#)).

## LED"s on the N1 controller

Tab. 1 LED,s on the N1 controller

| LED / color | Normal status | Description for WCS (water-air)            |
|-------------|---------------|--------------------------------------------|
| L1 - green  | blinking      | Program active                             |
| L2 - yellow | off           | Service active                             |
| L3 - n.a.   | n.a.          | not mounted                                |
| L4 - n.a.   | n.a.          | not mounted                                |
| L5 - green  | on, if -->    | Communication active<br>(telegram traffic) |
| L6 - n.a    | n.a.          | not mounted                                |
| L7 - yellow | off           | on, if --> malfunction                     |
| L8 - green  | on, if -->    | Condensation pressure2 is reached          |
| L9 - green  | on, if -->    | Condensation pressure1 is reached          |

| LED / color | Normal status | Description for WCS (water-air)                             |
|-------------|---------------|-------------------------------------------------------------|
| L10 - green | on            | off, if pressure too low<br>(pressure in coolant circuit)   |
| L11 - green | on            | off, if pressure too high<br>(pressure in coolant circuit)  |
| L12 - green | on, if -->    | 50Hz<br>off, if --> 60Hz                                    |
| L13 - green | n.a.          | not used                                                    |
| L14 - green | on            | Cooling pump protection F3 / OK<br>off, if --> over current |
| L15 - green | on            | Thermocontact: ventilator                                   |
| L16 - red   | off           | on, if --> MCU reset occurs                                 |
| L17 - green | on, if -->    | Compressor 1 active                                         |
| L18 - green | on, if -->    | Compressor 2 active                                         |
| L19 - green | on            | Cooling pump running                                        |
| L20 - green | n.a.          | not used                                                    |
| L21 - green | on, if -->    | Ventilator is working in Y                                  |
| L22 - green | on, if -->    | Ventilator is working in D                                  |
| L23 - green | on, if -->    | Cooling liquid regulation: colder                           |
| L24 - green | on, if -->    | Cooling liquid regulation: warmer                           |

#### The system control table

1. From user interface, select **Options -> Local Service**.
2. In the authentication page, enter the **Service Password** and Click **OK**.
3. In the Home menu, click the **Control** button.

4. Click **Table load/modify**:

⇒ The following page appears.




Fig. 10: Table load / modify page

## 5. To access the control system

- Select **T00** for TABLE NAME
- Select **CONTROL SYSTEM** for TABLE SOURCE

⇒ Click **GO**, and the following page appears




Fig. 11: Example of Parameter Set 1

6. Select **Dec**, and click **Read**

⇒ The corresponding table will be activated.

**NOTE**

Only use the <Read> button in the functions <Table load/modify>. All other functions are not necessary for service, and are normally for development staff only.

=> Improper handling can result in system malfunctions.

Tab. 2 T00 table contents

|                                | Normal condition                                      |
|--------------------------------|-------------------------------------------------------|
| gantry temperature             | current temperature inside gantry                     |
| gantry humidity                | current humidity inside gantry                        |
| WCS water temperature          | current temperature in secondary cooling circuit (R1) |
| compressor 1 active            | On, if compressor 1 active                            |
| compressor 2 active            | On, if compressor 2 active                            |
| ventilator Y/D active          | On                                                    |
| cooling pump active            | On                                                    |
| cooling med out of tolerance   | Off                                                   |
| cooling water out of tolerance | Off                                                   |
| cooling system out of tol.     | Off                                                   |
| cooling system defective       | Off                                                   |
| no water flow                  | Off                                                   |
| motor protector pump           | Off                                                   |
| motor protector ventilator     | Off                                                   |
| primary cooling media temp     | Off / Only for W/W                                    |

**NOTE**

**This table T00 include the current temperature value. Please check the tolerance as described in [\(Temperature value / p. 7\)](#).**

## Other components" parts in electrical circuit




Fig. 12: Other components parts

Tab. 3 Fuse, Measurement points and LED's

|    | Description                                               |
|----|-----------------------------------------------------------|
| K1 | Fan and compressor 1 contactor                            |
| K2 | Compressor 1 contactor                                    |
| K3 | Cooling pump contactor                                    |
| K4 | Fan star connection contactor                             |
| K5 | Fan delta connection contactor                            |
| F1 | Pressure control high protection<br>(refrigerant circuit) |
| F3 | Cooling pump over current protection                      |
| F4 | Fan motor over current protection                         |
| F5 | Pressure control low protection<br>(refrigerant circuit)  |
| F6 | Fan star connection (slow speed)                          |
| F7 | Fan delta connection (fast speed)                         |
| F8 | Circuit breaker 3A                                        |
| F9 | Circuit breaker 1A                                        |
| T1 | Transformer                                               |



## Water circuits




Fig. 13: WCS Water / Air Circuit Overview




Fig. 14: WCS Water / Air Refrigerant Circuit Overview

## Draining and filling the cooling system

Testing for leaks requires increasing or decreasing the water pressure.



### Draining process to decrease the water pressure

1. Switch off the WCS via breaker F14 in the PDC.
2. Connect the water hoses to the filling and draining valves (4/Fig. 15 / p. 18).
3. Route the water hoses to a syphon or a bucket.
4. The red shut-off valve for the gantry cooling system must be open.
  - ⇒ The lever is lined up with the connection pipe (1/Fig. 15 / p. 18).




Fig. 15: Shut-off valve at the WCS "Gantry water inlet"

- |        |                            |
|--------|----------------------------|
| Pos. 1 | Open                       |
| Pos. 2 | Close                      |
| Pos. 3 | Shut-off valve             |
| Pos. 4 | filling and draining valve |

5. Connect the manometer at the WCS "gantry water outlet".
6. Open the filling and draining valves (4/Fig. 15 / p. 18) by turning them 90° counterclockwise using the valve cap.
  - ⇒ The valve is open when the notch on the valve is lined up with the connection pipe.
7. Drain until the water pressure reaches  $2.0 \pm 0.2$  bar, and close the filling and draining valves (4/Fig. 15 / p. 18) by turning them 90° clockwise using the valve cap.
8. Switch off the WCS via breaker F14 in the PDC.

#### NOTE

**After the WCS has been running for a while, please check that the water pressure is  $2.0 \pm 0.2$  bar, with WCS switched off.**

### Filling the cooling system

Please always refer to the latest version of

Installation Instructions CT02-023.814.01.xx.xx.: Filling the cooling system.

## Reasons for a drop in pressure

### After installation

Drop in pressure may occur once or twice after start-up. After filling the system a maximum of 3 times, pressure should remain constant. The pressure should not drop by more than 1.0 bar in 6 months.

Possible causes of pressure drop may include:

- Excessive amounts of air in the cooling system.
- Warm water was used to fill the system (loss in volume when cooled to 18°C).
- The hose for filling the system was not filled with water when connected to the WCS.
- Widening of water hoses within the first weeks after system installation.

### After the system has been filled several times

After 6 months the drop in pressure is >1.0 bar, and the minimum value of the static pressure is 1.0 bar.

Possible causes of pressure drop may include:

- Screw connections are not water-tight.
- Leaks in the affected system components (e.g., in the water pump, 3-way valve, Spirovent, expansion reservoir, heat exchanger, connection WCS <--> gantry, within gantry). Continue with [\(Testing for leaks with WCS switched off / p. 19\)](#).

## Testing for leaks with WCS switched off

|             |
|-------------|
| <b>NOTE</b> |
|-------------|

**Do not operate the system in customer mode at the increased water pressure as used for the test. Always return the water pressure to the nominal operating pressure of 2.0 ± 0.2 bar, with the WCS switched off.**

**=> Please refer to [\(Draining process to decrease the water pressure / p. 18\)](#).**

|             |
|-------------|
| <b>NOTE</b> |
|-------------|

**Do not perform scans at this time.**



1. Switch off the WCS via breaker F14 in the PDC.
2. Connect the filling hose to the cooling system.
  - ⇒ Fill the hose with water prior to connecting it to the system to avoid introducing air into the system. Air in the filling hose will lead later to a loss of pressure.
3. Fill it until the pressure reaches a maximum (< 7.5 bar) as described in [\(Filling the cooling system / p. 18\)](#).

**NOTE**

The maximum pressure depends on the water supplier and the location of the site. In case the pressure is not able to reach 7.5 bar, test will be performed with your maximum.

4. Record the water pressure in the WCS:
  - ⇒ If the pressure drops by more than 10% within 1 hour, there is a leak in the system.
5. Inspect the WCS, the gantry, and the full length of all the water hoses for possible leaks (use blotting paper).
6. If leaks are found, replace the corresponding parts. Please refer to the respective replacement instructions:
  - Replacement of the Parts: Gantry CT02-023.841.01.xx.xx
  - Replacement of the Parts: WCS CT-023.841.06.xx.xx
7. After replacing the parts, make sure that leaks are not found using this test.
8. If leaks are not found with this test, individual components can be tested.

**NOTE**

Always return the water pressure to the nominal operating pressure of  $2.0 \pm 0.2$  bar, with the WCS switched off. Please refer to [\(Draining process to decrease the water pressure / p. 18\)](#).

## Testing for leaks in individual WCS components

**NOTE**

Do not operate the system in customer mode at the increased water pressure as used for the test. Always return the water pressure to the nominal operating pressure of  $2.0 \pm 0.2$  bar, with the WCS switched off.

=> Please refer to [\(Draining process to decrease the water pressure / p. 18\)](#).

**NOTE**

Do not perform scans at this time.

If leaks are not found with the previous test using increased pressure, individual components can be tested as follows.

### 3-way valve in gantry right stand and WCS

Overview: location of the 3-way valves



Fig. 16: 3-Way Valve in Gantry stand

Proceed as follows:

1. Remove the black insulation covering the 3-way valve.
2. Place absorbent cloth under the 3-way valve, because water may escape during test.
3. Open and close both 3-way valves approx. 10 times. A valve can be opened and closed by turning the white screw of the adjustment drive CW (open) or CCW (closed) with a 5 mm Allen key ([1/Fig. 17 / p. 21](#)).




Fig. 17: Example of Adjustment drive

4. Inspect the 3-way valve for leaks (use blotting paper).
5. If leaks are found, replace the corresponding 3-way valve. Please refer to the respective replacement instructions:  
3-Way Valve in Gantry Stand CT02-023.841.01.xx.xx  
3-Way Valve in WCS CT-023.841.06.xx.xx
6. After replacing the 3-way valve, make sure that leaks are not found using this test.
7. Drain until the water pressure reaches  $2.0 \pm 0.2$  bar as described in ([Draining process to decrease the water pressure / p. 18](#)).
8. Install the black insulation covering the 3-way valve.

**NOTE**

Always return the water pressure to the nominal operating pressure of  $2.0 \pm 0.2$  bar, with the WCS switched off. Please refer to [\(Draining process to decrease the water pressure / p. 18\)](#).

**Cooling pump**

Overview: location of the devices

00-04/



Fig. 18: Cooling pump




Fig. 19: Shut-off valve

**NOTE**

Do not perform the following test with increased water pressure. Even a short-term test with increased water pressure will destroy the water pump. Ensure that the water pressure is the nominal operating pressure of  $2.0 \pm 0.2$  bar, with the WCS switched off.

=> Please refer to [\(Draining process to decrease the water pressure / p. 18\)](#).

Proceed as follows:

1. Slide the black insulation sleeves away from the nuts.
2. For approx. 5 minutes, close the shut-off valve ([2/ Fig. 19 / p. 22](#)) of the cooling pump being tested. This forces the pump to work under increased pressure.
3. Inspect the pump for leaks (use blotting paper).
4. If leaks are found, replace the cooling pump as described in the replacement instructions:  
Replacement of Parts CT-023.841.06.xx.xx
5. After replacing the cooling pump, make sure the leaks are not found with this test.

## Additional hints

### Pressure Value

- The nominal water pressure should be  $2.0 \pm 0.2$  bar, with the WCS switched OFF.
- When the WCS is switched off (pump off), pressure levels at the water inlet (in front of the pump) and at the water outlet (behind the pump) are equal.
- When the WCS is switched on (pump on), a relative pressure difference of approx. 2.5 bar can be measured in front and behind the pump.
- The absolute pressure measured via a manometer depends on the length of the hose and the type of operation (50/60 Hz).

### Anti-freeze agent

- It is prohibited to add an anti-freeze agent to the secondary (gantry) cooling circuit.

### Others

- Excessive amounts of air in the cooling system (after filling it) cause the pump to over-heat, resulting in damage to the motor.
- To operate an additional external water pump or valve, a potential-free contact (max. 1 Amp, 220 Volts) is available in the PDC - PDS.

| Chapter | Change               | Reason        |
|---------|----------------------|---------------|
| FMP     | SOMATOM 40 was added | 136660 (P30G) |